A

Siemens Healthcare Diagnostics Sàrl, 8047 Zürich

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Dimension Vista is a registered trademark of Siemens Healthcare Diagnostics.

Urgent Field Safety Notice

VC-18-01.A.OUS

November 2017

Dimension Vista®

LOCI 8 CAL Lot 7CD076

Calibration Failures or QC/ Patient Sample Shifts on the Estradiol Assay

Our records indicate that your facility may have received the following product:

Table 1. Dimension Vista® Affected Product

Calibrator	Catalog Number	Siemens Material Number (SMN)	Lot Number	Expiration Date	1 st Distribution Date
Dimension Vista® LOCI 8 CAL	KC646	10488192	7CD076	2018-04-01	2017-05-10

Reason for Correction

The purpose of this communication is to inform you of an issue with the product indicated in Table 1 above and provide recommended actions for your laboratory.

Siemens Healthcare Diagnostics has confirmed an increase in the number of failed calibration events or negative bias with Quality Control (QC) and patient samples when using Dimension Vista® LOCI 8 Calibrator lot 7CD076 to calibrate Estradiol (E2) Assay.

As a result of the issue on LOCI 8 calibrator lot 7CD076 the following two scenarios may occur:

1. Estradiol (E2) calibration failures due to slope and percent deviation being out of acceptable range (slope recovering >1.05 and percent deviation recovering >5%).
2. In the event an acceptable calibration is achieved with the LOCI 8 Calibrator lot 7CD076, customers may observe an overall negative bias for Estradiol with Quality Control (QC) and patient samples. Depending on the target ranges used, QC results may fall outside of acceptable ranges.

A summary of bias results is provided in Table 2.

Table 2: Estradiol (E2) Serum Sample % Bias

E2 Concentrations	Average bias	Potential Range of Bias
< 80 pg/mL (293.7 pmol/L)	-32%	-33% at a concentration of 51.3 pg/mL (188.3 pmol/L) -30% at a concentration of 74.9 pg/mL (275.0 pmol/L)
80 pg/mL (293.7 pmol/L) to 200 pg/mL (734.2 pmol/L)	-23%	-26% at a concentration of 105.5 pg/mL (387.3 pmol/L) -19% at a concentration of 185.1 pg/mL (679.5 pmol/L)
201 pg/mL (737.9 pmol/L) to 500 pg/mL (1835.5 pmol/L)	-10%	-15% at a concentration of 242.5 pg/mL (890.2 pmol/L) -6% at a concentration of 465.9 pg/mL (1710.3 pmol/L)
>500 pg/mL (1835.5 pmol/L)	6%	+7% at a concentration of 693.5 pg/mL (2545.8 pmol/L) -1% at a concentration of 1343.3 pg/mL (4931.3 pmol/L)

The LOCI 8 Calibrator is used to calibrate three other Dimension Vista assays; Follicle Stimulating Hormone (FSH), Luteinizing Hormone (LH), and Prolactin (PRL). These other assays continue to perform without experiencing this issue.

Risk to Health

The probability of misinterpretation of estradiol results due to the observed negative bias is remote and limited to the low end of the measuring range (e.g. <80 pg/mL [293.7 pmol/L]). Estradiol is not used in isolation but is correlated with clinical history and presentation, as well as additional diagnostic laboratory testing. Siemens is not recommending a review of previously generated results.

Actions to be Taken by the Customer

- Calibrate the E2 assay using an alternate lot of LOCI 8 CAL.
- Discontinue use of and discard LOCI 8 CAL lot 7CD076.
- Review your inventory of this product to determine your laboratory's replacement needs and to provide information to Siemens for reporting to the authorities.
- Complete and return the Field Correction Effectiveness Check Form attached to this letter within 30 days. Please review this letter with your Medical Director.
- If you have received any complaints of illness or adverse events associated with the products listed in Table 1, immediately contact your local Siemens Customer Care Center or your local Siemens Technical Support Representative.

Please retain this letter with your laboratory records, and forward this letter to those who may have received this product. We apologize for the inconvenience this situation may cause. If you have any questions, please contact your Siemens Customer Care Center or your local Siemens Technical Support Representative.

Dimension Vista is a registered trademark of Siemens Healthcare Diagnostics