# RANDOX

## Urgent Field Safety Notice

**Date:** 16 Fev 2018

**Référence:** 311

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## Urgent Field Safety Notice

**Contact Reference:**

Radox Technical Services  
Radox Laboratories Ltd,  
55 Diamond Road,  
Crumlin,  
United Kingdom,  
BT29 4QY

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
## RX Instrument Carryover Avoidance

<b>Technical Bulletin No.:</b>	RXTB-0083	<b>Issue Date:</b>	16 Feb 2018
<b>Affected Analysers:</b>	ALL RX Analysers		
<b>Importance:</b>	MEDIUM		
<b>Items Required:</b>	N/A		
<b>Specialist Tools/Software:</b>	N/A		
<b>Software Update Required:</b>	N/A		

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Column 1	Column 2
Albumin	Iron
Glucose GODPAP	Phosphate
Uric acid	Phosphate
ALT	Phosphate
AST	Phosphate
LD	Phosphate
CK	Phosphate
ALP (AMP)	Magnesium
ALP (DEA)	Magnesium
Glucose Hexokinase	Magnesium
CKMB	Magnesium
Triglycerides	Magnesium
Cholesterol	Magnesium

Uric Acid	Magnesium
Potassium	Sodium
Direct Bilirubin	Sodium
Transferrin	Sodium
Total Protein	Sodium
Total Protein	Potassium
Total Protein	Copper
Transferrin	Chloride
Direct Bilirubin	Chloride
Cholesterol	Lipase
Triglycerides	Lipase
Potassium	LDH
Potassium	GLDH
Urea	GLDH
Creatinine*	CRP
Fructosamine	Bile acids
Amylase	Magnesium
Pancreatic amylase	Magnesium
Liquid CO2	Magnesium
Direct LDL	Lipase
Cholesterol	Calcium
TIBC	Iron
Myoglobin	HFABP
Cystatin C	HFABP
Adiponectin	HFABP
Micro-albumin	Calcium

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- C1 wash – RX8143
- Acid Wash – WS8397

Method 1	Method e2	Réactif Pipette	Solution Lavage
Albumin	Iron	R1>R1	C1 Wash
Glucose Oxidase	Inorganic phosphorus	R1>R1	C1 Wash
Uric Acid	Inorganic phosphorus	R1>R1 R2>R2	C1 Wash
Cholesterol	Calcium	R1>R1	Acid Wash
Albumin	Calcium	R1>R1	C1 Wash
TIBC	Iron	N/A	Amend running order/ test separately
CKMB	Magnesium	R1>R1	C1 Wash
Triglycerides	Magnesium	R1>R1	C1 Wash
CK	Magnesium	R1>R1	Acid Wash
ALP AMP	Magnesium	R1>R1	Acid Wash
ALP DEA	Magnesium	R1>R1	Acid Wash
Glucose Hexokinase	Magnesium	R1>R1	Acid Wash
Amylase	Magnesium	R1>R1	Acid Wash
LCO <sub>2</sub>	Magnesium	R1>R1	Acid Wash
Micro-albumin	Calcium	R1>R1	Acid Wash

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**Randox Customer Services Action Centre: +44 (0) 28 9445 1070**

**[technical.services@randox.com](mailto:technical.services@randox.com)**