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Avis de produit (Field Safety Customer Notification SYC18-02.A.OUS)

Syva® Emit® II Plus 6-Acetylmorphine (6-AM) Method with Incorrect Calibrator Level Listed in the Application Sheet

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Robert Schlatter
Head Quality & Technology

Patricia Vaith
Manager Product Marketing

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Syva[®] Emit[®] II Plus 6-Acetylmorphine (6-AM) Method with Incorrect Calibrator Level Listed in the Application Sheet

Our records indicate that your facility may have received the following product:

Table 1. Dimension[®] Systems Affected Product(s): (Not lot specific)

Assay	Catalog Number	Siemens Material Number (SMN)
Syva [®] Emit [®] II Plus 6-Acetylmorphine (6-AM)	9R129UL (1000mL)	10481483
Syva [®] Emit [®] II Plus 6-Acetylmorphine (6-AM)	9R039UL (28mL)	10470440

Reason for Urgent Field Safety Notification

The purpose of this communication is to inform you of an issue with the Dimension[®] Systems Application Sheet for the products indicated in Table 1 above and provide instructions on actions that your laboratory must take.

Siemens Healthcare Diagnostics has confirmed that the Application Sheet (Rev.D 6AM.1, 2012-08-08 msw:967q) for the Syva[®] Emit[®] II Plus 6-Acetylmorphine (6-AM) method has the incorrect calibrator level listed in the Qualitative Calibration Steps described on page 4 for the Dimension[®] EXL[™] 200, Dimension[®] EXL[™] with LM, Dimension[®] RxL Max[™], Dimension[®] Xpand[®] Plus, Dimension[®] Xpand[®] Systems.

In Step 2, cups 2 and 3 incorrectly list Calibrator Level “3”. The correct cut-off Calibrator Level to use for this application is Level 2 as indicated on Page 2 of the Application Sheet and in the Instructions for Use for the 6-AM method.

The use of Calibrator Level 3 when performing a qualitative calibration on the Dimension analyzers may result in the reporting of patient results as falsely negative for samples between 10 and 15 ng/mL. Quality Control (QC) material may detect this issue if you are using QC material that is near the cutoff (10 ng/mL).

The Qualitative Calibration section (page 4, step 2) of the Application Sheet will be corrected to indicate calibrator level 2 must be used in cups 2 and 3.

Risk to Health

The risk to health due to this issue is negligible. This assay provides only a preliminary analytical test result. Siemens is not recommending a review of previously generated results.

Actions to be Taken by the Customer

- Please review this letter with your Medical Director.
- Ensure that you are using Calibrator Level 2 during qualitative calibration of the Syva® Emit® II Plus 6-Acetylmorphine (6-AM) method on the Dimension® EXL™ 200, Dimension® EXL™ with LM, Dimension® RxL Max™, Dimension® Xpand® Plus, Dimension® Xpand® Systems.
- Complete and return the Field Correction Effectiveness Check Form attached to this letter within 30 days.

Please retain this letter with your laboratory records, and forward this letter to those who may have received this product.

We apologize for the inconvenience this situation may cause. If you have any questions, please contact your Siemens Customer Care Center or your local Siemens technical support representative.

Syva Emit and Dimension are registered trademarks of Siemens Healthcare Diagnostics.