
Urgent Field Safety Notice

THERAKOS™ CELLEX™ Photopheresis System

FSCA-identifier: 2018-08-03 CELLEX Labeling V34584

Type of action: Advice given by MANUFACTURER regarding the use of the device and/or the follow up of patients, users, or others

Date: 09-AUG-2018

Attention: Photopheresis Department

Details on affected devices:

The CELLEX™ Photopheresis System Operator's Manual is affected by this Field Safety Notice.

The Operator's Manual numbers for the CELLEX™ instrument are as identified below:

| Language | Software Version | Material Number |
|--------------------|-------------------------|------------------------|
| English - US | 3.0 | 1460415 |
| English - US | 5.1 | 1470096 |
| Italian | 4.1 | 1460440 |
| Italian | 5.1 | 1470103 |
| German | 4.1 | 1460438 |
| German | 5.1 | 1470100 |
| French | 4.1 | 1460437 |
| French | 5.1 | 1470106 |
| Turkish | 4.1 | 1460547 |
| Turkish | 5.1 | 1470154 |
| Swedish | 4.1 | 1460540 |
| Swedish | 5.1 | 1470136 |
| Portuguese | 4.1 | 1460454 |
| Portuguese | 5.1 | 1470163 |
| Norwegian | 4.1 | 1460539 |
| Norwegian | 5.1 | 1470145 |
| Greek | 4.1 | 1460529 |
| Greek | 5.1 | 1470181 |
| Finnish | 4.1 | 1460537 |
| Finnish | 5.1 | 1470172 |
| Dutch | 4.1 | 1460517 |
| Dutch | 5.1 | 1470190 |
| Danish | 4.1 | 1460526 |
| Danish | 5.1 | 1470127 |
| English – UK | 4.1 | 1460436 |
| English – UK | 5.1 | 1470097 |
| Spanish | 4.1 | 1460439 |
| Spanish | 5.1 | 1470109 |
| Hungarian | 5.1 | 1470407 |
| Croatian | 5.1 | 1470437 |
| Canadian – English | 3.0 | 1460451 |
| Canadian – French | 3.0 | 1460452 |



Description of the problem:

Mallinckrodt has received reports of thromboembolic events associated with the use of the THERAKOS™ CELLEX™ Photopheresis System in the treatment of Graft versus Host Disease (GvHD). Patients with GvHD have an increased risk of thromboembolic events. The Caution Statement will be updated in the device Operator's Manual section, titled "Anticoagulation."

Advise on action to be taken by the user:

A Technical Bulletin will be issued as follows:

The following CAUTION statement will be added to "Anticoagulation" section of the device Operator's Manual:

"Special attention to adequate anticoagulation is advised when treating patients with Graft versus Host Disease (GvHD), a condition associated with an increased risk of thromboembolic events.

Thromboembolic events, including pulmonary embolism and deep vein thrombosis, have been reported with the use of the THERAKOS™ CELLEX™ Photopheresis System in the treatment of GvHD, an indication not approved in some countries, including the United States and Canada."

Once received, please add this Technical Bulletin to your Operator's Manual. The Technical Bulletin will be released as follows:

- English Language: issued by 01-OCT-2018
- Translations/Non-English Language: issued by 01-DEC-2018

The Technical Bulletin will be posted to mytherakos.com once approved and released.

No product will need to be returned.

Transmission of this Field Safety Notice:

This Field Safety Notice should be communicated to all personnel within your organisation who in any way are engaged in the administration of the THERAKOS™ CELLEX™ Photopheresis System.

Contact reference person:

Megan Vernak, Director, Product Monitoring – Specialty Brands
Mallinckrodt Pharmaceuticals
1425 US Route 206
Bedminster, NJ 07921

Customer Care (North America): 877-566-9466

Customer Care (Rest-of-World)

By language:

English: +441548600009

French: +33182880867

German: +4932221093619

Spanish: +34932202094

Italian: +39 051 042 0666

ecphelp@therakos.com

The undersigned confirms that this notice has been notified to the appropriate Regulatory Agency.

Megan Vernak
Director, Product Monitoring – Specialty Brands