Medline International France SAS
5 rue Charles Lindbergh
44110 Châteaubriant - France
+33 (0)2 40 81 57 62
Fax +33 (0)2 40 81 56 34
www.medline.com/fr

URGENT: FIELD SAFETY NOTICE
Medical Device Safety Advisory Notice

Châteaubriant, 10th October 2018

For the attention of: the Pharmacist responsible for medical device vigilance and the Biomedical Engineering Department.

SECURITY INFORMATION of Medline Sterile Procedure Trays

Medline reference: **FSN-18/12**

Description: Medline Sterile Procedure Trays containing a Kocher clamp.

Product Codes and lots concerned: See details in the **table** of the Acknowledgment receipt

Dear Customer,

This letter is to inform you of an urgent safety notification concerning Medline sterile procedure trays including "kocher" clamps. Its objective is to remind customers on the recommendations of use and to inform you that there were no information for use included in the trays manufactured from June 2018 to October 2018.

Recommendations for Use:

The "kocher" clamps used for maternity procedures of newborns require visual inspection during use if they are used to clamp the umbilical cords (see Appendix – Operating instructions).
Other devices such as Barr clamps have a locking system that does not require monitoring

Actions required:

- Thank you for communicating the information for use to all user services of the sterile procedure trays listed in the table of the acknowledgement form BELOW.
- Thank you for providing users with the information for use (attached), missing from the kits manufactured between 2018-06 and 2018-10

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Instruction for Use will be included in the new productions.

Please return the acknowledgement form by fax or by email no later than October 31, 2018

Yours Sincerely,
Quality and Regulatory Affairs Dept.

Medline International France SAS
5 rue Charles Lindbergh
44110 Châteaubriant - France
+33 (0)2 40 81 57 62
Fax +33 (0)2 40 81 56 34
www.medline.com/fr

**Acknowledgement receipt to be faxed to the following Fax: XXXXXXXXXX or
E-mail: XXXXXXXXXX**

Please complete this form and return to Medline by either fax or email as soon as possible, but **not later than October 31st, 2018**.

Reference: FSN-18/12

Table:

For your facility, the references of Sterile Procedure Trays concerned by this notification are:

Reference
XXXXXXXX
XXXXXXXX
XXXXXXXX

I have read and understand the instructions provided by Medline and acknowledge receipt of the FSN-18/12 regarding the Medline Sterile Procedure Trays by signing below.

I also agree to further distribute and communicate this important information within my facility as required.

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