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23 November 2018

Ansprechpartner/in

## IMPORTANT SAFETY INFORMATION FOR USERS

### *Notification of a product recall*

**Affected product:** MP Prosthesis stem

**Item-Number:** 172-916/12

**Lot-Numbers:** All after manufacturing date: 2013-10 

**LINK® Reference-Number:** R-2018-03

Dear Sir or Madame,

a risk has been identified for the use of the MP prosthesis stem. For safety reasons WALDEMAR LINK GmbH & Co. KG will recall these products. It is noticed that the smallest available stem size is not compatible with two corresponding instruments.

As you were supplied with the affected products, please read carefully the information on the following pages.

Please return the fax reply to us in any case until **12 December 2018**, even if you have none of the listed products in stock or if these products do not exhibit the defect in question.


In conclusion, we apologize for any inconveniences caused and thank you for your understanding in this matter.

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For the smallest available stem size (Art. No. 172-916/12) the reaming guide and the guide rod can't be screwed into the intended depth. The reason for this is a reduced depth of the core hole in the stem.

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A reaming guide that is not screwed in completely may lead to a proximal preparation in the depth that is not sufficient and to a neck that can't be positioned correctly.

Corrective reaming can lead to a surgery time extension.

The correct fit of the neck can't be verified by means of the guide rod.

This failure can only occur with the smallest available stem size (Art. No. 172-916/12).

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Safety information and recall of the affected articles from the market

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-----

- Should you have any of the Prosthesis stems 172-916/12 in your inventory, please send the products back to Waldemar Link GmbH & Co. KG or confirm implantation.
- An exchange will not take place because this article is taken out of the product portfolio temporarily.
- Affected of this recall are all lot and serial numbers after manufacturing date: 10-2013.
- Please return the fax reply to us in any event until the **12 December 2018** as documentation of the recall. This applies even if you have none of the listed products in stock or if these products do not exhibit the defect in question. Please ensure that all users of the above products within your organization and other relevant persons have been notified of this **safety notice**. If you have transferred the products to third parties, please pass on a copy of this information or notify the contact person indicated below.

The responsible Competent Authorities have been notified of this Field Safety Corrective Action.

**Kontakte / Contacts Waldemar Link GmbH & Co. KG:****Reklamationsmanagement / Complaint Management**

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