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Closing Confirmation Sheet

DRI OCT Triton, DRI OCT Triton plus Software, Ver.1.1.4

Customer name	Installation date / / /	
Address	Agency/Supplier	Business office
	The body model: The body S/N:	
TEL () -		

● Version Up check procedure of a software

#	Procedure item	Check point	PASS	FAIL															
1	It be Ver1.1.3 in front of Version Up. (Select "information" from the "SETUP" menu.) <div style="border: 1px solid black; padding: 5px; width: fit-content;"> <table style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 15%;">Basic</td> <td>Serial No</td> <td></td> </tr> <tr> <td></td> <td>Region Code</td> <td></td> </tr> <tr> <td></td> <td>System</td> <td style="border: 2px solid red; text-align: center;">1.1.3</td> </tr> <tr> <td></td> <td>OCT Capture Count</td> <td></td> </tr> <tr> <td></td> <td>Fundus Capture Count</td> <td></td> </tr> </table> </div>	Basic	Serial No			Region Code			System	1.1.3		OCT Capture Count			Fundus Capture Count		It's confirmed that it's a repair target.	<input type="checkbox"/>	<input type="checkbox"/>
Basic	Serial No																		
	Region Code																		
	System	1.1.3																	
	OCT Capture Count																		
	Fundus Capture Count																		
2	Version Up, error confirmation	An error message don't be indicated.	<input type="checkbox"/>	<input type="checkbox"/>															
3	Version behind version Up be Ver.1.1.4. (Select "information" from the "SETUP" menu.) <div style="border: 1px solid black; padding: 5px; width: fit-content;"> <table style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 15%;">Basic</td> <td>Serial No</td> <td></td> </tr> <tr> <td></td> <td>Region Code</td> <td></td> </tr> <tr> <td></td> <td>System</td> <td style="border: 2px solid red; text-align: center;">1.1.4</td> </tr> <tr> <td></td> <td>OCT Capture Count</td> <td></td> </tr> <tr> <td></td> <td>Fundus Capture Count</td> <td></td> </tr> </table> </div>	Basic	Serial No			Region Code			System	1.1.4		OCT Capture Count			Fundus Capture Count		Make sure that the version information is as follows. <div style="border: 2px solid red; display: inline-block; padding: 2px;"> System 1.1.4 </div>	<input type="checkbox"/>	<input type="checkbox"/>
Basic	Serial No																		
	Region Code																		
	System	1.1.4																	
	OCT Capture Count																		
	Fundus Capture Count																		
4	Checking existing data	Check existing data by the tool software "PowerSaveLRChecker". If wrong data exist make necessary correction. <table style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 60%;">Nuber of checked images</td> <td style="width: 40%;"></td> </tr> <tr> <td>Hours taken for the check</td> <td></td> </tr> </table>	Nuber of checked images		Hours taken for the check		<input type="checkbox"/>	<input type="checkbox"/>											
Nuber of checked images																			
Hours taken for the check																			

Comment

Name	Signature	Initials	Date
TESTING PERFORMED BY			
CUSTOMER			