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## SureSigns VS3/VS4 Vital Signs Monitors Pulse Rate Software Installation Instructions

### Contents

<b>Contents</b> .....	1
<b>Overview</b> .....	1
<b>Philips InCenter</b> .....	2
<b>Prior to Software Installation</b> .....	2
Retrieving Documentation.....	2
Checking System Serial Number.....	2
Inspecting the Battery .....	4
Retrieving Software.....	4
Clearing Patient Data.....	4
Other Recommendations.....	7
<b>Upgrading Software</b> .....	7
<b>Verification of Software Installation</b> .....	10
Power-on Self Test.....	10
System Settings Test.....	11
Battery Management Test.....	11
System Serial Number & Software Version Test.....	12
<b>Installation Completion</b> .....	13
Email Address.....	13

### Overview

These Installation Instructions are intended for biomedical engineers, technicians, or personnel responsible for troubleshooting, repairing, and maintaining Philips patient monitoring systems.

These Installation Instructions are strictly to update the device with the software to monitor the health of the battery. These instructions require access to password-protected functions on the device.

These Installation Instructions are to be used in conjunction with the monitor's Service Guide for the following devices with software releases up to and including A.07.32.

Product	Description
863069	SureSigns VS3 NBP
863070	SureSigns VS3 NBP, Temp
863071	SureSigns VS3 NBP, SpO <sub>2</sub>

863072	SureSigns VS3 NBP, SpO <sub>2</sub> , Rec
863073	SureSigns VS3 NBP, SpO <sub>2</sub> , Temp
863074	SureSigns VS3 NBP, SpO <sub>2</sub> , Temp, Rec
863283	SureSigns VS4 NBP, SpO <sub>2</sub>
863286	SureSigns VS4 Government Bundle

## Philips InCenter

You will need access to the Philips InCenter, which requires an active registration and password.

Before you register, obtain the serial number of the monitor you plan to include in your registration. The serial number is located on the product identification label on the rear of case of the monitor.

To register:

Step	
1.	Access the Philips InCenter website at <b><a href="http://incenter.medical.philips.com">incenter.medical.philips.com</a></b> .
2.	Click <b>Need Help?</b> .
3.	Under <b>Software Updates</b> , Click <b>Click here for SureSigns patient monitor account registration</b> . The <b>SureSigns InCenter Registration Form</b> appears.
4.	Enter your personal information and answer the questions, and then click <b>Submit</b> .  After your information is processed, an email with temporary login information is sent to the email address you entered. It may take 3 days to receive the email.

## Prior to Software Installation

### Retrieving Documentation




The following documentations should be readily available when performing the installation. These documents can be retrieved from Philips InCenter.

- Installation and Configuration Guide (453564649971/453564650041)
- Service Guide (453564649981/453564650051) and Addendum (453564865971)
- Instruction for Use (453564303651/453564643641) and Addendum (453564862721)

### Checking System Serial Number

The serial number displayed on the monitor must match the serial number on the back of the monitor label.

You can verify the monitor's serial number by following the steps below.

Step	
1.	<p>From the main screen, press the <b>System Menu</b> key  .</p> <p>The <b>System Menu</b> appears. See figure 1.</p>  <p><b>Figure 1 System Menu</b></p>
2.	Locate the <b>Serial Number</b> (halfway down the <b>System Menu</b> Screen).
3.	<p>Locate the serial number on the back of the monitor.</p> <p>The serial number appears in a label (marked <b>SN</b> ) on the back of the monitor. See figure 2.</p>  <p><b>Figure 2 Back View</b></p>
4.	If the serial number matches, proceed to the <i>Inspecting the Battery</i> section.
5.	<p>If the serial number does not match, the serial number will need to be reset per the steps outlined in the Service Guide Section, section “Resetting the serial number”.</p> <p>The Service Guide can be downloaded from Philips InCenter.</p>

---

## Inspecting the Battery

It is very important that a Battery Maintenance, as described in Service Guide Addendum is performed. This document can be downloaded from Philips InCenter.

Performing Battery Maintenance requires access to the **Battery Info** screen. (password required)

---

**NOTE** - A replacement battery may need to be ordered and available before continuing with the software installation.

---

The battery in the device must meet the expected SOH and cycle count (charge/discharge).

## Retrieving Software

The system software will need to be downloaded from the Philips InCenter to a USB Flash Drive.

InCenter Login: [incenter.medical.philips.com](http://incenter.medical.philips.com)

Location of Software:

Service > Software > Software Downloads > Patient Monitoring > Bedside Monitors > SureSigns > SureSigns VS3 or SureSigns VS4 > A.07.33

---

## NOTES –

- Philips recommends using a SanDisk® or Kingston® USB flash drive for software upgrades.
  - The software upgrade folder must be located in the top directory of the USB flash drive.
- 

## Clearing Patient Data

All patient data **MUST** be transferred to your EMR before performing the following steps.

---

**WARNING** - These steps will delete all patient data on the monitor, make sure you have transferred the data to your EMR.

---

Two methods are available for deleting patient data:

- Clearing data by using the **Clear Data** option (password required)
- Deleting patient records from the **Display Mode** button on the main screen

To clear the data by using the **Clear Data** option:



---


**NOTE** - This method will clear all patient data in memory. It does not clear calibration and hardware configuration data.

---

Step	
1.	From the main screen, press the <b>System Menu</b> key  . The <b>System Menu</b> appears. See figure 1.
2.	Select the <b>System Admin</b> button, a <b>Password</b> window appears.

<p>3.</p>	<p>Enter the Administrator password. Click on the <b>OK</b> button.</p> <p><b>NOTE</b> – <i>The default Administor password is 2-1-5. For VS4, it is recommended to change the default Administrator password to a password defined by your institution.</i></p> <p>The <b>System Admin Menu</b> appears. See figure 3.</p>  <p><b>Figure 3 System Admin Menu</b></p>
<p>4.</p>	<p>Select the <b>Diagnostics</b> button.</p> <p>The <b>System Diagnostics</b> menu appears. See figure 4.</p>  <p><b>Figure 4 System Diagnostics</b></p>
<p>5.</p>	<p>Select the <b>Maintenance&gt;&gt;</b> button. A <b>Password</b> window appears.</p>
<p>6.</p>	<p>Enter the Maintenance password (1-2-9). Click on the <b>OK</b> button.</p> <p>The complete <b>System Diagnostics</b> menu appears, including all <b>Maintenance</b> functions. See figure 5.</p>


**Figure 5 System Diagnostics with Maintenance functions**

7.	<p>Select the <b>Clear Data</b> Button.</p> <p>A confirmation window appears. <b>Clear Data and restart system?</b></p> <p>Select <b>Yes</b>.</p> <p>The monitor clears the patient data and error log, and then restarts.</p>
----	--

To delete patient records from the **Display Mode** button on the main screen:

**NOTE** - This method will clear all patient records in memory. It does not clear calibration and hardware configuration data.

Step	
1.	From the main screen, select the <b>Display Mode</b> button.
2.	<p>Select <b>Display Patient Records</b>.</p>  <p>A list of <b>Patient Records</b> appears.</p>
3.	Select the <b>Delete</b> button.
4.	<p>Select <b>Delete All Records</b>.</p> <p>A confirmation window appears. <b>Are you sure?</b></p> <p>Select <b>Yes</b>.</p> <p>The monitor clears the patient data.</p>

## Other Recommendations

Before you upgrade the software, it is recommended that you back up the system settings by exporting the current configuration settings or by recording them on the worksheets provided in the *Installation and Configuration Guide* for your monitor.

When you upgrade the software:

- Charge the battery before upgrading the software.
- Never perform a software upgrade with the monitor connected to a patient.
- Disconnect any USB peripherals.
- Do not upgrade software through a USB hub.
- If the USB port has a clamp in place, remove the clamp to ensure that the USB flash drive can be inserted completely into the USB port.

After the upgrade starts:

- Do not unplug the monitor.
- Do not remove the USB flash drive.
- Do not press any key.

## Upgrading Software

The following are the steps to upgrade the system software with a USB flash drive.

When you upgrade the software, the monitor saves all of the system settings, including the System Admin Menu settings, and the Error log.

Before you place an upgraded monitor back into service, Philips recommends that your facility establish a protocol to train users about any changes in the functionality of the monitor resulting from the software upgrade.

---

**NOTE** – You will need the Administrator password to perform this task.



---

To perform a software upgrade:

Step	
1.	Connect the monitor to AC power and power up the monitor. <b>NOTE</b> - <i>Your monitor must be connected to AC power and have a fully charged battery before upgrading the software.</i>
2.	Verify <b>Serial Number</b> displayed on the device matches the serial number on the back of the device label. Refer to <i>Checking System Serial Number</i> section for complete details.
3.	Verify Battery Maintenance has been performed. Refer to <i>Inspecting the Battery</i> section for complete details.
4.	Verify Patient Data has been deleted. Refer to <i>Clearing Patient Data</i> section for complete details. <b>NOTE</b> - <i>All Patient data must be deleted before performing the upgrade.</i>
5.	Download the system software from Philips InCenter onto a USB flash drive. <b>NOTES</b> –

	<p>The software upgrade folder must be located in the top directory of the USB flash drive.</p> <p>The USB flash drive must not be write-protected, as the upgrade code will need to be exported to the USB after the upgrade has been completed.</p> <p>Refer to <i>Retrieving software</i> section for complete details.</p>
6.	<p>Insert the USB flash drive with the software upgrade into the USB port on the back of the monitor.</p> <p><b>NOTE</b> - The USB flash drive must not be write-protected, as the upgrade code will need to be exported to the USB after the upgrade has been completed.</p>
7.	<p>From the main screen, select <b>System Menu</b> key  .</p> <p>The <b>System Menu</b> appears. See figure 1.</p>
8.	<p>Select the <b>System Admin</b> button, a <b>Password</b> window appears.</p>
9.	<p>Enter the Administrator password. Click on the <b>OK</b> button.</p> <p><b>NOTE</b> – The default Administrator password is 2-1-5. For VS4, it is recommended to change the default Administrator password to a password defined by your institution.</p> <p>The <b>System Admin Menu</b> appears. See figure 3.</p>
10.	<p>Select <b>Upgrade Software</b>.</p>
11.	<p>The monitor searches for a valid software image on the USB flash drive, and then displays the software image information in the <b>Upgrade Software</b> window. See figure 6.</p> <div data-bbox="418 1171 1247 1749" data-label="Image"> </div> <p><b>Figure 6 Upgrade Software</b></p> <p><b>NOTE</b> - If the monitor cannot find a valid software image, the screen returns to the <b>System Admin Menu</b>.</p>

	<p><b>Caution: Do not downgrade the software to an earlier version, because it may cause hardware incompatibility and loss of system settings and patient records.</b></p>
12.	<p>Select the <b>Upgrade</b> button to start the upgrade.</p> <p><b>Caution:</b></p> <p><b>After the upgrade starts:</b></p> <ul style="list-style-type: none"> <li>• Do not unplug the monitor.</li> <li>• Do not remove the USB flash drive.</li> <li>• Do not press any keys.</li> </ul> <p><b>If the upgrade is inadvertently interrupted and the main board data is lost, the main board will need to be replaced. Refer to the Service Guide for more information.</b></p>
13.	<p>The <b>Upgrade in Progress</b> indicator increments while the upgrade is in progress.</p> <p>Keep the USB flash drive in the monitor until all steps are completed.</p> <p>When the upgrade is complete, the <b>Checking Memory CRC</b> and <b>Upgrade Successful</b> messages appear.</p> <p>After the software upgrade, the monitor automatically shuts down and restarts.</p> <p><b>NOTE</b> - <i>If the upgrade fails for any reason the monitor will initiate an error message and stop the upgrade. Close the upgrade software menu. The monitor will continue running with the original revision software.</i></p>
14.	<p>After the monitor restarts, the main screen is displayed.</p>
15.	<p>From the main screen, select <b>System Menu</b> key  .</p> <p>The <b>System Menu</b> appears. See figure 1.</p>
16.	<p>Locate the <b>Export Upgrade Code</b>, at the bottom of the screen.</p> <p>The <b>Export Upgrade Code</b> button will create or update the <i>UpgradeCode.csv</i> file on the USB flash drive with the Serial number and Upgrade Code of the device.</p> <p><b>NOTE</b> - <i>When the same USB flash drive is used to update several units, the <b>UpgradeCode.csv</b> file is also updated with each successful Export Upgrade Code.</i></p> <p><i>For example, if the USB flash drive was used to update 5 devices, the UpgradeCode.csv file will contain the confirmation information for each of the 5 devices.</i></p>
17.	<p>After the system has successfully exported the upgrade code, a message will appear and you can remove the USB Flash Drive.</p>

	 <p>If no USB Flash drive, the device is unable to export the file, the following message appears.</p>  <p><b>NOTE</b> - Please contact your local Service Response Center for additional assistance if unable to export the upgrade code.</p>
18.	Select <b>OK</b> to return to the main screen.
19.	<p>Remove the USB Flash Drive from the back of the monitor.</p> <p><b>NOTES –</b></p> <p><i>The USB Flash Drive will be needed for the Installation Completion steps. Before placing the monitor back in service, it is recommended to perform the following steps detailed in the Verification of Software Installation:</i></p> <ul style="list-style-type: none"> <li>• <i>Power-on Self test</i></li> <li>• <i>System settings test</i></li> <li>• <i>Battery Management test</i></li> <li>• <i>System serial number and upgrade code test</i></li> </ul>

## Verification of Software Installation

Philips requires that you complete the following tests after upgrading the software.

### Power-on Self Test

Step	
1.	Connect the monitor to an AC power source.
2.	Press the <b>On/Standby</b> key to power up the monitor.
3.	<p>Make sure that the monitor powers up successfully as described in the following sequence:</p> <ul style="list-style-type: none"> <li>• The screen displays color bars for about five seconds.</li> <li>• The LCD turns off for three seconds, and the charging LED lights.</li> <li>• The Philips splash screen appears for one second, and a startup tone sounds.</li> <li>• It may be possible, that the Date/Time Menu appears for a new time setting. Enter Date &amp; Time.</li> </ul>

	<ul style="list-style-type: none"> <li>The main screen appears.</li> </ul>
<b>Verification Results</b>	
Pass	Expected result is that the monitor starts up and displays the main screen.
Fail	<p>Refer to the device Service Guide for details, if any of the below occurs:</p> <ul style="list-style-type: none"> <li>If the start-up and power up sequences are not as expected.</li> <li>If the LEDs do not function as expected.</li> <li>If the display does not function as expected.</li> <li>If you do not hear a startup tone, or the monitor displays the <b>Speaker Malfunc</b> error message.</li> </ul> <p><b>NOTE</b> - Please contact your local Service Response Center for additional assistance if unable to obtain expected results after the software upgrade.</p>

### System Settings Test

<b>Step</b>	
1.	From the main screen, verify that your system settings are preserved and as expected.
<b>Verification Results</b>	
Pass	Expected result is that the system settings are preserved and as expected.
Fail	<p>Refer to the device Service Guide for details, if any of the below occurs:</p> <ul style="list-style-type: none"> <li>If the system settings are not preserved and as expected.</li> </ul> <p><b>NOTE</b> - Please contact your local Service Response Center for additional assistance if unable to obtain expected results after the software upgrade.</p>

### Battery Management Test

<b>Step</b>	
1.	<p>From the main screen, verify the bottom left hand corner does not display any battery warning or alarm messages.</p> <p>If a battery warning or alarm is displayed, refer to the provided Service Guide addendum for actions to resolve the battery failure.</p>
2.	<p>From the main screen, elect the Battery Status Icon .</p> <p>Verify the following displayed values are not in RED.</p> <ul style="list-style-type: none"> <li><b>State of Health (SOH)</b></li> <li><b>Cycle Count</b></li> </ul> <p>If the above parameters and values are in RED, refer to the Service Guide addendum for actions to resolve the displayed RED values.</p>

Verification Results	
Pass	<p>Expected results are:</p> <ul style="list-style-type: none"> <li>• Monitor main screen displays <b>no</b> battery error code, battery warning or alarm messages.</li> <li>• Battery Status displays <b>no</b> SOH or Cycle Count parameters or values in RED.</li> </ul>
Fail	<p>Refer to the device Service Guide for details if any of the below occurs:</p> <ul style="list-style-type: none"> <li>• Battery error codes, Battery warning or alarms are displayed on the lower left hand corner.</li> <li>• Battery status displays SOH or Cycle Count parameter and value in RED.</li> </ul> <p><b>NOTE</b> - Please contact your local Service Response Center for additional assistance if unable to obtain expected results after the software upgrade.</p>

### System Serial Number & Software Version Test

Step	
1.	<p>From the main screen, select <b>System Menu</b> key  .</p> <p>The <b>System Menu</b> appears. See figure 1.</p>
2.	<p>Locate the <b>Serial Number</b> (halfway down the <b>System Menu</b> Screen). Verify the Serial Number matches the Serial Number on the back of the monitor.</p>
3.	<p>Locate the Software Version (just below the Serial Number). Verify the software version is A.07.33.</p>
Verification Results	
Pass	<p>Expected results are:</p> <ul style="list-style-type: none"> <li>• Serial Number matches Serial Number on the back of the monitor.</li> <li>• Software Version matches the Software Version A.07.33.</li> </ul>
Fail	<p>Refer to the device Service Guide for details if any of the below occurs:</p> <ul style="list-style-type: none"> <li>• Serial Number <b>does not</b> match Serial Number on the back of the monitor.</li> <li>• Software Version <b>does not</b> match the Software Version A.07.33.</li> </ul> <p><b>NOTE</b> - Please contact your local Service Response Center for additional assistance if unable to obtain expected results after the software upgrade.</p>

## Installation Completion

After completing the above verifications, please provide confirmation that the monitors in your institution have been updated with the Battery Management Software.

The confirmation for each upgraded monitor should have been successfully exported to the *UpgradeCode.csv* file on the USB Flash drive. The .csv file contains the system serial number and upgrade code for each upgraded monitor.

Philips has provided an easy method for sending the device confirmation back to Philips.

### Email Address

Step	
1.	Insert the USB Flash Drive used to perform the software upgrade into the USB port on the PC.
2.	Attach the <b><i>UpgradeCode.csv</i></b> file to an email.
3.	Send the email to the following address: <b><u>recall.response2@philips.com</u></b>
4.	An email will be sent confirming the receipt of the email. <b>NOTE</b> - <i>You may want to keep a copy of the email and UpgradeCode.csv file, in case you are contacted for clarification at a later date.</i>

Philips Medical Systems  
3000 Minuteman Road  
Andover, MA 01810 USA  
©Copyright 2018 Koninklijke Philips N.V.  
All rights are reserved.

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This Addendum contains updated information for the SureSigns VS3, VS4 Release A.07 Service Guide. Please store this Addendum with your SureSigns VS3, VS4 documentation for future reference.

The topic of this Addendum is about changes in the Battery replacement, test, software upgrade and error codes. To optimize the safe use of your monitor and its performance, adherence to the instructions included in this addendum is required.

### Maintaining the Battery

#### About the Battery

The rechargeable lithium ion battery used in the VS series vital signs monitors is a smart battery with built-in circuitry that communicates battery status information to the monitor. Battery power lasts a minimum of four (4) hours of continuous monitoring with no printing and one NBP measurement every 15 minutes.

To properly maintain the battery and prevent damage to the monitor, observe these guidelines:

- If a battery shows damage or signs of leakage, replace it **immediately**.
- **Never** use a faulty battery in a monitor.
- **Never** dispose of the battery in a normal waste container.
- **Never** leave a battery inside the monitor if the monitor is not being used for a long period of time.
- **Never** store a battery that is charged to more than 50% capacity.

---

#### NOTES -


- Images shown are from a VS4 vital signs monitor. The VS3 may appear slightly different.
  - For information about the battery status indicator, please see the Instructions for Use provided with your monitor.
- 

### Viewing Battery Information

As a battery ages, its capacity decreases, and the battery status indicator becomes increasingly less accurate as the total number of charges and discharges increase. You can view the **Battery Info** screen to display information about the battery. Ensure that the monitor is connected to AC power before attempting to review battery information.

**NOTE** - If a VS3 monitor with a **Hardware ID** (displayed on the **System Menu**) of **1-1-A1** is disconnected from AC power, the **Battery Info** window displays the following message: **No data from battery. Please see Service Guide.** To view the information about the battery, reconnect the monitor to AC power.

**To view information about the battery:**

Step	
1.	On the main screen, touch <b>System</b> .  <b>NOTE</b> - <i>The VS4 vital signs monitor utilizes a touch screen. On a VS3 monitor, you must rotate the control wheel to highlight the desired menu option, then press the wheel to select it.</i>
2.	Touch <b>System Admin</b> and enter the administrator password (default is 215). The System Admin Menu appears.
3.	Touch <b>Diagnostics</b> to open the System Diagnostics menu.   <p>Figure 1 System Admin Menu</p>
4.	Touch <b>Battery Info</b> to open the Battery Info screen.

System Diagnostics

Monitoring Suspended

NBP Cycle Count: 8  
LCD Usage Hours: 166.9  
Errors: 500      Error Log

Self Test      Display Test  
Recorder Test      Audio Test  
**Battery Info**      LED Test  
Keys Test      Touch Screen Calibration

Network Test  
Maintenance >>

Return

Figure 2 System Diagnostics

5. To view the list of results, select the list to activate scrolling.

Battery Info

Parameter	Value
SOH	94% (Full Capacity/Design Capacity)
Manufacturer	EONEMOLI
Chemistry	LIION
Serial Number	#35592
Manufact. Date	04/23/2017 (mm/dd/yyyy)
Cycle Count	1
Max Error	2%
Relative Charge	74%
Absolute Charge	76%
Battery Name	E013R3
Voltage	11404 (mV)
Current	-795 (mA)
Temperature	300.60 (K) 27.60(C)
Full Capacity	7403 (mAh)
Design Voltage	11100 (mV)

Recondition      Return

Figure 3 Battery Info Screen

The **Battery Info** screen provides detailed information about battery capacity and charging status, including:

- **SOH:** State of Health, the Full Capacity/Design Capacity.
- **Cycle Count:** The number of full charge and discharge cycles calculated by the battery.

- **Max Error:** The expected margin of error in the state of the charge calculation. The Max Error value is the difference between the Relative Charge value and the Absolute Charge value.
- **Relative Charge:** The predicted remaining battery capacity, expressed as a percentage of Full Capacity. The value in the Relative Charge decreases as the battery ages. The Battery Status icon is a graphic representation of the relative charge.
- **Absolute Charge:** The predicted remaining battery capacity, expressed as a percentage of Design Capacity.
- **Full Capacity:** The predicted capacity of the battery when it is fully charged. The value in the Full Capacity field decreases as the battery ages. The difference between the value in the Full Capacity field and the value in the Design Capacity field is an indication of battery condition.
- **Design Capacity:** The capacity of a new battery.

If the battery **SOH** is  $\leq 80\%$ , or battery cycle count exceeds the recommended limit of **300 cycles**, the battery will need to be replaced (see Figure 3).

---

**NOTE** - If the message, “No data from battery. Please see Service Guide.” appears, you must reset the battery. Refer to your SureSigns VS3, VS4 Service Guide for more information.

---

---

## WARNING

The risk of battery failure increases when a battery remains in use past 300 charge-discharge cycles or when the SOH is  $\leq 80\%$ . Such failures can result in overheating that in rare cases can cause the battery to ignite.

---

## Reconditioning the Battery

Reconditioning the battery reduces the Max Error value, in turn, increases the accuracy of the Relative Charge. Philips recommends that you condition the battery by fully discharging and recharging it when the Max Error is 10% or greater.

To recondition the battery:

Step	
1.	Open the Battery Info window (see <b>Viewing Battery Information</b> ).
2.	Disconnect the monitor from AC power.
3.	Rotate the wheel to select <b>Recondition</b> , then press the wheel. The <b>Relative Charge</b> percentage will decrease to 0%.

4.	When the monitor shuts down, reconnect the monitor to AC power and allow the battery to recharge to 100%.
5.	Repeat steps 1 - 4. <b>NOTE</b> - <i>If the battery does not recharge after four reconditioning cycles, replace it.</i>

## Replacing the Battery

Replace the battery if the following conditions occur:

- After reconditioning, if the monitor operates for less than one hour on a fully charged battery before the low battery (**Low Batt**) alarm occurs, or
- The **Max Error** cannot be brought  $\leq$  8% after several recondition cycles.
- The battery **SOH** is  $\leq$  80%.
- The battery **Cycle Count** is greater than 300<sup>1</sup>.
- The monitor reports “Replace Battery” alert.
- The monitor reports “Replace Battery Immediately” alert.


---

### WARNING

**Dispose of aged batteries in an environmentally responsible manner. Do not dispose of the battery in normal waste containers. Consult your hospital administrator to find out about local arrangements.**


---

To replace the battery:

Step	
1.	Shut down the monitor.
2.	Disconnect the AC power cord from the rear of the monitor.
3.	Tip the monitor and insert a flathead screwdriver into the slot in the bottom of the case under the battery cover. 

---

<sup>1</sup> This applies to the Lithium-ion, Smart Battery. The only supported batteries in the SureSigns VS Series Monitors are the Lithium-ion, Smart Battery (Part Number 989803194541 or 989803144631).

	<p><b>CAUTION</b></p> <p><b>Do not try to remove the battery cover by inserting the screwdriver into the vents in the cover. This can damage the battery case cover.</b></p>
4.	Twist the screwdriver slightly to pop the battery cover off of the case.
5.	<p>Press down on the battery retaining clips and pull the old battery out of the monitor using the ribbon.</p> 
6.	Orient the replacement battery so that the contacts and ribbon are on the right.
7.	Insert the replacement battery into the empty compartment, and push until the battery snaps into place behind the metal retaining clips.
8.	Replace the battery cover by inserting it into the battery opening with the two tabs facing up.
9.	Snap the cover into place by pressing it firmly against the monitor case.

## Battery Test

Philips requires that you complete Battery Test after:

- Upgrading the software
- Opening the monitor for any reason
- Replacing any internal part (except NBP parts, SpO<sub>2</sub> board)

**To test the battery:**

Step	
1.	<p>Make sure there are no battery error codes, battery warnings or alarm messages in the lower left corner of the main menu.</p> <p><b>NOTE</b> - <i>If a battery error code, warning or alarm message is displayed, refer to the Error Codes below for resolution.</i></p>
2.	<p>Select the Battery Status Icon . Make sure there are no battery parameters in RED.</p> <p><b>NOTE</b> - <i>If a battery value is RED, refer to the Error Codes below for resolution.</i></p>

## Ordering a Replacement Battery

To order a replacement lithium ion battery, contact your Philips Customer Care center. For more information on how to contact Philips Customer Care for your country, go to <http://www.healthcare.philips.com>. Select your country and language, then navigate to the Customer Care page.

Product	Description	Battery PN*	Battery PN**
863069	SureSigns VS3 NBP	989803194541 (11.1V 7800 mAh, ME202EK)	989803144631 (11.1V 7200 mAh, ME202CJ)
863070	SureSigns VS3 NBP/Temp		
863071	SureSigns VS3 NBP/SPO2		
863072	SureSigns VS3 NBP/SPO2/Rec		
863073	SureSigns VS3 NBP/SPO2/Temp		
863074	SureSigns VS3 NBP/SPO2/Temp/Rec		
863283	SureSigns VS4 NBP/SPO2		
863286	SureSigns VS4 Government Bundle		

\* World-wide, except China

\*\* China only

## Battery Messages and Alarms

The condition of the battery is reported by technical alarms and error codes.

### Technical Alarms

The following battery technical alarms appear in the monitor's message area:

- **Low Batt:** remaining battery power is less than 30%.
- **Extreme Low Batt:** remaining battery power is less than 21%.
- **Unknown Battery, slow chrg** or **Unrecognized Battery, slow charge:** The battery plugged in is not supported.
- **Battery Charge Failure:** Battery charge power failure or battery malfunction.
- **Battery Overcharged:** Battery is overcharged.
- **HI Battery Temp, chrg stop** or **High Battery Temp, charge stopped:** The battery temperature is  $\geq 50^{\circ}\text{C}$  or the battery issues an "Over Temperature" status.
- **Batt Overhear, shut down** or **Battery Overheat, shutting down:** Battery Temperature is  $\geq 60^{\circ}\text{C}$  or the battery issues an "Over Temperature" status.

The battery will stop charging immediately and the monitor will shut down automatically after you confirm the error message or if the error message is displayed for 25 seconds.

- **Replace Battery:** The battery cycle count is > 300 or the battery SOH is  $\leq$  80%.
- **Replace Battery Immediately:** The battery cycle count is > 400 or the battery SOH is  $\leq$  75%.

## Error Codes

An error code (for example, **257 System Error**, indicating a battery charger power failure) appears in the Error Log. To view the Error Log, see “Viewing, Printing, and Exporting the Error Log” in the Service Guide. For a complete list of error codes and the actions to take, see Chapter “Troubleshooting,” in the Service Guide.

The following error codes have been added.

Code	Description/Cause	Error Message	Action
257	Battery Charger Power Failure	System Error	Replace battery. If the problem persists, replace the main board. If the problem persists, replace the AC power module.
261	Battery Over Charge	System Error	Run monitor without AC power. If the problem persists, replace the battery. If the problem persists, replace the main board.
264 <sub>temp</sub> value	Battery Over Temperature	System Error	Perform hard shutdown. If the problem persists, replace the battery. If the problem persists, replace the main board.
265	Battery cycle count is > 300 or Battery SOH is $\leq$ 80%.	System Error	Replace battery.
265	Battery cycle count is > 400 or Battery SOH is $\leq$ 75%.	System Error	Replace battery.
266	Unrecognized Battery	System Error	Replace battery.

Philips Medical Systems  
3000 Minuteman Road  
Andover, MA 01810 USA  
©Copyright 2019 Koninklijke Philips N.V.  
All rights are reserved.

453564865971 Revision A