

Healthcare facility
Address

To the attention of the vigilance Safety Officer
and orthopedic surgery departments


Valence, June 27th 2019,

Ref. AMPLITUDE: COMP-1061 (EN)

Object: **BATCH RECALL**
Instrument – Broach handle - AMPLITUDE

Reason for recall

Following feedbacks from healthcare facilities, a risk of breakage of Nitinol cables assembled on curved broach handles used for surgeries of Amplitude hip prostheses has been identified. The broach handles on the market exist in 2 versions: with Nitinol cables or braided cables (see below pictures). This recall concerns only the broach handles with Nitinol cables.



Picture of concerned broach handle



Picture of braided cable
Not concerned by this recall




Picture of Nitinol cable (smooth cable)
Concerned by this recall

Circumstances and risks for the user and/or the patient

No patient risk associated to this breakage was identified. In case of breakage of the cable, the broach will remain attached to the handle. The consequence could be an increase in surgical time to disassembly the broach from handle. This disassembly can be done by moving the mobile axis. The use of headless pin available in the instrument set (reference 2-0114000) is recommended to pull the mobile axis and therefore release the broach.

Headless pin 2-0111400

Hole of mobile axis



Concerned devices

The traceability data indicate that you have received devices of one or several of the involved batch numbers:

Reference	Designation	Batch numbers
2-0103100	Broach handle	1600624
		1800012
		1800013
		1800014
		1800015
		1800780
		1800781
		1800782
		1800783

The batches are partially involved (only broaches with Nitinol cables are concerned by this recall)



What you must do

Please circulate this notice to the appropriate individuals.

Your local representative will contact you to organize the replacement of the devices.

We remind you that any adverse event experienced using these devices must be declared to the competent authority and your local representative.

Other information

Swissmedic and the French competent authority have been advised about this recall.

We apologize for inconvenience and thank you for your comprehension.

Mireille LEMERY

Director, Quality and Regulatory Affairs – vigilance correspondant

vigilance@amplitude-ortho.com