**Follow up Urgent Field Safety Notice  
Resolution**

**IMC19-07.D.OUS**

**October 2022**

**IMMULITE®  
IMMULITE® 1000  
IMMULITE® 2000  
IMMULITE® 2000 XPi**

**Resolution of the IMMULITE Systems Low Discordant Progesterone Results on Some Patient Samples**

---

**Table 1. IMMULITE Systems Product(s)**

<b>Assay</b>	<b>Siemens Material Number (SMN)/Catalog Number</b>	<b>Unique Device Identification (UDI)</b>	<b>Lot Number</b>
PRG	10381128/LKPW1	00630414985732	312 and above
PRG	10381181/L2KPW2 10381170/L2KPW6	00630414961972 00630414961989	593 and above

**Reason for Communication**

In July 2019 and January 2020, Siemens Healthcare Diagnostics Inc. issued an Urgent Field Safety Notice (UFSN) IMC19-07.A.OUS and IMC19-07.B.OUS respectively to inform all customers of low discordant results for some patient samples when using the progesterone assay on the IMMULITE systems.

Investigation findings have identified that the presence of anti-progesterone antibodies and/or alkaline phosphatase antibodies in human serum caused this issue. Siemens has determined that increasing the progesterone antibody titer and reformulating the reagent wedge buffer mitigates the low discordant progesterone results.

This has been implemented beginning with kit lot 312 on the IMMULITE/IMMULITE 1000 and beginning with kit lot 593 on the IMMULITE 2000/2000 XPi.

*Resolution of the IMMULITE Systems Low Discordant Progesterone Results on Some Patient Samples*

**Actions to be Taken by the Customer**

- Please review this letter with your Medical Director.

Please retain this letter with your laboratory records and forward this letter to those who may have received this product.

IMMULITE is a trademark of Siemens Healthcare Diagnostics Inc.

*Resolution of the IMMULITE Systems Low Discordant Progesterone Results on Some Patient Samples*

**FIELD CORRECTION EFFECTIVENESS CHECK**

**Resolution of the IMMULITE Systems Low Discordant Progesterone Results on Some Patient Samples**

This response form is to confirm receipt of the enclosed Siemens Healthcare Diagnostics Urgent Field Safety Notice (UFSN) IMC19-07.D.OUS dated October 2022 regarding Resolution of the IMMULITE Systems Low Discordant Progesterone Results on Some Patient Samples. Please read each question and indicate the appropriate answer.

Return this completed form to Siemens Healthcare Diagnostics as per the instructions provided at the bottom of this page.

1. I have read and understood the UFSN instructions provided in this letter. Yes  No

---

Name of person completing questionnaire: \_\_\_\_\_

Title: \_\_\_\_\_

Institution: \_\_\_\_\_ Instrument Serial Number: \_\_\_\_\_

Street: \_\_\_\_\_

City: \_\_\_\_\_ State: \_\_\_\_\_

Phone: \_\_\_\_\_ Country: \_\_\_\_\_

Customer Sold To #: \_\_\_\_\_ Customer Ship To #: \_\_\_\_\_

Please send a scanned copy of the completed form via email to XXXX@XXXX

Or to fax this completed form to the Customer Care Center at XXXXXX

If you have any questions, contact your local Siemens Healthineers technical support representative.

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Par courrier: Siemens Healthineers International AG Zweigniederlassung Zürich,