Urgent Field Safety Notice

Article-No.: 10-806-12, 10-807-12, 10-808-12, 10-809-12, 10-810-12, 10-811-12, 10-812-12,
10-815-12, 10-867-12

Product Description: Reusable surgical panetti and acoustic neurinoma instruments for surgical-invasive procedures in the ENT area or for head and neck surgery excluding the areas of the central circulatory system (CCS) and the central nervous system (CNS).

Enclosed you will find details of a Field Safety Notice initiated by SPIGGLE & THEIS Medizintechnik GmbH for the above listed products. This notice is addressed to users and is intended to provide them with specific instructions for the proper use of the below listed products. Our records indicate that you have received at least one of the products in question and are therefore affected by this measure.

Please read the attached Field Safety Notice and then return the signed confirmation form stating that you have received a copy of the Field Safety Notice and that you understand the instructions for the proper use of the products in question.

You can continue to use the products safely in accordance with the information and precautions in the instructions for use.

Please complete the confirmation form even if you no longer have the products in question in use.

We also use this form to update our records and to avoid further unnecessary communications.

As part of this corrective action, we would like to ask you to complete the enclosed confirmation form and return it to us as follows:

Mrs. Martina Zell (Safty Officer)

or

Mr. Thomas Nüsse (QMB)

Phone: +49 (0) 2206 9081 - 41

Phone: +49 (0) 2206 9081 - 26

Mobile: +49 (0) 162 275 0907

Mobile: +49 (0) 1724300458

m.zell@spiggle-theis.com

t.nuesse@spiggle-theis.com

We hereby confirm that this safety notice has been brought to the attention of the responsible German supervisory authority, the BfArM (Federal Institute for Drugs and Medical Devices), the responsible Swiss supervisory authority, Swissmedic and our Notified Body DQS (0297).

We would like to take this opportunity to thank you for your cooperation in this matter!

NAME OF ESTABLISHMENT:

DEPARTMENT:

NAME:

POSITION:

DATE:

SIGNATURE:

Urgent Field Safety Notice

This letter provides important information regarding the potential for the panetti and acoustic neurinoma instruments to fracture when they are overstrained by excessive force application.

A corresponding warning will be added to the Instruction for Use (IFU).

Specifically twisting and levering may result in a fracture of the fine microsurgical blades of the panetti and acoustic neurinoma instruments.

There have been no deaths or serious adverse events. To ensure patient safety, please comply with the instructions and warnings included in the IFU.

- **The instruments are intended exclusively for microsurgical procedures on soft tissue types such as fine tissue incisions and for the suction process, but not for rinsing. The instruments are not suitable for bone surgery.**
- **Attention: The instruments are sensitive to impacts, throwing, dropping, the use of wire brushes/steel wool and abrasives as well as excessive use of force. It is therefore imperative to take appropriate care when handling them. The instruments are not allowed to be overstressed by twisting or levering. This may result in damage to or breakage of the product or injury to the patient.**
- **Caution: If an instrument breakage occurs due to excessive force on the microsurgical instrument and metal pieces enter the operating field, they should be removed from the operating field as safely, and in a manner that is as minimally invasive, as possible e.g. by flushing or other appropriate intervention. Metal pieces may be localized by X-ray examination in order to remove them completely. In some cases the physician may decide against removal.**

Thank you for reviewing this information. Additional details are included in the following pages. Please follow information indicated in the newly revised IFU for the benefit of your patients. If you have any questions regarding this letter, or would like additional information, please feel free to contact us.

All instruments must be inspected for damage and deterioration prior to clinical use. Particular attention must be paid to fractures, cracks and bends. Damaged instruments must be sorted out. As far as possible, all instruments must be subjected to a functional test.

Immediately inform SPIGGLE & THEIS Medizintechnik GmbH of any adverse events concerning the use of the subject devices. Comply with any local regulations concerning the reporting of adverse events to local Competent Authorities.

As an immediate part of this corrective action an additional warning on the instructions for use, wording "Instrument for microsurgical intervention on soft tissues. Not suitable for bone surgery!" will be added.

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Instructions for use

of reusable surgical Panetti and Acoustic Neurinoma
Instruments

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RxOnly

Caution: Federal law restricts this
device to sale by or on the order of
a physician.
For USA and Canada only.

Hersteller
Manufacturer
Fabricant
Produttore
Fabricante
Výrobce
Κατασκευαστής
制造商

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Reusable Panetti and Acoustic Neurinoma Instruments

Suction dissector for ear drum

10-800-12	30-1822	Dissection and suction of the cholesteatoma around the ear drum
10-801-12	30-1824	
10-813-12	30-1874	
10-814-12	30-1876	
15-800-12	30-1848	
15-801-12	30-1850	
15-813-12		

Suction dissector for sinus tympani

10-802-12	Dissection and suction of the cholesteatoma in the sinus tympani
10-803-12	
15-802-12	
15-803-12	
30-1826	
30-1828	

Suction dissector for atticus

10-804-12	Dissection and suction of the cholesteatoma in the atticus
10-805-12	
15-804-12	
15-805-12	
30-1830	
30-1832	

Suction dissector, blunt

10-817-12	Removal of cholesteatoma in the antrum
10-818-12	
30-1856	
30-1858	

Suction finder, blunt

10-819-12	Exploring the antrum
-----------	----------------------


Suction micro hook

10-806-12	Cutting through the incudo-stapedial-joint; lifting of the footplate
10-867-12	
10-807-12	
30-1834	
30-1835	
30-1836	

Suction separator

10-808-12	Dissolving and removing of the ear drum
15-808-12	
10-809-12	
10-811-12	
30-1838	
30-1840	
30-1844	

Suction knife

10-810-12	Incision of the external meatus; preparation of the meatal flap
10-812-12	
30-1842	
30-1846	

Suction sickle knife

10-815-12	Cutting tissue
30-1852	
30-1870	

Suction curette

10-816-12	Removal of bone structure
30-1854	
30-1872	

Suction raspatory

20-000-15	Preparation and suction of mucosa from the bone at the septum
-----------	---

Septum knife

20-001-15	Incision, preparation and suction of mucosa from the bone at the septum
20-002-15	
20-003-15	
20-000-15	

Suction dissector for acoustic neuroma surgery

10-832-17	Dissection and suction of acoustic neuroma
10-833-17	

Suction sickle knife, curved for acoustic neuroma surgery

10-830-17	Cutting tissue
-----------	----------------

Suction curette, curved for acoustic neuroma surgery

10-831-17	Removal of bone structure
-----------	---------------------------

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ENGLISH

1. Intended use and indication

The products described herein are reusable surgical Panetti and Acoustic Neurinoma Instruments for use in the field of ENT or head and neck surgery excluding areas of the central circulatory system (CCS) and the central nervous system (CNS). The instruments are intended exclusively for microsurgical procedures on soft tissue types such as fine tissue incisions and for the suction process, but not for flushing. Panetti and acoustic neuroma instruments with fine blades and micro hooks are not suitable for bone surgery.

2. Contraindication

Allergic reactions to the materials used.
No other contraindications are known if applied properly.

3. Side effects

No side effects are known if applied properly.

4. Directions for use

These instructions for use describe the standard use of the reusable surgical Panetti and Acoustic Neurinoma Instruments. Special instructions for use with regard to application and/or processing are included in the instructions for use of the 3mach® System Screen Basket for Panetti ear and nose suction instruments (GA_109) which must also be observed when using the instruments.

4.1. General directions for use

Prior to use, examine the instruments for:

- Accurate functioning
- Sufficient sharpness
- Tears
- Surface damage
- Loose parts
- Other damage

Do not use the instruments in the event that any damage is observed.

4.2. Special directions for use

The product should be used by skilled and experienced medical experts only (e.g. qualified medical specialist). The information enclosed (e.g. instructions for use) cannot replace relevant medical and technical basic knowledge. The user must obtain such knowledge in special professional further education courses. The acquisition of special medical knowledge and the associated diagnostic and therapeutic consequences are the sole responsibility of the product user.

- The surgeon has the responsibility for selecting the correct instruments during surgery. The surgeon must also make sure that surgical staff receive reasonable information and training with regard to the instruments and are familiar with the handling of the instruments (for instance with regard to the risk of glare that potentially results from the combined use of instruments and light-emitting devices).
- SPIGGLE & THEIS Medizintechnik GmbH assumes no responsibility for any use of the instruments in an area of application other than the intended one.
- Connect the proximal end to a suction pump via suction adapter and suction tubing. Verify that all parts are securely connected before use.
- The instruments can be connected to a commercially available syringe and used to gently remove secretions.
- Using the instruments for any purpose other than that intended may cause damage or destruction of the instruments, and we disclaim all liability for any damage caused by such abuse.
- If there is a blockage in the instrument or to avoid blockage of the instrument, the instrument should be flushed with a syringe with the suction control closed after each application during the operation.
- Spare products should be readily available in order to prevent delays in use in the event of a blockage.
- Some instruments may be sharp-pointed and/or sharp-edged. Such instruments pose a risk of injury. To avoid bodily injury to third parties, make sure the instruments are stored in protective covers.
- The instruments must be cleaned and sterilised prior to use. We recommend cleaning the instruments twice prior to first use.
For instructions on cleaning and sterilisation, please refer to the chapter below entitled "6. Instructions for processing and sterilising of instruments".
- Although the instruments are made of stainless steel, prolonged storage on surfaces made of base metals (e.g. chromium-plated surfaces) or in the vicinity of chemicals should be avoided since the transmission of rust film cannot be fully excluded.
- Due to the material composition, the instruments are not suitable for use with other magnetically sensitive medical products, (MRI) devices, or instruments.
- The instruments are electrically conductive. To avoid electric shock, make sure the instruments do not come between the patient and a source of electric power.
- The instruments are not suitable for coagulation (e.g. by connection to HF instrument).

Caution: The instruments are sensitive to impacts, throwing, dropping, the use of wire brushes/steel wool and abrasives as well as excessive use of force. It is therefore imperative to take appropriate care when handling them. The instruments are not allowed to be overstressed by twisting or levering. This may result in damage to or breakage of the product or injury to the patient.

Caution: If an instrument breakage occurs due to excessive force on the microsurgical instrument and metal pieces enter the operating field, they should be removed from the operating field as safely, and in a manner that is as minimally invasive, as possible e.g. by flushing or other appropriate intervention. Metal pieces may be localized by X-ray examination in order to remove them completely. In some cases the physician may decide against removal.

5. Storage

The products are delivered in non-sterile condition. Any stock of the products should be kept in suitable containers until initial sterilisation in order to avoid damage.

6. Instructions for processing and sterilising of instruments

- Before first use (of brand-new instruments) and before every reuse, make sure to clean, disinfect, and sterilise instrument thoroughly. Correct sterilisation is only possible with cleaned and disinfected instruments!
- Handle all instruments with the utmost care at all times when transporting, cleaning, maintaining, sterilising and storing the instruments. Fine tips require extra careful handling.
- Always process instruments in homogenous batches sorted by material composition. In case of base materials, a rust film may develop, which could be transferred to the instrument.
- Dismountable instruments shall be dismantled in proper form for preparation

and treatment purposes. In this context, particular attention shall be paid to detachable components such as plastic inlays and/or protective devices made of silicone.

- Where protective covers are used for the protection of sensitive components against mechanical stress during transport and sterilisation, these covers shall be removed in principle for cleaning purposes, since successful cleaning cannot be guaranteed otherwise. The instructions for use and processing of the respective manufacturer of the protective sleeves must be observed with regard to processing.
- For cleaning and disinfection, make sure to use standard agents listed by DGHM/VAH. Please refer to the specifications provided by the manufacturers of the relevant products for instructions for use, exposure times, and suitability of disinfectants and cleaning agents for certain types of instruments.
- We recommend the use of cleaning agents containing a low portion of silicate to avoid spots (silication) on stainless steel instruments.
- Carefully follow the instructions for use of any devices, chemicals etc. used in pre-processing.
- Strictly follow dosing instructions for chemicals and instructions for exposure times and exposure temperatures applicable in cleaning and disinfecting as well as in passivation (see chapter 7.4. Maintenance of instruments).
- Do not overload dishwashers and ultrasonic equipment. Avoid rinsing and acoustical shadows. Make sure items are safely fastened while in the machine.
- Please observe the applicable national provisions on the processing of medical devices.

Depending on the design of the product and the materials used, the maximum number of processing cycles may vary. Our medical products can be reprocessed as long as they continue to function properly. That means that the instruments may be reused as long as they are clean, undamaged, and adequately functional. The user assumes full responsibility if using damaged or soiled instruments. We deny any liability in case of non-compliance with the above.

Caution: When treating patients suffering from Creutzfeldt-Jakob disease (CJD) or its variants (vCJD), disposable instruments should be used since we cannot assume liability for reuse of such instruments. Used instruments should be destroyed in order to prevent cross-infection of other patients and operators.

6.1. Preparation of processing at usage site:

- After the operation, immediately remove any residue of blood, secretions, tissue and bones; do not allow residue to dry! Dried residue will cause corrosion.
- Immediately after use on the patient, immerse the instruments in a suitable cleaning agent/disinfectant. Do not immerse in solution of sodium chloride (NaCl) as chlorides attack the passive layer of instruments and induce pittings. Ensure that the air can escape from the cavities so that all areas are completely wetted. Immersion prevents residue from drying (protein fixation). It is recommended that instruments are processed within one hour after usage.
- Rinse the instruments thoroughly with a suitable cleaner/disinfectant. The patency can be checked additionally by sucking up the cleaner, keeping the suction control closed.
- Do not use fixating agents or hot water (> 40 °C) since both may cause fixation of residues and impair the cleaning effect.
- Store the instruments in a moist environment.

6.2. Transport

Store and transport the instruments to the processing site in a closed container in order to avoid damage to instruments and contamination of the environment.

6.3. Cleaning, disinfecting, and drying

- Ensure that there are no remnants in the lumina.

Pre-cleaning

- Wipe off all visible residues using a moist cloth.
- Immerse instruments in 0.5 % cleaning solution for 10 minutes. Make sure that all instrument parts are covered with cleaning solution.
- Brush the instrument surfaces under running water until all visible contamination has been removed. We expressly advise against the use of wire brushes/steel wool!
- Flush the cavities for 20 seconds using a water gun (3.8 bar static pressure). Keep the suction interrupter closed while flushing, if applicable.
- Clean instruments in the ultrasound bath for 15 minutes at 40 °C (with 0.5 % cleaning solution); if applicable, put the instruments into a metal basket. The channels should first be filled with the cleaning solution from the ultrasound bath, e.g. with a syringe.
- The cavities of the instruments are rinsed with pulsating cold water jets for at least 20 seconds using a water gun (3.8 bar static pressure). Keep the suction interrupter closed while flushing, if applicable.

As a general rule, the instruments should be flushed with a water gun in the water bath to exclude the possibility of contaminating the surrounding area and staff!

Machine cleaning

- Place the instruments in an appropriate Screen Basket and connect them correspondingly as described in the instructions for use of the 3mach® System Screen Basket for the Panetti ear and nose suction instruments (GA_109).
- Do not overload dishwashers. Make sure items are safely fastened while inside the machine.

For cleaning, use the Vario-TD programme of the cleaning and disinfecting device. For the cleaning process, the use of demineralised water is recommended.

Vario-TD programme:

- 4 minutes of pre-cleaning with cold water
- Emptying

- Dosage of cleaner at 40 °C
- 5 minutes of cleaning at 55 °C with a 0.5 % alkaline cleaner or at 45 °C with a 0.5 % enzymatic cleaner
- Emptying
- 3 minutes of neutralisation
- Emptying
- 2 minutes of intermediate flushing
- Emptying

Machine disinfection

The cleaning/disinfection device is provided with a thermal disinfection programme that follows the cleaning process. Thermal machine disinfection should be carried out in compliance with the national requirements for the A0 value (see DIN EN ISO 15883-1). Disinfection should be carried out using demineralised water. We recommend an A0 value of 3000 (e.g. A0 3000 = 90 °C and 5 minutes of application time) for all our reusable Panetti and Acoustic Neurinoma Instruments.

Warning

Insufficient flushing or excessive exposure to disinfectants or cleaning agents may cause corrosion of the instruments. For exposure times, please refer to the package insert of the cleaning agent or disinfectant employed.

Machine drying

Dry the outside of the instruments by means of the drying cycle of the cleaning/disinfection device. If necessary, instruments may in addition be dried manually with a lint-free cloth. Use sterile compressed air to dry the cavities of instruments. Every cleaning and disinfecting device must be provided by the manufacturer with a suitable drying method (see ISO 15883-1). Please refer to the directions and operating instructions provided by the manufacturer of the cleaning/disinfection device.

6.4. Function test, maintenance

Visually inspect instruments for cleanliness. In particular, please pay attention to parts, which are difficult to access, and to uneven surfaces and/or soldering joints. If necessary, repeat the reprocessing cycle until the instrument is visibly clean. The end of the product's life is determined by wear and damage caused by usage. All instruments shall be checked for damage and wear and tear. In particular, attention shall be paid to cracks, fissures and distortions. Damaged instruments must be discarded. All instruments must be subjected to a function test.

Maintenance of instruments

Prior to sterilisation, the instruments must be treated with the care product in accordance with the specifications provided by the manufacturer of the care product (water vapour-permeable paraffin oils respectively medicinal white oils). It must be ensured in any case that the cleaning agent used does not impair the results of sterilisation. Instruments made of stainless steel should be passivated at regular intervals and/or when necessary for maintenance and to be protected against corrosion (pitting corrosion!). Please observe the manufacturer's information for the passivation agent used.

6.5. Packaging

The instruments must be placed in the sterilising device in a standard sterilisation packaging in accordance with DIN EN ISO 11607-1. The packaging must have a clearly visible label which shows all relevant information regarding its contents, i.e. type and date of sterilisation, sterility expiration date, batch number, and sterilisation device number, if any. Individual packaging: Standard packaging materials may be used. The bag must be big enough for the instrument so that no pressure is applied to the sealing, potentially causing tearing of the bag.

6.6. Sterilisation

- Sterilisation of products by means of a fractioned pre-vacuum vapour sterilisation procedure (in acc. with DIN EN 556-1/DIN EN ISO 17665-1) taking into consideration the national requirements in each case:
 - Minimum requirement:
 - Sterilisation temperature: minimum of 132 °C.
 - Exposure time: minimum of 4 minutes (full cycle).
 - Drying: minimum of 10 minutes.

Steam must be free of any substances in order to prevent staining and corrosion. The recommended limit values for substances contained in feed water and steam condensate are defined in EN 285.

When sterilising several instruments at the same time, do not exceed the maximum capacity of the sterilizer (see manufacturer's specifications). Finally, use the relevant indications to make sure the sterilisation result is correct. According to the guidelines of the Robert Koch Institute, the preparation of instruments is completed upon the documented release of the medical device for use.

Storage

Storage must protect from dust, moisture and contamination. Avoid direct sunlight during storage.

Handling of sterile packaging

Check sterile packaging for integrity prior to removing the instruments from the packaging. Please ensure that the relevant aseptic instructions are complied with when removing the instruments from the packaging.

Information on reprocessing validation

The processing method described above has been verified by a validated procedure. The following materials and machines were used in the context of the validations:

- Alkaline cleaner: Neodisher® MediClean; Dr. Weigert; Hamburg.
- Enzymatic cleaner: Neodisher® MediZyme; Dr. Weigert; Hamburg.
- Cleaning/disinfecting device: Miele G 7836 CD.
- Insertion trolley: Miele E 429
- Steam steriliser: Selectomat 666-1 (MMM).
- Sterile packaging: Steriking® Klarsichtrollen (Wipak).

Other chemicals and machines than those mentioned may also be used. In that case, ask the manufacturers and/or suppliers if their products offer the same level of performance as the products used for validation of the procedure. If you choose a processing method different from the method described above, you are responsible for assessing the suitability of such alternative method.

Please also observe the legal provisions applicable in your country and the hygiene rules at the surgery or hospital. This is true in particular to different provisions on effective prion inactivation.

For further information on the preliminary treatment of instruments visit the website of the instrument preparation working group at <http://www.a-k-i.org> or the website of the Robert Koch Institute.

7. Disposal

When disposing the instruments, take all steps necessary to avoid risk of injury and infection. Collect and enclose sharp-edged and sharp-pointed instruments in semi-solid and unbreakable containers. Protect against unauthorised access during storage.

Contaminated instruments must be disposed of as hazardous waste and handled so as to avoid contamination of third parties.

8. Maintenance and repairs

Maintenance and repairs may be performed by the manufacturer only. Maintenance work or repairs performed on the instruments by the user or by an unauthorised third party will result in loss of all guarantee and warranty claims. Instruments returned to the manufacturer must undergo full processing and sterilisation cycle first. Sharp-edged and sharp-pointed instruments must be properly protected prior to shipping.

For services, repairs, or spare parts outside of Germany please contact the dealer or sales representative of your country.

9. Declaration of Conformity

We declare under sole responsibility the conformity of the products in their supplied version with the requirements of MDD 93/42/EEC and/or MDR 2017/745. The products comply with the relevant provisions of the MDD 93/42/EEC and/or MDR 2017/745 and are labelled with the CE-mark.

With reservation of changes and typographical errors!
© SPIGGLE & THEIS Medizintechnik GmbH, 201907_01.

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