

### **URGENT - Field Safety Notice**

### **Medical Device Correction**

**Forte, Forte JETStream, Forte JETStream AZ, Forte JETStream AZ upgrade,  
Forte JETStream upgrade, Diamond Select Forte, Diamond Select Forte  
JETStream**

***Forte Detector Unimpeded Motion***

## **Discontinue system use until further notice**

Dear Customer,

A problem has been detected in the Philips Forte product line that, if it were to re-occur, could pose a risk for patients or users. This Field Safety Notice is intended to inform you about:

- what the problem is and under what circumstances it can occur
- the actions that should be taken by the customer / user in order to prevent risks for patients or users
- the actions planned by Philips to correct the problem.

**This document contains important information for the continued safe and proper use of your equipment**

Please review the following information with all members of your staff who need to be aware of the contents of this communication. It is important to understand the implications of this communication.

Please retain a copy with the equipment Instruction for Use.

If you need any further information or support concerning this issue, please contact your local Philips representative:

**0800 80 3000**

This notice has been reported to the appropriate Regulatory Agency.

Sincerely,



### URGENT - Field Safety Notice

### Medical Device Correction

**Forte, Forte JETStream, Forte JETStream AZ, Forte JETStream AZ upgrade, Forte JETStream upgrade, Diamond Select Forte, Diamond Select Forte JETStream**

#### *Forte Detector Unimpeded Motion*

### Discontinue system use until further notice

<b>AFFECTED PRODUCTS</b>	882020 Forte, 882290 Forte JETStream, 882291 Forte JETStream upgrade, 882320 Forte JETStream AZ, 882321 Forte JETStream AZ upgrade, 889456 Diamond Select Forte, 889471 Diamond Select Forte JETStream
<b>PROBLEM DESCRIPTION</b>	As a result of a customer reported problem during detector radius movement, Philips identified an issue affecting the Forte family of gamma cameras that could result in either detector 1 or detector 2 falling unimpeded vertically to the end stops of its travel limit.
<b>HAZARD INVOLVED</b>	There is a possibility that a detector may fall unimpeded vertically to the end stops of its travel limit, due to a mechanical failure, possibly making contact with a patient that could result in entrapment, serious injury, or death.
<b>HOW TO IDENTIFY AFFECTED PRODUCTS</b>	This issue applies to all Forte gamma cameras:  Forte, Forte JETStream, Forte JETStream upgrade, Forte JETStream AZ, Forte JETStream AZ upgrade, Diamond Select Forte, Diamond Select Forte JETStream  The product type can be determined by finding the name of the system of the product label. Refer to the representative photos below.



### URGENT - Field Safety Notice

### Medical Device Correction

**Forte, Forte JETStream, Forte JETStream AZ, Forte JETStream AZ upgrade, Forte JETStream upgrade, Diamond Select Forte, Diamond Select Forte JETStream**

*Forte Detector Unimpeded Motion*

**Discontinue system use until further notice**

	
<p><b>ACTION TO BE TAKEN BY CUSTOMER / USER</b></p>	<ul style="list-style-type: none"> <li>➤ <b>Discontinue use of the system until further notice.</b></li> <li>➤ Inform those who need to be aware within your organization or any organization where the potentially affected devices have been transferred (If appropriate).</li> <li>➤ Maintain this notice with your system Instructions for Use (IFU) until the correction is made to the system.</li> </ul>



## **URGENT - Field Safety Notice**

### **Medical Device Correction**

**Forte, Forte JETStream, Forte JETStream AZ, Forte JETStream AZ upgrade,  
Forte JETStream upgrade, Diamond Select Forte, Diamond Select Forte  
JETStream**

***Forte Detector Unimpeded Motion***

### **Discontinue system use until further notice**

<b>ACTIONS PLANNED BY PHILIPS</b>	Philips Healthcare is distributing this FSN to all affected customers/users and will deploy a solution addressing the issue upon completion of the investigation.
<b>FURTHER INFORMATION AND SUPPORT</b>	If you need any further information or support concerning this issue, please contact your local Philips representative:  <b>0800 80 3000</b>



### Customer Response Form

**INSTRUCTIONS:** Please complete this form with the customer upon delivery of the Field Safety Notice and return to the BIU. **Please email completed form to [customercare.ch@philips.com](mailto:customercare.ch@philips.com).**

Customer Name:	System Code	System Serial No.
Address:		

Our records indicate that your firm has received affected systems. By signing this form, you acknowledge having received, read, and understood the content of this letter and have taken appropriate actions.

**Signature:** \_\_\_\_\_ **Date:** \_\_\_\_\_

**Print Name:** \_\_\_\_\_ **Email:** \_\_\_\_\_

