Urgent Safety Information

regarding the
NADAL® hCG Pregnancy Test
Product Code: 141000
LOT: 169757

Important customer information

Moers, 29.11.2019

For all users and distributors of the following *in-vitro* diagnostic test

Title: NADAL® hCG Pregnancy Test
Product Code: 141000
LOT / Batch: 169757
Expiry date: 30.11.2020

Dear customer,

This urgent safety information concerns the NADAL® hCG pregnancy test.

Issue description:

In some cases in the batch listed above, the pregnancy test strips are enclosed in primary packaging (protective pouch) which bears, on one side, the nal von minden Drug-Screen designation and visual instructions intended for the evaluation of drug test results. These visual instructions are not intended for use in evaluating the results of a pregnancy test, as in drug tests 2 lines is a negative result and 1 line is a positive result. The opposite is true in pregnancy tests.

Implications:

Evaluating the pregnancy test according to the visual instructions on the incorrect pouch will produce an incorrect result, as positive and negative results are reversed.

Action to be taken:

In carrying out the test and evaluating the test results of this batch of pregnancy tests, refer only to the instructions for use included in the enclosed package insert. **Ignore the visual instructions printed on the protective pouch.**

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Laborservice

Beratung & Service

For questions, please contact your nal von minden contact person or use the contact details listed below.

Please ensure that all users of the aforementioned product and other relevant staff members within your organisation are aware of this **urgent safety information**. If this product has been passed along to third parties, please forward a copy of this information or inform the contact person listed below. Please retain this information as long as these safety measures are in effect. The relevant regulatory body has received a copy of this urgent safety information.

Please fill out and send back the enclosed acknowledgement of receipt at your earliest convenience.

We apologise for any inconvenience caused.

Yours sincerely,

Sandra von Minden
Safety Officer for Medical Devices

nal von minden GmbH
Carl-Zeiss-Street 12
47445 Moers
Deutschland
Tel.: +49 2841 99 820 – 0
Fax: +49 2841 99 820 – 1
E-Mail: s.vonminden@nal-vonminden.com

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Laborservice

Beratung & Service

Customer acknowledgement of receipt

With your signature, you acknowledge that you have taken note of the urgent safety information of 29.11.2019 concerning

Title: NADAL® hCG Pregnancy Test
Product Code: 141000
LOT / Batch: 169757
Expiry date: 30.11.2020

Complete address (stamp):

Date: _____

Signature: _____

Name: _____

At your earliest convenience, please return this acknowledgement of receipt to:

nal von minden GmbH
Fax: +49 2841 99 820 – 1
E-Mail: info@nal-vonminden.com

or to your contact person at nal von minden GmbH

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14. Referenzen

- Ed. Lothar Thomas, Labor und Diagnose (8th Edition, 2012) TH Books, Chapter 38.3, 1859 ff
- Butler SA et al. (2001) Detection of Early Pregnancy Forms of Human Chorionic Gonadotropin by Home Pregnancy Test Devices, Clinical Chemistry 47(12): 2131-2136
- Lenton EA et al. (1982) Plasma concentrations of human chorionic gonadotropin from the time of implantation until the second week of pregnancy, Fertility and Sterility, 37(6): 773-778
- Fritz MA et al. (1987) Doubling time of human chorionic gonadotropin (hCG) in early normal pregnancy: relationship to hCG concentration and gestational age, Fertility and Sterility, 47(4): 584-589
- Stenman U-H et al. (2006) The classification, functions and clinical use of different isoforms of hCG, Human Reproduction Update 12(6): 769-84
- <http://americanpregnancy.org/while-pregnant/hcg-levels/> from 15-8-19, 2:08 pm
- Gronowski AM et al. (2009) False-Negative Results in Point-of-Care Qualitative Human Chorionic Gonadotropin(hCG) Devices Due to Excess hCGβ Core Fragment, Clinical Chemistry, 55 (7): 1389-1394

Rev. 2, 2019-09-18 OM/JB/KB

This package insert contains information on the following NADAL® hCG Pregnancy Tests:

Ref. Nr.:	Cut-off	Packaging
141000	10 mIU/ml in urine (Early Test)	Test strip, individually packed
141000N-10	10 mIU/ml in urine (Early Test)	Test strip, individually packed
143003	10 mIU/ml in urine (Early Test)	Container with single test strips
141002	25 mIU/ml in urine	Test strip, individually packed
143004	25 mIU/ml in urine	Container with single test strips
151003	10 mIU/ml in urine or serum (Early Test)	Test strip, individually packed
151002	25 mIU/ml in urine or serum	Test strip, individually packed
153003	25 mIU/ml in urine or serum	Container with single test strips

1. Intended Use

The NADAL® hCG Pregnancy Test (141002/ 143004/ 151002/ 153003) or the NADAL® hCG Pregnancy Early Test (141000/ 141000N-10/143003/151003) is a rapid, visual immunoassay for the qualitative, presumptive detection of human chorionic gonadotropin (hCG) in human urine or serum specimens. The test is intended for use as an aid in early detection of pregnancy and designed for professional use only.

2. Introduction and Clinical Significance

Human chorionic gonadotropin (hCG) is a glycoprotein hormone produced by the developing placenta shortly after fertilisation and can already be detected in the early stage of pregnancy in urine and serum.

In normal pregnancies, hCG can be detected in both serum and urine as early as 7 days after conception. The hCG level continues to rise very rapidly, frequently exceeding 50 mIU/ml by the first missed menstrual period, and culminates in the 10th-12th weeks of pregnancy in the range of 30,000-290,000 mIU/ml.

The early appearance of hCG in urine and blood soon after conception and its subsequent rapid rise in concentration makes it an ideal marker for the early detection of pregnancy.

Elevated hCG levels, comparable to the early stages of pregnancy may be associated with a hydatidiform mole or choriocarcinoma. For this reason, these possibilities should be ruled out before a pregnancy is diagnosed.


The NADAL® hCG Pregnancy Test is a rapid, visual immunoassay for the detection of hCG at the concentration of 10 mIU/ml or 25 mIU/ml in urine (141000/141000N-10/ 143003/141002/143004) and/or serum (151002/151003/ 153003).

The test uses anti-hCG antibodies to detect hCG in the urine (or serum). The specificity of this immune reaction rules out cross-reactivity with structurally-related hormones (hFSH, hLH and hTSH) at above the maximum anticipated physiological levels.

3. Test Principle

The NADAL® hCG Pregnancy Test detects human chorionic gonadotropin through visual interpretation of colour development on the test strip. Anti-hCG antibodies are immobilised in the test line region of the membrane. Anti-mouse antibodies are immobilised in the control line region. During testing, the specimen reacts with anti-hCG antibodies conjugated to coloured particles and pre-coated onto the conjugate pad of the membrane. The mixture then migrates along the membrane by capillary action and interacts with reagents on the membrane. If there is sufficient hCG in the specimen, a coloured line will appear in the test line region of the membrane. The presence of this coloured line indicates a positive result, while its absence indicates a negative result.

The appearance of a coloured line in the control line region (C) serves as a procedural control, indicating that the proper volume of specimen has been added and membrane wicking has occurred.

4. Reagents and Materials Supplied

141000N-10 for urine	10 single test strips individually packed
141000/141002 for urine	50 single test strips individually packed
143003/ 143004 for urine	25 single test strips in container
151003 for urine or serum	50 single test strips individually packed
151002 for urine or serum	50 single test strips individually packed
153003 for urine or serum	25 single test strips in container

- 1 package insert

5. Additional Materials Required

- Specimen collection container
- Timer

6. Storage & Stability

Test kits should be stored at 2-30°C until the expiry date printed on the packaging. The test strip should remain in the sealed foil pouch or the sealed tube until use. Reseal the tube immediately after removing a test strip. This will ensure that

the remaining test strips remain dry. Once the tube has been opened, the test strips can be used for up to 3 months. Note the date when the tube is first opened and ensure that the tests are used within this period of time. Do not freeze tests.

Care should be taken to protect the components of the test kit from contamination. Do not use test strips if there is evidence of microbial contamination or precipitation. Biological contamination of dispensing equipment, specimen collection containers or reagents can lead to inaccurate results.

7. Warnings and Precautions

- For professional *in-vitro* diagnostic use only.
- Carefully read through the test procedure prior to testing.
- Do not use the test beyond the expiration date indicated on the package.
- Do not use the test if the foil pouch is damaged.
- Do not reuse tests.
- The test strip should remain in the sealed foil pouch until use.
- Do not dip the test strip beyond the maximum immersion mark.
- Do not add samples to the reaction area (result area).
- In order to avoid contamination, do not touch the reaction area (result area) or the immersion area.
- Only urine samples (and/or serum samples for 151002/151003/153003) and no other liquids should be used for testing.
- Bring samples to room temperature (15-30°C) prior to testing.
- Avoid cross-contamination of specimens by using a new specimen collection container for each specimen obtained.
- Do not eat, drink or smoke in the area where specimens and test kits are handled.
- Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are being assayed.
- Handle all specimens as if they contain infectious agents. Observe established precautions for microbiological risks throughout all procedures and standard guidelines for the appropriate disposal of specimens.
- The test kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not completely guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious, and handled in accordance with usual safety precautions (e.g., do not ingest or inhale).
- Humidity and high temperature can adversely affect test results.
- Used testing materials should be discarded according to local regulations.

8. Specimen Collection and Preparation

Urine samples:

- Collect a urine specimen in a clean and dry plastic or glass specimen collection container which is free from any preservatives. From the first day on of the missed period the hCG levels are normally so high that urine samples

collected at any time of the day may be used. However, excessive fluid intake should be avoided (dilution effect).

- If the test is used at an earlier point in time, the first morning urine should be used, as it generally contains the highest hCG concentration. For tests with a cut-off level of 25 mIU/ml, testing can be carried out starting from 2 days before the expected period. For tests with a cut-off level of 10 mIU/ml, the earliest point in time for testing is 4 days before the expected period. However, the risk of false negative results increases the earlier the testing is carried out (see also point 12 'Limitations').
- Specimens and tests should be brought to room temperature for the correct execution of the test procedure. Urine samples exhibiting visible precipitates should be centrifuged, filtered or allowed to settle to obtain clear supernatant for testing. Only clear supernatant should be used for testing.
- If testing is not to be performed immediately, urine specimens should be stored at 2-8°C for up to 48 hours. For long-term storage, specimens should be kept below -20°C.
- Prior to testing, frozen specimens must be completely thawed and mixed well. Specimens should not be frozen and thawed repeatedly.


Serum sample (only for 151002, 151003 and 153003):

- Collect a blood specimen aseptically by venipuncture.
- To avoid haemolysis, separate serum from cruro as soon as possible. Use only clear, non-haemolysed specimens. Specimens containing particulate matter should be clarified by centrifugation prior to assay.
- Icteric, lipemic, haemolysed and contaminated sera may lead to incorrect results.
- Ideally, testing should be performed immediately after specimen collection. Do not leave specimens at room temperature for prolonged periods of time. Serum specimens should be stored at 2-8°C for up to 3 days. For long-term storage, specimens should be kept below -20°C.
- Bring specimens to room temperature prior to testing. Prior to testing, frozen specimens must be completely thawed and mixed well. Specimens should not be frozen and thawed repeatedly.

9. Test Procedure

Bring tests and specimens to room temperature (15-30°C) before testing.


1. Remove a test strip from the sealed foil pouch or from the sealed tube. **Note:** If a strip is removed from the tube, close it immediately afterwards. Label the strip with a patient or control identification. For best results the assay should be performed without delay.
2. Hold the test strip vertically at the marked end and immerse it into the specimen for at least 10-15 seconds. In order to avoid contamination, do not touch the membrane of the test strip, and make sure not to immerse the test strip beyond the maximum immersion mark (MAX).

or

As the test begins to run, you will observe a coloured liquid migrate along the membrane.

- Remove the test strip from the specimen and place it on a dry, clean, non-absorbent surface.
- Start the timer.
- Wait for the coloured line(s) to appear. Read the test result after 5 minutes, though samples with high hCG concentrations may show positive results in as little as 3 minutes. Do not interpret the result after more than 10 minutes.

Note: With a longer reading time, the risk of false positive results increases, as increasingly lower hCG concentrations are detected with increasing incubation time. As a result, in rare cases, basal concentrations in non-pregnant women may also be detected – especially in menopausal women.

10. Result Interpretation

Positive:

Two coloured lines appear on the membrane. One line appears in the control line region (C) and the other line appears in the test line region (T).

Negative:

Only one coloured line appears in the control line region (C). No coloured line appears in the test line region (T).

Invalid:

The control line (C) fails to appear. Results from any test which has not produced a control line at the specified reading time must be discarded. Please review the procedure and repeat the test with a new test strip. If the problem persists, discontinue using the test kit immediately and contact your distributor.

Note:

- The colour intensity in the test line region (T) may vary depending on the concentration of the analyte present in the specimen. Therefore, any shade of colour in the test line region should be considered positive. Note that this is a qualitative test only and it cannot determine the concentration of the analyte in the specimen.
- Insufficient specimen volume, incorrect operating procedure or expired tests are the most likely reasons for the control line failure.
- In the case of a strong test line (T), the control line (C) may be very pale. The test result should nevertheless be considered valid.
- In the case of a negative test result (possibly with strongly diluted samples) where pregnancy is nonetheless suspected, the test should be repeated after 48-72 hours or, in addition, a quantitative test should be performed,

since the hCG concentration may not be high enough. The first morning urine contains the highest hCG concentration.

11. Quality Control

An internal procedural control is included in the test strip:

A coloured line appearing in the control line region (C) is considered an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique.

Good laboratory practice (GLP) recommends the use of control materials to ensure proper test kit performance.

12. Limitations

- The NADAL® hCG Pregnancy Tests are for professional *in-vitro* diagnostic use and should only be used for the qualitative detection of human chorionic gonadotropin (hCG).
- Very diluted urine specimens, exhibiting low specific gravity, may not contain representative levels of hCG. If pregnancy is still suspected after a negative test result, a first morning urine sample should be obtained 48-72 hours later and tested.
- If tests are used earlier than before the first day of a missed period, there is a risk of false negative results due to individual fluctuations of hCG concentrations. According to the literature, if 10 mIU/ml tests are used 4 days before the expected period, only about 53-74% of pregnancies are initially detected when using urine samples. If pregnancy is suspected, testing should therefore be carried out again with morning urine at a later point in time. Alternatively, the hCG concentration can be determined using a more sensitive quantitative method.
- In healthy men and healthy, fertile, non-pregnant women, no hCG can be normally detected using NADAL® hCG Pregnancy Tests. Basal hCG values increase in the menopause, so that in rare cases, especially with 10 mIU/ml tests, false positive results may occur.
- Very low levels of hCG (<50 mIU/mL) are present in urine or serum shortly after the implantation of a fertilised ovum in the womb. However, since a significant number of first trimester pregnancies miscarry for natural reasons, a test result that is weakly positive should be confirmed using other clinical methods. If necessary, the test should be carried out again with the first morning urine at a later point in time. The result of the retesting should be clearly positive.
- If, after a weak positive result, the result of the retesting using the first morning urine at a later point in time is negative, a spontaneous miscarriage might be the cause. Approximately 22% of clinically-undetected pregnancies end with an unnoticed miscarriage and are perceived as a slightly delayed period.
- A number of conditions other than pregnancy, including trophoblastic diseases, such as molar pregnancy, choriocarcinoma, testicular tumors, prostate cancer and lung cancer as well as medication containing hCG cause elevated levels of hCG which can be detected as positive by the test.

- Genito-urinary infections can lead to a false-positive hCG test result. An increased number of leucocytes, a basic pH, an increased amount of nitrite, erythrocytes and proteins in urine among others are indicative of genito-urinary infections. In case of suspicion, please verify this by means of urinalysis strips or specific swab tests. Very suspicious urine samples should be verified by laboratory testing. In the case of a positive result, no pregnancy should be diagnosed until the aforementioned possibilities are excluded.
- The tests are calibrated for intact hCG. The proportions of intact hCG and other hCG forms in serum and urine may differ according to the pregnancy stage. Especially at the beginning of pregnancy, when hCG values are still low, this can lead, in rare cases, to deviating results for both sample materials. In rare cases, at the beginning of pregnancy, deviating results may also be observed with tests based on the detection of beta-hCG.
- For urine samples: In rare cases, at a more advanced stage of pregnancy (approximately week 5-8 of pregnancy), increased amounts of hCG-beta-core fragments may be present in the urine hindering the detection reaction and leading to false negative results in spite of high hCG concentration. [7]
- For serum samples: As with all assays that utilise anti-mouse antibodies, there is a possibility of interference with anti-mouse antibodies (HAMA) in the sample. Samples from patients who have received preparations of monoclonal antibodies for diagnosis or therapy may contain HAMA (immunologically interfering substances). Such samples may cause false positive or false negative test results. Heterophile antibodies in serum may also cause erroneous results.
- As with all diagnostic tests, a definitive clinical diagnosis should not be based on the results of a single test, but must only be made by the physician after all clinical and laboratory findings have been evaluated. It is therefore recommended to confirm the obtained results with another testing method or to interpret them in conjunction with clinical findings.

Expected values:

The reference range (serum, 97.5 percentile, beta hCG) of healthy, fertile, non-pregnant women is <3 mIU/ml. In postmenopausal women this value increases to <5 mIU/ml. However, some isolated higher values can be found in other literature. After fertilisation there is a sharp rise so that, approximately 10 days after fertilisation (4 days before the missed period), the values of ≥ 10 mIU/ml are often already reached. The doubling time varies during the course of the pregnancy and is approximately 2.5 days at the beginning of the pregnancy. Due to variations in liquid absorption there is a high individual variability for hCG concentrations in urine. In the course of the first trimester the concentrations of 30,000-290,000 are determined, which then decrease in the 2nd and 3rd trimester to values of 45,000-100,000 mIU/ml. After birth hCG serum values decrease and are no longer detectable after approximately 11-17 days. [1], [6]

13. Performance Characteristics

Detection limit and measuring range

The lower detection limit of the NADAL® hCG Pregnancy tests is set at 10 mIU/ml (for 141000/141000N-10/143003/151003) or 25 mIU/ml (for 141002/143004/151002/153003) (Calibration with the 4th IRP of hCG). NADAL® hCG Pregnancy tests show no prozone effect up to a maximum concentration of 600,000 mIU/ml.

Analytical specificity

Cross-reactivity

In a cross-reactivity study with the structurally related hormones hLH (300 mIU/ml), hFSH (1000 mIU/ml) and hTSH (1 mIU/ml) the NADAL® hCG Pregnancy tests showed no cross-reactivity. The concentrations used are higher than the physiologically maximum expected values of healthy women.

Interference

hCG-free urine samples and urine samples containing hCG concentrations at the detection limit were spiked with numerous compounds. None of the compounds interfered with result formation.

Substance	Concentration
Acetaminophen	0.2 mg/ml
Acetylsalicylic acid	0.2 mg/ml
Ascorbic acid	0.2 mg/ml
Atropine	0.2 mg/ml
Bilirubin	0.2 mg/ml
EDTA	0.2 mg/ml
Ethanol	1.0%
Caffeine	0.2 mg/ml
Gentisic acid	0.2 mg/ml
Glucose	20 mg/ml
Haemoglobin	10 µg/ml
Urea	20 mg/ml
Uric acid	0.2 mg/ml
Creatinine	0.2 mg/ml
Oxalic acid	0.4 mg/ml
THC	0.1 mg/ml

Diagnostic performance characteristics

In clinical studies, NADAL® hCG Pregnancy tests were tested with clinical sample material in comparison with another rapid test with the identical cut-off or enzyme immunoassay (EIA). Relative diagnostic performance characteristics were calculated on the basis of the results. The results are presented in the following tables.

Ref. 141002/143004 (urine test 25 mIU/ml), sample: urine

NADAL® hCG Pregnancy Test	Other hCG rapid test		
	Positive	Negative	Total
	Positive	66	0
Negative	0	142	142
Total	66	142	208

Relative diagnostic sensitivity: >99.9%

Relative diagnostic specificity: >99.9%

Overall agreement: >99.9%

Ref. 141000/141000N-10/143003 (urine test 10 mIU/ml), sample: urine

NADAL® hCG Pregnancy Test	Other hCG rapid test		
	Positive	Negative	Total
	Positive	79	0
Negative	0	131	131
Total	79	131	210

Relative diagnostic sensitivity: >99.9%

Relative diagnostic specificity: >99.9%

Overall agreement: >99.9%

Ref. 151002/153003 (urine/serum test 25 mIU/ml), sample: urine

		EIA		
		Positive	Negative	Total
NADAL® hCG Pregnancy Test	Positive	130	0	130
	Negative	0	178	178
	Total	130	178	308

Relative diagnostic sensitivity: >99.9%

Relative diagnostic specificity: >99.9%

Overall agreement: >99.9%

Ref. 151002/153003 (urine/serum test 25 mIU/ml), sample: serum

		EIA		
		Positive	Negative	Total
NADAL® hCG Pregnancy Test	Positive	169	0	169
	Negative	0	250	250
	Total	169	250	419

Relative diagnostic sensitivity: >99.9%

Relative diagnostic specificity: >99.9%

Overall agreement: >99.9%

Ref. 151003, (urine/serum test 10 mIU/ml), sample: urine/serum

		EIA		
		Positive	Negative	Total
NADAL® hCG Pregnancy Test	Positive	94	0	94
	Negative	0	106	106
	Total	94	106	200

Relative diagnostic sensitivity: >99.9%

Relative diagnostic specificity: >99.9%

Overall agreement: >99.9%

14. References

1. Ed. Lothar Thomas, Labor und Diagnose (8th Edition, 2012) TH Books, Chapter 38.3, 1859 ff
2. Butler SA et al. (2001) Detection of Early Pregnancy Forms of Human Chorionic Gonadotropin by Home Pregnancy Test Devices, Clinical Chemistry 47(12): 2131-2136
3. Lenton EA et al. (1982) Plasma concentrations of human chorionic gonadotropin from the time of implantation until the second week of pregnancy, Fertility and Sterility, 37(6): 773-778
4. Fritz MA et al. (1987) Doubling time of human chorionic gonadotropin (hCG) in early normal pregnancy: relationship to hCG concentration and gestational age, Fertility and Sterility, 47(4): 584-589
5. Stenman U-H et al. (2006) The classification, functions and clinical use of different isoforms of HCG, Human Reproduction Update 12(6): 769-84
6. <http://americanpregnancy.org/while-pregnant/hcg-levels/from-15-8-19,2:08pm>
7. Gronowski AM et al. (2009) False-Negative Results in Point-of-Care Qualitative Human Chorionic Gonadotropin(hCG) Devices Due to Excess hCGβ Core Fragment, Clinical Chemistry, 55 (7): 1389-1394

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Concordance générale : >99,9%

14. Bibliographie

- Ed. Lothar Thomas, Labor und Diagnose (8th Edition, 2012) TH Books, Chapter 38.3, 1859 ff
- Butler SA et al. (2001) Detection of Early Pregnancy Forms of Human Chorionic Gonadotropin by Home Pregnancy Test Devices, Clinical Chemistry 47(12): 2131-2136
- Lenton EA et al. (1982) Plasma concentrations of human chorionic gonadotropin from the time of implantation until the second week of pregnancy, Fertility and Sterility, 37(6): 773-778
- Fritz MA et al. (1987) Doubling time of human chorionic gonadotropin (hCG) in early normal pregnancy: relationship to hCG concentration and gestational age, Fertility and Sterility, 47(4): 584-589
- Stenman U-H et al. (2006) The classification, functions and clinical use of different isoforms of HCG, Human Reproduction Update 12(6): 769-84
- <http://americanpregnancy.org/while-pregnant/hcg-levels/from-15-8-19,2:08-pm>
- Gronowski AM et al. (2009) False-Negative Results in Point-of-Care Qualitative Human Chorionic Gonadotropin(hCG) Devices Due to Excess hCGβ Core Fragment, Clinical Chemistry, 55 (7): 1389-1394

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Concordancia general: >99,9%

14. Referencias

- Ed. Lothar Thomas, Labor und Diagnose (8th Edition, 2012) TH Books, Chapter 38.3, 1859 ff
- Butler SA et al. (2001) Detection of Early Pregnancy Forms of Human Chorionic Gonadotropin by Home Pregnancy Test Devices, Clinical Chemistry 47(12): 2131-2136
- Lenton EA et al. (1982) Plasma concentrations of human chorionic gonadotropin from the time of implantation until the second week of pregnancy, Fertility and Sterility, 37(6): 773-778
- Fritz MA et al. (1987) Doubling time of human chorionic gonadotropin (hCG) in early normal pregnancy: relationship to hCG concentration and gestational age, Fertility and Sterility, 47(4): 584-589
- Stenman U-H et al. (2006) The classification, functions and clinical use of different isoforms of hCG, Human Reproduction Update 12(6): 769-84
- <http://americanpregnancy.org/while-pregnant/hcg-levels/from-15-8-19>, 2:08 pm
- Gronowski AM et al. (2009) False-Negative Results in Point-of-Care Qualitative Human Chorionic Gonadotropin(hCG) Devices Due to Excess hCGβ Core Fragment, Clinical Chemistry, 55 (7): 1389-1394

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14. Bibliografia

1. Ed. Lothar Thomas, Labor und Diagnose (8th Edition, 2012) TH Books, Chapter 38.3, 1859 ff
2. Butler SA et al. (2001) Detection of Early Pregnancy Forms of Human Chorionic Gonadotropin by Home Pregnancy Test Devices, Clinical Chemistry 47(12): 2131-2136
3. Lenton EA et al. (1982) Plasma concentrations of human chorionic gonadotropin from the time of implantation until the second week of pregnancy, Fertility and Sterility, 37(6): 773-778
4. Fritz MA et al. (1987) Doubling time of human chorionic gonadotropin (hCG) in early normal pregnancy: relationship to hCG concentration and gestational age, Fertility and Sterility, 47(4): 584-589
5. Stenman U-H et al. (2006) The classification, functions and clinical use of different isoforms of hCG, Human Reproduction Update 12(6): 769-84
6. <http://americanpregnancy.org/while-pregnant/hcg-levels/from-15-8-19-2:08-pm>
7. Gronowski AM et al. (2009) False-Negative Results in Point-of-Care Qualitative Human Chorionic Gonadotropin(hCG) Devices Due to Excess hCGβ Core Fragment, Clinical Chemistry, 55 (7): 1389-1394

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14. Referências

- Ed. Lothar Thomas, Labor und Diagnostik (8th Edition, 2012) TH Books, Chapter 38.3, 1859 ff
- Butler SA et al. (2001) Detection of Early Pregnancy Forms of Human Chorionic Gonadotropin by Home Pregnancy Test Devices, Clinical Chemistry 47(12): 2131-2136
- Lenton EA et al. (1982) Plasma concentrations of human chorionic gonadotropin from the time of implantation until the second week of pregnancy, Fertility and Sterility, 37(6): 773-778
- Fritz MA et al. (1987) Doubling time of human chorionic gonadotropin (hCG) in early normal pregnancy: relationship to hCG concentration and gestational age, Fertility and Sterility, 47(4): 584-589
- Stenman U-H et al. (2006) The classification, functions and clinical use of different isoforms of hCG, Human Reproduction Update 12(6): 769-84
- <http://americanpregnancy.org/while-pregnant/hcg-levels/from-15-8-19-2:08-pm>
- Gronowski AM et al. (2009) False-Negative Results in Point-of-Care Qualitative Human Chorionic Gonadotropin(hCG) Devices Due to Excess hCGβ Core Fragment, Clinical Chemistry, 55 (7): 1389-1394

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Celková shoda: > 99,9 %

14. Reference

- Ed. Lothar Thomas, Labor und Diagnose (8th Edition, 2012) TH Books, Chapter 38.3, 1859 ff
- Butler SA et al. (2001) Detection of Early Pregnancy Forms of Human Chorionic Gonadotropin by Home Pregnancy Test Devices, Clinical Chemistry 47(12): 2131-2136
- Lenton EA et al. (1982) Plasma concentrations of human chorionic gonadotropin from the time of implantation until the second week of pregnancy, Fertility and Sterility, 37(6): 773-778
- Fritz MA et al. (1987) Doubling time of human chorionic gonadotropin (hCG) in early normal pregnancy: relationship to hCG concentration and gestational age, Fertility and Sterility, 47(4): 584-589
- Stenman U-H et al. (2006) The classification, functions and clinical use of different isoforms of HCG, Human Reproduction Update 12(6): 769-84
- <http://americanpregnancy.org/while-pregnant/hcg-levels/from-15-8-19,2:08pm>
- Gronowski AM et al. (2009) False-Negative Results in Point-of-Care Qualitative Human Chorionic Gonadotropin(hCG) Devices Due to Excess hCGβ Core Fragment, Clinical Chemistry, 55 (7): 1389-1394

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Samlet aftale: >99,9%

14. Referencer

- Ed. Lothar Thomas, Labor und Diagnose (8th Edition, 2012) TH Books, Chapter 38.3, 1859 ff
- Butler SA et al. (2001) Detection of Early Pregnancy Forms of Human Chorionic Gonadotropin by Home Pregnancy Test Devices, Clinical Chemistry 47(12): 2131-2136
- Lenton EA et al. (1982) Plasma concentrations of human chorionic gonadotropin from the time of implantation until the second week of pregnancy, Fertility and Sterility, 37(6): 773-778
- Fritz MA et al. (1987) Doubling time of human chorionic gonadotropin (hCG) in early normal pregnancy: relationship to hCG concentration and gestational age, Fertility and Sterility, 47(4): 584-589
- Stenman U-H et al. (2006) The classification, functions and clinical use of different isoforms of HCG, Human Reproduction Update 12(6): 769-84
- <http://americanpregnancy.org/while-pregnant/hcg-levels/from-15-8-19,2:03pm>
- Gronowski AM et al. (2009) False-Negative Results in Point-of-Care Qualitative Human Chorionic Gonadotropin(HCG) Devices Due to Excess hCGβ Core Fragment, Clinical Chemistry, 55 (7): 1389-1394

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Algehele overeenkomst: >99,9%

14. Referenties

- Ed. Lothar Thomas, Labor und Diagnose (8th Edition, 2012) TH Books, Chapter 38.3, 1859 ff
- Butler SA et al. (2001) Detection of Early Pregnancy Forms of Human Chorionic Gonadotropin by Home Pregnancy Test Devices, Clinical Chemistry 47(12): 2131-2136
- Lenton EA et al. (1982) Plasma concentrations of human chorionic gonadotropin from the time of implantation until the second week of pregnancy, Fertility and Sterility, 37(6): 773-778
- Fritz MA et al. (1987) Doubling time of human chorionic gonadotropin (hCG) in early normal pregnancy: relationship to hCG concentration and gestational age, Fertility and Sterility, 47(4): 584-589
- Stenman U-H et al. (2006) The classification, functions and clinical use of different isoforms of hCG, Human Reproduction Update 12(6): 769-84
- http://americanpregnancy.org/while-pregnant/hcg-levels/from 15-8-19, 2:08 pm
- Gronowski AM et al. (2009) False-Negative Results in Point-of-Care Qualitative Human Chorionic Gonadotropin(hCG) Devices Due to Excess hCGβ Core Fragment, Clinical Chemistry, 55 (7): 1389-1394

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