Rev 1: September 2018

FSN Ref: FSN-AG191103

FSCA Ref: FSCA-AG191103

Date: 15/01/2020

**Urgent Field Safety Notice**  
**STERNAL RETRACTOR IMA – REF DC30200-00**

For Attention of\*: [REDACTED]

|  |
|--|
| Contact details of local representative (name, e-mail, telephone, address etc.)* |
|--|

|  |
|--|
| <b>MCM Medsys AG – Cassandra Somma - Glutz-Blotzheim-Str.3 – 4500 Solothurn<br/>Switzerland.</b> |
|--|

**Urgent Field Safety Notice (FSN)**  
**STERNAL RETRACTOR IMA – REF DC30200-00**

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| <b>2 Reason for Field Safety Corrective Action (FSCA)*</b> |  |
|--|--|
| 2  | <p style="color: red;"><b>1. Description of the product problem*</b></p> <p>Damage of the device and burr formation during use due to a lack of lubrication. / Endommagement du dispositif et formation de bavures lors de l'utilisation par manque de lubrification</p> |
| 2  | <p style="color: red;"><b>2. Hazard giving rise to the FSCA*</b></p> <p>Risk for the patient : In worst case, parts could fall unnoticed into the patient during</p>   |

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| <b>3. Type of Action to mitigate the risk*</b> |   |
|--|---|
| <b>3.</b>                                      | <p><b>1. Action To Be Taken by the User*</b></p> <p> <input checked="" type="checkbox"/> Identify Device    <input type="checkbox"/> Quarantine Device    <input type="checkbox"/> Return Device    <input type="checkbox"/> Destroy Device </p> <p> <input type="checkbox"/> On-site device modification/inspection<br/> <input type="checkbox"/> Follow patient management recommendations<br/> <input checked="" type="checkbox"/> Take note of amendment/reinforcement of Instructions For Use (IFU)<br/> <input type="checkbox"/> Other                      <input type="checkbox"/> None </p> <p>Information already transmitted on November 12<sup>th</sup>, 2019 to the distributor as part of the expert report and the compliance of the incriminated MD.<br/>Information déjà transmise le 12 novembre 2019 au distributeur dans le cadre du rapport d'expertise et de la mise en conformité du DM incriminé.</p> |
| <b>3.</b>                                      | <p><b>2. By when should the action be completed?</b></p> <p>Laser marking on the new devices already applied since May 25th, 2018. Update of the instructions for use I AQ 008-20191220-FR-EN-DE-IT-v4 realised on December 12th, 2019. / Marquage sur les dispositifs neufs déjà en place depuis le 25/05/2018. Mise à jour de la notice IAQ 008-20191220-FR-EN-DE-IT-v4 réalisée le 20/12/2019.</p>   |

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| <b>4. General Information*</b> |   |
|--------------------------------|---|
| 4.                             | 1. FSN Type*<br>New   |
| 4.                             | 2. For updated FSN, reference number and date of previous FSN<br>NA   |
| 4.                             | 3. For Updated FSN, key new information as follows:<br>NA   |
| 4.                             | 4. Further advice or information already expected in follow-up FSN? *<br>No   |
| 4                              | 5. If follow-up FSN expected, what is the further advice expected to relate to:<br>NA   |
| 4                              | 6. Anticipated timescale for follow-up FSN<br>NA  |
| 4.                             | 7. Manufacturer information<br>(For contact details of local representative refer to page 1 of this FSN)  |
|                                | a. Company Name<br>LANDANGER  |
|                                | b. Address<br>ZI LA VENDUE – 52906 CHAUMONT - FRANCE  |
|                                | c. Website address<br>landanger@landanger.com –<br>aguillaumot@landanger.com  |
| 4.                             | 8. The Competent (Regulatory) Authority of your country has been informed about this communication to customers. *                              |
| 4.                             | 9. List of attachments/appendices:<br>NA  |
| 4.                             | 10. Name/Signature<br>Alain GUILLAUMOT, quality manager<br> |

| <b>Transmission of this Field Safety Notice</b> |  |
|---|--|
|   | <p>This notice needs to be passed on all those who need to be aware within your organisation or to any organisation where the potentially affected devices have been transferred. (As appropriate)</p> <p>Please transfer this notice to other organisations on which this action has an impact. (As appropriate)</p> <p>Please maintain awareness on this notice and resulting action for an appropriate period to ensure effectiveness of the corrective action.</p> <p>Please report all device-related incidents to the manufacturer, distributor or local representative, and the national Competent Authority if appropriate, as this provides important feedback..*</p> |

Note: Fields indicated by \* are considered necessary for all FSNs. Others are optional.