

# Important Notice

**Product code and description:** IE700 - Optilite Analyser  
**Date:** 19 June 2020  
**Subject:** New Software Version: Version 7.0.2

**Dear Customer,**

You were recently provided a Field Safety Notice (GLB132169, issued in April 2020) in relation to a potential software calibration curve issue when using Optilite software version 7.0 or 7.0.1. We are pleased to inform you that the Optilite software update to version 7.0.2 has now been released by The Binding Site. This software update will correct the calibration curve issue experienced with software versions 7.0 and 7.0.1..

You can expect to be contacted by a Binding Site representative within 3 months to schedule an upgrade for the Optilite Analyser(s) at your site. Once the Optilite software has been upgraded to version 7.0.2 the Binding Site representative will complete an acknowledgement form which will require a signature for our records.

In the meantime, please continue to follow the instructions issued in the earlier Field Safety Notice, which allow continued safe operation of the analyser and have been included below for your reference.

**Actions that should have been taken by the user:**

The Binding Site requested that the following instructions were followed to ensure that the software issue affecting calibration curve selection in Software Version 7.0 and 7.0.1 cannot occur:

1. Upon the receipt of the Field Safety Notice (GLB132169) perform a full re-calibration of all assays. Only one kit lot of each assay should be calibrated.  
*(Note: this step does not need to be performed again, if all steps were taken and continued to be followed after receipt of the Field Safety Notice)*
2. Ensure only one kit lot for each assay is in use for testing.
3. To re-calibrate assays, even if not prompted by the software, when any parameter changes are made including minor parameter changes in the test definitions screen, or the QC triggers (see list of relevant parameter changes in the Field Safety Notice).

These actions must be followed until Software Version 7.0.2 is provided. In the meantime, should you have any additional questions on this, please contact your local Binding Site technical support.

Kind regards,



Colin Booth  
*Director of Quality and Regulatory*

**Transmission of this Important information:**

This notice needs to be passed on to all those who need to be aware within your organisation or to any organisation where the potentially affected devices have been transferred.

Please transfer this notice to other organisations on which this information has an impact.

Please maintain awareness on this notice and resulting action until Software Version 7.0.2 is provided to ensure effectiveness of the corrective action.

Should you require any further information please contact

Your local Binding Site Representative

or

Technical Support Group

UK: +44(0) 1214569696

[technical.support@bindingsite.com](mailto:technical.support@bindingsite.com)