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**Produit : TaqPath™ COVID-19 CE-IVD RT-PCR kit (Catalog #A48067)  
Applied Biosystems™ COVID-19 CE-IVD Interpretive Software v1.3 and v2.3**

5 août 2020

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**Urgent Field Safety Notice**

**Mandatory Upgrade to TaqPath™ COVID-19 CE-IVD Interpretive Software v1.3 and v2.3**

Date:

Institution:	<<Customer_Name>>	<<Address_Line_2>>
Address:	<<Address_Line_4>>	<<Address_Line_3>>
	<<Address_Line_1>>	<<City>>, <<State>> <<Post Code>> <<Country>>

Dear Valued Customer,

This letter is to notify you of a mandatory upgrade for the TaqPath™ COVID-19 CE-IVD Interpretive Software.

- If you are using TaqPath COVID-19 CE-IVD Interpretive Software v1.2 (used with the 7500 series Real-Time PCR instruments) you need to upgrade to software v1.3.
- If you are using TaqPath COVID-19 CE-IVD Interpretive Software v2.1, or v2.2 (used with the QuantStudio Real-Time PCR platforms), you need to upgrade to software v2.3.

An internal review of a cross-section of customer data and feedback from using the TaqPath™ COVID-19 CE-IVD RT-PCR Kit resulted in two observed issues:

1. The previous interpretive software versions may not recognize suboptimal extractions (poorly extracted samples). The development team initially determined that if the extraction failed, the MS2 Ct would be >37. However, a review of customer data has shown that suboptimal extractions can result in lower MS2 Ct values than expected. This issue could potentially cause a weakly positive sample that is prepared by suboptimal extraction to be falsely called a valid negative sample, thereby constituting a false negative.
2. The MS2 assay, which detects the Internal Positive Control (IPC), may be erroneously called amplified in a small percentage of Positive Control (PC) samples. This issue could cause a plate to be designated invalid, thereby requiring unnecessary retesting of an entire batch of samples.

**Software Algorithm Improvements**

Versions 1.3 and 2.3 address both issues by applying an optimized Ct cutoff for MS2. Our investigation indicates that most users will not experience these errors. However, given that errors could occur, we have optimized the MS2 settings to reduce the risk of false test outcomes and to decrease the retest burden on customers. No other analysis settings have been changed.

**Instructions for Use Clarifications**

In addition to these software algorithm improvements, we have updated the Instructions for Use to highlight the importance of mixing the RT-PCR reaction plates to further mitigate potential false positives. Insufficient mixing of qPCR plates could produce a “waterfall” effect in which decreasing signal in the baseline region ultimately produces erroneous Ct values, potentially leading to an incorrect positive call for one or more targets.

**Actions Required:**

- Notify all COVID-19 CE-IVD Interpretive Software users in your facility of this communication.
- Stop using COVID-19 CE-IVD Interpretive Software v1.2, v2.1, and v2.2.
- Go to [www.thermofisher.com/COVID19EUATraining-ISWv13-33](http://www.thermofisher.com/COVID19EUATraining-ISWv13-33) to complete the training and upgrade the COVID-19 CE-IVD Interpretive Software. Re-registration is required for this upgrade.
- Sign and return acknowledgement of Receipt – page **3** of this communication.

We apologize for any inconvenience. Thermo Fisher Scientific is committed to supplying innovative, high-quality products. If you have any technical questions or concerns, please refer to the Frequently Asked Questions document provided with this communication. If you still have questions or concerns, please contact Technical Support at [\[800-955-6288 option 4 or techsupport@thermofisher.com\]](tel:800-955-6288) *add your country TS info here*.

Sincerely,

Genetic Sciences Division  
Thermo Fisher Scientific

DRAFT

**CUSTOMER RESPONSE SHEET**  
**Acknowledgement and Receipt Form**  
**(Customer response is required)**

Product Name	SKU Number	Number of Instances Upgraded <b>(Customer to populate)</b>
COVID19 INTRPTV CEIVD V1.2	100093771	
COVID19 INTPRTV CEIVD V2.1	100094143	
COVID19 INTPRTV CEIVD V2.2	100094318	

Customer Name and Address:

---

---

I have read and understand the instructions provided in the letter dated <enter date>.  YES  NO

I have notified all software users in my facility:	<input type="checkbox"/> YES <input type="checkbox"/> NO
If "No," please explain:	

Customer Signature of Receipt: \_\_\_\_\_ Date: \_\_\_\_\_

Please return response sheet by email to [xxx.xxxx@thermofisher.com](mailto:xxx.xxxx@thermofisher.com).

For any questions or comments please contact Technical Support at 1-XXX-XXX-XXXX.

Thermo Fisher Scientific Reference Number: PR XXXXX