

**URGENT**

Ortho Clinical Diagnostics

Month ##, 2021

**URGENT FIELD SAFETY NOTIFICATION****Potential Intermittent False Positives on ORTHO VISION and ORTHO VISION Max BioVue Analyzers When Testing High Titer Samples**

Dear Customer,

Ortho Clinical Diagnostics previously issued a letter (Ref. CL2020-206a) regarding the potential for intermittent false positive results being generated on ORTHO VISION® or ORTHO VISION® Max analyzers when processing ORTHO Sera anti-D (IAT) and Indirect Antiglobulin Test (IAT) crossmatch (XM) tests after pipetting plasma or serum samples with high titer ABO antibodies. An investigation was conducted to determine if any other tests might be impacted. Along with the previously affected assays in the table below, the newly identified assays that were impacted are noted with \* in the table below.

Affected Product	Product Code Unique Device Identifier No.
ORTHO VISION® Analyzer for ORTHO® BioVue Cassettes	<b>6904579</b> (10758750012831)
ORTHO VISION® Max Analyzer for ORTHO® BioVue Cassettes	<b>6904578</b> (10758750012848)

Impacted Test Type	Associated Product Codes
ORTHO™ Sera Antigen Typing <ul style="list-style-type: none"><li>• ORTHO™ Sera Anti-Fya</li><li>• ORTHO™ Sera Anti-Fyb</li><li>• ORTHO™ Sera Anti-S</li><li>• ORTHO™ Sera Anti-s</li><li>• ORTHO™ Sera Anti-D (IAT)</li></ul>	<ul style="list-style-type: none"><li>• 6904486*</li><li>• 6904487*</li><li>• 6904490*</li><li>• 6904491*</li><li>• 6904493</li></ul>
IAT Crossmatch and IAT Autocontrol* performed on <ul style="list-style-type: none"><li>• AHG Anti-IgG Ortho BioVue® System cassettes</li><li>• Anti-Human Globulin Anti-IgG, -C3d; polyspecific Ortho BioVue System® cassettes</li></ul>	<ul style="list-style-type: none"><li>• 707400/707450</li><li>• 707300/707350</li></ul>
IAT Dilution Series* performed on <ul style="list-style-type: none"><li>• AHG Anti-IgG Ortho BioVue® System cassettes</li><li>• Anti-Human Globulin Anti-IgG, -C3d; polyspecific Ortho BioVue System® cassettes</li></ul> <p>In conjunction with type A<sub>1</sub> or B (reagent) red blood cells</p>	<ul style="list-style-type: none"><li>• 707400/707450</li><li>• 707300/707350</li></ul>

## Issue Description

---

Ortho Clinical Diagnostics received complaints of discordant positive reactions for healthy donor and patient samples.

## Impact to Results

---

As is described in the ORTHO VISION / ORTHO VISION Max Reference Guides, plasma samples with antibody titers >1024 may intermittently cause carryover in subsequent test columns.

Ortho's investigation of low frequency intermittent false positive test results has determined that, for the specific test types listed in the table above, type O plasma samples with ABO antibody titers  $\geq 1024$  may cause false positive test results due to carryover.

Unexpected positive results may lead the blood bank laboratory to perform additional testing to confirm the sample results.

If a high titer ABO antibody sample ( $\geq 1024$ ) is suspected of being processed it is recommended to review column results for all plasma and/or cell dispenses that occurred after the high titer sample. Discuss any concerns you may have regarding previously reported results with your Laboratory Medical Director to determine the appropriate course of action.

## Carryover Claims

---

The Carryover Claims were established using an Anti-D antibody and state: *Testing of the ORTHO VISION® Analyzer indicated that a sample with a high-titered antibody (>1024) when tested may intermittently cause carryover in the subsequent sample test columns. The number of columns affected is dependent on the magnitude of the high titer sample. Testing also indicated that carryover was not observed in samples with antibody titers of 512 or 1024.*

During an investigation on low frequency intermittent false positive results it was concluded that the potential for intermittent carryover exists for Anti-A at an antibody titer  $\geq 1024$  and may generate a false positive result with red blood cells that express higher copy number of A antigen, while having no observable effect on the expected negative results with red blood cells expressing lower copy number of A antigen.

## Investigation

---

A limited number of customers, approximately 1.42% of the ORTHO VISION Platform's install base, have notified ORTHO about experiencing this issue.

Our investigation is focused on two potential causes of antibody carryover on ORTHO VISION systems and its impact on results of the test types listed in the table above.

Acute carryover – A false positive result caused by carryover of ABO antibodies from a previously pipetted high titer ( $\geq 1024$ ) samples into subsequently pipetted tests on the ORTHO VISION analyzers. This would occur when antibodies from the previously pipetted plasma or serum samples are not completely washed from the pipettor and contaminates subsequently pipetted test fluids.

Chronic carryover - The probability of a false positive test result due to carryover of ABO antibodies can slowly increase over the course of the lifetime of the pipettor assembly due to repeated exposure to antibody containing samples during normal pipetting processes. Daily probe maintenance cleans the internal surface of the pipettor assembly and applies a protective coating of Bovine Serum Albumin (BSA) to prevent this build up effect. As is described below, Ortho has validated enhanced daily maintenance procedures to more effectively limit potential impact of the chronic carryover.

This investigation has determined that generation of false positives due to ABO antibody carryover on ORTHO VISION and ORTHO VISION Max systems, is dependent on the following factors.

- Impacted test type(s) is run on ORTHO VISION / ORTHO VISION Max analyzers (impacted test types listed in above table).
- Titer of ABO antibodies in the previously pipetted serum or plasma sample ( $\geq 1024$ ).
- Aspiration volume of previously pipetted serum or plasma sample, and of fluids subsequently pipetted for the impacted test.
- The amount of time that pipetted test fluids reside in the pipettor between aspiration and dispense.
- A and/or B antigen copy number of the Red Blood Cells (RBCs) in the impacted test.

This issue is not reagent-related per se and, therefore, a replacement of reagents will not resolve it.

### **Resolution**

---

Ortho recommends following the mitigations defined in the ORTHO VISION/ VISION Max reference guides, if a plasma sample with ABO antibody titer equal to or greater than 1024 has been processed.

Ortho has also validated an optional enhanced daily maintenance procedure using 0.5 M NaOH (instead of 0.1 M NaOH that is currently used) when chronic carryover is suspected. Refer to the enclosed Technical Bulletin (Pub No. J68774) for its use. This optional enhanced daily maintenance procedure is meant to serve as an incremental improvement to reduce effects of chronic carryover in the interim.

Ortho is actively working on further improvements to mitigate carryover of ABO antibodies into subsequently pipetted fluids which will be introduced in a future software release.

---

### **REQUIRED ACTIONS**

---

- Please refer to ORTHO VISION/ORTHOR VISION Max reference guide recommendations (Pub No. J55655/Pub. No. J55657) if you suspect carryover has occurred. Specifically, If a high titer antibody sample ( $\geq 1024$ ) is suspected of being processed it is recommended to review column results for all plasma and/or cells dispenses that occurred after the high titer sample and to perform the Daily Probe Maintenance procedure.
- Assess your laboratory's needs for whether it would be beneficial to implement the enhanced daily maintenance procedure in the enclosed Technical Bulletin.
- Complete the enclosed Confirmation of Receipt form no later than **Month XX, 2020**.
- Please forward this notification if the product was relocated outside of your facility.
- Save this notification with your user documentation.

### **Contact Information**

---

If you have questions, please contact Ortho Care™ Technical Solutions Center at **insert number**.

**Insert signatory if appropriate in your region.**

Enclosure:

Technical Bulletin (Pub No. J68774)

## Confirmation of Receipt – Response Required

Communication ID: CL2021-110a\_EU

Date of Issue: 2021-MM-DD

### URGENT FIELD SAFETY NOTICE

## Potential Intermittent False Positives on ORTHO VISION and ORTHO VISION Max BioVue Analyzers When Testing High Titer Samples

Please return this completed form by **fax or scan to PDF** and email so that we can complete our records no later than:

**DD-MM-YYYY**Send to: **Name**e-Mail Address: **email address**Fax: **Fax Number**

### Your Name and Address

Verify your name and mailing address:

Please complete this section if any of this information has changed

Institution/

Contact Name: \_\_\_\_\_

Address: \_\_\_\_\_

City: \_\_\_\_\_

State/Prov: \_\_\_\_\_

Zip/Postal Code: \_\_\_\_\_

Phone: \_\_\_\_\_

Fax: \_\_\_\_\_

e-Mail: \_\_\_\_\_

### Please Confirm

I received the Urgent Field Safety Notice (Ref. CL2021-110a\_EU) regarding the potential for intermittent false positive results being generated on ORTHO VISION® or ORTHO VISION® Max BioVue analyzers with test types listed in the customer letter, after pipetting high titer plasma or serum samples.

I understand that I am advised to refer to ORTHO VISION/ORTHOTH VISION Max reference guide recommendations (Pub No. J55655/Pub. No. J55657) if I suspect carryover has occurred. Specifically, If a high titer antibody sample (>1024) is suspected of being processed it is recommended to review column results for all plasma and/or cells dispenses that occurred after the high titer sample and to perform the Daily Probe Maintenance procedure.

I have assessed my laboratory's needs for whether it would be beneficial to implement the enhanced daily maintenance procedure in the enclosed Technical Bulletin.

I understand this Urgent Field Safety issue is not reagent-related per se and, therefore, a replacement of reagents will not resolve it.

Signature: \_\_\_\_\_

Required  
Your signature confirms  
that you have received  
and understand this  
communication

Print Name: \_\_\_\_\_

Phone Number: \_\_\_\_\_

Date: \_\_\_\_\_

Your Comments: \_\_\_\_\_

# Technical Bulletin

ORTHO VISION® Analyzer for ORTHO BioVue® Cassettes  
ORTHO VISION® Max Analyzer for ORTHO BioVue® Cassettes

Pub. No.: J68774EN  
Issued: 2021-03-12

## Optional use of 0.5 M NaOH in Daily Probe Maintenance

### 1 Purpose

This Technical Bulletin provides important information regarding the optional usage of 0.5 M NaOH in daily maintenance to mitigate the possibility of antibody carryover.

### 2 Applicability

This Technical Bulletin is applicable for both the ORTHO VISION® and ORTHO VISION® Max Analyzers for ORTHO BioVue® Cassettes.

### 3 Changes to User's Guides

The next version of the Reference Guide and Self-Service Customer Procedure Guide may be revised to include the information within this bulletin:

- ORTHO VISION® Analyzer – ORTHO BioVue® Cassettes Reference Guide (J55655)
- ORTHO VISION® Max Analyzer – ORTHO BioVue® Cassettes Reference Guide (J55657)
- ORTHO VISION® Analyzer– ORTHO BioVue® Cassettes Self-Service Customer Procedure Guide (J55658)
- ORTHO VISION® Max Analyzer– ORTHO BioVue® Cassettes Self-Service Customer Procedure Guide (J55660)

#### 3.1 ORTHO VISION® Analyzer – ORTHO BioVue® Cassettes Reference Guide (J55655) and ORTHO VISION® Max Analyzer – ORTHO BioVue® Cassettes Reference Guide (J55657)

##### 3.1.1 Chapter 8 Resource Categories: Section: Materials Required but Not Supplied

The materials listed below are required for use with the system but are not supplied:

- 70% Isopropyl alcohol
- ORTHO™ 7% BSA (Bovine Serum Albumin) (5ml) (REF: 6844285)

**Note:** ORTHO™ 7% BSA is required for the Daily Probe Maintenance task, to condition the probe after decontamination.

- Buffered saline
- Deionized or distilled water
- ORTHO VISION® Dilution Trays
- Mild detergent
- ORTHO VISION® Evaporation Caps
- ORTHO Red Cell Diluent
- ORTHO BLISS
- Reagent Red Blood Cells

# Ortho Clinical Diagnostics

- ORTHO BioVue Cassettes
- ORTHO Control Samples
- 0.1 M NaOH or 0.5 M NaOH

**Note:** Some ORTHO™ Sera tests require NaOH (0.1 M or 0.5 M) to decontaminate the probe.


**Note:** 0.1 M NaOH is recommended for Daily Probe Maintenance. If antibody carryover is suspected, you may use 0.5 M NaOH for Daily Probe Maintenance.

## 3.2 ORTHO VISION® Analyzer – ORTHO BioVue® Cassettes Self-Service Customer Procedure Guide (J55658) and ORTHO VISION® Max Analyzer – ORTHO BioVue® Cassettes Self-Service Customer Procedure Guide (J55660)

### 3.2.1 Chapter 2 Maintenance: Section: Daily Probe Maintenance

The next release of the guide, Step 2 may be modified for ORTHO VISION® and ORTHO VISION® Max Analyzers Self-Service Customer Procedure Guide.

1. When prompted, open the Load Station Door.



2. Add 5 ml of 0.1 M NaOH or 0.5 M NaOH to a 10 ml vial with a supported barcode **1**.  
Place the vial into position 3 **2** of a Diluent Rack.

**Note:** 0.1 M NaOH is recommended for Daily Probe Maintenance. If antibody carryover is suspected, you may use 0.5 M NaOH for Daily Probe Maintenance.

## 4 Storing this Bulletin

Place this bulletin in the front section of your ORTHO VISION / ORTHO VISION Max Analyzer ORTHO BioVue Cassettes Reference Guide and Self-Service Customer Procedure Guide.

# Ortho Clinical Diagnostics

This page intentionally left blank.

# Ortho Clinical Diagnostics



IVD

EC REP

Ortho-Clinical Diagnostics  
1500 Boulevard Sébastien Brant  
B.P. 30335  
67411 Illkirch  
CEDEX, France



Ortho-Clinical Diagnostics  
Felindre Meadows  
Pencoed  
Bridgend  
CF35 5PZ  
United Kingdom

Pub. No.: J68774EN

All trademarks are registered to  
Ortho Clinical Diagnostics

© Ortho Clinical Diagnostics 2021

