

14 September 2020

URGENT - FIELD SAFETY NOTICE

Type of Action	Recall
Teleflex Reference:	EIF-000439
Commercial Name	Rusch TracFlex Plus Phonation Set, Cuffed Ruschcare TracFlex Plus Phonation Set, Cuffed
Product Code / Lot Number	See Appendix 2 for a list of product codes and lots in scope.

Dear Customer,

Teleflex Medical has voluntarily issued a recall for the product codes and lot numbers listed above. All affected lot numbers are listed within **Appendix 2**.

Description of the problem & immediate actions required

Teleflex Medical is recalling the above-mentioned products following receipt of two complaints reporting a leak in the cuff. The issue was identified prior to use during the pre-use integrity check as outlined in the IFU therefore there was no patient involvement. In the unlikely event that a user does not perform the pre-use check per IFU, and a cuff leak is present, it may lead to an inadequate seal and a risk of a break in ventilation requiring medical intervention.

Our records indicate you have received products that are subject to this recall.

Depending on your device location please adhere to the following Action list:

Device location	Action List Number
Medical facilities	1
Distributors	2

Action list number 1 – Medical facilities

1. We request that you check your inventory for product within the scope of this FSCA. Users should cease use and distribution of impacted product and quarantine immediately.
2. If you do have stock in scope of this FSCA, mark the according checkbox on the Acknowledgement Form (Appendix 1) and contact customer service by calling the phone number mentioned below. Customer service will issue you with a return number. Write the return number into the respective field in the Acknowledgement Form and return this form immediately to Customer Service.
3. If you do not have stock in scope of this FSCA mark the according checkbox on the Acknowledgement Form (Appendix 1) and return the form to the fax number or e-Mail address mentioned below.
4. Teleflex (or your local dealer) will issue a credit note upon receipt of the returned affected product.

Action list number 2 – Distributors

1. Provide this field safety notice to all customers who have received product in scope of this FSCA. Your customer is then required to complete the acknowledgement form and return to you.

2. We request that you check your inventory for product within the scope of this FSCA. Cease use and distribution of impacted product and quarantine immediately. You may then return all product in scope to Teleflex.
3. As a distributor, you are then required to confirm to Teleflex that you have completed the field activity outlined above. Upon completion of your actions, please forward the completed Acknowledgement Form to Customer Service.
4. Please be aware that all European Economic Area/Switzerland (EEA/CH) and Turkey Member State Competent Authorities in which Teleflex distribute directly will be notified by Teleflex.
5. If you have further distributed product outside of your country, please notify Teleflex by return email to the e-Mail address below.
6. If you are a distributor and/or have a reporting responsibility within or outside the EEA/CH/TR region, please notify your local Competent Authority of this action. Please forward the notification and all communication with your local competent authority to Teleflex.

Teleflex

Teleflex informs all customers, employees of Teleflex and distributors of this Field Safety Corrective Action.

Transmission of this Field Safety Notice

This notice should be passed on to all persons who need to be aware within your organisation or to any organisation where the potentially affected devices have been transferred. Please consider end users, clinicians, risk managers, supply chain/distribution centres etc. in the circulation of this notice. Maintain awareness of this notice until all required actions have been completed in your organisation.

Contact reference person

Should you require any further information or support concerning this issue, please contact:

Customer Service:

Contact: Shane Kenny

FAX: +353 0 14370773

Telephone: +353 0 90 64608769

Email: recalls.intl@teleflex.com

Please be advised that all Economic Area/Switzerland (EEA/CH) and Turkey Member State Competent Authorities to which Teleflex distribute directly will be notified by Teleflex. Teleflex is committed to providing high quality, safe and effective products. We sincerely apologise for any inconvenience this action may cause your operations. If you have any other questions, feel free to contact your local sales representative or Customer Service.

For and on behalf of Teleflex,

Padraig Hegarty

Padraig Hegarty VP, QA (Manufacturing)

FIELD SAFETY CORRECTIVE ACTION
ACKNOWLEDGEMENT FORM

PRODUCT FIELD ACTION BY TELEFLEX - IMMEDIATE ATTENTION REQUIRED

Ref. EIF-000439

RETURN COMPLETED FORM BY IMMEDIATELY TO:

FAX: +353 (0) 1 4370773

Email: Recalls.Intl@teleflex.com

<input type="checkbox"/> We confirm receipt of this FSN and completed the required actions contained therein. We confirm that our inventory does NOT include products affected by this Field Action.	<input type="checkbox"/> We confirm receipt of this FSN and completed the required actions contained therein. We confirm our inventory DOES include products affected by this Field Action. The use and further distribution of the affected products is stopped. All products are put on hold and the amount below will be returned. Return Authorisation No _____
---	---

PLEASE PRINT PRODUCT QUANTITY NUMBERS CLEARLY.

COMMERCIAL NAME OF AFFECTED PRODUCTS:		
PRODUCT NUMBER	LOT NUMBER	QUANTITY (Returning)

- Include a copy of the **completed Acknowledgement Form** in the returns package with the returned units
- Ensure the **RAN number is clearly visible** on the returns package.
- Please label returns as **"Field Action Returns"**

Complete this Acknowledgement form and return immediately by using the fax number or e-mail address above.

INSTITUTION NAME (EG NAME OF HOSPITAL, HEALTH CARE ORGANISATION)	
INSTITUTION ADDRESS	Phone / Fax
FORM COMPLETED BY:	Stamp
PRINT NAME: _____	
SIGNATURE: _____	
DATE	

Appendix 2

Product Codes	Lots in Scope							
121902-000070	17AT21	17CT20	17JT06	19KT07	19LT20			
121902-000080	17BT03	17GT28	17JT18	17JT25	19BT11	19JT67	19KT42	19LT20
121902-000090	17BT03	17ET20	17GT28	17JT13	17JT25	19DT26	19ET82	19GT64
	19KT10	19KT17	19LT23	19LT39				
121902-000100	17BT03	17BT09	17GT28	17JT12	19BT11	19DT26	19FT51	19JT67
	19KT26	19LT20						
121902-000110	17BT03	17BT16	17CT20	17GT28	19JT67			
858002-000070	17AT21	17JT06	19JT67					
858002-000080	17BT03	19BT12	19ET19	19FT51	19HT30	19JT67	19KT18	19LT23
858002-000090	17AT21	17BT03	18FT46	19BT08	19ET82	19HT17	19KT10	19LT38
858002-000100	17BT03	19ET56	19HT17	19JT67	19LT23			
858002-000110	17BT03	19LT20						