Date: 11.09.2020

**Dringend: Sicherheitshinweis**  
**E-ventus BX Peripheral Stent Graft System**

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**Clinical experience with the use of E-ventus BX Peripheral Stent Graft System as bridging stent:**

JOTEC is conducting three observational, prospective, non-randomised, multicenter studies in which peripheral covered stents are used as bridging stents. In all three studies the treatment of the patient is at the discretion of the physician and the patient. Participating physicians are asked to provide their observations collected during routine care by transferring the clinical data from the patient file into the eCRF.

**Information about the studies:**

**PLIANT:**

**ClinicalTrials.gov Identifier:** NCT02209194

**Title: PLIANT** - A post-market registry in **P**atients with **iL**iac **A**neurysm undergoing endovascular stenting with a **N**ew generation of low profile E-iliac **S**Tent Graft System

**Objective:** To evaluate clinical and technical success as well as safety and feasibility of the E-iliac Stent Graft System used in endovascular treatment of uni- or bilateral aorto-iliac or iliac aneurysm. Main study target is the exclusion of aneurysm with primary patency of the arteria iliaca interna and the arteria iliaca externa on iliac implantation side.

**Patient Population:** A total number of 40 male and female patients with unilateral or bilateral common iliac aneurysm, treated with the E-iliac Stent Graft System.

**Study status:**

- Study start: July 2014
- End of study: July 2019
- Participating centers: 11
- Patients enrolled: 45
- Follow-up: 3 years

**Peripheral covered stents used to bridge the side branch of the E-iliac Stent Graft with the internal iliac artery:**

- **Total number:** **62**
  - E-ventus BX (JOTEC) 46
  - Advanta V12 (Getinge) 13
  - Lifestream (Bard) 2
  - B-Graft (Bently) 1

**Observed incidents related to the peripheral covered stents:**

- None

**PLIANTII:**

**ClinicalTrials.gov Identifier:**

**Title: PLIANTII** - Prospective multicenter registry to examine the real-world performance of the E-liac Stent Graft System for treatment of uni- or bilateral aorto-iliac or iliac aneurysms

**Objective:** The PLIANT II registry is undertaken to examine the real-world outcome after treatment of consecutive patients with uni- or bilateral aorto-iliac or iliac aneurysms using the E-liac Stent Graft System.

**Patient Population:** A total number of 400 - 500 male and female patients with uni- or bilateral aorto-iliac or iliac aneurysm, treated with the E-liac Stent Graft System. Recruitment period is estimated to last 5 years and is conducted in up to 60 clinical European centers.

**Study status:**

- Start of enrollment: July 2018
- Patients enrolled: 90 (recruitment is ongoing)
- Participating centers: 26

**Peripheral covered stents used to bridge the side branch of the E-liac Stent Graft with the internal iliac artery:**

- **Total number:** **136**
  - E-ventus BX (JOTEC) 57
  - Advanta V12 (Getinge) 31
  - Viabahn VBX (Gore) 30
  - B-Graft (Bently) 10
  - Viabahn Endoprosthesis (Gore) 3
  - Fluency (Bard) 2
  - Covera (Bard) 2
  - Unknown 1

**Observed incidents related to the peripheral covered stents:**

- One patient with endoleak type IIIId due to PTFE detachment (E-ventus BX)

**CONNECT:**

**ClinicalTrials.gov Identifier:** NCT03295682

**Title:** **CONNECT** - Registry in Patients with Thoracoabdominal Aneurysms Treated with Multi-branch Stent Grafts Tailored to Their Individual Anatomies

**Objective:** The CONNECT registry is undertaken to evaluate clinical and technical success as well as safety and feasibility of the in-house manufactured multi-branch stent graft systems used in endovascular treatment of thoracoabdominal aortic aneurysms that cannot be treated with commercially available devices. Primary endpoint of this study is the rate of patients with stable or decreasing aneurysm size at 12 months follow-up.

**Patient Population:** A total number of 40 male and female patients with asymptomatic thoracoabdominal aortic aneurysm, treated with a contract manufactured multi-branch stent graft system.

**Study status:**

- Start of enrollment: September 2017
- Patients enrolled: 37 (recruitment is ongoing)
- Participating centers: 7

**Peripheral covered stents used to bridge the side branches of the multi-branch stent graft with the visceral arteries (celiac trunk, superior mesenteric artery) and renal arteries:**

- **Total number:** **161**
  - E-ventus BX (JOTEC) 89
  - Viabahn VBX / Endoprosthesis(Gore) 44
  - Lifestream (Bard) 13
  - Advanta V12 (Getinge) 9
  - B-Graft (Bently) 2
  - Unknown 4

**Observed incidents related to the peripheral covered stents:**

- None