

Teleflex Medical  
 IDA Business & Technology Park,  
 Dublin Road, Athlone  
 Westmeath, Ireland

15 October 2020

**URGENT – FIELD SAFETY NOTICE**

<b>Type of Action</b>	<b>Recall</b>
<b>Teleflex Reference</b>	EIF-000427-01
<b>Commercial Name</b>	Lasertube (Rubber) Laser resistant tracheal tube, cuffed; Endotracheal tube for laser surgery
<b>Product Code/Lot Number</b>	Refer to Appendix 2

Dear Customer,

**Details of affected devices**

Teleflex has initiated a voluntary Field Safety Corrective Action (FSCA) for the above listed product, refer to Appendix 2 for a list of product codes and lots impacted. This recall is an amendment to a recent recall reference EIF-000427. The expansion is necessary to include additional lots based on the outcome of the root cause investigation carried out. If you did not receive a copy of the original recall letter, it is because we have identified that you had only received products from the expanded scope of products involved, and did not receive any affected lots of the Rüsich Lasertubes that were included in the original recall scope. Regardless, if you had originally responded to the initial notification, or are just receiving notification for the first time, Teleflex requests that you complete and sign the attached response form confirming that you had been notified of this action and took the required action.

**Description of the problem & immediate actions required**

Teleflex is initiating a voluntary recall for the above-mentioned products due to reports indicating that the laser guard foil partially separated and/or slightly detached at the edges in the presence of moisture. If the defect is present and is not recognized prior to use, adverse health consequences may result from the use of the device during laser therapy in the trachea and larynx including potential for mucosal cell trauma/bleeding, scarring, infection and pain. No patient injuries have been reported.

Depending on your device location please adhere to the following Action list:

Device location	Action List Number
Medical facilities (hospitals, medical staff, etc.)	<b>1</b>
Distributors	<b>2</b>

**Action list number 1 – Medical facilities**

1. We request that you check your inventory for product within the scope of this FSCA. Users should cease use and distribution of impacted product and quarantine immediately.
2. If you do not have stock in scope of this FSCA mark the according checkbox on the Acknowledgement Form (Appendix 1) and return the form to the fax number or e-Mail address mentioned below.

3. If you do have stock in scope of this FSCA, mark the according checkbox on the Acknowledgement Form (Appendix 1) and contact customer service by calling the phone number mentioned below. Customer service will issue you with a return number. Write the return number into the respective field in the Acknowledgement Form and return this form immediately to Customer Service.
4. Teleflex (or your local dealer) will issue a credit note upon receipt of the returned affected product.

#### **Action list number 2 – Distributors**

1. Provide this field safety notice to all customers who have received product in scope of this FSCA. Your customer is then required to complete the acknowledgement form and return to you.
2. We request that you check your inventory for product within the scope of this FSCA. Cease use and distribution of impacted product and quarantine immediately. You may then return all product in scope, refer to Appendix 2 for the list of impacted codes & lots, to Teleflex.
3. As a distributor, you are then required to confirm to Teleflex that you have completed the field activity outlined above. Upon completion of your actions, please forward the completed Acknowledgement Form to Customer Service.
4. Please be aware that all European Economic Area/Switzerland (EEA/CH) and Turkey Member State Competent Authorities in which Teleflex distribute directly will be notified by Teleflex.
5. If you have further distributed product outside of your country, please notify Teleflex by return email to the e-Mail address below.
6. If you are a distributor and/or have a reporting responsibility within or outside the EEA/CH/TR region, please notify your local Competent Authority of this action. Please forward the notification and all communication with your local competent authority to Teleflex.

#### **Teleflex**

Teleflex informs all customers, employees of Teleflex and distributors of this Field Safety Corrective Action.

#### **Transmission of this Field Safety Notice**

This notice should be passed on to all persons who need to be aware within your organisation or to any organisation where the potentially affected devices have been transferred. Please consider end users, clinicians, risk managers, supply chain/distribution centres etc. in the circulation of this notice. Maintain awareness of this notice until all required actions have been completed in your organisation.



**Contact reference person**

Should you require any further information or support concerning this issue, please contact:

**Customer Service:**

**Contact:** Nicole Morawiec

**FAX:** +41 (0) 31 818 40 93

**Telephone:** +41 (0) 31 818 40 90

**Email:** info.ch@teleflex.com

Please be advised that all Economic Area/Switzerland (EEA/CH) and Turkey Member State Competent Authorities to which Teleflex distribute directly will be notified by Teleflex. Teleflex is committed to providing high quality, safe and effective products. We sincerely apologise for any inconvenience this action may cause your operations. If you have any other questions, feel free to contact your local sales representative or Customer Service.

*For and on behalf of Teleflex,*

*Padraig Hegarty*

*Padraig Hegarty VP, Global QA (Manufacturing)*

**FIELD SAFETY CORRECTIVE ACTION**  
**ACKNOWLEDGEMENT FORM**

**PRODUCT FIELD ACTION BY TELEFLEX – IMMEDIATE ATTENTION REQUIRED**

Ref. EIF-000427-01

**RETURN COMPLETED FORM IMMEDIATELY TO:**

**FAX: +41 (0) 31 818 40 93**

**Email: info.ch@teleflex.com**

<input type="checkbox"/> We confirm receipt of this FSN and completed the required actions contained therein. We confirm that our inventory does <b>NOT</b> include products affected by this Field Action.	<input type="checkbox"/> We confirm receipt of this FSN and completed the required actions contained therein. We confirm our inventory <b>DOES</b> include products affected by this Field Action. The use and further distribution of the affected products is stopped. All products are put on hold and the amount below will be returned. <b>Return Authorisation No</b> _____
---	--

**PLEASE PRINT PRODUCT QUANTITY NUMBERS CLEARLY**

PRODUCT NUMBER	LOT NUMBER	QUANTITY (Returning)

- Include a copy of the **completed Acknowledgement Form** in the returns package with the returned units
- Ensure the **RAN number** is clearly visible on the returns package
- Please label returns as **"Field Safety Returns"**

**Complete this Acknowledgement form and return immediately by using fax number or e-Mail address above.**

<b>INSTITUTION NAME (EG NAME OF HOSPITAL, HEALTH CARE ORGANISATION)</b>	
<b>INSITIUTION ADDRESS</b>	<b>Phone/FAX</b>
<b>FORM COMPLETED BY:</b>	<b>Stamp</b>
<b>PRINT NAME:</b> _____	
<b>SIGNATURE:</b> _____	
<b>DATE</b>	

**Appendix 2: EIF-000427-01 Scope of Product and Lots**

<b>Product Code</b>	<b>Lot Number</b>	<b>Product Code</b>	<b>Lot Number</b>
<b>102004-000040</b>	19491	<b>102004-000050 cont.</b>	19041
	19501		19051
	19511		19061
	20021		19071
	20041		19081
	20091		19091
	18451		19101
	18501		19111
	19031		19121
	19041		19141
	19051		19151
	19071		19161
	19081		19171
	19091		19181
	19101		19201
	19111		19211
	19151		19231
	19181		19241
	19191		19251
	19231		19281
	19241		19291
	19281		19301
	19301		19311
	19371		19381
	19391		19391
	19401		19401
19431	19421		
<b>102004-000050</b>	19471	<b>102004-000060</b>	19441
	19501		19451
	19511		19481
	20021		19491
	20031		19501
	20041		20081
	20051		20091
	20061		20101
	20071		20111
	20101		20121
	20111		20131
	20121		18461
	18451		18491
	18461		18501
	18491		19021
	18501		19031
	18511		19041
	19021		19051
	19031		19061

Product Code	Lot Number	Product Code	Lot Number
<b>102004-000060</b> cont.	19071	<b>102004-000070</b> cont.	19091
	19081		19101
	19091		19111
	19101		19121
	19111		19141
	19121		19171
	19131		19191
	19141		19211
	19151		19231
	19161		19261
	19171		19271
	19191		19281
	19211		19291
	19221		19301
	19241		19311
	19251		19351
	19261		19361
	19271		19381
	19281		19391
	19291		19401
	19301		19411
	19311		19431
	19351		19441
	19381		19451
	19391		19481
	19401		19501
	19411		19511
	19421		20041
	19431		20111
	<b>102004-000070</b>		19451
19461		18461	
19471		19021	
19501		19031	
19511		19061	
20021		19101	
20031		19131	
20051		19161	
20091		19261	
20111		19291	
18461		19361	
18491		19401	
18511		19421	
19031			
19061			
19071			
19081			