23. März 2021

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Ricardo Escolá
Senior Manager Quality,
Regulatory & Compliance

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The screenshot displays the 'User Defined Protocol Overview' section of the Ortho Clinical Diagnostics software. The table lists various user-defined protocols (UDP) with columns for Test, Original Test, Results, Last Modified, and Last Modifier. The 'UDP_4 ABO' row is highlighted with a red circle. Below the table, the 'Modify' button in the bottom navigation bar is also circled in red.

| Test | Original Test | Results | Last Modified | Last Modifier |
|----------------------|--------------------------|-----------|-----------------------|---------------|
| UDP_Fya antitoxin | 4 RAS Fya IgG | Fya | 2/23/2021 11:32:20 AM | Admin SOL |
| UDP_Cw | 08 RAS Lea Rvs | Anti-Cw | 2/23/2021 4:34:46 PM | Admin SOL |
| UDP_08 RAS Fya IgG | 08 RAS Fya IgG | Fya | 2/23/2021 11:51:27 AM | Admin SOL |
| UDP_4 ABO(FWD/RVS)/R | 4 ABO(FWD/RVS)/Rh-00 | ABO Rh | 2/23/2021 11:54:26 AM | Admin SOL |
| UDP_DS BV IAT IgG 08 | DS BV IAT IgG 08 37 RRBC | DIISeries | 2/23/2021 4:24:19 PM | Admin SOL |
| UDP_08 RAS Lea Rvs | 08 RAS Lea Rvs | Lea | 3/1/2021 5:00:35 PM | Admin SOL |
| UDP_4 ABO | 4 ABO(FWD/RVS)/Rh-00 | ABO Rh | 3/3/2021 3:30:00 PM | Admin SOL |

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The screenshot shows the 'Ortho Clinical Diagnostics Admin SOL' interface. At the top, there is a navigation bar with icons for Home, Resources, Samples, Results, Errors, QC, Setup, Software, Maintenance, and Diagnostics. The user is logged in as 'Admin SOL' on '03/03/2021 15:32'. The main area is divided into a left sidebar, a central input area, and a right sidebar. The left sidebar contains 'Test Name*' (UDP_4 ABO), 'Assays*' (circled in red), 'Cassettes', and 'Result Interpretation'. The central input area has a text field with 'UDP_4 ABO' and a virtual keyboard below it. The right sidebar contains a menu with 'General', 'Users', 'Testing' (highlighted), 'Results', 'System', 'Interfaces', 'Maintenance', and 'User Defined Reports'. At the bottom, there are 'Cancel' and 'Save' buttons, and a navigation bar with 'Back', 'Modify', 'Create Test', 'Delete Selected', 'Search', 'Help', and 'Resume Processing'.

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The screenshot shows the 'Ortho Clinical Diagnostics Admin SOL' interface. The left sidebar now shows 'Assays*' with a list of assays: 'Column 1*' (Anti-A), 'Column 2*' (Anti-B), 'Column 3*' (Anti-D), 'Column 4*' (Ctrl), 'Column 5*' (AI-Cells, circled in red), and 'Column 6*' (B-Cells). The central input area is empty. The right sidebar and bottom navigation bar are the same as in the previous screenshot.

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

The screenshot shows the Admin SOL interface for Ortho Clinical Diagnostics. The top navigation bar includes Home, Resources, Samples, Results, Errors, QC, Setup, Software, Maintenance, and Diagnostics. The main content area displays a list of modules with their descriptions. The 'Testing' menu item in the right-hand sidebar is highlighted. At the bottom right, the 'Stop Processing' button is circled in red.

| Module | Description |
|------------------------|--|
| Profiles | Grouping Tests into Profiles |
| Test Settings | Global Settings for Automatic Acceptance of Test Results |
| Reagent Kits | Register User-Defined Reagent Kit and Link to an Ortho Reagent Kit |
| QC Settings | Configure Global QC Settings |
| User Defined Protocols | User Defined Protocol Management: Enabled |

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The screenshot shows the Admin SOL interface for Ortho Clinical Diagnostics. The top navigation bar is the same as in the previous screenshot. The main content area displays the same list of modules. The 'Test Settings' row is circled in red. At the bottom center, the 'Show Details' button is circled in red.

| Module | Description |
|------------------------|--|
| Profiles | Grouping Tests into Profiles |
| Test Settings | Global Settings for Automatic Acceptance of Test Results |
| Reagent Kits | Register User-Defined Reagent Kit and Link to an Ortho Reagent Kit |
| QC Settings | Configure Global QC Settings |
| User Defined Protocols | User Defined Protocol Management: Enabled |

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The screenshot shows the 'Automatic Result Acceptance' configuration page in the Admin SOL interface. The page is divided into two main sections: 'Automatic Result Acceptance' and 'Manual Result Acceptance'. The 'Automatic Result Acceptance' section contains a table of rules and their current settings. The 'Require Manual Review of UDP' rule is highlighted with a red circle. The 'Manual Result Acceptance' section contains a single rule, 'Enforce Different Accept User', which is set to 'Yes'. A sidebar on the right lists various system settings categories, with 'Testing' selected. At the bottom, the 'Edit Automatic' button is also highlighted with a red circle.

| Auto Acceptance/Manual Review Rules (Ignored if Enforce Manual Review is enabled) | |
|---|-----|
| All Maintenance Tasks Must Be Completed Successfully | No |
| All Column Reactions Must Be Without Errors | Yes |
| Valid Test Result Interpretation (includes nondiscrepant columns) | Yes |
| Above/Below Positive Reaction Threshold | No |
| Enforce System must be Completed | No |
| Require Manual Review of UDP | No |

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