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Ricardo Escolá
Senior Manager Quality,
Regulatory & Compliance

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The screenshot displays the 'User Defined Protocol Overview' table in the Ortho Clinical Diagnostics software. The table lists various UDP tests with columns for Test, Original Test, Results, Last Modified, and Last Modifier. The 'UDP_4 ABO' row is highlighted with a red circle. Below the table, the 'Modify' button in the bottom navigation bar is also circled in red.

Test	Original Test	Results	Last Modified	Last Modifier
UDP_Fya antitoxin	4 RAS Fya IgG	Fya	2/23/2021 11:32:20 AM	Admin SOL
UDP_Cw	08 RAS Lea Rvs	Anti-Cw	2/23/2021 4:34:46 PM	Admin SOL
UDP_08 RAS Fya IgG	08 RAS Fya IgG	Fya	2/23/2021 11:51:27 AM	Admin SOL
UDP_4 ABO(FWD/RVS)/R	4 ABO(FWD/RVS)/Rh-00	ABO Rh	2/23/2021 11:54:26 AM	Admin SOL
UDP_DS BV IAT IgG 08	DS BV IAT IgG 08 37 RRBC	DIISeries	2/23/2021 4:24:19 PM	Admin SOL
UDP_08 RAS Lea Rvs	08 RAS Lea Rvs	Lea	3/1/2021 5:00:35 PM	Admin SOL
UDP_4 ABO	4 ABO(FWD/RVS)/Rh-00	ABO Rh	3/3/2021 3:30:00 PM	Admin SOL

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The screenshot shows the 'Ortho Clinical Diagnostics Admin SOL' interface. The top navigation bar includes icons for Home, Resources, Samples, Results, Errors, QC, Setup, Software, Maintenance, and Diagnostics. The main area is divided into a left sidebar and a central workspace. The sidebar contains 'Test Name*' (UDP_4 ABO), 'Assays*' (circled in red), 'Cassettes', and 'Result Interpretation'. The central workspace shows a text input field with 'UDP_4 ABO' and a virtual keyboard. A right-hand menu lists 'General', 'Users', 'Testing' (highlighted), 'Results', 'System', 'Interfaces', 'Maintenance', and 'User Defined Reports'. At the bottom, there are 'Cancel' and 'Save' buttons, and a navigation bar with 'Back', 'Modify', 'Create Test', 'Delete Selected', 'Search', 'Help', and 'Resume Processing'.

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The screenshot shows the 'Ortho Clinical Diagnostics Admin SOL' interface with the 'Assays*' section expanded. The 'Assays*' section contains a list of columns: 'Column 1*' (Anti-A), 'Column 2*' (Anti-B), 'Column 3*' (Anti-D), 'Column 4*' (Ctrl), 'Column 5*' (A1-Cells, circled in red), and 'Column 6*' (B-Cells). The rest of the interface, including the top navigation bar, right-hand menu, and bottom navigation bar, is identical to the previous screenshot.

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The screenshot shows the Admin SOL interface for Ortho Clinical Diagnostics. The top navigation bar includes Home, Resources, Samples, Results, Errors, QC, Setup, Software, Maintenance, and Diagnostics. The 'Setup' menu is expanded, showing a list of modules: Profiles, Test Settings, Reagent Kits, QC Settings, and User Defined Protocols. The 'Test Settings' module is selected. On the right side, there is a sidebar menu with options: General, Users, Testing, Results, System, Interfaces, Maintenance, and User Defined Reports. At the bottom right, there are three buttons: Search, Help, and Stop Processing. The 'Stop Processing' button is circled in red.

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The screenshot shows the Admin SOL interface for Ortho Clinical Diagnostics. The top navigation bar is the same as in the previous screenshot. The 'Test Settings' module is selected, and the 'Show Details' button is circled in red. The sidebar menu is also visible on the right side. At the bottom right, there are three buttons: Search, Help, and Resume Processing.

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The screenshot shows the 'Automatic Result Acceptance' configuration page in the Admin SQL interface. The page is divided into two main sections: 'Automatic Result Acceptance' and 'Manual Result Acceptance'. The 'Automatic Result Acceptance' section contains a table of rules and their current settings. The 'Require Manual Review of UDP' rule is highlighted with a red circle. The 'Edit Automatic Result' button at the bottom of the page is also circled in red.

Auto Acceptance/Manual Review Rules (ignored if Enforce Manual Review is enabled)	
All Maintenance Tasks Must Be Completed Successfully	No
All Column Reactions Must Be Without Errors	Yes
Valid Test Result Interpretation (includes nondiscrepant columns)	Yes
Above/Below Positive Reaction Threshold	No
Entire Order must be Completed	No
Require Manual Review of UDP	No

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