23 Marzo 2021

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Ricardo Escolá
Senior Manager Quality,
Regulatory & Compliance

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The screenshot displays the 'User Defined Protocol Overview' section of the Ortho Clinical Diagnostics software. The interface includes a top navigation bar with icons for Home, Resources, Samples, Results, Errors, QC, Setup, Software, Maintenance, and Diagnostics. The user is logged in as 'Admin SOL' on '03/03/2021 15:30'. The main table lists various user-defined protocols with columns for Test, Original Test, Results, Last Modified, and Last Modifier. The 'UDP_4 ABO' row is highlighted with a red circle. Below the table, a bottom navigation bar contains icons for Back, Modify (circled in red), Create Test, and Delete Selected. A right-hand sidebar menu lists options: General, Users, Testing (selected), Results, System, Interfaces, Maintenance, and User Defined Reports. Search, Help, and Resume Processing buttons are also present at the bottom right.

Test	Original Test	Results	Last Modified	Last Modifier
UDP_Fya antitoxin	4 RAS Fya IgG	Fya	2/23/2021 11:32:20 AM	Admin SOL
UDP_Cw	08 RAS Lea Rvs	Anti-Cw	2/23/2021 4:34:46 PM	Admin SOL
UDP_08 RAS Fya IgG	08 RAS Fya IgG	Fya	2/23/2021 11:51:27 AM	Admin SOL
UDP_4 ABO(FWD/RVS)/R	4 ABO(FWD/RVS)/Rh-00	ABO Rh	2/23/2021 11:54:26 AM	Admin SOL
UDP_DS BV IAT IgG 08	DS BV IAT IgG 08 37 RRBC	DilSeries	2/23/2021 4:24:19 PM	Admin SOL
UDP_08 RAS Lea Rvs	08 RAS Lea Rvs	Lea	3/1/2021 5:00:35 PM	Admin SOL
UDP_4 ABO	4 ABO(FWD/RVS)/Rh-00	ABO Rh	3/3/2021 3:30:00 PM	Admin SOL

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The screenshot shows the 'Admin SOL' interface for 'Ortho Clinical Diagnostics'. The top navigation bar includes icons for Home, Resources, Samples, Results, Errors, QC, Setup, Software, Maintenance, and Diagnostics. The main content area is divided into three sections: 'Test Name*' (UDP_4 ABO), 'Assays*', and 'Cassettes'. The 'Assays*' section is highlighted with a red circle. A virtual keyboard is overlaid on the screen, and a 'Cancel' and 'Save' button are visible at the bottom right. The right sidebar contains a menu with options: General, Users, Testing (highlighted), Results, System, Interfaces, Maintenance, and User Defined Reports.

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The screenshot shows the 'Admin SOL' interface for 'Ortho Clinical Diagnostics'. The top navigation bar is the same as in the previous screenshot. The main content area shows the 'Test Name*' (UDP_4 ABO) and the 'Assays*' section expanded to show a list of columns. The 'Column 5*' row is highlighted with a red circle. The 'Cassettes' and 'Result Interpretation' sections are also visible. The right sidebar and bottom navigation bar are the same as in the previous screenshot.

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The screenshot shows the Admin SOL interface for Ortho Clinical Diagnostics. The top navigation bar includes Home, Resources, Samples, Results, Errors, QC, Setup, Software, Maintenance, and Diagnostics. The main content area displays a table of modules with the following data:

Module	Description
Profiles	Grouping Tests into Profiles
Test Settings	Global Settings for Automatic Acceptance of Test Results
Reagent Kits	Register User-Defined Reagent Kit and Link to an Ortho Reagent Kit
QC Settings	Configure Global QC Settings
User Defined Protocols	User Defined Protocol Management: Enabled

On the right side, there is a vertical menu with options: General, Users, Testing, Results, System, Interfaces, Maintenance, and User Defined Reports. At the bottom right, there are buttons for Search, Help, and Stop Processing (highlighted with a red circle).

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The screenshot shows the Admin SOL interface with the 'Test Settings' row selected in the table. The 'Show Details' button at the bottom is highlighted with a red circle.

Module	Description
Profiles	Grouping Tests into Profiles
Test Settings	Global Settings for Automatic Acceptance of Test Results
Reagent Kits	Register User-Defined Reagent Kit and Link to an Ortho Reagent Kit
QC Settings	Configure Global QC Settings
User Defined Protocols	User Defined Protocol Management: Enabled

At the bottom, there are buttons for Disable User-Defined Pr..., Show Details (highlighted with a red circle), Search, Help, and Resume Processing.

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