

	Field Safety Notice FSN C.G.M. Divisione Medicale Meta Ref. no. 2021_001	DATE 03-05-2021 REV. 00 PAG. 1 di 7
---	---	--

Urgent Field Safety Notice
Device Commercial Names as provided in Appendix 1

To the kind attention of:

List of will be part of the FSN in the different destination countries

- Customer Name,
- Address
- Postal code, City name
- e-mail
- Telephone

	Field Safety Notice FSN C.G.M. Divisione Medica Meta Ref. no. 2021_001	DATE 03-05-2021 REV. 00 PAG. 2 di 7
---	---	---

Urgent Field Safety Notice (FSN)
Device Names as provided in Appendix 1

This letter contains important information which require your **immediate attention**.

1. Information on Affected Devices*	
1	1. Device Type(s)*
.	See Appendix 01
1	2. Commercial name(s)
.	See Appendix 01
1	3. Unique Device Identifier(s) (UDI-DI)
.	Not available
1	4. Primary clinical purpose of device(s)*
.	See Appendix 01
1	5. Device Model/Catalogue/part number(s)*
.	See Appendix 01
1	6. Software version
.	Not relevant
1	7. Affected serial or lot number range
.	See Appendix 01
1	8. Associated devices
.	Unknown.

2 Reason for Field Safety Corrective Action (FSCA)*	
2.	1. Description of the product problem C.G.M. Divisione Medica Meta is the legal manufacturer of the following devices: - sterile scraper for use as a collecting bone flakes in oral surgical operations. - Set for Uterine Suction with tube and canula - membrane fixation tacks for oral surgical operations - Umbilical Cord Clamp - Closed Circuit Urine Bag - Amniotic Membrane Perforator - Magnetic Mat for Surgical Instrument Those products are supplied to the market in sterile status, following the Etylene Oxide sterilization process performed overtime by Steril Milano Srl, one of the largest EO sterilization service providers in Italy. C.G.M. Divisione Medica Meta has become aware of sterilization issues notified by the contract sterilizer Steril Milano, with potential impact on efficacy of the Ethylene Oxide (EtO) sterilization processes at the contract sterilizer Steril Milano and sterile status of the devices placed on the market. According to our investigation, we have identified certain batches for which we are unable to guarantee the primary sterility, even though, for the time being, based on our



Field Safety Notice
FSN C.G.M. Divisione Medicale
Meta Ref. no. 2021_001

DATE 03-05-2021
REV. 00
PAG. 3 di 7

	test results, we have no evidence of non-sterile status of the goods. Those batches are listed in the attached Appendix 1 “List of Impacted Batches” .
2.	2. Hazard giving rise to the FSCA The falsification of relevant data, especially linked to the preconditioning cycle and the sterilization cycle, could play a crucial role in the functionality and effectiveness of the devices' sterilisation processes. As specified in the risk analysis of the technical files, the ineffective sterilization of the devices listed above, could have consequences for patient's health with potential side effects linked to unsterile products, patient infection and worsening of their health conditions. C.G.M. Medical Division META didn't receive any notification of adverse events or serious patient harm associated with this safety corrective action. Even in the past years, our company didn't receive claims for adverse events referred to the millions of devices sold. Based on the following reasons, no specific patient follow-up activities are required for the product used: 1) a preventive antibiotic therapy is prescribed before surgery procedures, 2) low level of microbial contamination of the products - detected by periodic Bioburden Tests - guarantee good disinfection of devices, 3) no adverse events occurred for over 2 Million products sold in over 20 years, 4) Several Sterility Tests performed on devices sterilised with batches affected by this FSN resulted "sterile". 5) The sterilization colour change indicators are always checked during incoming controls and no deviation was never detected. All the products identified as potentially not sterile delivered to your Company are listed in Appendix 1 “List of Impacted Batches of the present FSN”.
2.	3. Probability of problem arising All analysis performed in the past shown that the products were correctly sterile. Right now, further analysis is ongoing. Therefore we can't define a percentage, yet.
2.	4. Predicted risk to patient/users From the Health Hazard Evaluation of our devices, exposure to microbiological contamination could lead to bacterial infection and worsening of the patient health conditions.
2.	5. Further information to help characterise the problem NA
2.	6. Background on Issue NA
2.	7. Other information relevant to FSCA NA



Field Safety Notice
FSN C.G.M. Divisione Medica
Meta Ref. no. 2021_001

DATE 03-05-2021
 REV. 00
 PAG. 4 di 7

3. Type of Action to mitigate the risk											
3.	1. Action To Be Taken by the User ☒ Identify Device ☒ Quarantine Device ☒ Return Device, when requested by C.G.M. Divisione Medica Meta ☒ Destroy Device, when requested by C.G.M. Divisione Medica Meta ☐ On-site device modification/inspection ☐ Follow patient management recommendations ☐ Take note of amendment/reinforcement of Instructions For Use (IFU) ☐ Other ☐ None Once received this official notification, in order to prevent potential impact of the medical procedure, each user shall: 1) Identify and segregate all items listed in Appendix 01, still available at their premises, 2) Translate FSCA and Acknowledgment letter for Healthcare Facilities, provided in Annex 03, in your national languages, 3) Fill in the acknowledgment letter provided in the Appendix 02, including the number of segregate devices and returned devices, 4) Within 5 working days from receiving the official notification, return to C.G.M. Divisione Medica Meta premises, E.Villa n.7, I-42124 Reggio Emilia (RE) – Italy, or destroy all the segregated devices, according to instruction provided by META, As required, we have provided this notification to the relevant Regulatory Agencies of the countries where the devices have been distributed. Please refer to your local sales agent for any further information you may need or, in alternative, contact directly C.G.M. Divisione Medica Meta customer service at telephone number +39 0522 502311 or mail helpdesk@metahosp.com										
3.	2. By when should the action be completed? Within 5 (five) calendar days from the issue date <table border="1" style="width: 100%; border-collapse: collapse; margin-top: 10px;"> <thead> <tr> <th style="width: 10%;">ID#</th> <th style="width: 60%;">Actions description</th> <th style="width: 30%;">By when</th> </tr> </thead> <tbody> <tr> <td style="text-align: center;">1</td> <td>Identify and segregate all items listed in Appendix 01, still available at your premises</td> <td>Immediately or within 1 calendar day</td> </tr> <tr> <td style="text-align: center;">2</td> <td>Translate FSCA and Acknowledgment letter for Healthcare Facilities, provided in Annex 03, in your national languages</td> <td>Immediately or within 3 calendar day</td> </tr> </tbody> </table>		ID#	Actions description	By when	1	Identify and segregate all items listed in Appendix 01, still available at your premises	Immediately or within 1 calendar day	2	Translate FSCA and Acknowledgment letter for Healthcare Facilities, provided in Annex 03, in your national languages	Immediately or within 3 calendar day
ID#	Actions description	By when									
1	Identify and segregate all items listed in Appendix 01, still available at your premises	Immediately or within 1 calendar day									
2	Translate FSCA and Acknowledgment letter for Healthcare Facilities, provided in Annex 03, in your national languages	Immediately or within 3 calendar day									

 META	Field Safety Notice FSN C.G.M. Divisione Medica Meta Ref. no. 2021_001	DATE 03-05-2021 REV. 00 PAG. 5 di 7
--	---	--

	<table border="1"> <tr> <td style="text-align: center;">3</td> <td>Fill the Acknowledgment Letter provided in the Appendix 02, including the number of received devices, used or sold devices, remaining and segregated devices.</td> <td>Within 7 calendar days from the receipt of the present communication</td> </tr> <tr> <td style="text-align: center;">4</td> <td>Return to C.G.M. Divisione Medica Meta premises, Via E.Villa n.7, I-42124 Reggio Emilia (RE) – Italy, or destroy all the segregated devices, according to instruction provided by META</td> <td>Within 30 calendar days from receiving the official notification</td> </tr> </table>	3	Fill the Acknowledgment Letter provided in the Appendix 02, including the number of received devices, used or sold devices, remaining and segregated devices.	Within 7 calendar days from the receipt of the present communication	4	Return to C.G.M. Divisione Medica Meta premises, Via E.Villa n.7, I-42124 Reggio Emilia (RE) – Italy, or destroy all the segregated devices, according to instruction provided by META	Within 30 calendar days from receiving the official notification
3	Fill the Acknowledgment Letter provided in the Appendix 02, including the number of received devices, used or sold devices, remaining and segregated devices.	Within 7 calendar days from the receipt of the present communication					
4	Return to C.G.M. Divisione Medica Meta premises, Via E.Villa n.7, I-42124 Reggio Emilia (RE) – Italy, or destroy all the segregated devices, according to instruction provided by META	Within 30 calendar days from receiving the official notification					
3.	3. Particular considerations for: N/A						
3.	4. Is customer Reply Required? See Acknowledgment Letter in Appendix 02, to be returned within 7 calendar days from the issue date.						
3.	5. Action Being Taken by the Manufacturer ☒ Product Removal ☐ Software upgrade ☐ Other Device re-working ☐ On-site device modification/inspection ☐ IFU or labelling change ☐ None Based on the evaluation and sterility test performed, as conservative approach and a protective measure to maintain patient health, we decided to replace the devices listed in Appendix 01. C.G.M. Divisione Medica Meta has sent a Field Safety Notice to all affected customers. The Field Safety Notice identifies the problem, the affected products, the risk factors and the actions that must be taken by the users and distributors.						
3	6. By when should the action be completed?	Before 30 calendar days from the issue date					
3.	7. Is the FSN required to be communicated to the patient /lay user?	No					
3	8. If yes, has manufacturer provided additional information suitable for the patient/lay user in a patient/lay or non-professional user information letter/sheet? No Not appended to this FSN						

	4. General Information*	
4.	1. FSN Type*	New

	Field Safety Notice FSN C.G.M. Divisione Medica Meta Ref. no. 2021_001	DATE 03-05-2021 REV. 00 PAG. 6 di 7
---	---	---

4.	2. For updated FSN, reference number and date of previous FSN	NA
4.	3. For Updated FSN, key new information as follows:	
	NA	
4.	4. Further advice or information already expected in follow-up FSN?	No
4	5. If follow-up FSN expected, what is the further advice expected to relate to:	
	NA	
4	6. Anticipated timescale for follow-up FSN	NA
4.	7. Manufacturer information (For contact details of local representative refer to page 1 of this FSN)	
	a. Company Name	C.G.M. Divisione Medica Meta S.p.A.
	b. Address	Via E.Villa n.7, I-42124 Reggio Emilia (RE) – Italy
	c. Website address	http://www.metahosp.com/
4.	8. The Competent (Regulatory) Authority of your country has been informed about this communication to customers. Yes	
4.	9. List of attachments/appendices:	1. <u>Appendix 01: List of affected devices</u> 2. <u>Appendix 02: Acknowledgment letter for Distributor</u> 3. <u>Appendix 03: FSCA and Acknowledgment letter for Healthcare Facilities</u>
4.	4. Name/Signature	Insert Name and Title here and signature below

	Transmission of this Field Safety Notice
	This notice needs to be passed on all those who need to be aware within your organisation or to any organisation where the potentially affected devices have been transferred. (As appropriate) Please transfer this notice to other organisations and to all users on which this action has an impact. (As appropriate) Please maintain awareness on this notice and resulting action for an appropriate period to ensure effectiveness of the corrective action.

	Field Safety Notice FSN C.G.M. Divisione Medicale Meta Ref. no. 2021_001	DATE 03-05-2021 REV. 00 PAG. 7 di 7
---	---	--

	Please report all device-related incidents to the manufacturer, distributor or local representative, and the national Competent Authority if appropriate, as this provides important feedback.
--	---



Urgent Field Safety Notice FSN-01-2021

ANNEX 1 rev 1

BATCHES LIST OF MEDICAL DEVICES INVOLVED IN SWITZERLAND -MARKET WITH EXPIRY DATE STILL VALID

REF	Class	DEVICE	BATCH	CUSTOMER	q.ty	Exp
3987	IIa	SAFESCRAPER TWIST CURVE	21-26620	HEICO DENT GmbH	-60	2023-08
3987	IIa	SAFESCRAPER TWIST CURVE	10-10719	HEICO DENT GmbH	-60	2022-03
3987	IIa	SAFESCRAPER TWIST CURVE	7-06519	HEICO DENT GmbH	-60	2022-02
3987	IIa	SAFESCRAPER TWIST CURVE	17-19818	HEICO DENT GmbH	-60	2021-06
3987	IIa	SAFESCRAPER TWIST CURVE	14-15018	HEICO DENT GmbH	-60	2021-05

REF	Class	DEVICE	BATCH	CUSTOMER	q.ty	Exp
4049	IIa	MICROSS	11-10720	HEICO DENT GmbH	-80	2023-04
4049	IIa	MICROSS	10-10619	HEICO DENT GmbH	-80	2022-03

REF	Class	DEVICE	BATCH	CUSTOMER	q.ty	Exp
3598	IIa	SAFESCRAPER TWIST	10-10719	Karr Dental Ag	-60	2022-03

REF	Class	DEVICE	BATCH	CUSTOMER	q.ty	Exp
4049	IIa	MICROSS	1-00820	Karr Dental Ag	-40	2022-12

REF	Class	DEVICE	BATCH	CUSTOMER	q.ty	Exp
4049	IIa	MICROSS	17-24819	BDS Dental AG	-12	08/2022
4049	IIa	MICROSS	10-10619	BDS Dental AG	-12	03/2022
4049	IIa	MICROSS	24-29118	BDS Dental AG	-10	09/2021

HEICO DENT GmbH
Strahlholz 13
9056, Gais
SWITZERLAND
VAT NUMBER: CHE464.430.084
Tel. 0041717939000
info@heicodent.ch

Karr Dental Ag
Verenastrasse 4b
8832, Wollerau
SWITZERLAND
Tel. +41 44 727 40 07
Sandro Lardo (s.lardo@karrdental.ch)



C.G.M. S.p.A. – Divisione Medica META

Via E. Villa n.7 - 42124 Reggio Emilia (RE) – Italy ☎ phone +39 0522 502305 website: www.metahosp.com
info@metahosp.com

Sede Legale: Via Modena 22/24 – 42015 Correggio (RE) – Italy - P.IVA 00678290354

BDS Dental AG
Römerstrasse 149
Winterthur - 8404
SWITZERLAND
Tel +41 52 397 3020
Jasminka Janik (info@bds-dental.ch)

**C.G.M. S.p.A. – Divisione Medica META**

Via E. Villa n.7 - 42124 Reggio Emilia (RE) – Italy ☎ phone +39 0522 502305 website: www.metahosp.com
info@metahosp.com

Sede Legale: Via Modena 22/24 – 42015 Correggio (RE) – Italy - P.IVA 00678290354

Appendix 02 - Acknowledgment Letter for Distributors

Please read in conjunction with FIELD SAFETY NOTICE FSCA-01-2021 and return completed and signed as soon as possible or within 5 days from its receipt to hepldesk@metahosp.com

Tick all that apply		
☐	I confirm this notice has been read, understood and that all recommended actions have been implemented as required.	Customer/Distributor/Importer to complete, sign or enter N/A
☐	I have checked my internal stock and our clients stock and quarantined all inventories	Customer/Distributor/Importer to complete, sign or enter N/A
☐	I have identified all Healthcare organization and all end users where the devices listed in Annex 1 have been shipped and on which this action has an impact,	Customer/Distributor/Importer to complete, sign or enter N/A
☐	I have informed the identified Healthcare organization and all end users of this FSN	Date of communication:
☐	I have received confirmation of reply from all identified Healthcare organization and all end users	Date of receiving last communication:
☐	I have filled-in the Table 1 , with the number of remaining, segregated and returned/destroy devices to your premise.	Add quantity, Lot/REF/Date Returned (same information as requested by the Customer Reply form)
☐	Our organization has none of the affected devices in inventory	
☐	Our Healthcare clients and end users has none of the affected devices in inventory	

Comments, if any

--

Name of Trust / Organisation :			
Address :			
Postcode :		Country :	
Telephone number :		E-mail address :	
Name of your supplier for this product			
Name, title and signature of person completing this form:			

It is important that your organisation takes the actions detailed in the FSN and confirms that you have received the FSN.

Your organisation's reply is the evidence we need to monitor the progress of the corrective action

[illegible]

HEADED PAPER

Attachment 02 - Acknowledgment Letter for Healthcare Facilities

Please read in conjunction with FIELD SAFETY CORRECTIVE ACTION FSCA-01-2021 and return completed and signed as soon as possible or within 5 days from its receipt to **(EMAIL ADDRESS OF DISTRIBUTOR)**

Tick all that apply		
☐	I confirm this notice has been read, understood and that all recommended actions have been implemented as required.	Customer/Distributor/Importer to complete, sign or enter N/A
☐	I have checked my internal stock and quarantined all inventories	Customer/Distributor/Importer to complete, sign or enter N/A
☐	I have filled-in the Table 1 , with the number of remaining, segregated and returned devices to your premises.	Add quantity, Lot/REF/Date Returned (same information as requested by the Customer Reply form)
☐	Our Healthcare organization has none of the affected devices in inventory	

Comments, if any

--

Name of Trust / Organisation :			
Address :			
Postcode :		Country :	
Telephone number :		E-mail address :	
Name of your supplier for this product			
Name, title and signature of person completing this form:			

It is important that your organisation takes the actions detailed in the FSN and confirms that you have received the FSN.
Your organisation's reply is the evidence we need to monitor the progress of the corrective action

HEADED PAPER

Table 1: Mapping devices in inventory[illegible]