**Urgent Field Safety Notice**

**StimLoc™ Burr Hole Cover**

**Part Numbers 3387S-40, 3389S-40, 3391S-40, 3550S-01, 924256, DB-5000**

Excess Nylon Flash on Specific Lots

Recall

June 2021

Medtronic Reference: FA980

Dear Healthcare Professional,

Medtronic is voluntarily recalling specific lots of the StimLoc™ Burr Hole Cover due to a manufacturing issue resulting in excess nylon flash on the clip. Please review the information contained in this letter, quarantine any affected product (see Attachment A for list of affected lot numbers) in your inventory for return to Medtronic.

**Issue Description:**

Medtronic has identified that specific lots of the StimLoc clips may have excess nylon molding flash. This issue can prevent the clip from snapping into the StimLoc base, which could allow for lead movement or lead damage. This could potentially result in delivery of therapy to the wrong location, inadequate delivery of therapy, or intracranial hemorrhage. As of May 19, 2021, Medtronic has received 26 complaints related to this issue, two of which reported surgical replacement or revision of the lead or Burr Hole Cover. There have been no reports of intracranial hemorrhage.

*Excess flash (red circle)*

**Patient Management Recommendations:**

Previous lead placements with the affected lots where there was successful assembly and fixation with the StimLoc base, clip, and cap do not require additional patient follow-up. If you were able to complete placement, position the lead in the base exit slot, and snap the cap onto the base; the lead will be properly secured as described in the instructions for use (IFU).

As described in the IFU, during any lead placement procedure, avoid significant contact between the lead body and the support clip to prevent changing the position of the lead. As also described in the IFU, ensure the lead body is positioned in the base exit slot before installing the cap to prevent damage to lead. At any time during an implant procedure, if

there is any difficulty with anchoring the lead with kit components, do not use those items, and use a new lead and StimLoc Burr Hole Cover assembly.

**Required Actions:**

1. Please inform your materials department to identify, segregate, and quarantine unused affected units within their inventory. The list of affected lot numbers is included as Attachment A with this letter.
2. Contact Medtronic to return affected units and schedule replacement(s).

**Additional Information:**

Medtronic is communicating this information to the appropriate regulatory agencies.

We regret any inconvenience this may cause. We are committed to patient safety and appreciate your prompt attention to this matter. If you have any questions, please contact your Medtronic Representative.

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