COOK MEDICAL EUROPE LTD.  
O'HALLORAN ROAD  
NATIONAL TECHNOLOGY PARK  
LIMERICK, V94 N8X2, IRELAND  
TEL: +353 61 334440 FAX: +353 61 334441  
WWW.COOKMEDICAL.EU

N. rif. FSN: QCR-110/2021FA0008

FSCA Rif: QCR-110/2021FA0008

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## Customer Service Centres by Country:

Country	Telephone	Fax	E-mail
Austria	+43 179567121	+43 179567126	oe.orders@cookmedical.com
Belgium/French	+32 27001633	+32 27001634	be.orders@cookmedical.com
Belgium/Flemish	+32 27001633	+32 27001634	nl.orders@cookmedical.com
Denmark	+45 38487607	+45 38487608	da.orders@cookmedical.com
Finland	+358 972519996	+358 972519997	fi.orders@cookmedical.com
France	+33 171230269	+33 171230270	fr.orders@cookmedical.com
Germany	+49 6950072804	+49 6950072805	de.orders@cookmedical.com
Hungary	+36 17779199	+36 17779200	hu.orders@cookmedical.com
Ireland	+353 61239252	+353 61239253	ie.orders@cookmedical.com
Italy	+39 0269682853	+39 0269682854	it.orders@cookmedical.com
Netherlands	+31 202013367	+31 202013368	nl.orders@cookmedical.com
Norway	+47 23162968	+47 23162969	no.orders@cookmedical.com
Poland	+48 223060159	+48 223060160	pl.orders@cookmedical.com
Spain	+34 912702691	+34 912702692	es.orders@cookmedical.com
Sweden	+46 858769468	+46 858769467	se.orders@cookmedical.com
Switzerland/French	+41 448009609	+41 448009610	fr.orders@cookmedical.com
Switzerland/Italian	+41 448009609	+41 448009610	it.orders@cookmedical.com
Switzerland/German	+41 448009609	+41 448009610	de.orders@cookmedical.com
United Kingdom	+44 2073654183	+44 2073654184	uk.orders@cookmedical.com

Customer service lines are available (excluding national holidays) Monday through Thursday between the hours of 08:30 CET and 17:30 CET, and from 08:30 CET to 16:00 CET on Fridays.

**Please note: for all other countries, contact the nearest Customer Service centre or your local Cook Medical representative.**

<b>EN</b> <b>2</b>	<b>Small Gauge Ova-Stiff™ Ovum Aspiration Needles, Ova-Stiff™ EchoTip® Single Lumen Ovum Aspiration Needle</b> <b>Ovum Aspiration Needle Single Lumen, Single Lumen Ovum Aspiration Needle</b> Guidance for use, warnings and precautions
<b>BG</b> <b>4</b>	<b>Small Gauge Ova-Stiff™ Ovum Aspiration Needles, Ova-Stiff™ EchoTip® Single Lumen Ovum Aspiration Needle</b> <b>Ovum Aspiration Needle Single Lumen, Single Lumen Ovum Aspiration Needle</b> Указания за употреба, предупреждения и предпазни мерки
<b>CS</b> <b>6</b>	<b>Small Gauge Ova-Stiff™ Ovum Aspiration Needles, Ova-Stiff™ EchoTip® Single Lumen Ovum Aspiration Needle</b> <b>Ovum Aspiration Needle Single Lumen, Single Lumen Ovum Aspiration Needle</b> Návod k použití, varování a bezpečnostní opatření
<b>DA</b> <b>8</b>	<b>Small Gauge Ova-Stiff™ Ovum Aspiration Needles, Ova-Stiff™ EchoTip® Single Lumen Ovum Aspiration Needle</b> <b>Ovum Aspiration Needle Single Lumen, Single Lumen Ovum Aspiration Needle</b> Brugsvejledning, advarsler og forholdsregler
<b>DE</b> <b>10</b>	<b>Small Gauge Ova-Stiff™ Ovum Aspiration Needles, Ova-Stiff™ EchoTip® Single Lumen Ovum Aspiration Needle</b> <b>Ovum Aspiration Needle Single Lumen, Single Lumen Ovum Aspiration Needle</b> Anwendungsanleitung, Warnhinweise und Vorsichtsmaßnahmen
<b>EL</b> <b>12</b>	<b>Small Gauge Ova-Stiff™ Ovum Aspiration Needles, Ova-Stiff™ EchoTip® Single Lumen Ovum Aspiration Needle</b> <b>Ovum Aspiration Needle Single Lumen, Single Lumen Ovum Aspiration Needle</b> Οδηγίες χρήσης, προειδοποιήσεις και προφυλάξεις
<b>ES</b> <b>14</b>	<b>Small Gauge Ova-Stiff™ Ovum Aspiration Needles, Ova-Stiff™ EchoTip® Single Lumen Ovum Aspiration Needle</b> <b>Ovum Aspiration Needle Single Lumen, Single Lumen Ovum Aspiration Needle</b> Guía de uso, advertencias y precauciones
<b>ET</b> <b>16</b>	<b>Small Gauge Ova-Stiff™ Ovum Aspiration Needles, Ova-Stiff™ EchoTip® Single Lumen Ovum Aspiration Needle</b> <b>Ovum Aspiration Needle Single Lumen, Single Lumen Ovum Aspiration Needle</b> Kasutusjuhised, hoiatused ja ettevaatusabinõud
<b>FI</b> <b>18</b>	<b>Small Gauge Ova-Stiff™ Ovum Aspiration Needles, Ova-Stiff™ EchoTip® Single Lumen Ovum Aspiration Needle</b> <b>Ovum Aspiration Needle Single Lumen, Single Lumen Ovum Aspiration Needle</b> Käyttöohje, varoitukset ja varotoimet
<b>FR</b> <b>20</b>	<b>Small Gauge Ova-Stiff™ Ovum Aspiration Needles, Ova-Stiff™ EchoTip® Single Lumen Ovum Aspiration Needle</b> <b>Ovum Aspiration Needle Single Lumen, Single Lumen Ovum Aspiration Needle</b> Directives d'utilisation, avertissements et mises en garde
<b>HR</b> <b>22</b>	<b>Small Gauge Ova-Stiff™ Ovum Aspiration Needles, Ova-Stiff™ EchoTip® Single Lumen Ovum Aspiration Needle</b> <b>Ovum Aspiration Needle Single Lumen, Single Lumen Ovum Aspiration Needle</b> Priručnik za uporabu, upozorenja i mjere opreza
<b>HU</b> <b>24</b>	<b>Small Gauge Ova-Stiff™ Ovum Aspiration Needles, Ova-Stiff™ EchoTip® Single Lumen Ovum Aspiration Needle</b> <b>Ovum Aspiration Needle Single Lumen, Single Lumen Ovum Aspiration Needle</b> Útmutató a használatához, a figyelemzetetésekhez és az óvintézkedésekhez
<b>IT</b> <b>26</b>	<b>Small Gauge Ova-Stiff™ Ovum Aspiration Needles, Ova-Stiff™ EchoTip® Single Lumen Ovum Aspiration Needle</b> <b>Ovum Aspiration Needle Single Lumen, Single Lumen Ovum Aspiration Needle</b> Linee guida per l'uso, avvertenze e precauzioni
<b>LT</b> <b>28</b>	<b>Small Gauge Ova-Stiff™ Ovum Aspiration Needles, Ova-Stiff™ EchoTip® Single Lumen Ovum Aspiration Needle</b> <b>Ovum Aspiration Needle Single Lumen, Single Lumen Ovum Aspiration Needle</b> Naudojimo nurodymai, įspėjimai ir atsargumo priemonės
<b>LV</b> <b>30</b>	<b>Small Gauge Ova-Stiff™ Ovum Aspiration Needles, Ova-Stiff™ EchoTip® Single Lumen Ovum Aspiration Needle</b> <b>Ovum Aspiration Needle Single Lumen, Single Lumen Ovum Aspiration Needle</b> Lietošanas norādījumi, brīdinājumi un piesardzības pasākumi
<b>NL</b> <b>32</b>	<b>Small Gauge Ova-Stiff™ Ovum Aspiration Needles, Ova-Stiff™ EchoTip® Single Lumen Ovum Aspiration Needle</b> <b>Ovum Aspiration Needle Single Lumen, Single Lumen Ovum Aspiration Needle</b> Aanwijzingen voor het gebruik, waarschuwingen en voorzorgsmaatregelen
<b>NO</b> <b>34</b>	<b>Small Gauge Ova-Stiff™ Ovum Aspiration Needles, Ova-Stiff™ EchoTip® Single Lumen Ovum Aspiration Needle</b> <b>Ovum Aspiration Needle Single Lumen, Single Lumen Ovum Aspiration Needle</b> Veiledning for bruk, advarsler og forholdsregler
<b>PL</b> <b>36</b>	<b>Small Gauge Ova-Stiff™ Ovum Aspiration Needles, Ova-Stiff™ EchoTip® Single Lumen Ovum Aspiration Needle</b> <b>Ovum Aspiration Needle Single Lumen, Single Lumen Ovum Aspiration Needle</b> Wskazówki dotyczące użycia, ostrzeżenia i środki ostrożności
<b>PT</b> <b>38</b>	<b>Small Gauge Ova-Stiff™ Ovum Aspiration Needles, Ova-Stiff™ EchoTip® Single Lumen Ovum Aspiration Needle</b> <b>Ovum Aspiration Needle Single Lumen, Single Lumen Ovum Aspiration Needle</b> Orientações de utilização, advertências e precauções
<b>RO</b> <b>40</b>	<b>Small Gauge Ova-Stiff™ Ovum Aspiration Needles, Ova-Stiff™ EchoTip® Single Lumen Ovum Aspiration Needle</b> <b>Ovum Aspiration Needle Single Lumen, Single Lumen Ovum Aspiration Needle</b> Instrucțiuni de utilizare, avertismente și precauții
<b>RU</b> <b>42</b>	<b>Small Gauge Ova-Stiff™ Ovum Aspiration Needles, Ova-Stiff™ EchoTip® Single Lumen Ovum Aspiration Needle</b> <b>Ovum Aspiration Needle Single Lumen, Single Lumen Ovum Aspiration Needle</b> Руководство по эксплуатации, предупреждения и меры предосторожности
<b>SK</b> <b>44</b>	<b>Small Gauge Ova-Stiff™ Ovum Aspiration Needles, Ova-Stiff™ EchoTip® Single Lumen Ovum Aspiration Needle</b> <b>Ovum Aspiration Needle Single Lumen, Single Lumen Ovum Aspiration Needle</b> Pokyny na použitie, varovania a bezpečnostné opatrenia
<b>SL</b> <b>46</b>	<b>Small Gauge Ova-Stiff™ Ovum Aspiration Needles, Ova-Stiff™ EchoTip® Single Lumen Ovum Aspiration Needle</b> <b>Ovum Aspiration Needle Single Lumen, Single Lumen Ovum Aspiration Needle</b> Smernice za uporabo, opozorila in previdnostni ukrepi
<b>SV</b> <b>48</b>	<b>Small Gauge Ova-Stiff™ Ovum Aspiration Needles, Ova-Stiff™ EchoTip® Single Lumen Ovum Aspiration Needle</b> <b>Ovum Aspiration Needle Single Lumen, Single Lumen Ovum Aspiration Needle</b> Bruksanvisning, varningar och försiktighetsåtgärder
<b>TR</b> <b>50</b>	<b>Small Gauge Ova-Stiff™ Ovum Aspiration Needles, Ova-Stiff™ EchoTip® Single Lumen Ovum Aspiration Needle</b> <b>Ovum Aspiration Needle Single Lumen, Single Lumen Ovum Aspiration Needle</b> Kullanım, uyarılar ve önlemler için kılavuz

## SMALL GAUGE OVA-STIFF™ OVUM ASPIRATION NEEDLES OVA-STIFF™ ECHOTIP® SINGLE LUMEN OVUM ASPIRATION NEEDLE OVUM ASPIRATION NEEDLE SINGLE LUMEN SINGLE LUMEN OVUM ASPIRATION NEEDLE

### 1 DEVICE DESCRIPTION

All models of the Single Lumen Ovum Pick-Up Needles consist of a stainless-steel single lumen needle with a bevelled tip for puncture, an echogenic tip for enhanced visualization under ultrasound and a manipulation handle.

The Small Gauge Ova-Stiff™ Ovum Aspiration Needles (K-OPS model) also includes an aspiration line and a silicone stopper with Luer fitting to allow connection to a vacuum source.

The Ova-Stiff™ EchoTip® Single Lumen Ovum Aspiration Needle (K-OSN model) also includes an aspiration line and a silicone stopper with vacuum line and a Luer fitting to allow connection to a vacuum source.

The Ovum Aspiration Needle Single Lumen (K-UCI model) also includes a flushing line with a Luer fitting to allow connection to a Luer syringe, an aspiration line and a silicone stopper with Luer fitting to allow connection to a vacuum source.

The Single Lumen Ovum Aspiration Needle (K-OPAA model) also includes a Luer fitting on the handle for connection to a Luer lock syringe.

### 2 INTENDED USE

The Single Lumen Ovum Pick-up Needles are intended to be used for ultrasound guided transvaginal aspiration and flushing of ovarian follicles during ovum pick-up procedures.

### 3 CONTRAINDICATIONS

The Single Lumen Ovum Pick-Up Needles are not for use on patients:

- With an active vaginal or intrauterine infection
- With a sexually transmitted disease
- With a recent uterine perforation
- With a recent cesarean section
- After recent pelvic inflammatory disease
- Who are currently pregnant

### 4 WARNINGS

- The Single Lumen Ovum Pick-Up Needles should only be used by persons or under the guidance of qualified persons having adequate training and familiarity with ovum pick-up procedures.
- There is not sufficient clinical evidence of safety and/or efficacy to support the use of Single Lumen Ovum Pick-Up Needles for transabdominal ovum aspiration procedures.
- There is not sufficient clinical evidence of safety and/or efficacy to support the use of Single Lumen Ovum Pick-Up Needles for independent cystic aspiration procedures.
- Cases of significant hemoperitoneum post-operatively have rarely occurred hours to days after oocyte retrieval. All patients should be screened pre-operatively for increased risk of bleeding. At completion of the procedure a final ultrasound review of the ovaries and the recto-uterine pouch should ensure no active bleeding is seen. All patients with increased pain post-operatively should be reviewed.
- If an endometrioma or a hemorrhagic follicle is inadvertently punctured, the needle and collective tubing should be immediately withdrawn and flushed, and the collection tube should be changed. Treat patients in accordance with the higher risk of infection.
- The enlarged ovary following the controlled hyperstimulation as part of infertility treatments is at risk of twisting and restricting the ovarian bloody supply. Tangential needle insertion may increase the risk of that ovarian torsion. Firm pressure on the ovary with the ultrasound transducer may assist to stabilize the ovary. To further stabilize a mobile ovary in one place supra-pubic pressure can be applied (with the help of an assistant).
- This single use device is not designed for re-use. Attempts to reprocess (re-sterilize) and / or to re-use may lead to device failure and / or infection.

### 5 PRECAUTIONS

- Hematuria may occur due to the needle penetrating a filled bladder during transvaginal ultrasound aspiration. This complication typically resolves spontaneously within a day.
- Extravasation of urine may occur within the abdominal cavity if a needle puncture traverses the bladder. Patients should be monitored for evidence of this known complication; however, there is typically no prolonged associated discomfort or adverse sequelae.
- Infection may be introduced via needle puncture and result in urinary tract infection (UTI), pelvic inflammatory disease (PID), peritonitis, uterine infection or cystitis. Recommendations to minimize occurrence include use of only biocompatible materials, flushing the needle prior to use (and any other accessories used) with compatible sterile culture media, and closely adhering to sterile techniques.
- Vaginal/cervical bleeding has been reported to be associated with the transvaginal route for oocyte retrieval via needle aspiration. Bleeding is typically easily controlled with direct pressure or by a suture.
- Pelvic hematomata have been reported in the vaginal vault within an ovary, under the bladder and in the broad ligament. Patients should be monitored for evidence of this known complication; however, there is typically no further intervention required.
- Ultrasound equipment must be of sufficient resolution to visualize the needle tip and pelvic anatomy. If the needle tip and pelvic anatomy cannot be visualized the procedure cannot be successful.

### 6 HOW SUPPLIED

STERILE – DO NOT RESTERILIZE – SINGLE USE ONLY

Single Lumen Ovum Pick-Up Needles are supplied sterile in peel-open packages. They have been sterilised using Ethylene oxide and are intended for single use only.

Do not resterilize the device. Do not use after the expiration date printed on the label. Store in a dark, dry, cool place. Avoid extended exposure to light and extremes of temperature and humidity.

### 7 INSPECTION OF DEVICE

1. Check that the device packaging is not damaged and that the 'use by' date has not expired.  
**CAUTION: Do not use if the sterile packaging is damaged or unintentionally opened before use.**  
**CAUTION: Do not use after the expiration date printed on the label.**
2. Carefully remove the needle from the packaging maintaining the sterility of the product.
3. The sterile needle should be inspected for irregularities in needle tip profile or kinking of any supplied tubing.

### 8 DEVICE PREPARATION FOR K-OSN, K-OPS AND K-UCI MODELS

1. For **K-OSN Models**, connect the Luer fitting on the vacuum line of the needle to a new disposable vacuum line with hydrophobic filter attached to vacuum pump.  
For **K-OPS** and **K-UCI Models**, connect the Luer fitting on the silicone stopper to a new disposable vacuum line with hydrophobic filter attached to vacuum pump.
2. Ensure the vacuum pump has been prepared as per respective manufacturer's instructions.
3. Fit the silicone stopper onto the collection tube (designed to fit 14 mL Falcon tubes).
4. Place the needle tip in a spare collection tube containing approximately 5 mL of sterile culture medium and apply vacuum. Ensure culture medium is aspirated into collection tube.  
**CAUTION: Ensure all connections are tight and leak proof. Failure to confirm connections and system patency may result in failure to complete the procedure.**  
**NOTE:** The vacuum pump pressure should be calibrated to give the desired flow rate for the specific type of needle being used as per the vacuum pump manufacturer's instructions.  
**NOTE:** If blockages occur, remove silicone stopper from collection tube and connect a luer syringe filled with sterile culture medium into the aspiration line port. Inject the culture medium and attempt to force the blockage out. This ejected fluid should be checked for oocytes. Refer to the vacuum pump manufacturer's instructions for needle blockage troubleshooting. If blockage persists, remove blocked device and replace with new device for procedure completion.
5. Discard the used collection tube and replace silicone stopper into collection tube primed with sterile culture medium.
6. Ensure that the ultrasound equipment, ultrasound gel, transducer protective sheath and needle guide have been prepared as per respective manufacturer's instructions.

### 9 DEVICE OPERATION FOR K-OSN, K-OPS AND K-UCI MODELS

1. Position the patient in the lithotomy position on operating table. Local or general anaesthetic may be administered as necessary.
2. Introduce an ultrasound transducer into the vaginal fornix to visualize the ovary and follicles. Identify the follicles to be aspirated. Check for the presence of blood vessels in and around the ovary and determine a direct path into the ovarian follicles to be aspirated.  
**NOTE:** For mobile ovaries, firm pressure with the ultrasound transducer may assist to stabilize the ovary.
3. Check the needle guide is appropriately placed on the ultrasound transducer with no gaps and that the transducer protective sheath does not hinder path of the needle projection.  
**CAUTION: In passage of the needle between personnel, be aware of needle tip location to avoid needle stick injury.**
4. Insert the needle into the needle guide, ensuring there is smooth movement along its length.
5. Visualize the ovary and follicles to be aspirated and advance the needle tip to puncture the vaginal wall.  
**WARNING: The needle tip must be visible at all times as it is maneuvered to the ovaries and into each follicle. If visibility is lost, do not advance the device further. Manipulate the ultrasound transducer to obtain different views. If still unable to visualize the needle tip, withdraw the needle from the patient.**  
**NOTE:** The combination of the needle bevel and EchoTip® enhances visualization of the position of the needle tip. The handle indent indicates bevel orientation as well as facilitating grip.
6. Line up the target follicle using the needle guide on the ultrasound monitor and advance the needle tip into the centre of an ovarian follicle via a rapid, stabbing motion.  
**WARNING: To further stabilize a mobile ovary in one place, supra-pubic pressure can be applied (with the help of an assistant).**  
**WARNING: Be aware that big iliac vessels may resemble a follicle on an ultrasound image when in a cross-sectional plane. Rotate the ultrasound transducer to visualize anatomy to avoid inadvertent puncture.**
7. Apply vacuum to aspirate the follicular contents into the collection tube. As the follicle collapses, rotate the needle tip within the follicle to ensure complete emptying of the follicle.  
**NOTE:** To avoid aspiration difficulty, ensure tubing does not become kinked.  
**NOTE:** The vacuum pump pressure should be calibrated to give the desired flow rate for the specific type of needle being used as per the vacuum pump manufacturer's instructions.  
**NOTE:** If blockages occur, stop aspiration, remove the needle from the patient, remove silicone stopper from collection tube and connect a luer syringe filled with sterile culture medium into the aspiration line port. Inject the culture medium and attempt to force the blockage out. This ejected fluid should be checked for oocytes. Refer to the vacuum pump manufacturer's instructions for needle blockage troubleshooting. If blockage persists, remove blocked device and replace with new device for procedure completion.
8. Move needle tip to the adjacent follicle and repeat **Step 7**, maintaining vacuum aspiration. Where possible, keep the needle tip within the stroma or follicles of that ovary until follicular aspiration is complete.  
**CAUTION: Avoid recurrent penetration of the vaginal vault and ovarian capsule to minimize risks of infection and bleeding.**
9. Replace the silicone stopper into primed collection tubes at appropriate aspiration intervals. To avoid aspiration difficulty, ensure aspirated contents do not reach the lip of silicone stopper in collection tube.
10. If required the follicle can be flushed, as described below.  
**Follicle Flushing for K-OSN and K-OPS Models:**
  - a. Stop aspiration and remove the silicone stopper from the collection tube.
  - b. Use a Luer syringe filled with a follicle-flushing buffer. Insert the syringe into the aspiration line port.
  - c. With the needle tip in the collapsed follicle, slowly inject (1–2 mL per second) the flushing buffer to refill the follicle.
  - d. Replace the stopper and aspirate the follicular contents.**Follicle Flushing for K-UCI Models:**
  - a. Stop aspiration.
  - b. Use a Luer syringe filled with a sterile follicle-flushing buffer. Insert the syringe into the Luer fitting on the flushing line.
  - c. With the needle tip in the collapsed follicle, slowly inject (1–2 mL per second) the flushing buffer to refill the follicle.
  - d. Aspirate the follicular contents.
11. Confirm all follicles in the ovary have been aspirated. Stop aspiration before withdrawing the needle through the vaginal vault.
12. Place the needle tip in a spare collection tube containing approximately 5 mL of sterile culture medium and apply vacuum. Ensure culture medium is aspirated into the collection tube.  
**NOTE:** This aspirated content should be checked for oocytes.

13. Reposition the ultrasound transducer to visualize the remaining ovary. Repeat **Steps 5 to 10**.
14. At the completion of the aspiration procedure, stop aspiration and remove the needle from the patient.
15. Place the needle tip in a spare collection tube containing approximately 5 mL of sterile culture medium and apply vacuum. Ensure culture medium is aspirated into collection tube.  
**NOTE:** This aspirated content should be checked for oocytes.
16. Before the ultrasound transducer is removed, check for any evidence of continuing bleeding in the ovaries or peritoneal cavity.
17. Check the vagina for signs of bleeding. If vaginal bleeding is experienced post procedure, apply steady pressure. Consider inserting a suture to control continuing vaginal vault or cervical bleeding.

#### **10 DEVICE PREPARATION FOR K-OPAA MODELS:**

1. Fit an empty Luer lock syringe to the Luer fitting on the handle.
2. Place the needle tip in a spare collection tube containing approximately 5 mL of sterile culture medium.
3. Withdraw the culture medium into the empty syringe.
4. Discard the used syringe and used culture medium.
5. Replace with a Luer lock syringe primed with culture medium for aspiration.  
**CAUTION: Ensure all connections are tight and leak proof. Failure to confirm connections and system patency may result in failure to complete the procedure.**
6. Ensure that the ultrasound equipment, ultrasound gel, transducer protective sheath and needle guide have been prepared as per respective manufacturer's instructions.

#### **11 DEVICE OPERATION FOR K-OPAA MODELS:**

1. Position the patient in the lithotomy position on operating table. Local or general anaesthetic may be administered as necessary.
2. Introduce an ultrasound transducer into the vaginal fornix to visualize the ovary and follicles. Identify the follicles to be aspirated. Check for the presence of blood vessels in and around the ovary and determine a direct path into the ovarian follicles to be aspirated.  
**NOTE:** For mobile ovaries, firm pressure with the ultrasound transducer may assist to stabilise the ovary.
3. Check the needle guide is appropriately placed on the ultrasound transducer with no gaps and that the transducer protective sheath does not hinder path of the needle projection.  
**CAUTION: In passage of the needle between personnel, be aware of needle tip location to avoid needle stick injury.**
4. Insert the needle into the needle guide, ensuring there is smooth movement along its length.
5. Visualise the ovary and follicles to be aspirated and advance the needle tip to puncture the vaginal wall.  
**WARNING: The needle tip must be visible at all times as it is manoeuvred to the ovaries and into each follicle. If visibility is lost, do not advance the device further. Manipulate the ultrasound transducer to obtain different views. If still unable to visualize the needle tip, withdraw the needle from the patient.**  
**NOTE:** The combination of the needle bevel and EchoTip® enhances visualization of the position of the needle tip. The handle indent indicates bevel orientation as well as facilitating grip.
6. Line up the target follicle using the needle guide on the ultrasound monitor and advance the needle tip into the centre of an ovarian follicle via a rapid, stabbing motion.  
**WARNING: To further stabilize a mobile ovary in one place, supra-pubic pressure can be applied (with the help of an assistant).**  
**WARNING: Be aware that big iliac vessels may resemble a follicle on an ultrasound image when in a cross-sectional plane. Rotate the ultrasound transducer to visualize anatomy to avoid inadvertent puncture.**
7. Use the connected syringe to withdraw the follicular contents. As the follicle collapses, rotate the needle tip within the follicle to ensure complete emptying of the follicle.
8. Move needle tip to the adjacent follicle and repeat **Step 7**.  
**CAUTION: Avoid recurrent penetration of the vaginal vault and ovarian capsule to minimise risks of infection and bleeding.**
9. If required the follicle can be flushed, as described below.  
**Follicle Flushing for K-OPAA Models:**
  - a. Remove the Luer lock syringe being used for aspiration from handle.
  - b. Connect a Luer lock syringe filled with sterile follicle-flushing buffer to the handle.
  - c. With the needle tip in the collapsed follicle, slowly inject (1–2 mL per second) the flushing buffer to refill the follicle.
  - d. Remove the flushing syringe from the handle and replace with the aspiration syringe and continue to aspirate the follicular contents.
10. Confirm all follicles in the ovary have been aspirated. Remove the needle from the patient.
11. Place the needle tip in a spare collection tube containing approximately 5 mL of sterile culture medium and aspirate with syringe. Ensure culture medium is aspirated into the syringe.  
**NOTE:** This aspirated content should be checked for oocytes.
12. Reposition the ultrasound transducer to visualize the remaining ovary. Repeat **Steps 5 to 9**.
13. At the completion of the aspiration procedure, stop aspiration and remove the needle from the patient.
14. Place the needle tip in a spare collection tube containing approximately 5 mL of sterile culture medium and aspirate with syringe. Ensure culture medium is aspirated into the syringe.  
**NOTE:** This aspirated content should be checked for oocytes.
15. Before the ultrasound transducer is removed, check for any evidence of continuing bleeding in the ovaries or peritoneal cavity.
16. Check the vagina for signs of bleeding. If vaginal bleeding is experienced post procedure, apply steady pressure. Consider inserting a suture to control continuing vaginal vault or cervical bleeding.

#### **12 DISPOSAL OF DEVICES**

1. Dispose of aspiration needle in appropriate sharps container and any consumables in accordance with institutional guidelines.  
**WARNING:** This single use device is not designed for re-use. Attempts to reprocess (re-sterilize) and/or to re-use may lead to device failure and /or infection.

## SMALL GAUGE OVA-STIFF™ OVUM ASPIRATION NEEDLES OVA-STIFF™ ECHOTIP® SINGLE LUMEN OVUM ASPIRATION NEEDLE OVUM ASPIRATION NEEDLE SINGLE LUMEN SINGLE LUMEN OVUM ASPIRATION NEEDLE

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## COOK® ECHOTIP® DOUBLE LUMEN ASPIRATION NEEDLE

### 1 DEVICE DESCRIPTION

The Cook® EchoTip® Double Lumen Aspiration Needle consists of: a stainless-steel double lumen needle with a bevelled tip for puncture and an echogenic tip for enhanced visualisation under ultrasound, a manipulation handle, aspiration line, a silicone stopper with Luer fitting to allow connection to a vacuum source and flushing line.

### 2 INTENDED USE

The Double Lumen Ovum Pick-up Needles are intended to be used for ultrasound guided transvaginal simultaneous or independent aspiration and flushing of ovarian follicles during ovum pick-up procedures.

### 3 CONTRAINDICATIONS

The Double Lumen Ovum Pick-Up Needles are not for use on patients:

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### 4 WARNINGS

- The Double Lumen Ovum Pick-Up Needles should only be used by persons or under the guidance of qualified persons having adequate training and familiarity with ovum pick-up procedures.
- There is not sufficient clinical evidence of safety and/or efficacy to support the use of Double Lumen Ovum Pick-Up Needles for transabdominal ovum aspiration procedures.
- There is not sufficient clinical evidence of safety and/or efficacy to support the use of Double Lumen Ovum Pick-Up Needles for independent cystic aspiration procedures.
- Cases of significant hemoperitoneum post-operatively have rarely occurred hours to days after oocyte retrieval. All patients should be screened pre-operatively for increased risk of bleeding. At completion of the procedure a final ultrasound review of the ovaries and the recto-uterine pouch should ensure no active bleeding is seen. All patients with increased pain post-operatively should be reviewed.
- If an endometrioma or a hemorrhagic follicle is inadvertently punctured, the needle and collective tubing should be immediately withdrawn and flushed, and the collection tube should be changed. Treat patients in accordance with the higher risk of infection.
- The enlarged ovary following the controlled hyperstimulation as part of infertility treatments is at risk of twisting and restricting the ovarian blood supply. Tangential needle insertion may increase the risk of that ovarian torsion. Firm pressure on the ovary with the ultrasound transducer may assist to stabilise the ovary. To further stabilise a mobile ovary in one place supra-pubic pressure can be applied (with the help of an assistant).
- This single use device is not designed for re-use. Attempts to reprocess (re-sterilize) and/or to re-use may lead to device failure and/or infection.

### 5 PRECAUTIONS

- Hematuria may occur due to the needle penetrating a filled bladder during transvaginal ultrasound aspiration. This complication typically resolves spontaneously within a day.
- Extravasation of urine may occur within the abdominal cavity if a needle puncture traverses the bladder. Patients should be monitored for evidence of this known complication; however, there is typically no prolonged associated discomfort or adverse sequelae.
- Infection may be introduced via needle puncture and result in urinary tract infection (UTI), pelvic inflammatory disease (PID), peritonitis, uterine infection or cystitis. Recommendations to minimise occurrence include use of only biocompatible materials, flushing the needle prior to use (and any other accessories used) with compatible sterile culture media, and closely adhering to sterile techniques.
- Vaginal/cervical bleeding has been reported to be associated with the transvaginal route for oocyte retrieval via needle aspiration. Bleeding is typically easily controlled with direct pressure or by a suture.
- Pelvic hematomata have been reported in the vaginal vault within an ovary, under the bladder and in the broad ligament. Patients should be monitored for evidence of this known complication; however, there is typically no further intervention required.
- Due to the larger gauge of double lumen cannula, additional intra-operative pain management may be required.
- Ultrasound equipment must be of sufficient resolution to visualise the needle tip and pelvic anatomy. If the needle tip and pelvic anatomy cannot be visualised the procedure cannot be successful.

### 6 HOW SUPPLIED

STERILE – DO NOT RESTERILIZE – SINGLE USE ONLY

Double Lumen Ovum Pick-Up Needles are supplied sterile in peel-open packages. They have been sterilised using Ethylene oxide and are intended for single use only.

Do not resterilise the device. Do not use after the expiration date printed on the label. Store in a dark, dry, cool place. Avoid extended exposure to light and extremes of temperature and humidity.

### 7 INSPECTION OF DEVICE

1. Check that the device packaging is not damaged and that the 'use by' date has not expired.
  - CAUTION: Do not use if the sterile packaging is damaged or unintentionally opened before use.**
  - CAUTION: Do not use after the expiration date printed on the label.**
2. Carefully remove the needle from the packaging maintaining the sterility of the product.
3. The sterile needle should be inspected for irregularities in needle tip profile or kinking of any supplied tubing.

### 8 DEVICE PREPARATION

1. Connect the Luer fitting on the silicone stopper to a new disposable vacuum line with hydrophobic filter attached to vacuum pump.
2. Ensure the vacuum pump has been prepared as per respective manufacturer's instructions.
3. Fit the silicone stopper onto the collection tube (designed to fit 14 mL Falcon tubes).
4. Place the needle tip in a spare collection tube containing approximately 5 mL of sterile culture medium and apply vacuum. Ensure culture medium is aspirated into collection tube.

**CAUTION: Ensure all connections are tight and leak proof. Failure to confirm connections and system patency may result in failure to complete the procedure.**

**NOTE:** The vacuum pump pressure should be calibrated to give the desired flow rate for the specific type of needle being used as per the vacuum pump manufacturer's instructions.

**NOTE:** If blockages occur, remove silicone stopper from collection tube and connect a luer syringe filled with sterile culture medium into the aspiration line port. Inject the culture medium and attempt to force the blockage out. Refer to the vacuum pump manufacturer's instructions for additional blockage troubleshooting. If blockage persists, remove blocked device and replace with new device for procedure completion.

5. Discard the used collection tube and replace silicone stopper into collection tube primed with sterile culture medium.
6. Ensure that the ultrasound equipment, ultrasound gel, transducer protective sheath and needle guide have been prepared as per respective manufacturer's instructions.

### 9 DEVICE OPERATION

1. Position the patient in the lithotomy position on operating table. Local or general anaesthetic may be administered as necessary.
2. Introduce an ultrasound transducer into the vaginal fornix to visualize the ovary and follicles. Identify the follicles to be aspirated. Check for the presence of blood vessels in and around the ovary and determine a direct path into the ovarian follicles to be aspirated.
  - NOTE:** For mobile ovaries, firm pressure with the ultrasound transducer may assist to stabilize the ovary.
3. Check the needle guide is appropriately placed on the ultrasound transducer with no gaps and that the transducer protective sheath does not hinder path of the needle projection.
  - CAUTION: In passage of the needle between personnel, be aware of needle tip location to avoid needle stick injury.**
4. Insert the needle into the needle guide, ensuring there is smooth movement along its length.
5. Visualize the ovary and follicles to be aspirated and advance the needle tip to puncture the vaginal wall.
  - WARNING: The needle tip must be visible at all times as it is maneuvered to the ovaries and into each follicle. If visibility is lost, do not advance the device further. Manipulate the ultrasound transducer to obtain different views. If still unable to visualize the needle tip, withdraw the needle from the patient.**
  - NOTE:** The combination of the needle bevel and EchoTip® enhances visualization of the position of the needle tip. The handle indent indicates bevel orientation as well as facilitating grip.
6. Line up the target follicle using the needle guide on the ultrasound monitor and advance the needle tip into the centre of an ovarian follicle via a rapid, stabbing motion.
  - WARNING: To further stabilize a mobile ovary in one place, supra-pubic pressure can be applied (with the help of an assistant).**
  - WARNING: Be aware that big iliac vessels may resemble a follicle on an ultrasound image when in a cross-sectional plane. Rotate the ultrasound transducer to visualize anatomy to avoid inadvertent puncture.**
7. Apply vacuum to aspirate the follicular contents into the collection tube. As the follicle collapses, rotate the needle tip within the follicle to ensure complete emptying of the follicle.
  - NOTE:** To avoid aspiration difficulty, ensure tubing does not become kinked.
  - NOTE:** The vacuum pump pressure should be calibrated to give the desired flow rate for the specific type of needle being used as per the vacuum pump manufacturer's instructions.
  - NOTE:** If blockages occur, stop aspiration, remove the needle from the patient, remove silicone stopper from collection tube and connect a luer syringe filled with sterile culture medium into the aspiration line port. Inject the culture medium and attempt to force the blockage out. This ejected fluid should be checked for oocytes. Refer to the vacuum pump manufacturer's instructions for needle blockage troubleshooting. If blockage persists, remove blocked device and replace with new device for procedure completion.
8. Move needle tip to the adjacent follicle and repeat **Step 7**, maintaining vacuum aspiration. Where possible, keep the needle tip within the stroma or follicles of that ovary until follicular aspiration is complete.
  - CAUTION: Avoid recurrent penetration of the vaginal vault and ovarian capsule to minimise risks of infection and bleeding.**
9. Replace the silicone stopper into primed collection tubes at appropriate aspiration intervals. To avoid aspiration difficulty, ensure aspirated contents do not reach the lip of silicone stopper in collection tube.
10. If required the follicle can be flushed, as described below.
  - a. Stop aspiration.
  - b. Use a Luer syringe filled with a sterile follicle-flushing buffer. Insert the syringe into the Luer fitting on the flushing line.
  - c. With the needle tip in the collapsed follicle, slowly inject (1-2 mL per second) the flushing buffer to refill the follicle.
  - d. Aspirate the follicular contents.
11. Confirm all follicles in the ovary have been aspirated. Stop aspiration before withdrawing the needle through the vaginal vault.
12. Place the needle tip in a spare collection tube containing approximately 5 mL of sterile culture medium and apply vacuum. Ensure culture medium is aspirated into the collection tube.
  - NOTE:** This aspirated content should be checked for oocytes.
13. Reposition the ultrasound transducer to visualize the remaining ovary. Repeat **Steps 5 to 10**.
14. At the completion of the aspiration procedure, stop aspiration and remove the needle from the patient.
15. Place the needle tip in a spare collection tube containing approximately 5 mL of sterile culture medium and apply vacuum. Ensure culture medium is aspirated into collection tube.
  - NOTE:** This aspirated content should be checked for oocytes.
16. Before the ultrasound transducer is removed, check for any evidence of continuing bleeding in the ovaries or peritoneal cavity.
17. Check the vagina for signs of bleeding. If vaginal bleeding is experienced post procedure, apply steady pressure. Consider inserting a suture to control continuing vaginal vault or cervical bleeding.

### 10 DISPOSAL OF DEVICES

1. Dispose of aspiration needle in appropriate sharps container and any consumables in accordance with institutional guidelines.
  - WARNING:** This single use device is not designed for re-use. Attempts to reprocess (re-sterilize) and/or to re-use may lead to device failure and/or infection.

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