



Wednesday, October 27, 2021

Quality Notification
URGENT FIELD SAFETY NOTICE

Illumina is contacting you regarding an issue affecting certain lots of the VeriSeq™ NIPT Sample Prep Kit. Our records indicate that you have received material from these lots. This notification outlines the issue summary, Illumina actions, and customer actions required.

Table 1: Affected Products

| Product Name | Part Number | Component Name | Catalog Number | Catalog Lot Numbers |
|---|-------------|-------------------------------------|----------------|-------------------------------------|
| VeriSeq™ NIPT Sample Prep Kit (24 Sample) | 20025895 | VeriSeq™ NIPT Library Prep Box (24) | 20026030 | A162813-2 A164108-2 A163501-2 |
| VeriSeq™ NIPT Sample Prep Kit (48 Sample) | 15066801 | VeriSeq™ NIPT Library Prep Box (48) | 15066809 | A162738-2 A163577-2 |

Please see specific impacted index adapter plate barcode numbers in Appendix

Issue Summary

Illumina has received reports of samples in column 11 having low Non-Excluded Sites (NES) resulting in iFACT failures leading to no results. This occurs due to a shift in the library pool indices in column 11 that the sequencing run does not expect causing a sample failure. Illumina has identified 154 out of 704 plates are impacted by this issue. The root cause investigation is still underway.

The Appendix below has the index adapter plate barcodes of the 154 plates affected by this issue.

For any prior runs performed with impacted plates, if all samples have passed, there is no impact and no follow up is required. If you experienced sample failure in a previously completed run with impacted plates, please contact tech support.

The impacted index adapter plates have been kitted in both 24-sample VeriSeq™ NIPT Sample Prep Kit and 48-sample VeriSeq™ NIPT Sample Prep Kits (see lot numbers in the table above).

Technical Support:

techsupport@illumina.com

Customer Care:

customercare@illumina.com

© 2021 Illumina, Inc. All rights reserved. All trademarks are the property of Illumina, Inc. or their respective owners. For specific trademark information, see www.illumina.com/company/legal.html.

FSN2021-1197 (M-AMR-00354)



Recall

24-Sample Batches

During the 24-sample run, samples get associated with incorrect indices in column 11 due to a sample shift in the index adapter plate. These indices are not expected by the analysis software, causing samples 21 and 22 in wells E3 and F3 of the sample plate to fail. This issue results in iFACT sample failure and leads to a delay in patient results. No incorrect results are generated.

48-Sample Batches

During the 48-sample run, iFACT sample failure will not occur if the kit is used per VeriSeq™NIPT Solution v2 Package Insert instructions because indices in column 11 are not used. The 48-sample kits are included in this notification as an awareness because certain affected adapter plates (listed in Appendix below) were kitted in 48-sample kits (see Table 1 for lot numbers).

Illumina Actions

Illumina has stopped shipment of all impacted kits and is working to implement new quality checks to prevent the occurrence of this issue in the future.

The appropriate Competent Authorities are being notified of this issue.

Required Customer Actions

For 24-sample and 48-sample kits, please use the Appendix below to inspect your current inventory for potentially affected index adapter plates. If you have a plate listed in the Appendix, discontinue use of the impacted plate(s) and destroy per your local procedure. Please contact Technical Support (techsupport@illumina.com) for replacement.

If you have questions about past completed runs, please contact Technical Support or your local FAS.

We ask that you complete the attached FSN2021-1197 Verification Form to confirm receipt of this notification and return to Technical Support at techsupport@illumina.com within 5 business days.

We understand the impact that issues affecting test result delivery can cause. We greatly appreciate your continued business and regret any inconvenience this may have caused.

Technical Support:

techsupport@illumina.com

Customer Care:

customercare@illumina.com

© 2021 Illumina, Inc. All rights reserved. All trademarks are the property of Illumina, Inc. or their respective owners. For specific trademark information, see www.illumina.com/company/legal.html.

FSN2021-1197 (M-AMR-00354)



Recall

illumina®

Sincerely,

Gary Workman

VP, Quality

Electronically
signed by: Gary
Workman
Reason:
Approver
Date: Oct 26,
2021 11:06 EDT

Karen Gutekunst

VP, Regulatory Affairs

Electronically
signed by:
Karen
Gutekunst
Reason:
Approver
Date: Oct 26,
2021 09:22 PDT

Why You're Receiving This Notification

We are sending this notification to you because our records indicate that you are one of the appropriate contacts for your organization. We occasionally need to inform our customers of product changes, product obsolescence, or quality issues.

Accordingly, please note that these notifications contain important information about our products and are not marketing communications. You may, therefore, receive these notifications even if you have opted out of receiving marketing material from Illumina. If you are not the appropriate individual in your organization to receive these types of notifications, please email customernotifications@illumina.com with the appropriate contact. For more information, please see our [Privacy Policy](#).

Technical Support:

techsupport@illumina.com

Customer Care:

customercare@illumina.com

© 2021 Illumina, Inc. All rights reserved. All trademarks are the property of Illumina, Inc. or their respective owners. For specific trademark information, see www.illumina.com/company/legal.html.

FSN2021-1197 (M-AMR-00354)

Page 3 of 6

For In Vitro Diagnostic Use. Not available in all countries or regions.

**Appendix:****Affected Plate Barcodes**

| | | | | | |
|----------|----------|----------|----------|----------|----------|
| 14514076 | 14514235 | 14514249 | 14514100 | 14514386 | 14514344 |
| 14514075 | 14514234 | 14514248 | 14514085 | 14514376 | 14514343 |
| 14514074 | 14514233 | 14514247 | 14514084 | 14514375 | 14514342 |
| 14514073 | 14514232 | 14514246 | 14514083 | 14514374 | 14514341 |
| 14514072 | 14514231 | 14514245 | 14514082 | 14514373 | 14514340 |
| 14514071 | 14514230 | 14514244 | 14514081 | 14514372 | 14514339 |
| 14514070 | 14514229 | 14514243 | 14514077 | 14514361 | 14514338 |
| 14514069 | 14514228 | 14514242 | 14514723 | 14514266 | 14514336 |
| 14514067 | 14514227 | 14514241 | 14514722 | 14514265 | 14514335 |
| 14514066 | 14514226 | 14514240 | 14514721 | 14514264 | 14514399 |
| 14514065 | 14514225 | 14514239 | 14514720 | 14514263 | 14514390 |
| 14514064 | 14514224 | 14514238 | 14514719 | 14514262 | 14514389 |
| 14514063 | 14514223 | 14514237 | 14514718 | 14514261 | 14514388 |
| 14514062 | 14514211 | 14514236 | 14514716 | 14514260 | 14514334 |
| 14514061 | 14514210 | 14514360 | 14514715 | 14514259 | 14514044 |
| 14514060 | 14514209 | 14514359 | 14514714 | 14514257 | 14514043 |
| 14514059 | 14514208 | 14514358 | 14514713 | 14514256 | 14514042 |
| 14514058 | 14514207 | 14514357 | 14514712 | 14514255 | 14514040 |
| 14514056 | 14514206 | 14514356 | 14514708 | 14514254 | 14514039 |
| 14514055 | 14514104 | 14514355 | 14514707 | 14514253 | 14514038 |
| 14514054 | 14514103 | 14514354 | 14514702 | 14514252 | 14514037 |
| 14514053 | 14514102 | 14514353 | 14514400 | 14514251 | 14514036 |
| 14514052 | 14514101 | 14514352 | 14514398 | 14514347 | 14514035 |
| 14514051 | 14514048 | 14514351 | 14514397 | 14514346 | 14514349 |
| 14514050 | 14514047 | 14514350 | 14514395 | 14514345 | 14514348 |
| 14514049 | 14514046 | 14514045 | 14514391 | | |

Technical Support:techsupport@illumina.com**Customer Care:**customercare@illumina.com

© 2021 Illumina, Inc. All rights reserved. All trademarks are the property of Illumina, Inc. or their respective owners. For specific trademark information, see www.illumina.com/company/legal.html.



FSN2021-1197 Verification Form

Dear Customer,

On <DATE>, Illumina sent you an Urgent Field Safety Notice FSN2021-1197 regarding VeriSeq™NIPT Solution v2 Sample Prep Kit Column 11 Index Adapter Plate issue.

Please complete the table below to confirm that you received the notice, as well as the disposition of the recalled product. We ask that you kindly return the completed form to Illumina within 5 business days using one of the following methods:

- Scan the completed, signed form and email it to TechSupport@illumina.com
- Fax to 858-736-8426, **Attn: Illumina QA, FSN2021-1197 Notice**
- Return to Illumina at the address and “attn to:” below:
Attn: Illumina QA, FSN2021-1197 Notice
5200 Illumina Way, San Diego, CA 92122

| Verification Form | | | | |
|---|----------------|-------------|--|------------------|
| Company Name | | | | |
| Product/Device Name | Catalog Number | Part Number | Lot Number | Quantity Shipped |
| | | | | |
| Your Information | | | | |
| Print full name: | | | | |
| Print title of person completing form: | | | | |
| Customer Responses | | | | |
| I confirm receipt of FSN2021-1197 customer communication and that I read and understood its content. | | | Yes <input type="checkbox"/> No <input type="checkbox"/> | |
| I have performed all relevant actions requested by this communication. | | | Yes <input type="checkbox"/> No <input type="checkbox"/> | |
| The information and required actions have been brought to the attention of all relevant users and actions have been executed. | | | Yes <input type="checkbox"/> No <input type="checkbox"/> | |
| Product has been consumed. No affected kits are available for return, destruction, or inspection. | | | <input type="checkbox"/> | |

Technical Support:

techsupport@illumina.com

Customer Care:

customercare@illumina.com

© 2021 Illumina, Inc. All rights reserved. All trademarks are the property of Illumina, Inc. or their respective owners. For specific trademark information, see www.illumina.com/company/legal.html.

FSN2021-1197 (M-AMR-00354)

Page 5 of 6



Recall

| | |
|---|--|
| I have completed the inspection and kits are appropriate for use. | <input type="checkbox"/> |
| I have destroyed affected devices | <input type="checkbox"/> |
| Other actions (if any): | |
| Quantity Consumed: | Date Complete (DD/MM/YY) |
| Lot/Serial #(s) and Quantity: | |
| Comments: | |
| Distributor/Importer Responses | Not applicable <input type="checkbox"/> |
| I have identified customers that received or may have received the product. | Yes <input type="checkbox"/> No <input type="checkbox"/> |
| I have attached customer list, or I have listed the customers below. | Yes <input type="checkbox"/> No <input type="checkbox"/> |
| I have received confirmation from all identified customers | Yes <input type="checkbox"/> No <input type="checkbox"/> |
| I have informed the identified customers of this recall | Date (MM/DD/YY) |

Signature of Person Completing Form

Date

Technical Support:

techsupport@illumina.com

Customer Care:

customercare@illumina.com

© 2021 Illumina, Inc. All rights reserved. All trademarks are the property of Illumina, Inc. or their respective owners. For specific trademark information, see www.illumina.com/company/legal.html.

FSN2021-1197 (M-AMR-00354)