

## **Urgent Field Safety Notice** – risk of device unexpected shutdown

### **AquaVENT® FD140i Gas Flow Driver for CPAP and HFOT**

Product codes AMFD140i-UK, AMFD140i-EU, AMFD140i-UK2, AMFD140i-EU2

Please pass this Field Safety Notice (FSN) to all persons in your organisation who need to be aware of it.

|                          |   |
|--------------------------|---|
| <b>Type of Action:</b>   | To communicate an identified issue which may result in unexpected shutdown of device during clinical use  |
| <b>Device:</b>           | AquaVENT® FD140i Gas Flow Driver  |
| <b>Manufacturer:</b>     | Armstrong Medical Limited (Coleraine, Northern Ireland)   |
| <b>Date of Issue:</b>    | 11 January 2022   |
| <b>For Attention of:</b> | Nursing and medical staff (caregivers) working in critical care, emergency care and high-dependency areas of hospitals and all others to whom potentially affected devices have been transferred, including distributors. |
| <b>Scope of Action:</b>  | Manufacturing serial number specific corrective action  |
| <b>Keywords:</b>         | Breathing System, Flow Driver, HFOT, CPAP   |

#### **Summary**

AquaVENT® FD140i Gas Flow Driver is a respiratory therapy device. It assists respiration using continuous positive airway pressure (CPAP) and high flow oxygen therapy (HFOT) in patients in a hospital setting.

We have been made aware of a number of incidents during clinical use of the AquaVENT® FD140i Gas Flow Driver. Investigation of the reported incidents suggests that the devices involved were in use on battery power, rather than on mains power.

**During the reported incidents, the devices underwent unexpected power shutdown, without appropriate alarms. This required intervention by caregivers to avert patient harm.** The device defect has been traced to a characteristic of software version 1.01, which runs on all devices in the field at this time.



The unwanted behaviour relates to a loss of communication between the FLASH memory stored on the main printed circuit board and the back-up battery when the device is powered OFF when receiving a mains power charging current.

In software version 1.01, the device requires connection to mains power when the device is powered ON as a means to update the FLASH memory with the accurate remaining % capacity within the battery.


The loss of communication potentially allows storage, within the FLASH memory, of an erroneous battery % charge level that could be higher or lower than the actual remaining % capacity within the battery. If a device, in that condition, were to be powered ON without mains power supply, the device – which is then running on battery power – may shutdown (power OFF) unexpectedly during clinical therapy, without warning and without the expected audible and visual alarms associated with battery capacity depletion.

#### **Affected devices:**

See appendices A and B. AquaVENT® FD140i Gas Flow Drivers sold under product codes AMFD140i-UK, AMFD140i-UK2, AMFD140i-EU and AMFD140i-EU2 and in the serial number ranges 10200006 to 10200284, 1021109 to 1021140 and 20200187 to 20210046.

#### **Action in the field:**

All devices having software version 1.01 will be updated with new software release version 1.02. Our representatives will make contact with you soon to arrange for the new release software to be installed. The current software running on your device is displayed during device switch ON at conclusion of the Self Check. Please see graphic, *right*.



### Advice for continued safe use (until replacement software has been installed):

The device should be used on mains power. The device can be used on battery power ONLY where mains power supply is not available or convenient to use, such as during a short time-period during patient transportation. If use on battery power is required, the device must first have been powered ON whilst receiving mains charging. Prior to clinical use on battery power, the device should be fully charged (that is, displaying 100% remaining battery charge accompanied by a solid green battery indicator LED) before disconnecting from mains power supply and running the device on battery power. **Failure to follow this advice for continued safe use could result in patient injury secondary to oxygen desaturation and / or alveolar collapse.**



### Other actions:

The device user manual will be revised and re-issued as version 06. On page 08, the following text will be added:

*Clinical use of FD140i should be on mains power. Use on battery power should be restricted to short periods where mains power is not available or convenient to use, such as during patient transportation. In the event that FD140i is to be used on battery power, ensure that it is charged prior to use - preferably showing 100% remaining battery capacity on the display. Full battery charge level is associated with a solid green indicator LED. This LED is located on the bottom right of the front panel, above the 'ON/OFF' button.*

## Appendix A

### Serial number format

#### Example 1: 10190001

10190001 First "10" represents device with paramagnetic oxygen sensor option

10190001 "19" represents 2019 (year of manufacture)


10190001 Last "1" represents the 1st device manufactured

#### Example 2: 20200099


20200099 First "20" represents device with oxygen fuel cell option

20200099 "20" represents 2020 (year of manufacture)

20200099 Represents the 99th device manufactured



Example 1



Example 2

## Appendix B

### List of affected serial numbers by product code

| Product Code: AMFD140i-UK |          |          |          |
|---------------------------|----------|----------|----------|
| 10200006                  | 10200047 | 10200091 | 10200148 |
| 10200007                  | 10200048 | 10200092 | 10200149 |
| 10200008                  | 10200049 | 10200093 | 10200150 |
| 10200009                  | 10200050 | 10200094 | 10200152 |
| 10200010                  | 10200053 | 10200095 | 10200153 |
| 10200011                  | 10200054 | 10200096 | 10200197 |
| 10200012                  | 10200055 | 10200097 | 10200198 |
| 10200013                  | 10200056 | 10200098 | 10200199 |
| 10200014                  | 10200057 | 10200099 | 10200200 |
| 10200015                  | 10200058 | 10200100 | 10200201 |
| 10200016                  | 10200059 | 10200101 | 10200202 |
| 10200017                  | 10200060 | 10200102 | 10200204 |
| 10200018                  | 10200061 | 10200103 | 10200205 |
| 10200019                  | 10200062 | 10200104 | 10200206 |
| 10200020                  | 10200063 | 10200105 | 10200207 |
| 10200021                  | 10200064 | 10200106 | 10200208 |
| 10200022                  | 10200065 | 10200107 | 10200209 |
| 10200023                  | 10200066 | 10200110 | 10200227 |
| 10200024                  | 10200067 | 10200111 | 10200228 |
| 10200025                  | 10200068 | 10200113 | 10200234 |
| 10200026                  | 10200069 | 10200114 | 10200235 |
| 10200027                  | 10200070 | 10200115 | 10200236 |
| 10200028                  | 10200071 | 10200116 | 10200239 |
| 10200029                  | 10200072 | 10200118 | 10200241 |
| 10200030                  | 10200073 | 10200119 | 10200242 |
| 10200031                  | 10200074 | 10200120 | 10200243 |
| 10200032                  | 10200075 | 10200122 | 10200244 |
| 10200033                  | 10200076 | 10200123 | 10200245 |
| 10200034                  | 10200077 | 10200124 | 10200246 |
| 10200035                  | 10200078 | 10200125 | 10200248 |
| 10200036                  | 10200079 | 10200127 | 10200250 |
| 10200037                  | 10200080 | 10200128 | 10200253 |
| 10200038                  | 10200081 | 10200129 | 10200254 |
| 10200039                  | 10200082 | 10200132 | 10200261 |
| 10200040                  | 10200083 | 10200134 | 10200280 |
| 10200041                  | 10200084 | 10200136 | 10200281 |
| 10200042                  | 10200085 | 10200142 | 10200283 |
| 10200043                  | 10200086 | 10200143 | 10200284 |
| 10200044                  | 10200087 | 10200145 | 10210139 |
| 10200045                  | 10200089 | 10200146 | 10210140 |
| 10200046                  | 10200090 | 10200147 |          |

| Product Code: AMFD140i-UK2 |          |          |          |
|----------------------------|----------|----------|----------|
| 20200187                   | 20200192 | 20200216 | 20200226 |
| 20200188                   | 20200211 | 20200217 | 20210044 |
| 20200189                   | 20200212 | 20200218 | 20210045 |
| 20200190                   | 20200214 | 20200222 | 20210046 |
| 20200191                   | 20200215 | 20200224 |          |

| Product Code: AMFD140i-EU |          |          |          |
|---------------------------|----------|----------|----------|
| 10200157                  | 10200171 | 10200258 | 10200271 |
| 10200158                  | 10200172 | 10200259 | 10200272 |
| 10200159                  | 10200173 | 10200260 | 10200273 |
| 10200161                  | 10200175 | 10200262 | 10200274 |
| 10200162                  | 10200176 | 10200263 | 10200275 |
| 10200163                  | 10200177 | 10200264 | 10200276 |
| 10200165                  | 10200178 | 10200265 | 10200277 |
| 10200166                  | 10200179 | 10200266 | 10200278 |
| 10200167                  | 10200180 | 10200267 | 10200279 |
| 10200168                  | 10200255 | 10200268 | 10210109 |
| 10200169                  | 10200256 | 10200269 |          |
| 10200170                  | 10200257 | 10200270 |          |

| Product Code: AMFD140i-EU2 |          |          |          |
|----------------------------|----------|----------|----------|
| 10200186                   | 20200225 | 20200299 | 20200306 |
| 20200213                   | 20200231 | 20200300 | 20200307 |
| 20200219                   | 20200232 | 20200301 | 20200308 |
| 20200220                   | 20200233 | 20200302 | 20200309 |
| 20200221                   | 20200238 | 20200304 |          |
| 20200223                   | 20200252 | 20200305 |          |

Armstrong Medical Limited confirms that this Field Safety Notice has been notified to the UK Competent Authority - Medicines and Healthcare Products Regulatory Agency (MHRA).

Armstrong Medical Limited confirms that this Field Safety Notice has been notified to all Competent Authorities, in jurisdictions where the device is made available on the market.

## Field Safety Notice Response Form

FSN Reference: SI21-68 Date: 11 January 2022

Please complete the information below and return to [quality@armstrongmedical.net](mailto:quality@armstrongmedical.net).  
Alternatively, please telephone Armstrong Medical on 00 44 (0)28 70356029 and ask for the Sales Department.

Hospital or Delivery Location Name: \_\_\_\_\_

Hospital or Delivery Location Address: \_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

---

☐ We confirm that we have received this FSN and have distributed it to relevant individuals or departments within our organisation.

☐ **Armstrong Medical Distributors Only:** We confirm that we have received this FSN and have distributed it to all customers that have been supplied with the products listed above.

### Form Completed by:

Name: \_\_\_\_\_

Department or Position: \_\_\_\_\_

e-mail Address: \_\_\_\_\_

Contact telephone number: \_\_\_\_\_

Date: \_\_\_\_\_