Rappel / recall
FSCA-03-2023-1**TURP / POST OPERATIVE Irrigation/aspiration set (WITH PUMP SEGMENT) – For Endoflow® II Two ways**

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*This letter contains important information which require your **immediate attention**.*

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
Dear Customer,

ROCAMED France, manufacturer of the SOUR9200ST12 device reference, is conducting a corrective safety action regarding lot R006497 of TURP / POST OPERATIVE Irrigation/aspiration set (WITH PUMP SEGMENT) – For Endoflow® II Two ways.

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
ROCAMED SAM has received complaints about Urology: TURP / POST OPERATIVE Irrigation/aspiration set (WITH PUMP SEGMENT) – For Endoflow® II Two ways of reference SOUR9200ST12 of the batch R006497.

The fault is on the suction pipe. The part that fits in the suction pump head is missing as well as the poka yoke allowing the tubing to be inserted in the right direction.

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After a review of the batch files and investigations, we were able to confirm a drift in the production process on this batch R006497.

The product specification for reference SOUR9200ST12 uses ONG499206V2 (Suction Urology tubing with pump body), while the production site used ONG499206 (without pump body).

However, as a preventive measure, this FSN will be shared with all customers who have received the impacted batch R006497.

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Any establishment having in stock of reference SOUR9200ST12, of the batch R006497, are requested to return them to us.

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Rocamed has compliant lots available to replace the devices affected by this recall.

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Please then contact your usual contact to arrange for the turn.

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If you have any question about this Field Safety Notice, please contact us:

- at telephone number +377 97 98 41 93 or email mr@promepla.com.
- at telephone number +377 97 98 41 36 or email za@promepla.com
- at telephone number 04 94 90 54 00 or email ay@promepla.com

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ANSM, Swissmedic, SAHPRA, Ministry of Health, FAMHP - Federal Agency for Medicines and Health Products and Market Surveillance & Vigilance, Dutch Health and Youth Care Inspectorate (IGJ).

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In issuing this FSCA, we are aware that we have taken a conservative approach to the issue and we are counting on your understanding and full support to ensure the prompt and effective follow-up of our security action.

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We apologize for any inconvenience this may cause you and your customers and remain at your disposal for any other recommendations or requests.

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Actions to be implemented by distributor, user and clients:

- *Identify and isolate all devices from this batch, wherever they are located, and return them to ROCAMED;*
- *Fill out the acknowledgment form – in Appendix of the FSN - and return this form to ROCAMED within 5 days;*
- *Inform all users and individuals who may have been provided with these devices about this FSN.*

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