

Medtronic (Schweiz) AG

Talstrasse 9
3053 Münchenbuchsee
www.medtronic.com

Tel. 031 868 01 00
Fax 031 868 01 99
E-Mail swisscontact@medtronic.com

Urgent Field Safety Notice

Mahurkar™ Acute Triple Lumen Catheters and Mahurkar™ Acute High Pressure Triple Lumen Catheters

Recall

June 2023

Medtronic Reference: FA1333

EU Manufacturer Single Registration Number (SRN): US-MF-000028763

Dear Risk Manager/Healthcare Professional:

The purpose of this letter is to advise you that Medtronic is voluntarily initiating a recall for specific lots of **Mahurkar™ Acute Triple Lumen Catheters** and **Mahurkar™ Acute High Pressure Triple Lumen Catheters**.

Please note: This recall does **not** include any Mahurkar **Elite** Catheters.

You are receiving this letter as Medtronic records indicate your facility may have at least one of the Mahurkar TM Acute Triple Lumen Catheters and/or at least one of the Mahurkar TM Acute High Pressure Triple Lumen Catheters outlined in Attachment B. Medtronic initiated this action to prevent the use of potentially affected product that may impact patients.

Issue Description:

During manufacturing related testing, the catheter center lumen was found to have an occlusion in the tip of the catheter. Investigation identified the source of the occlusion as excessive MDX, a silicone-based lubricant which coats the catheter tip. As of 10 June 2023, there have been zero (0) confirmed complaints. Additionally, there have been zero (0) reported adverse events and there have been zero (0) reported deaths.

Risk to Health:

An incorrect application of MDX to catheters may result in the hazardous situation whereby the catheter is occluded, partially or fully, and/or uncured or excessive MDX may dislodge from the catheter. An occurrence of the hazardous situation may lead to potential harms identified as full catheter obstruction resulting in delay to treatment and partial obstruction resulting in reduced flow or particulate dislodgement that may result in delay to treatment, hemolysis, embolism/embolus or thrombosis/thrombus.

Patient Recommendation:

MahurkarTM Acute Triple Lumen Catheters and MahurkarTM Acute High Pressure Triple Lumen Catheters are intended for short term use of up to 29 days. For patients with affected lot(s) currently in place, a replacement procedure is recommended. If a patient is found to have a catheter from an affected lot, the patient's medical team should assess the overall patient risk when considering the timing of a replacement. Clinicians should continue to follow current product Instructions For Use (IFU) along with facility specific policies and procedures for routine assessment of the hemodialysis access device for patency, function, and efficacy.

Required Actions:

1. Immediately quarantine and discontinue use of all unused Mahurkar TM Acute Triple Lumen Catheters and Mahurkar TM Acute High Pressure Triple Lumen Catheters referenced in Attachment B - List of affected Lot numbers (see Attachment A for guidance to identify impacted product).

To help you identify if you have affected product, please visit our website www.Medtronic.com/Mahurkar-Triple-Lumen-Catheter-Recall. Here you will find a tool to help you determine if the product you have is affected by this recall.

Please note: This recall does not include any Mahurkar **Elite** Catheters.

2. Return all unused affected product(s) to Medtronic. Your Medtronic Sales Representative can assist in returning any affected product.
3. Please complete the enclosed Customer Acknowledgement Form and email to **rs.dusregulatory@medtronic.com**.
4. This notice should be passed on to all those who need to be aware within your organization or to any organization including but not limited to Nephrologists, Intensivists, physicians, renal nurses, critical care nurses, or other dialysis staff where the potentially affected devices have been transferred. Please maintain a copy of this notice in your records.

Additional Information:

Medtronic has notified the Competent Authority of your country of this action.

We regret any inconvenience this may cause. We are committed to patient safety and appreciate your prompt attention to this matter. If you have any questions regarding this communication, please contact your Medtronic Representative.

Sincerely,

Medtronic (Schweiz) AG

Enclosures:

Attachment A: Identifying Affected Product

Attachment B: List of Affected Lot Numbers

Attachment A: Identifying Affected Product

The **MAHURKAR™* Triple Lumen Catheter** is a radiopaque, polyurethane tube with two clear silicone catheter extensions and three internal lumina distinguished by color-coded adapters on the extensions:

- Red adapter (arterial)
- Blue adapter (venous)
- Clear adapter (medial)

The 12 Fr/Ch catheter is available in various implantable lengths as shown on the box and device labels. The **MAHURKAR™* Triple Lumen Catheter** is intended for short term central venous access for hemodialysis, apheresis, and infusion.

The **MAHURKAR™* Acute High Pressure Triple Lumen Catheter** is a radiopaque, polyurethane tube with two clear silicone catheter extensions and one clear polyurethane Infusion lumen. The three internal lumina can be distinguished by the color-coded luer-lock adapters on the silicone rubber extensions:

- Red adapter (arterial)
- Blue adapter (venous)
- Clear adapter (medial)

The proximal lumen provides "arterial" outflow from the patient; the distal lumen provides "venous" return; the medial lumen is for infusion of fluid, blood products, medications, blood sampling, pressure injection of contrast media and central venous pressure monitoring.



Attachment A: Identifying Affected Product

COVIDIEN™ REF 8888345611

MAHURKAR™

Acute Triple Lumen Catheter Kit

Curved Extensions
12 Fr/Ch (4.0 mm) x 16 cm
Kit de cathéter d'urgence à triple lumière
Extensions courbées

Kit, drieluimiger Katheter für die akute Versorgung, Gefäßmente Extensionen
Kit can catetere a triple lume per trattamento acuto, Prolunga curve
Equipo de catéter de triple lumen para enfermos agudos, Extensiones curvas
Katheterkit med trippelumen för akutvård, Böjda förlängningslängor
Set met catheter met drielubbel lumen voor acute zorg, Gebogen verlengslangen
Kit de cateter de lumen triple para cuidados intensivos, Extensiones curvas
Akutthöjdon tarkoitettu kolmiluiminen katetrissaaja, Kaarettu laajennusosat
Katheterset med tre lumen til akutbehandling, Buede forlængerslanger
Kit katētrā s trojtrijm lūmenem na krīzēkodēbā stāvstīvēst, Zakrīvēnē priedēdņa
Kit katētrā s trojtrijm lūmen, pro akuti pēl, Zahvātē produživaci dily
Souprava katētru s trojtrijm lūmen, pro akuti pēl, Zahvātē produživaci dily
Akut ellātāshoz szolgālő, hármas lūmenő katéterkészlet, ívelt hosszabbító
Комплект трёхпросветного катетера краткосрочного сосудистого доступа, Изогнутые удлинительные трубки
Zestaw cewnika trójkanalowego do krótkotrwałego stosowania, Zakrzywione przewody przedłużające
Akut Bakım Üç Lümenli Kateter Kiti, Eğri Uzatmalar
Kathetersett med tre lumen til akutbehandling, Buede forlængerslanger
Súprava katētra s trojtrijm lūmenem na krīzēkodēbā stāvstīvēst, Zakrīvēnē priedēdņa
Trusi cateter cu trei lumene pentru îngrijirea pacienților cu boli acute, Extensii curbate
Комплект катетер с тремя лумена за краткотрайна употреба, Изогнати удължители

MAHURKAR is a trademark of Saktharam D. Mahurkar, M.D., used under license.
© 2011 Covidien. Made in Costa Rica.
Covidien Inc, 15 Hampshire Street, Mansfield, MA 02048 USA.
Covidien Ireland Limited, IDA Business & Technology Park, Tuilamore.
www.covidien.com
CR8626124-R1

LOT XXXXXXXXXX Use by YYYY-MM-DD
FPO
(01)20884521006789(17)YYMMDD(10)XXXXXXX

STERILE EO

Not made with natural rubber latex

Single use

Rx ONLY

Not made with DEHP

Caution, consult accompanying documents

Non-pyrogenic
Apyrogenen
Pyrogenfrei
Apyrogeno
Pyrogenfrei
Niet-pyrogen
Apyrogenico
Pyrogeniscum
Pyrogenfrei
Mij naprotivova
Apyrogenini
Nem pirogen
Свободные от пирогенов

Keep away from sunlight

CE 0123

5

COVIDIEN™ REF 8888345611HP

MAHURKAR™

Acute High Pressure Triple Lumen Catheter Kit

Curved Extensions
12 Fr/Ch (4.0 mm) x 16 cm
Kit de cathéter d'hémodialyse aiguë triple lumière haute pression
Extensions courbées

Kit mit Hochdruck-Dreilumen-Katheter für die akute Versorgung, Gebogene Extensionen
Kit can catetere tripleme per alta presión, per tratamiento acuto, Prolunga curve
Kit de cateter temporal de triple lumen para altas flujos, Extensiones curvas
Högflödeskatheterkit med trippelumen för akutvård, Böjda förlängningslängor
Kit met hogedruktriple lumen voor acute zorg, Gebogen verlengslangen
Kit de cateter de triple lumen para alta presión em cuidados intensivos, Extensiones curvas
Akutthöjdon tarkoitettu suuren paineen kolmiluiminen katetrissaaja, Käytät jatkopäät
Sart med akut dialyskateter til højt flow med triple lumen, Buede forlængerslanger
Kit katētrā s trojtrijm lūmen, pro akuti pēl, Zahvātē produživaci dily
Souprava vysokotlakého katētru s trojtrijm lūmen, pro akuti pēl, Zahvātē produživaci dily
Akut ellātāshoz szolgālő, hármas lūmenes katéterkészlet, ívelt hosszabbító
Комплект катетера высокого давления для краткосрочного сосудистого доступа, Строгие просветные, Изогнутые удлинительные трубки

MAHURKAR is a trademark of Saktharam D. Mahurkar, M.D., used under license.
© 2011 Covidien. Made in Costa Rica.
Covidien Inc, 15 Hampshire Street, Mansfield, MA 02048 USA.
Covidien Ireland Limited, IDA Business & Technology Park, Tuilamore.
www.covidien.com
CR8637149-01

LOT XXXXXXXXXX Use by YYYY-MM-DD
FPO
(01)20884521128146(17)YYMMDD(10)XXXXXXX

STERILE EO

Not made with natural rubber latex

Single use

Rx ONLY

Not made with DEHP

Caution, consult accompanying documents

Non-pyrogenic
Apyrogenen
Pyrogenfrei
Apyrogeno
Pyrogenfrei
Niet-pyrogen
Apyrogenico
Pyrogeniscum
Pyrogenfrei
Mij naprotivova
Apyrogenini
Nem pirogen
Свободные от пирогенов

Keep away from sunlight

CE 0123

5

COVIDIEN™ REF 8888345611

MAHURKAR™

Acute Triple Lumen Catheter Kit

Curved Extensions
12 Fr/Ch (4.0 mm) x 16 cm

Kit, drieluimiger Katheter für die akute Versorgung, Gefäßmente Extensionen
Kit can catetere a triple lume per trattamento acuto, Prolunga curve
Equipo de catéter de triple lumen para enfermos agudos, Extensiones curvas
Katheterkit med trippelumen för akutvård, Böjda förlängningslängor
Set met catheter met drielubbel lumen voor acute zorg, Gebogen verlengslangen
Kit de cateter de lumen triple para cuidados intensivos, Extensiones curvas
Akutthöjdon tarkoitettu kolmiluiminen katetrissaaja, Kaarettu laajennusosat
Katheterset med tre lumen til akutbehandling, Buede forlængerslanger
Kit katētrā s trojtrijm lūmenem na krīzēkodēbā stāvstīvēst, Zakrīvēnē priedēdņa
Kit katētrā s trojtrijm lūmen, pro akuti pēl, Zahvātē produživaci dily
Souprava katētru s trojtrijm lūmen, pro akuti pēl, Zahvātē produživaci dily
Akut ellātāshoz szolgālő, hármas lūmenő katéterkészlet, ívelt hosszabbító
Комплект трёхпросветного катетера краткосрочного сосудистого доступа, Изогнутые удлинительные трубки
Zestaw cewnika trójkanalowego do krótkotrwałego stosowania, Zakrzywione przewody przedłużające
Akut Bakım Üç Lümenli Kateter Kiti, Eğri Uzatmalar
Kathetersett med tre lumen til akutbehandling, Buede forlængerslanger
Súprava katētra s trojtrijm lūmenem na krīzēkodēbā stāvstīvēst, Zakrīvēnē priedēdņa
Trusi cateter cu trei lumene pentru îngrijirea pacienților cu boli acute, Extensii curbate
Комплект катетер с тремя лумена за краткотрайна употреба, Изогнати удължители

MAHURKAR is a trademark of Saktharam D. Mahurkar, M.D., used under license.
© 2011 Covidien. Made in Costa Rica.
Covidien Inc, 15 Hampshire Street, Mansfield, MA 02048 USA.
Covidien Ireland Limited, IDA Business & Technology Park, Tuilamore.
www.covidien.com
CR8626124-R1

LOT XXXXXXXXXX Use by YYYY-MM-DD
FPO
(01)20884521006789(17)YYMMDD(10)XXXXXXX

STERILE EO

Not made with natural rubber latex

Single use

Rx ONLY

Not made with DEHP

Caution, consult accompanying documents

Non-pyrogenic
Apyrogenen
Pyrogenfrei
Apyrogeno
Pyrogenfrei
Niet-pyrogen
Apyrogenico
Pyrogeniscum
Pyrogenfrei
Mij naprotivova
Apyrogenini
Nem pirogen
Свободные от пирогенов

Keep away from sunlight

CE 0123

5

COVIDIEN™ REF 8888345611HP

MAHURKAR™

Acute High Pressure Triple Lumen Catheter Kit

Curved Extensions
12 Fr/Ch (4.0 mm) x 16 cm

Kit mit Hochdruck-Dreilumen-Katheter für die akute Versorgung, Gebogene Extensionen
Kit can catetere tripleme per alta presión, per tratamiento acuto, Prolunga curve
Kit de cateter temporal de triple lumen para altas flujos, Extensiones curvas
Högflödeskatheterkit med trippelumen för akutvård, Böjda förlängningslängor
Kit met hogedruktriple lumen voor acute zorg, Gebogen verlengslangen
Kit de cateter de triple lumen para alta presión em cuidados intensivos, Extensiones curvas
Akutthöjdon tarkoitettu suuren paineen kolmiluiminen katetrissaaja, Käytät jatkopäät
Sart med akut dialyskateter til højt flow med triple lumen, Buede forlængerslanger
Kit katētrā s trojtrijm lūmen, pro akuti pēl, Zahvātē produživaci dily
Souprava vysokotlakého katētru s trojtrijm lūmen, pro akuti pēl, Zahvātē produživaci dily
Akut ellātāshoz szolgālő, hármas lūmenes katéterkészlet, ívelt hosszabbító
Комплект катетера высокого давления для краткосрочного сосудистого доступа, Строгие просветные, Изогнутые удлинительные трубки

MAHURKAR is a trademark of Saktharam D. Mahurkar, M.D., used under license.
© 2011 Covidien. Made in Costa Rica.
Covidien Inc, 15 Hampshire Street, Mansfield, MA 02048 USA.
Covidien Ireland Limited, IDA Business & Technology Park, Tuilamore.
www.covidien.com
CR8637149-01

LOT XXXXXXXXXX Use by YYYY-MM-DD
FPO
(01)20884521128146(17)YYMMDD(10)XXXXXXX

STERILE EO

Not made with natural rubber latex

Single use

Rx ONLY

Not made with DEHP

Caution, consult accompanying documents

Non-pyrogenic
Apyrogenen
Pyrogenfrei
Apyrogeno
Pyrogenfrei
Niet-pyrogen
Apyrogenico
Pyrogeniscum
Pyrogenfrei
Mij naprotivova
Apyrogenini
Nem pirogen
Свободные от пирогенов

Keep away from sunlight

CE 0123

5

Attachment B: List of Affected Lot Numbers

Product Description	CFN	GTIN	Lot Number			
8888340629 12FR 20CM MAHKA ACUTE 3LCATH	8888340629	20884521057019	1921300069			
8888345611HP 12FR 16CM LUM HP CURV KIT	8888345611HP	20884521128146	2113300288	2203300175	2021000109	
8888345637HP 12FR 24CM LUM HP CURV KIT	8888345637HP	20884521128184	2102600094			