23 novembre 2023

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
www.bd.com/PI-23-4912.

The diagram shows a vertical rectangular label for a Bard Mission biopsy kit. At the top, it says "BARD MISSION" and "Descriptor". Below this, there are three dashed boxes representing size options: "XXg x XXcm", "XXg x XX.Xcm", and "XXg x XX.Xcm".

Two callout boxes on the left point to specific fields on the label:

- Codice prodotto (RIF)** points to the **REF** (Catalogue number) field, which contains "XXXXXXXX".
- Numero di lotto** points to the **LOT** (Lot number) field, which contains "XXXXXXXXXX".

Other fields on the label include:

- Use-By Date**: "YYYY-MM-DD"
- PD** (Penetration depth): "10mm and 20mm"
- Two barcode fields: "01 XXXXXXXXXXXXXXX" and "17 YYYYYY 10 XXXXXXXXXXXXX"
- Description Translations** (empty dashed box)
- Usage instructions: "Consult instructions for Use", "Rx Only", "Not Made with Natural Rubber Latex", "Single Use", "Do Not Use if the Product Sterile Barrier System or its Packaging is Compromised", "Do Not Re-sterilize", "STERILIZED Sterilized Using Ethylene Oxide", "Non-Pyrogenic".
- Manufacturer information: "BD, the BD Logo, Bard and Mission are trademarks of Becton, Dickinson and Company or its affiliates. ©2020 BD. All rights reserved." and "© 2020 BD. All rights reserved."
- Origin: "Made in Dominican Republic"
- Logos: "BARD Biopsy", "CE 2197", and "Authorized Representative in the European Community: BD Switzerland Sàrl, Terre Bonne Park – A4, Route de Crassier 17, 1262 Eysins – Switzerland"
- Manufacturer contact: "Manufacturer: Bard Peripheral Vascular, Inc., 925 West 3rd Street, Tempe, AZ 85281 USA, Tel: 1-800-894-9515, 1-800-321-4254, Fax: 1-800-966-7052, 1-800-640-5376, www.cbard.com/biopsy"
- Reference: "37269579 Rev. 1 12/00"

12th April 2023

PRODUCT NOTIFICATION – PI-23-4640
Bard® Mission™ Disposable Core Biopsy Instrument Kit
REF: See Table 1 Lot Numbers: See Table 1
Type of Action: Product Removal

For the Attention of: Medical Director, Risk Manager, Medical Device Safety Officer, Lab Manager, Clinical Engineering.

Dear customer,

BD is issuing this product notification to remove specific lots of **Bard® Mission™ Disposable Core Biopsy Instrument Kit**. According to our distribution records your organisation may have received the impacted product in Table 1. Product was distributed between August 2022 and September 2022.

Manufacturer's SRN: US-MF-000017556

Product Code (REF)	Lot Number	Expiry Date	UDI
2016MSK (20G x 16CM)	0001458844	28 Feb 2025	(01)00801741097096(17)250228(10)0001458844

Table 1: Impacted product

This product removal is limited to the product code and lot numbers listed in Table 1. No other product codes or lot numbers are affected.

Description of the Problem

BD has received complaints regarding mismatch between the coaxial and the needle in Mission Kit devices, where the coaxial cannula internal diameter and length exceed the external diameter and length of the biopsy needle, prohibiting the biopsy needle from properly fitting into the coaxial cannula and accessing the target tissue.

Potential Risk

If these cannulas are used to perform biopsies, there may be the need to obtain an additional device, the procedure and overall patient care may be prolonged, and there may be insufficient sample acquisition requiring a repeat procedure. Patients may experience tissue injury and discomfort due to need to obtain new access. Providers can mitigate these risks by verifying that the biopsy needle properly fits the coaxial cannula prior to initiating the procedure.

If complications are encountered while using these products, the health consequences are expected to be non-life threatening and providers should apply the standard of care for patient treatment. In the

case of a failed biopsy procedure, any repeat procedures should be performed using a different, clinically appropriate product.

Actions taken by BD

BD is investigating and will implement appropriate measures to prevent recurrence of this product issue.

Customer Actions:

- Cease use of any unused affected **Bard® Mission™ Disposable Core Biopsy Instrument Kit**.
- Identify and quarantine all unused affected **Bard® Mission™ Disposable Core Biopsy Instrument Kit**.
- Make a note of the lot number and destroy all unused affected units.
- Complete and return the Customer Response Form **even if you no longer have any inventory remaining in your facility by 28th April 2023**.
- Circulate this notice to all those who need to be aware within your organization or to any organization where the potentially affected products have been transferred.
- If you experience any issues with **Bard® Mission™ Disposable Core Biopsy Instrument Kit**, please report as a complaint as per your normal process.

Distributor Actions:

- Cease distribution.
- Identify, quarantine, making a note of the lot number then destroy all unused affected **Bard® Mission™ Disposable Core Biopsy Instrument Kit**.
- Identify the facilities where you have distributed affected product and notify them immediately of this notice.
 - Have your customers complete and return the Customer Response form to your organisation for reconciliation purposes by **28th April 2023**.
- Complete and return the Customer Response Form following completion of your reconciliation activities.
- If you experience any issues with **Bard® Mission™ Disposable Core Biopsy Instrument Kit**, please report as a complaint as per your normal process.

	End User with Inventory	End User with ZERO inventory	Where to send completed form
Purchased directly from BD	Complete the form in its entirety. Upon receipt, BD will process the response, and you will receive replacements for unused product	Complete form and check the box indicating “no inventory”	<<insert BD email address>>
Purchased from a distributor/3rd party	Complete all fields on the form and contact your distributor to arrange for replacements as available or credit.	Complete form and check the box indicating “no inventory”	Return the form to your distributor


BD Switzerland Sàrl
Terre Bonne Park – A4
Route de Crassier 17
1262 Eysins – Switzerland
Tél: +41 21 556 30 00
Fax: +41 21 556 30 99
www.BD.com

Contact Reference Person

If you have any questions about this, please contact your local BD representative or the local BD office on <<insert telephone details here>> or e-mail <<insert contact email address here>>.

BD is committed to advancing the world of health. Our primary objectives are patient safety and user safety and providing you with quality products. We apologise for the inconvenience this situation may cause you and thank you in advance for helping BD to resolve this matter as quickly and effectively as possible.

Sincerely,

Lorna Darrock
Associate Director, Post Market Quality
EMEA Quality

Customer Response Form - PI-23-4640

Bard® Mission™ Disposable Core Biopsy Instrument Kit

REF: See Table 1 Lot Numbers: See Table 1

Return to <insert fax/email address here> as soon as possible or **no later than the 28th April 2023**

- I confirm this notice has been read, understood and that all recommended actions have been implemented as required.

Tick the appropriate box below

We do not have any of the affected product as listed in **Table 1** in our facility. Affected product has been used.

All product that is not available for destruction will be considered as dispositioned at your location and therefore physically unavailable unless otherwise specified.

OR

We have the following units of the affected product as listed in **Table 1** in our possession and I confirm that the units have been destroyed *(Please complete the table below with the lot number and the number of units destroyed. Replacement product will only be sent on completion and return of this form).*

REF:	Lot Number/s:	Units destroyed/returned <i>(insert quantity below)</i>

Account/Organisation Name:	
Department <i>(if applicable):</i>	
Address:	
Postcode:	City:
Contact Name:	
Job Title:	
Contact Telephone Number:	Contact E-mail Address:
Name of your supplier for this product <i>(if not direct from BD)</i>	
Signature:	Date:

This form must be returned to BD before this action can be considered closed for your account.

**If you were forwarded this Product notification via a distributor/3rd party, please return your completed form to that organisation for reconciliation purposes.*

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