B. Braun Medical AG

Seesatz 17

CH-6204 Sempach

www.bbraun.ch

cliente

Sempach, 28. novembre 2023

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B. Braun Medical AG

Dr. Karin Birbaum

Stv. FvP & Safety Officer Swiss Sales Divisions

Himjat Bajra

Product Manager Neurosurgery & Interventional Vascular Therapy

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GB Warning

Re-use of single-use devices creates a potential risk of patient or user. It may lead to contamination and/or impairment of functional capability. Contamination and/or limited functionality of the device may lead to injury, illness or death of the patient. Multidirectional stopcocks and manifolds may tend to fissure formation when influenced as follows:

- getting into contact with infusion solutions and drugs
- alcoholic disinfectants
- sudden changes in temperature necessitated by therapy
- residues of blood or thrombocytes as well as precipitates of substances physicochemically incompatible at the junction points.

These fissures may lead to fluid efflux or air entry.

We recommend to examine the entire infusion line regarding these hazards.

In case of doubt, we recommend, especially as to the following drugs, to administer them via a separate stopcock (Discofix® C), which is resistant against drugs:

- Ciclosporin
- Cyclophosphamid
- Etoposid
- Lipids
- Nimodipin
- Phenytoin-sodium
- Propofol

Please observe absolutely the instructions for use and special information of the respective producer.

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