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Date : 12-01-2024

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Email address*	Confirm email address *
Password *	- Confirm Password *
	کة ××××××××××××××××××××
Name *	Last name *
2 Name	😫 Last name
Telephone number.*	Cardio protected area *
	🔂 City / Institution / Area
Serial No. *	Key*
♥ 2007XXXXX	©₽ XXXXXXXXXXXXXXXXXXXXXX
L have read and acc	ept the privacy policy

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Вехеп саявно	Program. Todos los dispositivos 👻								⊕ en → <b>1</b>
Coordination	All the devices								
<ul> <li>Programs</li> <li>Groups</li> </ul>	Devices								Additevice
Contents	Search Search								< +
Devices	Status Name	Serial No.	Model	Software	Usage	Manufacturer	Location	Connected	Description

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----

	2	ſ		3
	Verify the information			
Device data				
Name *	Manufacturer * Bexen Cardio		Reanibex 100	
Serial No. *				
Select the initial location The connected Reanibex 100 loc				
Where is the device?* Location details*				
Latitude * Longitude *				
	Device data Name * Serial No. * Select the initial locatio The connected Reanibex 100 loc Where is the device? *	Fill the require Device data Name * Serial No, * Select the initial location The connected Reanibex 100 location will be automatically Where is the device? * Location details *	Fill the required fields         Device data         Name*         Bexen Cardio         Serial No.*         Select the initial location         The connected Reanibex 100 location will be automatically updated         Where is the device?*         Location details *	Fill the required fields         Device data         Name*       Manufacture**         Bexen Cardio       Reanibex 100         Serial No,*       Rey*         Select the initial location       Key*         The connected Reanibex 100 location will be automatically updated         Where is the device?*       Location details*

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×

DEVICE STATUS · Revision re	quired	
A Electrodes not connected during tes	t	
Last synchronization · 2023-10-19	11:37:55	
USB SYNC		15
Test Type · Power On Test		4.
Accessories		
Battery	Pads	
94% Battery Charge	Not Connected	Qp
		ACCEPT

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Bexen CARDIO			() EN 🛩 🚺
Condination	Résource center		
<ul> <li>Programs</li> <li>Groups</li> </ul>	Bexen Cardio Reso	ources	
Contents	Reanibex USB SYNC	Reanibex Alert	Reanibex training
Co Devices	PC application to synchronize the Reanibex 100 defibrillator. Through this synchronization.	Mobile application to check the status of the defibrillator and its accessories without the	Web application for the training of volunteer communities.
Lt Users	changes in configuration, updates, or data transfer occur.	need for a PC connection: valid only for Android devices.	
Swpport	Download App	Download App	View Web

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¢Sy	nchronization process of Reanibex 100
	1. Plug the REANIBEX 100 to the computer using the USB cable as shown in the image. Check device

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Δ	The synchronization proccess may take up to 30 minutes. Please, do not disconnect the device until all the steps are compl	eted.	
1) Synchronizing Episodes			0
2 Synchronizing Tests			0
3 Synchronizing Settings	Synchronization Status		0
Checking pending Updates	The Synchronization of files has finished successfully. Please, press the "Start Update" button in order to finish the update process.		Ø
Downloading Software	Start Update	Not recessary-	E
Downloading Settings			0
Downloading Audios		Not necessary	E
Verifying Files			ø

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If the device plays an audio requesting to connect the paids plassic connect the paids or the update will not be performed. If the device plays an audio requesting to replace the paids plassic connect new pads or the update will not be performed.
Once in test mode, please, click the "Next" button below.
The update proccess may take several minutes. Do not disconnect the REANIBEK 100 battery nor turn the device off during the update proccess or the device may become unusable.
(s) "Update In Progress"
"Update In Progress"  Extended remaining time: 08:57      The the test is completed, a message will be issued inducting the result. If the device has communications available, the test rend will also be sent (the may take usered mondes).
"Update In Progress" Extracted remaining time: 08.57      The first the sense will be insued indicating the result. If the device has communications available, the test     ment will also be sent (this may take several minutes).     The force will start the update process. This process will take several minutes. After the update, the device will issue a     message indicating the result of the update (Device many or update out performed).
Source the test is completed, a message will be issued indicating the result. If the device has communications available, the test read will also be sent (this may take saveral minutes).     The device will start the update process. This process will take serveral minutes. After the update, the device will issue a message indicating the result of the update (Device mady or update end performed).     The update is successfuld, the device will automatically perform a last and sum this result of the:
So      S

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connect the REANIBEX 1	00 to the computer using the USB cable, and click on the "Verify Update" b start the validation.	utton t
	(?) Verify Update	

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