

Unimax Medical Systems Inc.
8F-2, No. 127, Lane 235, Pao Chiao Road, Xindian Dist,
New Taipei City 231, Taiwan
Tel: +886-2-8919-1698
www.unimaxmeds.com

Rev 2: February 2020
FSN Ref: FSN-EP24002

FSCA Ref: FCA-EP24002

Date: 2024.02.05

Field Safety Notice
Detachable Endo Retrieval Pouch

For Attention of*: Theatre Manager

Contact details of local representative (name, e-mail, telephone, address etc.)*
Mölnlycke Health Care AG Contact's Name: Natalia Kasperowicz Email: info.ch@mölnlycke.com Phone: +41 (0) 44 744 54 00 Country: Switzerland Address: Brandstrasse 24, 8952 Schlieren – Switzerland Post Code: 8952

Rev 2: February 2020
FSN Ref: FSN-EP24002

FSCA Ref: FCA-EP24002

Field Safety Notice (FSN) **Detachable Endo Retrieval Pouch**

1. Information on Affected Devices*																																																																								
1.	1. Device Type(s)*																																																																							
	Detachable Endo Retrieval Pouch Small (250-300ml) / 10mm introducer diameter Medium /Large (500-700ml)/ 10mm introducer diameter Extra Large (1150-1500ml)/12mm and 15 mm introducer diameter																																																																							
1.	2. Commercial name(s)*																																																																							
	Detachable Endo Retrieval Pouch																																																																							
1.	3. Unique Device Identifier(s) (UDI-DI)																																																																							
	07323190272792 (model 899102) 07323190272808 (model 899103) 07323190272815 (model 899104) 07323190272907 (model 899112)																																																																							
1.	4. Primary clinical purpose of device(s)*																																																																							
	The detachable endo pocket is a device that is used to collect and extract specimens during laparoscopic surgery.																																																																							
1.	5. Device Model/Catalogue/part number(s)*																																																																							
	899102; 899103; 899104; 899112																																																																							
1.	6. Affected serial or lot number range																																																																							
	<table><tr><td colspan="6">Model no# 899102</td></tr><tr><td>6252207126</td><td>6252208039</td><td>6252208096</td><td>6252208124</td><td>6252209038</td><td>6252209071</td></tr><tr><td>6252209263</td><td>6252207098</td><td>6252210099</td><td>6252212169</td><td></td><td></td></tr><tr><td colspan="6">Model no# 899103</td></tr><tr><td>6252207125</td><td>6252208121</td><td>6252208125</td><td>6252209039</td><td>6252209072</td><td>6252209264</td></tr><tr><td>6252211001</td><td>6252211014</td><td>6252212024</td><td>6252212170</td><td>6252301132</td><td></td></tr><tr><td colspan="6">Model no# 899104</td></tr><tr><td>6252209041</td><td></td><td></td><td></td><td></td><td></td></tr><tr><td colspan="6">Model no# 899105</td></tr><tr><td>6252207131</td><td>6252208101</td><td>6252208132</td><td>6252209053</td><td>6252209080</td><td>6252211005</td></tr><tr><td>6252212026</td><td>6252212178</td><td></td><td></td><td></td><td></td></tr></table>						Model no# 899102						6252207126	6252208039	6252208096	6252208124	6252209038	6252209071	6252209263	6252207098	6252210099	6252212169			Model no# 899103						6252207125	6252208121	6252208125	6252209039	6252209072	6252209264	6252211001	6252211014	6252212024	6252212170	6252301132		Model no# 899104						6252209041						Model no# 899105						6252207131	6252208101	6252208132	6252209053	6252209080	6252211005	6252212026	6252212178				
Model no# 899102																																																																								
6252207126	6252208039	6252208096	6252208124	6252209038	6252209071																																																																			
6252209263	6252207098	6252210099	6252212169																																																																					
Model no# 899103																																																																								
6252207125	6252208121	6252208125	6252209039	6252209072	6252209264																																																																			
6252211001	6252211014	6252212024	6252212170	6252301132																																																																				
Model no# 899104																																																																								
6252209041																																																																								
Model no# 899105																																																																								
6252207131	6252208101	6252208132	6252209053	6252209080	6252211005																																																																			
6252212026	6252212178																																																																							

2. Reason for Field Safety Corrective Action (FSCA)*	
2.	1. Description of the product problem*
	The mechanism of the listed article number operates in a way that the tube within detaches during the removal process. If the tube is not precisely fixed, part of the tube may stretch out from the opening after detachment and fall into the abdomen of the patient.
2.	2. Hazard giving rise to the FSCA*
	The reported incidence is potentially serious to patients as the extending part may fall into the cavity
2.	3. Probability of problem arising
	Overall occurrence rate: within 0.0001
2.	4. Predicted risk to patient/users
	Prolonged surgery or surgical intervention

Rev 2: February 2020

FSN Ref: FSN-EP24002

FSCA Ref: FCA-EP24002


2.	5. Background on Issue
	The device is used to contain and remove specimen removed during laparoscopic surgery. The mechanism of the listed article number operates in a way that the tube within detaches during the removal process. If the tube is not precisely fixed, part of the tube may stretch out from the opening after detachment and fall into the abdomen of the patient. It was thus decided to proceed with a field safety corrective action to replace the current version with an improved design variant thus reducing the potential for the tube stretching out / falling into the patients abdomen.

3. Type of Action to mitigate the risk*		
3.	1. Action To Be Taken by the User* <input checked="" type="checkbox"/> Identify Device <input checked="" type="checkbox"/> Quarantine Device <input type="checkbox"/> Return Device <input checked="" type="checkbox"/> Destroy Device <input type="checkbox"/> On-site device modification / inspection <input type="checkbox"/> Follow patient management recommendations <input type="checkbox"/> Take note of amendment / reinforcement of Instructions For Use (IFU) <input type="checkbox"/> Other <input type="checkbox"/> None	
3.	2. By when should the action be completed?	Estimated within 6 months
3.	3. Is customer Reply Required? * (If yes, form attached specifying deadline for return)	Yes
3.	4. Action Being Taken by the Manufacturer* <div style="display: flex; justify-content: space-between;"> <div> <input checked="" type="checkbox"/> Product Removal <input type="checkbox"/> Software upgrade <input type="checkbox"/> Other </div> <div> <input type="checkbox"/> On-site device modification/inspection <input type="checkbox"/> IFU or labelling change <input type="checkbox"/> None </div> </div>	
3.	5. By when should the action be completed?	- Distribution of the FSN scheduled for 1-3 weeks. - Expected return of customer reply form in 1-2 months - Replacement of affected devices scheduled for 3-6 months depends on the quantity.
3.	6. Is the FSN required to be communicated to the patient /lay user?	No

Rev 2: February 2020

FSN Ref: FSN-EP24002

FSCA Ref: FCA-EP24002

4. General Information*		
4.	1. FSN Type*	New
4.	2. Further advice or information already expected in follow-up FSN? *	Not planned yet
4.	3. Manufacturer information (For contact details of local representative refer to page 1 of this FSN)	
	a. Company Name	Unimax Medical Systems Inc.
	b. Address	8F-2, No. 127, Lane 235, Pao Chiao Road, Xindian Dist., New Taipei City, Taiwan
	c. Website address	http://www.unimaxmeds.com
4.	4. The Competent (Regulatory) Authority of your country has been informed about this communication to customers. *	
	Yes.	
4.	5. List of attachments/appendices:	If extensive consider providing web-link instead.
4.	6. Name/Signature	Magellan C.L. Yu / Regulatory Affairs Specialist
		

Transmission of this Field Safety Notice	
	<p>This notice needs to be passed on all those who need to be aware within your organisation or to any organisation where the potentially affected devices have been transferred. (As appropriate)</p> <p>Please transfer this notice to other organisations on which this action has an impact. (As appropriate)</p> <p>Please maintain awareness on this notice and resulting action for an appropriate period to ensure effectiveness of the corrective action.</p> <p>Please report all device-related incidents to the manufacturer, distributor or local representative, and the national Competent Authority if appropriate, as this provides important feedback.*</p>

Note: Fields indicated by * are considered necessary for all FSNs. Others are optional.



Rev 1: July 2018

Unimax Medical Systems Inc.

8F-2, No. 127, Lane 235, Pao Chiao Road, Xindian Dist,

New Taipei City 231, Taiwan

Tel: +886-2-8919-1698

www.unimaxmeds.com

Customer Reply Form

1. Field Safety Notice (FSN) information	
FSN Reference number*	FSN-EP24002
FSN Date*	2024.02.05
Product/ Device name*	Detachable Endo Retrieval Pouch
Product Code(s)	1. 899102 2. 899103 3. 899104 4. 899112
Batch/Serial Number (s)	1. 899102: 6252207126/6252208039/6252208096/6252208124/ 6252209038/6252209071/6252209263/6252207098/ 6252210099/6252212169. 2. 899103: 6252207125/6252208121/6252208125/6252209039/ 6252209072/6252209264/6252211001/6252211014/ 6252212024/6252212170/6252301132 3. 899104: 6252209041 4. 899114: 6252207131/6252208101/6252208132/6252209053/ 6252209080/6252211005/6252212026/6252212178

2. Customer Details	
Account Number	
Healthcare Organisation Name*	
Organisation Address*	
Department/Unit	
Shipping address if different to above	
Contact Name*	
Title or Function	
Telephone number*	
Email*	

3. Customer action undertaken on behalf of Healthcare Organisation				
<input type="checkbox"/>	I confirm receipt of the Field Safety Notice and that I read and understood its content. I do not have any affected devices.			
<input type="checkbox"/>	I confirm receipt of the Field Safety Notice and that I read and understood its content. I have destroyed the affected single packed devices. I have completed the table with the details of affected devices quantity, its article and lot/batch number.	Quantity	Article/Material Number	Lot/Batch Number
		N/A	Comments	
Print Name*		Customer print name here		
Signature*		Customer sign here		
Date*				

4. Return acknowledgement to sender	
Email	ra@unimaxmeds.com
Customer Helpline	+886-2-8919-1698
Postal Address	Unimax Medical Systems Inc. 8F-2, No. 127, Lane 235, Pao Chiao Rd., Xindian District, New Taipei City, Taiwan
Web Portal	http://www.unimaxmeds.com/
Fax	+886-2-8919-1528
Deadline for returning the customer reply form*	Within 1-2 months

Mandatory fields are marked with *

It is important that your organisation takes the actions detailed in the FSN and confirms that you have received the FSN.

Your organisation's reply is the evidence we need to monitor the progress of the corrective actions.