



COOK MEDICAL EUROPE LTD.  
O'HALLORAN ROAD  
NATIONAL TECHNOLOGY PARK  
LIMERICK, V94 N8X2, IRELAND  
TEL: +353 61 334440 FAX: +353 61 334441  
WWW.COOKMEDICAL.EU

FSN & FSCA Ref: 2024FA0001

Date: 28 February 2024

**Urgent Field Safety Notice**  
**Dilators - Coons Taper**  
**Entuit™ Secure Adjustable Gastrointestinal Suture Anchor Sets**  
**Chiba Biopsy Needles**

For Attention of: Chief Executive / Risk Management / Purchasing

**Contact details of local representative (name, e-mail, telephone, address etc.)**

Cook Medical Europe Ltd.  
O'Halloran Road  
National Technology Park  
Limerick, Ireland  
E-mail: [European.FieldAction@CookMedical.com](mailto:European.FieldAction@CookMedical.com)  
Phone: Please refer to the attached Country Contacts List

For any further information or support concerning the information within this FSN,  
please contact your local Cook Medical Sales Representative or Cook Medical Europe  
Ltd.



COOK MEDICAL EUROPE LTD.  
O'HALLORAN ROAD  
NATIONAL TECHNOLOGY PARK  
LIMERICK, V94 N8X2, IRELAND  
TEL: +353 61 334440 FAX: +353 61 334441  
WWW.COOKMEDICAL.EU

FSN & FSCA Ref: 2024FA0001

## **Urgent Field Safety Notice**

### **Dilators - Coons Taper Entuit™ Secure Adjustable Gastrointestinal Suture Anchor Sets Chiba Biopsy Needles**

#### **Risk Addressed by FSN**

<b>1. Information on Affected Devices</b>	
1.	<b>1. Device Type(s)</b> Dilators are designed with gradual tapers, may include hydrophilic coating, and are available in multiple inner and outer diameter and length combinations.  The Entuit™ Secure Adjustable Gastrointestinal Suture Anchor Set consists of two or three introducer needles with preloaded anchors and a .035 inch wire guide.  The Chiba Biopsy Needle is ideal for initial access and biopsy. The match-ground bevelled tip provides steerability. The Chiba needle that features EchoTip® enhances visibility under ultrasound.
	<b>2. Commercial name(s)</b> Dilators - Coons Taper Entuit™ Secure Adjustable Gastrointestinal Suture Anchor Sets Chiba Biopsy Needles
	<b>3. Primary clinical purpose of device(s)</b> Dilators - Coons Taper are intended to be used for dilating puncture sites or catheter tracts for percutaneous placement of devices for vascular and non-vascular applications such as in the venous, arterial, biliary and renal systems.  Entuit™ Secure Adjustable Gastrointestinal Suture Anchor Sets are intended for anchoring the anterior wall of the stomach to the abdominal wall prior to introduction of interventional catheters.  Chiba Biopsy Needles are intended for used for aspiration biopsy.
1.	<b>4. Device Model/Catalogue/Part Number(s)</b> <u>Dilators - Coons Taper</u> Reference Part Number (RPN): JCD22.0-38-20-COONS Order Number (GPN): G04443  <u>Entuit™ Secure Adjustable Gastrointestinal Suture Anchor Sets</u> Reference Part Number (RPN): GIAS-SRM-ADJ-2 Order Number (GPN): G35562  <u>Chiba Biopsy Needles</u> Reference Part Numbers (RPNs): DCHN-22-15.0, DCHN-22-15.0-U Order Numbers (GPNs): G00012, G03314
	<b>5. Affected serial or lot number range</b>
	Dilators - Coons Taper: 15767086 Entuit™ Secure Adjustable Gastrointestinal Suture Anchor Sets: 15791906 Chiba Biopsy Needles: 15785990, 15784638, 15784722



COOK MEDICAL EUROPE LTD.  
O'HALLORAN ROAD  
NATIONAL TECHNOLOGY PARK  
LIMERICK, V94 N8X2, IRELAND  
TEL: +353 61 334440 FAX: +353 61 334441  
WWW.COOKMEDICAL.EU

FSN & FSCA Ref: 2024FA0001

2. Reason for Field Safety Corrective Action (FSCA)		
2.	<b>1. Description of the product problem</b> Cook Medical identified that devices from the affected device lots may have packaging seals that do not meet peel strength specifications. This creates a potential of a sterile barrier failure that could result in compromised device sterility.  You are receiving this letter as Cook Medical records indicate that impacted products were shipped to your facility.	
2.	<b>2. Hazard giving rise to the FSCA</b> The affected devices may be non-sterile or contaminated with microorganisms. Potential adverse events that may occur if an affected product is used include increased procedural time and infection, potentially being life-threatening and/or requiring medical/surgical intervention.  To date, Cook Medical has not received any customer complaints related to the adverse patient effects listed above for the affected lots. However, please be advised that compromised device sterility may go undetected by the user.	
3. Type of Action to Mitigate the Risk		
3.	<b>1. Actions To Be Taken by the User</b> <input checked="" type="checkbox"/> Identify Device(s) <input checked="" type="checkbox"/> Quarantine Device(s) <input checked="" type="checkbox"/> Return Device(s) to Cook Medical <input checked="" type="checkbox"/> Other  Please complete the enclosed Customer Reply Form. Where devices are indicated as being returned, our Customer Services department will contact you to organize the return and issue you with the relevant Returns Authorization number. Please include contact details on the Customer Reply form.  Returned Device(s) should be addressed to: Cook Medical EUDC Robert-Koch-Straße, 2 52499 Baesweiler GERMANY  Credit will be provided for the returned affected device(s) where applicable.	
3.	<b>2. Is Customer Reply Required?</b> Form is attached specifying deadline for return.	Yes
3.	<b>3. Action Being Taken by the Manufacturer</b> <input checked="" type="checkbox"/> Product Removal	
3.	<b>4. Patient Management</b> Physicians should follow their institution's protocols/guidelines for the standard of patient care after the procedure for identification and treatment(s) of infection or any complications.	



COOK MEDICAL EUROPE LTD.  
O'HALLORAN ROAD  
NATIONAL TECHNOLOGY PARK  
LIMERICK, V94 N8X2, IRELAND  
TEL: +353 61 334440 FAX: +353 61 334441  
WWW.COOKMEDICAL.EU

FSN & FSCA Ref: 2024FA0001

4. General Information		
4.	1. FSN Type	New
4.	2. Further advice or information already expected in follow-up FSN?	No
4.	3. Manufacturer information Refer to page 1 of this FSN for contact details of local representative.	
	a. Company Name	Cook Incorporated
	b. Address	750 Daniels Way Bloomington, IN 47402, United States
4.	4. The Competent (Regulatory) Authority of your country has been informed about this communication to customers.	
4.	5. Name/Signature	
		Larry D. Pool Director, Post Market Cook Incorporated

Transmission of this Field Safety Notice
<p>This notice needs to be passed on all those who need to be aware within your organisation or to any organisation where the potentially affected devices have been transferred.</p> <p>Please transfer this notice to other organisations on which this action has an impact.</p> <p>Please maintain awareness on this notice and resulting action for an appropriate period to ensure effectiveness of the corrective action.</p> <p>Please report all device-related incidents to the manufacturer, distributor or local representative, and the national Competent Authority if appropriate, as this provides important feedback.</p>



COOK MEDICAL EUROPE LTD.  
O'HALLORAN ROAD  
NATIONAL TECHNOLOGY PARK  
LIMERICK, V94 N8X2, IRELAND  
TEL: +353 61 334440 FAX: +353 61 334441  
WWW.COOKMEDICAL.EU

## Field Action Customer Reply Form

1. Field Safety Notice (FSN) information	
FSN Reference number	2024FA0001
FSN Date	28 February 2024
Product/Device name	Dilator - Coons Taper Entuit™ Secure Adjustable Gastrointestinal Suture Anchor Set Chiba Biopsy Needle
Product Part Number(s)	<u>Dilator - Coons Taper</u> Reference Part Number (RPN): JCD22.0-38-20-COONS  <u>Entuit™ Secure Adjustable Gastrointestinal Suture Anchor Set</u> Reference Part Number (RPN): GIAS-SRM-ADJ-2  <u>Chiba Biopsy Needle</u> Reference Part Numbers (RPNs): DCHN-22-15.0, DCHN-22-15.0-U
Batch/Serial Number(s)	Dilators - Coons Taper: 15767086 Entuit™ Secure Adjustable Gastrointestinal Suture Anchor Sets: 15791906 Chiba Biopsy Needles: 15785990, 15784638, 15784722

  

2. Customer Details	
Account Number	
Healthcare Organisation Name	
Organisation Address	
Contact Name	
Title or Function	
Telephone number	
Email	



COOK MEDICAL EUROPE LTD.  
O'HALLORAN ROAD  
NATIONAL TECHNOLOGY PARK  
LIMERICK, V94 N8X2, IRELAND  
TEL: +353 61 334440 FAX: +353 61 334441  
WWW.COOKMEDICAL.EU

3. Customer action undertaken on behalf of Healthcare Organisation		
Please mark boxes below to indicate actions have been completed.		
<input type="checkbox"/>	I confirm receipt of the Field Safety Notice and that I read and understood its content.	
<input type="checkbox"/>	The information and required actions have been brought to the attention of all relevant users and executed.	
<input type="checkbox"/>	I have affected devices to return - enter Lot number and quantities in table below.	
<input type="checkbox"/>	No affected devices remain in our organisation's inventory	
Print Name		
Signature		
Date		

4. Return acknowledgement to sender	
Email	European.FieldAction@CookMedical.com
Fax	+ 353 61 239294
Deadline for returning the customer reply form	Please return this form within 5 business days of receipt, even if you do not have any of the affected devices.
Customer Helpline	Please refer to the attached Country Contacts List.

If you are returning any affected devices, please indicate the part number, lot number and quantity:

Product Part Number	Product Lot Number	Quantity



COOK MEDICAL EUROPE LTD.  
O'HALLORAN ROAD  
NATIONAL TECHNOLOGY PARK  
LIMERICK, V94 N8X2, IRELAND  
TEL: +353 61 334440 FAX: +353 61 334441  
WWW.COOKMEDICAL.EU

It is important that your organisation takes the actions detailed in the FSN and confirms that you have received the FSN.

Your organisation's reply is the evidence we need to monitor the progress of the corrective actions.