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FSCA/FSN Referenz: 11.03.2024

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Dear Valued Customer,

by choosing MectaCem-X HV with gentamicin and MectaCem-X LV with gentamicin you are applying a top-quality product that is characterized by high-quality and reliable processing properties.

This safety notice is to inform you of the incorrect information present in the Instruction for Uses (IFU) of some batches of bone cements produced between March 2023 and January 2024 and listed in the above table.

Problem description

During the review of the information contained in the IFU an incorrect information was discovered.

In the antibiotic device IFU there is an error in the list of contraindications which led to the creation of a contraindication formatted from two contraindications, reducing the cases in which MectaCem-X with gentamicin devices are contraindicated.

WRONG

MectaCem-X with gentamicin are contraindicated:

- where the loss of musculature or neuromuscular compromise in the affected limb would render the surgical procedure unjustifiable in patients who are allergic or sensitive to any of its components.

CORRECT

MectaCem-X with gentamicin are contraindicated:

- where the loss of musculature or neuromuscular compromise in the affected limb would render the surgical procedure unjustifiable
- in patients who are allergic or sensitive to any of its components.

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Risk evaluation

At the end of the evaluation of the possible impact on the patient, it was concluded that the problem encountered could have harmful effects on the patient where the loss of musculature or neuromuscular compromise in the affected limb would render the surgical procedure unjustifiable in case the patient is allergic or hypersensitive to one of the components of the bone cement.

Therefore, it was evaluated to notify the customer/hospital of the error on the leaflet. Other actions were not considered necessary because the problems of allergy or hypersensitivity to bone cement components are known.

Actions requested from the Customer/Hospital

- Note the contraindication error in leaflet.
- Share this communication with anyone who needs to know within your organization and any organization to which the affected devices have been delivered.
- Send to the sender (TECRES / MEDACTA) confirmation of receipt and FSN application by 19/04/2024.
- File a copy for your records.

Further information

TECRES has already updated the IFU so that future batches do not have this problem.

The Italian Competent Authority has been informed of the contents of this safety notice.

TECRES's priority is patient safety.

TECRES thanks you for the attention dedicated to the contents of this communication. For any further information or clarification, please contact quality@tecres.com.

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soggetta ad attività di direzione e coordinamento di Demetra Holding S.p.A.