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Anhang 2: Antwortformular (BRF)

URGENT: FIELD SAFETY NOTICE CUSTOMER NOTIFICATION
VARIPULSE Bi-Directional Ablation Catheters

Product name	Product code	LOT
VARIPULSE Bi-Directional Ablation Catheters	D141201	ALL LOTS

Attachment 1: Instructions for Use (IFU) Modification Details

Section: WARNINGS AND PRECAUTIONS

New or Modified Warnings:

- The safety and efficacy of the VARIPULSE™ Catheter when used with a VARIPULSE™ Generator has been determined in the AdmiRE trial in the treatment of subjects with symptomatic paroxysmal atrial fibrillation (PAF). The safety and efficacy of their use in the treatment of other arrhythmias have not been established.
- The safety and effective use of this device outside of the pulmonary veins for the treatment of atrial fibrillation has not been clinically established and may increase the risk of patient injury.
- The safety and effectiveness of catheter ablation for the treatment of atrial fibrillation in patients with significant left ventricular dysfunction, advanced heart failure, substantial left atrial enlargement, and structural heart disease has not been established.
- Implantable pacemakers, implantable cardioverter/defibrillators (ICDs), and metal stents can disrupt the electric field produced by PF ablation.
 - Pay close attention when applying Pulse Field (PF) ablation in patients with implantable pacemakers and implantable cardioverter/defibrillator (ICDs). Implantable pacemakers and ICDs may be adversely affected by PF energy. Safety and Effectiveness have not been evaluated for patients with pacemakers and ICDs. It is important to have temporary external sources of pacing and defibrillation available during ablation and to consider temporarily reprogramming the pacing system to minimum output or OFF mode to minimize the risk of inappropriate pacing. Exercise extreme caution during ablation when in close proximity to atrial or ventricular permanent leads; program the ICD to the OFF mode during the ablation procedure; perform complete analysis after ablation.
- Carefully follow the catheter directions for use. In clinical studies (AdmiRE and InspiRE), the median number of ablations delivered was 23 (69 applications) and 16 (49 applications), respectively. Safety and effectiveness were established in all treated subjects, 75% of whom had fewer than 28 ablations (84 applications) delivered in the AdmiRE clinical study and 17 ablations (52 applications) in the InspiRE.
- The risk of patient injury may increase if the number of ablations performed exceeds those evaluated in clinical testing and device performance may be impacted.

- Maintain appropriate levels of recommended anticoagulation therapy before and after the procedure. The safety of discontinuing anticoagulation therapy following catheter ablation of atrial fibrillation has not been established. Anticoagulation therapy in such patients should be administered in accordance with the 2023 ACC/AHA/ACCP/HRS Guideline for the Diagnosis and Management of Atrial Fibrillation.
- If there is uncertainty regarding the patient's anticoagulation status or rhythm prior to the procedure, evaluate the need to perform a Transesophageal Echocardiogram (TEE) prior to the procedure to confirm absence of mural thrombus and/or thrombus in the left atrial appendage.
- Confirm ACT \geq 350 seconds prior to ablation with the catheter and check every 15-30 minutes while the catheter is in the left atrium. Failure to maintain appropriate anticoagulation levels may increase the risk for thromboembolic events.
- Use both direct imaging guidance (such as fluoroscopy and/or intracardiac echocardiography) and electrogram data to monitor catheter advancement and reduce risk of tissue injury.
 - When using the catheter with conventional systems (such as fluoroscopy or intracardiac echocardiography), with the CARTO™ 3 System, careful catheter manipulation must be performed in order to avoid cardiac damage, perforation, or tamponade. Catheter advancement should be done under direct imaging guidance.
 - Do not use excessive force to advance or withdraw the catheter when resistance is encountered. The firmness of the braided tip dictates that care must be taken to prevent perforation of the heart.
- Unsheathe the catheter from the sheath in the middle of the chamber to prevent entrapment or perforation.
- When using the catheter in or around the atrioventricular valve, take precautions to prevent the catheter from slipping into the ventricles and becoming entangled with the valve. Entanglement may lead to damage to the valve, entrapment of the device, and/or need for surgical intervention.

Section: DIRECTIONS FOR USE

New or Modified Directions for Use

- **Performing PF Ablation with the Catheter**
 - One ablation sequence consists of three applications. Each application has a ten second interval delay. This ablation sequence is programmed in the TRUPULSE™ Generator.
 - Reposition the catheter between each complete ablation (up to 3 consecutive applications) to avoid stacking of ablations in the same location.

- **Workflow for Pulmonary Vein Isolation**
 - For each pulmonary vein, create a concentric ring of lesion using 4 ablations (2 ostial with the contracted loop and 2 antral with an open loop) (12 applications).
 - The catheter position should be co-axial to the vein.
 - Place an electrode catheter in the right ventricle or a ventricular branch of the coronary sinus to be ready for ventricular pacing in the event of vagal response.
 - Perform 1 ostial ablation (3 applications) with the contracted loop at the same catheter position without moving the catheter.
 - Within the ostial area (with the contracted loop), reposition the catheter and perform an additional ablation (3 applications). The new catheter position should cover any potential gaps from the first ostial catheter position.
 - Ensure that all applications are complete (1 complete ablation) (Generator screen reads 100%).
 - If any ablation is not complete (less than 100%), then repeat the incomplete ablation in the same location.
 - After completing the 2 ostial ablations (6 applications with the contracted loop), move the catheter into an antral side (with an open loop) and follow a similar workflow completing 1 ablation (3 applications) at the same catheter position and repositioning to complete the final antral ablation (a total of 6 applications with an open loop).
 - Deliver 4 ablations (12 applications) per vein.
 - In anatomy with a common PV, deliver 8 ablations (24 applications) per common.
 - If isolation is not obtained, apply additional applications of PF ablation where needed.
 - It is recommended to check all rotations of the CARTO™ map to ensure full circumferential co-axial lesions were achieved on all map surfaces.
 - In clinical studies (admIRE and inspIRE), the median number of ablations delivered was 23 (69 applications) and 16 (49 applications), respectively. Safety and effectiveness were established in all treated subjects, 75% of whom had fewer than 28 ablations (84 applications) delivered in the admIRE clinical study and 17 ablations (52 applications) in the InspIRE clinical study.

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Alle Inhalte © Johnson & Johnson AG 2024