Urgent Field Safety Notice

PPC25-02.A.OUS

BN II System

BN ProSpec-System

Atellica NEPH 630 System

Atellica CH Analysator

Atellica CI-Analysator

Titel N Latex FLC LAMBDA Lot 473298/473298A — Falsche Haltbarkeitsinformationen

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Produkte

Assay	Siemens-Materialnummer/UDI	Chargennummer	Datum der Herstellung	Verfallsdatum
N Latex FLC LAMBDA	10482438/ 00842768031427	473298	08-Mai-2024	07-November-2025
		473298A		

Problembeschreibung

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Customer
To whom it may concern

Certificate of Analysis

P.O.	
Material:	Commission-No.
10482438	
Kit-Lot:	Customer-Order-No.
473298	
Total lot size:	Customer-No.
791 PC	
Delivery amount:	Storage at:
	+ 02°C to 08°C
Release date:	
2024-11-08	
Manufacturing date:	Expiration date:
2024-05-08	2025-11-07
Legal Material Number:	
OPJB035	
Released by:	
Stefan Barth	

Material: 10482438
N Latex FLC lambda
Inspection ID-No.: 400000743152

Kit components:

Level	Material-No. Material-Name	Lot	Release date
0010	11523525 N Latex FLC lambda	473298	2024-11-07

Material: 10482438
N Latex FLC lambda
Inspection ID-No.: 400000743152

Inspection results:

Material-No.	Lot
11523525 N Latex FLC lambda	473298

Sensitivity

The difference of absorbance between a specific protein concentration (see remark) must exceed the required lower limit, if it is compared with the blank value.

	Result Unit	Assigned Value	Upper Limit	Lower Limit
	154 Bit	-	-	70 Bit
Inspection remark	Protein concentration: 1,8 mg/L			

Specificity

Recoveries of known concentration samples should be within the range of the given upper and lower limit.

N FLC Control SL1

	Result Unit	Assigned Value	Upper Limit	Lower Limit
	12.7 mg/l	12.1 mg/l	13.9 mg/l	10.3 mg/l
Inspection remark	lot 473580			

Recovery from target value

	Result Unit	Assigned Value	Upper Limit	Lower Limit
	105 %	-	120 %	80 %
Inspection remark	lot 473580			

N FLC Control SL2

	Result Unit	Assigned Value	Upper Limit	Lower Limit
	32.1 mg/l	33.5 mg/l	38.5 mg/l	28.5 mg/l
Inspection remark	lot 473680			

Recovery from target value

Material: 10482438
N Latex FLC lambda
Inspection ID-No.: 400000743152

Material-No.	Lot
11523525 N Latex FLC lambda	473298

	Result Unit	Assigned Value	Upper Limit	Lower Limit
	96 %	-	120 %	80 %
Inspection remark	lot 473680			

Internal QC Control

	Result Unit	Assigned Value	Upper Limit	Lower Limit
	30.2 mg/l	30.6 mg/l	35.2 mg/l	26.0 mg/l
Inspection remark	lot 473475			

Recovery from target value

	Result Unit	Assigned Value	Upper Limit	Lower Limit
	99 %	-	120 %	80 %
Inspection remark	lot 473475			

Reproducibility

When known concentration samples are assayed 3 times within run, CV has to be within the range of the given upper limit.

N FLC Control SL1

	Result Unit	Assigned Value	Upper Limit	Lower Limit
	0.80 %	-	6.00 %	-
Inspection remark	lot 473580			

N FLC Control SL2

	Result Unit	Assigned Value	Upper Limit	Lower Limit
	1.40 %	-	6.00 %	-
Inspection remark	lot 473680			

Internal QC Control

	Result Unit	Assigned Value	Upper Limit	Lower Limit
	0.60 %	-	6.00 %	-
Inspection remark	lot 473475			

This document was generated by means of electronic data system, which was designed and validated to comply with the FDA21 CFR Part 11 electronic records and signature. This document indicates an electronic signature for batch release.
This lot has been tested and approved for release and shipment by the responsible legal manufacturer Quality department.

COA Verification
Stefan Barth / Quality Assurance

Siemens Healthcare Diagnostics Products GmbH
Emil-von-Behring-Str. 76
35041 Marburg
Germany

Customer
To whom it may concern

Certificate of Analysis

P.O.	
Material:	Commission-No.
10482438	
Kit-Lot:	Customer-Order-No.
473298A	
Total lot size:	Customer-No.
865 PC	
Delivery amount:	Storage at:
	+ 02°C to 08°C
Release date:	
2024-11-22	
Manufacturing date:	Expiration date:
2024-05-08	2025-11-07
Legal Material Number:	
OPJB035	
Released by:	
Stefan Barth	

Material: 10482438
N Latex FLC lambda
Inspection ID-No.: 400000745044

Kit components:

Level	Material-No. Material-Name	Lot	Release date
0010	11523525 N Latex FLC lambda	473298	2024-11-07

Material: 10482438
N Latex FLC lambda
Inspection ID-No.: 400000745044

Inspection results:

Material-No.	Lot
11523525 N Latex FLC lambda	473298

Sensitivity

The difference of absorbance between a specific protein concentration (see remark) must exceed the required lower limit, if it is compared with the blank value.

	Result Unit	Assigned Value	Upper Limit	Lower Limit
	154 Bit	-	-	70 Bit
Inspection remark	Protein concentration: 1,8 mg/L			

Specificity

Recoveries of known concentration samples should be within the range of the given upper and lower limit.

N FLC Control SL1

	Result Unit	Assigned Value	Upper Limit	Lower Limit
	12.7 mg/l	12.1 mg/l	13.9 mg/l	10.3 mg/l
Inspection remark	lot 473580			

Recovery from target value

	Result Unit	Assigned Value	Upper Limit	Lower Limit
	105 %	-	120 %	80 %
Inspection remark	lot 473580			

N FLC Control SL2

	Result Unit	Assigned Value	Upper Limit	Lower Limit
	32.1 mg/l	33.5 mg/l	38.5 mg/l	28.5 mg/l
Inspection remark	lot 473680			

Recovery from target value

Material: 10482438
N Latex FLC lambda
Inspection ID-No.: 400000745044

Material-No.	Lot
11523525 N Latex FLC lambda	473298

	Result Unit	Assigned Value	Upper Limit	Lower Limit
Inspection remark lot 473680	96 %	-	120 %	80 %

Internal QC Control

	Result Unit	Assigned Value	Upper Limit	Lower Limit
Inspection remark lot 473475	30.2 mg/l	30.6 mg/l	35.2 mg/l	26.0 mg/l

Recovery from target value

	Result Unit	Assigned Value	Upper Limit	Lower Limit
Inspection remark lot 473475	99 %	-	120 %	80 %

Reproducibility

When known concentration samples are assayed 3 times within run, CV has to be within the range of the given upper limit.

N FLC Control SL1

	Result Unit	Assigned Value	Upper Limit	Lower Limit
Inspection remark lot 473580	0.80 %	-	6.00 %	-

N FLC Control SL2

	Result Unit	Assigned Value	Upper Limit	Lower Limit
Inspection remark lot 473680	1.40 %	-	6.00 %	-

Internal QC Control

	Result Unit	Assigned Value	Upper Limit	Lower Limit
Inspection remark lot 473475	0.60 %	-	6.00 %	-

This document was generated by means of electronic data system, which was designed and validated to comply with the FDA21 CFR Part 11 electronic records and signature. This document indicates an electronic signature for batch release.
This lot has been tested and approved for release and shipment by the responsible legal manufacturer Quality department.

COA Verification
Stefan Barth / Quality Assurance