Together, improving life

August 2025

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Sarah Pastrnak
Product Pipeline Regulatory Leader
W. L. Gore & Associates, Inc.
1505 N Fourth Street
Flagstaff, AZ 86005, USA
T +1 928 864 4146
spastrna@wlgore.com

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IFU Summary of Changes for the GORE® ACUSEAL Vascular Graft

This document summarizes the changes made to the Instructions for Use (IFU) for the GORE® ACUSEAL Vascular Graft. Table 1 provides comparison of verbiage found in the current version and in the updated version.

Table 1: Summary of Changes

Section	Current IFU Text	Updated IFU Text	Reason for Change
Warning	“When using the GORE® ACUSEAL Vascular Graft, avoid using traumatic instruments, handling the graft with excessive force or high rates of force, cutting the graft incorrectly, or placing the graft in an undersized tissue tunnel. These actions may lead to separation of graft layers. Potential consequences include partial or complete occlusion due to hemodynamically significant stenosis or thrombosis and related serious harms, including additional interventions to resolve.”	“When using the GORE® ACUSEAL Vascular Graft, avoid using traumatic instruments, handling the graft with excessive force or high rates of force, cutting or suturing the graft incorrectly, or placing the graft in an undersized tissue tunnel. When performing reinterventions on the GORE® ACUSEAL Vascular Graft, avoid the use of oversized balloons. These actions may lead to separation of graft layers (delamination). Potential consequences include partial or complete occlusion due to hemodynamically significant stenosis or thrombosis and related serious harms, including additional interventions to resolve.”	Suturing incorrectly (ie – not passing suture through all layers of the graft wall) and the use of oversized balloons may lead to delamination.
Potential Device or Procedure Related Adverse Events	No reference to delamination	“ A possible complication which may occur with the use of the GORE® ACUSEAL Vascular Graft is delamination. Delamination may result in partial or complete occlusion due to hemodynamically significant stenosis or thrombosis and related serious harms, including additional interventions to resolve. ”	Delamination is a potential device related adverse event.

Together, improving life

RETURN ACKNOWLEDGMENT FORM

URGENT Medical Device Safety Correction / Field Safety Notice

GORE® ACUSEAL Vascular Graft

Please return this completed Return Acknowledgment Form via email or fax to Sedgwick using the information listed below **within two weeks of receipt of this letter.**

Attn: Event Number 3081: 2017233.08/18/2025.002-C

Fax:

Email Address:

Phone:

Please check and complete the following boxes to acknowledge receipt of this document sent in reference to the Medical Device Safety Correction / Field Safety Notice of the **GORE® ACUSEAL Vascular Graft.**

I have read and I understand the documents

Print Name and Title of Person Completing Form:	Facility/Business Name:
Date the notification was received:	Telephone Number: Email:
Signed*	Date:
<i>*Your signature provides confirmation that you have received and understood this Medical Device Safety Correction / Field Safety Notice.</i>	

It is important that your organization takes the actions detailed in this Medical Device Safety Correction / Field Safety Notice and that you confirm you have received the included documents. Your organization's reply is the evidence we need to monitor the progress of the corrective actions.

MD206058 Attachment 7

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© 2025 W. L. Gore & Associates, Inc.