



## **URGENT FIELD SAFETY NOTICE** **DISTRIBUTORS**

**Commercial name:** Rectal Tubes, Unoversal Drainage Systems, SimpaVac, Sterile Suction Connecting Tubes, Sterile Connecting Pieces, Suction Handles/Sets (FilterFlow™/Deltaflo), Oxygen Catheters, Sterile Nasal Oxygen Cannulas, Sterile Oxygen Connecting Tubes, and Sterile Forceps. (see Attachment 1 for full product details)

**Issue Date:** 04 Jan 2019

**Reference/REF No:** See Attachment 1

**FSCA ID:** 2018-005 Phase 3

**Type of action:** Recall / Product Disposal

**Recall Centre email:** [convatecproductrecall@stericycle.com](mailto:convatecproductrecall@stericycle.com)

**Please note that this action only applies to specific product codes and does not affect all product codes of the above-mentioned products.**

---

### **Description of the problem:**

ConvaTec has voluntarily initiated a recall of specific product codes of the above-mentioned products.

Internal assessment of product packaging integrity has confirmed that these devices are not meeting our expectations or those of our customers. Transportation testing conducted on the product packaging failed, confirming the potential for a breach in the sterile barrier. Using a non-sterile device on a patient may expose the patient to infectious agents increasing the patient risk of developing infection. ConvaTec has not received any reports of incidents related to the packaging seal issue.

### **Product Descriptions**

**Rectal Tubes** are single-use devices which are inserted via the rectum to drain faeces, perform lavage, administer medicine and/or enemas to the lower intestine, for colonoscopy exploration and colonic decompression.

**Unoversal Drainage Systems** are intended to be used for wound drainage by means of low vacuum where a medium or small volume of post-operative bleeding is expected. They are single use devices that are intended for use in short-term treatment.

**SimpaVac** is a canister intended for low vacuum post-operative wound drainage in conventional (open) surgery.

**Sterile Suction Connecting Tubes** are intended to connect two existing devices, (e.g. suction handles and canisters), or as an extension between an existing device and a vacuum source.

**Sterile Connecting Pieces** are intended to connect two existing devices such as suction handles and collection canisters.



**Suction Handle/Sets (FilterFlow™/Deltaflo)** is a device intended primarily for orthopaedic surgical procedures to remove and collect debris and bone fragments chips, cement, marrow and blood clots without blocking the suction device, using a combination of suctioning and filtering.

**Oxygen Catheters** are single patient use devices, requiring short-term oxygen therapy and administer oxygen (in low concentration) directly via the nasopharyngeal or nasal route.

**Sterile Nasal Oxygen Cannulas** are intended for direct administration of low concentrations of oxygen via the bi-nasal route.

**Sterile Oxygen Connecting Tubes** are intended for administration of low/medium concentration of oxygen in oxygen therapy. The connecting/extension tubes channel the oxygen from a stationary or mobile oxygen supply unit to the patient via a dispensing device. Dispensing devices can be e.g. oxygen catheters, nasal oxygen cannulas or facemasks. When using connecting/extension tubes between the patient dispensing device and the supply unit, the patient is allowed more freedom of movement without interruption of the oxygen flow.

**Sterile Forceps:** These are single use tweezers for multiple applications. Used in clinical procedures to grasp, manipulate, maneuver or assist with the removal of debris, particulates or wound tissue.


#### **Product Identification Procedure**

Only the identified product part codes within this notice may have a potential breach in the sterile barrier packaging.

For this reason and to address any potential risk of harm, all the affected products should **not be used**.

The only way to identify affected product is by comparing Product code/REF and LOT/Batch number to the recalled product list (see Attachment 1). There is no other discernable difference between affected and unaffected product.

See Attachment 2 for examples of package labelling that highlights the location of the product code and LOT/Batch number on the device labels. These are located on the primary packaging and the shipping carton. The product code (reference number) is preceded by the word 'REF' and the LOT/Batch number is preceded by

the symbol,  .



**Advice on action to be taken by distributor.**

Our records show that you have taken delivery of affected product. Please follow the steps below:

1. Please examine the enclosed questionnaires. Immediately stop distributing and quarantine all affected products.
2. Please forward copies of the 'Field Safety Notice END USERS' to your customers, asking them to return the affected products to you along with a completed 'Recall Response Form for END USERS'.
3. When the completed Response Form(s) have been returned to you, please complete the 'Recall Response Form for DISTRIBUTORS' and contact the dedicated Recall Centre (see page 4) for further instructions on how to dispose of affected product and to arrange credit. **Products will not be automatically replaced. For replacement product (s), please place new orders with ConvaTec using the usual process. ConvaTec will supply new orders as soon as possible.**
4. Please also return the completed 'Recall Response Form for DISTRIBUTORS', and all 'Recall Response Form for END USERS' to your 'Recall Centre' via Fax/ E-mail.

**PLEASE PROVIDE A COMPLETED RESPONSE AS SOON AS POSSIBLE.**

**Recall Centre:**  
 Email: [convatecproductrecall@stericycle.com](mailto:convatecproductrecall@stericycle.com)

Please continue to report any non-recall related information such as any adverse events involving these products to your distributor or the ConvaTec Customer Care Line (refer to the **Regional ConvaTec contact** list for details).

**Transmission of this Field Safety Notice:**


This notice should be sent to all others who have received the affected devices within your organization or to any organization where the affected devices have been transferred.

ConvaTec is committed to providing quality products and services to our customers and we sincerely apologize for any inconvenience this notice may cause.

If you have any questions relating to this recall, please contact the Recall Centre. For other support, please contact your distributor or local ConvaTec representative (see **Regional ConvaTec contact** list for details).

The relevant National Authorities have been advised about this Field Safety Corrective Action.

**Authorisation:**

<u>Name</u> Duncan Rowley	<u>Title</u> Director, Regulatory Affairs and Quality Assurance, EMEA	<u>Address</u> ConvaTec Limited, First Avenue, Deeside Industrial Park, Deeside, CH5 2NU, U.K.	
<u>Date</u>	04 Jan 2019	<u>Signature</u>	



**Recall Centre: Dedicated help desk for this Field Safety Corrective Action:**

Tel: +44-208 6101 608

Fax: +44-207 6601 568

Email: [convatecproductrecall@stericycle.com](mailto:convatecproductrecall@stericycle.com)

**'Regional  
Conva  
Tec**

**contact' for other support:**

**Austria, Belgium, Croatia, Cyprus, Estonia, France, Germany, Greece, Hungary, Israel, Liechtenstein, Luxemburg, Malta, Netherlands, Portugal, Romania, Switzerland, Ireland**

Tel: +41 (0) 52 630 54 01

Fax: +41 (0) 52 630 54 99

Email: [ccc.customerservice@convatec.com](mailto:ccc.customerservice@convatec.com)

**Denmark**

Tel: +45 4816 7030

Fax: +45 8025 3413

Email: [customerservicenordic@convatec.com](mailto:customerservicenordic@convatec.com)

**Finland**

Tel: +358 (0) 20 7659 630

Email: [mail.fi@convatec.com](mailto:mail.fi@convatec.com)

**Italy**

Tel: 800500190

Email: [clienti.convatec@convatec.com](mailto:clienti.convatec@convatec.com)

**Norway**

Tel: +47 22686095

Fax: +47 80019602

Email: [customerservicenordic@convatec.com](mailto:customerservicenordic@convatec.com)

**Sweden**

Tel: +46 (0)42 332010

Fax: +46 200887486

Email: [customerservicenordic@convatec.com](mailto:customerservicenordic@convatec.com)

**Turkey**

Tel: +90 216 416 52 00/ +90 312 219 83 61

Fax: +90 216 416 28 30 / +90 312 219 53 21

Email: [musterihizmetleri@convatec.com](mailto:musterihizmetleri@convatec.com)

**United Kingdom**

Tel: +44 (0) 1244 832206

Fax: 0800 279 9004

Email: [unomedical-uk.customerservice@convatec.com](mailto:unomedical-uk.customerservice@convatec.com)

**Czech Republic, Iceland, Poland, Slovakia, Spain**

Tel: +41 (0) 52 630 54 01

Fax: +41 (0) 52 630 54 99

Email: [ccc.customerservice@convatec.com](mailto:ccc.customerservice@convatec.com)



**RECALL RESPONSE FORM for DISTRIBUTORS**  
**URGENT FIELD SAFETY NOTICE**  
**PLEASE COMPLETE AND RETURN by Fax/Email**

Consignee of the device:

<b>Consignee Account No:</b>	<i>Column D</i>
<b>Consignee Name:</b>	<i>Column AA</i>
<b>Consignee Address:</b>	<i>Column AA, AF, AC, AE, AD, B</i>

The following products, Rectal Tubes, Unoversal Drainage Systems, SimpaVac, Sterile Suction Connecting Tubes, Sterile Connecting Pieces, Suction Handles/Sets (FilterFlow™/Deltaflo), Oxygen Catheters, Sterile Nasal Oxygen Cannulas, Sterile Oxygen Connecting Tubes, and Sterile Forceps have been distributed to your facility:

<b>Invoice #</b>	<b>Sales Order #</b>	<b>Product Code / REF No.</b>	<b>SAP Code</b>	<b>LOT No.</b>	<b>Quantity Delivered (pieces)</b>
<i>Column AL</i>	<i>Column M</i>	<i>Column R</i>	<i>Column O</i>	<i>Column S</i>	<i>Column X</i>

**Please answer each of the following.**

Have You Distributed the Product Further?      NO       YES\*

\*If YES, have you notified down to your customer?    NO\*       YES

\*If NO explain why not:

1.  We do NOT have any affected product(s).
2.  We have the following affected product(s):

**Record quantity (pieces) for each LOT to be disposed:**

<b>LOT No.</b>	<b>Units on Hand</b>	<b>LOT No.</b>	<b>Units on Hand</b>	<b>LOT No.</b>	<b>Units on Hand</b>	<b>LOT No.</b>	<b>Units on Hand</b>	<b>LOT No.</b>	<b>Units on Hand</b>



LOT No.	Units on Hand	LOT No.	Units on Hand	LOT No.	Units on Hand	LOT No.	Units on Hand	LOT No.	Units on Hand

Please provide details of affected Rectal Tubes, Unoversal Drainage Systems, SimpaVac, Sterile Suction Connecting Tubes, Sterile Connecting Pieces, Suction Handles/Sets (FilterFlow™/Deltaflo), Oxygen Catheters, Sterile Nasal Oxygen Cannulas, Sterile Oxygen Connecting Tubes, and Sterile Forceps products that were distributed to your customers:

Customer Name	Product Code / REF No.	SAP Code	LOT No.	Quantity (pieces)



FORM Completed and Returned From:

<b>Name (CAPITAL LETTERS):</b>	
<b>Position:</b>	
<b>Company Name:</b>	
<b>Address:</b>	
<b>Phone No:</b>	
<b>Signature:</b>	
<b>Date (dd/mmm/yyyy):</b>	

- As requested, I have also provided a copy of all 'Recall Response Forms for END USERS' returned from customers who received affected product (please tick to confirm).

**Recall Centre: Dedicated help desk for this Field Safety Corrective Action:**

Tel: +44-208 6101 608

Fax: +44-207 6601 568

Email: [convatecproductrecall@stericycle.com](mailto:convatecproductrecall@stericycle.com)



**Attachment 1: List showing affected Product Codes/REF and batches.**

REF/CODE	DESCRIPTION	LOT/BATCH Number
13001181	UNO MINIVAC SET REDON DRN CH6(50/200)INT	203693, 208419, 224350, 307673, 320642, 495175
13001182	MINI SET DE REDON CH06	153977, 159115, 160981, 166817, 181085, 182768, 185144, 188275, 201951, 221101, 223429, 226870, 234889, 237479, 300656, 302931, 305252, 310342, 311314, 312600, 314969, 315284, 316732, 318568, 320361
	UNO MINIVAC SET REDON DRN CH6(50/200)	306398, 314307
	UNO MINIVAC SET REDON DRN CH6(50/200)AUS	152327, 152545, 152687, 155400, 155630, 155911, 158207, 160981, 162465, 162489, 165661, 166258, 166358, 175453, 176571, 178879, 179958, 181085, 182768, 183810, 184728, 185144, 186373, 188275, 190373, 195242, 199290, 203008, 206419, 208562, 212316, 213262, 214604, 217367, 219084, 222168, 223429, 224377, 226870, 228221, 230187, 233575, 234889, 236098, 300656, 300738, 301428, 302259, 302931, 303266, 304730, 305252, 306030, 306280, 306398, 307006, 307301, 307795, 308069, 309106, 309349, 309534, 310342, 311314, 312600, 313045, 313428, 314307, 314969, 315284, 316732, 318568, 320361, 497717
13001183	UNO MINIVAC SET REDON DRN CH06(50/200)IT	160560, 164812, 193766, 198434, 211599, 227240, 229179, 233683, 235314, 300442, 303313, 304312, 305057, 310121, 312601, 313802, 314865, 316376, 318264, 318644, 320290, 321372
13001184	UNO MINIVAC SET REDON DRN CH6(50/200)INT	207658, 301585, 310609, 320915
13002181	MINIVAC 40ML M/ULMER DRAEN&TROCER CH08	164467, 197379, 201035, 304750, 305653, 307137, 310731, 314970, 315521, 321009, 321179
	Minivac setti 30ml + dreeni ch8 10kpl	197379, 201035, 214586, 305653, 307137, 310731
	UNO MINIVAC SET ULMER DRN CH8(50/200)INT	164467, 169614, 201035, 214586, 300452, 304750, 305653, 307137, 310731, 314970, 321179
13002182	MINI SET DE REDON CH08	150531, 153978, 154220, 201085, 201999, 207358, 211085, 214570, 216436, 229746, 234149, 237464, 301363, 302288, 302932, 309981, 312007, 312602, 314971, 315185, 320239
	UNO MINIVAC SET ULMER DRN CH8(50/200)AUS	150531, 153978, 154220, 170092, 197394, 198972, 201085, 201999, 207358, 208208, 210460, 211085, 211393, 212920, 213817, 214570, 216436, 217368, 220363, 222534, 225302, 228794, 229746, 234149, 236625, 237464, 300834, 301363, 302288, 302932, 303472, 304975, 309568, 309868, 309981, 311417, 312007, 312602, 320239
13002183	UNO MINIVAC SET ULMER DRN CH08(50/200)IT	153663, 156313, 157366, 160561, 161966, 164813, 167011, 181343, 188132, 188534, 195259, 200400, 238680, 302748, 306199, 312362, 312466, 314215, 314866, 315755, 316377, 318834, 321373
13002184	UNO MINIVAC SET ULMER DRN CH8(50/200)INT	151999, 152000, 173153, 203997, 207659, 211620, 237166, 301586, 320916



REF/CODE	DESCRIPTION	LOT/BATCH Number
13002185	UNO MINIVAC SET ULMER DRN CH8(50/200)INT	208882
13004181	Thyroideasetti 100ml+2xdreeni ch10 10kpl	13002185, 13004181
	UNO THYROIDEA SET (10/60) INT	204329
13004182	UNO THYROIDEA SET (10/60) AUS	150540, 154002, 159420, 300565, 300576, 300580, 300585, 305767, 312476, 312603, 498229
13004184	UNO THYROIDEA SET (10/60) INT	193787, 300236, 300238, 300259, 300278, 300288, 304976
14003181	REKTALSONDE CH 18/40 CM	157678, 174278, 184057, 200654, 227711, 229628, 240906, 247935, 7A00676, 7C04348, 7C04349, 7E00774, 7M02048
	UNO RECTAL TUBE CH18/40CM(100/800) INT	184057, 227711, 7A00676, 7B01025, 7C04348, 7C04349, 7E00774, 7M02048, 8F00459
14003182	Rektální rourka CH 18/40 cm	6K05088, 7A05626, 7B01328, 7C01496
	UNO RECTAL TUBE CH18/40CM (100/800) AUS	150326, 150687, 151669, 152641, 154319, 159650, 162041, 162772, 162938, 163257, 163796, 164676, 165268, 167943, 168513, 172556, 174978, 177810, 177921, 181572, 181702, 183713, 184650, 185359, 186978, 187841, 188416, 190304, 191210, 193641, 194580, 196712, 199775, 200279, 202361, 203908, 204650, 207571, 209153, 210098, 210785, 212199, 213186, 214930, 215386, 217631, 218362, 218994, 219581, 220593, 222132, 222876, 223671, 226312, 227137, 228094, 229098, 231069, 234459, 235207, 236045, 236998, 238547, 239669, 240551, 241962, 246823, 247488, 248275, 249623, 250476, 251351, 252088, 253684, 497919, 6K05088, 6K05743, 6L01515, 6L06087, 7A00674, 7B01328, 7C01496, 7C05505, 7E00773, 7E02113, 7F02484, 7G00079, 7H01918, 7H03195, 7J02311, 7K03634, 7L03177, 7M02047, 8A02901, 8A05957, 8B02381, 8C02338, 8D00542, 8E01926, 8E02698, 8F00458
14003183	SONDA RECTAL CH/18 40CM-100UDS	150049
	UNO RECTAL TUBE CH18/40CM(100/800) INT	150049, 172997
14003185	UNO RECTAL TBE CH18/40CM(100/800) INT	171461
14006181	REKTALSONDE CH 20/40 CM	227724, 229633, 6K06899, 8C03003
	UNO RECTAL TUBE CH20/40CM(100/600) INT	227724, 229633, 6K06899, 8C03003
14006182	Rektální rourka CH 20/40 cm	6K06894

REF/CODE	DESCRIPTION	LOT/BATCH Number
	UNO RECTAL TUBE CH20/40CM(100/600) AUS	150343, 150689, 152764, 154328, 155618, 156616, 157300, 158189, 158883, 159643, 162483, 162768, 163877, 164755, 165259, 172462, 174977, 177914, 181648, 183700, 185354, 186976, 188404, 191988, 192380, 194572, 196012, 197227, 197575, 199177, 199767, 200271, 208486, 209149, 213182, 216008, 218359, 220588, 222872, 223667, 225600, 227133, 228088, 231065, 232064, 238544, 239663, 240549, 245290, 246817, 249618, 251344, 253680, 497930, 499171, 6K06894, 6L01519, 6L06646, 6M02324, 7A04130, 7B06812, 7C01493, 7C05508, 7E03907, 7E04955, 7F02482, 7G00078, 7H01917, 7J03377, 7K03632, 7L03172, 8A03196, 8B02388, 8C02339, 8D01513
<b>14006183</b>	SONDA RECTAL CH/20 40CM-100UDS	172395, 176496, 254028, 7A00673, 7E04954, 7F00109, 7F02481
	UNO RECTAL TUBE CH20/40CM(100/600) INT	152995, 173237, 180384, 7F00109, 8D02810
<b>14008181</b>	REKTAALIPUTKI CH22/40CM PHTAL.TON 100KPL	154963, 157679, 165481, 169160, 177763, 179313, 180825, 185122, 190723, 192442, 193195, 201917, 208990, 211326, 211838, 216338, 217359, 221017, 221820, 224306, 226091, 226756, 227712, 236592, 241665, 246362, 247936, 6L06125, 7A04716, 7D01706, 7F00139, 7M00168, 8B01946
	REKTALSONDE CH 22/40 CM	152380, 154963, 157679, 165481, 169160, 170441, 174392, 176429, 177763, 179313, 180825, 183348, 185122, 190723, 192442, 193195, 194968, 197248, 198952, 201917, 208990, 216338, 217359, 221017, 221820, 224306, 226091, 226756, 227712, 236592, 241665, 246362, 247936, 250034, 6K06897, 6K06902, 6L01269, 6L04561, 7A04716, 7D01706, 7F00139, 7M00168, 8B01946
	UNO RECTAL TUBE CH22/40CM(100/600) INT	152380, 154963, 157679, 165481, 170441, 174392, 183348, 190723, 192442, 193195, 208990, 211326, 211838, 214486, 216338, 217359, 217984, 224306, 226091, 226756, 227712, 228716, 236592, 241665, 246362, 247936, 250034, 6L01269, 6L06125, 7A04716, 7D01706, 7F00139, 7M00168, 8B01946, 8F01629
<b>14008182</b>	Rektální rourka CH 22/40 cm	7B02810, 7C02293

REF/CODE	DESCRIPTION	LOT/BATCH Number
	UNO RECTAL TUBE CH22/40CM(100/600) AUS	150344, 150690, 151676, 151885, 152484, 152651, 154329, 155220, 155891, 156617, 158190, 158707, 159645, 162038, 162587, 162769, 162933, 164756, 166387, 166562, 167053, 167966, 168510, 170229, 172463, 173369, 176052, 177915, 178732, 179452, 181568, 183701, 184644, 185355, 186339, 187833, 188405, 190298, 191193, 191989, 193637, 194573, 195181, 196013, 197576, 199178, 199768, 200273, 205503, 208487, 209150, 210096, 210782, 212195, 213183, 214926, 216009, 217629, 218360, 218991, 219909, 221394, 222130, 222873, 223668, 224704, 225601, 226311, 227134, 228089, 228650, 229095, 230105, 231066, 232686, 233410, 234455, 234794, 235205, 236043, 238545, 240550, 241957, 245291, 246818, 247486, 248273, 250471, 251196, 252082, 253063, 253681, 497628, 6K06893, 6L04528, 6L04581, 6L06127, 6M00666, 6M01803, 7A01806, 7A04717, 7B02810, 7C03435, 7C04988, 7D01705, 7D04861, 7E00997, 7E02944, 7F00120, 7F04252, 7G01830, 7H00927, 7H03094, 7H05988, 7K01255, 7K02643, 7K04185, 7L05966, 7M02406, 8A04413, 8B02896, 8C01998, 8D03387, 8E02387, 8E02991, 8F01631
<b>14008183</b>	SONDA RECTAL CH/22 40CM- 100UDS	169151, 184033, 185119, 199202, 204247, 205148, 208916, 211237, 214468, 216334, 217974, 220995, 224229, 227697, 230603, 238050, 239181, 240885, 250008, 498358, 6K01000, 6L04555, 6M00670, 7A04715, 7G01831, 7J01188, 8B02895, 8D02036, 8E02386
	UNO RECTAL TUBE CH22/40CM(100/600) IT	152989, 153632, 156832, 169151, 172957, 176101, 180366, 192005, 199202, 206923, 209173, 240567, 498358, 6K01000, 8B02895, 8F01630
<b>14010182</b>	Rektální rourka CH 24/40 cm	207255, 234733, 6K06895
	UNO RECTAL TUBE CH24/40CM(100/600) INT	150686, 152757, 154895, 156578, 156958, 157216, 158780, 159584, 163244, 169401, 176018, 177889, 190274, 193611, 194545, 195164, 200161, 203898, 207255, 207551, 208175, 211489, 218970, 220571, 221384, 223653, 225585, 228064, 233384, 239645, 248266, 251322, 253042, 6K06895, 7B01057, 7C05994, 7D03500, 7E04516, 7F06008, 7H05991, 7K06060, 8A02033, 8D02018, 8D03385, 8E04228
<b>14010183</b>	SONDA RECTAL CH/24 40CM- 100UDS	169182, 254029, 7H04608
	UNO RECTAL TUBE CH24/40CM(100/600) INT	254029, 8C05204
<b>14012182</b>	Rektální rourka CH 25/40 cm	7B06456

REF/CODE	DESCRIPTION	LOT/BATCH Number
	UNO RECTAL TUBE CH25/40CM(100/600) AUS	151677, 152315, 152485, 152652, 152765, 154597, 155221, 155619, 156618, 156816, 157231, 158708, 158884, 159646, 160506, 162040, 162770, 162934, 163256, 164757, 165260, 166388, 166943, 167732, 167942, 168511, 171897, 172548, 173370, 174969, 176054, 177916, 179453, 181569, 183702, 184645, 184925, 185356, 186977, 187834, 188406, 190299, 191195, 193638, 194574, 195182, 196014, 196706, 197577, 199179, 200274, 202359, 203596, 203906, 204649, 205439, 205504, 207569, 209151, 210783, 212196, 213184, 214927, 215385, 216010, 218361, 218992, 220589, 222131, 222606, 222874, 223669, 225602, 227135, 228090, 229096, 230106, 231067, 232687, 233411, 234456, 235698, 236044, 236991, 238546, 239664, 241403, 245292, 246819, 248274, 249074, 249619, 250472, 251345, 252083, 253064, 253682, 488023, 497629, 497931, 497932, 499172, 6K05134, 6K06905, 6L01265, 6L04524, 6L04587, 6M00664, 6M01807, 6M03539, 7A04720, 7B02809, 7B06456, 7C03433, 7D01679, 7D01681, 7D02071, 7D04867, 7E02941, 7F00126, 7F01942, 7F04209, 7G01112, 7H01941, 7H04601, 7H05992, 7K01257, 7K04183, 7L03601, 7L05973, 7M02409, 8A03195, 8A05104, 8B02879, 8C05203, 8D02032, 8D04424, 8E02383, 8F02423
<b>14012183</b>	SONDA RECTAL CH/25 40CM-100UDS	153232, 188205, 201849, 205149, 206029, 208917, 216292, 218381, 220996, 224230, 227698, 230604, 237447, 238051, 240886, 498359, 6K01001, 6L04559, 6M00671, 7A04718, 7A04721, 7B06455, 7C05996, 7E02943, 7F02839, 7G01111, 7J01190, 7L03603, 7M00165, 8A00951, 8A05105, 8C00107, 8E02384, 8F01628
	UNO RECTAL TUBE CH25/40CM(100/600) IT	152990, 153232, 153633, 156032, 158371, 160519, 172958, 176102, 180367, 181715, 185376, 198381, 209174, 218381, 233438, 235224, 248294, 498359, 499311, 7A04718, 7B06455, 7C03432, 7C05996, 7K02646, 7M00165, 8C00107, 8D02033, 8E02384, 8F01628
<b>14012185</b>	UNO RECTAL TUBE CH25/40CM(100/600) INT	8B03348
<b>14014182</b>	UNO RECTAL TUBE CH26/40CM(100/600) INT	157518, 179653, 7H04609, 8E04246
<b>14014183</b>	SONDA RECTAL CH/26 40CM-100UDS	214497, 224263, 238105, 498410, 6K01014, 6L04568, 7A04724
	UNO RECTAL TUBE CH26/40CM(100/600) INT	162846, 166866, 173238, 6K01014, 7A04724, 7A04725
<b>14016181</b>	Rektaaliputki ch28/40cm phtal.ton 50kpl	160948, 180827, 190679, 193197, 200656, 211328, 213659, 215178, 216339, 217985, 225253, 226757, 228669, 230696, 240842, 249463, 250122, 251607, 252585, 6L06121, 6M01804, 7K04187
	REKTALSONDE CH 28/40 CM	157700, 159908, 160948, 184058, 190679, 193197, 200656, 211328, 216339, 217985, 226757, 227713, 228669, 240842, 249463, 250122, 251607, 252585, 6L04564, 6L06121, 6M01804, 7K04187

REF/CODE	DESCRIPTION	LOT/BATCH Number
	UNO RECTAL TUBE CH28/40CM(50/400) INT	157700, 184058, 216339, 217985, 249463, 250122, 252585, 6L04564, 6L06121, 6M01804, 7K04187, 8D02035
<b>14016182</b>	Rektální rourka CH 28/40 cm	7B02811
	UNO RECTAL TUBE CH28/40CM(50/400) AUS	150321, 151665, 152640, 154315, 154592, 155610, 156579, 156580, 157217, 159630, 160489, 163790, 164907, 165517, 166199, 167419, 167933, 171636, 171943, 172866, 177305, 177890, 181639, 184920, 186964, 187807, 189204, 190275, 191168, 191970, 192756, 194546, 195165, 195991, 197545, 198689, 203585, 204643, 205441, 208146, 208471, 210091, 211128, 211490, 212178, 216697, 218971, 219892, 220572, 222118, 222608, 222854, 229089, 230085, 231040, 233385, 234443, 236026, 239646, 241944, 245923, 247471, 250458, 253043, 497455, 497609, 497911, 6K00999, 6K06896, 6K06904, 6L01270, 6L04580, 6L06123, 7A01842, 7B02811, 7D01703, 7F00124, 7H00929, 7J01191, 7K06059, 8A02040, 8C02415, 8D02034, 8F01632
<b>14016183</b>	SONDA RECTAL CH/28 40CM- 50UDS	176492, 178592, 211239, 214469, 215138, 217975, 224231, 230605, 238052, 240887, 250009, 251277, 498360, 6K01012, 6L04560, 6M00667, 7A04722, 7C01105, 7J01202, 8A00949, 8E02385
	UNO RECTAL TUBE CH28/40CM(50/400) IT	153634, 156033, 156254, 156255, 172959, 176103, 180368, 218382, 228120, 233439, 251277, 6K01012, 6M00667, 7A04722, 7C01105, 7J01202, 8A00949, 8E02385
<b>14018181</b>	Rektaaliputki ch30/40cm phtal.ton 50kpl	177657, 184059, 200783, 216340, 234841, 246401, 6M00669, 7A01808, 7B06482, 7F06006, 7L03598, 8A00920
	REKTALSONDE CH 30/40 CM	157720, 184059, 200783, 216340, 234841, 246401, 6L01264, 6M00669, 7B06482, 7H04602, 7L03598, 8A00920
	UNO RECTAL TUBE CH30/40CM(50/400) INT	177657, 184059, 200783, 216340, 234841, 246401, 6L01264, 6M00669, 7A01808, 7B06482, 7F06006, 7H04602, 7L03598, 8A00920, 8D03383, 8F02428
<b>14018182</b>	UNO RECTAL TUBE CH30/40CM(50/400) IT	150346, 151682, 151888, 152489, 154336, 154598, 155233, 156635, 159560, 159659, 162044, 162538, 162801, 163264, 164683, 166213, 166955, 168675, 170757, 172565, 176104, 177942, 179462, 181657, 184931, 186986, 188433, 192006, 192788, 193655, 195205, 197599, 200300, 203618, 205440, 205529, 208510, 209175, 210801, 212220, 213207, 214941, 216019, 218383, 219009, 219918, 220607, 222138, 222607, 224717, 229122, 231085, 234473, 236059, 238576, 240568, 241972, 245306, 245959, 246836, 247503, 249080, 251202, 253090, 253699, 497300, 6K02736, 6K06906, 6L01271, 6L04554, 6L04586, 6M01805, 7A01810, 7A04728, 7B02812, 7C01102, 7C04989, 7D02092, 7E01690, 7F00121, 7F04253, 7G01828, 7H04603, 7K01253, 7L03599, 7L05968, 8A00917, 8C00105, 8D03382, 8D04425
<b>14018183</b>	SONDA RECTAL CH/30 40CM- 50UDS	250010, 498361, 6M00668, 7F00122, 7H00936, 7J01203, 7K02649, 7L05969, 8A05102, 8C02417

REF/CODE	DESCRIPTION	LOT/BATCH Number
	UNO RECTAL TUBE CH30/40CM(50/400) IT	162802, 172960, 498361, 7F00122, 7L03600, 7L05969, 8A05102, 8C02417
<b>14031022</b>	UNO RECTAL TUBE CH28/40CM (100/1000) IT	153191, 155227, 156625, 156626, 164765, 165267, 165442, 167739, 168091, 170403, 171639, 174110, 174368, 176067, 178739, 184654, 185110, 185363, 185837, 187446, 192779, 193257, 194149, 194469, 195189, 197233, 205521, 209164, 216390, 217341, 218595, 219516, 220600, 222470, 223314, 226816, 227770, 228697, 229604, 230118, 232077, 233426, 234562, 237008, 240562, 240957, 245952, 246339, 247500, 251259, 251361, 251564, 497324, 6K03729, 6K06907, 6L01267, 6L06116, 7C01104, 7C06000, 7F00125, 7H00930, 7J02388, 7K02647, 7K04186, 8A04411, 8C02416
<b>14032022</b>	UNO RECTAL TUBE CH30/40CM (100/1000) IT	169690, 171640, 176078, 183158, 194585, 203610, 208505, 210437, 211451, 212863, 219288, 219517, 223315, 226819, 228700, 229608, 251258, 6K00534, 6K02735, 6L01268, 6L04574, 6L06119, 7F00123, 7F04254, 7H03097, 7K01254, 7K04190, 7L05967, 8A00918
<b>14033022</b>	UNO RECTAL TUBE CH32/40CM (75/750) IT	200610, 210378
<b>16002181</b>	ILTFORBINDELSESSLANGE 2TRAGTE HVID 200CM	153983, 160970, 178884, 179880, 183590, 187550, 193324, 201033, 202014, 233060, 300833, 301378, 302958, 304428, 306286, 310344, 310733, 315683, 316382, 317993
	UNO OXY CONN TBE 2 FNL 200CM(50/300) INT	153983, 160970, 178884, 179880, 183590, 187550, 193324, 201033, 202014, 233060, 300833, 301378, 302958, 304428, 306286, 310344, 310733, 315683, 316382, 317993
<b>16003181</b>	ILTFORBINDELSESSLANGE 2TRAGTE HVID 300CM	167230, 179893, 184225, 187515, 192526, 193331, 201984, 208841, 210511, 211374, 221854, 226889, 227840, 302525, 302966, 304095, 305761, 306287, 306548, 308744, 309359, 309573, 310734, 311319, 312011, 313432, 313831, 314221, 316298, 495225
	UNO OXY CONN TBE 2 FNL 300CM(1/50) INT	167230, 179893, 183563, 184225, 187515, 193331, 201984, 211374, 214592, 219689, 221121, 221854, 226889, 230746, 302525, 302966, 304095, 304444, 305761, 306287, 307682, 308744, 309359, 309573, 310734, 311319, 313432, 313831, 314221, 316298, 495225
<b>16013182</b>	UNO CONN TUBE FINGERTIP 210CM(1/50) INT	155077, 155284, 168328, 168638, 168725, 179519, 181311, 184720, 187917, 192799, 196080, 201410, 203007, 206483, 211574, 213787, 214568, 214979, 216410, 216806, 218090, 221099, 224336, 228220, 230722, 300489, 301427, 302033, 304295, 305265, 306144, 306774, 307016, 308273, 309109, 492531, 497478, 499706
<b>16021181</b>	FORBINDELSESSLANGE 210CM, FINGERTIP/TRAG	170528, 211367, 211907, 214585, 215713, 224349, 227833, 229723, 233090, 234881, 236633, 237471, 300636, 304089, 304429, 305266, 307017, 307504, 308745, 309890, 309991, 310953, 312613, 313131, 313832, 314222



REF/CODE	DESCRIPTION	LOT/BATCH Number
	UNO CONN TBE FT FNL CH24/210CM(1/50) INT	151081, 152335, 162526, 166265, 169610, 170528, 181298, 195289, 199328, 203689, 203992, 211367, 211907, 212355, 213799, 224349, 226422, 228303, 230738, 234159, 235307, 236633, 300450, 300636, 303164, 304089, 304429, 305266, 305660, 306870, 307017, 307683, 309574, 309890, 310953, 312613, 313131, 313832, 314222
<b>16021182</b>	UNO CONN TBE FT FNL CH24/210CM(1/50) INT	151300, 152530, 152678, 152679, 156701, 159625, 162061, 162600, 163920, 164796, 165323, 167543, 168598, 170557, 174208, 181338, 186371, 188516, 189243, 191282, 192053, 195238, 196079, 197640, 200368, 203004, 206413, 208558, 210893, 211569, 213255, 214166, 214974, 216796, 221441, 223724, 226384, 227197, 228213, 228833, 229153, 230179, 233566, 234496, 235261, 236094, 300484, 300726, 300938, 301421, 301538, 301795, 302028, 302189, 302690, 303254, 303649, 304151, 304298, 304562, 305267, 305998, 306288, 306601, 307018, 307368, 307505, 308073, 308396, 308586, 308920, 309110, 309360, 309992, 310735, 311320, 312336, 313590, 314696, 315222, 316384, 317778, 497486, 497692, 499191, 499711
<b>16021184</b>	UNO CONN TBE FT FNL CH24/210CM(1/50) INT	152011, 152012, 157543, 173160, 180438, 193779, 200410, 204007, 207669, 211637, 217665, 233673, 237180, 300462, 301598, 304708, 306463, 306891, 308397, 310954, 313313, 318828, 490116
<b>16024182</b>	UNO CONN TBE FT SP FNL CH24/210(1/50)INT	151318, 152342, 155109, 155644, 158724, 159273, 159761, 162477, 165051, 166446, 168607, 168656, 173474, 173985, 179548, 181335, 183827, 184744, 187952, 189257, 191311, 194657, 197665, 199802, 201426, 203696, 205666, 207670, 210929, 212360, 215000, 216107, 217666, 219150, 222189, 224769, 227242, 228315, 231158, 233676, 234537, 237181, 300463, 300797, 301466, 302064, 303173, 303334, 303959, 304620, 305042, 305268, 305661, 306033, 306145, 306289, 306404, 307019, 307313, 307506, 308074, 308348, 308587, 308746, 308921, 309361, 309575, 309993, 310604, 310955, 312012, 312614, 313504, 314223, 314597, 315522, 317349, 318136, 318866, 493893, 494639, 497709, 497710, 499733, 499734
<b>16031182</b>	UNO CONN TBE FNL CH24/350CM (1/50) INT	151295, 159705, 166726, 168639, 172895, 184723, 191285, 194631, 202435, 206420, 214980, 216076, 216808, 225681, 228224, 235267, 300365, 302035, 302824, 306290, 307507, 308075, 309994, 310814, 493789, 497479, 497480, 497481, 499707
<b>16039181</b>	ILTFORBINDELSESSLANGE, 300 CM	153998, 157759, 169475, 181116, 183559, 184251, 193325, 194250, 202015, 204328, 205232, 207373, 208872, 226886, 227834, 230739, 300617, 300649, 301379, 301745, 302959, 305662, 305835, 306291, 307369, 307684, 308398, 308747, 309576, 309995, 310345, 310736, 495065
	UNO COMM TBE CH30/300CM(1/50) INT	202015, 205232, 300649



REF/CODE	DESCRIPTION	LOT/BATCH Number
16039182	UNO COMM TBE CH30/300CM(1/50) INT	151319, 161736, 168608, 171892, 181344, 187953, 195278, 200411, 210145, 211638, 215001, 219151, 224770, 228316, 233677, 300398, 303336, 303960, 304796, 306146, 306524, 307020, 308588, 309111, 309362, 311159, 311572, 497236
16042182	UNO CONN TBE 2 FNL CH24/400CM(1/40) INT	151303, 152531, 155082, 159627, 164329, 166730, 167546, 168647, 169384, 172896, 173973, 179524, 181319, 182984, 185531, 187926, 189249, 190366, 191290, 192802, 193750, 197644, 199281, 202436, 203968, 206414, 210121, 210894, 213256, 214168, 215412, 216799, 219945, 220635, 222955, 228215, 230180, 232144, 234498, 235262, 237811, 300356, 300729, 301540, 302029, 303257, 303643, 304146, 304568, 304770, 305269, 305663, 306292, 306405, 306775, 307021, 308076, 308399, 308589, 308836, 308922, 309363, 309577, 309997, 311531, 493814, 497490, 497491, 497952, 499713
16049181	FORBINDELSESSLANGE 210CM, 2 TRAGTE CH 24	164468, 309364
	UNO CONN TBE 2 FNL CH24/210CM(1/50) INT	151086, 152343, 164468, 166269, 169618, 195279, 199314, 203697, 204008, 206452, 211368, 212361, 214587, 226430, 228317, 300464, 300981, 301382, 302522, 302860, 302962, 303174, 304431, 306167, 306776, 307022, 307508, 308748, 309364, 309998, 310584
16049182	UNO CONN TBE 2 FNL CH24/210CM(1/50) INT	150391, 151980, 156700, 157254, 157905, 158902, 159616, 167122, 168588, 168724, 169381, 170556, 171882, 173972, 174207, 178820, 180425, 181310, 182105, 182981, 184718, 190363, 191281, 193747, 196077, 199276, 200367, 201409, 203666, 206417, 208561, 211573, 213260, 214978, 216804, 219948, 221443, 222166, 225679, 226388, 228219, 230185, 232148, 236097, 300735, 301424, 302573, 302694, 303949, 303979, 305270, 305664, 306147, 306293, 306406, 306777, 307023, 307509, 308077, 308749, 308923, 309999, 310737, 311029, 312253, 312864, 314309, 314565, 315886, 316386, 317460, 497690, 499704
16049183	UNO CONN TBE 2 FNL CH24/210CM (1/50) IT	153656, 156289, 198438, 199308, 209249, 210906, 211591, 224767, 228310, 229173, 235309, 236127, 300395, 301818, 302740, 308078, 308293, 309112, 310000, 310738, 311532, 311728, 312365, 313050, 313857, 314226, 314869, 316029, 316784, 317779, 499329
16049184	UNO CONN TBE 2 FNL CH24/210CM(1/50) INT	152015, 152016, 157544, 157545, 166068, 166069, 185556, 192081, 193780, 200412, 205667, 207671, 226431, 233678, 237182, 301599, 303659, 306464, 308079, 308400, 309440, 309578, 311321, 312337, 312502, 313314, 314416, 315181, 318753, 494973, 494974
16050022	UNO CONN TBE 2 FNL CH24/300CM (1/50) IT	153975, 155436, 158089, 160059, 160984, 161154, 161352, 166814, 176573, 178894, 179878, 181135, 181327, 182985, 192815, 194651, 199293, 203037, 206952, 207364, 208214, 216420, 217370, 232740, 303307, 308828, 309113, 309365, 309579, 309926, 310124, 310739, 310956, 492952, 495141, 497696, 497961
16050181	FORBINDELSESSLANGE 300CM, 2 TRAGTE CH 24	153970, 160971, 181117, 183560, 197380, 208873, 229725, 301494, 302135, 302537, 304760, 305665, 312887, 495147





REF/CODE	DESCRIPTION	LOT/BATCH Number
	UNO CONN TBE 2 FNL CH24/300CM (1/50) INT	183560, 302135, 304760
<b>16050182</b>	UNO CONN TBE 2 FNL CH24/300CM (1/50) IT	151314, 152332, 152561, 152692, 156717, 157351, 158620, 158910, 159644, 162607, 163050, 163930, 165337, 166428, 167559, 168594, 168651, 169388, 170569, 174217, 176574, 180433, 181342, 183815, 184737, 185541, 186382, 187943, 188526, 189255, 191299, 192069, 194653, 195254, 196099, 196800, 199294, 199799, 200394, 202439, 205647, 206435, 207647, 208585, 209260, 212345, 213286, 214997, 216098, 216838, 219959, 220658, 222186, 222990, 224760, 226414, 227226, 228285, 229184, 230201, 232741, 233624, 235289, 236123, 237826, 300438, 300518, 300769, 300968, 302055, 302226, 302427, 302828, 303309, 303655, 304159, 304308, 304601, 305043, 305271, 305772, 306034, 306407, 307024, 307314, 307510, 307685, 308349, 308750, 308924, 309114, 309580, 310001, 310957, 311684, 312338, 313433, 314310, 315148, 316674, 317351, 319267, 497697, 497698, 497699, 499196, 499725
<b>16050183</b>	UNO CONN TBE 2 FNL CH24/300CM (1/50) IT	153657, 156060, 156291, 156292, 156293, 156294, 160546, 164808, 192077, 195269, 200398, 205669, 209250, 221476, 228311, 229174, 235310, 237160, 300396, 301819, 302741, 304165, 308985, 311533, 312366, 312616, 313803, 314227, 314975, 316946, 317705, 319368, 498218
<b>16050184</b>	UNO CONN TBE 2 FNL CH24/300CM (1/50) INT	152017, 157547, 165906, 166447, 184746, 192082, 193782, 196809, 306892, 310346, 311384, 312503, 314417, 316580, 317862, 494975
<b>16056182</b>	UNO CONN TBE 2 FNL CH30/400CM(1/40) INT	151981, 155083, 158614, 159720, 162462, 164797, 165035, 165896, 168648, 173141, 179525, 181320, 185532, 191291, 193751, 200372, 206484, 207633, 209220, 211570, 213816, 214569, 216411, 216800, 218091, 221100, 221442, 224376, 228216, 230181, 233567, 235263, 237101, 300485, 300604, 301541, 302030, 303258, 303644, 304569, 305773, 306168, 306602, 306779, 307686, 308080, 308401, 308751, 309115, 309366, 310002, 310816, 490090, 494955, 497550, 499714
<b>16059182</b>	UNO CONN TBE FNL CH24/210CM(1/50) INT	151296, 152529, 152677, 159706, 163922, 164326, 165032, 167544, 168640, 168727, 170558, 178821, 183808, 186372, 187919, 188517, 200369, 201411, 203667, 206421, 209229, 214170, 214981, 219086, 224741, 228225, 231131, 237814, 300493, 300742, 301431, 302036, 303152, 303248, 304150, 304296, 305272, 306148, 307025, 307389, 308752, 310003, 310347, 311685, 312339, 312865, 313593, 314786, 315150, 316675, 317246, 318515, 497482, 497691, 499708

REF/CODE	DESCRIPTION	LOT/BATCH Number
16064182	UNO CONN TBE FNL CH30/300CM(1/50) INT	151304, 155084, 155931, 159250, 159716, 163651, 164327, 168643, 169383, 173493, 176811, 179521, 181293, 183809, 187921, 189246, 191286, 201412, 203005, 206938, 209224, 211159, 213258, 219082, 222963, 223725, 227199, 230183, 232146, 235265, 300417, 300487, 302031, 302193, 302257, 303153, 303262, 303980, 304300, 304565, 305999, 306149, 306295, 306780, 307026, 307390, 307687, 308753, 309581, 310004, 310348, 497492, 497493, 498881, 499715
16064184	UNO CONN TBE FNL CH30/300CM(1/50) INT	152059, 152060, 173158, 185560, 192108, 193815, 196851, 223000, 226424, 233672, 237167, 300456, 301588, 302743, 303518, 304192, 306465, 306893, 308081, 310349, 495006
16069181	ILTFORBINDELSESSLANGE M/1 TRAGT, 210 CM	202016, 227835, 230773, 238141, 300917, 302960, 308673, 310005, 314297
	UNO CONN TBE FNL CONN CH24/210(1/50) INT	151082, 152336, 152521, 159570, 162528, 166266, 169611, 179545, 190395, 195290, 199330, 202016, 203690, 203994, 205233, 206448, 208839, 212356, 226423, 227835, 228304, 229192, 230773, 233661, 234903, 235308, 235727, 236667, 238141, 300451, 300672, 300917, 301779, 302960, 303165, 304984, 306169, 306781, 307027, 307688, 308673, 308754, 309201, 309582, 310005, 314297, 497703
16069182	UNO CONN TBE FNL CONN CH24/210(1/50) INT	150451, 156704, 162603, 163925, 167547, 200375, 222966, 301432, 305856, 306158, 312254, 315190
16069183	UNO CONN TBE FNL CONN CH24/210CM(1/50)IT	153659, 156299, 162146, 167226, 211594, 220665, 235312, 306525, 312367, 315481, 318993
16069184	UNO CONN TBE FNL CONN CH24/210(1/50) INT	166448, 173163, 193783, 200413, 211639, 217667, 223001, 226432, 233679, 301600, 303491, 308082, 308402, 311386, 313434, 494976
16074182	UNO CONN TBE CARIFLEX CH24/400(1/40) INT	152583, 154360, 155921, 164492, 165663, 166449, 167955, 171613, 173283, 179549, 182119, 186388, 188541, 192083, 195281, 196106, 197666, 200414, 202454, 203045, 203698, 205668, 210146, 212362, 213293, 219152, 219968, 223002, 224771, 225709, 227243, 228859, 232755, 233680, 234538, 236131, 300399, 300798, 301601, 302755, 302861, 303394, 304985, 306180, 306408, 306894, 307391, 308986, 309202, 309441, 309891, 311799
16086242	UNO OXY CONN TBE CH16/150CM(50/300) INT	151302, 159715, 181314, 189245, 194632, 208049, 211158, 216803, 219947, 222164, 222962, 226386, 231129, 232145, 300359, 301542, 303261, 304299, 305274, 305666, 309583, 311161, 313594, 316387, 318867, 497489
16087242	UNO OXY CONN TBE CH16/210CM(50/250) INT	155076, 159623, 168726, 180426, 187918, 191283, 199792, 203972, 209227, 213263, 216807, 222167, 227201, 229158, 237106, 300490, 301429, 302574, 303151, 303648, 304574, 305275, 306150, 306409, 307029, 307511, 308590, 309368, 310006, 310960, 311534, 312618, 313302, 313435, 314051, 315151, 317186, 318754, 497475, 499705
16089182	UNO OXY CONN TBE CH16/150CM(50/300) INT	151998, 152566, 155638, 162065, 165662, 166419, 172172, 172476, 176577, 179534, 183814, 185538, 187014

REF/CODE	DESCRIPTION	LOT/BATCH Number
16092242	UNO OXY CONN TBE CH16/180CM(50/250) INT	150449, 151297, 155079, 159707, 165895, 168641, 181313, 187920, 193748, 197641, 202413, 205599, 210125, 214171, 219949, 227203, 229159, 300311, 300743, 301801, 302037, 302572, 303249, 304575, 305667, 306410, 307030, 308350, 308591, 309116, 309584, 310146, 312619, 313788, 314870, 317461, 318139, 319615, 497483, 497484, 499709
16100181	UNO CONN TUBE CH24/210CM (1/50) INT	152337, 154354, 159571, 162529, 166267, 169612, 185548, 190396, 195291, 203691, 209284, 212357, 300977, 302856, 303166, 304986, 306170, 306782, 307392, 307689, 309369, 310411, 310961, 497704
16100182	UNO CONN TUBE CH24/210CM (1/50) IT	151317, 152571, 152697, 156727, 158622, 159757, 162609, 163815, 163932, 166737, 167563, 168746, 170570, 171058, 173984, 174220, 176580, 178824, 182988, 184738, 187950, 190387, 191306, 192078, 193774, 195270, 196806, 197661, 198439, 202441, 203043, 206430, 209252, 210907, 213281, 214192, 216106, 216849, 219145, 220353, 221477, 223755, 224768, 227238, 228312, 231157, 232159, 233663, 234156, 236129, 237161, 237823, 300667, 300914, 301464, 301582, 301822, 302437, 303302, 303385, 304310, 304403, 304576, 304783, 305276, 306035, 307031, 307393, 307690, 308755, 308925, 309203, 309930, 310007, 310351, 310585, 310962, 497706, 498597, 499732
16100183	UNO CONN TUBE CH24/210CM (1/50) IT	155293, 156846, 159574, 161967, 164809, 171655, 174221, 176581, 182989, 188135, 191307, 193776, 196100, 203044, 205670, 206947, 209253, 216850, 219956, 222188, 223756, 225708, 228313, 236130, 302592, 304311, 499064
16100184	UNO CONN TUBE CH24/210CM (1/50) INT	152019, 152020, 157548, 166071, 166072, 173164, 196810, 200415, 204009, 211640, 223003, 226433, 233091, 233681, 237183, 301602, 302066, 303513, 304710, 305853, 306895, 308083, 308403, 309442, 310352, 310740, 494977, 494978
16102181	Imuletku OP-Flex N/N ch 24/350cm 50 kpl	201973, 213047, 229724, 302620, 303004, 303550, 306181, 308756, 309117, 310125, 310412, 495145
	OP-Flex, 350 cm forbindelsesslange	190804, 201034, 201973, 213047, 229724, 302620, 302783, 303004, 303550, 305774, 306181, 306466, 307394, 308404, 308756, 309370, 310125, 310412
	UNO CONN TUBE CH24/350CM (1/50) INT	159572, 190397, 190804, 201034, 201973, 213047, 229724, 233662, 300978, 302620, 302783, 302816, 302857, 303004, 303167, 303550, 305277, 306181, 306296, 306466, 307394, 309117, 309370, 309585, 310125, 310412, 495145, 497705
16102182	TUBULURE RACCORDEMENT OP- FLEX 3.5M	150530, 236656, 300655, 306603, 307032, 307395, 308757, 310008, 310611



REF/CODE	DESCRIPTION	LOT/BATCH Number
	UNO CONN TUBE CH24/350CM (1/50) INT	150530, 151306, 152323, 152532, 152680, 155087, 155286, 155628, 156702, 158715, 158904, 159719, 162460, 162601, 163043, 163923, 165034, 165327, 166253, 166727, 167008, 167545, 168599, 169601, 170559, 170766, 174209, 176566, 179522, 181084, 181294, 181885, 181921, 184725, 185143, 185527, 187922, 189247, 189888, 191287, 191740, 192054, 193355, 194633, 195009, 195239, 196788, 198427, 199278, 199793, 200370, 202431, 203469, 203967, 204339, 205207, 205597, 206412, 207388, 208226, 208557, 208855, 209218, 210892, 212313, 212950, 213254, 214973, 216073, 216795, 219078, 222162, 224740, 227196, 228212, 230178, 231126, 233565, 234148, 234495, 236656, 237478, 237807, 300483, 300655, 300725, 300935, 301419, 301794, 302187, 303250, 304141, 304286, 304566, 304768, 305278, 305668, 306000, 306411, 306603, 307032, 307395, 308592, 308757, 308926, 309118, 309371, 309540, 309586, 309931, 310008, 310611, 497495, 497693, 498583, 499193, 499716
16102183	UNO CONN TUBE CH24/350CM (1/50) IT	150629, 153658, 156061, 156296, 156297, 156298, 160548, 160549, 160550, 187023, 188537, 195271, 211593, 237162, 300433, 302742, 492915
16102184	UNO CONN TUBE CH24/350CM (1/50) INT	152022, 152023, 157549, 165908, 166450, 173166, 185557, 193785, 200416, 204010, 207672, 211641, 217668, 223004, 226434, 233682, 237184, 300465, 301603, 302067, 303339, 306896, 308084, 308405, 308927, 309443, 494979
16151022	UNO OXY CONN TBE CH14/200CM(50/300) INT	202419, 304694, 316076, 497813
16172182	UNO CONN TBE FT 2FNL CH24/300CM(1/50)INT	157771, 160989, 164484, 169465, 170508, 170509, 172045, 174499, 175445, 178888, 178918, 179176, 179897, 181092, 181898, 181924, 182783, 183605, 184262, 184274, 190835, 192531, 193373, 194284, 195651, 197416, 198025, 198999, 199563, 199721, 201069, 204348, 207015, 208219, 208387, 208880, 209916, 210514, 211099, 211908, 212956, 213836, 214621, 217390, 218647, 219692, 220372, 221150, 221857, 222550, 223480, 224402, 226925, 229280, 229756, 230759, 231636, 232469, 233687, 234161, 234916, 235730, 236671, 236770, 300598, 300679, 300805, 301496, 301911, 302299, 303010, 303694, 303782, 304463, 305279, 306526, 306604, 306783, 307033, 307315, 307512, 308085, 308274, 308351, 308593, 309372, 309587, 310009, 491520, 495054, 495152
16173182	UNO CONN TBE FT 2FNL CH24/500CM(1/40)INT	157772, 163690, 164474, 172046, 194263, 201070, 207399, 208220, 213837, 219693, 219837, 226926, 233074, 234162, 300923, 302156, 302300, 302985, 304464, 305280, 306182, 307396, 309373
19005182	UNO OXY CATH GREEN CH8/40CM(100/1000)INT	152316, 172635, 178761, 194599, 211153, 228155, 247513, 6M00339, 7D04805, 7H04072, 8B04615
19006182	UNO OXY CATH GREEN CH10/40(100/1000) IT	151889, 152494, 154334, 157309, 159901, 162798, 165652, 166212, 172564, 174979, 180325, 181581, 184667, 192004, 201385, 205527, 208508, 239701, 497299, 497937, 7J05966, 8C05605
19006183	UNO OXY CATH GREEN CH10/40(100/1000) IT	153631, 156250

REF/CODE	DESCRIPTION	LOT/BATCH Number
	UNO OXY CATH GREEN CH10/40CM(1X100PK) IT	218380, 250498, 498357, 6L06019, 7D01371, 7G05467, 8B04926
<b>19007182</b>	UNO OXY CATH GREEN CH12/40(100/1000) INT	172636, 228156, 7D04808
<b>19008182</b>	UNO OXY CATH GREEN CH14/40(100/1000) INT	154341, 206928, 212241
<b>19017182</b>	UNO OXY GREEN FOAM CH10/40CM(100/600)INT	150353, 152505, 152665, 154342, 155242, 155898, 156035, 156647, 157242, 158893, 159665, 159917, 162058, 162558, 162836, 162959, 163886, 163887, 165297, 165825, 167771, 169704, 170466, 172469, 173391, 174846, 176448, 176557, 177728, 177984, 178762, 179471, 180998, 181494, 182690, 183757, 184104, 184686, 185474, 186356, 187003, 188226, 188904, 190732, 191259, 193680, 194185, 194600, 196755, 198962, 199236, 200984, 202411, 202781, 204302, 205551, 208205, 208528, 209007, 211071, 211843, 213242, 213738, 215140, 216037, 217331, 219033, 220342, 222913, 223342, 225633, 226768, 228157, 229645, 230145, 231098, 233468, 234125, 234482, 234861, 237427, 239119, 240582, 241987, 245969, 246864, 247128, 249470, 249653, 250060, 251372, 252113, 252609, 497218, 497941, 499182, 6K02048, 6K04860, 6K04861, 6L01078, 6L02707, 6L04906, 6M00323, 7A01930, 7B03375, 7B06788, 7C05475, 7E00344, 7E03853, 7F02103, 7G00535, 7H02463, 7J03839, 7L00934, 7M02335, 8A02881, 8B01985, 8D00284, 8D03435, 8E01100, 8F00409
<b>19019181</b>	Happikatetri ch 14 superlonkärki 100 kpl	163588, 200617, 222490
	UNO OXY GREEN FOAM CH14/40CM(100/600)INT	200617, 222490
<b>19019182</b>	UNO OXY GREEN FOAM CH14/40CM(100/600) IT	150350, 151887, 152488, 152656, 154335, 155232, 155894, 156634, 157232, 162043, 162945, 164682, 165279, 166954, 168521, 170743, 172468, 176100, 179461, 180326, 181582, 183741, 184668, 185375, 191231, 193654, 196725, 199201, 202366, 205528, 208509, 213205, 216018, 219008, 222891, 225617, 228119, 230129, 231084, 234470, 236058, 238575, 240566, 245957, 246833, 249634, 251363, 497185, 6K02075, 6K04863, 6L03927, 6M01173, 7A05387, 7B04091, 7C02886, 7C05249, 7E00343, 7E04878, 7G00965, 7H02756, 7L01742, 8A02878, 8B04299, 8C02278, 8D02623, 8F00081
<b>19020182</b>	UNO OXY GREEN CH10/40CM (100/1000) INT	150354, 155899, 184687
<b>19024182</b>	UNO OXY ADJ FOAM CH10/40CM(100/600)INT	154344, 155900, 163807, 172592, 181598
<b>19030182</b>	SONDE OXYGENE AVEC COMPRESSE 40CM CH10	158961, 161621, 162132, 234755, 6K04073, 7A02135, 7C02888, 7F00831, 7G05466, 7L01411, 7L06810, 8B00814, 8B04925
	UNO OXY ADJ FOAM CH10/40CM(100/600)INT	8B04925, 158809, 190983, 8B00814, 8E02430

REF/CODE	DESCRIPTION	LOT/BATCH Number
19034182	SONDE OXYGENE SANS COMPRESSE 40CM CH08	161639, 193154, 196510, 211292, 220983, 223369, 225274, 227756, 234787, 250106, 6K06844, 6L01664, 7B05103, 8A04758
	UNO OXY CATH BLUE CH8/40CM(100/1000) AUS	161639, 225274, 7B05103, 8A04758
19034242	UNO OXY CATH BLUE CH8/40CM(100/1000) INT	168958, 214949
19035182	SONDE OXYGENE SANS COMPRESSE 40CM CH10	158984, 165437, 167477, 167814, 169136, 170463, 174323, 176402, 176463, 178623, 180963, 181851, 185076, 188171, 188873, 198901, 202768, 205133, 206746, 208966, 209144, 211293, 214507, 218614, 219576, 221048, 224276, 225275, 226797, 227757, 230592, 233021, 234788, 236503, 239152, 240946, 247971, 249989, 252541, 498299, 6L01464, 6L06594, 7A02134, 7D03814, 7H06064, 7M00379, 8A01457
	UNO OXY CATH BLUE CH10/40CM(100/1000)AUS	158984, 167814, 8A01457, 8E00901
19035242	UNO OXY CATH BLUE CH10/40CM(100/1000)INT	8C01512
19035908	OXY CATH BLU CH10/40CMTHRMFSHR(100/1000)	193299
19037182	SONDE OXYGENE SANS COMPRESSE 40CM CH12	150085, 167516, 167815, 169102, 176536, 178624, 181852, 200727, 201889, 202735, 207257, 208967, 209145, 211294, 218029, 224277, 227758, 229667, 232411, 234789, 240866, 247909, 249990, 498300, 6L02344, 6M00333, 7A02137, 7H03884, 8C04561
	UNO OXY CATH CH12/40CM (100/1000) AUS	167516, 7H03884, 8C04561
19038182	SONDE OXYGENE SANS COMPRESSE 40CM CH14	150005, 167108, 179771, 209146, 238031, 6L05092
	UNO OXY CATH CH14/40CM (100/1000) AUS	6L05092, 8A05947
19040052	UNO OXY BLUE FOAM CH10/40CM(100/600) INT	151178, 154979, 159679, 163311, 165017, 166710, 168701, 172833, 179480, 181531, 185331, 187769, 189186, 191143, 195137, 195973, 199123, 200222, 203561, 206331, 209118, 210756, 213136, 214887, 215964, 218940, 221357, 222831, 226276, 227088, 228012, 232025, 233351, 238512, 239592, 245902, 247448, 250413, 252030, 497012, 497013, 499520, 6K01705, 6K01706, 6L00243, 6L01459, 6L06016, 7A01929, 7A04285, 7A07021, 7B05241, 7C00989, 7C04315, 7D01369, 7D03813, 7E01096, 7E04074, 7F00829, 7H00707, 7J03838, 7K01811, 7K03180, 7K06464, 7L04766, 8A01332, 8B03350, 8B05301, 8C04371, 8D03638, 8E01903, 8E03759, 8F03546

REF/CODE	DESCRIPTION	LOT/BATCH Number
19041052	UNO OXY BLUE FOAM CH12/40CM(100/600) INT	151193, 159670, 160261, 168491, 168925, 172850, 174868, 181545, 184614, 187785, 190260, 193599, 195150, 197534, 200231, 202393, 203572, 206339, 207536, 210081, 213144, 214899, 215978, 222108, 223631, 225573, 227098, 228035, 231025, 235191, 236008, 237730, 239612, 241380, 248247, 250429, 253008, 497049, 499557, 6K01707, 6K06673, 6L01077, 6L06017, 6M02973, 7A05311, 7B02834, 7B05100, 7C03161, 7D02295, 7E02029, 7E05384, 7F02821, 7H05577, 7J05964, 7K05782, 7L01741, 7L03868, 8A01460, 8A03138, 8A04354, 8B03605, 8C01670, 8C04555, 8D03918, 8E03491, 8F01694
19042052	UNO OXY BLUE FOAM CH14/40CM(100/600) INT	155004, 159619, 168926, 172852, 181634, 185341, 189198, 191160, 197539, 199754, 202919, 206893, 210760, 213146, 215983, 216670, 222111, 224683, 227100, 229073, 232031, 236012, 238519, 241382, 247458, 250431, 253010, 497053, 499560, 6K04859, 6K07201, 6L01659, 6M00326, 7A05352, 7B05099, 7C01671, 7D02950, 7E04877, 7F03913, 7H04186, 7K03179, 7L06809, 8A02879, 8B04300, 8C02277, 8E00900
19046022	UNO OXY CATH GREEN CH08/40(100/1000) INT	219771, 231058, 245930
19047022	UNO OXY CATH GREEN CH10/40(100/1000) INT	153100, 154130, 155222, 156817, 161093, 164758, 166390, 166854, 172883, 178733, 181570, 184646, 190300, 192768, 199771, 201376, 206906, 212197, 216710, 222875, 224781, 227136, 235206, 236992, 245295, 251346, 497933, 6K04970, 6K04972, 6M00623, 6M02981, 7B05102, 7C00498, 7C01480, 7C06393, 7D01647, 7F00830, 7F03121, 7J00284, 7J01365, 7J05965, 7K03589, 7L04767, 8A01458, 8C00396, 8C01896, 8E01102, 8E02431, 8F00759
19049022	UNO OXY CATH GREEN CH14/40CM (1X100) IT	240261, 6K00314, 7A06711
	UNO OXY CATH GREEN CH14/40CM(100/1000)IT	6K00314, 7J03169
19055022	UNO OXY GRN FOAM CH10/40CM(100/1000) INT	151500, 154135, 158076, 158512, 160664, 167765, 171647, 172585, 174131, 176119, 178756, 179463, 181669, 182954, 185452, 188480, 194606, 197616, 199234, 205540, 206369, 214118, 215401, 216024, 218393, 219023, 221418, 223689, 232088, 233451, 238588, 241976, 246844, 250510, 253098, 497301, 6K04950, 6K04951, 6M00311, 6M02650, 7A02128, 7B01022, 7F01895, 7F06051, 7G02616, 7J02350, 7L03565, 7M02930, 8A06298, 8D02801
19056022	UNO OXY GRN FOAM CH12/40CM(100/1000) INT	154136, 157310, 167766, 174836, 179464, 185453, 190334, 202377, 215402, 218394, 220612, 221419, 232089, 237024, 237756, 240574, 6K02077, 6K05744, 7D00476, 7E05383, 7K03036, 8A05944, 8C01893, 8D02500, 8E03492

REF/CODE	DESCRIPTION	LOT/BATCH Number
19057022	UNO OXY GREEN FOAM CH14/40CM(100/1000)IT	154134, 156636, 160661, 162440, 167743, 179458, 181656, 188424, 190315, 191223, 195202, 202365, 206360, 207582, 217638, 218389, 220609, 225621, 226211, 235229, 237018, 241434, 250507, 6K05478, 6L03368, 6M00325, 7D04824, 7F00082, 7F04200, 7J01854, 7L02113, 8A03736, 8C03623, 8D02622
19068022	UNO OXY ADJ FOAM CH12/40CM(100/1000) INT	219772
19069022	UNO OXY ADJ FOAM CH14/40CM(100/1000) INT	6M03680
19074052	UNO OXY BLUE FOAM CH10/40CM(100/600) INT	155060, 164193, 168897, 169336, 181721, 187883, 194608, 199785, 206370, 213224, 216743, 227157, 228130, 237757, 247505, 253101, 499645, 6K03598, 6M00622, 7B03376, 7C06394, 7F00066, 7H03169, 7K03185, 7M02363, 8B01986, 8D02625, 8F00757
29001182	UNO NASAL OXY CANN STR(50/300) IT	150458, 152578, 155920, 156063, 156685, 158915, 159125, 162068, 162532, 162871, 163079, 163934, 165341, 167014, 168606, 169369, 169389, 169609, 179567, 180459, 181334, 184743, 187026, 190386, 191305, 194656, 196807, 199313, 202448, 205617, 208576, 230191, 233595, 303657, 304319, 304996, 306426, 306903, 308113, 308604, 308944, 310363, 311585, 312031, 312473, 314299, 315204, 316106, 317867, 318764, 499197
29002182	LUNETTE OXY STE CONNECT M-F TUB2M ADULTE	211389, 214626, 218138, 220375, 221136, 224375, 226146, 226935, 227854, 301768, 302535, 303404, 304124, 304475, 305697, 307070, 307722, 308771, 309209, 309611, 310037, 312651, 313021, 314988, 315528, 316403
29003242	UNO NASAL OXY CANN STR FOAM (50/300) INT	187948, 214191, 223753, 304610, 310038, 314989
29006242	UNO NASAL OXY CANN STR (50/300) INT	151301, 155081, 158613, 160147, 168729, 180458, 185530, 189248, 199796, 203663, 207632, 214167, 214975, 216798, 227198, 232143, 237100, 237810, 300728, 303148, 303256, 304769, 305336, 306537, 306808, 308605, 309612, 311030, 312652, 313445, 314058, 315372, 316404, 318142, 497488, 499712
29008181	HAPPIVIIKSET, STERIILI, VIHREÄT 50 KPL	184221, 211361, 221109, 224391, 301374, 304118, 304415, 305698
	ILTBRILLE LIGE CH10/200CM STER (GROEN)	157782, 161438, 167241, 170546, 181114, 183575, 184221, 197376, 201960, 204323, 205225, 211361, 214580, 218095, 219655, 221109, 224391, 227824, 230731, 300830, 301374, 301740, 302951, 303767, 304415, 305698
	UNO NASAL OXY CANN STR FNL (50/300) INT	157782, 167241, 181114, 184221, 197376, 204323, 218095, 221109, 224391, 226134, 301374, 304118, 304415, 305698, 307723, 309613
29008182	UNO NASAL OXY CANN STR FNL (50/300) INT	150398, 150453, 154350, 162063, 163841, 174228, 176583, 178028, 179565, 181297, 182111, 185537, 186374













REF/CODE	DESCRIPTION	LOT/BATCH Number
29010242	UNO NASAL OXY CANN STR FOAM (50/300) INT	151308, 155092, 157906, 159723, 160728, 166728, 168644, 178826, 181317, 185528, 187924, 191289, 193749, 199279, 202432, 206936, 209219, 211568, 222954, 223723, 226385, 228214, 231127, 234497, 237808, 300937, 301796, 303253, 303463, 304567, 305337, 305775, 306008, 306621, 307071, 308114, 308772, 309614, 310039, 499719
29011242	UNO NASAL OXY CANN STR FNL (50/300) INT	166729, 168645, 171883, 181318, 185529, 187925, 190365, 194634, 196789, 199794, 202433, 206937, 211157, 216074, 216797, 219944, 225676, 231128, 232142, 237809, 300727, 301420, 302401, 303640, 305338, 305776, 306622, 306809, 307534, 308773, 309390, 315161, 499194, 499720
29014181	ILTBRILLE BUET CH10/200CM, STER (GROEN)	150543, 152420, 154013, 157783, 161439, 170547, 172063, 174534, 175473, 179183, 179922, 181115, 183589, 184365, 187544, 189895, 193315, 202009, 205226, 208837, 211408, 214614, 224345, 227825, 229743, 230768, 300611, 300903, 301777, 302553, 302952, 303768, 304119, 304416, 310749
	UNO NASAL OXY CANN CRV GREEN(50/300) INT	150543, 152420, 154013, 157783, 161439, 170547, 174534, 179183, 179922, 181115, 183589, 184365, 189895, 202009, 205226, 211408, 224345, 227825, 229743, 230768, 301777, 304119
29015181	ILTBRILLE BUET CH10/210CM, STER (KLAR)	226891
29015182	UNO NASAL OXY CANN CRV CLR (50/300) INT	167020, 169606, 191329, 196845, 202425, 208589, 216104, 225705, 230210, 234535, 302236, 302852, 304191, 304350, 304997, 305844, 306427, 306904, 307375, 308115, 308774, 308945, 309615, 310365, 312653, 314300, 314427, 315205, 317254, 317892, 318766
29017181	ILTBRILLE LIGE CH10/200CM STER (KLAR)	226135
29018182	UNO NASAL OXY CANN CRV GREEN(50/300) INT	150541, 152569, 152695, 156724, 161437, 163928, 167557, 170542, 170567, 172062, 174229, 181341, 187502, 188285, 190796, 196096, 227208, 233053, 233589, 234153, 234896, 236627, 237483, 300661, 301366, 301439, 301487, 302127, 302296, 304153, 305339, 309456, 310040, 310366, 310559, 310599, 310750, 310984, 311172, 312486, 312654, 312871, 313022, 313306, 313838, 314059, 314428, 314707, 314789, 315032, 315276, 315757, 316284, 317108, 318280, 318439
29022022	UNO NOC STR NOSE PC FNL (50/300) INT	156694, 159937, 164806, 167503, 176266, 192079, 192816, 201422, 205655, 210919, 212350, 222992, 226417, 227235, 230208, 235304, 300771, 301577, 302436, 302739, 309146
33001181	Oeresug, buk	479440
	UNO PRCN HNDL SGL BEND 10CM(100/1000)NEU	479440, 481829
33001182	UNO PRCN HNDL SGL BEND 10CM(100/1000)EUR	152634, 163866, 170754, 181691, 195974, 202913, 216665, 219879, 227090, 230063, 234416, 236000, 246767, 249583, 479250, 481802, 485393, 486189, 492546, 497121, 498074, 6L06114, 6M03532, 6M03533, 7C02295, 7D04863, 7F00138, 7H01944, 7K04189, 8A03851, 8B04375, 8D04423

REF/CODE	DESCRIPTION	LOT/BATCH Number
33002182	UNO PRCN HNDL DBL BEND 10CM(100/1000)EUR	161719, 162580, 190269, 200235, 224673, 229060, 246753, 478852, 488004, 498073, 6L04518, 7F04256, 7L03605, 8A03855, 8C00111, 8C03634
33002183	UNO PRCN HNDL DBL BEND 10CM(100/1000)SEU	156251, 167748, 237013, 8F02433
	UNO PRECISION HNDL DBL BEND(100/1000) IT	237013, 496293, 6L01221, 7M00167, 8C02421
33004182	CANULE ORL COUDEE CH06	167139, 208948, 7C04976
	UNO PRCN HNDL SNG BEND 13CM(100/1000)EUR	7C04976
33005182	CANULE ORL BICOUDEE CH06	184113, 480633, 6M03538
	UNO PRCN HNDL DBL BEND 13CM(100/1000)EUR	480633, 8F02432
33094182	UNO FLTFLOW HNDL WO TBE CH24 (20/40) INT	162618, 163729, 164993, 184302, 185169, 188326, 190851, 192564, 194308, 195041, 195669, 197438, 198041, 205261, 206479, 208250, 209034, 210533, 211111, 213862, 214207, 215209, 216459, 217401, 220396, 230228, 233103, 300574, 300579, 300603, 304473, 304659, 305360, 306820, 307085, 307541, 307869, 312263, 312266, 312827, 316170, 496475, 498188
33094183	UNO FLTFLOW HNDL WO TBE CH24 (20/40) INT	164836, 178039, 187040, 200474, 312267, 497711, 497712, 498124
33094822	UNO DELTAFLO YANK WO TBE (20/40) INT	150399, 158617, 172878, 172889, 172892, 179580, 184765, 191355, 194691, 199823, 206469, 215226, 215227, 300845, 300950, 303236, 303237, 303238, 308647, 310906, 311342, 312894, 490151, 490152
33095181	OP-FLEX FILTER FLOW HANDLE (20/40kpl)	159515, 174636, 309214
	OP-FLEX HAANDSTYKKE M/FILTERFLOW	155445, 170648, 175494, 181175, 201122, 202036, 301770, 302992, 304126, 308832, 309214, 312261, 312264, 312268, 312833, 495352, 498189
	UNO OP-FLEX FILTERFLOW HNDL(20/40) INT	154521, 155445, 159515, 169312, 170648, 174636, 175494, 179237, 181175, 201122, 202036, 202308, 302992, 304126, 309214, 312261, 312268, 312833, 312834, 495352, 498189
33095182	CANULE OP-FLEX FILTER FLOW SS TUB CH24	170637, 170795, 174625, 175480, 203462, 206058, 207369, 212930, 213831, 300257, 302151, 308129, 312265, 312269, 498198
	UNO OP-FLEX FILTERFLOW HNDL(20/40) IT	150578, 151320, 154847, 154919, 155110, 157373, 159758, 163350, 166910, 167566, 168652, 170637, 170795, 174625, 175480, 179941, 202309, 207343, 207369, 208393, 208577, 211170, 211405, 212930, 216085, 218120, 219110, 223739, 229175, 300257, 300846, 302151, 302947, 303985, 307542, 307732, 308129, 308678, 308782, 312265, 312269, 497748, 498104, 498198

REF/CODE	DESCRIPTION	LOT/BATCH Number
40071181	SIMPAVAC, SUGEFLASKE 200 ML	160979, 164473, 164978, 178917, 179888, 181128, 183604, 184238, 184272, 197415, 198998, 201118, 201993, 207381, 208850, 211385, 214599, 219671, 221134, 221843, 224372, 226137, 226897, 230756, 300632, 300910, 301762, 302532, 302983, 304462, 305408, 305736, 306345, 308806, 310785, 311011, 311366, 312058, 314071, 314270, 314768, 316197, 316768, 316985, 495226
	UNO SIMPAVAC 200ML (1/50) INT	160979, 164473, 201118, 300910, 311366, 312058
40071182	UNO SIMPAVAC 200ML (1/50) INT	302105, 309562, 309664, 311210, 313499, 314130, 316198, 319555
40159181	TRAPPEKONNEKTOR MED HUN-LUER	150527, 153987, 155437, 157766, 159410, 160975, 164470, 164977, 169477, 169478, 170526, 170777, 181124, 183579, 184234, 194261, 198994, 201059, 201991, 211378, 211898, 213807, 214596, 230752, 304455, 305424, 305744, 307761
	UNO STEP CONN W/FEMALE LUER(35/140) INT	150527
42001181	KATETERSTUDS M/FLASKEKANYLE	160988, 174526, 178891, 179895, 181127, 182771, 183601, 184260, 196523, 197388, 198997, 201066, 202028, 205244, 208847, 211392, 211899, 214598, 218106, 219669, 221131, 224369, 227849, 230780, 234914, 236650, 300631, 301761, 302531, 302982, 304105, 304461, 304762, 305432, 305751, 309684, 310106, 311199, 311373, 312724, 313036, 313477, 314083, 314283, 315021, 315404, 316436, 318589
	UNO CONN FOR BOTTLE NDL(25/250) INT	160988, 178891, 183601, 184260, 201066, 211392, 218106, 227849, 234914, 236650, 304461, 304762, 305432, 312724, 315404, 318589
84011022	UNO FORCEPS BLUE 125MM (100/1200) IT	151542, 154165, 155440, 157353, 160670, 161548, 165044, 168744, 169488, 170635, 170794, 171055, 174708, 181414, 184298, 185622, 187575, 190419, 194693, 197688, 197706, 202460, 205693, 208602, 211199, 213851, 214203, 216455, 217680, 219724, 221500, 225723, 227888, 228338, 228820, 229765, 232765, 235334, 301810, 302432, 304695, 305780, 306458, 307128, 308183, 308977, 309190, 309565, 310113, 311569, 311701, 312077, 312387, 313480, 314086, 315555, 316574, 316575, 316576, 317128, 319471, 495142, 497700
84011181	ATULA SININEN 100 kpl	156450, 159524, 172084, 174632, 175483, 181165, 184303, 192568, 194304, 197445, 206079, 209027, 211443, 215216, 231643, 300643, 301769, 302536, 302871, 303039, 304125, 307141, 307775, 308667, 308823, 309237, 309692, 310114, 315763, 316100

REF/CODE	DESCRIPTION	LOT/BATCH Number
	PINCET, BLAA, STERIL	154040, 156450, 157795, 159524, 162011, 164994, 169490, 170647, 170800, 172084, 174632, 175483, 179236, 181165, 181935, 182806, 184303, 188933, 192568, 194304, 197445, 198047, 201155, 202070, 206079, 209027, 210531, 211443, 215216, 225329, 229770, 231643, 234935, 300643, 300927, 301769, 302536, 302800, 303016, 303039, 304125, 305039, 305439, 305757, 307141, 307361, 307775, 308667, 308823, 309237, 309435, 309692, 309924, 309962, 310114, 312747, 314287, 315405, 315763, 316100, 316212, 316997, 318170
	UNO FORCEPS BLUE (100/1200) NEU	151101, 152358, 169624, 172190, 175483, 181165, 181378, 185628, 190423, 197445, 199028, 199372, 203712, 209318, 211443, 212418, 216868, 219176, 226452, 228341, 231643, 233743, 235338, 300643, 302536, 302871, 303039, 303180, 304125, 305757, 307141, 308823, 309237, 309692, 309924, 312747, 315405, 315763, 316100, 318170, 497715, 498137
<b>84011182</b>	Pinzeta modrá	300936, 302287, 304563, 307129, 308166, 308633, 308978, 316050
	UNO FORCEPS BLUE (100/1200) EUR	150554, 150555, 151292, 152525, 152674, 152835, 159668, 160138, 162599, 163341, 168288, 176327, 181409, 182217, 183859, 184764, 187974, 188584, 189326, 190410, 191349, 193828, 196127, 197680, 198461, 199820, 200458, 202465, 203704, 204019, 204662, 207677, 210944, 211433, 212401, 213848, 214198, 215008, 215429, 216862, 217675, 218427, 221486, 224785, 225715, 226149, 226441, 227248, 227886, 228328, 229771, 230218, 231165, 231247, 231638, 232757, 233104, 233713, 234545, 234928, 235731, 300305, 300325, 300349, 300475, 300936, 301499, 301537, 301771, 302123, 302178, 302188, 302287, 302539, 302994, 303252, 303413, 303639, 303678, 304110, 304143, 304374, 304563, 304767, 305017, 305440, 305758, 306058, 306359, 306459, 306552, 306887, 307129, 307578, 307776, 308166, 308189, 308633, 308978, 309462, 311570, 312206, 313614, 316049, 316050, 316051, 317017, 317261, 318273, 495805, 497472, 497688, 498582, 499190
<b>84011183</b>	UNO FORCEPS BLUE (100/1200) IT	178037, 183867, 192119, 204022, 206474, 218444, 227256, 300308, 300337, 302747, 303389, 306222, 308167, 312388, 312729, 314288, 316438, 316998, 318662
	UNO FORCEPS BLUE (100/1200) SEU	156308, 156309
<b>84011184</b>	UNO FORCEPS BLUE (100/1200) PL	152063, 157585, 157586, 165918, 166473, 173196, 188598, 192126, 193842, 196878, 199825, 211658, 217683, 223017, 226450, 233732, 300332, 301570, 302431, 302730, 303500, 304717, 306888, 308168, 308450, 309925, 310399, 311410, 312102, 312518, 313347, 317884, 318779, 495004

Attachment 2: Example of Package Labelling

For the Affected product	
The shipper label	<div style="border: 1px solid black; padding: 5px;"> <p><b>REF 29001182</b> ← <b>Product code</b></p> <p>GB Twin nasal oxygen cannula  D Sauerstoff-Brille  F Lunette a oxygene  NL Zuurstofbril  I Cannula per ossigenoterapia</p> <p>Made in Belarus  <b>CH 10 210CM</b></p> <p style="text-align: center; font-size: 24px;"><b>6 X 50</b></p> <div style="display: flex; justify-content: space-between;"> <div style="text-align: center;">  </div> <div style="text-align: center;">  </div> <div style="font-size: 8px;"> ConvaTec Limited  First Avenue  Deeside Industrial Park  Deeside, Flintshire, CH5 2NU, UK </div> </div> <hr/> <div style="display: flex; justify-content: space-between;"> <div style="text-align: center;">  <b>2012-04</b>   <b>2017-03</b> </div> <div style="text-align: center;"> <b>LOT 630904</b> ← <b>Lot/Batch Number</b> </div> <div style="font-size: 8px;"> EAN No. 05705243273444 </div> </div> <hr/> <div style="text-align: center;">  <p style="font-size: 8px;">(01)05705243273444(17)170331(10)630904 v.05</p> </div> </div>
The market unit label	<div style="border: 1px solid black; padding: 5px;"> <p><b>50 REF 29001182</b> ← <b>Product Code</b> <span style="float: right;"><b>CH 10 210CM</b></span></p> <p>GB Twin nasal oxygen cannula  D Sauerstoff-Brille  F Lunette a oxygene  NL Zuurstofbril  I Cannula per ossigenoterapia</p> <div style="display: flex; justify-content: space-between;"> <div style="font-size: 8px;"> EAN No. 05705243273437 </div> <div style="text-align: center;">  </div> <div style="text-align: center;"> <b>LOT 630904</b> ← <b>Lot/Batch Number</b> </div> </div> <div style="display: flex; justify-content: space-between; margin-top: 10px;"> <div style="text-align: center;">  <b>2012-04</b>   <b>2017-03</b> </div> <div style="text-align: center;">  </div> </div> <hr/> <div style="text-align: center;">  <p style="font-size: 8px;">(01)05705243273437(17)170331(10)630904 v.05</p> </div> </div>