

URGENT FIELD SAFETY NOTICE

Commercial Name: Convex Two-Piece Skin Barriers i.e. Natura / Surfit / Combihesive Wafers
(see Attachment 3 for commercial name and full details)

Issue Date: 20 September 2019

X Original Notice ___ Revised Notice

Revision No.: Rev. 1 (for EU)

Reference/REF No: See Attachment 3: Product List

FSCA Ref.: 2019-006MA

Type of Action: Customer Notification

FSCA Helpdesk Email: Refer to the “ ‘Regional ConvaTec Contact’ for other support:” section

Please note that this action only applies to specific product codes and lot numbers of the above-mentioned products.

DESCRIPTION OF THE PROBLEM	<p>ConvaTec has voluntarily initiated a Field Safety Corrective Action in the form of a Customer Notification for specific product codes and lots of the above-mentioned Ostomy products.</p> <p>ConvaTec has received complaints associated with use of Convex Two-Piece Skin Barriers where the starter hole (stoma hole) is off-centre on product manufactured from February 2017 to September 2018. Specifically, we have had three instances globally where use of an off-centre Convex Two-Piece Skin Barrier led to a stomal cut. These cuts did not require hospitalization and the end-user was able to resolve the issue. The overall risk of injury is low. However, since this risk can be visually detected, we are taking precautionary measures to notify customers of this matter so that the end-user can thoroughly check their product to see if it: a) was manufactured during this period and b) is visually off-centre.</p>
INFORMATION ON AFFECTED DEVICES	<ul style="list-style-type: none"> • SurFit™/ Combihesive™ Wafers: • System 92 Secure Moldable Convex Durahesive™ Wafer • System 92 Secure Durahesive Wafer with Convex

PRODUCT IDENTIFICATION

How to tell the difference between affected product and unaffected product	<p><u>Confirmation of Specific Product Codes and Lot Numbers:</u></p> <p>This issue is limited to specific product codes and lot numbers, identified in Attachment 3.</p> <p>All customers may refer to Attachment 3 to verify if their product is included in the scope of this Field Safety Notice.</p> <p>Upon request, a customer scope listing based on your receipts will be provided. Please contact your ConvaTec Customer Services Team Member for assistance.</p> <p>For product in scope of this Field Safety Notice, End Users may refer to Attachment 1 for examples of package labeling highlighting the location of the LOT No. and Product Code / REF No.</p> <p><u>Ostomy Product End User Visual Verification of Potentially Affected Product:</u></p> <p>Ostomy Product End Users may refer to Attachment 2 for additional information on how to visually verify potentially affected product and / or contact the FSCA Help Desk should they have additional questions.</p>
---	---

DISTRIBUTOR ACTIONS

1	<p>Verify if you have affected product in your inventory and identify which of your customers have received potentially affected product using Attachment 3: Product List showing affected Product Codes/REF No. and LOT Numbers.</p>
2	<p>If you have distributed this product:</p> <ul style="list-style-type: none"> • To other Wholesalers: <ul style="list-style-type: none"> ○ Please forward this Field Safety Notice, including attachments, to them and ask that they follow these Distributor Actions. • To End Users <ul style="list-style-type: none"> ○ Please forward this Field Safety Notice, including attachments, to your end-users and ask that they follow the End User Actions Section (below).
3	<p>If you have received returned product from your customers that they have identified as off-centre or if you have any stock on hand which you have verified as included in the scope of this Field Safety Notice (as described in Attachment 3: Product List showing affected Product Codes/REF No. and LOT Numbers), you do not need to open or inspect the product, ConvaTec will accept product returns and provide full credit to customers. Replacement product can be ordered via the normal ordering process.</p> <p>Contact the FSCA Help Desk for further instructions on how to return affected product and to arrange credit.</p>
4	<p>Action: Please complete the Field Safety Notice Distributor Response Form and return to the FSCA Help Desk by fax or email as soon as possible. Include details of product to be returned, if applicable.</p> <p>Action: Please send a complete list of all consignees to the FSCA Help Desk by fax or email as soon as possible. Please indicate if you have sent this Field Safety Notice, including attachments, to your consignees. This information is not required but important to allow ConvaTec to perform effectiveness checks.</p> <p>Please contact the FSCA Help Desk if you have any additional questions.</p>

PLEASE PROVIDE A COMPLETED RESPONSE AS SOON AS POSSIBLE

RETAILER ACTIONS

1	<p>For Retailers with a physical location (e.g. a “brick and mortar” store)</p> <ul style="list-style-type: none"> Please post Attachment 2: Ostomy Product User Letter of this Field Safety Notice in a conspicuous location in your store. <p>For Retailers which have a digital location (e.g. an eRetailer)</p> <ul style="list-style-type: none"> Please forward Attachment 2: Ostomy Product User Letter and Attachment 3: Product List showing affected Product Codes/REF No. and LOT Numbers via physical mail or e-mail.
2	<p>If you have received returned product from your customers that they have identified as off-centre or if you have any stock on hand which you have verified as included in the scope of this Field Safety Notice (as described in Attachment 3: Product List showing affected Product Codes/REF No. and LOT Numbers), you do not need to open or inspect the product, ConvaTec will accept product returns and provide full credit to customers. Replacement product can be ordered via the normal ordering process.</p> <p>Contact the FSCA Help Desk for further instructions on how to return affected product and to arrange credit.</p>
3	<p>Action: Please complete the Field Safety Notice Distributor Response Form and return to the FSCA Help Desk by fax or email as soon as possible. Include details of product to be returned, if applicable.</p> <p>Action: Please send a complete list of all consignees to the FSCA Help Desk by fax or email as soon as possible. Please indicate if you have sent this Field Safety Notice, including attachments, to your consignees. This information is not required but important to allow ConvaTec to perform effectiveness checks.</p> <p>Please contact the FSCA Help Desk if you have any additional questions.</p>

PLEASE PROVIDE A COMPLETED RESPONSE AS SOON AS POSSIBLE

END USER ACTIONS (CLINICIANS, HOSPITALS, PATIENTS, OTHERS)

1	<p>Based on our records, the product you have received may potentially have an off-centre starter hole (stoma hole) which is visibly off-centre. Example images are provided in Attachment 1: Photographic Examples of Labeling and Affected Product.</p>
2	<p>This issue is limited to specific product codes and lots. Please verify that your lot is associated with the product included in the scope of this Field Safety Notice utilising Attachment 3: Product List showing affected Product Codes/REF No. and LOT Numbers.</p> <p>If you do not know where the product code and/ or lot number are located, please refer to Attachment 1: Photographic Examples of Labeling and Affected Product.</p>
3	<p>Additional supporting information is included for you in Attachment 2: Ostomy Product User Letter. Please provide this to patients via physical mail or e-mail, as applicable.</p>
3	<p>If you have any product you have identified as off-centre in your possession, please contact the FSCA Help Desk.</p> <p>ConvaTec will accept product returns and replace product free of charge. You do not need to open or inspect this product.</p>
4	<p>Action: Please complete the Field Safety Notice End User Response Form and return to the FSCA Help Desk by fax or email as soon as possible. Include details of product to be returned, if applicable.</p> <p>This information is required to allow ConvaTec to perform effectiveness checks.</p> <p>Please contact the FSCA Help Desk if you have any additional questions.</p>


PLEASE PROVIDE A COMPLETED RESPONSE AS SOON AS POSSIBLE

TRANSMISSION OF THIS FIELD SAFETY NOTICE	<p>This notice needs to be passed on all those who need to be aware within your organisation or to any organisation where the potentially affected devices have been transferred. (As appropriate)</p> <p>Please transfer this notice to other organisations on which this action has an impact. (As appropriate)</p> <p>Please maintain awareness on this notice and resulting action for an appropriate period to ensure effectiveness of the corrective action.</p> <p>Please report all device-related incidents to the manufacturer, distributor or local representative, and the national Competent Authority if appropriate, as this provides important feedback.</p>
---	--

ConvaTec is committed to providing quality products and services to our customers and we sincerely apologise for any inconvenience this notice may cause.

The relevant National Authorities have been advised about this Field Safety Corrective Action.

AUTHORISATION

<u>Name</u>	<u>Title</u>	<u>Address</u>	
David Pudwill	Senior Director, Ostomy Franchise Regulatory Affairs and Quality	ConvaTec Limited, First Avenue, Deeside Industrial Park, Deeside, CH5 2NU, U.K.	
<u>Date</u>	20 September 2019	<u>Signature</u>	

List of Attachments

Attachment 1: Photographic Examples of Labelling and Affected Product

Attachment 2: Ostomy Product User Letter

Attachment 3: Product List showing affected Product Codes/REF No. and LOT Numbers

'Regional ConvaTec Contact' for other support:

Country	Distributors, Retailers, Wholesalers should contact	End Users should contact
Austria	Stericycle Phone: 0800 281598 E-Mail: Convatec2p@stericycle.com	CIC Team Austria Phone: +43 800 21 63 39 E-Mail: convatec.austria@convatec.com
Belgium	Stericycle Phone: 0800 71 804 E-Mail: Convatec2p@stericycle.com	Francoise Engelen Phone: 0032-476461740 E-Mail: info.convatec@convatec.com
Bulgaria	Stericycle Phone: 44 20 8089 3767 E-Mail: Convatec2p@stericycle.com	Please contact your local distributor
Croatia	Stericycle Phone: 44 20 8089 3767 E-Mail: Convatec2p@stericycle.com	Please contact your local distributor
Czech Republic	ConvaTec Czech Republic Phone: 800321222 E-Mail: info@convatec.cz	ConvaTec Czech Republic Phone: 800122111 E-Mail: info@convatec.cz
Denmark	Stericycle Phone: 80 82 00 54 E-Mail: Convatec2p@stericycle.com	CIC Team Denmark Phone: +45 48167475 E-Mail: convatec.danmark@convatec.com
Estonia	Stericycle Phone: 44 20 8089 3767 E-Mail: Convatec2p@stericycle.com	Please contact your local distributor
Finland	Stericycle Phone: 0800 413608 E-Mail: Convatec2p@stericycle.com	ConvaTec Oy Phone: +358 20 7659630 E-Mail: mail.fi@convatec.com
France	Stericycle Phone: 0 805 08 02 60 E-Mail: Convatec2p@stericycle.com	ConvaTec Customer Service Team Phone: +33800 358 480 E-Mail: fr.serviceclient@convatec.com
Germany	Stericycle Phone: 0800-0008875 E-Mail: Convatec2p@stericycle.com	CIC Team Germany Phone: +49 800 766 200 E-Mail: convatec.servicede@convatec.com
Greece	ConvaTec Customer Support Team Phone: 30-210-6564740 E-Mail: Lila.Karipidi@convatec.com	ConvaTec Customer Support Team Phone: 30-210-6564740 E-Mail: Lila.Karipidi@convatec.com
Hungary	Stericycle Phone: 44 20 8089 3767 E-Mail: Convatec2p@stericycle.com	Please contact your local distributor
Iceland	Stericycle Phone: 44 20 8089 3767 E-Mail: Convatec2p@stericycle.com	Medor ehf Phone: 665 7011 E-Mail: signyf@medor.is
Ireland	Stericycle Phone: 44 20 8089 3767	Tracey Fairclough Phone: 07442188256

Country	Distributors, Retailers, Wholesalers should contact	End Users should contact
	E-Mail: Convatec2p@stericycle.com	E-Mail: tracey.fairclough@amcaregroup.co.uk
Israel	Stericycle Phone: 44 20 8089 3767 E-Mail: Convatec2p@stericycle.com	Please contact your local distributor
Italy	Stericycle Phone: 800-789926 E-Mail: Convatec2p@stericycle.com	ConvaTel Italia Phone: 800930930 E-Mail: convatel.italia@convatec.com
Luxembourg	Stericycle Phone: 44 20 8089 3767 E-Mail: Convatec2p@stericycle.com	Francoise Engelen Phone: 0032-476461740 E-Mail: info.convatec@convatec.com
Malta	Stericycle Phone: 44 20 8089 3767 E-Mail: Convatec2p@stericycle.com	Ian Pace Phone: +356-2133-4412 E-Mail: ian@ipl.com.mt
Netherlands	Stericycle Phone: 0800 0233049 E-Mail: Convatec2p@stericycle.com	CIC Woerden Phone: 0031-348436941 E-Mail: info.convatec@convatec.com
Norway	Stericycle Phone: 800 24 720 E-Mail: Convatec2p@stericycle.com	CIC Norway Phone: 800 30 995 E-Mail: convatec.kundeservice@convatec.com
Poland	Dział Pomocy ConvaTec Phone: 800120093 E-Mail: dzial.pomocy@convatec.com	Dział Pomocy ConvaTec Phone: 800120093 E-Mail: dzial.pomocy@convatec.com
Portugal	Customer Service Iberia Phone: 93 602 3700 E-Mail: CustomerServiceIberia@convatec.com	ConvaTel Phone: 800201678 E-Mail: convatel.pt@convatec.com
Romania	Stericycle Phone: 44 20 8089 3767 E-Mail: Convatec2p@stericycle.com	Please contact your local distributor
Slovakia	Stericycle Phone: 44 20 8089 3767 E-Mail: Convatec2p@stericycle.com	Zákaznícke centrum ConvaTec Nitra Slovakia Phone: +421 800 800 111 E-Mail: stomalinka@convatec.com
Slovenia	Stericycle Phone: 44 20 8089 3767 E-Mail: Convatec2p@stericycle.com	Please contact your local distributor
Spain	Customer Service Iberia Phone: 93 602 3700 E-Mail: CustomerServiceIberia@convatec.com	InfoTec Phone: 936023737 E-Mail: infotec@convatec.com
Sweden	Stericycle Phone: 020-10 93 35 E-Mail: Convatec2p@stericycle.com	Nordic Customer Service Phone: +46(0)42332010 E-Mail: customerservicenordic@convatec.com
Switzerland	Stericycle Phone: 0800-848035 E-Mail: Convatec2p@stericycle.com	ConvaTec Customer Service Phone: +41 800 55 11 10 E-Mail: convatec.switzerland@convatec.com

Country	Distributors, Retailers, Wholesalers should contact	End Users should contact
Turkey	Doğa Soysal, Fulya Tüzüner, or İrem Taspınar Phone: <ul style="list-style-type: none"> • +90 2164165200 • +90 3122198361 E-Mail: <ul style="list-style-type: none"> • doga.soysal@convatec.com • fulya.tuzuner@convatec.com • irem.taspinar@convatec.com 	Doğa Soysal, Fulya Tüzüner, or İrem Taspınar Phone: <ul style="list-style-type: none"> • +90 2164165200 • +90 3122198361 E-Mail: <ul style="list-style-type: none"> • doga.soysal@convatec.com • fulya.tuzuner@convatec.com • irem.taspinar@convatec.com
United Kingdom	Stericycle Phone: 0800 069 8202 E-Mail: Convatec2p@stericycle.com	Tracey Fairclough Phone: 07442188256 E-Mail: tracey.fairclough@amcaregroup.co.uk

**FIELD SAFETY NOTICE
DISTRIBUTOR RESPONSE FORM
PLEASE COMPLETE AND RETURN by Email**

Consignee of the device:

Consignee Account No:	
Consignee Name:	
Consignee Address:	

The following product, Convex Two-Piece Skin Barriers have been distributed to your facility:

Invoice #	Sales Order #	Product Code / REF No.	SAP Code	LOT No.	Quantity Delivered (pieces)

Distributors (Tick all that apply and give details, where applicable)	
<input type="checkbox"/>	I confirm the receipt, the reading and understanding of the Field Safety Notice.
<input type="checkbox"/>	I have checked my stock and quarantined inventory
<input type="checkbox"/>	I have identified customers that received or may have received this device
<input type="checkbox"/>	I have attached customer list
<input type="checkbox"/>	I have informed the identified customers of this FSN Date sent:
<input type="checkbox"/>	I have received confirmation of reply from all identified customers
<input type="checkbox"/>	Neither I nor any of my customers has any affected devices in inventory

It is important that your organisation takes the actions detailed in the FSN and confirms that you have received the FSN.
Your organisation's reply is the evidence we need to monitor the progress of this action.



Table 1. Quarantined Inventory: Record quantity (pieces) for each LOT to be returned.

LOT No.	Units on Hand	LOT No.	Units on Hand	LOT No.	Units on Hand	LOT No.	Units on Hand

Table 2. Customer List: Please provide details of affected Convex Two-Piece Skin Barrier products that were distributed to your customers.

Customer Name	Product Code / REF No.	SAP Code	LOT No.	Quantity (pieces)

FORM Completed and Returned From:	
Name (CAPITAL LETTERS):	
Position:	
Company Name:	
Address:	
Phone No:	
Signature:	
Date (dd/mmm/yyyy):	

**FIELD SAFETY NOTICE
CUSTOMER RESPONSE FORM
PLEASE COMPLETE AND RETURN by Email**

Consignee of the device:

Consignee Account No:	
Consignee Name:	
Consignee Address:	

The following products, Convex Two-Piece Skin Barriers have been distributed to your facility:

Invoice #	Sales Order #	Product Code / REF No.	SAP Code	LOT No.	Quantity Delivered (pieces)

Customer action undertaken on behalf of Healthcare Organisation (Tick all that apply)		
<input type="checkbox"/>	I confirm receipt of the Field Safety Notice and that I read and understand its content.	
<input type="checkbox"/>	I performed all actions requested by the FSN.	
<input type="checkbox"/>	The information and required actions have been brought to the attention of all relevant users and executed.	
<input type="checkbox"/>	I have checked my stock and quarantined inventory	
<input type="checkbox"/>	No affected devices are available for return	




It is important that your organisation takes the actions detailed in the FSN and confirms that you have received the FSN.
Your organisation's reply is the evidence we need to monitor the progress of this action.

Table 1. Quarantined Inventory: Record quantity (pieces) for each LOT to be returned.

LOT No.	Units on Hand	LOT No.	Units on Hand	LOT No.	Units on Hand	LOT No.	Units on Hand

FORM Completed and Returned From:	
Name (CAPITAL LETTERS):	
Position:	
Company Name:	
Address:	
Phone No:	
Signature:	
Date (dd/mmm/yyyy):	

Attachment 1: Photographic Examples of Labelling and Affected Product

For the AFFECTED product	
<p>...the shipper and market unit labeling will include the product code and lot number.</p> <p>See example photos to the right.</p> <p>Product Code (REF/SKU)</p> <p>Lot Number</p>	
<p>... the physical product will present as a skin barrier with a starter hole (stoma hole), as shown to the right.</p>	<div style="display: flex; justify-content: space-around;"> <div style="text-align: center;"> <p>Off Center Starter Hole (Stoma Hole) This product may be out of specification</p>  <p>Figure 1</p> </div> <div style="text-align: center;"> <p>The starter hole should be "centered" in relation to the flange.</p> <p>The red lines in the photo on the left should be more equal.</p> <p>The green lines in the photo on the right are more equal.</p> </div> <div style="text-align: center;"> <p>Centered Starter Hole (Stoma Hole) This product should be within specification</p>  <p>Figure 2</p> </div> </div>

**Urgent
Medical Device Field Safety Notice**

Dear Ostomy Product User,


We are writing to notify you of a visually-detectable problem affecting convex two-piece skin barriers manufactured from February 2017 to September 2018, specifically where the starter hole (stoma hole) is off-center. (For a list of affected product, see Appendix A of this letter.) The risk of injury is low and is contained to specific lots manufactured over a specific timeframe in the past. Not all product or boxes belonging to the specific lots have this problem. However, since this risk can be visually-detected, we are taking precautionary measures to notify you of this matter so that you can thoroughly check your product to see if it:

- a) was manufactured during this period and
- b) is visually off-center.

Our records indicate you may be in receipt of this product. Example images are provided below:

Off Center Starter Hole (Stoma Hole)

This product may be out of specification




The starter hole should be "centered" in relation to the flange.

The red lines in the photo on the left should be more equal.

The green lines in the photo on the right are more equal.

Centered Starter Hole (Stoma Hole)

This product should be within specification





We are sending this communication to make you aware of the issue and outlining the steps you can take if you identify off-center product at point of use.

- Please examine your inventory. If you have received product that resembles **Figure 1 on the left**, isolate the affected product and contact the ConvaTec Customer Service Team as listed in Appendix B to arrange for return and replacement of the affected product.
- If you have received product that resembles **Figure 2 on the right**, no additional action is needed.
- If you are unsure, please contact the ConvaTec Customer Service Team for assistance.

The ConvaTec Customer Service Teams are a group of highly trained product specialists and nurses to help you get the care and information you need throughout your journey as an ostomate.

ConvaTec strives to provide the highest quality products and services to every customer. Our goal is to manufacture products that offer comfort, confidence and security. We want to reassure you that we have addressed the issue on the manufacturing line and corrected product has been distributed into the market.

If you have any further questions or concerns about this notice, please contact the applicable country-specific point of contact as outlined in Appendix B.

Thank you for being a valued ostomy product user. ConvaTec is passionate about improving people's lives and we put people at the center of everything we do.

Sincerely,

A handwritten signature in black ink, appearing to read "David Pudwill". The signature is stylized and cursive.

David Pudwill
Senior Director, Ostomy Franchise Regulatory Affairs and Quality Assurance
ConvaTec

'Regional ConvaTec Contact' for other support:

Country	Distributors, Retailers, Wholesalers should contact	End Users should contact
Austria	Stericycle Phone: 0800 281598 E-Mail: Convatec2p@stericycle.com	CIC Team Austria Phone: +43 800 21 63 39 E-Mail: convatec.austria@convatec.com
Belgium	Stericycle Phone: 0800 71 804 E-Mail: Convatec2p@stericycle.com	Francoise Engelen Phone: 0032-476461740 E-Mail: info.convatec@convatec.com
Bulgaria	Stericycle Phone: 44 20 8089 3767 E-Mail: Convatec2p@stericycle.com	Please contact your local distributor
Croatia	Stericycle Phone: 44 20 8089 3767 E-Mail: Convatec2p@stericycle.com	Please contact your local distributor
Czech Republic	ConvaTec Czech Republic Phone: 800321222 E-Mail: info@convatec.cz	ConvaTec Czech Republic Phone: 800122111 E-Mail: info@convatec.cz
Denmark	Stericycle Phone: 80 82 00 54 E-Mail: Convatec2p@stericycle.com	CIC Team Denmark Phone: +45 48167475 E-Mail: convatec.denmark@convatec.com
Estonia	Stericycle Phone: 44 20 8089 3767 E-Mail: Convatec2p@stericycle.com	Please contact your local distributor
Finland	Stericycle Phone: 0800 413608 E-Mail: Convatec2p@stericycle.com	ConvaTec Oy Phone: +358 20 7659630 E-Mail: mail.fi@convatec.com
France	Stericycle Phone: 0 805 08 02 60 E-Mail: Convatec2p@stericycle.com	ConvaTec Customer Service Team Phone: +33800 358 480 E-Mail: fr.serviceclient@convatec.com
Germany	Stericycle Phone: 0800-0008875 E-Mail: Convatec2p@stericycle.com	CIC Team Germany Phone: +49 800 766 200 E-Mail: convatec.servicede@convatec.com
Greece	ConvaTec Customer Support Team Phone: 30-210-6564740 E-Mail: Lila.Karipidi@convatec.com	ConvaTec Customer Support Team Phone: 30-210-6564740 E-Mail: Lila.Karipidi@convatec.com
Hungary	Stericycle Phone: 44 20 8089 3767 E-Mail: Convatec2p@stericycle.com	Please contact your local distributor
Iceland	Stericycle Phone: 44 20 8089 3767 E-Mail: Convatec2p@stericycle.com	Medor ehf Phone: 665 7011 E-Mail: signyf@medor.is
Ireland	Stericycle	Tracey Fairclough

Country	Distributors, Retailers, Wholesalers should contact	End Users should contact
	Phone: 44 20 8089 3767 E-Mail: Convatec2p@stericycle.com	Phone: 07442188256 E-Mail: tracey.fairclough@amcaregroup.co.uk
Israel	Stericycle Phone: 44 20 8089 3767 E-Mail: Convatec2p@stericycle.com	Please contact your local distributor
Italy	Stericycle Phone: 800-789926 E-Mail: Convatec2p@stericycle.com	ConvaTel Italia Phone: 800930930 E-Mail: convatel.italia@convatec.com
Luxembourg	Stericycle Phone: 44 20 8089 3767 E-Mail: Convatec2p@stericycle.com	Francoise Engelen Phone: 0032-476461740 E-Mail: info.convatec@convatec.com
Malta	Stericycle Phone: 44 20 8089 3767 E-Mail: Convatec2p@stericycle.com	Ian Pace Phone: +356-2133-4412 E-Mail: ian@ipl.com.mt
Netherlands	Stericycle Phone: 0800 0233049 E-Mail: Convatec2p@stericycle.com	CIC Woerden Phone: 0031-348436941 E-Mail: info.convatec@convatec.com
Norway	Stericycle Phone: 800 24 720 E-Mail: Convatec2p@stericycle.com	CIC Norway Phone: 800 30 995 E-Mail: convatec.kundeservice@convatec.com
Poland	Dział Pomocy ConvaTec Phone: 800120093 E-Mail: dzial.pomocy@convatec.com	Dział Pomocy ConvaTec Phone: 800120093 E-Mail: dzial.pomocy@convatec.com
Portugal	Customer Service Iberia Phone: 93 602 3700 E-Mail: CustomerServiceIberia@convatec.com	ConvaTel Phone: 800201678 E-Mail: convatel.pt@convatec.com
Romania	Stericycle Phone: 44 20 8089 3767 E-Mail: Convatec2p@stericycle.com	Please contact your local distributor
Slovakia	Stericycle Phone: 44 20 8089 3767 E-Mail: Convatec2p@stericycle.com	Zákaznícke centrum ConvaTec Nitra Slovakia Phone: +421 800 800 111 E-Mail: stomalinka@convatec.com
Slovenia	Stericycle Phone: 44 20 8089 3767 E-Mail: Convatec2p@stericycle.com	Please contact your local distributor
Spain	Customer Service Iberia Phone: 93 602 3700 E-Mail: CustomerServiceIberia@convatec.com	InfoTec Phone: 936023737 E-Mail: infotec@convatec.com
Sweden	Stericycle Phone: 020-10 93 35 E-Mail: Convatec2p@stericycle.com	Nordic Customer Service Phone: +46(0)42332010 E-Mail: customerservicenordic@convatec.com
Switzerland	Stericycle Phone: 0800-848035	ConvaTec Customer Service Phone: +41 800 55 11 10

Country	Distributors, Retailers, Wholesalers should contact	End Users should contact
	E-Mail: Convatec2p@stericycle.com	E-Mail: convatec.switzerland@convatec.com
Turkey	Doğa Soysal, Fulya Tüzüner, or İrem Taspınar Phone: <ul style="list-style-type: none"> • +90 2164165200 • +90 3122198361 E-Mail: <ul style="list-style-type: none"> • doga.soysal@convatec.com • fulya.tuzuner@convatec.com • irem.taspınar@convatec.com 	Doğa Soysal, Fulya Tüzüner, or İrem Taspınar Phone: <ul style="list-style-type: none"> • +90 2164165200 • +90 3122198361 E-Mail: <ul style="list-style-type: none"> • doga.soysal@convatec.com • fulya.tuzuner@convatec.com • irem.taspınar@convatec.com
United Kingdom	Stericycle Phone: 0800 069 8202 E-Mail: Convatec2p@stericycle.com	Tracey Fairclough Phone: 07442188256 E-Mail: tracey.fairclough@amcaregroup.co.uk

ATTACHMENT 3: Global List of Product Codes and Lot Numbers

PURPOSE

This attachment is a listing of all global products included in the scope of this Ostomy Field Safety Notice, which has been organized by Product Code. Within each product code, the lot numbers are organized numerically, smallest to largest.

Upon request, a custom scope listing based on your receipts will be provided. Please contact your ConvaTec Customer Service Team Member for assistance.

Column Header	Details	What Level of Packaging includes this information?
Product Code	The Product Code is also known as the International Commodity Code (ICC), the Catalog Number, or the SKU.	This number is included on the Shipper Label and the Market Unit Label.
Material	The Material is also known as the SAP Material Number. This is a supplemental number often utilized in distribution channels.	This number is located on the Shipper Label only and is not part of the Market Unit Labeling.
Material Description	<p>The Material Description is a description of the finished product from the internal ConvaTec system, SAP. It is a brief, coded description of the product. Utilizing the first line item example, Product Code 125026, Lot 8A06575, SURFIT PLUS WFR CVX 13X45MM (1X5PK):</p> <ul style="list-style-type: none"> • WFR: Wafer (Skin Barrier) • CVX: Convex • 13x45mm: 13mm stoma hole with a 45mm flange • 1x5pk: one market unit (1) includes five (5) skin barriers 	A description of the material is on the Shipper Label and on the Market Unit Labeling; however, the description on these labels may not exactly match (verbatim) the Material Description as listed in this attachment.
Lot Number	A traceable number which is specific to the product.	This number is included on the Shipper Label and the Market Unit Label.

Product Code	Material	Material Description	Lot Number
125026	1012910	SURFIT PLUS WFR CVX 13X45MM (1X5PK) FR	8A06575
125026	1053605	SURFIT PLUS WFR CVX 13X45MM (1X5PK) INT	8A06576
125027	1012911	SURFIT PLUS WFR CVX 16X45MM (1X5PK) FR	7G00162
125027	1053606	SURFIT PLUS WFR CVX 16X45MM (1X5PK) INT	7G00167
125027	1013422	COMBIHESIVE 2S WFR CVX 16/45MM(1X5PK)CEU	7G00168
125027	1012913	SURFIT PLUS WFR CVX 16X45MM (1X5PK) IT	7G00169
125027	1002346	COMBIHESIVE II WFR CVX 16X45MM(1X5PK)SCN	7J00668
125027	1012913	SURFIT PLUS WFR CVX 16X45MM (1X5PK) IT	7J00669
125027	1013422	COMBIHESIVE 2S WFR CVX 16/45MM(1X5PK)CEU	7J00670
125027	1002346	COMBIHESIVE II WFR CVX 16X45MM(1X5PK)SCN	8C03534
125027	1013422	COMBIHESIVE 2S WFR CVX 16/45MM(1X5PK)CEU	8C03535
125028	1012914	SURFIT PLUS WFR CVX 19X45MM (1X5PK) IT	7D00609
125028	1002348	SURFIT PLUS WFR CVX 19X45MM (1X5PK) INT	7D00610
125028	1012353	COMBIHESIVE WFR IIS D/CVX 19MM (1X5PK)BE	7D00611
125028	1012903	SURFIT PLUS WFR CVX 19X45MM (1X5PK) FR	7D00643
125028	1012914	SURFIT PLUS WFR CVX 19X45MM (1X5PK) IT	7L02963
125028	1012903	SURFIT PLUS WFR CVX 19X45MM (1X5PK) FR	7L02964
125028	1002347	COMBIHESIVE II WFR CVX 19X45MM(1X5PK)SCN	7L03102
125028	1002348	SURFIT PLUS WFR CVX 19X45MM (1X5PK) INT	8A00728
125028	1002348	SURFIT PLUS WFR CVX 19X45MM (1X5PK) INT	8B04995
125028	1012353	COMBIHESIVE WFR IIS D/CVX 19MM (1X5PK)BE	8B04996
125028	1012914	SURFIT PLUS WFR CVX 19X45MM (1X5PK) IT	8B04997
125028	1002348	SURFIT PLUS WFR CVX 19X45MM (1X5PK) INT	8D00465
125028	1012903	SURFIT PLUS WFR CVX 19X45MM (1X5PK) FR	8D00466
125028	1013423	COMBIHESIVE 2S WFR CVX 19/45MM(1X5PK)CEU	8D00467
125028	1013423	COMBIHESIVE 2S WFR CVX 19/45MM(1X5PK)CEU	8F03145
125028	1002348	SURFIT PLUS WFR CVX 19X45MM (1X5PK) INT	8F03149
125028	1012353	COMBIHESIVE WFR IIS D/CVX 19MM (1X5PK)BE	8F03150
125028	1012914	SURFIT PLUS WFR CVX 19X45MM (1X5PK) IT	8F03151
125028	1013423	COMBIHESIVE 2S WFR CVX 19/45MM(1X5PK)CEU	8G03301
125028	1012914	SURFIT PLUS WFR CVX 19X45MM (1X5PK) IT	8G03332
125029	1002349	COMBIHESIVE II WFR CVX 22X45MM(1X5PK)SCN	7C05870
125029	1012915	SURFIT PLUS WFR CVX 22X45MM (1X5PK) IT	7C05871
125029	1002350	SURFIT PLUS WFR CVX 22X45MM (1X5PK) INT	7G04445
125029	1002349	COMBIHESIVE II WFR CVX 22X45MM(1X5PK)SCN	7G04446
125029	1012915	SURFIT PLUS WFR CVX 22X45MM (1X5PK) IT	7G04448
125029	1012354	COMBIHESIVE WFR IIS D/CVX 22MM (1X5PK)BE	7G04449
125029	1012904	SURFIT PLUS WFR CVX 22X45MM (1X5PK) FR	7J01269
125029	1012915	SURFIT PLUS WFR CVX 22X45MM (1X5PK) IT	7J01270
125029	1002349	COMBIHESIVE II WFR CVX 22X45MM(1X5PK)SCN	7J01271

Product Code	Material	Material Description	Lot Number
125029	1002350	SURFIT PLUS WFR CVX 22X45MM (1X5PK) INT	7J01292
125029	1002349	COMBIHESIVE II WFR CVX 22X45MM(1X5PK)SCN	8B01186
125029	1012354	COMBIHESIVE WFR IIS D/CVX 22MM (1X5PK)BE	8D00470
125029	1012915	SURFIT PLUS WFR CVX 22X45MM (1X5PK) IT	8D00471
125029	1002350	SURFIT PLUS WFR CVX 22X45MM (1X5PK) INT	8D00472
125029	1002350	SURFIT PLUS WFR CVX 22X45MM (1X5PK) INT	8E00347
125029	1002349	COMBIHESIVE II WFR CVX 22X45MM(1X5PK)SCN	8E00349
125029	1012354	COMBIHESIVE WFR IIS D/CVX 22MM (1X5PK)BE	8E00350
125029	1012904	SURFIT PLUS WFR CVX 22X45MM (1X5PK) FR	8E00992
125029	1013424	COMBIHESIVE 2S WFR CVX 22/45MM(1X5PK)CEU	8E01071
125030	1002352	COMBIHESIVE II WFR CVX 25X45MM(1X5PK)SCN	7C05894
125030	1012355	COMBIHESIVE WFR IIS D/CVX 25MM (1X5PK)BE	7G03585
125030	1002353	COMBIHESIVE 2 WFR CVX 25X45MM (1X5PK) IT	7G03586
125030	1002353	COMBIHESIVE 2 WFR CVX 25X45MM (1X5PK) IT	8A02705
125030	1013425	COMBIHESIVE 2S WFR CVX 25/45MM(1X5PK)CEU	8A02706
125030	1002354	SURFIT PLUS WFR CVX 25X45MM (1X5PK) INT	8C04244
125030	1002352	COMBIHESIVE II WFR CVX 25X45MM(1X5PK)SCN	8F02589
125030	1002354	SURFIT PLUS WFR CVX 25X45MM (1X5PK) INT	8F02590
125030	1012355	COMBIHESIVE WFR IIS D/CVX 25MM (1X5PK)BE	8F04927
125030	1002354	SURFIT PLUS WFR CVX 25X45MM (1X5PK) INT	8F04928
125030	1013425	COMBIHESIVE 2S WFR CVX 25/45MM(1X5PK)CEU	8F04929
125030	1002351	COMBIHESIVE 2 WFR CVX 25X45MM(1X5PK) FR	8F04932
125031	1002358	SURFIT PLUS WFR CVX 28X45MM (1X5PK) INT	7C04211
125031	1002358	SURFIT PLUS WFR CVX 28X45MM (1X5PK) INT	7D01175
125031	1013426	COMBIHESIVE 2S WFR CVX 28/45MM(1X5PK)CEU	7D03659
125031	1013426	COMBIHESIVE 2S WFR CVX 28/45MM(1X5PK)CEU	7D03959
125031	1002356	COMBIHESIVE II WFR CVX 28X45MM(1X5PK)SCN	7F00603
125031	1012356	COMBIHESIVE WFR IIS D/CVX 28MM (1X5PK)BE	7F00604
125031	1002357	COMBIHESIVE 2 WFR CVX 28X45MM (1X5PK) IT	7F00605
125031	1002358	SURFIT PLUS WFR CVX 28X45MM (1X5PK) INT	7F00606
125031	1002356	COMBIHESIVE II WFR CVX 28X45MM(1X5PK)SCN	8B01383
125031	1013426	COMBIHESIVE 2S WFR CVX 28/45MM(1X5PK)CEU	8B01384
125031	1002357	COMBIHESIVE 2 WFR CVX 28X45MM (1X5PK) IT	8B01385
125031	1002358	SURFIT PLUS WFR CVX 28X45MM (1X5PK) INT	8B03229
125031	1012356	COMBIHESIVE WFR IIS D/CVX 28MM (1X5PK)BE	8B03230
125031	1002355	COMBIHESIVE 2 WFR CVX 28X45MM (1X5PK) FR	8B03231
125031	1002358	SURFIT PLUS WFR CVX 28X45MM (1X5PK) INT	8C00877
125031	1013426	COMBIHESIVE 2S WFR CVX 28/45MM(1X5PK)CEU	8C00879
125031	1002356	COMBIHESIVE II WFR CVX 28X45MM(1X5PK)SCN	8D02693
125031	1013426	COMBIHESIVE 2S WFR CVX 28/45MM(1X5PK)CEU	8D02694
125031	1012356	COMBIHESIVE WFR IIS D/CVX 28MM (1X5PK)BE	8D05081

Product Code	Material	Material Description	Lot Number
125031	1002356	COMBIHESIVE II WFR CVX 28X45MM(1X5PK)SCN	8D05092
125031	1002355	COMBIHESIVE 2 WFR CVX 28X45MM (1X5PK) FR	8D05093
125031	1013426	COMBIHESIVE 2S WFR CVX 28/45MM(1X5PK)CEU	8E01174
125031	1002358	SURFIT PLUS WFR CVX 28X45MM (1X5PK) INT	8E01175
125032	1002362	SURFIT PLUS WFR CVX 32X45MM (1X5PK) INT	7C04224
125032	1012357	COMBIHESIVE WFR IIS D/CVX 32MM (1X5PK)BE	7E01287
125032	1002360	COMBIHESIVE II WFR CVX 32X45MM(1X5PK)SCN	7E01288
125032	1002362	SURFIT PLUS WFR CVX 32X45MM (1X5PK) INT	7E01289
125032	1002359	COMBIHESIVE 2 WFR CVX 32X45MM (1X5PK) FR	7E01290
125032	1002361	COMBIHESIVE 2 WFR CVX 32X45MM (1X5PK) IT	7F00608
125033	1012358	COMBIHESIVE WFR IIS D/CVX 35MM (1X5PK)BE	7E01292
125033	1013428	COMBIHESIVE 2S WFR CVX 35/45MM(1X5PK)CEU	7E01293
125033	1002364	COMBIHESIVE II WFR CVX 35X45MM(1X5PK)SCN	7E01294
125033	1002363	COMBIHESIVE 2 WFR CVX 35X45MM (1X5PK) FR	7E01295
125033	1002366	SURFIT PLUS WFR CVX 35X45MM (1X5PK) INT	7E03364
125033	1002366	SURFIT PLUS WFR CVX 35X45MM (1X5PK) INT	8B04219
125033	1013428	COMBIHESIVE 2S WFR CVX 35/45MM(1X5PK)CEU	8B04220
125033	1002364	COMBIHESIVE II WFR CVX 35X45MM(1X5PK)SCN	8C03753
125034	1002367	COMBIHESIVE II WFR CVX 38X57MM(1X5PK)SCN	7E02412
125034	1013429	COMBIHESIVE 2S WFR CVX 38/57MM(1X5PK)CEU	7E02413
125034	1002367	COMBIHESIVE II WFR CVX 38X57MM(1X5PK)SCN	7H04904
125034	1012359	COMBIHESIVE WFR IIS D/CVX 38MM (1X5PK)BE	7H04905
125034	1013429	COMBIHESIVE 2S WFR CVX 38/57MM(1X5PK)CEU	7H04906
125034	1002368	SURFIT PLUS WFR CVX 38X57MM (1X5PK) INT	7H05074
125034	1012916	SURFIT PLUS WFR CVX 38X57MM (1X5PK) IT	8A04956
125034	1013429	COMBIHESIVE 2S WFR CVX 38/57MM(1X5PK)CEU	8A04957
125034	1002367	COMBIHESIVE II WFR CVX 38X57MM(1X5PK)SCN	8A04958
125034	1012359	COMBIHESIVE WFR IIS D/CVX 38MM (1X5PK)BE	8A04959
125034	1013429	COMBIHESIVE 2S WFR CVX 38/57MM(1X5PK)CEU	8B02249
125034	1002368	SURFIT PLUS WFR CVX 38X57MM (1X5PK) INT	8F02600
125034	1012359	COMBIHESIVE WFR IIS D/CVX 38MM (1X5PK)BE	8F02613
125034	1012916	SURFIT PLUS WFR CVX 38X57MM (1X5PK) IT	8F02614
125035	1013430	COMBIHESIVE 2S WFR CVX 41/57MM(1X5PK)CEU	7H04907
125035	1002370	SURFIT PLUS WFR CVX 41X57MM (1X5PK) INT	7H04908
125035	1002369	COMBIHESIVE II WFR CVX 41X57MM(1X5PK)SCN	7H04909
125035	1012376	COMBIHESIVE WFR IIS D/CVX 41MM (1X5PK)BE	7H04910
125035	1002370	SURFIT PLUS WFR CVX 41X57MM (1X5PK) INT	7K04883
125035	1013430	COMBIHESIVE 2S WFR CVX 41/57MM(1X5PK)CEU	7K04884
125035	1012917	SURFIT PLUS WFR CVX 41X57MM (1X5PK) IT	7K04885
125035	1012376	COMBIHESIVE WFR IIS D/CVX 41MM (1X5PK)BE	7K04886
125035	1002369	COMBIHESIVE II WFR CVX 41X57MM(1X5PK)SCN	8B01729

Product Code	Material	Material Description	Lot Number
125035	1013430	COMBIHESIVE 2S WFR CVX 41/57MM(1X5PK)CEU	8B01730
125036	1002372	SURFIT PLUS WFR CVX 45X57MM (1X5PK) INT	7J03062
125036	1013431	COMBIHESIVE 2S WFR CVX 45/57MM(1X5PK)CEU	7J03064
125036	1014051	SURFIT PLUS WFR CVX 45X57MM (1X5PK) IT	7J03065
125036	1002372	SURFIT PLUS WFR CVX 45X57MM (1X5PK) INT	8C02094
125036	1002371	COMBIHESIVE II WFR CVX 45X57MM(1X5PK)SCN	8C02095
125036	1013431	COMBIHESIVE 2S WFR CVX 45/57MM(1X5PK)CEU	8G02097
125037	1012919	SURFIT PLUS WFR CVX 50X57MM (1X5PK) IT	7L03278
125037	1013432	COMBIHESIVE 2S WFR CVX 50/57MM(1X5PK)CEU	7L03279
125037	1002374	SURFIT PLUS WFR CVX 50X57MM (1X5PK) INT	7L03280
125037	1002373	COMBIHESIVE II WFR CVX 50X57MM(1X5PK)SCN	7L03281
125037	1012919	SURFIT PLUS WFR CVX 50X57MM (1X5PK) IT	7M02850
125278	1051431	VARICARE S4S CVX WFR 13/45MM (1X5PK) JP	7J00006
125279	1051432	VARICARE S4S CVX WFR 16/45MM (1X5PK) JP	7F04003
125279	1051432	VARICARE S4S CVX WFR 16/45MM (1X5PK) JP	7J00008
125279	1051432	VARICARE S4S CVX WFR 16/45MM (1X5PK) JP	8F00683
125280	1051433	VARICARE S4S CVX WFR 19/45MM (1X5PK) JP	7M03972
125280	1051433	VARICARE S4S CVX WFR 19/45MM (1X5PK) JP	8A00733
125280	1051433	VARICARE S4S CVX WFR 19/45MM (1X5PK) JP	8E04639
125281	1051434	VARICARE S4S CVX WFR 22/45MM (1X5PK) JP	7E01836
125281	1051434	VARICARE S4S CVX WFR 22/45MM (1X5PK) JP	7K00524
125282	1051435	VARICARE S4S CVX WFR 25/45MM (1X5PK) JP	7J00001
125282	1051435	VARICARE S4S CVX WFR 25/45MM (1X5PK) JP	8E03843
125283	1051436	VARICARE S4S CVX WFR 28/45MM (1X5PK) JP	7F03154
125283	1051436	VARICARE S4S CVX WFR 28/45MM (1X5PK) JP	7M03580
125283	1051436	VARICARE S4S CVX WFR 28/45MM (1X5PK) JP	8F01049
125284	1051437	VARICARE S4S CVX WFR 32/45MM (1X5PK) JP	7H02190
125285	1051438	VARICARE S4S CVX WFR 35/45MM (1X5PK) JP	7J04969
125285	1051438	VARICARE S4S CVX WFR 35/45MM (1X5PK) JP	8G03496
402200	1052415	SURFIT NATURA S4S WFR CVX13/45MM1X5NL/GB	7J01294
402200	1052403	SURFIT NATURA S4S WFR CVX 13/45MM(1X5)DE	7J01295
402201	1052416	SURFIT NATURA S4S WFR CVX16/45MM1X5NL/GB	7F03877
402201	1052404	SURFIT NATURA S4S WFR CVX 16/45MM(1X5)DE	7F03878
402201	1052416	SURFIT NATURA S4S WFR CVX16/45MM1X5NL/GB	7H04125
402201	1052404	SURFIT NATURA S4S WFR CVX 16/45MM(1X5)DE	8A06015
402201	1052416	SURFIT NATURA S4S WFR CVX16/45MM1X5NL/GB	8D01629
402202	1052417	SURFIT NATURA S4S WFR CVX19/45MM1X5NL/GB	7H01478
402202	1052405	SURFIT NATURA S4S WFR CVX 19/45MM(1X5)DE	7H01479
402202	1052405	SURFIT NATURA S4S WFR CVX 19/45MM(1X5)DE	7K05366
402202	1052417	SURFIT NATURA S4S WFR CVX19/45MM1X5NL/GB	7K05368
402202	1052417	SURFIT NATURA S4S WFR CVX19/45MM1X5NL/GB	8D05099

Product Code	Material	Material Description	Lot Number
402202	1052417	SURFIT NATURA S4S WFR CVX19/45MM1X5NL/GB	8G04010
402202	1052405	SURFIT NATURA S4S WFR CVX 19/45MM(1X5)DE	8G04012
402203	1052418	SURFIT NATURA S4S WFR CVX22/45MM1X5NL/GB	7J00665
402203	1052406	SURFIT NATURA S4S WFR CVX 22/45MM(1X5)DE	7K03422
402203	1052418	SURFIT NATURA S4S WFR CVX22/45MM1X5NL/GB	7K03423
402203	1052406	SURFIT NATURA S4S WFR CVX 22/45MM(1X5)DE	7L01060
402203	1052418	SURFIT NATURA S4S WFR CVX22/45MM1X5NL/GB	7L01062
402203	1052406	SURFIT NATURA S4S WFR CVX 22/45MM(1X5)DE	8C00960
402203	1052418	SURFIT NATURA S4S WFR CVX22/45MM1X5NL/GB	8C03915
402204	1052407	SURFIT NATURA S4S WFR CVX 25/45MM(1X5)DE	7F00955
402204	1052419	SURFIT NATURA S4S WFR CVX25/45MM1X5NL/GB	7F00956
402204	1052419	SURFIT NATURA S4S WFR CVX25/45MM1X5NL/GB	7H01474
402204	1052407	SURFIT NATURA S4S WFR CVX 25/45MM(1X5)DE	7H01475
402204	1052419	SURFIT NATURA S4S WFR CVX25/45MM1X5NL/GB	7K02141
402204	1052407	SURFIT NATURA S4S WFR CVX 25/45MM(1X5)DE	7K02153
402204	1052419	SURFIT NATURA S4S WFR CVX25/45MM1X5NL/GB	8D05096
402204	1052407	SURFIT NATURA S4S WFR CVX 25/45MM(1X5)DE	8G01612
402204	1052419	SURFIT NATURA S4S WFR CVX25/45MM1X5NL/GB	8G01613
402204	1052419	SURFIT NATURA S4S WFR CVX25/45MM1X5NL/GB	8H04519
402204	1052407	SURFIT NATURA S4S WFR CVX 25/45MM(1X5)DE	8H04521
402205	1052408	SURFIT NATURA S4S WFR CVX 28/45MM(1X5)DE	7F06139
402205	1052408	SURFIT NATURA S4S WFR CVX 28/45MM(1X5)DE	7G03644
402205	1052420	SURFIT NATURA S4S WFR CVX28/45MM1X5NL/GB	7G03645
402205	1052420	SURFIT NATURA S4S WFR CVX28/45MM1X5NL/GB	7M03576
402205	1052408	SURFIT NATURA S4S WFR CVX 28/45MM(1X5)DE	8D00474
402205	1052420	SURFIT NATURA S4S WFR CVX28/45MM1X5NL/GB	8D00476
402205	1052408	SURFIT NATURA S4S WFR CVX 28/45MM(1X5)DE	8G04015
402205	1052420	SURFIT NATURA S4S WFR CVX28/45MM1X5NL/GB	8J00540
402206	1052421	SURFIT NATURA S4S WFR CVX32/45MM1X5NL/GB	7G03647
402206	1052409	SURFIT NATURA S4S WFR CVX 32/45MM(1X5)DE	7G03648
402206	1052409	SURFIT NATURA S4S WFR CVX 32/45MM(1X5)DE	7L01338
402206	1052421	SURFIT NATURA S4S WFR CVX32/45MM1X5NL/GB	7L01339
402206	1052421	SURFIT NATURA S4S WFR CVX32/45MM1X5NL/GB	7M03584
402206	1052421	SURFIT NATURA S4S WFR CVX32/45MM1X5NL/GB	8C00953
402206	1052421	SURFIT NATURA S4S WFR CVX32/45MM1X5NL/GB	8C03490
402207	1052410	SURFIT NATURA S4S WFR CVX 35/45MM(1X5)DE	7K02047
402207	1052422	SURFIT NATURA S4S WFR CVX35/45MM1X5NL/GB	7K02048
402207	1052422	SURFIT NATURA S4S WFR CVX35/45MM1X5NL/GB	8A00737
402207	1052410	SURFIT NATURA S4S WFR CVX 35/45MM(1X5)DE	8A00738
402207	1052422	SURFIT NATURA S4S WFR CVX35/45MM1X5NL/GB	8F00089
402207	1052422	SURFIT NATURA S4S WFR CVX35/45MM1X5NL/GB	8G03493

Product Code	Material	Material Description	Lot Number
402208	1052423	SURFIT NATURA S4S WFR CVX38/57MM1X5NL/GB	7H04929
402208	1052411	SURFIT NATURA S4S WFR CVX 38/57MM(1X5)DE	7H04930
402208	1052423	SURFIT NATURA S4S WFR CVX38/57MM1X5NL/GB	7J04231
402208	1052411	SURFIT NATURA S4S WFR CVX 38/57MM(1X5)DE	7J04233
402208	1052423	SURFIT NATURA S4S WFR CVX38/57MM1X5NL/GB	8H00327
402209	1052424	SURFIT NATURA S4S WFR CVX41/57MM1X5NL/GB	7H03495
402209	1052424	SURFIT NATURA S4S WFR CVX41/57MM1X5NL/GB	7H03496
402209	1052412	SURFIT NATURA S4S WFR CVX 41/57MM(1X5)DE	7H03497
402210	1052425	SURFIT NATURA S4S WFR CVX45/57MM1X5NL/GB	7H04917
402210	1052413	SURFIT NATURA S4S WFR CVX 45/57MM(1X5)DE	7H04918
402211	1052414	SURFIT NATURA S4S WFR CVX 50/57MM(1X5)DE	7H05736
402211	1052426	SURFIT NATURA S4S WFR CVX50/57MM1X5NL/GB	7H05737
125028	1210693	SURFIT PLUS WFR CVX 19X45MM (1X5PK) ES	7L02965
125028	1210693	SURFIT PLUS WFR CVX 19X45MM (1X5PK) ES	8A00729
125028	1210693	SURFIT PLUS WFR CVX 19X45MM (1X5PK) ES	8D00468
125028	1210693	SURFIT PLUS WFR CVX 19X45MM (1X5PK) ES	8F03146
125029	1210694	SURFIT PLUS WFR CVX 22X45MM (1X5PK) ES	8B01190
125030	1210695	COMBIHESIVE 2 WFR CVX 25X45MM (1X5PK) ES	7G03584
125030	1210695	COMBIHESIVE 2 WFR CVX 25X45MM (1X5PK) ES	8A02704
125030	1210695	COMBIHESIVE 2 WFR CVX 25X45MM (1X5PK) ES	8C04243
125030	1210695	COMBIHESIVE 2 WFR CVX 25X45MM (1X5PK) ES	8F02591
125030	1210695	COMBIHESIVE 2 WFR CVX 25X45MM (1X5PK) ES	8F04930
125031	1210696	COMBIHESIVE 2 WFR CVX 28X45MM (1X5PK) ES	8B03228
125032	1210697	COMBIHESIVE 2 WFR CVX 32X45MM (1X5PK) ES	7E03289
125033	1210698	COMBIHESIVE 2 WFR CVX 34X45MM (1X5PK) ES	8B04221
125033	1210698	COMBIHESIVE 2 WFR CVX 34X45MM (1X5PK) ES	8C03752
125278	1704525	NATURA WFR FLEX CVX 13/45MM (1X5PK) EUR	7L03067
125278	1135962	C/HESIVE NAT WFR FLEX 13/45MM (1X5)GB/NL	7L03070
125278	1136074	C/HESIVE NAT WFR FLEX 13/45MM (1X5PK) DE	7L03072
125279	1704527	NATURA WFR FLEX CVX 16/45MM (1X5PK) EUR	7E02269
125279	1136075	C/HESIVE NAT WFR FLEX 16/45MM (1X5PK) DE	7E02312
125279	1136075	C/HESIVE NAT WFR FLEX 16/45MM (1X5PK) DE	7K05185
125279	1704527	NATURA WFR FLEX CVX 16/45MM (1X5PK) EUR	7K05186
125279	1135963	C/HESIVE NAT WFR FLEX 16/45MM (1X5)GB/NL	7K05187
125279	1135963	C/HESIVE NAT WFR FLEX 16/45MM (1X5)GB/NL	8F00670
125279	1704527	NATURA WFR FLEX CVX 16/45MM (1X5PK) EUR	8F00671
125280	1136076	C/HESIVE NAT WFR FLEX 19/45MM (1X5PK) DE	7D01171
125280	1135964	C/HESIVE NAT WFR FLEX 19/45MM (1X5)GB/NL	7D01173
125280	1136076	C/HESIVE NAT WFR FLEX 19/45MM (1X5PK) DE	7D04717
125280	1135964	C/HESIVE NAT WFR FLEX 19/45MM (1X5)GB/NL	8A00732
125280	1136076	C/HESIVE NAT WFR FLEX 19/45MM (1X5PK) DE	8A00735

Product Code	Material	Material Description	Lot Number
125280	1135964	C/HESIVE NAT WFR FLEX 19/45MM (1X5)GB/NL	8E04640
125281	1136077	C/HESIVE NAT WFR FLEX 22/45MM (1X5PK) DE	7C05397
125281	1135965	C/HESIVE NAT WFR FLEX 22/45MM (1X5)GB/NL	7C05398
125281	1704529	NATURA WFR FLEX CVX 22/45MM (1X5PK) EUR	7E01837
125281	1136077	C/HESIVE NAT WFR FLEX 22/45MM (1X5PK) DE	7E01838
125281	1135965	C/HESIVE NAT WFR FLEX 22/45MM (1X5)GB/NL	7K00519
125281	1136077	C/HESIVE NAT WFR FLEX 22/45MM (1X5PK) DE	7K00520
125281	1135965	C/HESIVE NAT WFR FLEX 22/45MM (1X5)GB/NL	7K05908
125281	1704529	NATURA WFR FLEX CVX 22/45MM (1X5PK) EUR	7K05911
125281	1704529	NATURA WFR FLEX CVX 22/45MM (1X5PK) EUR	7L01342
125281	1136077	C/HESIVE NAT WFR FLEX 22/45MM (1X5PK) DE	7L01343
125281	1135965	C/HESIVE NAT WFR FLEX 22/45MM (1X5)GB/NL	8E01678
125282	1704530	NATURA WFR FLEX CVX 25/45MM (1X5PK) EUR	7E01711
125282	1704530	NATURA WFR FLEX CVX 25/45MM (1X5PK) EUR	7F01779
125282	1704530	NATURA WFR FLEX CVX 25/45MM (1X5PK) EUR	7G00165
125282	1136078	C/HESIVE NAT WFR FLEX 25/45MM (1X5PK) DE	7G02896
125282	1136078	C/HESIVE NAT WFR FLEX 25/45MM (1X5PK) DE	7M03979
125282	1704530	NATURA WFR FLEX CVX 25/45MM (1X5PK) EUR	8B00892
125282	1704530	NATURA WFR FLEX CVX 25/45MM (1X5PK) EUR	8E03842
125282	1135966	C/HESIVE NAT WFR FLEX 25/45MM (1X5)GB/NL	8H01556
125282	1136078	C/HESIVE NAT WFR FLEX 25/45MM (1X5PK) DE	8H01557
125283	1135967	C/HESIVE NAT WFR FLEX 28/45MM (1X5)GB/NL	7C05414
125283	1704531	NATURA WFR FLEX CVX 28/45MM (1X5PK) EUR	7F03153
125283	1136079	C/HESIVE NAT WFR FLEX 28/45MM (1X5PK) DE	7H04122
125283	1704531	NATURA WFR FLEX CVX 28/45MM (1X5PK) EUR	7L01043
125283	1136079	C/HESIVE NAT WFR FLEX 28/45MM (1X5PK) DE	7L01045
125283	1704531	NATURA WFR FLEX CVX 28/45MM (1X5PK) EUR	7M03582
125283	1135967	C/HESIVE NAT WFR FLEX 28/45MM (1X5)GB/NL	7M03583
125283	1136079	C/HESIVE NAT WFR FLEX 28/45MM (1X5PK) DE	8J00539
125284	1135968	C/HESIVE NAT WFR FLEX 32/45MM (1X5)GB/NL	7H02189
125284	1704532	NATURA WFR FLEX CVX 32/45MM (1X5PK) EUR	8C03529
125285	1136081	C/HESIVE NAT WFR FLEX 35/45MM (1X5PK) DE	7D04718
125285	1704533	NATURA WFR FLEX CVX 35/45MM (1X5PK) EUR	7F01762
125285	1135969	C/HESIVE NAT WFR FLEX 35/45MM (1X5)GB/NL	7J02058
125285	1704533	NATURA WFR FLEX CVX 35/45MM (1X5PK) EUR	8G03495
125286	1136082	C/HESIVE NAT WFR FLEX 38/57MM (1X5PK) DE	7F05915
125286	1136082	C/HESIVE NAT WFR FLEX 38/57MM (1X5PK) DE	7H06312
125286	1135970	C/HESIVE NAT WFR FLEX 38/57MM (1X5)GB/NL	7H06313
125286	1704534	NATURA WFR FLEX CVX 38/57MM (1X5PK) EUR	7H06316
125286	1136082	C/HESIVE NAT WFR FLEX 38/57MM (1X5PK) DE	7J04227
125286	1704534	NATURA WFR FLEX CVX 38/57MM (1X5PK) EUR	7J04228

Product Code	Material	Material Description	Lot Number
125286	1135970	C/HESIVE NAT WFR FLEX 38/57MM (1X5)GB/NL	8A05991
125286	1704534	NATURA WFR FLEX CVX 38/57MM (1X5PK) EUR	8E01296
125287	1136071	C/HESIVE NAT WFR FLEX 41/57MM (1X5)GB/NL	8A06017
125287	1136083	C/HESIVE NAT WFR FLEX 41/57MM (1X5PK) DE	8A06018
125287	1704535	NATURA WFR FLEX CVX 41/57MM (1X5PK) EUR	8A06019
125287	1136083	C/HESIVE NAT WFR FLEX 41/57MM (1X5PK) DE	8D00709
125287	1136071	C/HESIVE NAT WFR FLEX 41/57MM (1X5)GB/NL	8D05101
125288	1704536	NATURA WFR FLEX CVX 45/57MM (1X5PK) EUR	7H04927
125288	1704536	NATURA WFR FLEX CVX 45/57MM (1X5PK) EUR	7M00468
125288	1136084	C/HESIVE NAT WFR FLEX 45/57MM (1X5PK) DE	7M00469
125289	1700992	NATURA WFR FLEX CVX 50/57MM (1X5PK) EUR	7G01363
125289	1136073	C/HESIVE NAT WFR FLEX 50/57MM (1X5)GB/NL	7G01364
125289	1136085	C/HESIVE NAT WFR FLEX 50/57MM (1X5PK) DE	7G01366
125289	1700992	NATURA WFR FLEX CVX 50/57MM (1X5PK) EUR	7J00630
125289	1136085	C/HESIVE NAT WFR FLEX 50/57MM (1X5PK) DE	7K07046
125289	1136073	C/HESIVE NAT WFR FLEX 50/57MM (1X5)GB/NL	7K07047
404592	1156500	NATURA WFR MOLDABLE CVX 13/45MM (1X10)US	7C05400
404592	1156500	NATURA WFR MOLDABLE CVX 13/45MM (1X10)US	7C05401
404592	1175793	C/HESIVE NAT WFR MLD CVX 13/22MM(1X10)JP	7F06140
404592	1175760	C/HESIVE NAT WFR MLD CVX 13/22MM(1X10)DE	7F06141
404592	1175793	C/HESIVE NAT WFR MLD CVX 13/22MM(1X10)JP	7K00745
404592	1701043	NATURA MLDC CVX 13/22MM (1X10PK) EUR	7K00746
404592	1174445	C/HESIVE NAT MLDC CVX13/22MM(1X10) GB/NL	7M00630
404592	1175760	C/HESIVE NAT WFR MLD CVX 13/22MM(1X10)DE	7M00631
404592	1156500	NATURA WFR MOLDABLE CVX 13/45MM (1X10)US	7M00632
404592	1701043	NATURA MLDC CVX 13/22MM (1X10PK) EUR	7M03042
404592	1156500	NATURA WFR MOLDABLE CVX 13/45MM (1X10)US	7M03044
404592	1156500	NATURA WFR MOLDABLE CVX 13/45MM (1X10)US	8A00115
404592	1156500	NATURA WFR MOLDABLE CVX 13/45MM (1X10)US	8A03532
404592	1701043	NATURA MLDC CVX 13/22MM (1X10PK) EUR	8A04344
404592	1156500	NATURA WFR MOLDABLE CVX 13/45MM (1X10)US	8A04345
404592	1156500	NATURA WFR MOLDABLE CVX 13/45MM (1X10)US	8A04652
404592	1156500	NATURA WFR MOLDABLE CVX 13/45MM (1X10)US	8A06020
404592	1175793	C/HESIVE NAT WFR MLD CVX 13/22MM(1X10)JP	8A06807
404592	1174445	C/HESIVE NAT MLDC CVX13/22MM(1X10) GB/NL	8A06808
404592	1175760	C/HESIVE NAT WFR MLD CVX 13/22MM(1X10)DE	8A06810
404592	1156500	NATURA WFR MOLDABLE CVX 13/45MM (1X10)US	8C00510
404592	1156500	NATURA WFR MOLDABLE CVX 13/45MM (1X10)US	8C02727
404592	1156500	NATURA WFR MOLDABLE CVX 13/45MM (1X10)US	8C04760
404592	1174445	C/HESIVE NAT MLDC CVX13/22MM(1X10) GB/NL	8C04761
404592	1156500	NATURA WFR MOLDABLE CVX 13/45MM (1X10)US	8D00831

Product Code	Material	Material Description	Lot Number
404592	1156500	NATURA WFR MOLDABLE CVX 13/45MM (1X10)US	8E03229
404592	1175760	C/HESIVE NAT WFR MLD CVX 13/22MM(1X10)DE	8E03230
404592	1701043	NATURA MLDC CVX 13/22MM (1X10PK) EUR	8E03231
404592	1156500	NATURA WFR MOLDABLE CVX 13/45MM (1X10)US	8E05049
404592	1175793	C/HESIVE NAT WFR MLD CVX 13/22MM(1X10)JP	8G02863
404593	1156501	NATURA WFR MOLDABLE CVX 22/45MM (1X10)US	7D00656
404593	1156501	NATURA WFR MOLDABLE CVX 22/45MM (1X10)US	7D04454
404593	1701044	NATURA WFR MLD CVX 22/33MM (1X10PK) EUR	7E02404
404593	1175761	C/HESIVE NAT WFR MLD CVX 22/33MM(1X10)DE	7E02405
404593	1175761	C/HESIVE NAT WFR MLD CVX 22/33MM(1X10)DE	7F01769
404593	1156501	NATURA WFR MOLDABLE CVX 22/45MM (1X10)US	7G02078
404593	1174446	C/HESIVE NAT MLDC CVX22/33MM(1X10) GB/NL	7G02852
404593	1175761	C/HESIVE NAT WFR MLD CVX 22/33MM(1X10)DE	7G02853
404593	1175794	C/HESIVE NAT WFR MLD CVX 22/33MM(1X10)JP	7G02854
404593	1701044	NATURA WFR MLD CVX 22/33MM (1X10PK) EUR	7G04435
404593	1156501	NATURA WFR MOLDABLE CVX 22/45MM (1X10)US	7G04441
404593	1156501	NATURA WFR MOLDABLE CVX 22/45MM (1X10)US	7H02557
404593	1156501	NATURA WFR MOLDABLE CVX 22/45MM (1X10)US	7J01301
404593	1175761	C/HESIVE NAT WFR MLD CVX 22/33MM(1X10)DE	7K00530
404593	1156501	NATURA WFR MOLDABLE CVX 22/45MM (1X10)US	7K04909
404593	1701044	NATURA WFR MLD CVX 22/33MM (1X10PK) EUR	7K04911
404593	1174446	C/HESIVE NAT MLDC CVX22/33MM(1X10) GB/NL	7L00138
404593	1175794	C/HESIVE NAT WFR MLD CVX 22/33MM(1X10)JP	7L00139
404593	1156501	NATURA WFR MOLDABLE CVX 22/45MM (1X10)US	7M01936
404593	1175794	C/HESIVE NAT WFR MLD CVX 22/33MM(1X10)JP	7M02490
404593	1156501	NATURA WFR MOLDABLE CVX 22/45MM (1X10)US	8A00074
404593	1156501	NATURA WFR MOLDABLE CVX 22/45MM (1X10)US	8B00644
404593	1156501	NATURA WFR MOLDABLE CVX 22/45MM (1X10)US	8B02183
404593	1156501	NATURA WFR MOLDABLE CVX 22/45MM (1X10)US	8B03434
404593	1156501	NATURA WFR MOLDABLE CVX 22/45MM (1X10)US	8B04381
404593	1174446	C/HESIVE NAT MLDC CVX22/33MM(1X10) GB/NL	8B05257
404593	1175794	C/HESIVE NAT WFR MLD CVX 22/33MM(1X10)JP	8B05258
404593	1156501	NATURA WFR MOLDABLE CVX 22/45MM (1X10)US	8C00509
404593	1156501	NATURA WFR MOLDABLE CVX 22/45MM (1X10)US	8C01171
404593	1156501	NATURA WFR MOLDABLE CVX 22/45MM (1X10)US	8C02729
404593	1156501	NATURA WFR MOLDABLE CVX 22/45MM (1X10)US	8C03723
404593	1156501	NATURA WFR MOLDABLE CVX 22/45MM (1X10)US	8D02090
404593	1156501	NATURA WFR MOLDABLE CVX 22/45MM (1X10)US	8D03559
404593	1174446	C/HESIVE NAT MLDC CVX22/33MM(1X10) GB/NL	8E01304
404593	1701044	NATURA WFR MLD CVX 22/33MM (1X10PK) EUR	8E01305
404593	1156501	NATURA WFR MOLDABLE CVX 22/45MM (1X10)US	8E01307

Product Code	Material	Material Description	Lot Number
404593	1156501	NATURA WFR MOLDABLE CVX 22/45MM (1X10)US	8E02086
404593	1174446	C/HESIVE NAT MLDC CVX22/33MM(1X10) GB/NL	8E02335
404593	1156501	NATURA WFR MOLDABLE CVX 22/45MM (1X10)US	8E03031
404593	1156501	NATURA WFR MOLDABLE CVX 22/45MM (1X10)US	8E04637
404593	1156501	NATURA WFR MOLDABLE CVX 22/45MM (1X10)US	8F02984
404593	1156501	NATURA WFR MOLDABLE CVX 22/45MM (1X10)US	8F03648
404593	1156501	NATURA WFR MOLDABLE CVX 22/45MM (1X10)US	8F04417
404593	1156501	NATURA WFR MOLDABLE CVX 22/45MM (1X10)US	8G00535
404593	1156501	NATURA WFR MOLDABLE CVX 22/45MM (1X10)US	8G03022
404593	1701044	NATURA WFR MLD CVX 22/33MM (1X10PK) EUR	8G03023
404593	1156501	NATURA WFR MOLDABLE CVX 22/45MM (1X10)US	8G04263
404593	1156501	NATURA WFR MOLDABLE CVX 22/45MM (1X10)US	8H01875
404593	1156501	NATURA WFR MOLDABLE CVX 22/45MM (1X10)US	8H01908
404593	1156501	NATURA WFR MOLDABLE CVX 22/45MM (1X10)US	8H04937
404594	1156502	NATURA WFR MOLDABLE CVX 33/57MM (1X10)US	7E01633
404594	1156502	NATURA WFR MOLDABLE CVX 33/57MM (1X10)US	7E04176
404594	1156502	NATURA WFR MOLDABLE CVX 33/57MM (1X10)US	7E04177
404594	1156502	NATURA WFR MOLDABLE CVX 33/57MM (1X10)US	7E04777
404594	1701045	NATURA MLDC CVX 33/57MM (1X10PK) EUR	7F00003
404594	1156502	NATURA WFR MOLDABLE CVX 33/57MM (1X10)US	7F02063
404594	1156502	NATURA WFR MOLDABLE CVX 33/57MM (1X10)US	7F05062
404594	1175795	C/HESIVE NAT WFR MLD CVX 33/57MM(1X10)JP	7F05064
404594	1701045	NATURA MLDC CVX 33/57MM (1X10PK) EUR	7G00578
404594	1174447	C/HESIVE NAT MLDC CVX33/57MM(1X10) GB/NL	7G05697
404594	1175762	C/HESIVE NAT WFR MLD CVX 33/57MM(1X10)DE	7G05698
404594	1156502	NATURA WFR MOLDABLE CVX 33/57MM (1X10)US	7H04130
404594	1156502	NATURA WFR MOLDABLE CVX 33/57MM (1X10)US	7H04920
404594	1174447	C/HESIVE NAT MLDC CVX33/57MM(1X10) GB/NL	7H04921
404594	1156502	NATURA WFR MOLDABLE CVX 33/57MM (1X10)US	7J02079
404594	1175795	C/HESIVE NAT WFR MLD CVX 33/57MM(1X10)JP	7L01034
404594	1701045	NATURA MLDC CVX 33/57MM (1X10PK) EUR	7L01035
404594	1175762	C/HESIVE NAT WFR MLD CVX 33/57MM(1X10)DE	7L04697
404594	1156502	NATURA WFR MOLDABLE CVX 33/57MM (1X10)US	7L04700
404594	1156502	NATURA WFR MOLDABLE CVX 33/57MM (1X10)US	7L05271
404594	1156502	NATURA WFR MOLDABLE CVX 33/57MM (1X10)US	7M00428
404594	1156502	NATURA WFR MOLDABLE CVX 33/57MM (1X10)US	7M01598
404594	1156502	NATURA WFR MOLDABLE CVX 33/57MM (1X10)US	7M01995
404594	1156502	NATURA WFR MOLDABLE CVX 33/57MM (1X10)US	7M03355
404594	1701045	NATURA MLDC CVX 33/57MM (1X10PK) EUR	7M03432
404594	1175795	C/HESIVE NAT WFR MLD CVX 33/57MM(1X10)JP	7M03434
404594	1156502	NATURA WFR MOLDABLE CVX 33/57MM (1X10)US	8A00740

Product Code	Material	Material Description	Lot Number
404594	1175762	C/HESIVE NAT WFR MLD CVX 33/57MM(1X10)DE	8A00741
404594	1174447	C/HESIVE NAT MLDC CVX33/57MM(1X10) GB/NL	8A00742
404594	1701045	NATURA MLDC CVX 33/57MM (1X10PK) EUR	8B03275
404594	1156502	NATURA WFR MOLDABLE CVX 33/57MM (1X10)US	8D00463
404594	1174447	C/HESIVE NAT MLDC CVX33/57MM(1X10) GB/NL	8D00595
404594	1156502	NATURA WFR MOLDABLE CVX 33/57MM (1X10)US	8D02102
404594	1156502	NATURA WFR MOLDABLE CVX 33/57MM (1X10)US	8D02814
404594	1156502	NATURA WFR MOLDABLE CVX 33/57MM (1X10)US	8D03557
404594	1175795	C/HESIVE NAT WFR MLD CVX 33/57MM(1X10)JP	8D05072
404594	1701045	NATURA MLDC CVX 33/57MM (1X10PK) EUR	8D05074
404594	1156502	NATURA WFR MOLDABLE CVX 33/57MM (1X10)US	8E00637
404594	1156502	NATURA WFR MOLDABLE CVX 33/57MM (1X10)US	8E01067
404594	1156502	NATURA WFR MOLDABLE CVX 33/57MM (1X10)US	8E01627
404594	1174447	C/HESIVE NAT MLDC CVX33/57MM(1X10) GB/NL	8E03028
404594	1156502	NATURA WFR MOLDABLE CVX 33/57MM (1X10)US	8E03029
404594	1175795	C/HESIVE NAT WFR MLD CVX 33/57MM(1X10)JP	8F03571
404594	1156502	NATURA WFR MOLDABLE CVX 33/57MM (1X10)US	8F03634
404594	1156502	NATURA WFR MOLDABLE CVX 33/57MM (1X10)US	8F04925
404594	1701045	NATURA MLDC CVX 33/57MM (1X10PK) EUR	8G03296
404594	1156502	NATURA WFR MOLDABLE CVX 33/57MM (1X10)US	8G03730
404594	1174447	C/HESIVE NAT MLDC CVX33/57MM(1X10) GB/NL	8G03932
404594	1156502	NATURA WFR MOLDABLE CVX 33/57MM (1X10)US	8H01554
404594	1175795	C/HESIVE NAT WFR MLD CVX 33/57MM(1X10)JP	8H02437
404594	1156502	NATURA WFR MOLDABLE CVX 33/57MM (1X10)US	8H02438
404594	1156502	NATURA WFR MOLDABLE CVX 33/57MM (1X10)US	8H04517
404594	1701045	NATURA MLDC CVX 33/57MM (1X10PK) EUR	8J00972
411449	1189481	C/HESIVE S92 MOLDABLE CVX13/45MM1X10PKFR	7D05763
411449	1191965	SURFIT PLUS MOLDABLE CVX13/45MM1X10PKINT	7D05766
411449	1189481	C/HESIVE S92 MOLDABLE CVX13/45MM1X10PKFR	7J01296
411449	1191965	SURFIT PLUS MOLDABLE CVX13/45MM1X10PKINT	7K03430
411449	1189481	C/HESIVE S92 MOLDABLE CVX13/45MM1X10PKFR	8E01630
411450	1189482	C/HESIVE S92 MOLDABLE CVX22/45MM1X10PKFR	7D02980
411450	1191966	SURFIT PLUS MOLDABLE CVX22/45MM1X10PKINT	7D03850
411450	1189482	C/HESIVE S92 MOLDABLE CVX22/45MM1X10PKFR	7E00892
411450	1191966	SURFIT PLUS MOLDABLE CVX22/45MM1X10PKINT	7G01988
411450	1189482	C/HESIVE S92 MOLDABLE CVX22/45MM1X10PKFR	7L01056
411450	1191966	SURFIT PLUS MOLDABLE CVX22/45MM1X10PKINT	7L07677
411450	1191966	SURFIT PLUS MOLDABLE CVX22/45MM1X10PKINT	8J02466
411451	1189483	C/HESIVE S92 MOLDABLE CVX33/57MM1X10PKFR	7G03101
411451	1191968	SURFIT S92MOLDABLE CVX33/57MM1X10PKINT	7G03115
411451	1191968	SURFIT S92MOLDABLE CVX33/57MM1X10PKINT	7K01486

Product Code	Material	Material Description	Lot Number
411451	1189483	C/HESIVE S92 MOLDABLE CVX33/57MM1X10PKFR	7K01698
411451	1189483	C/HESIVE S92 MOLDABLE CVX33/57MM1X10PKFR	7L06598
411451	1191968	SURFIT S92MOLDABLE CVX33/57MM1X10PKINT	8A00745
411451	1191968	SURFIT S92MOLDABLE CVX33/57MM1X10PKINT	8E03026
411451	1189483	C/HESIVE S92 MOLDABLE CVX33/57MM1X10PKFR	8F00894
411451	1191968	SURFIT S92MOLDABLE CVX33/57MM1X10PKINT	8F00895
411451	1191968	SURFIT S92MOLDABLE CVX33/57MM1X10PKINT	8F04686
411451	1191968	SURFIT S92MOLDABLE CVX33/57MM1X10PKINT	8J01450
411451	1189483	C/HESIVE S92 MOLDABLE CVX33/57MM1X10PKFR	8J01492
411452	1187547	C/HESIVE S92 MOLDABLE CVX13/45MM1X5PK IT	7D05767
411452	1187547	C/HESIVE S92 MOLDABLE CVX13/45MM1X5PK IT	7J01297
411452	1190777	C/HESIVE S92 MOLDABLE CVX13/45MM1X5PK BE	7K03428
411452	1187547	C/HESIVE S92 MOLDABLE CVX13/45MM1X5PK IT	7K03432
411452	1187541	C/HESIVE S92 MOLDABLE CVX13/45MM1X5PKSCN	8E01628
411452	1190777	C/HESIVE S92 MOLDABLE CVX13/45MM1X5PK BE	8E01631
411453	1187548	C/HESIVE S92 MOLDABLE CVX22/45MM1X5PK IT	7D02278
411453	1187542	C/HESIVE S92 MOLDABLE CVX22/45MM1X5PKSCN	7D03851
411453	1190778	C/HESIVE S92 MOLDABLE CVX22/45MM1X5PK BE	7G01703
411453	1187548	C/HESIVE S92 MOLDABLE CVX22/45MM1X5PK IT	7G01705
411453	1187542	C/HESIVE S92 MOLDABLE CVX22/45MM1X5PKSCN	7G01989
411453	1187542	C/HESIVE S92 MOLDABLE CVX22/45MM1X5PKSCN	7L01053
411453	1187548	C/HESIVE S92 MOLDABLE CVX22/45MM1X5PK IT	7L01054
411453	1190778	C/HESIVE S92 MOLDABLE CVX22/45MM1X5PK BE	7L07676
411453	1187542	C/HESIVE S92 MOLDABLE CVX22/45MM1X5PKSCN	7L07678
411453	1187548	C/HESIVE S92 MOLDABLE CVX22/45MM1X5PK IT	7L07679
411453	1187545	C/HESIVE S92 MOLDABLE CVX22/45MM1X5PKCEU	7L07681
411453	1187542	C/HESIVE S92 MOLDABLE CVX22/45MM1X5PKSCN	8E00276
411453	1190778	C/HESIVE S92 MOLDABLE CVX22/45MM1X5PK BE	8E00277
411453	1187545	C/HESIVE S92 MOLDABLE CVX22/45MM1X5PKCEU	8J02468
411454	1187543	C/HESIVE S92 MOLDABLE CVX33/57MM1X5PKSCN	7G03112
411454	1187549	C/HESIVE S92 MOLDABLE CVX33/57MM1X5PK IT	7G03113
411454	1190779	C/HESIVE S92 MOLDABLE CVX33/57MM1X5PK BE	7G03114
411454	1187549	C/HESIVE S92 MOLDABLE CVX33/57MM1X5PK IT	7J00624
411454	1190779	C/HESIVE S92 MOLDABLE CVX33/57MM1X5PK BE	7J00626
411454	1187543	C/HESIVE S92 MOLDABLE CVX33/57MM1X5PKSCN	7J00627
411454	1187546	C/HESIVE S92 MOLDABLE CVX33/57MM1X5PKCEU	7J02059
411454	1187546	C/HESIVE S92 MOLDABLE CVX33/57MM1X5PKCEU	7K01699
411454	1190779	C/HESIVE S92 MOLDABLE CVX33/57MM1X5PK BE	7K01700
411454	1187549	C/HESIVE S92 MOLDABLE CVX33/57MM1X5PK IT	7L06595
411454	1187543	C/HESIVE S92 MOLDABLE CVX33/57MM1X5PKSCN	7L06597
411454	1190779	C/HESIVE S92 MOLDABLE CVX33/57MM1X5PK BE	7L06599

Product Code	Material	Material Description	Lot Number
411454	1187546	C/HESIVE S92 MOLDABLE CVX33/57MM1X5PKCEU	8A00744
411454	1190779	C/HESIVE S92 MOLDABLE CVX33/57MM1X5PK BE	8E03024
411454	1187546	C/HESIVE S92 MOLDABLE CVX33/57MM1X5PKCEU	8E03025
411454	1187543	C/HESIVE S92 MOLDABLE CVX33/57MM1X5PKSCN	8F00891
411454	1187549	C/HESIVE S92 MOLDABLE CVX33/57MM1X5PK IT	8F00893
411454	1187546	C/HESIVE S92 MOLDABLE CVX33/57MM1X5PKCEU	8F00900
411454	1187546	C/HESIVE S92 MOLDABLE CVX33/57MM1X5PKCEU	8F04688
411454	1187549	C/HESIVE S92 MOLDABLE CVX33/57MM1X5PK IT	8J01449
411454	1190779	C/HESIVE S92 MOLDABLE CVX33/57MM1X5PK BE	8J01494
413177	1161266	SURFIT NATURA WFR CVX 13/45MM (1X10PK)US	7C05884
413177	1701082	NATURA WFR CVX 13/45MM (1X10PK) EUR	7L03068
413177	1161266	SURFIT NATURA WFR CVX 13/45MM (1X10PK)US	7L03069
413177	1161266	SURFIT NATURA WFR CVX 13/45MM (1X10PK)US	8C03530
413178	1161267	SURFIT NATURA WFR CVX 16/45MM (1X10PK)US	7C05886
413178	1701083	NATURA WFR CVX 16/45MM (1X10PK) EUR	7E02270
413178	1161267	SURFIT NATURA WFR CVX 16/45MM (1X10PK)US	7E02314
413178	1161267	SURFIT NATURA WFR CVX 16/45MM (1X10PK)US	7F04002
413178	1161267	SURFIT NATURA WFR CVX 16/45MM (1X10PK)US	7M03586
413178	1701083	NATURA WFR CVX 16/45MM (1X10PK) EUR	8A00726
413178	1161267	SURFIT NATURA WFR CVX 16/45MM (1X10PK)US	8A00727
413178	1161267	SURFIT NATURA WFR CVX 16/45MM (1X10PK)US	8B03741
413178	1161267	SURFIT NATURA WFR CVX 16/45MM (1X10PK)US	8F00684
413178	1161267	SURFIT NATURA WFR CVX 16/45MM (1X10PK)US	8H04702
413179	1161268	SURFIT NATURA WFR CVX 19/45MM (1X10PK)US	7D01169
413179	1161268	SURFIT NATURA WFR CVX 19/45MM (1X10PK)US	8A00731
413179	1161268	SURFIT NATURA WFR CVX 19/45MM (1X10PK)US	8E04641
413180	1161269	SURFIT NATURA WFR CVX 22/45MM (1X10PK)US	7C01423
413180	1161269	SURFIT NATURA WFR CVX 22/45MM (1X10PK)US	7C02570
413180	1161269	SURFIT NATURA WFR CVX 22/45MM (1X10PK)US	7J01730
413180	1701085	NATURA WFR CVX 22/45MM (1X10PK) EUR	7K00521
413180	1161269	SURFIT NATURA WFR CVX 22/45MM (1X10PK)US	7K00522
413180	1701085	NATURA WFR CVX 22/45MM (1X10PK) EUR	7K05907
413180	1161269	SURFIT NATURA WFR CVX 22/45MM (1X10PK)US	7K05910
413180	1161269	SURFIT NATURA WFR CVX 22/45MM (1X10PK)US	7M02816
413180	1161269	SURFIT NATURA WFR CVX 22/45MM (1X10PK)US	8C05312
413180	1161269	SURFIT NATURA WFR CVX 22/45MM (1X10PK)US	8E01680
413180	1161269	SURFIT NATURA WFR CVX 22/45MM (1X10PK)US	8F00696
413180	1701085	NATURA WFR CVX 22/45MM (1X10PK) EUR	8F00698
413180	1161269	SURFIT NATURA WFR CVX 22/45MM (1X10PK)US	8F00901
413181	1161270	SURFIT NATURA WFR CVX 25/45MM (1X10PK)US	7B01217
413181	1161270	SURFIT NATURA WFR CVX 25/45MM (1X10PK)US	7B04035

Product Code	Material	Material Description	Lot Number
413181	1161270	SURFIT NATURA WFR CVX 25/45MM (1X10PK)US	7C01401
413181	1161270	SURFIT NATURA WFR CVX 25/45MM (1X10PK)US	7C04233
413181	1701086	NATURA WFR CVX 25/45MM (1X10PK) EUR	7E01710
413181	1161270	SURFIT NATURA WFR CVX 25/45MM (1X10PK)US	7E01772
413181	1161270	SURFIT NATURA WFR CVX 25/45MM (1X10PK)US	7E02668
413181	1161270	SURFIT NATURA WFR CVX 25/45MM (1X10PK)US	7E04174
413181	1161270	SURFIT NATURA WFR CVX 25/45MM (1X10PK)US	7F01780
413181	1701086	NATURA WFR CVX 25/45MM (1X10PK) EUR	7G00166
413181	1161270	SURFIT NATURA WFR CVX 25/45MM (1X10PK)US	7G04454
413181	1161270	SURFIT NATURA WFR CVX 25/45MM (1X10PK)US	7K03820
413181	1161270	SURFIT NATURA WFR CVX 25/45MM (1X10PK)US	7M03980
413181	1701086	NATURA WFR CVX 25/45MM (1X10PK) EUR	7M03981
413181	1161270	SURFIT NATURA WFR CVX 25/45MM (1X10PK)US	8B00893
413181	1161270	SURFIT NATURA WFR CVX 25/45MM (1X10PK)US	8E03844
413181	1701086	NATURA WFR CVX 25/45MM (1X10PK) EUR	8G00927
413181	1161270	SURFIT NATURA WFR CVX 25/45MM (1X10PK)US	8G00928
413181	1161270	SURFIT NATURA WFR CVX 25/45MM (1X10PK)US	8H01558
413181	1161270	SURFIT NATURA WFR CVX 25/45MM (1X10PK)US	8H03658
413181	1161270	SURFIT NATURA WFR CVX 25/45MM (1X10PK)US	8H04303
413181	1701086	NATURA WFR CVX 25/45MM (1X10PK) EUR	8H04356
413182	1161271	SURFIT NATURA WFR CVX 28/45MM (1X10PK)US	7C05415
413182	1161271	SURFIT NATURA WFR CVX 28/45MM (1X10PK)US	7D01141
413182	1161271	SURFIT NATURA WFR CVX 28/45MM (1X10PK)US	7F03152
413182	1701087	NATURA WFR CVX 28/45MM (1X10PK) EUR	7F03155
413182	1161271	SURFIT NATURA WFR CVX 28/45MM (1X10PK)US	7H04051
413182	1701087	NATURA WFR CVX 28/45MM (1X10PK) EUR	7L01044
413182	1161271	SURFIT NATURA WFR CVX 28/45MM (1X10PK)US	7L06593
413182	1701087	NATURA WFR CVX 28/45MM (1X10PK) EUR	8B03860
413182	1701087	NATURA WFR CVX 28/45MM (1X10PK) EUR	8F01047
413182	1161271	SURFIT NATURA WFR CVX 28/45MM (1X10PK)US	8F01050
413182	1161271	SURFIT NATURA WFR CVX 28/45MM (1X10PK)US	8J00537
413183	1161272	SURFIT NATURA WFR CVX 32/45MM (1X10PK)US	7B04470
413183	1161272	SURFIT NATURA WFR CVX 32/45MM (1X10PK)US	7F01265
413183	1161272	SURFIT NATURA WFR CVX 32/45MM (1X10PK)US	7H02188
413183	1161272	SURFIT NATURA WFR CVX 32/45MM (1X10PK)US	8A02999
413183	1161272	SURFIT NATURA WFR CVX 32/45MM (1X10PK)US	8B03435
413183	1161272	SURFIT NATURA WFR CVX 32/45MM (1X10PK)US	8B05002
413183	1161272	SURFIT NATURA WFR CVX 32/45MM (1X10PK)US	8C03533
413183	1161272	SURFIT NATURA WFR CVX 32/45MM (1X10PK)US	8D05105
413183	1161272	SURFIT NATURA WFR CVX 32/45MM (1X10PK)US	8F03646
413184	1161273	SURFIT NATURA WFR CVX 35/45MM (1X10PK)US	7C05888

Product Code	Material	Material Description	Lot Number
413184	1161273	SURFIT NATURA WFR CVX 35/45MM (1X10PK)US	7C06118
413184	1161273	SURFIT NATURA WFR CVX 35/45MM (1X10PK)US	7F01763
413184	1701089	NATURA WFR CVX 35/45MM (1X10PK) EUR	7J04970
413184	1161273	SURFIT NATURA WFR CVX 35/45MM (1X10PK)US	7J05042
413184	1161273	SURFIT NATURA WFR CVX 35/45MM (1X10PK)US	8G03497
413184	1701089	NATURA WFR CVX 35/45MM (1X10PK) EUR	8G03499
413185	1161274	SURFIT NATURA WFR CVX 38/57MM (1X10PK)US	7C05891
413185	1161274	SURFIT NATURA WFR CVX 38/57MM (1X10PK)US	7C06116
413185	1161274	SURFIT NATURA WFR CVX 38/57MM (1X10PK)US	7D04719
413185	1161274	SURFIT NATURA WFR CVX 38/57MM (1X10PK)US	7F05916
413185	1701090	NATURA WFR CVX 38/57MM (1X10PK) EUR	7H06314
413185	1161274	SURFIT NATURA WFR CVX 38/57MM (1X10PK)US	7H06315
413185	1701090	NATURA WFR CVX 38/57MM (1X10PK) EUR	7J04230
413185	1161274	SURFIT NATURA WFR CVX 38/57MM (1X10PK)US	8A04962
413185	1161274	SURFIT NATURA WFR CVX 38/57MM (1X10PK)US	8A06014
413185	1161274	SURFIT NATURA WFR CVX 38/57MM (1X10PK)US	8C03676
413185	1161274	SURFIT NATURA WFR CVX 38/57MM (1X10PK)US	8E01297
413185	1161274	SURFIT NATURA WFR CVX 38/57MM (1X10PK)US	8F04684
413185	1161274	SURFIT NATURA WFR CVX 38/57MM (1X10PK)US	8G02137
413186	1161275	SURFIT NATURA WFR CVX 41/57MM (1X10PK)US	7F05917
413186	1161275	SURFIT NATURA WFR CVX 41/57MM (1X10PK)US	8A01385
413186	1161275	SURFIT NATURA WFR CVX 41/57MM (1X10PK)US	8D00499
413186	1701091	NATURA WFR CVX 41/57MM (1X10PK) EUR	8D05102
413186	1161275	SURFIT NATURA WFR CVX 41/57MM (1X10PK)US	8D05103
413187	1161276	SURFIT NATURA WFR CVX 45/57MM (1X10PK)US	7F00550
413187	1161276	SURFIT NATURA WFR CVX 45/57MM (1X10PK)US	7H04926
413187	1161276	SURFIT NATURA WFR CVX 45/57MM (1X10PK)US	7H05499
413187	1701092	NATURA WFR CVX 45/57MM (1X10PK) EUR	7M00470
413187	1161276	SURFIT NATURA WFR CVX 45/57MM (1X10PK)US	7M00471
413187	1161276	SURFIT NATURA WFR CVX 45/57MM (1X10PK)US	8B02277
413187	1161276	SURFIT NATURA WFR CVX 45/57MM (1X10PK)US	8B03309
413187	1161276	SURFIT NATURA WFR CVX 45/57MM (1X10PK)US	8D02868
413187	1161276	SURFIT NATURA WFR CVX 45/57MM (1X10PK)US	8F01515
413187	1161276	SURFIT NATURA WFR CVX 45/57MM (1X10PK)US	8J02475
413188	1161277	SURFIT NATURA WFR CVX 50/57MM (1X10PK)US	7E02407
413188	1161277	SURFIT NATURA WFR CVX 50/57MM (1X10PK)US	7G01362
413188	1701093	NATURA WFR CVX 50/57MM (1X10PK) EUR	7G01525
413188	1161277	SURFIT NATURA WFR CVX 50/57MM (1X10PK)US	7J00631
413188	1701093	NATURA WFR CVX 50/57MM (1X10PK) EUR	7K07045
413188	1161277	SURFIT NATURA WFR CVX 50/57MM (1X10PK)US	8B01726
413188	1161277	SURFIT NATURA WFR CVX 50/57MM (1X10PK)US	8D00097

Product Code	Material	Material Description	Lot Number
413188	1161277	SURFIT NATURA WFR CVX 50/57MM (1X10PK)US	8H00324
413188	1161277	SURFIT NATURA WFR CVX 50/57MM (1X10PK)US	8J00420
420233	1701098	NATURA WFR MLD CVX 13/45MM (1X5PK) EUR	7J04280
420233	1701098	NATURA WFR MLD CVX 13/45MM (1X5PK) EUR	7J06070
420233	1701098	NATURA WFR MLD CVX 13/45MM (1X5PK) EUR	7M03041
420233	1701098	NATURA WFR MLD CVX 13/45MM (1X5PK) EUR	8C04758
420233	1701098	NATURA WFR MLD CVX 13/45MM (1X5PK) EUR	8E03322
420233	1701098	NATURA WFR MLD CVX 13/45MM (1X5PK) EUR	8E05050
420234	1701099	NATURA WFR MLD CVX 22/45MM (1X5PK) EUR	7F01770
420234	1701099	NATURA WFR MLD CVX 22/45MM (1X5PK) EUR	7G03384
420234	1701099	NATURA WFR MLD CVX 22/45MM (1X5PK) EUR	7G04437
420234	1701099	NATURA WFR MLD CVX 22/45MM (1X5PK) EUR	7K00529
420234	1701099	NATURA WFR MLD CVX 22/45MM (1X5PK) EUR	7L00137
420234	1701099	NATURA WFR MLD CVX 22/45MM (1X5PK) EUR	7M01882
420235	1701100	NATURA WFR MLD CVX 33/57MM (1X5PK) EUR	7F00002
420235	1701100	NATURA WFR MLD CVX 33/57MM (1X5PK) EUR	7F05063
420235	1701100	NATURA WFR MLD CVX 33/57MM (1X5PK) EUR	7J02080
420235	1701100	NATURA WFR MLD CVX 33/57MM (1X5PK) EUR	7K02855
420235	1701100	NATURA WFR MLD CVX 33/57MM (1X5PK) EUR	7L04696
420235	1701100	NATURA WFR MLD CVX 33/57MM (1X5PK) EUR	7M03433
420235	1701100	NATURA WFR MLD CVX 33/57MM (1X5PK) EUR	8A00739
420235	1701100	NATURA WFR MLD CVX 33/57MM (1X5PK) EUR	8B04155
420235	1701100	NATURA WFR MLD CVX 33/57MM (1X5PK) EUR	8C00709
420235	1701100	NATURA WFR MLD CVX 33/57MM (1X5PK) EUR	8D02104
420235	1701100	NATURA WFR MLD CVX 33/57MM (1X5PK) EUR	8D05073
420235	1701100	NATURA WFR MLD CVX 33/57MM (1X5PK) EUR	8F03632
420235	1701100	NATURA WFR MLD CVX 33/57MM (1X5PK) EUR	8G03297
420235	1701100	NATURA WFR MLD CVX 33/57MM (1X5PK) EUR	8J00971