

Month xx, 2020

URGENT FIELD SAFETY NOTICE
Datascope/Getinge IAB Potential Endotoxin Contamination

Linear 7.5Fr 25cc IAB

Linear 7.5Fr 40cc IAB

Linear 7.5Fr 34cc IAB

Sensation 7Fr 34cc IAB

Sensation 7Fr 40cc IAB

Sensation Plus 7.5 Fr 40cc IAB

Sensation Plus 8Fr 50cc IAB

MEGA 7.5Fr 30cc IAB

MEGA 7.5Fr 40cc IAB

MEGA 8Fr 50cc IAB

KIT Part Numbers and LOT numbers	The IAB Kit LOTS that were shipped to you and may potentially contain affected serial numbers are listed in Table 1
Affected IAB Serial Numbers:	Affected serial numbers to screen for are listed in Table 2
IAB Manufacturing Dates:	[February 3, 2017] through [February 21, 2020]
KIT Distribution Dates:	[February 9, 2017] through [May 21, 2020]

PLEASE FORWARD THIS INFORMATION TO ALL CURRENT AND POTENTIAL DATASCOPE/GETINGE IAB USERS WITHIN YOUR HOSPITAL / FACILITY. IF YOU HAVE FURTHER DISTRIBUTED ANY OF THE AFFECTED PRODUCTS, FORWARD THIS INFORMATION TO THE RECIPIENT.

IF YOU HAVE FURTHER DISTRIBUTED ANY OF THE AFFECTED PRODUCTS, FORWARD THIS INFORMATION TO THE RECIPIENT.

Dear Risk Manager,

Datascope/Getinge is initiating a voluntary Recall-Removal involving certain Intra-Aortic Balloon Catheters (IABs) that may not meet the requirement for endotoxin per AAMI ST72. Datascope/Getinge performs functional testing on a small number of units from every lot prior to sterilization and these functionally-tested units may pose an elevated risk of endotoxin contamination compared to normal production IABs. These functionally-tested units can be identified by serial number and represent less than 1% of total IABs distributed in this timeframe.

Identification of the issue

The issue is not detectable by the end user. Datascope/Getinge has determined that the affected IABs may not meet the requirement for endotoxin per industry Standard AAMI ST72. All affected IABs can be identified using the Product Labeling and serial number on the Y-fitting of the IAB, quarantined and returned to Datascope/Getinge.

Risk to Health

Physicians would not be able to detect elevated levels of endotoxin before a device is used on a patient. Moreover, since patients in need of this therapy are at higher risk for a systemic inflammatory response including fever, the sole presence of such signs and symptoms usually would not allow identification of a pyrogenic device as root cause. Several drugs administered during and after


procedures or during intensive care stay may even limit the extent of activation of endotoxin-associated cascades.

To date, Getinge/Maquet has not received any complaints or adverse events regarding this issue.


Actions to be taken:

Our records indicate that you received potentially affected IAB Serial Numbers identified in [Table 2]. Please complete the steps below:

- Monitor patients for pyrogens reaction/humoral immune response/ coagulation and complement cascades/ inflammation. Monitor and treat any signs of inflammation according to your facility's protocols and clinical judgment.
- Please examine your inventory immediately, remove and quarantine any unexpired affected IAB following the steps below:
 - Identify any unexpired IAB Kits referencing the *part and lot numbers* listed in [Table 1].
The Kit part and lot numbers can be found on the Outer Shelf Carton as circled in RED in the *example* below. REF denotes part number and LOT denotes lot/batch number.



- Using steps below, Identify any IAB Y-fittings with *serial numbers* listed in [Table 2].
 - Remove tape on **one** side of the Shelf Carton.
 - Hinge the box open, keeping the other side taped.
 - Do not remove any items from the carton. Instead, Lift the Accessories up and locate the Y-Fitting through the clear Mylar overwrap.
 - Read the serial number off the IAB Y-Fitting. (See below for example. IAB SN will appear is the area circled in RED.)



- If you find an IAB with a serial number listed in [**Table 2**], immediately place the entire Kit in quarantine for return and replacement or credit.
- If the affected serial number is not found, close the carton and re-seal. Cut out and affix one of the enclosed labels to designate the kit has been checked for affected recall Serial Numbers and may be placed back in stock. Color copies may be made and affixed with tape.
- **If you have unexpired affected product to return, please contact Getinge Customer Service at X-XXX-XXX-XXXX insert local SSU information to request a return authorization number (RMA) and shipping instructions.**
- **Whether you have affected product or not, please complete and return the form (including all pages) to acknowledge this recall by fax or email to insert SSU information here.**

This voluntary recall only affects specific IAB serial numbers manufactured between February 3, 2017 and February 21, 2020. No other products are affected by this voluntary recall.

Sincerely,

Insert SSU information here

URGENT FIELD SAFETY NOTICE
Datascope/Getinge IAB Endotoxin Issue Recall
Customer Response Form

Linear 7.5Fr 25cc IAB
 Linear 7.5Fr 40cc IAB
 Linear 7.5Fr 34cc IAB

Sensation 7Fr 34cc IAB
 Sensation 7Fr 40cc IAB
 Sensation Plus 7.5 Fr 40cc IAB
 Sensation Plus 8Fr 50cc IAB

MEGA 7.5Fr 30cc IAB
 MEGA 7.5Fr 40cc IAB
 MEGA 8Fr 50cc IAB

Facility Name

Address

Address

☐

If you have **NO AFFECTED UNEXPIRED**
SERIAL NUMBERS at your facility,

please check box and complete this form!

If you have any un-expired affected SERIAL NUMBERS for return and replacement or credit, please complete the form below by circling any serial numbers if found. If you have affected product to return, please contact Getinge Customer Service at [enter SSU information here](#).

Finished Goods Kit Part Number	Finished Goods Batch	Serial Number	Quantity Returned	Getinge RMA No.

ACKNOWLEDGMENT:

By signing below, I acknowledge that I have read and understand this Medical Device Recall Notice for the recalled Maquet/Getinge Intra-Aortic Balloon Catheter Endotoxin recall. I confirm that all IAB users at this facility have been notified accordingly.

Confirm the Facility Name you are responding for	
Address	
Signature	Date
Name	Phone
Title	Department

FAX BACK TO: [SSU info](#) OR EMAIL TO: [SSU info](#)

Kit Description	Kit PNs	Kit LOT numbers
MEGA 7.5Fr 30cc IAB	0684-00-0294-01	3000047044
MEGA 7.5Fr 30cc IAB	0684-00-0294-01	3000049910
MEGA 7.5Fr 30cc IAB	0684-00-0294-01	3000055308
MEGA 7.5Fr 30cc IAB	0684-00-0294-01	3000059274
MEGA 7.5Fr 30cc IAB	0684-00-0294-01	3000064252
MEGA 7.5Fr 30cc IAB	0684-00-0294-01	3000067462
MEGA 7.5Fr 30cc IAB	0684-00-0294-01	3000071220
MEGA 7.5Fr 30cc IAB	0684-00-0294-01	3000077343
MEGA 7.5Fr 30cc IAB	0684-00-0294-01	3000084725
MEGA 7.5Fr 30cc IAB	0684-00-0294-01	3000086853
MEGA 7.5Fr 30cc IAB	0684-00-0294-01	3000087475
MEGA 7.5Fr 30cc IAB	0684-00-0294-01	3000087971
MEGA 7.5Fr 30cc IAB	0684-00-0294-01	3000096716
MEGA 7.5Fr 30cc IAB	0684-00-0294-01	3000097647
MEGA 7.5Fr 30cc IAB	0684-00-0294-01	3000100873
MEGA 7.5Fr 30cc IAB	0684-00-0294-01	3000102114
MEGA 7.5Fr 30cc IAB	0684-00-0294-01	3000109816
MEGA 7.5Fr 30cc IAB	0684-00-0294-01U	3000047042
MEGA 7.5Fr 30cc IAB	0684-00-0294-01U	3000047931
MEGA 7.5Fr 30cc IAB	0684-00-0294-01U	3000049913
MEGA 7.5Fr 30cc IAB	0684-00-0294-01U	3000055307
MEGA 7.5Fr 30cc IAB	0684-00-0294-01U	3000056433
MEGA 7.5Fr 30cc IAB	0684-00-0294-01U	3000059275
MEGA 7.5Fr 30cc IAB	0684-00-0294-01U	3000060230
MEGA 7.5Fr 30cc IAB	0684-00-0294-01U	3000062405
MEGA 7.5Fr 30cc IAB	0684-00-0294-01U	3000064251
MEGA 7.5Fr 30cc IAB	0684-00-0294-01U	3000067461
MEGA 7.5Fr 30cc IAB	0684-00-0294-01U	3000067855
MEGA 7.5Fr 30cc IAB	0684-00-0294-01U	3000071981
MEGA 7.5Fr 30cc IAB	0684-00-0294-01U	3000071219
MEGA 7.5Fr 30cc IAB	0684-00-0294-01U	3000071982
MEGA 7.5Fr 30cc IAB	0684-00-0294-01U	3000077342
MEGA 7.5Fr 30cc IAB	0684-00-0294-01U	3000080542
MEGA 7.5Fr 30cc IAB	0684-00-0294-01U	3000081080
MEGA 7.5Fr 30cc IAB	0684-00-0294-01U	3000084724
MEGA 7.5Fr 30cc IAB	0684-00-0294-01U	3000086852
MEGA 7.5Fr 30cc IAB	0684-00-0294-01U	3000087972
MEGA 7.5Fr 30cc IAB	0684-00-0294-01U	3000096715
MEGA 7.5Fr 30cc IAB	0684-00-0294-01U	3000090374
MEGA 7.5Fr 30cc IAB	0684-00-0294-01U	3000092556
MEGA 7.5Fr 30cc IAB	0684-00-0294-01U	3000097646
MEGA 7.5Fr 30cc IAB	0684-00-0294-01U	3000102113
MEGA 7.5Fr 30cc IAB	0684-00-0294-01U	3000106149
MEGA 7.5Fr 30cc IAB	0684-00-0294-01U	3000112151
MEGA 7.5Fr 30cc IAB	0684-00-0294-01U	3000114260
MEGA 7.5Fr 30cc IAB	0684-00-0294-02	3000047043
MEGA 7.5Fr 30cc IAB	0684-00-0294-02	3000049911
MEGA 7.5Fr 30cc IAB	0684-00-0294-02	3000064746
MEGA 7.5Fr 30cc IAB	0684-00-0294-02	3000067463
MEGA 7.5Fr 30cc IAB	0684-00-0294-02	3000084473

MEGA 7.5Fr 30cc IAB	0684-00-0294-02	3000097648
MEGA 7.5Fr 30cc IAB	0684-00-0294-02	3000112153
MEGA 7.5Fr 30cc IAB	0684-00-0294-02U	3000047041
MEGA 7.5Fr 30cc IAB	0684-00-0294-02U	3000047932
MEGA 7.5Fr 30cc IAB	0684-00-0294-02U	3000049912
MEGA 7.5Fr 30cc IAB	0684-00-0294-02U	3000067856
MEGA 7.5Fr 30cc IAB	0684-00-0294-02U	3000071983
MEGA 7.5Fr 30cc IAB	0684-00-0294-02U	3000087476
MEGA 7.5Fr 30cc IAB	0684-00-0294-02U	3000096793
MEGA 7.5Fr 30cc IAB	0684-00-0294-02U	3000097649
MEGA 7.5Fr 30cc IAB	0684-00-0294-02U	3000112152
MEGA 7.5Fr 40cc IAB	0684-00-0295-01	3000047433
MEGA 7.5Fr 40cc IAB	0684-00-0295-01	3000048714
MEGA 7.5Fr 40cc IAB	0684-00-0295-01	3000049793
MEGA 7.5Fr 40cc IAB	0684-00-0295-01	3000053296
MEGA 7.5Fr 40cc IAB	0684-00-0295-01	3000054596
MEGA 7.5Fr 40cc IAB	0684-00-0295-01	3000055303
MEGA 7.5Fr 40cc IAB	0684-00-0295-01	3000056631
MEGA 7.5Fr 40cc IAB	0684-00-0295-01	3000060227
MEGA 7.5Fr 40cc IAB	0684-00-0295-01	3000064257
MEGA 7.5Fr 40cc IAB	0684-00-0295-01	3000070321
MEGA 7.5Fr 40cc IAB	0684-00-0295-01	3000071835
MEGA 7.5Fr 40cc IAB	0684-00-0295-01	3000072581
MEGA 7.5Fr 40cc IAB	0684-00-0295-01	3000072582
MEGA 7.5Fr 40cc IAB	0684-00-0295-01	3000076925
MEGA 7.5Fr 40cc IAB	0684-00-0295-01	3000078405
MEGA 7.5Fr 40cc IAB	0684-00-0295-01	3000079518
MEGA 7.5Fr 40cc IAB	0684-00-0295-01	3000080107
MEGA 7.5Fr 40cc IAB	0684-00-0295-01	3000083511
MEGA 7.5Fr 40cc IAB	0684-00-0295-01	3000084721
MEGA 7.5Fr 40cc IAB	0684-00-0295-01	3000087689
MEGA 7.5Fr 40cc IAB	0684-00-0295-01	3000092131
MEGA 7.5Fr 40cc IAB	0684-00-0295-01	3000092215
MEGA 7.5Fr 40cc IAB	0684-00-0295-01	3000095792
MEGA 7.5Fr 40cc IAB	0684-00-0295-01	3000096364
MEGA 7.5Fr 40cc IAB	0684-00-0295-01	3000097534
MEGA 7.5Fr 40cc IAB	0684-00-0295-01	3000100872
MEGA 7.5Fr 40cc IAB	0684-00-0295-01	3000102116
MEGA 7.5Fr 40cc IAB	0684-00-0295-01	3000106750
MEGA 7.5Fr 40cc IAB	0684-00-0295-01	3000107623
MEGA 7.5Fr 40cc IAB	0684-00-0295-01	3000112279
MEGA 7.5Fr 40cc IAB	0684-00-0295-01U	3000047432
MEGA 7.5Fr 40cc IAB	0684-00-0295-01U	3000047933
MEGA 7.5Fr 40cc IAB	0684-00-0295-01U	3000049794
MEGA 7.5Fr 40cc IAB	0684-00-0295-01U	3000052255
MEGA 7.5Fr 40cc IAB	0684-00-0295-01U	3000055306
MEGA 7.5Fr 40cc IAB	0684-00-0295-01U	3000056632
MEGA 7.5Fr 40cc IAB	0684-00-0295-01U	3000061197
MEGA 7.5Fr 40cc IAB	0684-00-0295-01U	3000064258
MEGA 7.5Fr 40cc IAB	0684-00-0295-01U	3000065635
MEGA 7.5Fr 40cc IAB	0684-00-0295-01U	3000066086

MEGA 7.5Fr 40cc IAB	0684-00-0295-01U	3000066265
MEGA 7.5Fr 40cc IAB	0684-00-0295-01U	3000067264
MEGA 7.5Fr 40cc IAB	0684-00-0295-01U	3000067841
MEGA 7.5Fr 40cc IAB	0684-00-0295-01U	3000068601
MEGA 7.5Fr 40cc IAB	0684-00-0295-01U	3000068602
MEGA 7.5Fr 40cc IAB	0684-00-0295-01U	3000074729
MEGA 7.5Fr 40cc IAB	0684-00-0295-01U	3000074850
MEGA 7.5Fr 40cc IAB	0684-00-0295-01U	3000079519
MEGA 7.5Fr 40cc IAB	0684-00-0295-01U	3000081079
MEGA 7.5Fr 40cc IAB	0684-00-0295-01U	3000080541
MEGA 7.5Fr 40cc IAB	0684-00-0295-01U	3000083512
MEGA 7.5Fr 40cc IAB	0684-00-0295-01U	3000084722
MEGA 7.5Fr 40cc IAB	0684-00-0295-01U	3000087474
MEGA 7.5Fr 40cc IAB	0684-00-0295-01U	3000087973
MEGA 7.5Fr 40cc IAB	0684-00-0295-01U	3000089282
MEGA 7.5Fr 40cc IAB	0684-00-0295-01U	3000092129
MEGA 7.5Fr 40cc IAB	0684-00-0295-01U	3000092130
MEGA 7.5Fr 40cc IAB	0684-00-0295-01U	3000095793
MEGA 7.5Fr 40cc IAB	0684-00-0295-01U	3000097645
MEGA 7.5Fr 40cc IAB	0684-00-0295-01U	3000096762
MEGA 7.5Fr 40cc IAB	0684-00-0295-01U	3000097644
MEGA 7.5Fr 40cc IAB	0684-00-0295-01U	3000097643
MEGA 7.5Fr 40cc IAB	0684-00-0295-01U	3000103157
MEGA 7.5Fr 40cc IAB	0684-00-0295-01U	3000104738
MEGA 7.5Fr 40cc IAB	0684-00-0295-01U	3000106751
MEGA 7.5Fr 40cc IAB	0684-00-0295-01U	3000107624
MEGA 7.5Fr 40cc IAB	0684-00-0295-01U	3000112280
MEGA 7.5Fr 40cc IAB	0684-00-0295-02	3000047434
MEGA 7.5Fr 40cc IAB	0684-00-0295-02	3000047523
MEGA 7.5Fr 40cc IAB	0684-00-0295-02	3000048024
MEGA 7.5Fr 40cc IAB	0684-00-0295-02	3000048025
MEGA 7.5Fr 40cc IAB	0684-00-0295-02	3000049795
MEGA 7.5Fr 40cc IAB	0684-00-0295-02	3000052254
MEGA 7.5Fr 40cc IAB	0684-00-0295-02	3000053295
MEGA 7.5Fr 40cc IAB	0684-00-0295-02	3000054426
MEGA 7.5Fr 40cc IAB	0684-00-0295-02	3000054425
MEGA 7.5Fr 40cc IAB	0684-00-0295-02	3000055305
MEGA 7.5Fr 40cc IAB	0684-00-0295-02	3000059063
MEGA 7.5Fr 40cc IAB	0684-00-0295-02	3000060228
MEGA 7.5Fr 40cc IAB	0684-00-0295-02	3000061341
MEGA 7.5Fr 40cc IAB	0684-00-0295-02	3000064844
MEGA 7.5Fr 40cc IAB	0684-00-0295-02	3000072583
MEGA 7.5Fr 40cc IAB	0684-00-0295-02	3000081622
MEGA 7.5Fr 40cc IAB	0684-00-0295-02	3000082248
MEGA 7.5Fr 40cc IAB	0684-00-0295-02	3000082188
MEGA 7.5Fr 40cc IAB	0684-00-0295-02	3000083990
MEGA 7.5Fr 40cc IAB	0684-00-0295-02	3000084723
MEGA 7.5Fr 40cc IAB	0684-00-0295-02	3000099218
MEGA 7.5Fr 40cc IAB	0684-00-0295-02	3000104736
MEGA 7.5Fr 40cc IAB	0684-00-0295-02	3000104737
MEGA 7.5Fr 40cc IAB	0684-00-0295-02	3000106752

MEGA 7.5Fr 40cc IAB	0684-00-0295-02U	3000047309
MEGA 7.5Fr 40cc IAB	0684-00-0295-02U	3000047934
MEGA 7.5Fr 40cc IAB	0684-00-0295-02U	3000049796
MEGA 7.5Fr 40cc IAB	0684-00-0295-02U	3000066797
MEGA 7.5Fr 40cc IAB	0684-00-0295-02U	3000074731
MEGA 7.5Fr 40cc IAB	0684-00-0295-02U	3000083513
MEGA 7.5Fr 40cc IAB	0684-00-0295-02U	3000079520
MEGA 7.5Fr 40cc IAB	0684-00-0295-02U	3000083514
MEGA 7.5Fr 40cc IAB	0684-00-0295-02U	3000095794
MEGA 7.5Fr 40cc IAB	0684-00-0295-02U	3000096366
MEGA 7.5Fr 40cc IAB	0684-00-0295-02U	3000102115
MEGA 7.5Fr 40cc IAB	0684-00-0295-02U	3000106753
MEGA 7.5Fr 40cc IAB	0684-00-0295-02U	3000112278
MEGA 7.5Fr 40cc IAB	0684-00-0295-05	3000047355
MEGA 7.5Fr 40cc IAB	0684-00-0295-05	3000048228
MEGA 7.5Fr 40cc IAB	0684-00-0295-05	3000050257
MEGA 7.5Fr 40cc IAB	0684-00-0295-05	3000052253
MEGA 7.5Fr 40cc IAB	0684-00-0295-05	3000054424
MEGA 7.5Fr 40cc IAB	0684-00-0295-05	3000055304
MEGA 7.5Fr 40cc IAB	0684-00-0295-05	3000056630
MEGA 7.5Fr 40cc IAB	0684-00-0295-05	3000059062
MEGA 7.5Fr 40cc IAB	0684-00-0295-05	3000060229
MEGA 7.5Fr 40cc IAB	0684-00-0295-05	3000061195
MEGA 7.5Fr 40cc IAB	0684-00-0295-05	3000061196
MEGA 7.5Fr 40cc IAB	0684-00-0295-05	3000064259
MEGA 7.5Fr 40cc IAB	0684-00-0295-05	3000064355
MEGA 7.5Fr 40cc IAB	0684-00-0295-05	3000066796
MEGA 7.5Fr 40cc IAB	0684-00-0295-05	3000069443
MEGA 7.5Fr 40cc IAB	0684-00-0295-05	3000070319
MEGA 7.5Fr 40cc IAB	0684-00-0295-05	3000080729
MEGA 7.5Fr 40cc IAB	0684-00-0295-05	3000077201
MEGA 7.5Fr 40cc IAB	0684-00-0295-05	3000078053
MEGA 7.5Fr 40cc IAB	0684-00-0295-05	3000078919
MEGA 7.5Fr 40cc IAB	0684-00-0295-05	3000080419
MEGA 7.5Fr 40cc IAB	0684-00-0295-05	3000083599
MEGA 7.5Fr 40cc IAB	0684-00-0295-05	3000084691
MEGA 7.5Fr 40cc IAB	0684-00-0295-05	3000085499
MEGA 7.5Fr 40cc IAB	0684-00-0295-05	3000086793
MEGA 7.5Fr 40cc IAB	0684-00-0295-05	3000086794
MEGA 7.5Fr 40cc IAB	0684-00-0295-05	3000089346
MEGA 7.5Fr 40cc IAB	0684-00-0295-05	3000091874
MEGA 7.5Fr 40cc IAB	0684-00-0295-05	3000090851
MEGA 7.5Fr 40cc IAB	0684-00-0295-05	3000091873
MEGA 7.5Fr 40cc IAB	0684-00-0295-05	3000093117
MEGA 7.5Fr 40cc IAB	0684-00-0295-05	3000093118
MEGA 7.5Fr 40cc IAB	0684-00-0295-05	3000093119
MEGA 7.5Fr 40cc IAB	0684-00-0295-05	3000094614
MEGA 7.5Fr 40cc IAB	0684-00-0295-05	3000094613
MEGA 7.5Fr 40cc IAB	0684-00-0295-05	3000095319
MEGA 7.5Fr 40cc IAB	0684-00-0295-05	3000095854
MEGA 7.5Fr 40cc IAB	0684-00-0295-05	3000096761

MEGA 7.5Fr 40cc IAB	0684-00-0295-05	3000097533
MEGA 7.5Fr 40cc IAB	0684-00-0295-05	3000099217
MEGA 7.5Fr 40cc IAB	0684-00-0295-05	3000100871
MEGA 7.5Fr 40cc IAB	0684-00-0295-05	3000102117
MEGA 7.5Fr 40cc IAB	0684-00-0295-05	3000103158
MEGA 7.5Fr 40cc IAB	0684-00-0295-05	3000104739
MEGA 7.5Fr 40cc IAB	0684-00-0295-05	3000106898
MEGA 7.5Fr 40cc IAB	0684-00-0295-05	3000107622
MEGA 7.5Fr 40cc IAB	0684-00-0295-05	3000112281
MEGA 7.5Fr 40cc IAB	0684-00-0295-05	3000114254
MEGA 7.5Fr 40cc IAB	0684-00-0295-07	3000071786
MEGA 7.5Fr 40cc IAB	0684-00-0295-07	3000073203
MEGA 7.5Fr 40cc IAB	0684-00-0295-09	3000096365
MEGA 7.5Fr 40cc IAB	0684-00-0295-10	3000103156
MEGA 8Fr 50cc IAB	0684-00-0296-01	3000046414
MEGA 8Fr 50cc IAB	0684-00-0296-01	3000047431
MEGA 8Fr 50cc IAB	0684-00-0296-01	3000048243
MEGA 8Fr 50cc IAB	0684-00-0296-01	3000049080
MEGA 8Fr 50cc IAB	0684-00-0296-01	3000049569
MEGA 8Fr 50cc IAB	0684-00-0296-01	3000050495
MEGA 8Fr 50cc IAB	0684-00-0296-01	3000051154
MEGA 8Fr 50cc IAB	0684-00-0296-01	3000051751
MEGA 8Fr 50cc IAB	0684-00-0296-01	3000053311
MEGA 8Fr 50cc IAB	0684-00-0296-01	3000054591
MEGA 8Fr 50cc IAB	0684-00-0296-01	3000054592
MEGA 8Fr 50cc IAB	0684-00-0296-01	3000056356
MEGA 8Fr 50cc IAB	0684-00-0296-01	3000056629
MEGA 8Fr 50cc IAB	0684-00-0296-01	3000058277
MEGA 8Fr 50cc IAB	0684-00-0296-01	3000060542
MEGA 8Fr 50cc IAB	0684-00-0296-01	3000064249
MEGA 8Fr 50cc IAB	0684-00-0296-01	3000066800
MEGA 8Fr 50cc IAB	0684-00-0296-01	3000065563
MEGA 8Fr 50cc IAB	0684-00-0296-01	3000067963
MEGA 8Fr 50cc IAB	0684-00-0296-01	3000070668
MEGA 8Fr 50cc IAB	0684-00-0296-01	3000071291
MEGA 8Fr 50cc IAB	0684-00-0296-01	3000071979
MEGA 8Fr 50cc IAB	0684-00-0296-01	3000073516
MEGA 8Fr 50cc IAB	0684-00-0296-01	3000074146
MEGA 8Fr 50cc IAB	0684-00-0296-01	3000075999
MEGA 8Fr 50cc IAB	0684-00-0296-01	3000078921
MEGA 8Fr 50cc IAB	0684-00-0296-01	3000080138
MEGA 8Fr 50cc IAB	0684-00-0296-01	3000080137
MEGA 8Fr 50cc IAB	0684-00-0296-01	3000081624
MEGA 8Fr 50cc IAB	0684-00-0296-01	3000082197
MEGA 8Fr 50cc IAB	0684-00-0296-01	3000082848
MEGA 8Fr 50cc IAB	0684-00-0296-01	3000083430
MEGA 8Fr 50cc IAB	0684-00-0296-01	3000084785
MEGA 8Fr 50cc IAB	0684-00-0296-01	3000086995
MEGA 8Fr 50cc IAB	0684-00-0296-01	3000087888
MEGA 8Fr 50cc IAB	0684-00-0296-01	3000089898
MEGA 8Fr 50cc IAB	0684-00-0296-01	3000091753

MEGA 8Fr 50cc IAB	0684-00-0296-01	3000092239
MEGA 8Fr 50cc IAB	0684-00-0296-01	3000092240
MEGA 8Fr 50cc IAB	0684-00-0296-01	3000094450
MEGA 8Fr 50cc IAB	0684-00-0296-01	3000094451
MEGA 8Fr 50cc IAB	0684-00-0296-01	3000095604
MEGA 8Fr 50cc IAB	0684-00-0296-01	3000096970
MEGA 8Fr 50cc IAB	0684-00-0296-01	3000097957
MEGA 8Fr 50cc IAB	0684-00-0296-01	3000100722
MEGA 8Fr 50cc IAB	0684-00-0296-01	3000105063
MEGA 8Fr 50cc IAB	0684-00-0296-01	3000102628
MEGA 8Fr 50cc IAB	0684-00-0296-01	3000105064
MEGA 8Fr 50cc IAB	0684-00-0296-01	3000107101
MEGA 8Fr 50cc IAB	0684-00-0296-01	3000108962
MEGA 8Fr 50cc IAB	0684-00-0296-01	3000109818
MEGA 8Fr 50cc IAB	0684-00-0296-01	3000112782
MEGA 8Fr 50cc IAB	0684-00-0296-01U	3000046415
MEGA 8Fr 50cc IAB	0684-00-0296-01U	3000047430
MEGA 8Fr 50cc IAB	0684-00-0296-01U	3000048245
MEGA 8Fr 50cc IAB	0684-00-0296-01U	3000049085
MEGA 8Fr 50cc IAB	0684-00-0296-01U	3000049570
MEGA 8Fr 50cc IAB	0684-00-0296-01U	3000051153
MEGA 8Fr 50cc IAB	0684-00-0296-01U	3000051881
MEGA 8Fr 50cc IAB	0684-00-0296-01U	3000053310
MEGA 8Fr 50cc IAB	0684-00-0296-01U	3000054729
MEGA 8Fr 50cc IAB	0684-00-0296-01U	3000056538
MEGA 8Fr 50cc IAB	0684-00-0296-01U	3000059168
MEGA 8Fr 50cc IAB	0684-00-0296-01U	3000059618
MEGA 8Fr 50cc IAB	0684-00-0296-01U	3000060730
MEGA 8Fr 50cc IAB	0684-00-0296-01U	3000060451
MEGA 8Fr 50cc IAB	0684-00-0296-01U	3000064248
MEGA 8Fr 50cc IAB	0684-00-0296-01U	3000066798
MEGA 8Fr 50cc IAB	0684-00-0296-01U	3000065636
MEGA 8Fr 50cc IAB	0684-00-0296-01U	3000066799
MEGA 8Fr 50cc IAB	0684-00-0296-01U	3000067964
MEGA 8Fr 50cc IAB	0684-00-0296-01U	3000070669
MEGA 8Fr 50cc IAB	0684-00-0296-01U	3000071292
MEGA 8Fr 50cc IAB	0684-00-0296-01U	3000071980
MEGA 8Fr 50cc IAB	0684-00-0296-01U	3000073517
MEGA 8Fr 50cc IAB	0684-00-0296-01U	3000074144
MEGA 8Fr 50cc IAB	0684-00-0296-01U	3000074145
MEGA 8Fr 50cc IAB	0684-00-0296-01U	3000076000
MEGA 8Fr 50cc IAB	0684-00-0296-01U	3000078922
MEGA 8Fr 50cc IAB	0684-00-0296-01U	3000080140
MEGA 8Fr 50cc IAB	0684-00-0296-01U	3000082847
MEGA 8Fr 50cc IAB	0684-00-0296-01U	3000083431
MEGA 8Fr 50cc IAB	0684-00-0296-01U	3000083989
MEGA 8Fr 50cc IAB	0684-00-0296-01U	3000084730
MEGA 8Fr 50cc IAB	0684-00-0296-01U	3000087112
MEGA 8Fr 50cc IAB	0684-00-0296-01U	3000088376
MEGA 8Fr 50cc IAB	0684-00-0296-01U	3000089899
MEGA 8Fr 50cc IAB	0684-00-0296-01U	3000091754

MEGA 8Fr 50cc IAB	0684-00-0296-01U	3000092134
MEGA 8Fr 50cc IAB	0684-00-0296-01U	3000094449
MEGA 8Fr 50cc IAB	0684-00-0296-01U	3000095692
MEGA 8Fr 50cc IAB	0684-00-0296-01U	3000096971
MEGA 8Fr 50cc IAB	0684-00-0296-01U	3000097952
MEGA 8Fr 50cc IAB	0684-00-0296-01U	3000097953
MEGA 8Fr 50cc IAB	0684-00-0296-01U	3000100721
MEGA 8Fr 50cc IAB	0684-00-0296-01U	3000102697
MEGA 8Fr 50cc IAB	0684-00-0296-01U	3000105062
MEGA 8Fr 50cc IAB	0684-00-0296-01U	3000107100
MEGA 8Fr 50cc IAB	0684-00-0296-01U	3000109817
MEGA 8Fr 50cc IAB	0684-00-0296-01U	3000111681
MEGA 8Fr 50cc IAB	0684-00-0296-01U	3000113153
MEGA 8Fr 50cc IAB	0684-00-0296-02	3000048244
MEGA 8Fr 50cc IAB	0684-00-0296-02	3000049567
MEGA 8Fr 50cc IAB	0684-00-0296-02	3000050496
MEGA 8Fr 50cc IAB	0684-00-0296-02	3000051155
MEGA 8Fr 50cc IAB	0684-00-0296-02	3000053312
MEGA 8Fr 50cc IAB	0684-00-0296-02	3000054730
MEGA 8Fr 50cc IAB	0684-00-0296-02	3000059169
MEGA 8Fr 50cc IAB	0684-00-0296-02	3000060545
MEGA 8Fr 50cc IAB	0684-00-0296-02	3000060543
MEGA 8Fr 50cc IAB	0684-00-0296-02	3000066063
MEGA 8Fr 50cc IAB	0684-00-0296-02	3000064842
MEGA 8Fr 50cc IAB	0684-00-0296-02	3000064843
MEGA 8Fr 50cc IAB	0684-00-0296-02	3000066801
MEGA 8Fr 50cc IAB	0684-00-0296-02	3000066802
MEGA 8Fr 50cc IAB	0684-00-0296-02	3000066803
MEGA 8Fr 50cc IAB	0684-00-0296-02	3000076001
MEGA 8Fr 50cc IAB	0684-00-0296-02	3000080543
MEGA 8Fr 50cc IAB	0684-00-0296-02	3000080544
MEGA 8Fr 50cc IAB	0684-00-0296-02	3000081621
MEGA 8Fr 50cc IAB	0684-00-0296-02	3000082189
MEGA 8Fr 50cc IAB	0684-00-0296-02	3000084786
MEGA 8Fr 50cc IAB	0684-00-0296-02	3000086996
MEGA 8Fr 50cc IAB	0684-00-0296-02	3000091755
MEGA 8Fr 50cc IAB	0684-00-0296-02	3000092238
MEGA 8Fr 50cc IAB	0684-00-0296-02	3000097956
MEGA 8Fr 50cc IAB	0684-00-0296-02	3000102698
MEGA 8Fr 50cc IAB	0684-00-0296-02	3000105060
MEGA 8Fr 50cc IAB	0684-00-0296-02	3000106965
MEGA 8Fr 50cc IAB	0684-00-0296-02	3000112155
MEGA 8Fr 50cc IAB	0684-00-0296-02	3000112156
MEGA 8Fr 50cc IAB	0684-00-0296-02U	3000046416
MEGA 8Fr 50cc IAB	0684-00-0296-02U	3000047435
MEGA 8Fr 50cc IAB	0684-00-0296-02U	3000049568
MEGA 8Fr 50cc IAB	0684-00-0296-02U	3000056539
MEGA 8Fr 50cc IAB	0684-00-0296-02U	3000059619
MEGA 8Fr 50cc IAB	0684-00-0296-02U	3000064250
MEGA 8Fr 50cc IAB	0684-00-0296-02U	3000066804
MEGA 8Fr 50cc IAB	0684-00-0296-02U	3000066805

MEGA 8Fr 50cc IAB	0684-00-0296-02U	3000067857
MEGA 8Fr 50cc IAB	0684-00-0296-02U	3000070670
MEGA 8Fr 50cc IAB	0684-00-0296-02U	3000073518
MEGA 8Fr 50cc IAB	0684-00-0296-02U	3000074162
MEGA 8Fr 50cc IAB	0684-00-0296-02U	3000078923
MEGA 8Fr 50cc IAB	0684-00-0296-02U	3000080139
MEGA 8Fr 50cc IAB	0684-00-0296-02U	3000084731
MEGA 8Fr 50cc IAB	0684-00-0296-02U	3000097954
MEGA 8Fr 50cc IAB	0684-00-0296-02U	3000102699
MEGA 8Fr 50cc IAB	0684-00-0296-02U	3000108432
MEGA 8Fr 50cc IAB	0684-00-0296-02U	3000112154
MEGA 8Fr 50cc IAB	0684-00-0296-10	3000103022
Sensation 7Fr 34cc IAB	0684-00-0469-01	3000045657
Sensation 7Fr 34cc IAB	0684-00-0469-01	3000047164
Sensation 7Fr 34cc IAB	0684-00-0469-01	3000048708
Sensation 7Fr 34cc IAB	0684-00-0469-01	3000049477
Sensation 7Fr 34cc IAB	0684-00-0469-01	3000050927
Sensation 7Fr 34cc IAB	0684-00-0469-01	3000051753
Sensation 7Fr 34cc IAB	0684-00-0469-01	3000051754
Sensation 7Fr 34cc IAB	0684-00-0469-01	3000052790
Sensation 7Fr 34cc IAB	0684-00-0469-01	3000053815
Sensation 7Fr 34cc IAB	0684-00-0469-01	3000053816
Sensation 7Fr 34cc IAB	0684-00-0469-01	3000056436
Sensation 7Fr 34cc IAB	0684-00-0469-01	3000058893
Sensation 7Fr 34cc IAB	0684-00-0469-01	3000059906
Sensation 7Fr 34cc IAB	0684-00-0469-01	3000063263
Sensation 7Fr 34cc IAB	0684-00-0469-01	3000064521
Sensation 7Fr 34cc IAB	0684-00-0469-01	3000065233
Sensation 7Fr 34cc IAB	0684-00-0469-01	3000066677
Sensation 7Fr 34cc IAB	0684-00-0469-01	3000067689
Sensation 7Fr 34cc IAB	0684-00-0469-01	3000071413
Sensation 7Fr 34cc IAB	0684-00-0469-01	3000072399
Sensation 7Fr 34cc IAB	0684-00-0469-01	3000074269
Sensation 7Fr 34cc IAB	0684-00-0469-01	3000075922
Sensation 7Fr 34cc IAB	0684-00-0469-01	3000077813
Sensation 7Fr 34cc IAB	0684-00-0469-01	3000080874
Sensation 7Fr 34cc IAB	0684-00-0469-01	3000081078
Sensation 7Fr 34cc IAB	0684-00-0469-01	3000083096
Sensation 7Fr 34cc IAB	0684-00-0469-01	3000084472
Sensation 7Fr 34cc IAB	0684-00-0469-01	3000084787
Sensation 7Fr 34cc IAB	0684-00-0469-01	3000086850
Sensation 7Fr 34cc IAB	0684-00-0469-01	3000087366
Sensation 7Fr 34cc IAB	0684-00-0469-01	3000088618
Sensation 7Fr 34cc IAB	0684-00-0469-01	3000089991
Sensation 7Fr 34cc IAB	0684-00-0469-01	3000090768
Sensation 7Fr 34cc IAB	0684-00-0469-01	3000094337
Sensation 7Fr 34cc IAB	0684-00-0469-01	3000097300
Sensation 7Fr 34cc IAB	0684-00-0469-01	3000098247
Sensation 7Fr 34cc IAB	0684-00-0469-01	3000099929
Sensation 7Fr 34cc IAB	0684-00-0469-01	3000102977
Sensation 7Fr 34cc IAB	0684-00-0469-01	3000106755

Sensation 7Fr 34cc IAB	0684-00-0469-01	3000110208
Sensation 7Fr 34cc IAB	0684-00-0469-01	3000111136
Sensation 7Fr 34cc IAB	0684-00-0469-01C	3000058894
Sensation 7Fr 34cc IAB	0684-00-0469-01C	3000070665
Sensation 7Fr 34cc IAB	0684-00-0469-01C	3000073189
Sensation 7Fr 34cc IAB	0684-00-0469-01C	3000077203
Sensation 7Fr 34cc IAB	0684-00-0469-01C	3000080875
Sensation 7Fr 34cc IAB	0684-00-0469-01C	3000084690
Sensation 7Fr 34cc IAB	0684-00-0469-01C	3000086792
Sensation 7Fr 34cc IAB	0684-00-0469-01C	3000089348
Sensation 7Fr 34cc IAB	0684-00-0469-01C	3000094338
Sensation 7Fr 34cc IAB	0684-00-0469-01C	3000103213
Sensation 7Fr 34cc IAB	0684-00-0469-01C	3000110209
Sensation 7Fr 34cc IAB	0684-00-0469-01U	3000045656
Sensation 7Fr 34cc IAB	0684-00-0469-01U	3000047710
Sensation 7Fr 34cc IAB	0684-00-0469-01U	3000048709
Sensation 7Fr 34cc IAB	0684-00-0469-01U	3000049478
Sensation 7Fr 34cc IAB	0684-00-0469-01U	3000050926
Sensation 7Fr 34cc IAB	0684-00-0469-01U	3000051752
Sensation 7Fr 34cc IAB	0684-00-0469-01U	3000052791
Sensation 7Fr 34cc IAB	0684-00-0469-01U	3000053817
Sensation 7Fr 34cc IAB	0684-00-0469-01U	3000056437
Sensation 7Fr 34cc IAB	0684-00-0469-01U	3000059907
Sensation 7Fr 34cc IAB	0684-00-0469-01U	3000061290
Sensation 7Fr 34cc IAB	0684-00-0469-01U	3000061291
Sensation 7Fr 34cc IAB	0684-00-0469-01U	3000064520
Sensation 7Fr 34cc IAB	0684-00-0469-01U	3000065234
Sensation 7Fr 34cc IAB	0684-00-0469-01U	3000065987
Sensation 7Fr 34cc IAB	0684-00-0469-01U	3000066678
Sensation 7Fr 34cc IAB	0684-00-0469-01U	3000067690
Sensation 7Fr 34cc IAB	0684-00-0469-01U	3000068924
Sensation 7Fr 34cc IAB	0684-00-0469-01U	3000069876
Sensation 7Fr 34cc IAB	0684-00-0469-01U	3000071218
Sensation 7Fr 34cc IAB	0684-00-0469-01U	3000071977
Sensation 7Fr 34cc IAB	0684-00-0469-01U	3000071978
Sensation 7Fr 34cc IAB	0684-00-0469-01U	3000073188
Sensation 7Fr 34cc IAB	0684-00-0469-01U	3000074270
Sensation 7Fr 34cc IAB	0684-00-0469-01U	3000076924
Sensation 7Fr 34cc IAB	0684-00-0469-01U	3000077812
Sensation 7Fr 34cc IAB	0684-00-0469-01U	3000080877
Sensation 7Fr 34cc IAB	0684-00-0469-01U	3000081077
Sensation 7Fr 34cc IAB	0684-00-0469-01U	3000083095
Sensation 7Fr 34cc IAB	0684-00-0469-01U	3000086851
Sensation 7Fr 34cc IAB	0684-00-0469-01U	3000087367
Sensation 7Fr 34cc IAB	0684-00-0469-01U	3000088617
Sensation 7Fr 34cc IAB	0684-00-0469-01U	3000090767
Sensation 7Fr 34cc IAB	0684-00-0469-01U	3000092537
Sensation 7Fr 34cc IAB	0684-00-0469-01U	3000094336
Sensation 7Fr 34cc IAB	0684-00-0469-01U	3000097301
Sensation 7Fr 34cc IAB	0684-00-0469-01U	3000098248
Sensation 7Fr 34cc IAB	0684-00-0469-01U	3000099928

Sensation 7Fr 34cc IAB	0684-00-0469-01U	3000102978
Sensation 7Fr 34cc IAB	0684-00-0469-01U	3000106754
Sensation 7Fr 34cc IAB	0684-00-0469-01U	3000110207
Sensation 7Fr 34cc IAB	0684-00-0469-01U	3000111137
Sensation 7Fr 34cc IAB	0684-00-0469-07	3000046410
Sensation 7Fr 34cc IAB	0684-00-0469-07	3000083097
Sensation 7Fr 34cc IAB	0684-00-0469-07	3000111135
Sensation 7Fr 34cc IAB	0684-00-0469-09	3000103214
Sensation 7Fr 34cc IAB	0684-00-0469-09	3000110210
Sensation 7Fr 34cc IAB	0684-00-0469-09	3000111192
Sensation 7Fr 40cc IAB	0684-00-0470-01	3000046316
Sensation 7Fr 40cc IAB	0684-00-0470-01	3000047161
Sensation 7Fr 40cc IAB	0684-00-0470-01	3000047162
Sensation 7Fr 40cc IAB	0684-00-0470-01	3000048615
Sensation 7Fr 40cc IAB	0684-00-0470-01	3000049262
Sensation 7Fr 40cc IAB	0684-00-0470-01	3000049914
Sensation 7Fr 40cc IAB	0684-00-0470-01	3000050929
Sensation 7Fr 40cc IAB	0684-00-0470-01	3000051380
Sensation 7Fr 40cc IAB	0684-00-0470-01	3000052645
Sensation 7Fr 40cc IAB	0684-00-0470-01	3000054090
Sensation 7Fr 40cc IAB	0684-00-0470-01	3000055099
Sensation 7Fr 40cc IAB	0684-00-0470-01	3000056537
Sensation 7Fr 40cc IAB	0684-00-0470-01	3000059946
Sensation 7Fr 40cc IAB	0684-00-0470-01	3000061447
Sensation 7Fr 40cc IAB	0684-00-0470-01	3000062406
Sensation 7Fr 40cc IAB	0684-00-0470-01	3000065986
Sensation 7Fr 40cc IAB	0684-00-0470-01	3000066489
Sensation 7Fr 40cc IAB	0684-00-0470-01	3000066764
Sensation 7Fr 40cc IAB	0684-00-0470-01	3000066770
Sensation 7Fr 40cc IAB	0684-00-0470-01	3000068730
Sensation 7Fr 40cc IAB	0684-00-0470-01	3000069160
Sensation 7Fr 40cc IAB	0684-00-0470-01	3000069376
Sensation 7Fr 40cc IAB	0684-00-0470-01	3000071293
Sensation 7Fr 40cc IAB	0684-00-0470-01	3000071739
Sensation 7Fr 40cc IAB	0684-00-0470-01	3000073191
Sensation 7Fr 40cc IAB	0684-00-0470-01	3000073277
Sensation 7Fr 40cc IAB	0684-00-0470-01	3000074271
Sensation 7Fr 40cc IAB	0684-00-0470-01	3000075134
Sensation 7Fr 40cc IAB	0684-00-0470-01	3000075676
Sensation 7Fr 40cc IAB	0684-00-0470-01	3000076957
Sensation 7Fr 40cc IAB	0684-00-0470-01	3000079650
Sensation 7Fr 40cc IAB	0684-00-0470-01	3000081778
Sensation 7Fr 40cc IAB	0684-00-0470-01	3000084071
Sensation 7Fr 40cc IAB	0684-00-0470-01	3000084788
Sensation 7Fr 40cc IAB	0684-00-0470-01	3000084789
Sensation 7Fr 40cc IAB	0684-00-0470-01	3000086924
Sensation 7Fr 40cc IAB	0684-00-0470-01	3000085291
Sensation 7Fr 40cc IAB	0684-00-0470-01	3000087368
Sensation 7Fr 40cc IAB	0684-00-0470-01	3000088278
Sensation 7Fr 40cc IAB	0684-00-0470-01	3000088570
Sensation 7Fr 40cc IAB	0684-00-0470-01	3000090526

Sensation 7Fr 40cc IAB	0684-00-0470-01	3000090664
Sensation 7Fr 40cc IAB	0684-00-0470-01	3000092443
Sensation 7Fr 40cc IAB	0684-00-0470-01	3000096052
Sensation 7Fr 40cc IAB	0684-00-0470-01	3000097302
Sensation 7Fr 40cc IAB	0684-00-0470-01	3000098997
Sensation 7Fr 40cc IAB	0684-00-0470-01	3000101522
Sensation 7Fr 40cc IAB	0684-00-0470-01	3000103209
Sensation 7Fr 40cc IAB	0684-00-0470-01	3000106823
Sensation 7Fr 40cc IAB	0684-00-0470-01	3000110212
Sensation 7Fr 40cc IAB	0684-00-0470-01	3000111023
Sensation 7Fr 40cc IAB	0684-00-0470-01	3000112113
Sensation 7Fr 40cc IAB	0684-00-0470-01C	3000058895
Sensation 7Fr 40cc IAB	0684-00-0470-01C	3000064045
Sensation 7Fr 40cc IAB	0684-00-0470-01C	3000073190
Sensation 7Fr 40cc IAB	0684-00-0470-01C	3000077204
Sensation 7Fr 40cc IAB	0684-00-0470-01C	3000084689
Sensation 7Fr 40cc IAB	0684-00-0470-01C	3000086162
Sensation 7Fr 40cc IAB	0684-00-0470-01C	3000089349
Sensation 7Fr 40cc IAB	0684-00-0470-01C	3000094526
Sensation 7Fr 40cc IAB	0684-00-0470-01C	3000103212
Sensation 7Fr 40cc IAB	0684-00-0470-01C	3000110211
Sensation 7Fr 40cc IAB	0684-00-0470-01U	3000046315
Sensation 7Fr 40cc IAB	0684-00-0470-01U	3000047053
Sensation 7Fr 40cc IAB	0684-00-0470-01U	3000047709
Sensation 7Fr 40cc IAB	0684-00-0470-01U	3000048710
Sensation 7Fr 40cc IAB	0684-00-0470-01U	3000049479
Sensation 7Fr 40cc IAB	0684-00-0470-01U	3000050215
Sensation 7Fr 40cc IAB	0684-00-0470-01U	3000050930
Sensation 7Fr 40cc IAB	0684-00-0470-01U	3000051379
Sensation 7Fr 40cc IAB	0684-00-0470-01U	3000052646
Sensation 7Fr 40cc IAB	0684-00-0470-01U	3000054089
Sensation 7Fr 40cc IAB	0684-00-0470-01U	3000055098
Sensation 7Fr 40cc IAB	0684-00-0470-01U	3000056438
Sensation 7Fr 40cc IAB	0684-00-0470-01U	3000060258
Sensation 7Fr 40cc IAB	0684-00-0470-01U	3000060259
Sensation 7Fr 40cc IAB	0684-00-0470-01U	3000061292
Sensation 7Fr 40cc IAB	0684-00-0470-01U	3000063264
Sensation 7Fr 40cc IAB	0684-00-0470-01U	3000064522
Sensation 7Fr 40cc IAB	0684-00-0470-01U	3000066679
Sensation 7Fr 40cc IAB	0684-00-0470-01U	3000065125
Sensation 7Fr 40cc IAB	0684-00-0470-01U	3000065984
Sensation 7Fr 40cc IAB	0684-00-0470-01U	3000065985
Sensation 7Fr 40cc IAB	0684-00-0470-01U	3000066515
Sensation 7Fr 40cc IAB	0684-00-0470-01U	3000066680
Sensation 7Fr 40cc IAB	0684-00-0470-01U	3000067692
Sensation 7Fr 40cc IAB	0684-00-0470-01U	3000069161
Sensation 7Fr 40cc IAB	0684-00-0470-01U	3000069872
Sensation 7Fr 40cc IAB	0684-00-0470-01U	3000069873
Sensation 7Fr 40cc IAB	0684-00-0470-01U	3000071217
Sensation 7Fr 40cc IAB	0684-00-0470-01U	3000071891
Sensation 7Fr 40cc IAB	0684-00-0470-01U	3000073192

Sensation 7Fr 40cc IAB	0684-00-0470-01U	3000074272
Sensation 7Fr 40cc IAB	0684-00-0470-01U	3000075133
Sensation 7Fr 40cc IAB	0684-00-0470-01U	3000076154
Sensation 7Fr 40cc IAB	0684-00-0470-01U	3000076956
Sensation 7Fr 40cc IAB	0684-00-0470-01U	3000079648
Sensation 7Fr 40cc IAB	0684-00-0470-01U	3000079649
Sensation 7Fr 40cc IAB	0684-00-0470-01U	3000080226
Sensation 7Fr 40cc IAB	0684-00-0470-01U	3000081779
Sensation 7Fr 40cc IAB	0684-00-0470-01U	3000083059
Sensation 7Fr 40cc IAB	0684-00-0470-01U	3000084153
Sensation 7Fr 40cc IAB	0684-00-0470-01U	3000086058
Sensation 7Fr 40cc IAB	0684-00-0470-01U	3000087369
Sensation 7Fr 40cc IAB	0684-00-0470-01U	3000088277
Sensation 7Fr 40cc IAB	0684-00-0470-01U	3000089284
Sensation 7Fr 40cc IAB	0684-00-0470-01U	3000090665
Sensation 7Fr 40cc IAB	0684-00-0470-01U	3000092442
Sensation 7Fr 40cc IAB	0684-00-0470-01U	3000094528
Sensation 7Fr 40cc IAB	0684-00-0470-01U	3000094531
Sensation 7Fr 40cc IAB	0684-00-0470-01U	3000096053
Sensation 7Fr 40cc IAB	0684-00-0470-01U	3000096054
Sensation 7Fr 40cc IAB	0684-00-0470-01U	3000098998
Sensation 7Fr 40cc IAB	0684-00-0470-01U	3000098999
Sensation 7Fr 40cc IAB	0684-00-0470-01U	3000101523
Sensation 7Fr 40cc IAB	0684-00-0470-01U	3000103210
Sensation 7Fr 40cc IAB	0684-00-0470-01U	3000106825
Sensation 7Fr 40cc IAB	0684-00-0470-01U	3000107337
Sensation 7Fr 40cc IAB	0684-00-0470-01U	3000111022
Sensation 7Fr 40cc IAB	0684-00-0470-01U	3000112114
Sensation 7Fr 40cc IAB	0684-00-0470-07	3000046411
Sensation 7Fr 40cc IAB	0684-00-0470-07	3000069162
Sensation 7Fr 40cc IAB	0684-00-0470-07	3000069442
Sensation 7Fr 40cc IAB	0684-00-0470-07	3000071837
Sensation 7Fr 40cc IAB	0684-00-0470-07	3000094141
Sensation 7Fr 40cc IAB	0684-00-0470-07	3000111134
Sensation 7Fr 40cc IAB	0684-00-0470-09	3000094151
Sensation 7Fr 40cc IAB	0684-00-0470-09	3000094527
Sensation 7Fr 40cc IAB	0684-00-0470-09	3000096001
Sensation 7Fr 40cc IAB	0684-00-0470-09	3000099000
Sensation 7Fr 40cc IAB	0684-00-0470-09	3000103211
Sensation 7Fr 40cc IAB	0684-00-0470-09	3000106824
Sensation 7Fr 40cc IAB	0684-00-0470-09	3000107336
Sensation 7Fr 40cc IAB	0684-00-0470-09	3000111191
Sensation 7Fr 40cc IAB	0684-00-0470-09	3000112115
Linear 7.5Fr 25cc IAB	0684-00-0478-01	3000049797
Linear 7.5Fr 25cc IAB	0684-00-0478-01	3000053552
Linear 7.5Fr 25cc IAB	0684-00-0478-01	3000056193
Linear 7.5Fr 25cc IAB	0684-00-0478-01	3000056201
Linear 7.5Fr 25cc IAB	0684-00-0478-01	3000057749
Linear 7.5Fr 25cc IAB	0684-00-0478-01	3000066195
Linear 7.5Fr 25cc IAB	0684-00-0478-01	3000067965
Linear 7.5Fr 25cc IAB	0684-00-0478-01	3000069878

Linear 7.5Fr 25cc IAB	0684-00-0478-01	3000071222
Linear 7.5Fr 25cc IAB	0684-00-0478-01	3000072585
Linear 7.5Fr 25cc IAB	0684-00-0478-01	3000074287
Linear 7.5Fr 25cc IAB	0684-00-0478-01	3000078937
Linear 7.5Fr 25cc IAB	0684-00-0478-01	3000079446
Linear 7.5Fr 25cc IAB	0684-00-0478-01	3000080540
Linear 7.5Fr 25cc IAB	0684-00-0478-01	3000084470
Linear 7.5Fr 25cc IAB	0684-00-0478-01	3000090372
Linear 7.5Fr 25cc IAB	0684-00-0478-01	3000092553
Linear 7.5Fr 25cc IAB	0684-00-0478-01	3000095910
Linear 7.5Fr 25cc IAB	0684-00-0478-01	3000099930
Linear 7.5Fr 25cc IAB	0684-00-0478-01	3000103632
Linear 7.5Fr 25cc IAB	0684-00-0478-01	3000112148
Linear 7.5Fr 25cc IAB	0684-00-0478-01C	3000053297
Linear 7.5Fr 25cc IAB	0684-00-0478-01C	3000057125
Linear 7.5Fr 25cc IAB	0684-00-0478-01C	3000072586
Linear 7.5Fr 25cc IAB	0684-00-0478-01U	3000049798
Linear 7.5Fr 25cc IAB	0684-00-0478-01U	3000053553
Linear 7.5Fr 25cc IAB	0684-00-0478-01U	3000057748
Linear 7.5Fr 25cc IAB	0684-00-0478-01U	3000061193
Linear 7.5Fr 25cc IAB	0684-00-0478-01U	3000066264
Linear 7.5Fr 25cc IAB	0684-00-0478-01U	3000069270
Linear 7.5Fr 25cc IAB	0684-00-0478-01U	3000071221
Linear 7.5Fr 25cc IAB	0684-00-0478-01U	3000072584
Linear 7.5Fr 25cc IAB	0684-00-0478-01U	3000074288
Linear 7.5Fr 25cc IAB	0684-00-0478-01U	3000077341
Linear 7.5Fr 25cc IAB	0684-00-0478-01U	3000077814
Linear 7.5Fr 25cc IAB	0684-00-0478-01U	3000083060
Linear 7.5Fr 25cc IAB	0684-00-0478-01U	3000084471
Linear 7.5Fr 25cc IAB	0684-00-0478-01U	3000087477
Linear 7.5Fr 25cc IAB	0684-00-0478-01U	3000087478
Linear 7.5Fr 25cc IAB	0684-00-0478-01U	3000088619
Linear 7.5Fr 25cc IAB	0684-00-0478-01U	3000090373
Linear 7.5Fr 25cc IAB	0684-00-0478-01U	3000092552
Linear 7.5Fr 25cc IAB	0684-00-0478-01U	3000095911
Linear 7.5Fr 25cc IAB	0684-00-0478-01U	3000099932
Linear 7.5Fr 25cc IAB	0684-00-0478-01U	3000103631
Linear 7.5Fr 25cc IAB	0684-00-0478-02	3000053554
Linear 7.5Fr 25cc IAB	0684-00-0478-02	3000056540
Linear 7.5Fr 25cc IAB	0684-00-0478-02	3000061194
Linear 7.5Fr 25cc IAB	0684-00-0478-02	3000068728
Linear 7.5Fr 25cc IAB	0684-00-0478-02	3000092554
Linear 7.5Fr 25cc IAB	0684-00-0478-02	3000095912
Linear 7.5Fr 25cc IAB	0684-00-0478-02	3000099933
Linear 7.5Fr 25cc IAB	0684-00-0478-02	3000102246
Linear 7.5Fr 25cc IAB	0684-00-0478-02U	3000049800
Linear 7.5Fr 25cc IAB	0684-00-0478-02U	3000074289
Linear 7.5Fr 25cc IAB	0684-00-0478-02U	3000092555
Linear 7.5Fr 25cc IAB	0684-00-0478-02U	3000099931
Linear 7.5Fr 34cc IAB	0684-00-0479-01	3000047635
Linear 7.5Fr 34cc IAB	0684-00-0479-01	3000048712

Linear 7.5Fr 34cc IAB	0684-00-0479-01	3000048711
Linear 7.5Fr 34cc IAB	0684-00-0479-01	3000049311
Linear 7.5Fr 34cc IAB	0684-00-0479-01	3000051296
Linear 7.5Fr 34cc IAB	0684-00-0479-01	3000052326
Linear 7.5Fr 34cc IAB	0684-00-0479-01	3000053634
Linear 7.5Fr 34cc IAB	0684-00-0479-01	3000054821
Linear 7.5Fr 34cc IAB	0684-00-0479-01	3000056194
Linear 7.5Fr 34cc IAB	0684-00-0479-01	3000057204
Linear 7.5Fr 34cc IAB	0684-00-0479-01	3000058417
Linear 7.5Fr 34cc IAB	0684-00-0479-01	3000059272
Linear 7.5Fr 34cc IAB	0684-00-0479-01	3000059273
Linear 7.5Fr 34cc IAB	0684-00-0479-01	3000061448
Linear 7.5Fr 34cc IAB	0684-00-0479-01	3000062934
Linear 7.5Fr 34cc IAB	0684-00-0479-01	3000065133
Linear 7.5Fr 34cc IAB	0684-00-0479-01	3000065724
Linear 7.5Fr 34cc IAB	0684-00-0479-01	3000066487
Linear 7.5Fr 34cc IAB	0684-00-0479-01	3000068135
Linear 7.5Fr 34cc IAB	0684-00-0479-01	3000069877
Linear 7.5Fr 34cc IAB	0684-00-0479-01	3000071129
Linear 7.5Fr 34cc IAB	0684-00-0479-01	3000071890
Linear 7.5Fr 34cc IAB	0684-00-0479-01	3000073519
Linear 7.5Fr 34cc IAB	0684-00-0479-01	3000074021
Linear 7.5Fr 34cc IAB	0684-00-0479-01	3000075917
Linear 7.5Fr 34cc IAB	0684-00-0479-01	3000076955
Linear 7.5Fr 34cc IAB	0684-00-0479-01	3000079930
Linear 7.5Fr 34cc IAB	0684-00-0479-01	3000077811
Linear 7.5Fr 34cc IAB	0684-00-0479-01	3000079055
Linear 7.5Fr 34cc IAB	0684-00-0479-01	3000079517
Linear 7.5Fr 34cc IAB	0684-00-0479-01	3000080673
Linear 7.5Fr 34cc IAB	0684-00-0479-01	3000082253
Linear 7.5Fr 34cc IAB	0684-00-0479-01	3000083197
Linear 7.5Fr 34cc IAB	0684-00-0479-01	3000084068
Linear 7.5Fr 34cc IAB	0684-00-0479-01	3000084852
Linear 7.5Fr 34cc IAB	0684-00-0479-01	3000084853
Linear 7.5Fr 34cc IAB	0684-00-0479-01	3000085756
Linear 7.5Fr 34cc IAB	0684-00-0479-01	3000086925
Linear 7.5Fr 34cc IAB	0684-00-0479-01	3000088673
Linear 7.5Fr 34cc IAB	0684-00-0479-01	3000089610
Linear 7.5Fr 34cc IAB	0684-00-0479-01	3000091320
Linear 7.5Fr 34cc IAB	0684-00-0479-01	3000091698
Linear 7.5Fr 34cc IAB	0684-00-0479-01	3000093895
Linear 7.5Fr 34cc IAB	0684-00-0479-01	3000096739
Linear 7.5Fr 34cc IAB	0684-00-0479-01	3000100472
Linear 7.5Fr 34cc IAB	0684-00-0479-01	3000101934
Linear 7.5Fr 34cc IAB	0684-00-0479-01	3000105734
Linear 7.5Fr 34cc IAB	0684-00-0479-01	3000106532
Linear 7.5Fr 34cc IAB	0684-00-0479-01	3000106533
Linear 7.5Fr 34cc IAB	0684-00-0479-01	3000109449
Linear 7.5Fr 34cc IAB	0684-00-0479-01	3000113519
Linear 7.5Fr 34cc IAB	0684-00-0479-01	3000116780
Linear 7.5Fr 34cc IAB	0684-00-0479-01C	3000047354

Linear 7.5Fr 34cc IAB	0684-00-0479-01C	3000047707
Linear 7.5Fr 34cc IAB	0684-00-0479-01C	3000048027
Linear 7.5Fr 34cc IAB	0684-00-0479-01C	3000048963
Linear 7.5Fr 34cc IAB	0684-00-0479-01C	3000050259
Linear 7.5Fr 34cc IAB	0684-00-0479-01C	3000051645
Linear 7.5Fr 34cc IAB	0684-00-0479-01C	3000051646
Linear 7.5Fr 34cc IAB	0684-00-0479-01C	3000052792
Linear 7.5Fr 34cc IAB	0684-00-0479-01C	3000054317
Linear 7.5Fr 34cc IAB	0684-00-0479-01C	3000054873
Linear 7.5Fr 34cc IAB	0684-00-0479-01C	3000056073
Linear 7.5Fr 34cc IAB	0684-00-0479-01C	3000057383
Linear 7.5Fr 34cc IAB	0684-00-0479-01C	3000058898
Linear 7.5Fr 34cc IAB	0684-00-0479-01C	3000058899
Linear 7.5Fr 34cc IAB	0684-00-0479-01C	3000059609
Linear 7.5Fr 34cc IAB	0684-00-0479-01C	3000061524
Linear 7.5Fr 34cc IAB	0684-00-0479-01C	3000062026
Linear 7.5Fr 34cc IAB	0684-00-0479-01C	3000063084
Linear 7.5Fr 34cc IAB	0684-00-0479-01C	3000064055
Linear 7.5Fr 34cc IAB	0684-00-0479-01C	3000066140
Linear 7.5Fr 34cc IAB	0684-00-0479-01C	3000066513
Linear 7.5Fr 34cc IAB	0684-00-0479-01C	3000066514
Linear 7.5Fr 34cc IAB	0684-00-0479-01C	3000068859
Linear 7.5Fr 34cc IAB	0684-00-0479-01C	3000070407
Linear 7.5Fr 34cc IAB	0684-00-0479-01C	3000071128
Linear 7.5Fr 34cc IAB	0684-00-0479-01C	3000071415
Linear 7.5Fr 34cc IAB	0684-00-0479-01C	3000073402
Linear 7.5Fr 34cc IAB	0684-00-0479-01C	3000074029
Linear 7.5Fr 34cc IAB	0684-00-0479-01C	3000075865
Linear 7.5Fr 34cc IAB	0684-00-0479-01C	3000077274
Linear 7.5Fr 34cc IAB	0684-00-0479-01C	3000080105
Linear 7.5Fr 34cc IAB	0684-00-0479-01C	3000077560
Linear 7.5Fr 34cc IAB	0684-00-0479-01C	3000079568
Linear 7.5Fr 34cc IAB	0684-00-0479-01C	3000080106
Linear 7.5Fr 34cc IAB	0684-00-0479-01C	3000080444
Linear 7.5Fr 34cc IAB	0684-00-0479-01C	3000080732
Linear 7.5Fr 34cc IAB	0684-00-0479-01C	3000082660
Linear 7.5Fr 34cc IAB	0684-00-0479-01C	3000083200
Linear 7.5Fr 34cc IAB	0684-00-0479-01C	3000084477
Linear 7.5Fr 34cc IAB	0684-00-0479-01C	3000085500
Linear 7.5Fr 34cc IAB	0684-00-0479-01C	3000085845
Linear 7.5Fr 34cc IAB	0684-00-0479-01C	3000086491
Linear 7.5Fr 34cc IAB	0684-00-0479-01C	3000089458
Linear 7.5Fr 34cc IAB	0684-00-0479-01C	3000090850
Linear 7.5Fr 34cc IAB	0684-00-0479-01C	3000091740
Linear 7.5Fr 34cc IAB	0684-00-0479-01C	3000093114
Linear 7.5Fr 34cc IAB	0684-00-0479-01C	3000094658
Linear 7.5Fr 34cc IAB	0684-00-0479-01C	3000095998
Linear 7.5Fr 34cc IAB	0684-00-0479-01C	3000095999
Linear 7.5Fr 34cc IAB	0684-00-0479-01C	3000098650
Linear 7.5Fr 34cc IAB	0684-00-0479-01C	3000100381
Linear 7.5Fr 34cc IAB	0684-00-0479-01C	3000102537

Linear 7.5Fr 34cc IAB	0684-00-0479-01C	3000103918
Linear 7.5Fr 34cc IAB	0684-00-0479-01C	3000105531
Linear 7.5Fr 34cc IAB	0684-00-0479-01C	3000106633
Linear 7.5Fr 34cc IAB	0684-00-0479-01C	3000109450
Linear 7.5Fr 34cc IAB	0684-00-0479-01C	3000110736
Linear 7.5Fr 34cc IAB	0684-00-0479-01C	3000113520
Linear 7.5Fr 34cc IAB	0684-00-0479-01C	3000115115
Linear 7.5Fr 34cc IAB	0684-00-0479-01U	3000048818
Linear 7.5Fr 34cc IAB	0684-00-0479-01U	3000049422
Linear 7.5Fr 34cc IAB	0684-00-0479-01U	3000049791
Linear 7.5Fr 34cc IAB	0684-00-0479-01U	3000051295
Linear 7.5Fr 34cc IAB	0684-00-0479-01U	3000053635
Linear 7.5Fr 34cc IAB	0684-00-0479-01U	3000054874
Linear 7.5Fr 34cc IAB	0684-00-0479-01U	3000057267
Linear 7.5Fr 34cc IAB	0684-00-0479-01U	3000058418
Linear 7.5Fr 34cc IAB	0684-00-0479-01U	3000058848
Linear 7.5Fr 34cc IAB	0684-00-0479-01U	3000059458
Linear 7.5Fr 34cc IAB	0684-00-0479-01U	3000061477
Linear 7.5Fr 34cc IAB	0684-00-0479-01U	3000062027
Linear 7.5Fr 34cc IAB	0684-00-0479-01U	3000065130
Linear 7.5Fr 34cc IAB	0684-00-0479-01U	3000065722
Linear 7.5Fr 34cc IAB	0684-00-0479-01U	3000066766
Linear 7.5Fr 34cc IAB	0684-00-0479-01U	3000068923
Linear 7.5Fr 34cc IAB	0684-00-0479-01U	3000070426
Linear 7.5Fr 34cc IAB	0684-00-0479-01U	3000070663
Linear 7.5Fr 34cc IAB	0684-00-0479-01U	3000071289
Linear 7.5Fr 34cc IAB	0684-00-0479-01U	3000072588
Linear 7.5Fr 34cc IAB	0684-00-0479-01U	3000072587
Linear 7.5Fr 34cc IAB	0684-00-0479-01U	3000073520
Linear 7.5Fr 34cc IAB	0684-00-0479-01U	3000074020
Linear 7.5Fr 34cc IAB	0684-00-0479-01U	3000075914
Linear 7.5Fr 34cc IAB	0684-00-0479-01U	3000076953
Linear 7.5Fr 34cc IAB	0684-00-0479-01U	3000076954
Linear 7.5Fr 34cc IAB	0684-00-0479-01U	3000079932
Linear 7.5Fr 34cc IAB	0684-00-0479-01U	3000077810
Linear 7.5Fr 34cc IAB	0684-00-0479-01U	3000079516
Linear 7.5Fr 34cc IAB	0684-00-0479-01U	3000080539
Linear 7.5Fr 34cc IAB	0684-00-0479-01U	3000082534
Linear 7.5Fr 34cc IAB	0684-00-0479-01U	3000083198
Linear 7.5Fr 34cc IAB	0684-00-0479-01U	3000084476
Linear 7.5Fr 34cc IAB	0684-00-0479-01U	3000085814
Linear 7.5Fr 34cc IAB	0684-00-0479-01U	3000086855
Linear 7.5Fr 34cc IAB	0684-00-0479-01U	3000088103
Linear 7.5Fr 34cc IAB	0684-00-0479-01U	3000089732
Linear 7.5Fr 34cc IAB	0684-00-0479-01U	3000091321
Linear 7.5Fr 34cc IAB	0684-00-0479-01U	3000091322
Linear 7.5Fr 34cc IAB	0684-00-0479-01U	3000093116
Linear 7.5Fr 34cc IAB	0684-00-0479-01U	3000094007
Linear 7.5Fr 34cc IAB	0684-00-0479-01U	3000100469
Linear 7.5Fr 34cc IAB	0684-00-0479-01U	3000100470
Linear 7.5Fr 34cc IAB	0684-00-0479-01U	3000101888

Linear 7.5Fr 34cc IAB	0684-00-0479-01U	3000101932
Linear 7.5Fr 34cc IAB	0684-00-0479-01U	3000105387
Linear 7.5Fr 34cc IAB	0684-00-0479-01U	3000106634
Linear 7.5Fr 34cc IAB	0684-00-0479-01U	3000109452
Linear 7.5Fr 34cc IAB	0684-00-0479-02	3000048614
Linear 7.5Fr 34cc IAB	0684-00-0479-02	3000049263
Linear 7.5Fr 34cc IAB	0684-00-0479-02	3000049312
Linear 7.5Fr 34cc IAB	0684-00-0479-02	3000051297
Linear 7.5Fr 34cc IAB	0684-00-0479-02	3000051298
Linear 7.5Fr 34cc IAB	0684-00-0479-02	3000051299
Linear 7.5Fr 34cc IAB	0684-00-0479-02	3000053766
Linear 7.5Fr 34cc IAB	0684-00-0479-02	3000054338
Linear 7.5Fr 34cc IAB	0684-00-0479-02	3000054822
Linear 7.5Fr 34cc IAB	0684-00-0479-02	3000056937
Linear 7.5Fr 34cc IAB	0684-00-0479-02	3000059457
Linear 7.5Fr 34cc IAB	0684-00-0479-02	3000061523
Linear 7.5Fr 34cc IAB	0684-00-0479-02	3000065132
Linear 7.5Fr 34cc IAB	0684-00-0479-02	3000065723
Linear 7.5Fr 34cc IAB	0684-00-0479-02	3000066692
Linear 7.5Fr 34cc IAB	0684-00-0479-02	3000068136
Linear 7.5Fr 34cc IAB	0684-00-0479-02	3000068729
Linear 7.5Fr 34cc IAB	0684-00-0479-02	3000069123
Linear 7.5Fr 34cc IAB	0684-00-0479-02	3000070427
Linear 7.5Fr 34cc IAB	0684-00-0479-02	3000071290
Linear 7.5Fr 34cc IAB	0684-00-0479-02	3000071720
Linear 7.5Fr 34cc IAB	0684-00-0479-02	3000072523
Linear 7.5Fr 34cc IAB	0684-00-0479-02	3000074027
Linear 7.5Fr 34cc IAB	0684-00-0479-02	3000075916
Linear 7.5Fr 34cc IAB	0684-00-0479-02	3000075921
Linear 7.5Fr 34cc IAB	0684-00-0479-02	3000076853
Linear 7.5Fr 34cc IAB	0684-00-0479-02	3000079931
Linear 7.5Fr 34cc IAB	0684-00-0479-02	3000077809
Linear 7.5Fr 34cc IAB	0684-00-0479-02	3000079514
Linear 7.5Fr 34cc IAB	0684-00-0479-02	3000079515
Linear 7.5Fr 34cc IAB	0684-00-0479-02	3000080102
Linear 7.5Fr 34cc IAB	0684-00-0479-02	3000080731
Linear 7.5Fr 34cc IAB	0684-00-0479-02	3000083201
Linear 7.5Fr 34cc IAB	0684-00-0479-02	3000083202
Linear 7.5Fr 34cc IAB	0684-00-0479-02	3000085815
Linear 7.5Fr 34cc IAB	0684-00-0479-02	3000086856
Linear 7.5Fr 34cc IAB	0684-00-0479-02	3000086863
Linear 7.5Fr 34cc IAB	0684-00-0479-02	3000088104
Linear 7.5Fr 34cc IAB	0684-00-0479-02	3000088696
Linear 7.5Fr 34cc IAB	0684-00-0479-02	3000091478
Linear 7.5Fr 34cc IAB	0684-00-0479-02	3000091697
Linear 7.5Fr 34cc IAB	0684-00-0479-02	3000093943
Linear 7.5Fr 34cc IAB	0684-00-0479-02	3000097844
Linear 7.5Fr 34cc IAB	0684-00-0479-02	3000098672
Linear 7.5Fr 34cc IAB	0684-00-0479-02	3000101842
Linear 7.5Fr 34cc IAB	0684-00-0479-02	3000105735
Linear 7.5Fr 34cc IAB	0684-00-0479-02	3000108049

Linear 7.5Fr 34cc IAB	0684-00-0479-02	3000108048
Linear 7.5Fr 34cc IAB	0684-00-0479-02	3000113622
Linear 7.5Fr 34cc IAB	0684-00-0479-02	3000115918
Linear 7.5Fr 34cc IAB	0684-00-0479-02U	3000049792
Linear 7.5Fr 34cc IAB	0684-00-0479-02U	3000057344
Linear 7.5Fr 34cc IAB	0684-00-0479-02U	3000058579
Linear 7.5Fr 34cc IAB	0684-00-0479-02U	3000059459
Linear 7.5Fr 34cc IAB	0684-00-0479-02U	3000064514
Linear 7.5Fr 34cc IAB	0684-00-0479-02U	3000065131
Linear 7.5Fr 34cc IAB	0684-00-0479-02U	3000066767
Linear 7.5Fr 34cc IAB	0684-00-0479-02U	3000068137
Linear 7.5Fr 34cc IAB	0684-00-0479-02U	3000073521
Linear 7.5Fr 34cc IAB	0684-00-0479-02U	3000074028
Linear 7.5Fr 34cc IAB	0684-00-0479-02U	3000083199
Linear 7.5Fr 34cc IAB	0684-00-0479-02U	3000085816
Linear 7.5Fr 34cc IAB	0684-00-0479-02U	3000089283
Linear 7.5Fr 34cc IAB	0684-00-0479-02U	3000091479
Linear 7.5Fr 34cc IAB	0684-00-0479-02U	3000094938
Linear 7.5Fr 34cc IAB	0684-00-0479-02U	3000097938
Linear 7.5Fr 34cc IAB	0684-00-0479-02U	3000100471
Linear 7.5Fr 34cc IAB	0684-00-0479-02U	3000101844
Linear 7.5Fr 34cc IAB	0684-00-0479-02U	3000109451
Linear 7.5Fr 34cc IAB	0684-00-0479-07	3000066138
Linear 7.5Fr 34cc IAB	0684-00-0479-07	3000069441
Linear 7.5Fr 34cc IAB	0684-00-0479-07	3000069655
Linear 7.5Fr 34cc IAB	0684-00-0479-07	3000071127
Linear 7.5Fr 34cc IAB	0684-00-0479-07	3000072522
Linear 7.5Fr 34cc IAB	0684-00-0479-07	3000073968
Linear 7.5Fr 34cc IAB	0684-00-0479-07	3000079897
Linear 7.5Fr 34cc IAB	0684-00-0479-07	3000080464
Linear 7.5Fr 34cc IAB	0684-00-0479-07	3000082658
Linear 7.5Fr 34cc IAB	0684-00-0479-07	3000083203
Linear 7.5Fr 34cc IAB	0684-00-0479-07	3000084566
Linear 7.5Fr 34cc IAB	0684-00-0479-07	3000086902
Linear 7.5Fr 34cc IAB	0684-00-0479-07	3000101585
Linear 7.5Fr 34cc IAB	0684-00-0479-07	3000101586
Linear 7.5Fr 34cc IAB	0684-00-0479-07	3000103054
Linear 7.5Fr 34cc IAB	0684-00-0479-07	3000107908
Linear 7.5Fr 34cc IAB	0684-00-0479-07	3000110206
Linear 7.5Fr 34cc IAB	0684-00-0479-08	3000056261
Linear 7.5Fr 34cc IAB	0684-00-0479-08	3000056262
Linear 7.5Fr 34cc IAB	0684-00-0479-08	3000061585
Linear 7.5Fr 34cc IAB	0684-00-0479-08	3000076852
Linear 7.5Fr 34cc IAB	0684-00-0479-08	3000082252
Linear 7.5Fr 34cc IAB	0684-00-0479-08	3000100001
Linear 7.5Fr 34cc IAB	0684-00-0479-08	3000100147
Linear 7.5Fr 34cc IAB	0684-00-0479-08	3000101521
Linear 7.5Fr 34cc IAB	0684-00-0479-10	3000094153
Linear 7.5Fr 34cc IAB	0684-00-0479-10	3000094937
Linear 7.5Fr 34cc IAB	0684-00-0479-10	3000103155
Linear 7.5Fr 34cc IAB	0684-00-0479-10	3000103917

Linear 7.5Fr 34cc IAB	0684-00-0479-10	3000105386
Linear 7.5Fr 34cc IAB	0684-00-0479-10	3000106635
Linear 7.5Fr 40cc IAB	0684-00-0480-01	3000046733
Linear 7.5Fr 40cc IAB	0684-00-0480-01	3000047274
Linear 7.5Fr 40cc IAB	0684-00-0480-01	3000047604
Linear 7.5Fr 40cc IAB	0684-00-0480-01	3000048401
Linear 7.5Fr 40cc IAB	0684-00-0480-01	3000049310
Linear 7.5Fr 40cc IAB	0684-00-0480-01	3000051066
Linear 7.5Fr 40cc IAB	0684-00-0480-01	3000051995
Linear 7.5Fr 40cc IAB	0684-00-0480-01	3000053048
Linear 7.5Fr 40cc IAB	0684-00-0480-01	3000054067
Linear 7.5Fr 40cc IAB	0684-00-0480-01	3000054914
Linear 7.5Fr 40cc IAB	0684-00-0480-01	3000057090
Linear 7.5Fr 40cc IAB	0684-00-0480-01	3000058631
Linear 7.5Fr 40cc IAB	0684-00-0480-01	3000059460
Linear 7.5Fr 40cc IAB	0684-00-0480-01	3000060261
Linear 7.5Fr 40cc IAB	0684-00-0480-01	3000060731
Linear 7.5Fr 40cc IAB	0684-00-0480-01	3000060848
Linear 7.5Fr 40cc IAB	0684-00-0480-01	3000062743
Linear 7.5Fr 40cc IAB	0684-00-0480-01	3000064057
Linear 7.5Fr 40cc IAB	0684-00-0480-01	3000071409
Linear 7.5Fr 40cc IAB	0684-00-0480-01	3000065401
Linear 7.5Fr 40cc IAB	0684-00-0480-01	3000065942
Linear 7.5Fr 40cc IAB	0684-00-0480-01	3000066768
Linear 7.5Fr 40cc IAB	0684-00-0480-01	3000067500
Linear 7.5Fr 40cc IAB	0684-00-0480-01	3000068688
Linear 7.5Fr 40cc IAB	0684-00-0480-01	3000068689
Linear 7.5Fr 40cc IAB	0684-00-0480-01	3000070070
Linear 7.5Fr 40cc IAB	0684-00-0480-01	3000072315
Linear 7.5Fr 40cc IAB	0684-00-0480-01	3000072885
Linear 7.5Fr 40cc IAB	0684-00-0480-01	3000083860
Linear 7.5Fr 40cc IAB	0684-00-0480-01	3000075682
Linear 7.5Fr 40cc IAB	0684-00-0480-01	3000076775
Linear 7.5Fr 40cc IAB	0684-00-0480-01	3000079570
Linear 7.5Fr 40cc IAB	0684-00-0480-01	3000080053
Linear 7.5Fr 40cc IAB	0684-00-0480-01	3000077520
Linear 7.5Fr 40cc IAB	0684-00-0480-01	3000077941
Linear 7.5Fr 40cc IAB	0684-00-0480-01	3000078503
Linear 7.5Fr 40cc IAB	0684-00-0480-01	3000078600
Linear 7.5Fr 40cc IAB	0684-00-0480-01	3000078656
Linear 7.5Fr 40cc IAB	0684-00-0480-01	3000079242
Linear 7.5Fr 40cc IAB	0684-00-0480-01	3000080764
Linear 7.5Fr 40cc IAB	0684-00-0480-01	3000081888
Linear 7.5Fr 40cc IAB	0684-00-0480-01	3000082186
Linear 7.5Fr 40cc IAB	0684-00-0480-01	3000082885
Linear 7.5Fr 40cc IAB	0684-00-0480-01	3000083861
Linear 7.5Fr 40cc IAB	0684-00-0480-01	3000084565
Linear 7.5Fr 40cc IAB	0684-00-0480-01	3000084851
Linear 7.5Fr 40cc IAB	0684-00-0480-01	3000085193
Linear 7.5Fr 40cc IAB	0684-00-0480-01	3000087974
Linear 7.5Fr 40cc IAB	0684-00-0480-01	3000088594

Linear 7.5Fr 40cc IAB	0684-00-0480-01	3000088639
Linear 7.5Fr 40cc IAB	0684-00-0480-01	3000089730
Linear 7.5Fr 40cc IAB	0684-00-0480-01	3000090853
Linear 7.5Fr 40cc IAB	0684-00-0480-01	3000091126
Linear 7.5Fr 40cc IAB	0684-00-0480-01	3000091327
Linear 7.5Fr 40cc IAB	0684-00-0480-01	3000092503
Linear 7.5Fr 40cc IAB	0684-00-0480-01	3000093663
Linear 7.5Fr 40cc IAB	0684-00-0480-01	3000095290
Linear 7.5Fr 40cc IAB	0684-00-0480-01	3000095946
Linear 7.5Fr 40cc IAB	0684-00-0480-01	3000096738
Linear 7.5Fr 40cc IAB	0684-00-0480-01	3000097654
Linear 7.5Fr 40cc IAB	0684-00-0480-01	3000099417
Linear 7.5Fr 40cc IAB	0684-00-0480-01	3000100344
Linear 7.5Fr 40cc IAB	0684-00-0480-01	3000101784
Linear 7.5Fr 40cc IAB	0684-00-0480-01	3000103915
Linear 7.5Fr 40cc IAB	0684-00-0480-01	3000105898
Linear 7.5Fr 40cc IAB	0684-00-0480-01	3000106317
Linear 7.5Fr 40cc IAB	0684-00-0480-01	3000109026
Linear 7.5Fr 40cc IAB	0684-00-0480-01	3000111350
Linear 7.5Fr 40cc IAB	0684-00-0480-01	3000111973
Linear 7.5Fr 40cc IAB	0684-00-0480-01C	3000048959
Linear 7.5Fr 40cc IAB	0684-00-0480-01C	3000047353
Linear 7.5Fr 40cc IAB	0684-00-0480-01C	3000047706
Linear 7.5Fr 40cc IAB	0684-00-0480-01C	3000048026
Linear 7.5Fr 40cc IAB	0684-00-0480-01C	3000048961
Linear 7.5Fr 40cc IAB	0684-00-0480-01C	3000048960
Linear 7.5Fr 40cc IAB	0684-00-0480-01C	3000048962
Linear 7.5Fr 40cc IAB	0684-00-0480-01C	3000050258
Linear 7.5Fr 40cc IAB	0684-00-0480-01C	3000051647
Linear 7.5Fr 40cc IAB	0684-00-0480-01C	3000052948
Linear 7.5Fr 40cc IAB	0684-00-0480-01C	3000054290
Linear 7.5Fr 40cc IAB	0684-00-0480-01C	3000054291
Linear 7.5Fr 40cc IAB	0684-00-0480-01C	3000054967
Linear 7.5Fr 40cc IAB	0684-00-0480-01C	3000056074
Linear 7.5Fr 40cc IAB	0684-00-0480-01C	3000057384
Linear 7.5Fr 40cc IAB	0684-00-0480-01C	3000058896
Linear 7.5Fr 40cc IAB	0684-00-0480-01C	3000058897
Linear 7.5Fr 40cc IAB	0684-00-0480-01C	3000060262
Linear 7.5Fr 40cc IAB	0684-00-0480-01C	3000060733
Linear 7.5Fr 40cc IAB	0684-00-0480-01C	3000063082
Linear 7.5Fr 40cc IAB	0684-00-0480-01C	3000063083
Linear 7.5Fr 40cc IAB	0684-00-0480-01C	3000064046
Linear 7.5Fr 40cc IAB	0684-00-0480-01C	3000066141
Linear 7.5Fr 40cc IAB	0684-00-0480-01C	3000066142
Linear 7.5Fr 40cc IAB	0684-00-0480-01C	3000066284
Linear 7.5Fr 40cc IAB	0684-00-0480-01C	3000067688
Linear 7.5Fr 40cc IAB	0684-00-0480-01C	3000069656
Linear 7.5Fr 40cc IAB	0684-00-0480-01C	3000069657
Linear 7.5Fr 40cc IAB	0684-00-0480-01C	3000070408
Linear 7.5Fr 40cc IAB	0684-00-0480-01C	3000070664
Linear 7.5Fr 40cc IAB	0684-00-0480-01C	3000071058

Linear 7.5Fr 40cc IAB	0684-00-0480-01C	3000071414
Linear 7.5Fr 40cc IAB	0684-00-0480-01C	3000073065
Linear 7.5Fr 40cc IAB	0684-00-0480-01C	3000073403
Linear 7.5Fr 40cc IAB	0684-00-0480-01C	3000075684
Linear 7.5Fr 40cc IAB	0684-00-0480-01C	3000079569
Linear 7.5Fr 40cc IAB	0684-00-0480-01C	3000080103
Linear 7.5Fr 40cc IAB	0684-00-0480-01C	3000077273
Linear 7.5Fr 40cc IAB	0684-00-0480-01C	3000077519
Linear 7.5Fr 40cc IAB	0684-00-0480-01C	3000078683
Linear 7.5Fr 40cc IAB	0684-00-0480-01C	3000080104
Linear 7.5Fr 40cc IAB	0684-00-0480-01C	3000080443
Linear 7.5Fr 40cc IAB	0684-00-0480-01C	3000080782
Linear 7.5Fr 40cc IAB	0684-00-0480-01C	3000082659
Linear 7.5Fr 40cc IAB	0684-00-0480-01C	3000083600
Linear 7.5Fr 40cc IAB	0684-00-0480-01C	3000085381
Linear 7.5Fr 40cc IAB	0684-00-0480-01C	3000086267
Linear 7.5Fr 40cc IAB	0684-00-0480-01C	3000087315
Linear 7.5Fr 40cc IAB	0684-00-0480-01C	3000089459
Linear 7.5Fr 40cc IAB	0684-00-0480-01C	3000090849
Linear 7.5Fr 40cc IAB	0684-00-0480-01C	3000092396
Linear 7.5Fr 40cc IAB	0684-00-0480-01C	3000093113
Linear 7.5Fr 40cc IAB	0684-00-0480-01C	3000094005
Linear 7.5Fr 40cc IAB	0684-00-0480-01C	3000095038
Linear 7.5Fr 40cc IAB	0684-00-0480-01C	3000095945
Linear 7.5Fr 40cc IAB	0684-00-0480-01C	3000098249
Linear 7.5Fr 40cc IAB	0684-00-0480-01C	3000100378
Linear 7.5Fr 40cc IAB	0684-00-0480-01C	3000102536
Linear 7.5Fr 40cc IAB	0684-00-0480-01C	3000103919
Linear 7.5Fr 40cc IAB	0684-00-0480-01C	3000106447
Linear 7.5Fr 40cc IAB	0684-00-0480-01C	3000109023
Linear 7.5Fr 40cc IAB	0684-00-0480-01C	3000110738
Linear 7.5Fr 40cc IAB	0684-00-0480-01C	3000115116
Linear 7.5Fr 40cc IAB	0684-00-0480-01C	3000116665
Linear 7.5Fr 40cc IAB	0684-00-0480-01U	3000047198
Linear 7.5Fr 40cc IAB	0684-00-0480-01U	3000047708
Linear 7.5Fr 40cc IAB	0684-00-0480-01U	3000049208
Linear 7.5Fr 40cc IAB	0684-00-0480-01U	3000049723
Linear 7.5Fr 40cc IAB	0684-00-0480-01U	3000051122
Linear 7.5Fr 40cc IAB	0684-00-0480-01U	3000052154
Linear 7.5Fr 40cc IAB	0684-00-0480-01U	3000052155
Linear 7.5Fr 40cc IAB	0684-00-0480-01U	3000053152
Linear 7.5Fr 40cc IAB	0684-00-0480-01U	3000053792
Linear 7.5Fr 40cc IAB	0684-00-0480-01U	3000054880
Linear 7.5Fr 40cc IAB	0684-00-0480-01U	3000056191
Linear 7.5Fr 40cc IAB	0684-00-0480-01U	3000056936
Linear 7.5Fr 40cc IAB	0684-00-0480-01U	3000057345
Linear 7.5Fr 40cc IAB	0684-00-0480-01U	3000058847
Linear 7.5Fr 40cc IAB	0684-00-0480-01U	3000060260
Linear 7.5Fr 40cc IAB	0684-00-0480-01U	3000060587
Linear 7.5Fr 40cc IAB	0684-00-0480-01U	3000063258
Linear 7.5Fr 40cc IAB	0684-00-0480-01U	3000064097

Linear 7.5Fr 40cc IAB	0684-00-0480-01U	3000065030
Linear 7.5Fr 40cc IAB	0684-00-0480-01U	3000065853
Linear 7.5Fr 40cc IAB	0684-00-0480-01U	3000065943
Linear 7.5Fr 40cc IAB	0684-00-0480-01U	3000066282
Linear 7.5Fr 40cc IAB	0684-00-0480-01U	3000067497
Linear 7.5Fr 40cc IAB	0684-00-0480-01U	3000068770
Linear 7.5Fr 40cc IAB	0684-00-0480-01U	3000069317
Linear 7.5Fr 40cc IAB	0684-00-0480-01U	3000070409
Linear 7.5Fr 40cc IAB	0684-00-0480-01U	3000070410
Linear 7.5Fr 40cc IAB	0684-00-0480-01U	3000071410
Linear 7.5Fr 40cc IAB	0684-00-0480-01U	3000072886
Linear 7.5Fr 40cc IAB	0684-00-0480-01U	3000072916
Linear 7.5Fr 40cc IAB	0684-00-0480-01U	3000075681
Linear 7.5Fr 40cc IAB	0684-00-0480-01U	3000076776
Linear 7.5Fr 40cc IAB	0684-00-0480-01U	3000075131
Linear 7.5Fr 40cc IAB	0684-00-0480-01U	3000079466
Linear 7.5Fr 40cc IAB	0684-00-0480-01U	3000080054
Linear 7.5Fr 40cc IAB	0684-00-0480-01U	3000077522
Linear 7.5Fr 40cc IAB	0684-00-0480-01U	3000077523
Linear 7.5Fr 40cc IAB	0684-00-0480-01U	3000077524
Linear 7.5Fr 40cc IAB	0684-00-0480-01U	3000078051
Linear 7.5Fr 40cc IAB	0684-00-0480-01U	3000078684
Linear 7.5Fr 40cc IAB	0684-00-0480-01U	3000080537
Linear 7.5Fr 40cc IAB	0684-00-0480-01U	3000080783
Linear 7.5Fr 40cc IAB	0684-00-0480-01U	3000081157
Linear 7.5Fr 40cc IAB	0684-00-0480-01U	3000083094
Linear 7.5Fr 40cc IAB	0684-00-0480-01U	3000083092
Linear 7.5Fr 40cc IAB	0684-00-0480-01U	3000083093
Linear 7.5Fr 40cc IAB	0684-00-0480-01U	3000080538
Linear 7.5Fr 40cc IAB	0684-00-0480-01U	3000084692
Linear 7.5Fr 40cc IAB	0684-00-0480-01U	3000085501
Linear 7.5Fr 40cc IAB	0684-00-0480-01U	3000085502
Linear 7.5Fr 40cc IAB	0684-00-0480-01U	3000087316
Linear 7.5Fr 40cc IAB	0684-00-0480-01U	3000087317
Linear 7.5Fr 40cc IAB	0684-00-0480-01U	3000087975
Linear 7.5Fr 40cc IAB	0684-00-0480-01U	3000088620
Linear 7.5Fr 40cc IAB	0684-00-0480-01U	3000089731
Linear 7.5Fr 40cc IAB	0684-00-0480-01U	3000090848
Linear 7.5Fr 40cc IAB	0684-00-0480-01U	3000091184
Linear 7.5Fr 40cc IAB	0684-00-0480-01U	3000091328
Linear 7.5Fr 40cc IAB	0684-00-0480-01U	3000093664
Linear 7.5Fr 40cc IAB	0684-00-0480-01U	3000095289
Linear 7.5Fr 40cc IAB	0684-00-0480-01U	3000097650
Linear 7.5Fr 40cc IAB	0684-00-0480-01U	3000099419
Linear 7.5Fr 40cc IAB	0684-00-0480-01U	3000098615
Linear 7.5Fr 40cc IAB	0684-00-0480-01U	3000099420
Linear 7.5Fr 40cc IAB	0684-00-0480-01U	3000100266
Linear 7.5Fr 40cc IAB	0684-00-0480-01U	3000101785
Linear 7.5Fr 40cc IAB	0684-00-0480-01U	3000103920
Linear 7.5Fr 40cc IAB	0684-00-0480-01U	3000105895
Linear 7.5Fr 40cc IAB	0684-00-0480-01U	3000106449

Linear 7.5Fr 40cc IAB	0684-00-0480-01U	3000106450
Linear 7.5Fr 40cc IAB	0684-00-0480-01U	3000109025
Linear 7.5Fr 40cc IAB	0684-00-0480-01U	3000110414
Linear 7.5Fr 40cc IAB	0684-00-0480-01U	3000111358
Linear 7.5Fr 40cc IAB	0684-00-0480-01U	3000111359
Linear 7.5Fr 40cc IAB	0684-00-0480-01U	3000116581
Linear 7.5Fr 40cc IAB	0684-00-0480-02	3000046737
Linear 7.5Fr 40cc IAB	0684-00-0480-02	3000047275
Linear 7.5Fr 40cc IAB	0684-00-0480-02	3000047637
Linear 7.5Fr 40cc IAB	0684-00-0480-02	3000048223
Linear 7.5Fr 40cc IAB	0684-00-0480-02	3000048400
Linear 7.5Fr 40cc IAB	0684-00-0480-02	3000049421
Linear 7.5Fr 40cc IAB	0684-00-0480-02	3000049722
Linear 7.5Fr 40cc IAB	0684-00-0480-02	3000050669
Linear 7.5Fr 40cc IAB	0684-00-0480-02	3000051151
Linear 7.5Fr 40cc IAB	0684-00-0480-02	3000051152
Linear 7.5Fr 40cc IAB	0684-00-0480-02	3000052046
Linear 7.5Fr 40cc IAB	0684-00-0480-02	3000053049
Linear 7.5Fr 40cc IAB	0684-00-0480-02	3000054337
Linear 7.5Fr 40cc IAB	0684-00-0480-02	3000054915
Linear 7.5Fr 40cc IAB	0684-00-0480-02	3000057067
Linear 7.5Fr 40cc IAB	0684-00-0480-02	3000057385
Linear 7.5Fr 40cc IAB	0684-00-0480-02	3000057386
Linear 7.5Fr 40cc IAB	0684-00-0480-02	3000058630
Linear 7.5Fr 40cc IAB	0684-00-0480-02	3000060734
Linear 7.5Fr 40cc IAB	0684-00-0480-02	3000062028
Linear 7.5Fr 40cc IAB	0684-00-0480-02	3000062697
Linear 7.5Fr 40cc IAB	0684-00-0480-02	3000065235
Linear 7.5Fr 40cc IAB	0684-00-0480-02	3000065854
Linear 7.5Fr 40cc IAB	0684-00-0480-02	3000066769
Linear 7.5Fr 40cc IAB	0684-00-0480-02	3000067531
Linear 7.5Fr 40cc IAB	0684-00-0480-02	3000067532
Linear 7.5Fr 40cc IAB	0684-00-0480-02	3000067533
Linear 7.5Fr 40cc IAB	0684-00-0480-02	3000067891
Linear 7.5Fr 40cc IAB	0684-00-0480-02	3000068434
Linear 7.5Fr 40cc IAB	0684-00-0480-02	3000068731
Linear 7.5Fr 40cc IAB	0684-00-0480-02	3000069377
Linear 7.5Fr 40cc IAB	0684-00-0480-02	3000070489
Linear 7.5Fr 40cc IAB	0684-00-0480-02	3000071411
Linear 7.5Fr 40cc IAB	0684-00-0480-02	3000072230
Linear 7.5Fr 40cc IAB	0684-00-0480-02	3000072884
Linear 7.5Fr 40cc IAB	0684-00-0480-02	3000073275
Linear 7.5Fr 40cc IAB	0684-00-0480-02	3000075683
Linear 7.5Fr 40cc IAB	0684-00-0480-02	3000075675
Linear 7.5Fr 40cc IAB	0684-00-0480-02	3000075963
Linear 7.5Fr 40cc IAB	0684-00-0480-02	3000077521
Linear 7.5Fr 40cc IAB	0684-00-0480-02	3000077940
Linear 7.5Fr 40cc IAB	0684-00-0480-02	3000078502
Linear 7.5Fr 40cc IAB	0684-00-0480-02	3000078601
Linear 7.5Fr 40cc IAB	0684-00-0480-02	3000079315
Linear 7.5Fr 40cc IAB	0684-00-0480-02	3000079829

Linear 7.5Fr 40cc IAB	0684-00-0480-02	3000078934
Linear 7.5Fr 40cc IAB	0684-00-0480-02	3000080098
Linear 7.5Fr 40cc IAB	0684-00-0480-02	3000080780
Linear 7.5Fr 40cc IAB	0684-00-0480-02	3000082187
Linear 7.5Fr 40cc IAB	0684-00-0480-02	3000085290
Linear 7.5Fr 40cc IAB	0684-00-0480-02	3000087124
Linear 7.5Fr 40cc IAB	0684-00-0480-02	3000087976
Linear 7.5Fr 40cc IAB	0684-00-0480-02	3000088621
Linear 7.5Fr 40cc IAB	0684-00-0480-02	3000090854
Linear 7.5Fr 40cc IAB	0684-00-0480-02	3000091039
Linear 7.5Fr 40cc IAB	0684-00-0480-02	3000091040
Linear 7.5Fr 40cc IAB	0684-00-0480-02	3000091477
Linear 7.5Fr 40cc IAB	0684-00-0480-02	3000093115
Linear 7.5Fr 40cc IAB	0684-00-0480-02	3000095291
Linear 7.5Fr 40cc IAB	0684-00-0480-02	3000095796
Linear 7.5Fr 40cc IAB	0684-00-0480-02	3000097651
Linear 7.5Fr 40cc IAB	0684-00-0480-02	3000097652
Linear 7.5Fr 40cc IAB	0684-00-0480-02	3000098613
Linear 7.5Fr 40cc IAB	0684-00-0480-02	3000099418
Linear 7.5Fr 40cc IAB	0684-00-0480-02	3000100066
Linear 7.5Fr 40cc IAB	0684-00-0480-02	3000100346
Linear 7.5Fr 40cc IAB	0684-00-0480-02	3000101786
Linear 7.5Fr 40cc IAB	0684-00-0480-02	3000106316
Linear 7.5Fr 40cc IAB	0684-00-0480-02	3000108557
Linear 7.5Fr 40cc IAB	0684-00-0480-02	3000109401
Linear 7.5Fr 40cc IAB	0684-00-0480-02	3000110454
Linear 7.5Fr 40cc IAB	0684-00-0480-02	3000110735
Linear 7.5Fr 40cc IAB	0684-00-0480-02	3000115917
Linear 7.5Fr 40cc IAB	0684-00-0480-02C	3000097096
Linear 7.5Fr 40cc IAB	0684-00-0480-02U	3000046499
Linear 7.5Fr 40cc IAB	0684-00-0480-02U	3000047201
Linear 7.5Fr 40cc IAB	0684-00-0480-02U	3000049313
Linear 7.5Fr 40cc IAB	0684-00-0480-02U	3000054879
Linear 7.5Fr 40cc IAB	0684-00-0480-02U	3000056192
Linear 7.5Fr 40cc IAB	0684-00-0480-02U	3000058580
Linear 7.5Fr 40cc IAB	0684-00-0480-02U	3000062803
Linear 7.5Fr 40cc IAB	0684-00-0480-02U	3000065029
Linear 7.5Fr 40cc IAB	0684-00-0480-02U	3000065852
Linear 7.5Fr 40cc IAB	0684-00-0480-02U	3000066283
Linear 7.5Fr 40cc IAB	0684-00-0480-02U	3000069316
Linear 7.5Fr 40cc IAB	0684-00-0480-02U	3000070488
Linear 7.5Fr 40cc IAB	0684-00-0480-02U	3000072248
Linear 7.5Fr 40cc IAB	0684-00-0480-02U	3000076777
Linear 7.5Fr 40cc IAB	0684-00-0480-02U	3000075132
Linear 7.5Fr 40cc IAB	0684-00-0480-02U	3000079467
Linear 7.5Fr 40cc IAB	0684-00-0480-02U	3000079941
Linear 7.5Fr 40cc IAB	0684-00-0480-02U	3000078050
Linear 7.5Fr 40cc IAB	0684-00-0480-02U	3000080781
Linear 7.5Fr 40cc IAB	0684-00-0480-02U	3000085503
Linear 7.5Fr 40cc IAB	0684-00-0480-02U	3000088490
Linear 7.5Fr 40cc IAB	0684-00-0480-02U	3000089725

Linear 7.5Fr 40cc IAB	0684-00-0480-02U	3000090855
Linear 7.5Fr 40cc IAB	0684-00-0480-02U	3000093675
Linear 7.5Fr 40cc IAB	0684-00-0480-02U	3000097057
Linear 7.5Fr 40cc IAB	0684-00-0480-02U	3000100067
Linear 7.5Fr 40cc IAB	0684-00-0480-02U	3000101787
Linear 7.5Fr 40cc IAB	0684-00-0480-02U	3000105897
Linear 7.5Fr 40cc IAB	0684-00-0480-02U	3000109024
Linear 7.5Fr 40cc IAB	0684-00-0480-07	3000066139
Linear 7.5Fr 40cc IAB	0684-00-0480-07	3000067810
Linear 7.5Fr 40cc IAB	0684-00-0480-07	3000071126
Linear 7.5Fr 40cc IAB	0684-00-0480-07	3000072259
Linear 7.5Fr 40cc IAB	0684-00-0480-07	3000083859
Linear 7.5Fr 40cc IAB	0684-00-0480-07	3000077205
Linear 7.5Fr 40cc IAB	0684-00-0480-07	3000079828
Linear 7.5Fr 40cc IAB	0684-00-0480-07	3000082126
Linear 7.5Fr 40cc IAB	0684-00-0480-07	3000080463
Linear 7.5Fr 40cc IAB	0684-00-0480-07	3000084652
Linear 7.5Fr 40cc IAB	0684-00-0480-07	3000085504
Linear 7.5Fr 40cc IAB	0684-00-0480-07	3000087161
Linear 7.5Fr 40cc IAB	0684-00-0480-07	3000100856
Linear 7.5Fr 40cc IAB	0684-00-0480-07	3000103052
Linear 7.5Fr 40cc IAB	0684-00-0480-07	3000103053
Linear 7.5Fr 40cc IAB	0684-00-0480-07	3000107166
Linear 7.5Fr 40cc IAB	0684-00-0480-07	3000107281
Linear 7.5Fr 40cc IAB	0684-00-0480-07	3000107909
Linear 7.5Fr 40cc IAB	0684-00-0480-07	3000110205
Linear 7.5Fr 40cc IAB	0684-00-0480-08	3000046814
Linear 7.5Fr 40cc IAB	0684-00-0480-08	3000047059
Linear 7.5Fr 40cc IAB	0684-00-0480-08	3000056260
Linear 7.5Fr 40cc IAB	0684-00-0480-08	3000061584
Linear 7.5Fr 40cc IAB	0684-00-0480-08	3000076778
Linear 7.5Fr 40cc IAB	0684-00-0480-08	3000081156
Linear 7.5Fr 40cc IAB	0684-00-0480-08	3000100000
Linear 7.5Fr 40cc IAB	0684-00-0480-08	3000101129
Linear 7.5Fr 40cc IAB	0684-00-0480-08	3000101130
Linear 7.5Fr 40cc IAB	0684-00-0480-10	3000094152
Linear 7.5Fr 40cc IAB	0684-00-0480-10	3000095152
Linear 7.5Fr 40cc IAB	0684-00-0480-10	3000096000
Linear 7.5Fr 40cc IAB	0684-00-0480-10	3000097653
Linear 7.5Fr 40cc IAB	0684-00-0480-10	3000103055
Linear 7.5Fr 40cc IAB	0684-00-0480-10	3000103916
Linear 7.5Fr 40cc IAB	0684-00-0480-10	3000105384
Linear 7.5Fr 40cc IAB	0684-00-0480-10	3000106448
Linear 7.5Fr 40cc IAB	0684-00-0480-10	3000109400
Linear 7.5Fr 40cc IAB	0684-00-0480-10	3000116582
MEGA 8Fr 50cc IAB	0684-00-0498-01	3000047356
MEGA 8Fr 50cc IAB	0684-00-0498-01	3000048230
MEGA 8Fr 50cc IAB	0684-00-0498-01	3000049055
MEGA 8Fr 50cc IAB	0684-00-0498-01	3000051378
MEGA 8Fr 50cc IAB	0684-00-0498-01	3000054590
MEGA 8Fr 50cc IAB	0684-00-0498-01	3000056355

MEGA 8Fr 50cc IAB	0684-00-0498-01	3000056628
MEGA 8Fr 50cc IAB	0684-00-0498-01	3000059148
MEGA 8Fr 50cc IAB	0684-00-0498-01	3000060544
MEGA 8Fr 50cc IAB	0684-00-0498-01	3000068926
MEGA 8Fr 50cc IAB	0684-00-0498-01	3000068941
MEGA 8Fr 50cc IAB	0684-00-0498-01	3000073404
MEGA 8Fr 50cc IAB	0684-00-0498-01	3000074166
MEGA 8Fr 50cc IAB	0684-00-0498-01	3000074890
MEGA 8Fr 50cc IAB	0684-00-0498-01	3000076131
MEGA 8Fr 50cc IAB	0684-00-0498-01	3000076402
MEGA 8Fr 50cc IAB	0684-00-0498-01	3000078920
MEGA 8Fr 50cc IAB	0684-00-0498-01	3000080420
MEGA 8Fr 50cc IAB	0684-00-0498-01	3000081777
MEGA 8Fr 50cc IAB	0684-00-0498-01	3000083598
MEGA 8Fr 50cc IAB	0684-00-0498-01	3000085937
MEGA 8Fr 50cc IAB	0684-00-0498-01	3000090852
MEGA 8Fr 50cc IAB	0684-00-0498-01	3000089347
MEGA 8Fr 50cc IAB	0684-00-0498-01	3000091872
MEGA 8Fr 50cc IAB	0684-00-0498-01	3000094108
MEGA 8Fr 50cc IAB	0684-00-0498-01	3000094935
MEGA 8Fr 50cc IAB	0684-00-0498-01	3000096067
MEGA 8Fr 50cc IAB	0684-00-0498-01	3000095318
MEGA 8Fr 50cc IAB	0684-00-0498-01	3000096066
MEGA 8Fr 50cc IAB	0684-00-0498-01	3000097959
MEGA 8Fr 50cc IAB	0684-00-0498-01	3000100857
MEGA 8Fr 50cc IAB	0684-00-0498-01	3000100582
MEGA 8Fr 50cc IAB	0684-00-0498-01	3000102905
MEGA 8Fr 50cc IAB	0684-00-0498-01	3000105061
MEGA 8Fr 50cc IAB	0684-00-0498-01	3000107102
MEGA 8Fr 50cc IAB	0684-00-0498-01	3000108963
MEGA 8Fr 50cc IAB	0684-00-0498-01	3000111680
MEGA 8Fr 50cc IAB	0684-00-0498-01	3000113152
MEGA 8Fr 50cc IAB	0684-00-0498-07	3000068029
MEGA 8Fr 50cc IAB	0684-00-0498-07	3000071836
MEGA 8Fr 50cc IAB	0684-00-0498-07	3000074163
MEGA 8Fr 50cc IAB	0684-00-0498-07	3000103039
MEGA 8Fr 50cc IAB	0684-00-0498-07	3000106964
Trans-Ray 7Fr 34cc IAB	0684-00-0545-01	3000050506
Trans-Ray 7Fr 34cc IAB	0684-00-0545-01	3000050507
Trans-Ray 7Fr 34cc IAB	0684-00-0545-01	3000051382
Trans-Ray 7Fr 34cc IAB	0684-00-0545-01	3000053314
Trans-Ray 7Fr 34cc IAB	0684-00-0545-01	3000055312
Trans-Ray 7Fr 34cc IAB	0684-00-0545-01	3000062556
Trans-Ray 7Fr 34cc IAB	0684-00-0545-01	3000064523
Trans-Ray 7Fr 34cc IAB	0684-00-0545-01	3000066986
Trans-Ray 7Fr 34cc IAB	0684-00-0545-01	3000069732
Trans-Ray 7Fr 34cc IAB	0684-00-0545-01	3000077731
Trans-Ray 7Fr 34cc IAB	0684-00-0545-01	3000077732
Trans-Ray 7Fr 34cc IAB	0684-00-0545-01	3000080305
Trans-Ray 7Fr 34cc IAB	0684-00-0545-01	3000083332
Trans-Ray 7Fr 34cc IAB	0684-00-0545-01	3000083333

Trans-Ray 7Fr 34cc IAB	0684-00-0545-01	3000084968
Trans-Ray 7Fr 34cc IAB	0684-00-0545-01	3000087492
Trans-Ray 7Fr 34cc IAB	0684-00-0545-01	3000095469
Trans-Ray 7Fr 34cc IAB	0684-00-0545-01	3000096180
Trans-Ray 7Fr 34cc IAB	0684-00-0545-01	3000096181
Trans-Ray 7Fr 34cc IAB	0684-00-0545-01	3000096585
Trans-Ray 7Fr 34cc IAB	0684-00-0545-01	3000098052
Trans-Ray 7Fr 34cc IAB	0684-00-0545-01	3000103498
Trans-Ray 7Fr 34cc IAB	0684-00-0545-01	3000104616
Trans-Ray 7Fr 34cc IAB	0684-00-0545-01	3000104617
Trans-Ray 7Fr 34cc IAB	0684-00-0545-01	3000112275
Trans-Ray 7Fr 34cc IAB	0684-00-0545-01	3000113067
Trans-Ray 7Fr 40cc IAB	0684-00-0546-01	3000051381
Trans-Ray 7Fr 40cc IAB	0684-00-0546-01	3000053316
Trans-Ray 7Fr 40cc IAB	0684-00-0546-01	3000057969
Trans-Ray 7Fr 40cc IAB	0684-00-0546-01	3000062074
Trans-Ray 7Fr 40cc IAB	0684-00-0546-01	3000066987
Trans-Ray 7Fr 40cc IAB	0684-00-0546-01	3000062557
Trans-Ray 7Fr 40cc IAB	0684-00-0546-01	3000062558
Trans-Ray 7Fr 40cc IAB	0684-00-0546-01	3000066988
Trans-Ray 7Fr 40cc IAB	0684-00-0546-01	3000064524
Trans-Ray 7Fr 40cc IAB	0684-00-0546-01	3000066989
Trans-Ray 7Fr 40cc IAB	0684-00-0546-01	3000066990
Trans-Ray 7Fr 40cc IAB	0684-00-0546-01	3000067890
Trans-Ray 7Fr 40cc IAB	0684-00-0546-01	3000069733
Trans-Ray 7Fr 40cc IAB	0684-00-0546-01	3000074397
Trans-Ray 7Fr 40cc IAB	0684-00-0546-01	3000074398
Trans-Ray 7Fr 40cc IAB	0684-00-0546-01	3000077670
Trans-Ray 7Fr 40cc IAB	0684-00-0546-01	3000077671
Trans-Ray 7Fr 40cc IAB	0684-00-0546-01	3000078857
Trans-Ray 7Fr 40cc IAB	0684-00-0546-01	3000080306
Trans-Ray 7Fr 40cc IAB	0684-00-0546-01	3000080310
Trans-Ray 7Fr 40cc IAB	0684-00-0546-01	3000083334
Trans-Ray 7Fr 40cc IAB	0684-00-0546-01	3000083335
Trans-Ray 7Fr 40cc IAB	0684-00-0546-01	3000084272
Trans-Ray 7Fr 40cc IAB	0684-00-0546-01	3000084969
Trans-Ray 7Fr 40cc IAB	0684-00-0546-01	3000087490
Trans-Ray 7Fr 40cc IAB	0684-00-0546-01	3000087491
Trans-Ray 7Fr 40cc IAB	0684-00-0546-01	3000092779
Trans-Ray 7Fr 40cc IAB	0684-00-0546-01	3000092780
Trans-Ray 7Fr 40cc IAB	0684-00-0546-01	3000095470
Trans-Ray 7Fr 40cc IAB	0684-00-0546-01	3000096586
Trans-Ray 7Fr 40cc IAB	0684-00-0546-01	3000098053
Trans-Ray 7Fr 40cc IAB	0684-00-0546-01	3000098054
Trans-Ray 7Fr 40cc IAB	0684-00-0546-01	3000098055
Trans-Ray 7Fr 40cc IAB	0684-00-0546-01	3000098056
Trans-Ray 7Fr 40cc IAB	0684-00-0546-01	3000099590
Trans-Ray 7Fr 40cc IAB	0684-00-0546-01	3000099591
Trans-Ray 7Fr 40cc IAB	0684-00-0546-01	3000112273
Trans-Ray 7Fr 40cc IAB	0684-00-0546-01	3000107761
Trans-Ray 7Fr 40cc IAB	0684-00-0546-01	3000107762

Trans-Ray 7Fr 40cc IAB	0684-00-0546-01	3000112274
YAMATO PLUS-R 7.5Fr 30cc IAB	0684-00-0559-01	3000045349
YAMATO PLUS-R 7.5Fr 30cc IAB	0684-00-0559-01	3000051218
YAMATO PLUS-R 7.5Fr 30cc IAB	0684-00-0559-01	3000053318
YAMATO PLUS-R 7.5Fr 30cc IAB	0684-00-0559-01	3000055309
YAMATO PLUS-R 7.5Fr 30cc IAB	0684-00-0559-01	3000055310
YAMATO PLUS-R 7.5Fr 30cc IAB	0684-00-0559-01	3000055311
YAMATO PLUS-R 7.5Fr 30cc IAB	0684-00-0559-01	3000056703
YAMATO PLUS-R 7.5Fr 30cc IAB	0684-00-0559-01	3000057316
YAMATO PLUS-R 7.5Fr 30cc IAB	0684-00-0559-01	3000057317
YAMATO PLUS-R 7.5Fr 30cc IAB	0684-00-0559-01	3000057967
YAMATO PLUS-R 7.5Fr 30cc IAB	0684-00-0559-01	3000066213
YAMATO PLUS-R 7.5Fr 30cc IAB	0684-00-0559-01	3000066214
YAMATO PLUS-R 7.5Fr 30cc IAB	0684-00-0559-01	3000068109
YAMATO PLUS-R 7.5Fr 30cc IAB	0684-00-0559-01	3000068110
YAMATO PLUS-R 7.5Fr 30cc IAB	0684-00-0559-01	3000068111
YAMATO PLUS-R 7.5Fr 30cc IAB	0684-00-0559-01	3000069731
YAMATO PLUS-R 7.5Fr 30cc IAB	0684-00-0559-01	3000071412
YAMATO PLUS-R 7.5Fr 30cc IAB	0684-00-0559-01	3000072401
YAMATO PLUS-R 7.5Fr 30cc IAB	0684-00-0559-01	3000072402
YAMATO PLUS-R 7.5Fr 30cc IAB	0684-00-0559-01	3000073682
YAMATO PLUS-R 7.5Fr 30cc IAB	0684-00-0559-01	3000073816
YAMATO PLUS-R 7.5Fr 30cc IAB	0684-00-0559-01	3000073817
YAMATO PLUS-R 7.5Fr 30cc IAB	0684-00-0559-01	3000074887
YAMATO PLUS-R 7.5Fr 30cc IAB	0684-00-0559-01	3000083328
YAMATO PLUS-R 7.5Fr 30cc IAB	0684-00-0559-01	3000083329
YAMATO PLUS-R 7.5Fr 30cc IAB	0684-00-0559-01	3000087746
YAMATO PLUS-R 7.5Fr 30cc IAB	0684-00-0559-01	3000076462
YAMATO PLUS-R 7.5Fr 30cc IAB	0684-00-0559-01	3000078855
YAMATO PLUS-R 7.5Fr 30cc IAB	0684-00-0559-01	3000079319
YAMATO PLUS-R 7.5Fr 30cc IAB	0684-00-0559-01	3000079320
YAMATO PLUS-R 7.5Fr 30cc IAB	0684-00-0559-01	3000083326
YAMATO PLUS-R 7.5Fr 30cc IAB	0684-00-0559-01	3000080308
YAMATO PLUS-R 7.5Fr 30cc IAB	0684-00-0559-01	3000083327
YAMATO PLUS-R 7.5Fr 30cc IAB	0684-00-0559-01	3000087747
YAMATO PLUS-R 7.5Fr 30cc IAB	0684-00-0559-01	3000090530
YAMATO PLUS-R 7.5Fr 30cc IAB	0684-00-0559-01	3000090531
YAMATO PLUS-R 7.5Fr 30cc IAB	0684-00-0559-01	3000092773
YAMATO PLUS-R 7.5Fr 30cc IAB	0684-00-0559-01	3000096367
YAMATO PLUS-R 7.5Fr 30cc IAB	0684-00-0559-01	3000096368
YAMATO PLUS-R 7.5Fr 30cc IAB	0684-00-0559-01	3000095471
YAMATO PLUS-R 7.5Fr 30cc IAB	0684-00-0559-01	3000096369
YAMATO PLUS-R 7.5Fr 30cc IAB	0684-00-0559-01	3000096370
YAMATO PLUS-R 7.5Fr 30cc IAB	0684-00-0559-01	3000096371
YAMATO PLUS-R 7.5Fr 30cc IAB	0684-00-0559-01	3000096372
YAMATO PLUS-R 7.5Fr 30cc IAB	0684-00-0559-01	3000096373
YAMATO PLUS-R 7.5Fr 30cc IAB	0684-00-0559-01	3000096374
YAMATO PLUS-R 7.5Fr 30cc IAB	0684-00-0559-01	3000098749
YAMATO PLUS-R 7.5Fr 30cc IAB	0684-00-0559-01	3000098750
YAMATO PLUS-R 7.5Fr 30cc IAB	0684-00-0559-01	3000103503
YAMATO PLUS-R 7.5Fr 30cc IAB	0684-00-0559-01	3000103504

[illegible]

YAMATO PLUS-R 7.5Fr 35cc IAB	0684-00-0560-01	3000096378
YAMATO PLUS-R 7.5Fr 35cc IAB	0684-00-0560-01	3000098752
YAMATO PLUS-R 7.5Fr 35cc IAB	0684-00-0560-01	3000099588
YAMATO PLUS-R 7.5Fr 35cc IAB	0684-00-0560-01	3000099589
YAMATO PLUS-R 7.5Fr 40cc IAB	0684-00-0561-01	3000053317
YAMATO PLUS-R 7.5Fr 40cc IAB	0684-00-0561-01	3000056694
YAMATO PLUS-R 7.5Fr 40cc IAB	0684-00-0561-01	3000057318
YAMATO PLUS-R 7.5Fr 40cc IAB	0684-00-0561-01	3000057970
YAMATO PLUS-R 7.5Fr 40cc IAB	0684-00-0561-01	3000062073
YAMATO PLUS-R 7.5Fr 40cc IAB	0684-00-0561-01	3000062639
YAMATO PLUS-R 7.5Fr 40cc IAB	0684-00-0561-01	3000062640
YAMATO PLUS-R 7.5Fr 40cc IAB	0684-00-0561-01	3000068106
YAMATO PLUS-R 7.5Fr 40cc IAB	0684-00-0561-01	3000071287
YAMATO PLUS-R 7.5Fr 40cc IAB	0684-00-0561-01	3000072717
YAMATO PLUS-R 7.5Fr 40cc IAB	0684-00-0561-01	3000064723
YAMATO PLUS-R 7.5Fr 40cc IAB	0684-00-0561-01	3000072718
YAMATO PLUS-R 7.5Fr 40cc IAB	0684-00-0561-01	3000072719
YAMATO PLUS-R 7.5Fr 40cc IAB	0684-00-0561-01	3000074886
YAMATO PLUS-R 7.5Fr 40cc IAB	0684-00-0561-01	3000087744
YAMATO PLUS-R 7.5Fr 40cc IAB	0684-00-0561-01	3000090528
YAMATO PLUS-R 7.5Fr 40cc IAB	0684-00-0561-01	3000076632
YAMATO PLUS-R 7.5Fr 40cc IAB	0684-00-0561-01	3000078856
YAMATO PLUS-R 7.5Fr 40cc IAB	0684-00-0561-01	3000080307
YAMATO PLUS-R 7.5Fr 40cc IAB	0684-00-0561-01	3000083331
YAMATO PLUS-R 7.5Fr 40cc IAB	0684-00-0561-01	3000084962
YAMATO PLUS-R 7.5Fr 40cc IAB	0684-00-0561-01	3000092774
YAMATO PLUS-R 7.5Fr 40cc IAB	0684-00-0561-01	3000092775
YAMATO PLUS-R 7.5Fr 40cc IAB	0684-00-0561-01	3000095473
YAMATO PLUS-R 7.5Fr 40cc IAB	0684-00-0561-01	3000096376
YAMATO PLUS-R 7.5Fr 40cc IAB	0684-00-0561-01	3000096606
YAMATO PLUS-R 7.5Fr 40cc IAB	0684-00-0561-01	3000096379
YAMATO PLUS-R 7.5Fr 40cc IAB	0684-00-0561-01	3000098751
YAMATO PLUS-R 7.5Fr 40cc IAB	0684-00-0561-01	3000104521
YAMATO PLUS-R 7.5Fr 40cc IAB	0684-00-0561-01	3000115011
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01	3000047947
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01	3000047946
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01	3000047814
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01	3000048703
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01	3000049591
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01	3000050881
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01	3000051546
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01	3000051547
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01	3000053550
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01	3000054595
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01	3000054594
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01	3000055238
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01	3000056435
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01	3000057510
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01	3000057511
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01	3000058169
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01	3000058170

Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01	3000061118
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01	3000062404
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01	3000063262
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01	3000066059
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01	3000066060
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01	3000066061
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01	3000066062
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01	3000067530
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01	3000068522
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01	3000069870
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01	3000072066
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01	3000075678
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01	3000076244
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01	3000077031
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01	3000077765
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01	3000078404
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01	3000079679
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01	3000082472
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01	3000084069
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01	3000085231
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01	3000086994
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01	3000087688
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01	3000091205
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01	3000091206
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01	3000091828
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01	3000092976
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01	3000093557
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01	3000094739
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01	3000095037
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01	3000100585
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01	3000105260
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01	3000106962
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01	3000108882
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01	3000108883
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01	3000113809
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01	3000115337
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01U	3000046417
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01U	3000047163
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01U	3000047815
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01U	3000047818
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01U	3000049592
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01U	3000050882
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01U	3000051544
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01U	3000051545
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01U	3000053551
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01U	3000054593
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01U	3000055179
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01U	3000055239
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01U	3000056434
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01U	3000057382
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01U	3000057546

[illegible]

Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01U	3000094740
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01U	3000095150
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01U	3000095151
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01U	3000096363
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01U	3000097058
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01U	3000097059
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01U	3000098713
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01U	3000099522
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01U	3000100581
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01U	3000102215
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01U	3000102248
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01U	3000102249
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01U	3000103405
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01U	3000103406
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01U	3000105261
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01U	3000105262
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01U	3000106963
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01U	3000107331
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01U	3000107332
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01U	3000108589
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01U	3000108881
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01U	3000110520
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01U	3000111682
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-01U	3000113810
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-05	3000046748
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-09	3000102329
Sensation Plus 7.5 Fr 40cc IAB	0684-00-0568-09	3000103404
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01	3000048918
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01	3000050497
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01	3000050612
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01	3000051521
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01	3000052324
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01	3000052325
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01	3000053151
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01	3000054475
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01	3000054968
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01	3000057034
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01	3000057505
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01	3000057506
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01	3000057507
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01	3000058278
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01	3000061119
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01	3000063261
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01	3000065028
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01	3000065134
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01	3000065565
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01	3000066085
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01	3000066143
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01	3000067135
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01	3000067888
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01	3000067892

Sensation Plus 8Fr 50cc IAB	0684-00-0576-01	3000069378
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01	3000069364
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01	3000070068
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01	3000070921
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01	3000070922
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01	3000072316
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01	3000073278
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01	3000074484
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01	3000075427
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01	3000075677
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01	3000075964
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01	3000077095
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01	3000078049
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01	3000078764
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01	3000079621
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01	3000080191
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01	3000081446
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01	3000082347
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01	3000082996
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01	3000084154
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01	3000083510
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01	3000084917
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01	3000085293
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01	3000086993
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01	3000088593
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01	3000091236
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01	3000087370
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01	3000087687
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01	3000087686
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01	3000088899
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01	3000091235
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01	3000091986
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01	3000091234
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01	3000091987
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01	3000093061
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01	3000095036
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01	3000096972
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01	3000096973
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01	3000099219
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01	3000101990
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01	3000102705
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01	3000103207
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01	3000104387
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01	3000106150
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01	3000108350
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01	3000108591
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01	3000110954
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01	3000112931
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01	3000113623
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01U	3000047054
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01U	3000047058

[illegible]

[illegible]

Sensation Plus 8Fr 50cc IAB	0684-00-0576-01U	3000102703
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01U	3000102704
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01U	3000103208
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01U	3000104808
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01U	3000106152
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01U	3000108349
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01U	3000108590
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01U	3000108682
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01U	3000108683
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01U	3000109819
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01U	3000110955
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01U	3000111683
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01U	3000111684
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01U	3000112929
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01U	3000112930
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01U	3000113624
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01U	3000113625
Sensation Plus 8Fr 50cc IAB	0684-00-0576-01U	3000117020
Sensation Plus 8Fr 50cc IAB	0684-00-0576-05	3000046749
Sensation Plus 8Fr 50cc IAB	0684-00-0576-09	3000094529
Sensation Plus 8Fr 50cc IAB	0684-00-0576-09	3000094530
Sensation Plus 8Fr 50cc IAB	0684-00-0576-09	3000103206
Sensation Plus 8Fr 50cc IAB	0684-00-0576-09	3000104388
Sensation Plus 8Fr 50cc IAB	0684-00-0576-09	3000108351
Sensation Plus 8Fr 50cc IAB	0684-00-0576-09	3000112932
TRANS-RAY PLUS 7.5Fr 35cc IAB	0684-00-0607	3000047804
TRANS-RAY PLUS 7.5Fr 35cc IAB	0684-00-0607	3000053322
TRANS-RAY PLUS 7.5Fr 35cc IAB	0684-00-0607	3000053323
TRANS-RAY PLUS 7.5Fr 35cc IAB	0684-00-0607	3000060107
TRANS-RAY PLUS 7.5Fr 35cc IAB	0684-00-0607	3000060108
TRANS-RAY PLUS 7.5Fr 35cc IAB	0684-00-0607	3000062255
TRANS-RAY PLUS 7.5Fr 35cc IAB	0684-00-0607	3000067698
TRANS-RAY PLUS 7.5Fr 35cc IAB	0684-00-0607	3000069727
TRANS-RAY PLUS 7.5Fr 35cc IAB	0684-00-0607	3000071976
TRANS-RAY PLUS 7.5Fr 35cc IAB	0684-00-0607	3000073186
TRANS-RAY PLUS 7.5Fr 35cc IAB	0684-00-0607	3000073680
TRANS-RAY PLUS 7.5Fr 35cc IAB	0684-00-0607	3000073681
TRANS-RAY PLUS 7.5Fr 35cc IAB	0684-00-0607	3000074885
TRANS-RAY PLUS 7.5Fr 35cc IAB	0684-00-0607	3000076710
TRANS-RAY PLUS 7.5Fr 35cc IAB	0684-00-0607	3000076711
TRANS-RAY PLUS 7.5Fr 35cc IAB	0684-00-0607	3000077752
TRANS-RAY PLUS 7.5Fr 35cc IAB	0684-00-0607	3000077753
TRANS-RAY PLUS 7.5Fr 35cc IAB	0684-00-0607	3000079103
TRANS-RAY PLUS 7.5Fr 35cc IAB	0684-00-0607	3000079104
TRANS-RAY PLUS 7.5Fr 35cc IAB	0684-00-0607	3000080873
TRANS-RAY PLUS 7.5Fr 35cc IAB	0684-00-0607	3000080872
TRANS-RAY PLUS 7.5Fr 35cc IAB	0684-00-0607	3000081671
TRANS-RAY PLUS 7.5Fr 35cc IAB	0684-00-0607	3000084267
TRANS-RAY PLUS 7.5Fr 35cc IAB	0684-00-0607	3000084966
TRANS-RAY PLUS 7.5Fr 35cc IAB	0684-00-0607	3000084965
TRANS-RAY PLUS 7.5Fr 35cc IAB	0684-00-0607	3000086799

TRANS-RAY PLUS 7.5Fr 35cc IAB	0684-00-0607	3000088548
TRANS-RAY PLUS 7.5Fr 35cc IAB	0684-00-0607	3000088371
TRANS-RAY PLUS 7.5Fr 35cc IAB	0684-00-0607	3000090599
TRANS-RAY PLUS 7.5Fr 35cc IAB	0684-00-0607	3000090600
TRANS-RAY PLUS 7.5Fr 35cc IAB	0684-00-0607	3000090601
TRANS-RAY PLUS 7.5Fr 35cc IAB	0684-00-0607	3000090602
TRANS-RAY PLUS 7.5Fr 35cc IAB	0684-00-0607	3000090603
TRANS-RAY PLUS 7.5Fr 35cc IAB	0684-00-0607	3000092664
TRANS-RAY PLUS 7.5Fr 35cc IAB	0684-00-0607	3000092663
TRANS-RAY PLUS 7.5Fr 35cc IAB	0684-00-0607	3000092665
TRANS-RAY PLUS 7.5Fr 35cc IAB	0684-00-0607	3000096263
TRANS-RAY PLUS 7.5Fr 35cc IAB	0684-00-0607	3000095467
TRANS-RAY PLUS 7.5Fr 35cc IAB	0684-00-0607	3000096264
TRANS-RAY PLUS 7.5Fr 35cc IAB	0684-00-0607	3000096265
TRANS-RAY PLUS 7.5Fr 35cc IAB	0684-00-0607	3000096266
TRANS-RAY PLUS 7.5Fr 35cc IAB	0684-00-0607	3000098176
TRANS-RAY PLUS 7.5Fr 35cc IAB	0684-00-0607	3000098177
TRANS-RAY PLUS 7.5Fr 35cc IAB	0684-00-0607	3000102408
TRANS-RAY PLUS 7.5Fr 35cc IAB	0684-00-0607	3000102409
TRANS-RAY PLUS 7.5Fr 35cc IAB	0684-00-0607	3000103469
TRANS-RAY PLUS 7.5Fr 35cc IAB	0684-00-0607	3000104520
TRANS-RAY PLUS 7.5Fr 35cc IAB	0684-00-0607	3000104615
TRANS-RAY PLUS 7.5Fr 35cc IAB	0684-00-0607	3000106012
TRANS-RAY PLUS 7.5Fr 35cc IAB	0684-00-0607	3000107758
TRANS-RAY PLUS 7.5Fr 35cc IAB	0684-00-0607	3000107759
TRANS-RAY PLUS 7.5Fr 35cc IAB	0684-00-0607	3000110626
TRANS-RAY PLUS 7.5Fr 35cc IAB	0684-00-0607	3000113155
TRANS-RAY PLUS 7.5Fr 35cc IAB	0684-00-0607	3000112276
TRANS-RAY PLUS 7.5Fr 35cc IAB	0684-00-0607	3000113156
TRANS-RAY PLUS 7.5Fr 40cc IAB	0684-00-0608	3000045568
TRANS-RAY PLUS 7.5Fr 40cc IAB	0684-00-0608	3000047929
TRANS-RAY PLUS 7.5Fr 40cc IAB	0684-00-0608	3000047928
TRANS-RAY PLUS 7.5Fr 40cc IAB	0684-00-0608	3000047930
TRANS-RAY PLUS 7.5Fr 40cc IAB	0684-00-0608	3000053321
TRANS-RAY PLUS 7.5Fr 40cc IAB	0684-00-0608	3000060109
TRANS-RAY PLUS 7.5Fr 40cc IAB	0684-00-0608	3000060110
TRANS-RAY PLUS 7.5Fr 40cc IAB	0684-00-0608	3000067699
TRANS-RAY PLUS 7.5Fr 40cc IAB	0684-00-0608	3000062192
TRANS-RAY PLUS 7.5Fr 40cc IAB	0684-00-0608	3000067700
TRANS-RAY PLUS 7.5Fr 40cc IAB	0684-00-0608	3000067701
TRANS-RAY PLUS 7.5Fr 40cc IAB	0684-00-0608	3000069728
TRANS-RAY PLUS 7.5Fr 40cc IAB	0684-00-0608	3000069729
TRANS-RAY PLUS 7.5Fr 40cc IAB	0684-00-0608	3000069730
TRANS-RAY PLUS 7.5Fr 40cc IAB	0684-00-0608	3000071288
TRANS-RAY PLUS 7.5Fr 40cc IAB	0684-00-0608	3000073187
TRANS-RAY PLUS 7.5Fr 40cc IAB	0684-00-0608	3000076242
TRANS-RAY PLUS 7.5Fr 40cc IAB	0684-00-0608	3000077669
TRANS-RAY PLUS 7.5Fr 40cc IAB	0684-00-0608	3000079105
TRANS-RAY PLUS 7.5Fr 40cc IAB	0684-00-0608	3000079318
TRANS-RAY PLUS 7.5Fr 40cc IAB	0684-00-0608	3000080309
TRANS-RAY PLUS 7.5Fr 40cc IAB	0684-00-0608	3000080876

TRANS-RAY PLUS 7.5Fr 40cc IAB	0684-00-0608	3000084268
TRANS-RAY PLUS 7.5Fr 40cc IAB	0684-00-0608	3000084967
TRANS-RAY PLUS 7.5Fr 40cc IAB	0684-00-0608	3000084960
TRANS-RAY PLUS 7.5Fr 40cc IAB	0684-00-0608	3000084961
TRANS-RAY PLUS 7.5Fr 40cc IAB	0684-00-0608	3000086800
TRANS-RAY PLUS 7.5Fr 40cc IAB	0684-00-0608	3000088372
TRANS-RAY PLUS 7.5Fr 40cc IAB	0684-00-0608	3000090597
TRANS-RAY PLUS 7.5Fr 40cc IAB	0684-00-0608	3000090598
TRANS-RAY PLUS 7.5Fr 40cc IAB	0684-00-0608	3000092668
TRANS-RAY PLUS 7.5Fr 40cc IAB	0684-00-0608	3000092666
TRANS-RAY PLUS 7.5Fr 40cc IAB	0684-00-0608	3000092667
TRANS-RAY PLUS 7.5Fr 40cc IAB	0684-00-0608	3000092669
TRANS-RAY PLUS 7.5Fr 40cc IAB	0684-00-0608	3000096268
TRANS-RAY PLUS 7.5Fr 40cc IAB	0684-00-0608	3000095468
TRANS-RAY PLUS 7.5Fr 40cc IAB	0684-00-0608	3000096267
TRANS-RAY PLUS 7.5Fr 40cc IAB	0684-00-0608	3000096269
TRANS-RAY PLUS 7.5Fr 40cc IAB	0684-00-0608	3000096270
TRANS-RAY PLUS 7.5Fr 40cc IAB	0684-00-0608	3000096271
TRANS-RAY PLUS 7.5Fr 40cc IAB	0684-00-0608	3000096587
TRANS-RAY PLUS 7.5Fr 40cc IAB	0684-00-0608	3000098175
TRANS-RAY PLUS 7.5Fr 40cc IAB	0684-00-0608	3000098747
TRANS-RAY PLUS 7.5Fr 40cc IAB	0684-00-0608	3000098748
TRANS-RAY PLUS 7.5Fr 40cc IAB	0684-00-0608	3000102410
TRANS-RAY PLUS 7.5Fr 40cc IAB	0684-00-0608	3000102411
TRANS-RAY PLUS 7.5Fr 40cc IAB	0684-00-0608	3000103470
TRANS-RAY PLUS 7.5Fr 40cc IAB	0684-00-0608	3000104519
TRANS-RAY PLUS 7.5Fr 40cc IAB	0684-00-0608	3000106011
TRANS-RAY PLUS 7.5Fr 40cc IAB	0684-00-0608	3000107760
TRANS-RAY PLUS 7.5Fr 40cc IAB	0684-00-0608	3000110622
TRANS-RAY PLUS 7.5Fr 40cc IAB	0684-00-0608	3000110623
TRANS-RAY PLUS 7.5Fr 40cc IAB	0684-00-0608	3000112744
TRANS-RAY PLUS 7.5Fr 40cc IAB	0684-00-0608	3000113154

Master IAB Batch (Component level, for reference only)	Master IAB Serial Numbers- affected, to screen for	Finished Good Part Number (IAB Kit PN)	Finished Good Batch Number (IAB Kit lot number)
3000043852	30000438520131	D684-00-0561-01	3000053317
3000043852	30000438520131	D684-00-0561-01	3000056694
3000043852	30000438520130	D684-00-0561-01	3000053317
3000043852	30000438520130	D684-00-0561-01	3000056694
3000043852	30000438520134	D684-00-0561-01	3000053317
3000043852	30000438520134	D684-00-0561-01	3000056694
3000043852	30000438520132	D684-00-0561-01	3000053317
3000043852	30000438520132	D684-00-0561-01	3000056694
3000043852	30000438520138	D684-00-0561-01	3000053317
3000043852	30000438520138	D684-00-0561-01	3000056694
3000043969	30000439690418	D684-00-0608	3000045568
3000043969	30000439690418	D684-00-0608	3000047929
3000043969	30000439690422	D684-00-0608	3000045568
3000043969	30000439690422	D684-00-0608	3000047929
3000043969	30000439690420	D684-00-0608	3000045568
3000043969	30000439690420	D684-00-0608	3000047929
3000043969	30000439690426	D684-00-0608	3000045568
3000043969	30000439690426	D684-00-0608	3000047929
3000043969	30000439690424	D684-00-0608	3000045568
3000043969	30000439690424	D684-00-0608	3000047929
3000043970	30000439700268	D684-00-0560-01	3000051211
3000043970	30000439700268	D684-00-0560-01	3000051213
3000043970	30000439700268	D684-00-0560-01	3000051214
3000043970	30000439700268	D684-00-0560-01	3000051215
3000043970	30000439700268	D684-00-0560-01	3000051216
3000043970	30000439700268	D684-00-0560-01	3000051217
3000043970	30000439700268	D684-00-0560-01	3000055314
3000043970	30000439700268	D684-00-0560-01	3000056699
3000043970	30000439700263	D684-00-0560-01	3000051211
3000043970	30000439700263	D684-00-0560-01	3000051213
3000043970	30000439700263	D684-00-0560-01	3000051214
3000043970	30000439700263	D684-00-0560-01	3000051215
3000043970	30000439700263	D684-00-0560-01	3000051216
3000043970	30000439700263	D684-00-0560-01	3000051217
3000043970	30000439700263	D684-00-0560-01	3000055314
3000043970	30000439700263	D684-00-0560-01	3000056699
3000043970	30000439700270	D684-00-0560-01	3000051211
3000043970	30000439700270	D684-00-0560-01	3000051213
3000043970	30000439700270	D684-00-0560-01	3000051214
3000043970	30000439700270	D684-00-0560-01	3000051215
3000043970	30000439700270	D684-00-0560-01	3000051216
3000043970	30000439700270	D684-00-0560-01	3000051217
3000043970	30000439700270	D684-00-0560-01	3000055314
3000043970	30000439700270	D684-00-0560-01	3000056699
3000043970	30000439700264	D684-00-0560-01	3000051211
3000043970	30000439700264	D684-00-0560-01	3000051213
3000043970	30000439700264	D684-00-0560-01	3000051214
3000043970	30000439700264	D684-00-0560-01	3000051215
3000043970	30000439700264	D684-00-0560-01	3000051216
3000043970	30000439700264	D684-00-0560-01	3000051217
3000043970	30000439700264	D684-00-0560-01	3000055314
3000043970	30000439700264	D684-00-0560-01	3000056699
3000043970	30000439700267	D684-00-0560-01	3000051211
3000043970	30000439700267	D684-00-0560-01	3000051213
3000043970	30000439700267	D684-00-0560-01	3000051214

[illegible]

[illegible]

[illegible]

3000044130	30000441300110	D684-00-0568-01	3000047947
3000044130	30000441300111	D684-00-0568-01U	3000046417
3000044130	30000441300111	D684-00-0568-05	3000046748
3000044130	30000441300111	D684-00-0568-01U	3000047163
3000044130	30000441300111	D684-00-0568-01U	3000047163
3000044130	30000441300111	D684-00-0568-01	3000047947
3000044130	30000441300116	D684-00-0568-01U	3000046417
3000044130	30000441300116	D684-00-0568-05	3000046748
3000044130	30000441300116	D684-00-0568-01U	3000047163
3000044130	30000441300116	D684-00-0568-01U	3000047163
3000044130	30000441300116	D684-00-0568-01	3000047947
3000044130	30000441300117	D684-00-0568-01U	3000046417
3000044130	30000441300117	D684-00-0568-05	3000046748
3000044130	30000441300117	D684-00-0568-01U	3000047163
3000044130	30000441300117	D684-00-0568-01U	3000047163
3000044130	30000441300117	D684-00-0568-01	3000047947
3000044560	30000445600217	D684-00-0546-01	3000051381
3000044560	30000445600217	D684-00-0546-01	3000053316
3000044560	30000445600219	D684-00-0546-01	3000051381
3000044560	30000445600219	D684-00-0546-01	3000053316
3000044560	30000445600221	D684-00-0546-01	3000051381
3000044560	30000445600221	D684-00-0546-01	3000053316
3000044560	30000445600222	D684-00-0546-01	3000051381
3000044560	30000445600222	D684-00-0546-01	3000053316
3000044560	30000445600238	D684-00-0546-01	3000051381
3000044560	30000445600238	D684-00-0546-01	3000053316
3000044561	30000445610667	D684-00-0470-01U	3000047053
3000044561	30000445610667	D684-00-0470-01	3000047161
3000044561	30000445610667	D684-00-0470-01	3000047162
3000044561	30000445610667	D684-00-0470-01U	3000047709
3000044561	30000445610659	D684-00-0470-01U	3000047053
3000044561	30000445610659	D684-00-0470-01	3000047161
3000044561	30000445610659	D684-00-0470-01	3000047162
3000044561	30000445610659	D684-00-0470-01U	3000047709
3000044561	30000445610670	D684-00-0470-01U	3000047053
3000044561	30000445610670	D684-00-0470-01	3000047161
3000044561	30000445610670	D684-00-0470-01	3000047162
3000044561	30000445610670	D684-00-0470-01U	3000047709
3000044561	30000445610675	D684-00-0470-01U	3000047053
3000044561	30000445610675	D684-00-0470-01	3000047161
3000044561	30000445610675	D684-00-0470-01	3000047162
3000044561	30000445610675	D684-00-0470-01U	3000047709
3000044561	30000445610677	D684-00-0470-01U	3000047053
3000044561	30000445610677	D684-00-0470-01	3000047161
3000044561	30000445610677	D684-00-0470-01	3000047162
3000044561	30000445610677	D684-00-0470-01U	3000047709
3000044562	30000445620273	D684-00-0480-02	3000047275
3000044562	30000445620273	D684-00-0480-01C	3000047353
3000044562	30000445620273	D684-00-0480-01	3000047604
3000044562	30000445620273	D684-00-0480-01C	3000047706
3000044562	30000445620273	D684-00-0480-01C	3000048026
3000044562	30000445620273	D684-00-0480-01C	3000048961
3000044562	30000445620277	D684-00-0480-02	3000047275
3000044562	30000445620277	D684-00-0480-01C	3000047353
3000044562	30000445620277	D684-00-0480-01	3000047604
3000044562	30000445620277	D684-00-0480-01C	3000047706
3000044562	30000445620277	D684-00-0480-01C	3000048026
3000044562	30000445620277	D684-00-0480-01C	3000048961

3000044562	30000445620283	D684-00-0480-02	3000047275
3000044562	30000445620283	D684-00-0480-01C	3000047353
3000044562	30000445620283	D684-00-0480-01	3000047604
3000044562	30000445620283	D684-00-0480-01C	3000047706
3000044562	30000445620283	D684-00-0480-01C	3000048026
3000044562	30000445620283	D684-00-0480-01C	3000048961
3000044562	30000445620288	D684-00-0480-02	3000047275
3000044562	30000445620288	D684-00-0480-01C	3000047353
3000044562	30000445620288	D684-00-0480-01	3000047604
3000044562	30000445620288	D684-00-0480-01C	3000047706
3000044562	30000445620288	D684-00-0480-01C	3000048026
3000044562	30000445620288	D684-00-0480-01C	3000048961
3000044562	30000445620293	D684-00-0480-02	3000047275
3000044562	30000445620293	D684-00-0480-01C	3000047353
3000044562	30000445620293	D684-00-0480-01	3000047604
3000044562	30000445620293	D684-00-0480-01C	3000047706
3000044562	30000445620293	D684-00-0480-01C	3000048026
3000044562	30000445620293	D684-00-0480-01C	3000048961
3000044882	30000448820069	D684-00-0568-01	3000047946
3000044882	30000448820070	D684-00-0568-01	3000047946
3000044882	30000448820071	D684-00-0568-01	3000047946
3000044882	30000448820074	D684-00-0568-01	3000047946
3000044882	30000448820072	D684-00-0568-01	3000047946
3000044892	30000448920024	D684-00-0479-02	3000048614
3000044892	30000448920024	D684-00-0479-01	3000048711
3000044892	30000448920024	D684-00-0479-01U	3000048818
3000044892	30000448920024	D684-00-0479-01C	3000048963
3000044892	30000448920024	D684-00-0479-02	3000049263
3000044892	30000448920020	D684-00-0479-02	3000048614
3000044892	30000448920020	D684-00-0479-01	3000048711
3000044892	30000448920020	D684-00-0479-01U	3000048818
3000044892	30000448920020	D684-00-0479-01C	3000048963
3000044892	30000448920020	D684-00-0479-02	3000049263
3000044892	30000448920014	D684-00-0479-02	3000048614
3000044892	30000448920014	D684-00-0479-01	3000048711
3000044892	30000448920014	D684-00-0479-01U	3000048818
3000044892	30000448920014	D684-00-0479-01C	3000048963
3000044892	30000448920014	D684-00-0479-02	3000049263
3000044892	30000448920001	D684-00-0479-02	3000048614
3000044892	30000448920001	D684-00-0479-01	3000048711
3000044892	30000448920001	D684-00-0479-01U	3000048818
3000044892	30000448920001	D684-00-0479-01C	3000048963
3000044892	30000448920001	D684-00-0479-02	3000049263
3000044892	30000448920008	D684-00-0479-02	3000048614
3000044892	30000448920008	D684-00-0479-01	3000048711
3000044892	30000448920008	D684-00-0479-01U	3000048818
3000044892	30000448920008	D684-00-0479-01C	3000048963
3000044892	30000448920008	D684-00-0479-02	3000049263
3000044893	30000448930211	D684-00-0295-02U	3000047309
3000044893	30000448930211	D684-00-0295-05	3000047355
3000044893	30000448930211	D684-00-0295-01U	3000047432
3000044893	30000448930211	D684-00-0295-01	3000047433
3000044893	30000448930211	D684-00-0295-02	3000047434
3000044893	30000448930211	D684-00-0295-02	3000047434
3000044893	30000448930211	D684-00-0295-02	3000047523
3000044893	30000448930218	D684-00-0295-02U	3000047309
3000044893	30000448930218	D684-00-0295-05	3000047355
3000044893	30000448930218	D684-00-0295-01U	3000047432

[illegible]

[illegible]

3000045453	30000454530303	D684-00-0295-02	3000048025
3000045453	30000454530303	D684-00-0295-05	3000048228
3000045453	30000454530304	D684-00-0295-01U	3000047933
3000045453	30000454530304	D684-00-0295-02U	3000047934
3000045453	30000454530304	D684-00-0295-02	3000048024
3000045453	30000454530304	D684-00-0295-02	3000048024
3000045453	30000454530304	D684-00-0295-02	3000048025
3000045453	30000454530304	D684-00-0295-05	3000048228
3000045453	30000454530316	D684-00-0295-01U	3000047933
3000045453	30000454530316	D684-00-0295-02U	3000047934
3000045453	30000454530316	D684-00-0295-02	3000048024
3000045453	30000454530316	D684-00-0295-02	3000048024
3000045453	30000454530316	D684-00-0295-02	3000048025
3000045453	30000454530316	D684-00-0295-05	3000048228
3000045454	30000454540344	D684-00-0607	3000047804
3000045454	30000454540349	D684-00-0607	3000047804
3000045454	30000454540359	D684-00-0607	3000047804
3000045454	30000454540342	D684-00-0607	3000047804
3000045454	30000454540332	D684-00-0607	3000047804
3000045495	30000454950046	D684-00-0498-01	3000048230
3000045495	30000454950046	D684-00-0296-01	3000048243
3000045495	30000454950046	D684-00-0296-02	3000048244
3000045495	30000454950046	D684-00-0296-01U	3000048245
3000045495	30000454950046	D684-00-0498-01	3000049055
3000045495	30000454950046	D684-00-0296-01	3000049080
3000045495	30000454950046	D684-00-0296-01U	3000049085
3000045495	30000454950027	D684-00-0498-01	3000048230
3000045495	30000454950027	D684-00-0296-01	3000048243
3000045495	30000454950027	D684-00-0296-02	3000048244
3000045495	30000454950027	D684-00-0296-01U	3000048245
3000045495	30000454950027	D684-00-0498-01	3000049055
3000045495	30000454950027	D684-00-0296-01	3000049080
3000045495	30000454950027	D684-00-0296-01U	3000049085
3000045495	30000454950048	D684-00-0498-01	3000048230
3000045495	30000454950048	D684-00-0296-01	3000048243
3000045495	30000454950048	D684-00-0296-02	3000048244
3000045495	30000454950048	D684-00-0296-01U	3000048245
3000045495	30000454950048	D684-00-0498-01	3000049055
3000045495	30000454950048	D684-00-0296-01	3000049080
3000045495	30000454950048	D684-00-0296-01U	3000049085
3000045495	30000454950029	D684-00-0498-01	3000048230
3000045495	30000454950029	D684-00-0296-01	3000048243
3000045495	30000454950029	D684-00-0296-02	3000048244
3000045495	30000454950029	D684-00-0296-01U	3000048245
3000045495	30000454950029	D684-00-0498-01	3000049055
3000045495	30000454950029	D684-00-0296-01	3000049080
3000045495	30000454950029	D684-00-0296-01U	3000049085
3000045495	30000454950050	D684-00-0498-01	3000048230
3000045495	30000454950050	D684-00-0296-01	3000048243
3000045495	30000454950050	D684-00-0296-02	3000048244
3000045495	30000454950050	D684-00-0296-01U	3000048245
3000045495	30000454950050	D684-00-0498-01	3000049055
3000045495	30000454950050	D684-00-0296-01	3000049080
3000045495	30000454950050	D684-00-0296-01U	3000049085
3000045497	30000454970205	D684-00-0608	3000047928
3000045497	30000454970205	D684-00-0608	3000047928
3000045497	30000454970205	D684-00-0608	3000047928

[illegible]

[illegible]

[illegible]

[illegible]

3000046806	30000468060329	D684-00-0295-01U	3000049794
3000046806	30000468060329	D684-00-0295-02	3000049795
3000046806	30000468060329	D684-00-0295-02U	3000049796
3000046806	30000468060329	D684-00-0295-05	3000050257
3000046807	30000468070105	D684-00-0607	3000053322
3000046807	30000468070105	D684-00-0607	3000053323
3000046807	30000468070105	D684-00-0607	3000060107
3000046807	30000468070107	D684-00-0607	3000053322
3000046807	30000468070107	D684-00-0607	3000053323
3000046807	30000468070107	D684-00-0607	3000060107
3000046807	30000468070108	D684-00-0607	3000053322
3000046807	30000468070108	D684-00-0607	3000053323
3000046807	30000468070108	D684-00-0607	3000060107
3000046807	30000468070112	D684-00-0607	3000053322
3000046807	30000468070112	D684-00-0607	3000053323
3000046807	30000468070112	D684-00-0607	3000060107
3000046807	30000468070125	D684-00-0607	3000053322
3000046807	30000468070125	D684-00-0607	3000053323
3000046807	30000468070125	D684-00-0607	3000060107
3000046808	30000468080089	D684-00-0479-01C	3000050259
3000046808	30000468080089	D684-00-0479-01U	3000051295
3000046808	30000468080089	D684-00-0479-02	3000051298
3000046808	30000468080089	D684-00-0479-02	3000051298
3000046808	30000468080094	D684-00-0479-01C	3000050259
3000046808	30000468080094	D684-00-0479-01U	3000051295
3000046808	30000468080094	D684-00-0479-02	3000051298
3000046808	30000468080094	D684-00-0479-02	3000051298
3000046808	30000468080097	D684-00-0479-01C	3000050259
3000046808	30000468080097	D684-00-0479-01U	3000051295
3000046808	30000468080097	D684-00-0479-02	3000051298
3000046808	30000468080097	D684-00-0479-02	3000051298
3000046808	30000468080098	D684-00-0479-01C	3000050259
3000046808	30000468080098	D684-00-0479-01U	3000051295
3000046808	30000468080098	D684-00-0479-02	3000051298
3000046808	30000468080098	D684-00-0479-02	3000051298
3000046808	30000468080102	D684-00-0479-01C	3000050259
3000046808	30000468080102	D684-00-0479-01U	3000051295
3000046808	30000468080102	D684-00-0479-02	3000051298
3000046808	30000468080102	D684-00-0479-02	3000051298
3000047038	30000470380512	D684-00-0568-01	3000048703
3000047038	30000470380512	D684-00-0568-01	3000049591
3000047038	30000470380512	D684-00-0568-01U	3000049592
3000047038	30000470380520	D684-00-0568-01	3000048703
3000047038	30000470380520	D684-00-0568-01	3000049591
3000047038	30000470380520	D684-00-0568-01U	3000049592
3000047038	30000470380521	D684-00-0568-01	3000048703
3000047038	30000470380521	D684-00-0568-01	3000049591
3000047038	30000470380521	D684-00-0568-01U	3000049592
3000047038	30000470380529	D684-00-0568-01	3000048703
3000047038	30000470380529	D684-00-0568-01	3000049591
3000047038	30000470380529	D684-00-0568-01U	3000049592
3000047038	30000470380533	D684-00-0568-01	3000048703
3000047038	30000470380533	D684-00-0568-01	3000049591
3000047038	30000470380533	D684-00-0568-01U	3000049592
3000047040	30000470400240	D684-00-0608	3000053321
3000047040	30000470400240	D684-00-0608	3000060109
3000047040	30000470400245	D684-00-0608	3000053321
3000047040	30000470400245	D684-00-0608	3000060109

[illegible]

3000047057	30000470570232	D684-00-0296-01	3000051154
3000047057	30000470570232	D684-00-0296-01	3000051154
3000047057	30000470570232	D684-00-0296-02	3000051155
3000047421	30000474210289	D684-00-0480-01C	3000050258
3000047421	30000474210289	D684-00-0480-02	3000050669
3000047421	30000474210289	D684-00-0480-02	3000051151
3000047421	30000474210282	D684-00-0480-01C	3000050258
3000047421	30000474210282	D684-00-0480-02	3000050669
3000047421	30000474210282	D684-00-0480-02	3000051151
3000047421	30000474210278	D684-00-0480-01C	3000050258
3000047421	30000474210278	D684-00-0480-02	3000050669
3000047421	30000474210278	D684-00-0480-02	3000051151
3000047421	30000474210273	D684-00-0480-01C	3000050258
3000047421	30000474210273	D684-00-0480-02	3000050669
3000047421	30000474210273	D684-00-0480-02	3000051151
3000047421	30000474210275	D684-00-0480-01C	3000050258
3000047421	30000474210275	D684-00-0480-02	3000050669
3000047421	30000474210275	D684-00-0480-02	3000051151
3000047555	30000475550418	D684-00-0470-01	3000049914
3000047555	30000475550418	D684-00-0470-01U	3000050215
3000047555	30000475550418	D684-00-0470-01	3000050929
3000047555	30000475550418	D684-00-0470-01U	3000050930
3000047555	30000475550419	D684-00-0470-01	3000049914
3000047555	30000475550419	D684-00-0470-01U	3000050215
3000047555	30000475550419	D684-00-0470-01	3000050929
3000047555	30000475550419	D684-00-0470-01U	3000050930
3000047555	30000475550425	D684-00-0470-01	3000049914
3000047555	30000475550425	D684-00-0470-01U	3000050215
3000047555	30000475550425	D684-00-0470-01	3000050929
3000047555	30000475550425	D684-00-0470-01U	3000050930
3000047555	30000475550410	D684-00-0470-01	3000049914
3000047555	30000475550410	D684-00-0470-01U	3000050215
3000047555	30000475550410	D684-00-0470-01	3000050929
3000047555	30000475550410	D684-00-0470-01U	3000050930
3000047555	30000475550408	D684-00-0470-01	3000049914
3000047555	30000475550408	D684-00-0470-01U	3000050215
3000047555	30000475550408	D684-00-0470-01	3000050929
3000047555	30000475550408	D684-00-0470-01U	3000050930
3000047640	30000476400023	D684-00-0545-01	3000050506
3000047640	30000476400023	D684-00-0545-01	3000050507
3000047640	30000476400023	D684-00-0545-01	3000051382
3000047640	30000476400011	D684-00-0545-01	3000050506
3000047640	30000476400011	D684-00-0545-01	3000050507
3000047640	30000476400011	D684-00-0545-01	3000051382
3000047640	30000476400004	D684-00-0545-01	3000050506
3000047640	30000476400004	D684-00-0545-01	3000050507
3000047640	30000476400004	D684-00-0545-01	3000051382
3000047640	30000476400024	D684-00-0545-01	3000050506
3000047640	30000476400024	D684-00-0545-01	3000050507
3000047640	30000476400024	D684-00-0545-01	3000051382
3000047640	30000476400020	D684-00-0545-01	3000050506
3000047640	30000476400020	D684-00-0545-01	3000050507
3000047640	30000476400020	D684-00-0545-01	3000051382
3000047641	30000476410445	D684-00-0469-01U	3000050926
3000047641	30000476410445	D684-00-0469-01	3000050927
3000047641	30000476410453	D684-00-0469-01U	3000050926
3000047641	30000476410453	D684-00-0469-01	3000050927
3000047641	30000476410462	D684-00-0469-01U	3000050926

3000047641	30000476410462	D684-00-0469-01	3000050927
3000047641	30000476410467	D684-00-0469-01U	3000050926
3000047641	30000476410467	D684-00-0469-01	3000050927
3000047641	30000476410454	D684-00-0469-01U	3000050926
3000047641	30000476410454	D684-00-0469-01	3000050927
3000047646	30000476460133	D684-00-0480-01	3000051066
3000047646	30000476460133	D684-00-0480-01U	3000051122
3000047646	30000476460133	D684-00-0480-02	3000051152
3000047646	30000476460139	D684-00-0480-01	3000051066
3000047646	30000476460139	D684-00-0480-01U	3000051122
3000047646	30000476460139	D684-00-0480-02	3000051152
3000047646	30000476460131	D684-00-0480-01	3000051066
3000047646	30000476460131	D684-00-0480-01U	3000051122
3000047646	30000476460131	D684-00-0480-02	3000051152
3000047646	30000476460143	D684-00-0480-01	3000051066
3000047646	30000476460143	D684-00-0480-01U	3000051122
3000047646	30000476460143	D684-00-0480-02	3000051152
3000047646	30000476460135	D684-00-0480-01	3000051066
3000047646	30000476460135	D684-00-0480-01U	3000051122
3000047646	30000476460135	D684-00-0480-02	3000051152
3000047651	30000476510147	D684-00-0479-01	3000051296
3000047651	30000476510147	D684-00-0479-02	3000051299
3000047651	30000476510147	D684-00-0479-01C	3000051645
3000047651	30000476510147	D684-00-0479-01C	3000051646
3000047651	30000476510128	D684-00-0479-01	3000051296
3000047651	30000476510128	D684-00-0479-02	3000051299
3000047651	30000476510128	D684-00-0479-01C	3000051645
3000047651	30000476510128	D684-00-0479-01C	3000051646
3000047651	30000476510146	D684-00-0479-01	3000051296
3000047651	30000476510146	D684-00-0479-02	3000051299
3000047651	30000476510146	D684-00-0479-01C	3000051645
3000047651	30000476510146	D684-00-0479-01C	3000051646
3000047651	30000476510136	D684-00-0479-01	3000051296
3000047651	30000476510136	D684-00-0479-02	3000051299
3000047651	30000476510136	D684-00-0479-01C	3000051645
3000047651	30000476510136	D684-00-0479-01C	3000051646
3000047651	30000476510148	D684-00-0479-01	3000051296
3000047651	30000476510148	D684-00-0479-02	3000051299
3000047651	30000476510148	D684-00-0479-01C	3000051645
3000047651	30000476510148	D684-00-0479-01C	3000051646
3000047920	30000479200481	D684-00-0576-01	3000050497
3000047920	30000479200481	D684-00-0576-01	3000050612
3000047920	30000479200481	D684-00-0576-01U	3000050613
3000047920	30000479200483	D684-00-0576-01	3000050497
3000047920	30000479200483	D684-00-0576-01	3000050612
3000047920	30000479200483	D684-00-0576-01U	3000050613
3000047920	30000479200473	D684-00-0576-01	3000050497
3000047920	30000479200473	D684-00-0576-01	3000050612
3000047920	30000479200473	D684-00-0576-01U	3000050613
3000047920	30000479200474	D684-00-0576-01	3000050497
3000047920	30000479200474	D684-00-0576-01	3000050612
3000047920	30000479200474	D684-00-0576-01U	3000050613
3000047920	30000479200469	D684-00-0576-01	3000050497
3000047920	30000479200469	D684-00-0576-01	3000050612
3000047920	30000479200469	D684-00-0576-01U	3000050613
3000048319	30000483190037	D684-00-0568-01	3000050881
3000048319	30000483190037	D684-00-0568-01U	3000050882
3000048319	30000483190049	D684-00-0568-01	3000050881

3000048319	30000483190049	D684-00-0568-01U	3000050882
3000048319	30000483190039	D684-00-0568-01	3000050881
3000048319	30000483190039	D684-00-0568-01U	3000050882
3000048319	30000483190040	D684-00-0568-01	3000050881
3000048319	30000483190040	D684-00-0568-01U	3000050882
3000048319	30000483190027	D684-00-0568-01	3000050881
3000048319	30000483190027	D684-00-0568-01U	3000050882
3000048327	30000483270053	D684-00-0480-01C	3000051647
3000048327	30000483270053	D684-00-0480-01	3000051995
3000048327	30000483270070	D684-00-0480-01C	3000051647
3000048327	30000483270070	D684-00-0480-01	3000051995
3000048327	30000483270074	D684-00-0480-01C	3000051647
3000048327	30000483270074	D684-00-0480-01	3000051995
3000048327	30000483270066	D684-00-0480-01C	3000051647
3000048327	30000483270066	D684-00-0480-01	3000051995
3000048327	30000483270068	D684-00-0480-01C	3000051647
3000048327	30000483270068	D684-00-0480-01	3000051995
3000048410	30000484100152	D684-00-0470-01U	3000051379
3000048410	30000484100152	D684-00-0470-01	3000051380
3000048410	30000484100154	D684-00-0470-01U	3000051379
3000048410	30000484100154	D684-00-0470-01	3000051380
3000048410	30000484100155	D684-00-0470-01U	3000051379
3000048410	30000484100155	D684-00-0470-01	3000051380
3000048410	30000484100167	D684-00-0470-01U	3000051379
3000048410	30000484100167	D684-00-0470-01	3000051380
3000048410	30000484100173	D684-00-0470-01U	3000051379
3000048410	30000484100173	D684-00-0470-01	3000051380
3000048411	30000484110228	D684-00-0498-01	3000051378
3000048411	30000484110228	D684-00-0296-01	3000051751
3000048411	30000484110228	D684-00-0296-01U	3000051881
3000048411	30000484110234	D684-00-0498-01	3000051378
3000048411	30000484110234	D684-00-0296-01	3000051751
3000048411	30000484110234	D684-00-0296-01U	3000051881
3000048411	30000484110243	D684-00-0498-01	3000051378
3000048411	30000484110243	D684-00-0296-01	3000051751
3000048411	30000484110243	D684-00-0296-01U	3000051881
3000048411	30000484110242	D684-00-0498-01	3000051378
3000048411	30000484110242	D684-00-0296-01	3000051751
3000048411	30000484110242	D684-00-0296-01U	3000051881
3000048411	30000484110247	D684-00-0498-01	3000051378
3000048411	30000484110247	D684-00-0296-01	3000051751
3000048411	30000484110247	D684-00-0296-01U	3000051881
3000049084	30000490840495	D684-00-0469-01U	3000051752
3000049084	30000490840495	D684-00-0469-01	3000051753
3000049084	30000490840495	D684-00-0469-01	3000051754
3000049084	30000490840495	D684-00-0469-01	3000052790
3000049084	30000490840495	D684-00-0469-01U	3000052791
3000049084	30000490840495	D684-00-0469-01	3000053815
3000049084	30000490840500	D684-00-0469-01U	3000051752
3000049084	30000490840500	D684-00-0469-01	3000051753
3000049084	30000490840500	D684-00-0469-01	3000051754
3000049084	30000490840500	D684-00-0469-01	3000052790
3000049084	30000490840500	D684-00-0469-01U	3000052791
3000049084	30000490840500	D684-00-0469-01	3000053815
3000049084	30000490840501	D684-00-0469-01U	3000051752
3000049084	30000490840501	D684-00-0469-01	3000051753
3000049084	30000490840501	D684-00-0469-01	3000051754
3000049084	30000490840501	D684-00-0469-01	3000052790

[illegible]

3000049485	30000494850159	D684-00-0480-02	3000053049
3000049486	30000494860337	D684-00-0479-01	3000052326
3000049486	30000494860337	D684-00-0479-01C	3000052792
3000049486	30000494860336	D684-00-0479-01	3000052326
3000049486	30000494860336	D684-00-0479-01C	3000052792
3000049486	30000494860317	D684-00-0479-01	3000052326
3000049486	30000494860317	D684-00-0479-01C	3000052792
3000049486	30000494860324	D684-00-0479-01	3000052326
3000049486	30000494860324	D684-00-0479-01C	3000052792
3000049486	30000494860327	D684-00-0479-01	3000052326
3000049486	30000494860327	D684-00-0479-01C	3000052792
3000049920	30000499200064	D684-00-0568-01U	3000051544
3000049920	30000499200064	D684-00-0568-01U	3000051545
3000049920	30000499200064	D684-00-0568-01	3000051546
3000049920	30000499200064	D684-00-0568-01	3000051547
3000049920	30000499200051	D684-00-0568-01U	3000051544
3000049920	30000499200051	D684-00-0568-01U	3000051545
3000049920	30000499200051	D684-00-0568-01	3000051546
3000049920	30000499200051	D684-00-0568-01	3000051547
3000049920	30000499200072	D684-00-0568-01U	3000051544
3000049920	30000499200072	D684-00-0568-01U	3000051545
3000049920	30000499200072	D684-00-0568-01	3000051546
3000049920	30000499200072	D684-00-0568-01	3000051547
3000049920	30000499200056	D684-00-0568-01U	3000051544
3000049920	30000499200056	D684-00-0568-01U	3000051545
3000049920	30000499200056	D684-00-0568-01	3000051546
3000049920	30000499200056	D684-00-0568-01	3000051547
3000049920	30000499200059	D684-00-0568-01U	3000051544
3000049920	30000499200059	D684-00-0568-01U	3000051545
3000049920	30000499200059	D684-00-0568-01	3000051546
3000049920	30000499200059	D684-00-0568-01	3000051547
3000049924	30000499240495	D684-00-0576-01U	3000052323
3000049924	30000499240495	D684-00-0576-01	3000052324
3000049924	30000499240495	D684-00-0576-01	3000052325
3000049924	30000499240501	D684-00-0576-01U	3000052323
3000049924	30000499240501	D684-00-0576-01	3000052324
3000049924	30000499240501	D684-00-0576-01	3000052325
3000049924	30000499240507	D684-00-0576-01U	3000052323
3000049924	30000499240507	D684-00-0576-01	3000052324
3000049924	30000499240507	D684-00-0576-01	3000052325
3000049924	30000499240504	D684-00-0576-01U	3000052323
3000049924	30000499240504	D684-00-0576-01	3000052324
3000049924	30000499240504	D684-00-0576-01	3000052325
3000049924	30000499240510	D684-00-0576-01U	3000052323
3000049924	30000499240510	D684-00-0576-01	3000052324
3000049924	30000499240510	D684-00-0576-01	3000052325
3000049925	30000499250121	D684-00-0295-05	3000052253
3000049925	30000499250121	D684-00-0295-02	3000052254
3000049925	30000499250121	D684-00-0295-01U	3000052255
3000049925	30000499250129	D684-00-0295-05	3000052253
3000049925	30000499250129	D684-00-0295-02	3000052254
3000049925	30000499250129	D684-00-0295-01U	3000052255
3000049925	30000499250117	D684-00-0295-05	3000052253
3000049925	30000499250117	D684-00-0295-02	3000052254
3000049925	30000499250117	D684-00-0295-01U	3000052255
3000049925	30000499250127	D684-00-0295-05	3000052253
3000049925	30000499250127	D6	

3000049925	30000499250136	D684-00-0295-05	3000052253
3000049925	30000499250136	D684-00-0295-02	3000052254
3000049925	30000499250136	D684-00-0295-01U	3000052255
3000050491	30000504910788	D684-00-0480-01U	3000053152
3000050491	30000504910788	D684-00-0480-01U	3000053792
3000050491	30000504910788	D684-00-0480-01	3000054067
3000050491	30000504910788	D684-00-0480-01C	3000054290
3000050491	30000504910787	D684-00-0480-01U	3000053152
3000050491	30000504910787	D684-00-0480-01U	3000053792
3000050491	30000504910787	D684-00-0480-01	3000054067
3000050491	30000504910787	D684-00-0480-01C	3000054290
3000050491	30000504910782	D684-00-0480-01U	3000053152
3000050491	30000504910782	D684-00-0480-01U	3000053792
3000050491	30000504910782	D684-00-0480-01	3000054067
3000050491	30000504910782	D684-00-0480-01C	3000054290
3000050491	30000504910785	D684-00-0480-01U	3000053152
3000050491	30000504910785	D684-00-0480-01U	3000053792
3000050491	30000504910785	D684-00-0480-01	3000054067
3000050491	30000504910785	D684-00-0480-01C	3000054290
3000050491	30000504910779	D684-00-0480-01U	3000053152
3000050491	30000504910779	D684-00-0480-01U	3000053792
3000050491	30000504910779	D684-00-0480-01	3000054067
3000050491	30000504910779	D684-00-0480-01C	3000054290
3000050492	30000504920003	D684-00-0294-01U	3000055307
3000050492	30000504920003	D684-00-0294-01	3000055308
3000050492	30000504920003	D684-00-0294-01U	3000056433
3000050492	30000504920004	D684-00-0294-01U	3000055307
3000050492	30000504920004	D684-00-0294-01	3000055308
3000050492	30000504920004	D684-00-0294-01U	3000056433
3000050492	30000504920005	D684-00-0294-01U	3000055307
3000050492	30000504920005	D684-00-0294-01	3000055308
3000050492	30000504920005	D684-00-0294-01U	3000056433
3000050492	30000504920001	D684-00-0294-01U	3000055307
3000050492	30000504920001	D684-00-0294-01	3000055308
3000050492	30000504920001	D684-00-0294-01U	3000056433
3000050492	30000504920002	D684-00-0294-01U	3000055307
3000050492	30000504920002	D684-00-0294-01	3000055308
3000050492	30000504920002	D684-00-0294-01U	3000056433
3000050493	30000504930121	D684-00-0545-01	3000053314
3000050493	30000504930108	D684-00-0545-01	3000053314
3000050493	30000504930103	D684-00-0545-01	3000053314
3000050493	30000504930125	D684-00-0545-01	3000053314
3000050493	30000504930123	D684-00-0545-01	3000053314
3000050494	30000504940031	D684-00-0470-01	3000052645
3000050494	30000504940031	D684-00-0470-01U	3000052646
3000050494	30000504940040	D684-00-0470-01	3000052645
3000050494	30000504940040	D684-00-0470-01U	3000052646
3000050494	30000504940042	D684-00-0470-01	3000052645
3000050494	30000504940042	D684-00-0470-01U	3000052646
3000050494	30000504940032	D684-00-0470-01	3000052645
3000050494	30000504940032	D684-00-0470-01U	3000052646
3000050494	30000504940035	D684-00-0470-01	3000052645
3000050494	30000504940035	D684-00-0470-01U	3000052646
3000050639	30000506390079	D684-00-0296-01U	3000053310
3000050639	30000506390079	D684-00-0296-01	3000053311
3000050639	30000506390079	D684-00-0296-02	3000053312
3000050639	30000506390079	D684-00-0498-01	3000054590
3000050639	30000506390079	D684-00-0296-01	3000054591

[illegible]

3000051061	30000510610255	D684-00-0295-02	3000053295
3000051061	30000510610255	D684-00-0295-01	3000053296
3000051061	30000510610255	D684-00-0295-02	3000054426
3000051061	30000510610256	D684-00-0295-02	3000053295
3000051061	30000510610256	D684-00-0295-01	3000053296
3000051061	30000510610256	D684-00-0295-02	3000054426
3000051061	30000510610259	D684-00-0295-02	3000053295
3000051061	30000510610259	D684-00-0295-01	3000053296
3000051061	30000510610259	D684-00-0295-02	3000054426
3000051061	30000510610261	D684-00-0295-02	3000053295
3000051061	30000510610261	D684-00-0295-01	3000053296
3000051061	30000510610261	D684-00-0295-02	3000054426
3000051095	30000510950744	D684-00-0576-01U	3000053149
3000051095	30000510950744	D684-00-0576-01U	3000053150
3000051095	30000510950744	D684-00-0576-01	3000053151
3000051095	30000510950746	D684-00-0576-01U	3000053149
3000051095	30000510950746	D684-00-0576-01U	3000053150
3000051095	30000510950746	D684-00-0576-01	3000053151
3000051095	30000510950749	D684-00-0576-01U	3000053149
3000051095	30000510950749	D684-00-0576-01U	3000053150
3000051095	30000510950749	D684-00-0576-01	3000053151
3000051095	30000510950755	D684-00-0576-01U	3000053149
3000051095	30000510950755	D684-00-0576-01U	3000053150
3000051095	30000510950755	D684-00-0576-01	3000053151
3000051095	30000510950761	D684-00-0576-01U	3000053149
3000051095	30000510950761	D684-00-0576-01U	3000053150
3000051095	30000510950761	D684-00-0576-01	3000053151
3000051210	30000512100078	D684-00-0568-01	3000053550
3000051210	30000512100078	D684-00-0568-01U	3000053551
3000051210	30000512100078	D684-00-0568-01	3000054595
3000051210	30000512100083	D684-00-0568-01	3000053550
3000051210	30000512100083	D684-00-0568-01U	3000053551
3000051210	30000512100083	D684-00-0568-01	3000054595
3000051210	30000512100085	D684-00-0568-01	3000053550
3000051210	30000512100085	D684-00-0568-01U	3000053551
3000051210	30000512100085	D684-00-0568-01	3000054595
3000051210	30000512100090	D684-00-0568-01	3000053550
3000051210	30000512100090	D684-00-0568-01U	3000053551
3000051210	30000512100090	D684-00-0568-01	3000054595
3000051210	30000512100092	D684-00-0568-01	3000053550
3000051210	30000512100092	D684-00-0568-01U	3000053551
3000051210	30000512100092	D684-00-0568-01	3000054595
3000051212	30000512120152	D684-00-0469-01	3000053816
3000051212	30000512120152	D684-00-0469-01U	3000053817
3000051212	30000512120157	D684-00-0469-01	3000053816
3000051212	30000512120157	D684-00-0469-01U	3000053817
3000051212	30000512120158	D684-00-0469-01	3000053816
3000051212	30000512120158	D684-00-0469-01U	3000053817
3000051212	30000512120164	D684-00-0469-01	3000053816
3000051212	30000512120164	D684-00-0469-01U	3000053817
3000051212	30000512120167	D684-00-0469-01	3000053816
3000051212	30000512120167	D684-00-0469-01U	3000053817
3000051288	30000512880160	D684-00-0480-01C	3000054291
3000051288	30000512880160	D684-00-0480-02	3000054337
3000051288	30000512880160	D684-00-0480-02U	3000054879
3000051288	30000512880160	D684-00-0480-01U	3000054880
3000051288	30000512880160	D684-00-0480-01	3000054914
3000051288	30000512880160	D684-00-0480-02	3000054915

[illegible]

3000051289	30000512890307	D684-00-0479-01C	3000054317
3000051289	30000512890307	D684-00-0479-02	3000054338
3000051289	30000512890307	D684-00-0479-01	3000054821
3000051289	30000512890307	D684-00-0479-02	3000054822
3000051289	30000512890307	D684-00-0479-01C	3000054873
3000051289	30000512890307	D684-00-0479-01U	3000054874
3000051744	30000517440437	D684-00-0470-01U	3000054089
3000051744	30000517440437	D684-00-0470-01	3000054090
3000051744	30000517440437	D684-00-0470-01	3000054090
3000051744	30000517440437	D684-00-0470-01U	3000055098
3000051744	30000517440438	D684-00-0470-01U	3000054089
3000051744	30000517440438	D684-00-0470-01	3000054090
3000051744	30000517440438	D684-00-0470-01	3000054090
3000051744	30000517440438	D684-00-0470-01U	3000055098
3000051744	30000517440446	D684-00-0470-01U	3000054089
3000051744	30000517440446	D684-00-0470-01	3000054090
3000051744	30000517440446	D684-00-0470-01	3000054090
3000051744	30000517440446	D684-00-0470-01U	3000055098
3000051744	30000517440448	D684-00-0470-01U	3000054089
3000051744	30000517440448	D684-00-0470-01	3000054090
3000051744	30000517440448	D684-00-0470-01	3000054090
3000051744	30000517440448	D684-00-0470-01U	3000055098
3000051744	30000517440455	D684-00-0470-01U	3000054089
3000051744	30000517440455	D684-00-0470-01	3000054090
3000051744	30000517440455	D684-00-0470-01	3000054090
3000051744	30000517440455	D684-00-0470-01U	3000055098
3000051745	30000517450419	D684-00-0576-01	3000054475
3000051745	30000517450419	D684-00-0576-01U	3000054476
3000051745	30000517450419	D684-00-0576-01U	3000054477
3000051745	30000517450411	D684-00-0576-01	3000054475
3000051745	30000517450411	D684-00-0576-01U	3000054476
3000051745	30000517450411	D684-00-0576-01U	3000054477
3000051745	30000517450427	D684-00-0576-01	3000054475
3000051745	30000517450427	D684-00-0576-01U	3000054476
3000051745	30000517450427	D684-00-0576-01U	3000054477
3000051745	30000517450430	D684-00-0576-01	3000054475
3000051745	30000517450430	D684-00-0576-01U	3000054476
3000051745	30000517450430	D684-00-0576-01U	3000054477
3000051745	30000517450416	D684-00-0576-01	3000054475
3000051745	30000517450416	D684-00-0576-01U	3000054476
3000051745	30000517450416	D684-00-0576-01U	3000054477
3000051877	30000518770111	D684-00-0480-01C	3000054967
3000051877	30000518770111	D684-00-0480-01C	3000056074
3000051877	30000518770111	D684-00-0480-02U	3000056192
3000051877	30000518770111	D684-00-0480-08	3000056260
3000051877	30000518770111	D684-00-0480-01U	3000056936
3000051877	30000518770113	D684-00-0480-01C	3000054967
3000051877	30000518770113	D684-00-0480-01C	3000056074
3000051877	30000518770113	D684-00-0480-02U	3000056192
3000051877	30000518770113	D684-00-0480-08	3000056260
3000051877	30000518770113	D684-00-0480-01U	3000056936
3000051877	30000518770123	D684-00-0480-01C	3000054967
3000051877	30000518770123	D684-00-0480-01C	3000056074
3000051877	30000518770123	D684-00-0480-02U	3000056192
3000051877	30000518770123	D684-00-0480-08	3000056260
3000051877	30000518770123	D684-00-0480-01U	3000056936
3000051877	30000518		

3000051877	30000518770125	D684-00-0480-02U	3000056192
3000051877	30000518770125	D684-00-0480-08	3000056260
3000051877	30000518770125	D684-00-0480-01U	3000056936
3000051877	30000518770129	D684-00-0480-01C	3000054967
3000051877	30000518770129	D684-00-0480-01C	3000056074
3000051877	30000518770129	D684-00-0480-02U	3000056192
3000051877	30000518770129	D684-00-0480-08	3000056260
3000051877	30000518770129	D684-00-0480-01U	3000056936
3000052256	30000522560382	D684-00-0568-01U	3000054593
3000052256	30000522560382	D684-00-0568-01	3000054594
3000052256	30000522560389	D684-00-0568-01U	3000054593
3000052256	30000522560389	D684-00-0568-01	3000054594
3000052256	30000522560390	D684-00-0568-01U	3000054593
3000052256	30000522560390	D684-00-0568-01	3000054594
3000052256	30000522560392	D684-00-0568-01U	3000054593
3000052256	30000522560392	D684-00-0568-01	3000054594
3000052256	30000522560406	D684-00-0568-01U	3000054593
3000052256	30000522560406	D684-00-0568-01	3000054594
3000052334	30000523340118	D684-00-0295-05	3000054424
3000052334	30000523340118	D684-00-0295-02	3000054425
3000052334	30000523340121	D684-00-0295-05	3000054424
3000052334	30000523340121	D684-00-0295-02	3000054425
3000052334	30000523340125	D684-00-0295-05	3000054424
3000052334	30000523340125	D684-00-0295-02	3000054425
3000052334	30000523340108	D684-00-0295-05	3000054424
3000052334	30000523340108	D684-00-0295-02	3000054425
3000052334	30000523340113	D684-00-0295-05	3000054424
3000052334	30000523340113	D684-00-0295-02	3000054425
3000052336	30000523360239	D684-00-0296-01	3000054592
3000052336	30000523360239	D684-00-0296-01U	3000054729
3000052336	30000523360239	D684-00-0296-02	3000054730
3000052336	30000523360239	D684-00-0498-01	3000056355
3000052336	30000523360239	D684-00-0296-01	3000056356
3000052336	30000523360234	D684-00-0296-01	3000054592
3000052336	30000523360234	D684-00-0296-01U	3000054729
3000052336	30000523360234	D684-00-0296-02	3000054730
3000052336	30000523360234	D684-00-0498-01	3000056355
3000052336	30000523360234	D684-00-0296-01	3000056356
3000052336	30000523360231	D684-00-0296-01	3000054592
3000052336	30000523360231	D684-00-0296-01U	3000054729
3000052336	30000523360231	D684-00-0296-02	3000054730
3000052336	30000523360231	D684-00-0498-01	3000056355
3000052336	30000523360231	D684-00-0296-01	3000056356
3000052336	30000523360233	D684-00-0296-01	3000054592
3000052336	30000523360233	D684-00-0296-01U	3000054729
3000052336	30000523360233	D684-00-0296-02	3000054730
3000052336	30000523360233	D684-00-0498-01	3000056355
3000052336	30000523360233	D684-00-0296-01	3000056356
3000052336	30000523360249	D684-00-0296-01	3000054592
3000052336	30000523360249	D684-00-0296-01U	3000054729
3000052336	30000523360249	D684-00-0296-02	3000054730
3000052336	30000523360249	D684-00-0498-01	3000056355
3000052336	30000523360249	D684-00-0296-01	3000056356
3000053154	30000531540078	D684-00-0545-01	3000055312
3000053154	30000531540079	D684-00-0545-01	3000055312
3000053154	30000531540082	D684-00-0545-01	3000055312
3000053154	30000531540087	D684-00-0545-01	3000055312
3000053154	30000531540089	D684-00-0545-01	3000055312

[illegible]

3000053160	30000531600005	D684-00-0479-01C	3000056073
3000053160	30000531600005	D684-00-0479-01	3000056194
3000053160	30000531600005	D684-00-0479-08	3000056261
3000053160	30000531600005	D684-00-0479-08	3000056262
3000053160	30000531600008	D684-00-0479-01C	3000056073
3000053160	30000531600008	D684-00-0479-01	3000056194
3000053160	30000531600008	D684-00-0479-08	3000056261
3000053160	30000531600008	D684-00-0479-08	3000056262
3000053160	30000531600009	D684-00-0479-01C	3000056073
3000053160	30000531600009	D684-00-0479-01	3000056194
3000053160	30000531600009	D684-00-0479-08	3000056261
3000053160	30000531600009	D684-00-0479-08	3000056262
3000053476	30000534760054	D684-00-0568-01U	3000055179
3000053476	30000534760054	D684-00-0568-01	3000055238
3000053476	30000534760054	D684-00-0568-01U	3000055239
3000053476	30000534760057	D684-00-0568-01U	3000055179
3000053476	30000534760057	D684-00-0568-01	3000055238
3000053476	30000534760057	D684-00-0568-01U	3000055239
3000053476	30000534760059	D684-00-0568-01U	3000055179
3000053476	30000534760059	D684-00-0568-01	3000055238
3000053476	30000534760059	D684-00-0568-01U	3000055239
3000053476	30000534760062	D684-00-0568-01U	3000055179
3000053476	30000534760062	D684-00-0568-01	3000055238
3000053476	30000534760062	D684-00-0568-01U	3000055239
3000053476	30000534760075	D684-00-0568-01U	3000055179
3000053476	30000534760075	D684-00-0568-01	3000055238
3000053476	30000534760075	D684-00-0568-01U	3000055239
3000053480	30000534800228	D684-00-0469-01	3000056436
3000053480	30000534800228	D684-00-0469-01U	3000056437
3000053480	30000534800229	D684-00-0469-01	3000056436
3000053480	30000534800229	D684-00-0469-01U	3000056437
3000053480	30000534800242	D684-00-0469-01	3000056436
3000053480	30000534800242	D684-00-0469-01U	3000056437
3000053480	30000534800243	D684-00-0469-01	3000056436
3000053480	30000534800243	D684-00-0469-01U	3000056437
3000053480	30000534800251	D684-00-0469-01	3000056436
3000053480	30000534800251	D684-00-0469-01U	3000056437
3000053631	30000536310010	D684-00-0295-01	3000055303
3000053631	30000536310010	D684-00-0295-05	3000055304
3000053631	30000536310010	D684-00-0295-02	3000055305
3000053631	30000536310013	D684-00-0295-01	3000055303
3000053631	30000536310013	D684-00-0295-05	3000055304
3000053631	30000536310013	D684-00-0295-02	3000055305
3000053631	30000536310015	D684-00-0295-01	3000055303
3000053631	30000536310015	D684-00-0295-05	3000055304
3000053631	30000536310015	D684-00-0295-02	3000055305
3000053631	30000536310017	D684-00-0295-01	3000055303
3000053631	30000536310017	D684-00-0295-05	3000055304
3000053631	30000536310017	D684-00-0295-02	3000055305
3000053631	30000536310019	D684-00-0295-01	3000055303
3000053631	30000536310019	D684-00-0295-05	3000055304
3000053631	30000536310019	D684-00-0295-02	3000055305
3000053636	30000536360064	D684-00-0296-01U	3000056538
3000053636	30000536360064	D684-00-0296-02U	3000056539
3000053636	30000536360064	D684-00-0498-01	3000056628
3000053636	30000536360064	D684-00-0296-01	3000056629
3000053636	30000536360064	D684-00-0296-01	3000058277
3000053636	30000536360068	D684-00-0296-01U	3000056538

[illegible]

[illegible]

[illegible]

[illegible]

3000054904	30000549040579	D684-00-0576-01U	3000057035
3000054904	30000549040579	D684-00-0576-01U	3000057381
3000054904	30000549040579	D684-00-0576-01U	3000057592
3000054904	30000549040581	D684-00-0576-01	3000057034
3000054904	30000549040581	D684-00-0576-01U	3000057035
3000054904	30000549040581	D684-00-0576-01U	3000057381
3000054904	30000549040581	D684-00-0576-01U	3000057592
3000055038	30000550380145	D684-00-0568-01U	3000057382
3000055038	30000550380145	D684-DO-0567	0
3000055038	30000550380145	D684-00-0568-01	3000057510
3000055038	30000550380145	D684-00-0568-01	3000057511
3000055038	30000550380145	D684-00-0568-01U	3000057546
3000055038	30000550380132	D684-00-0568-01U	3000057382
3000055038	30000550380132	D684-DO-0567	0
3000055038	30000550380132	D684-00-0568-01	3000057510
3000055038	30000550380132	D684-00-0568-01	3000057511
3000055038	30000550380132	D684-00-0568-01U	3000057546
3000055038	30000550380148	D684-00-0568-01U	3000057382
3000055038	30000550380148	D684-DO-0567	0
3000055038	30000550380148	D684-00-0568-01	3000057510
3000055038	30000550380148	D684-00-0568-01	3000057511
3000055038	30000550380148	D684-00-0568-01U	3000057546
3000055038	30000550380141	D684-00-0568-01U	3000057382
3000055038	30000550380141	D684-DO-0567	0
3000055038	30000550380141	D684-00-0568-01	3000057510
3000055038	30000550380141	D684-00-0568-01	3000057511
3000055038	30000550380141	D684-00-0568-01U	3000057546
3000055038	30000550380139	D684-00-0568-01U	3000057382
3000055038	30000550380139	D684-DO-0567	0
3000055038	30000550380139	D684-00-0568-01	3000057510
3000055038	30000550380139	D684-00-0568-01	3000057511
3000055038	30000550380139	D684-00-0568-01U	3000057546
3000055039	30000550390301	D684-00-0560-01	3000057319
3000055039	30000550390301	D684-00-0560-01	3000057968
3000055039	30000550390301	D684-00-0560-01	3000058129
3000055039	30000550390301	D684-00-0560-01	3000062638
3000055039	30000550390303	D684-00-0560-01	3000057319
3000055039	30000550390303	D684-00-0560-01	3000057968
3000055039	30000550390303	D684-00-0560-01	3000058129
3000055039	30000550390303	D684-00-0560-01	3000062638
3000055039	30000550390307	D684-00-0560-01	3000057319
3000055039	30000550390307	D684-00-0560-01	3000057968
3000055039	30000550390307	D684-00-0560-01	3000058129
3000055039	30000550390307	D684-00-0560-01	3000062638
3000055039	30000550390309	D684-00-0560-01	3000057319
3000055039	30000550390309	D684-00-0560-01	3000057968
3000055039	30000550390309	D684-00-0560-01	3000058129
3000055039	30000550390309	D684-00-0560-01	3000062638
3000055039	30000550390310	D684-00-0560-01	3000057319
3000055039	30000550390310	D684-00-0560-01	3000057968
3000055039	30000550390310	D684-00-0560-01	3000058129
3000055039	30000550390310	D684-00-0560-01	3000062638
3000055042	30000550420193	D684-00-0561-01	3000057318
3000055042	30000550420193	D684-00-0561-01	3000057970
3000055042	30000550420193	D684-00-0561-01	3000062073
3000055042	30000550420193	D684-00-0561-01	3000062639
3000055042	30000550420194	D684-00-0561-01	3000057318
3000055042	30000550420194	D684-00-0561-01	3000057970

[illegible]

[illegible]

[illegible]

[illegible]

[illegible]

[illegible]

[illegible]

[illegible]

3000057874	30000578740174	D684-00-0296-02	3000060543
3000057874	30000578740174	D684-00-0498-01	3000060544
3000057874	30000578740174	D684-00-0296-02	3000060545
3000057874	30000578740174	D684-00-0296-02	3000066063
3000057877	30000578770228	D684-00-0608	3000060110
3000057877	30000578770228	D684-00-0608	3000067699
3000057877	30000578770230	D684-00-0608	3000060110
3000057877	30000578770230	D684-00-0608	3000067699
3000057877	30000578770237	D684-00-0608	3000060110
3000057877	30000578770237	D684-00-0608	3000067699
3000057877	30000578770243	D684-00-0608	3000060110
3000057877	30000578770243	D684-00-0608	3000067699
3000057877	30000578770246	D684-00-0608	3000060110
3000057877	30000578770246	D684-00-0608	3000067699
3000058108	30000581080177	D684-00-0480-01C	3000060733
3000058108	30000581080177	D684-00-0480-01C	3000060733
3000058108	30000581080177	D684-00-0480-02	3000060734
3000058108	30000581080177	D684-00-0480-01	3000060848
3000058108	30000581080177	D684-00-0480-08	3000061584
3000058108	30000581080177	D684-00-0480-02	3000062028
3000058108	30000581080177	D684-00-0480-01	3000062743
3000058108	30000581080178	D684-00-0480-01C	3000060733
3000058108	30000581080178	D684-00-0480-01C	3000060733
3000058108	30000581080178	D684-00-0480-02	3000060734
3000058108	30000581080178	D684-00-0480-01	3000060848
3000058108	30000581080178	D684-00-0480-08	3000061584
3000058108	30000581080178	D684-00-0480-02	3000062028
3000058108	30000581080178	D684-00-0480-01	3000062743
3000058108	30000581080180	D684-00-0480-01C	3000060733
3000058108	30000581080180	D684-00-0480-01C	3000060733
3000058108	30000581080180	D684-00-0480-02	3000060734
3000058108	30000581080180	D684-00-0480-01	3000060848
3000058108	30000581080180	D684-00-0480-08	3000061584
3000058108	30000581080180	D684-00-0480-02	3000062028
3000058108	30000581080180	D684-00-0480-01	3000062743
3000058108	30000581080181	D684-00-0480-01C	3000060733
3000058108	30000581080181	D684-00-0480-01C	3000060733
3000058108	30000581080181	D684-00-0480-02	3000060734
3000058108	30000581080181	D684-00-0480-01	3000060848
3000058108	30000581080181	D684-00-0480-08	3000061584
3000058108	30000581080181	D684-00-0480-02	3000062028
3000058108	30000581080181	D684-00-0480-01	3000062743
3000058108	30000581080185	D684-00-0480-01C	3000060733
3000058108	30000581080185	D684-00-0480-01C	3000060733
3000058108	30000581080185	D684-00-0480-02	3000060734
3000058108	30000581080185	D684-00-0480-01	3000060848
3000058108	30000581080185	D684-00-0480-08	3000061584
3000058108	30000581080185	D684-00-0480-02	3000062028
3000058108	30000581080185	D684-00-0480-01	3000062743
3000058109	30000581090082	D684-00-0576-01U	3000061115
3000058109	30000581090082	D684-00-0576-01	3000061119
3000058109	30000581090082	D684-00-0576-01U	3000062403
3000058109	30000581090084	D684-00-0576-01U	3000061115
3000058109	30000581090084	D684-00-0576-01	3000061119
3000058109	30000581090084	D684-00-0576-01U	3000062403
3000058109	30000581090088	D684-00-0576-01U	3000061115
3000058109	30000581090088	D684-00-0576-01	3000061119
3000058109	30000581090088	D684-00-0576-01U	3000062403

[illegible]

3000058291	30000582910279	D684-00-0470-01U	3000065985
3000058291	30000582910279	D684-00-0470-01	3000065986
3000058291	30000582910280	D684-00-0470-01U	3000065125
3000058291	30000582910280	D684-00-0470-01U	3000065984
3000058291	30000582910280	D684-00-0470-01U	3000065985
3000058291	30000582910280	D684-00-0470-01	3000065986
3000058291	30000582910285	D684-00-0470-01U	3000065125
3000058291	30000582910285	D684-00-0470-01U	3000065984
3000058291	30000582910285	D684-00-0470-01U	3000065985
3000058291	30000582910285	D684-00-0470-01	3000065986
3000058291	30000582910275	D684-00-0470-01U	3000065125
3000058291	30000582910275	D684-00-0470-01U	3000065984
3000058291	30000582910275	D684-00-0470-01U	3000065985
3000058291	30000582910275	D684-00-0470-01	3000065986
3000058291	30000582910286	D684-00-0470-01U	3000065125
3000058291	30000582910286	D684-00-0470-01U	3000065984
3000058291	30000582910286	D684-00-0470-01U	3000065985
3000058291	30000582910286	D684-00-0470-01	3000065986
3000058739	30000587390096	D684-00-0478-01U	3000061193
3000058739	30000587390096	D684-00-0478-02	3000061194
3000058739	30000587390097	D684-00-0478-01U	3000061193
3000058739	30000587390097	D684-00-0478-02	3000061194
3000058739	30000587390098	D684-00-0478-01U	3000061193
3000058739	30000587390098	D684-00-0478-02	3000061194
3000058739	30000587390099	D684-00-0478-01U	3000061193
3000058739	30000587390099	D684-00-0478-02	3000061194
3000058739	30000587390100	D684-00-0478-01U	3000061193
3000058739	30000587390100	D684-00-0478-02	3000061194
3000058740	30000587400208	D684-00-0295-05	3000061195
3000058740	30000587400208	D684-00-0295-05	3000061196
3000058740	30000587400208	D684-00-0295-01U	3000061197
3000058740	30000587400208	D684-00-0295-02	3000061341
3000058740	30000587400209	D684-00-0295-05	3000061195
3000058740	30000587400209	D684-00-0295-05	3000061196
3000058740	30000587400209	D684-00-0295-01U	3000061197
3000058740	30000587400209	D684-00-0295-02	3000061341
3000058740	30000587400210	D684-00-0295-05	3000061195
3000058740	30000587400210	D684-00-0295-05	3000061196
3000058740	30000587400210	D684-00-0295-01U	3000061197
3000058740	30000587400210	D684-00-0295-02	3000061341
3000058740	30000587400211	D684-00-0295-05	3000061195
3000058740	30000587400211	D684-00-0295-05	3000061196
3000058740	30000587400211	D684-00-0295-01U	3000061197
3000058740	30000587400211	D684-00-0295-02	3000061341
3000058740	30000587400212	D684-00-0295-05	3000061195
3000058740	30000587400212	D684-00-0295-05	3000061196
3000058740	30000587400212	D684-00-0295-01U	3000061197
3000058740	30000587400212	D684-00-0295-02	3000061341
3000058903	30000589030496	D684-00-0479-08	3000061585
3000058903	30000589030496	D684-00-0479-01C	3000062026
3000058903	30000589030496	D684-00-0479-01U	3000062027
3000058903	30000589030496	D684-00-0479-01	3000062934
3000058903	30000589030496	D684-00-0479-01C	3000063084
3000058903	30000589030515	D684-00-0479-08	3000061585
3000058903	30000589030515	D684-00-0479-01C	3000062026
3000058903	30000589030515	D684-00-0479-01U	3000062027
3000058903	30000589030515</		

3000058903	30000589030510	D684-00-0479-08	3000061585
3000058903	30000589030510	D684-00-0479-01C	3000062026
3000058903	30000589030510	D684-00-0479-01U	3000062027
3000058903	30000589030510	D684-00-0479-01	3000062934
3000058903	30000589030510	D684-00-0479-01C	3000063084
3000058903	30000589030500	D684-00-0479-08	3000061585
3000058903	30000589030500	D684-00-0479-01C	3000062026
3000058903	30000589030500	D684-00-0479-01U	3000062027
3000058903	30000589030500	D684-00-0479-01	3000062934
3000058903	30000589030500	D684-00-0479-01C	3000063084
3000058903	30000589030516	D684-00-0479-08	3000061585
3000058903	30000589030516	D684-00-0479-01C	3000062026
3000058903	30000589030516	D684-00-0479-01U	3000062027
3000058903	30000589030516	D684-00-0479-01	3000062934
3000058903	30000589030516	D684-00-0479-01C	3000063084
3000059179	30000591790309	D684-00-0608	3000062192
3000059179	30000591790309	D684-00-0608	3000067700
3000059179	30000591790312	D684-00-0608	3000062192
3000059179	30000591790312	D684-00-0608	3000067700
3000059179	30000591790317	D684-00-0608	3000062192
3000059179	30000591790317	D684-00-0608	3000067700
3000059179	30000591790320	D684-00-0608	3000062192
3000059179	30000591790320	D684-00-0608	3000067700
3000059179	30000591790324	D684-00-0608	3000062192
3000059179	30000591790324	D684-00-0608	3000067700
3000059185	30000591850050	D684-00-0561-01	3000062640
3000059185	30000591850050	D684-00-0561-01	3000068106
3000059185	30000591850050	D684-00-0561-01	3000071287
3000059185	30000591850050	D684-00-0561-01	3000072717
3000059185	30000591850051	D684-00-0561-01	3000062640
3000059185	30000591850051	D684-00-0561-01	3000068106
3000059185	30000591850051	D684-00-0561-01	3000071287
3000059185	30000591850051	D684-00-0561-01	3000072717
3000059185	30000591850052	D684-00-0561-01	3000062640
3000059185	30000591850052	D684-00-0561-01	3000068106
3000059185	30000591850052	D684-00-0561-01	3000071287
3000059185	30000591850052	D684-00-0561-01	3000072717
3000059185	30000591850053	D684-00-0561-01	3000062640
3000059185	30000591850053	D684-00-0561-01	3000068106
3000059185	30000591850053	D684-00-0561-01	3000071287
3000059185	30000591850053	D684-00-0561-01	3000072717
3000059185	30000591850054	D684-00-0561-01	3000062640
3000059185	30000591850054	D684-00-0561-01	3000068106
3000059185	30000591850054	D684-00-0561-01	3000071287
3000059185	30000591850054	D684-00-0561-01	3000072717
3000059283	30000592830520	D684-00-0480-02	3000062697
3000059283	30000592830520	D684-00-0480-02U	3000062803
3000059283	30000592830520	D684-00-0480-01C	3000063082
3000059283	30000592830520	D684-00-0480-01C	3000063083
3000059283	30000592830520	D684-00-0480-01U	3000063258
3000059283	30000592830520	D684-00-0480-01C	3000064046
3000059283	30000592830525	D684-00-0480-02	3000062697
3000059283	30000592830525	D684-00-0480-02U	3000062803
3000059283	30000592830525	D684-00-0480-01C	3000063082
3000059283	30000592830525	D684-00-0480-01C	3000063083
3000059283	30000592830525	D684-00-0480-01U	3000063258
3000059283	30000592830525	D684-00-0480-01C	3000064046
3000059283	30000592830530	D684-00-0480-02	3000062697

3000059283	30000592830530	D684-00-0480-02U	3000062803
3000059283	30000592830530	D684-00-0480-01C	3000063082
3000059283	30000592830530	D684-00-0480-01C	3000063083
3000059283	30000592830530	D684-00-0480-01U	3000063258
3000059283	30000592830530	D684-00-0480-01C	3000064046
3000059283	30000592830532	D684-00-0480-02	3000062697
3000059283	30000592830532	D684-00-0480-02U	3000062803
3000059283	30000592830532	D684-00-0480-01C	3000063082
3000059283	30000592830532	D684-00-0480-01C	3000063083
3000059283	30000592830532	D684-00-0480-01U	3000063258
3000059283	30000592830532	D684-00-0480-01C	3000064046
3000059283	30000592830535	D684-00-0480-02	3000062697
3000059283	30000592830535	D684-00-0480-02U	3000062803
3000059283	30000592830535	D684-00-0480-01C	3000063082
3000059283	30000592830535	D684-00-0480-01C	3000063083
3000059283	30000592830535	D684-00-0480-01U	3000063258
3000059283	30000592830535	D684-00-0480-01C	3000064046
3000059310	30000593100756	D684-00-0607	3000062255
3000059310	30000593100756	D684-00-0607	3000067698
3000059310	30000593100766	D684-00-0607	3000062255
3000059310	30000593100766	D684-00-0607	3000067698
3000059310	30000593100765	D684-00-0607	3000062255
3000059310	30000593100765	D684-00-0607	3000067698
3000059310	30000593100763	D684-00-0607	3000062255
3000059310	30000593100763	D684-00-0607	3000067698
3000059310	30000593100774	D684-00-0607	3000062255
3000059310	30000593100774	D684-00-0607	3000067698
3000059733	30000597330288	D684-00-0568-01U	3000063259
3000059733	30000597330290	D684-00-0568-01U	3000063259
3000059733	30000597330293	D684-00-0568-01U	3000063259
3000059733	30000597330297	D684-00-0568-01U	3000063259
3000059733	30000597330307	D684-00-0568-01U	3000063259
3000059918	30000599180515	D684-00-0576-01U	3000063260
3000059918	30000599180515	D684-00-0576-01	3000063261
3000059918	30000599180501	D684-00-0576-01U	3000063260
3000059918	30000599180501	D684-00-0576-01	3000063261
3000059918	30000599180510	D684-00-0576-01U	3000063260
3000059918	30000599180510	D684-00-0576-01	3000063261
3000059918	30000599180500	D684-00-0576-01U	3000063260
3000059918	30000599180500	D684-00-0576-01	3000063261
3000059918	30000599180497	D684-00-0576-01U	3000063260
3000059918	30000599180497	D684-00-0576-01	3000063261
3000059919	30000599190081	D684-00-0295-01	3000064257
3000059919	30000599190081	D684-00-0295-01U	3000064258
3000059919	30000599190083	D684-00-0295-01	3000064257
3000059919	30000599190083	D684-00-0295-01U	3000064258
3000059919	30000599190091	D684-00-0295-01	3000064257
3000059919	30000599190091	D684-00-0295-01U	3000064258
3000059919	30000599190092	D684-00-0295-01	3000064257
3000059919	30000599190092	D684-00-0295-01U	3000064258
3000059919	30000599190098	D684-00-0295-01	3000064257
3000059919	30000599190098	D684-00-0295-01U	3000064258
3000059924	30000599240103	D684-00-0480-01	3000064057
3000059924	30000599240103	D684-00-0480-01U	3000064097
3000059924	30000599240103	D684-00-0480-01C	3000066141
3000059924	30000599240103	D684-00-0480-01	3000071409
3000059924	30000599240103	D684-DO-0475	0
3000059924	30000599240109	D684-00-0480-01	3000064057

3000059924	30000599240109	D684-00-0480-01U	3000064097
3000059924	30000599240109	D684-00-0480-01C	3000066141
3000059924	30000599240109	D684-00-0480-01	3000071409
3000059924	30000599240109	D684-DO-0475	0
3000059924	30000599240113	D684-00-0480-01	3000064057
3000059924	30000599240113	D684-00-0480-01U	3000064097
3000059924	30000599240113	D684-00-0480-01C	3000066141
3000059924	30000599240113	D684-00-0480-01	3000071409
3000059924	30000599240113	D684-DO-0475	0
3000059924	30000599240118	D684-00-0480-01	3000064057
3000059924	30000599240118	D684-00-0480-01U	3000064097
3000059924	30000599240118	D684-00-0480-01C	3000066141
3000059924	30000599240118	D684-00-0480-01	3000071409
3000059924	30000599240118	D684-DO-0475	0
3000059924	30000599240124	D684-00-0480-01	3000064057
3000059924	30000599240124	D684-00-0480-01U	3000064097
3000059924	30000599240124	D684-00-0480-01C	3000066141
3000059924	30000599240124	D684-00-0480-01	3000071409
3000059924	30000599240124	D684-DO-0475	0
3000060579	30000605790141	D684-00-0295-02	3000064844
3000060579	30000605790150	D684-00-0295-02	3000064844
3000060579	30000605790156	D684-00-0295-02	3000064844
3000060581	30000605810345	D684-00-0568-01	3000063262
3000060581	30000605810345	D684-00-0568-01U	3000065031
3000060581	30000605810350	D684-00-0568-01	3000063262
3000060581	30000605810350	D684-00-0568-01U	3000065031
3000060581	30000605810355	D684-00-0568-01	3000063262
3000060581	30000605810355	D684-00-0568-01U	3000065031
3000060581	30000605810361	D684-00-0568-01	3000063262
3000060581	30000605810361	D684-00-0568-01U	3000065031
3000060581	30000605810363	D684-00-0568-01	3000063262
3000060581	30000605810363	D684-00-0568-01U	3000065031
3000060582	30000605820352	D684-00-0296-01U	3000064248
3000060582	30000605820352	D684-00-0296-01	3000064249
3000060582	30000605820352	D684-00-0296-02U	3000064250
3000060582	30000605820352	D684-00-0296-02	3000064842
3000060582	30000605820352	D684-00-0296-02	3000064843
3000060582	30000605820352	D684-00-0296-01U	3000066798
3000060582	30000605820352	D684-00-0296-01	3000066800
3000060582	30000605820352	D684-00-0296-02	3000066801
3000060582	30000605820352	D684-00-0296-02	3000066802
3000060582	30000605820352	D684-00-0296-02	3000066803
3000060582	30000605820354	D684-00-0296-01U	3000064248
3000060582	30000605820354	D684-00-0296-01	3000064249
3000060582	30000605820354	D684-00-0296-02U	3000064250
3000060582	30000605820354	D684-00-0296-02	3000064842
3000060582	30000605820354	D684-00-0296-02	3000064843
3000060582	30000605820354	D684-00-0296-01U	3000066798
3000060582	30000605820354	D684-00-0296-01	3000066800
3000060582	30000605820354	D684-00-0296-02	3000066801
3000060582	30000605820354	D684-00-0296-02	3000066802
3000060582	30000605820354	D684-00-0296-02	3000066803
3000060582	30000605820356	D684-00-0296-01U	3000064248
3000060582	30000605820356	D684-00-0296-01	3000064249
3000060582	30000605820356	D684-00-0296-02U	3000064250
3000060582	30000605820356	D684-00-0296-02	3000064842
3000060582	30000605820356	D684-00-0296-02	3000064843
3000060582	30000605820356	D684-00-0296-01U	3000066798

3000060582	30000605820356	D684-00-0296-01	3000066800
3000060582	30000605820356	D684-00-0296-02	3000066801
3000060582	30000605820356	D684-00-0296-02	3000066802
3000060582	30000605820356	D684-00-0296-02	3000066803
3000060582	30000605820357	D684-00-0296-01U	3000064248
3000060582	30000605820357	D684-00-0296-01	3000064249
3000060582	30000605820357	D684-00-0296-02U	3000064250
3000060582	30000605820357	D684-00-0296-02	3000064842
3000060582	30000605820357	D684-00-0296-02	3000064843
3000060582	30000605820357	D684-00-0296-01U	3000066798
3000060582	30000605820357	D684-00-0296-01	3000066800
3000060582	30000605820357	D684-00-0296-02	3000066801
3000060582	30000605820357	D684-00-0296-02	3000066802
3000060582	30000605820357	D684-00-0296-02	3000066803
3000060582	30000605820359	D684-00-0296-01U	3000064248
3000060582	30000605820359	D684-00-0296-01	3000064249
3000060582	30000605820359	D684-00-0296-02U	3000064250
3000060582	30000605820359	D684-00-0296-02	3000064842
3000060582	30000605820359	D684-00-0296-02	3000064843
3000060582	30000605820359	D684-00-0296-01U	3000066798
3000060582	30000605820359	D684-00-0296-01	3000066800
3000060582	30000605820359	D684-00-0296-02	3000066801
3000060582	30000605820359	D684-00-0296-02	3000066802
3000060582	30000605820359	D684-00-0296-02	3000066803
3000060738	30000607380150	D684-00-0608	3000067701
3000060738	30000607380150	D684-00-0608	3000069728
3000060738	30000607380150	D684-00-0608	3000069729
3000060738	30000607380150	D684-00-0608	3000069730
3000060738	30000607380153	D684-00-0608	3000067701
3000060738	30000607380153	D684-00-0608	3000069728
3000060738	30000607380153	D684-00-0608	3000069729
3000060738	30000607380153	D684-00-0608	3000069730
3000060738	30000607380158	D684-00-0608	3000067701
3000060738	30000607380158	D684-00-0608	3000069728
3000060738	30000607380158	D684-00-0608	3000069729
3000060738	30000607380158	D684-00-0608	3000069730
3000060738	30000607380161	D684-00-0608	3000067701
3000060738	30000607380161	D684-00-0608	3000069728
3000060738	30000607380161	D684-00-0608	3000069729
3000060738	30000607380161	D684-00-0608	3000069730
3000060738	30000607380170	D684-00-0608	3000067701
3000060738	30000607380170	D684-00-0608	3000069728
3000060738	30000607380170	D684-00-0608	3000069729
3000060738	30000607380170	D684-00-0608	3000069730
3000060742	30000607420401	D684-00-0479-01C	3000064055
3000060742	30000607420401	D684-00-0479-02U	3000064514
3000060742	30000607420401	D684-00-0479-01U	3000065130
3000060742	30000607420401	D684-00-0479-02U	3000065131
3000060742	30000607420401	D684-00-0479-02	3000065132
3000060742	30000607420401	D684-00-0479-01	3000065133
3000060742	30000607420401	D684-00-0479-01	3000065724
3000060742	30000607420404	D684-00-0479-01C	3000064055
3000060742	30000607420404	D684-00-0479-02U	3000064514
3000060742	30000607420404	D684-00-0479-01U	3000065130
3000060742	30000607420404	D684-00-0479-02U	3000065131
3000060742	30000607420404	D684-00-0479-02	3000065132
3000060742	30000607420404	D684-00-0479-01	3000065133
3000060742	30000607420404	D684-00-0479-01	3000065724

[illegible]

3000061114	30000611140225	D684-00-0560-01	3000068107
3000061114	30000611140231	D684-00-0560-01	3000064726
3000061114	30000611140231	D684-00-0560-01	3000068107
3000061316	30000613160153	D684-00-0545-01	3000064523
3000061316	30000613160153	D684-00-0545-01	3000066986
3000061316	30000613160159	D684-00-0545-01	3000064523
3000061316	30000613160159	D684-00-0545-01	3000066986
3000061316	30000613160163	D684-00-0545-01	3000064523
3000061316	30000613160163	D684-00-0545-01	3000066986
3000061316	30000613160172	D684-00-0545-01	3000064523
3000061316	30000613160172	D684-00-0545-01	3000066986
3000061316	30000613160173	D684-00-0545-01	3000064523
3000061316	30000613160173	D684-00-0545-01	3000066986
3000061317	30000613170194	D684-00-0546-01	3000064524
3000061317	30000613170194	D684-00-0546-01	3000066989
3000061317	30000613170194	D684-00-0546-01	3000066990
3000061317	30000613170194	D684-00-0546-01	3000067890
3000061317	30000613170187	D684-00-0546-01	3000064524
3000061317	30000613170187	D684-00-0546-01	3000066989
3000061317	30000613170187	D684-00-0546-01	3000066990
3000061317	30000613170187	D684-00-0546-01	3000067890
3000061317	30000613170188	D684-00-0546-01	3000064524
3000061317	30000613170188	D684-00-0546-01	3000066989
3000061317	30000613170188	D684-00-0546-01	3000066990
3000061317	30000613170188	D684-00-0546-01	3000067890
3000061317	30000613170196	D684-00-0546-01	3000064524
3000061317	30000613170196	D684-00-0546-01	3000066989
3000061317	30000613170196	D684-00-0546-01	3000066990
3000061317	30000613170196	D684-00-0546-01	3000067890
3000061317	30000613170185	D684-00-0546-01	3000064524
3000061317	30000613170185	D684-00-0546-01	3000066989
3000061317	30000613170185	D684-00-0546-01	3000066990
3000061317	30000613170185	D684-00-0546-01	3000067890
3000061318	30000613180212	D684-00-0561-01	3000064723
3000061318	30000613180212	D684-00-0561-01	3000072718
3000061318	30000613180221	D684-00-0561-01	3000064723
3000061318	30000613180221	D684-00-0561-01	3000072718
3000061318	30000613180228	D684-00-0561-01	3000064723
3000061318	30000613180228	D684-00-0561-01	3000072718
3000061318	30000613180234	D684-00-0561-01	3000064723
3000061318	30000613180234	D684-00-0561-01	3000072718
3000061318	30000613180235	D684-00-0561-01	3000064723
3000061318	30000613180235	D684-00-0561-01	3000072718
3000061337	30000613370143	D684-00-0478-01	3000066195
3000061337	30000613370143	D684-00-0478-01U	3000066264
3000061337	30000613370143	D684-00-0478-01	3000067965
3000061337	30000613370143	D684-00-0478-02	3000068728
3000061337	30000613370143	D684-00-0478-01U	3000069270
3000061337	30000613370143	D684-00-0478-01	3000069878
3000061337	30000613370149	D684-00-0478-01	3000066195
3000061337	30000613370149	D684-00-0478-01U	3000066264
3000061337	30000613370149	D684-00-0478-01	3000067965
3000061337	30000613370149	D684-00-0478-02	3000068728
3000061337	30000613370149	D684-00-0478-01U	3000069270
3000061337	30000613370149	D684-00-0478-01	3000069878
3000061337	30000613370142	D684-00-0478-01	3000066195
3000061337	30000613370142	D684-00-0478-01U	3000066264
3000061337	30000613370142	D684-00-0478-01	3000067965

[illegible]

3000061340	30000613400276	D684-00-0296-02U	3000066804
3000061340	30000613400286	D684-00-0296-01	3000065563
3000061340	30000613400286	D684-00-0296-01U	3000065636
3000061340	30000613400286	D684-DO-0497	0
3000061340	30000613400286	D684-00-0296-02U	3000066804
3000061340	30000613400290	D684-00-0296-01	3000065563
3000061340	30000613400290	D684-00-0296-01U	3000065636
3000061340	30000613400290	D684-DO-0497	0
3000061340	30000613400290	D684-00-0296-02U	3000066804
3000061340	30000613400293	D684-00-0296-01	3000065563
3000061340	30000613400293	D684-00-0296-01U	3000065636
3000061340	30000613400293	D684-DO-0497	0
3000061340	30000613400293	D684-00-0296-02U	3000066804
3000061340	30000613400298	D684-00-0296-01	3000065563
3000061340	30000613400298	D684-00-0296-01U	3000065636
3000061340	30000613400298	D684-DO-0497	0
3000061340	30000613400298	D684-00-0296-02U	3000066804
3000061342	30000613420031	D684-00-0480-02U	3000065029
3000061342	30000613420031	D684-00-0480-01U	3000065030
3000061342	30000613420031	D684-00-0480-02	3000065235
3000061342	30000613420036	D684-00-0480-02U	3000065029
3000061342	30000613420036	D684-00-0480-01U	3000065030
3000061342	30000613420036	D684-00-0480-02	3000065235
3000061342	30000613420038	D684-00-0480-02U	3000065029
3000061342	30000613420038	D684-00-0480-01U	3000065030
3000061342	30000613420038	D684-00-0480-02	3000065235
3000061342	30000613420044	D684-00-0480-02U	3000065029
3000061342	30000613420044	D684-00-0480-01U	3000065030
3000061342	30000613420044	D684-00-0480-02	3000065235
3000061342	30000613420048	D684-00-0480-02U	3000065029
3000061342	30000613420048	D684-00-0480-01U	3000065030
3000061342	30000613420048	D684-00-0480-02	3000065235
3000062043	30000620430270	D684-DC-0567	0
3000062043	30000620430270	D684-00-0568-01U	3000065564
3000062043	30000620430270	D684-00-0568-01U	3000065634
3000062043	30000620430284	D684-DC-0567	0
3000062043	30000620430284	D684-00-0568-01U	3000065564
3000062043	30000620430284	D684-00-0568-01U	3000065634
3000062043	30000620430266	D684-DC-0567	0
3000062043	30000620430266	D684-00-0568-01U	3000065564
3000062043	30000620430266	D684-00-0568-01U	3000065634
3000062043	30000620430272	D684-DC-0567	0
3000062043	30000620430272	D684-00-0568-01U	3000065564
3000062043	30000620430272	D684-00-0568-01U	3000065634
3000062043	30000620430274	D684-DC-0567	0
3000062043	30000620430274	D684-00-0568-01U	3000065564
3000062043	30000620430274	D684-00-0568-01U	3000065634
3000062257	30000622570106	D684-00-0576-01U	3000065402
3000062257	30000622570106	D684-00-0576-01U	3000065402
3000062257	30000622570106	D684-00-0576-01	3000065565
3000062257	30000622570106	D684-00-0576-01	3000066085
3000062257	30000622570107	D684-00-0576-01U	3000065402
3000062257	30000622570107	D684-00-0576-01U	3000065402
3000062257	30000622570107	D684-00-0576-01	3000065565
3000062257	30000622570107	D684-00-0576-01	3000066085
3000062257	30000622570108	D684-00-0576-01U	3000065402
3000062257	30000622570108	D684-00-0576-01U	3000065402
3000062257	30000622570108	D684-00-0576-01	3000065565

[illegible]

[illegible]

[illegible]

3000063955	30000639550258	D684-00-0479-02U	3000066767
3000063955	30000639550258	D684-00-0475	3000075861
3000063955	30000639550259	D684-00-0479-01	3000066487
3000063955	30000639550259	D684-00-0479-01C	3000066513
3000063955	30000639550259	D684-00-0479-01C	3000066514
3000063955	30000639550259	D684-00-0479-02	3000066692
3000063955	30000639550259	D684-00-0479-01U	3000066766
3000063955	30000639550259	D684-00-0479-02U	3000066767
3000063955	30000639550259	D684-00-0475	3000075861
3000063958	30000639581096	D684-00-0470-01U	3000067692
3000063958	30000639581096	D684-00-0470-01	3000068730
3000063958	30000639581096	D684-00-0470-01	3000069160
3000063958	30000639581096	D684-00-0470-01U	3000069161
3000063958	30000639581096	D684-00-0470-07	3000069162
3000063958	30000639581096	D684-00-0470-01	3000069376
3000063958	30000639581096	D684-00-0470-07	3000069442
3000063958	30000639581096	D684-00-0470-01U	3000069872
3000063958	30000639581096	D684-00-0470-01U	3000069873
3000063958	30000639581098	D684-00-0470-01U	3000067692
3000063958	30000639581098	D684-00-0470-01	3000068730
3000063958	30000639581098	D684-00-0470-01	3000069160
3000063958	30000639581098	D684-00-0470-01U	3000069161
3000063958	30000639581098	D684-00-0470-07	3000069162
3000063958	30000639581098	D684-00-0470-01	3000069376
3000063958	30000639581098	D684-00-0470-07	3000069442
3000063958	30000639581098	D684-00-0470-01U	3000069872
3000063958	30000639581098	D684-00-0470-01U	3000069873
3000063958	30000639581101	D684-00-0470-01U	3000067692
3000063958	30000639581101	D684-00-0470-01	3000068730
3000063958	30000639581101	D684-00-0470-01	3000069160
3000063958	30000639581101	D684-00-0470-01U	3000069161
3000063958	30000639581101	D684-00-0470-07	3000069162
3000063958	30000639581101	D684-00-0470-01	3000069376
3000063958	30000639581101	D684-00-0470-07	3000069442
3000063958	30000639581101	D684-00-0470-01U	3000069872
3000063958	30000639581101	D684-00-0470-01U	3000069873
3000063958	30000639581106	D684-00-0470-01U	3000067692
3000063958	30000639581106	D684-00-0470-01	3000068730
3000063958	30000639581106	D684-00-0470-01	3000069160
3000063958	30000639581106	D684-00-0470-01U	3000069161
3000063958	30000639581106	D684-00-0470-07	3000069162
3000063958	30000639581106	D684-00-0470-01	3000069376
3000063958	30000639581106	D684-00-0470-07	3000069442
3000063958	30000639581106	D684-00-0470-01U	3000069872
3000063958	30000639581106	D684-00-0470-01U	3000069873
3000063958	30000639581113	D684-00-0470-01U	3000067692
3000063958	30000639581113	D684-00-0470-01	3000068730
3000063958	30000639581113	D684-00-0470-01	3000069160
3000063958	30000639581113	D684-00-0470-01U	3000069161
3000063958	30000639581113	D684-00-0470-07	3000069162
3000063958	30000639581113	D684-00-0470-01	3000069376
3000063958	30000639581113	D684-00-0470-07	3000069442
3000063958	30000639581113	D684-00-0470-01U	3000069872
3000063958	30000639581113	D684-00-0470-01U	3000069873
3000064995	30000649950021	D684-00-0480-01U	3000067497
3000064995	30000649950021	D684-00-0480-01	3000067500
3000064995	30000649950021	D684-00-0480-02	3000067531
3000064995	30000649950021	D684-00-0480-02	3000067532

[illegible]

[illegible]

[illegible]

3000065939	30000659390049	D684-00-0567	3000096709
3000065939	30000659390050	D684-00-0294-01U	3000067461
3000065939	30000659390050	D684-00-0294-01	3000067462
3000065939	30000659390050	D684-00-0294-02	3000067463
3000065939	30000659390050	D684-00-0294-01U	3000067855
3000065939	30000659390050	D684-00-0294-02U	3000067856
3000065939	30000659390050	D684-00-0294-01U	3000071981
3000065939	30000659390050	D684-00-0549-01	3000089612
3000065939	30000659390050	D684-00-0567	3000096709
3000065939	30000659390051	D684-00-0294-01U	3000067461
3000065939	30000659390051	D684-00-0294-01	3000067462
3000065939	30000659390051	D684-00-0294-02	3000067463
3000065939	30000659390051	D684-00-0294-01U	3000067855
3000065939	30000659390051	D684-00-0294-02U	3000067856
3000065939	30000659390051	D684-00-0294-01U	3000071981
3000065939	30000659390051	D684-00-0549-01	3000089612
3000065939	30000659390051	D684-00-0567	3000096709
3000065941	30000659410562	D684-00-0568-01	3000068522
3000065941	30000659410562	D684-00-0568-01U	3000068523
3000065941	30000659410562	D684-00-0568-01U	3000069088
3000065941	30000659410563	D684-00-0568-01	3000068522
3000065941	30000659410563	D684-00-0568-01U	3000068523
3000065941	30000659410563	D684-00-0568-01U	3000069088
3000065941	30000659410565	D684-00-0568-01	3000068522
3000065941	30000659410565	D684-00-0568-01U	3000068523
3000065941	30000659410565	D684-00-0568-01U	3000069088
3000065941	30000659410567	D684-00-0568-01	3000068522
3000065941	30000659410567	D684-00-0568-01U	3000068523
3000065941	30000659410567	D684-00-0568-01U	3000069088
3000065941	30000659410568	D684-00-0568-01	3000068522
3000065941	30000659410568	D684-00-0568-01U	3000068523
3000065941	30000659410568	D684-00-0568-01U	3000069088
3000065944	30000659440257	D684-00-0480-01	3000068688
3000065944	30000659440257	D684-00-0480-01	3000068689
3000065944	30000659440257	D684-00-0480-01U	3000068770
3000065944	30000659440254	D684-00-0480-01	3000068688
3000065944	30000659440254	D684-00-0480-01	3000068689
3000065944	30000659440254	D684-00-0480-01U	3000068770
3000065944	30000659440256	D684-00-0480-01	3000068688
3000065944	30000659440256	D684-00-0480-01	3000068689
3000065944	30000659440256	D684-00-0480-01U	3000068770
3000065944	30000659440258	D684-00-0480-01	3000068688
3000065944	30000659440258	D684-00-0480-01	3000068689
3000065944	30000659440258	D684-00-0480-01U	3000068770
3000065944	30000659440259	D684-00-0480-01	3000068688
3000065944	30000659440259	D684-00-0480-01	3000068689
3000065944	30000659440259	D684-00-0480-01U	3000068770
3000065967	30000659670015	D684-00-0607	3000069727
3000065967	30000659670016	D684-00-0607	3000069727
3000065967	30000659670022	D684-00-0607	3000069727
3000065967	30000659670012	D684-00-0607	3000069727
3000065967	30000659670010	D684-00-0607	3000069727
3000065968	30000659680431	D684-00-0480-02	3000068731
3000065968	30000659680431	D684-00-0480-02U	3000069316
3000065968	30000659680431	D684-00-0480-01U	3000069317
3000065968	30000659680431	D684-00-0480-02	3000069377
3000065968	30000659680431	D684-00-0480-01C	3000069656
3000065968	30000659680428	D684-00-0480-02	3000068731

[illegible]

3000066273	30000662730110	D684-00-0576-01U	3000069363
3000066273	30000662730110	D684-00-0576-01U	3000069363
3000066273	30000662730110	D684-00-0576-01	3000069364
3000066273	30000662730112	D684-00-0576-01	3000069378
3000066273	30000662730112	D684-00-0576-01U	3000069363
3000066273	30000662730112	D684-00-0576-01U	3000069363
3000066273	30000662730112	D684-00-0576-01	3000069364
3000066273	30000662730125	D684-00-0576-01	3000069378
3000066273	30000662730125	D684-00-0576-01U	3000069363
3000066273	30000662730125	D684-00-0576-01U	3000069363
3000066273	30000662730125	D684-00-0576-01	3000069364
3000066502	30000665020359	D684-00-0479-02	3000068729
3000066502	30000665020359	D684-00-0479-01U	3000068923
3000066502	30000665020359	D684-00-0479-02	3000069123
3000066502	30000665020359	D684-00-0479-07	3000069441
3000066502	30000665020359	D684-00-0479-07	3000069655
3000066502	30000665020359	D684-00-0479-01	3000069877
3000066502	30000665020354	D684-00-0479-02	3000068729
3000066502	30000665020354	D684-00-0479-01U	3000068923
3000066502	30000665020354	D684-00-0479-02	3000069123
3000066502	30000665020354	D684-00-0479-07	3000069441
3000066502	30000665020354	D684-00-0479-07	3000069655
3000066502	30000665020354	D684-00-0479-01	3000069877
3000066502	30000665020375	D684-00-0479-02	3000068729
3000066502	30000665020375	D684-00-0479-01U	3000068923
3000066502	30000665020375	D684-00-0479-02	3000069123
3000066502	30000665020375	D684-00-0479-07	3000069441
3000066502	30000665020375	D684-00-0479-07	3000069655
3000066502	30000665020375	D684-00-0479-01	3000069877
3000066502	30000665020351	D684-00-0479-02	3000068729
3000066502	30000665020351	D684-00-0479-01U	3000068923
3000066502	30000665020351	D684-00-0479-02	3000069123
3000066502	30000665020351	D684-00-0479-07	3000069441
3000066502	30000665020351	D684-00-0479-07	3000069655
3000066502	30000665020351	D684-00-0479-01	3000069877
3000066502	30000665020367	D684-00-0479-02	3000068729
3000066502	30000665020367	D684-00-0479-01U	3000068923
3000066502	30000665020367	D684-00-0479-02	3000069123
3000066502	30000665020367	D684-00-0479-07	3000069441
3000066502	30000665020367	D684-00-0479-07	3000069655
3000066502	30000665020367	D684-00-0479-01	3000069877
3000066505	30000665050374	D684-00-0568-01	3000069870
3000066505	30000665050374	D684-00-0568-01U	3000069871
3000066505	30000665050372	D684-00-0568-01	3000069870
3000066505	30000665050372	D684-00-0568-01U	3000069871
3000066505	30000665050359	D684-00-0568-01	3000069870
3000066505	30000665050359	D684-00-0568-01U	3000069871
3000066505	30000665050365	D684-00-0568-01	3000069870
3000066505	30000665050365	D684-00-0568-01U	3000069871
3000066505	30000665050377	D684-00-0568-01	3000069870
3000066505	30000665050377	D684-00-0568-01U	3000069871
3000066506	30000665060036	D684-00-0479-01C	3000070407
3000066506	30000665060036	D684-00-0479-01U	3000070426
3000066506	30000665060036	D684-00-0479-02	3000070427
3000066506	30000665060036	D684-00-0479-01U	3000070663
3000066506	30000665060036	D684-00-0479-07	3000071127
3000066506	30000665060036	D684-0	

[illegible]

3000066985	30000669850334	D684-00-0480-01	3000070070
3000066985	30000669850334	D684-00-0480-01U	3000070409
3000066985	30000669850339	D684-00-0480-01C	3000069657
3000066985	30000669850339	D684-00-0480-01	3000070070
3000066985	30000669850339	D684-00-0480-01U	3000070409
3000067125	30000671250043	D684-00-0469-01U	3000069876
3000067125	30000671250043	D684-00-0469-01C	3000070665
3000067125	30000671250043	D684-00-0469-01U	3000071218
3000067125	30000671250043	D684-00-0469-01	3000071413
3000067125	30000671250043	D684-00-0469-01U	3000071977
3000067125	30000671250043	D684-00-0469-01U	3000071978
3000067125	30000671250030	D684-00-0469-01U	3000069876
3000067125	30000671250030	D684-00-0469-01C	3000070665
3000067125	30000671250030	D684-00-0469-01U	3000071218
3000067125	30000671250030	D684-00-0469-01	3000071413
3000067125	30000671250030	D684-00-0469-01U	3000071977
3000067125	30000671250030	D684-00-0469-01U	3000071978
3000067125	30000671250029	D684-00-0469-01U	3000069876
3000067125	30000671250029	D684-00-0469-01C	3000070665
3000067125	30000671250029	D684-00-0469-01U	3000071218
3000067125	30000671250029	D684-00-0469-01	3000071413
3000067125	30000671250029	D684-00-0469-01U	3000071977
3000067125	30000671250029	D684-00-0469-01U	3000071978
3000067125	30000671250041	D684-00-0469-01U	3000069876
3000067125	30000671250041	D684-00-0469-01C	3000070665
3000067125	30000671250041	D684-00-0469-01U	3000071218
3000067125	30000671250041	D684-00-0469-01	3000071413
3000067125	30000671250041	D684-00-0469-01U	3000071977
3000067125	30000671250041	D684-00-0469-01U	3000071978
3000067125	30000671250050	D684-00-0469-01U	3000069876
3000067125	30000671250050	D684-00-0469-01C	3000070665
3000067125	30000671250050	D684-00-0469-01U	3000071218
3000067125	30000671250050	D684-00-0469-01	3000071413
3000067125	30000671250050	D684-00-0469-01U	3000071977
3000067125	30000671250050	D684-00-0469-01U	3000071978
3000067128	30000671280156	D684-00-0576-01U	3000070067
3000067128	30000671280156	D684-00-0576-01	3000070068
3000067128	30000671280164	D684-00-0576-01U	3000070067
3000067128	30000671280164	D684-00-0576-01	3000070068
3000067128	30000671280168	D684-00-0576-01U	3000070067
3000067128	30000671280168	D684-00-0576-01	3000070068
3000067128	30000671280176	D684-00-0576-01U	3000070067
3000067128	30000671280176	D684-00-0576-01	3000070068
3000067128	30000671280177	D684-00-0576-01U	3000070067
3000067128	30000671280177	D684-00-0576-01	3000070068
3000067129	30000671290067	D684-00-0294-01U	3000071219
3000067129	30000671290067	D684-00-0294-01	3000071220
3000067129	30000671290067	D684-00-0294-01U	3000071982
3000067129	30000671290067	D684-00-0294-02U	3000071983
3000067129	30000671290052	D684-00-0294-01U	3000071219
3000067129	30000671290052	D684-00-0294-01	3000071220
3000067129	30000671290052	D684-00-0294-01U	3000071982
3000067129	30000671290052	D684-00-0294-02U	3000071983
3000067129	30000671290056	D684-00-0294-01U	3000071219
3000067129	30000671290056	D684-00-0294-01	3000071220
3000067129	30000671290056	D684-00-0294-01U	3000071982
3000067129	30		

3000067129	30000671290072	D684-00-0294-01	3000071220
3000067129	30000671290072	D684-00-0294-01U	3000071982
3000067129	30000671290072	D684-00-0294-02U	3000071983
3000067129	30000671290053	D684-00-0294-01U	3000071219
3000067129	30000671290053	D684-00-0294-01	3000071220
3000067129	30000671290053	D684-00-0294-01U	3000071982
3000067129	30000671290053	D684-00-0294-02U	3000071983
3000067130	30000671300089	D684-00-0295-05	3000069443
3000067130	30000671300094	D684-00-0295-05	3000069443
3000067130	30000671300080	D684-00-0295-05	3000069443
3000067130	30000671300079	D684-00-0295-05	3000069443
3000067130	30000671300096	D684-00-0295-05	3000069443
3000067516	30000675160171	D684-00-0296-01	3000070668
3000067516	30000675160171	D684-00-0296-01U	3000070669
3000067516	30000675160171	D684-00-0296-02U	3000070670
3000067516	30000675160171	D684-00-0296-01	3000071291
3000067516	30000675160171	D684-00-0296-01U	3000071292
3000067516	30000675160171	D684-00-0498-07	3000071836
3000067516	30000675160171	D684-00-0296-01	3000071979
3000067516	30000675160171	D684-00-0296-01U	3000071980
3000067516	30000675160175	D684-00-0296-01	3000070668
3000067516	30000675160175	D684-00-0296-01U	3000070669
3000067516	30000675160175	D684-00-0296-02U	3000070670
3000067516	30000675160175	D684-00-0296-01	3000071291
3000067516	30000675160175	D684-00-0296-01U	3000071292
3000067516	30000675160175	D684-00-0498-07	3000071836
3000067516	30000675160175	D684-00-0296-01	3000071979
3000067516	30000675160175	D684-00-0296-01U	3000071980
3000067516	30000675160177	D684-00-0296-01	3000070668
3000067516	30000675160177	D684-00-0296-01U	3000070669
3000067516	30000675160177	D684-00-0296-02U	3000070670
3000067516	30000675160177	D684-00-0296-01	3000071291
3000067516	30000675160177	D684-00-0296-01U	3000071292
3000067516	30000675160177	D684-00-0498-07	3000071836
3000067516	30000675160177	D684-00-0296-01	3000071979
3000067516	30000675160177	D684-00-0296-01U	3000071980
3000067516	30000675160178	D684-00-0296-01	3000070668
3000067516	30000675160178	D684-00-0296-01U	3000070669
3000067516	30000675160178	D684-00-0296-02U	3000070670
3000067516	30000675160178	D684-00-0296-01	3000071291
3000067516	30000675160178	D684-00-0296-01U	3000071292
3000067516	30000675160178	D684-00-0498-07	3000071836
3000067516	30000675160178	D684-00-0296-01	3000071979
3000067516	30000675160178	D684-00-0296-01U	3000071980
3000067516	30000675160179	D684-00-0296-01	3000070668
3000067516	30000675160179	D684-00-0296-01U	3000070669
3000067516	30000675160179	D684-00-0296-02U	3000070670
3000067516	30000675160179	D684-00-0296-01	3000071291
3000067516	30000675160179	D684-00-0296-01U	3000071292
3000067516	30000675160179	D684-00-0498-07	3000071836
3000067516	30000675160179	D684-00-0296-01	3000071979
3000067516	30000675160179	D684-00-0296-01U	3000071980
3000067524	30000675240051	D684-00-0608	3000071288
3000067524	30000675240051	D684-00-0608	3000073187
3000067524	30000675240053	D684-00-0608	3000071288
3000067524	30000675240053	D684-00-0608	3000073187
3000067524	30000675240056	D684-00-0608	

3000067524	30000675240066	D684-00-0608	3000071288
3000067524	30000675240066	D684-00-0608	3000073187
3000067524	30000675240068	D684-00-0608	3000071288
3000067524	30000675240068	D684-00-0608	3000073187
3000067525	30000675250489	D684-00-0568-01U	3000070817
3000067525	30000675250499	D684-00-0568-01U	3000070817
3000067525	30000675250508	D684-00-0568-01U	3000070817
3000067525	30000675250510	D684-00-0568-01U	3000070817
3000067525	30000675250518	D684-00-0568-01U	3000070817
3000067526	30000675260541	D684-00-0480-01C	3000070408
3000067526	30000675260541	D684-00-0480-01U	3000070410
3000067526	30000675260541	D684-00-0480-02U	3000070488
3000067526	30000675260541	D684-00-0480-02	3000070489
3000067526	30000675260541	D684-00-0480-01C	3000070664
3000067526	30000675260542	D684-00-0480-01C	3000070408
3000067526	30000675260542	D684-00-0480-01U	3000070410
3000067526	30000675260542	D684-00-0480-02U	3000070488
3000067526	30000675260542	D684-00-0480-02	3000070489
3000067526	30000675260542	D684-00-0480-01C	3000070664
3000067526	30000675260543	D684-00-0480-01C	3000070408
3000067526	30000675260543	D684-00-0480-01U	3000070410
3000067526	30000675260543	D684-00-0480-02U	3000070488
3000067526	30000675260543	D684-00-0480-02	3000070489
3000067526	30000675260543	D684-00-0480-01C	3000070664
3000067526	30000675260546	D684-00-0480-01C	3000070408
3000067526	30000675260546	D684-00-0480-01U	3000070410
3000067526	30000675260546	D684-00-0480-02U	3000070488
3000067526	30000675260546	D684-00-0480-02	3000070489
3000067526	30000675260546	D684-00-0480-01C	3000070664
3000067526	30000675260548	D684-00-0480-01C	3000070408
3000067526	30000675260548	D684-00-0480-01U	3000070410
3000067526	30000675260548	D684-00-0480-02U	3000070488
3000067526	30000675260548	D684-00-0480-02	3000070489
3000067526	30000675260548	D684-00-0480-01C	3000070664
3000067805	30000678050006	D684-00-0545-01	3000077732
3000067805	30000678050006	D684-00-0545-01	3000080305
3000067805	30000678050006	D684-00-0545-01	3000083332
3000067805	30000678050008	D684-00-0545-01	3000077732
3000067805	30000678050008	D684-00-0545-01	3000080305
3000067805	30000678050008	D684-00-0545-01	3000083332
3000067805	30000678050009	D684-00-0545-01	3000077732
3000067805	30000678050009	D684-00-0545-01	3000080305
3000067805	30000678050009	D684-00-0545-01	3000083332
3000067805	30000678050021	D684-00-0545-01	3000077732
3000067805	30000678050021	D684-00-0545-01	3000080305
3000067805	30000678050021	D684-00-0545-01	3000083332
3000067805	30000678050022	D684-00-0545-01	3000077732
3000067805	30000678050022	D684-00-0545-01	3000080305
3000067805	30000678050022	D684-00-0545-01	3000083332
3000067893	30000678930266	D684-00-0470-01U	3000071217
3000067893	30000678930266	D684-00-0470-01	3000071293
3000067893	30000678930266	D684-00-0470-01	3000071293
3000067893	30000678930266	D684-00-0470-01	3000071739
3000067893	30000678930266	D684-00-0470-07	3000071837
3000067893	30000678930266	D684-00-0470-01U	3000071891
3000067893	30000678930265	D684-00-0470-01U	3000071217
3000067893	30000678930265	D684-00-0470-01	3000071293
3000067893	30000678930265	D684-00-0470-01	3000071293

[illegible]

[illegible]

[illegible]

3000068051	30000680510067	D684-00-0469-01C	3000073189
3000068051	30000680510067	D684-00-0469-01	3000074269
3000068051	30000680510067	D684-00-0469-01U	3000074270
3000068051	30000680510067	D684-00-0469-01U	3000074270
3000068051	30000680510069	D684-00-0469-01	3000072399
3000068051	30000680510069	D684-00-0469-01U	3000073188
3000068051	30000680510069	D684-00-0469-01C	3000073189
3000068051	30000680510069	D684-00-0469-01	3000074269
3000068051	30000680510069	D684-00-0469-01U	3000074270
3000068051	30000680510069	D684-00-0469-01U	3000074270
3000068051	30000680510061	D684-00-0469-01	3000072399
3000068051	30000680510061	D684-00-0469-01U	3000073188
3000068051	30000680510061	D684-00-0469-01C	3000073189
3000068051	30000680510061	D684-00-0469-01	3000074269
3000068051	30000680510061	D684-00-0469-01U	3000074270
3000068051	30000680510061	D684-00-0469-01U	3000074270
3000068051	30000680510058	D684-00-0469-01	3000072399
3000068051	30000680510058	D684-00-0469-01U	3000073188
3000068051	30000680510058	D684-00-0469-01C	3000073189
3000068051	30000680510058	D684-00-0469-01	3000074269
3000068051	30000680510058	D684-00-0469-01U	3000074270
3000068051	30000680510058	D684-00-0469-01U	3000074270
3000068052	30000680520001	D684-00-0607	3000071976
3000068052	30000680520002	D684-00-0607	3000071976
3000068052	30000680520003	D684-00-0607	3000071976
3000068052	30000680520004	D684-00-0607	3000071976
3000068052	30000680520008	D684-00-0607	3000071976
3000068053	30000680530328	D684-00-0568-01	3000072066
3000068053	30000680530328	D684-00-0568-01U	3000072067
3000068053	30000680530328	D684-DC-0567	0
3000068053	30000680530328	D684-DO-0567	0
3000068053	30000680530328	D684-00-0568-01U	3000072404
3000068053	30000680530330	D684-00-0568-01	3000072066
3000068053	30000680530330	D684-00-0568-01U	3000072067
3000068053	30000680530330	D684-DC-0567	0
3000068053	30000680530330	D684-DO-0567	0
3000068053	30000680530330	D684-00-0568-01U	3000072404
3000068053	30000680530331	D684-00-0568-01	3000072066
3000068053	30000680530331	D684-00-0568-01U	3000072067
3000068053	30000680530331	D684-DC-0567	0
3000068053	30000680530331	D684-DO-0567	0
3000068053	30000680530331	D684-00-0568-01U	3000072404
3000068053	30000680530341	D684-00-0568-01	3000072066
3000068053	30000680530341	D684-00-0568-01U	3000072067
3000068053	30000680530341	D684-DC-0567	0
3000068053	30000680530341	D684-DO-0567	0
3000068053	30000680530341	D684-00-0568-01U	3000072404
3000068053	30000680530344	D684-00-0568-01	3000072066
3000068053	30000680530344	D684-00-0568-01U	3000072067
3000068053	30000680530344	D684-DC-0567	0
3000068053	30000680530344	D684-DO-0567	0
3000068053	30000680530344	D684-00-0568-01U	3000072404
3000068141	30000681410060	D684-00-0498-01	3000073404
3000068141	30000681410060	D684-00-0296-01	3000073516
3000068141	30000681410060	D684-00-0296-01U	3000073517
3000068141	30000681410060	D684-00-0296-02U	3000073518
3000068141	30000681410060	D684-00-0296-01U	3000074144
3000068141	30000681410060	D684-00-0296-01U	3000074145

[illegible]

[illegible]

3000069371	30000693710199	D684-00-0295-07	3000073203
3000069371	30000693710199	D684-00-0295-07	3000073203
3000069371	30000693710199	D684-00-0295-01U	3000074729
3000069371	30000693710199	D684-00-0295-01U	3000074729
3000069371	30000693710199	D684-00-0295-02U	3000074731
3000069371	30000693710199	D684-00-0295-01U	3000074850
3000069371	30000693710199	D684-DC-0293	0
3000069371	30000693710199	D684-00-0295-01	3000076925
3000069383	30000693830151	D684-00-0560-01	3000072716
3000069383	30000693830153	D684-00-0560-01	3000072716
3000069383	30000693830158	D684-00-0560-01	3000072716
3000069383	30000693830161	D684-00-0560-01	3000072716
3000069383	30000693830162	D684-00-0560-01	3000072716
3000069384	30000693840075	D684-00-0561-01	3000072719
3000069384	30000693840075	D684-00-0561-01	3000074886
3000069384	30000693840078	D684-00-0561-01	3000072719
3000069384	30000693840078	D684-00-0561-01	3000074886
3000069384	30000693840082	D684-00-0561-01	3000072719
3000069384	30000693840082	D684-00-0561-01	3000074886
3000069384	30000693840085	D684-00-0561-01	3000072719
3000069384	30000693840085	D684-00-0561-01	3000074886
3000069384	30000693840088	D684-00-0561-01	3000072719
3000069384	30000693840088	D684-00-0561-01	3000074886
3000069386	30000693860040	D684-00-0480-02	3000072884
3000069386	30000693860040	D684-00-0480-01	3000072885
3000069386	30000693860040	D684-00-0480-01U	3000072886
3000069386	30000693860041	D684-00-0480-02	3000072884
3000069386	30000693860041	D684-00-0480-01	3000072885
3000069386	30000693860041	D684-00-0480-01U	3000072886
3000069386	30000693860044	D684-00-0480-02	3000072884
3000069386	30000693860044	D684-00-0480-01	3000072885
3000069386	30000693860044	D684-00-0480-01U	3000072886
3000069386	30000693860043	D684-00-0480-02	3000072884
3000069386	30000693860043	D684-00-0480-01	3000072885
3000069386	30000693860043	D684-00-0480-01U	3000072886
3000069386	30000693860048	D684-00-0480-02	3000072884
3000069386	30000693860048	D684-00-0480-01	3000072885
3000069386	30000693860048	D684-00-0480-01U	3000072886
3000069503	30000695030618	D684-00-0576-01	3000072316
3000069503	30000695030618	D684-00-0576-01U	3000072830
3000069503	30000695030620	D684-00-0576-01	3000072316
3000069503	30000695030620	D684-00-0576-01U	3000072830
3000069503	30000695030626	D684-00-0576-01	3000072316
3000069503	30000695030626	D684-00-0576-01U	3000072830
3000069503	30000695030637	D684-00-0576-01	3000072316
3000069503	30000695030637	D684-00-0576-01U	3000072830
3000069503	30000695030639	D684-00-0576-01	3000072316
3000069503	30000695030639	D684-00-0576-01U	3000072830
3000069902	30000699020615	D684-00-0470-01C	3000073190
3000069902	30000699020615	D684-00-0470-01	3000073191
3000069902	30000699020615	D684-00-0470-01U	3000073192
3000069902	30000699020615	D684-00-0470-01	3000073277
3000069902	30000699020615	D684-00-0470-01	3000074271
3000069902	30000699020615	D684-00-0470-01U	3000074272
3000069902	30000699020622	D684-00-0470-01C	3000073190
3000069902	30000699020622	D684-00-0470-01	3000073191
3000069902	30000699020622	D684-00-0470-01U	3000073192
3000069902	30000699020622	D684-00-0470-01	3000073277

[illegible]

3000069903	30000699030195	D684-00-0479-02U	3000074028
3000069903	30000699030195	D684-00-0479-01C	3000074029
3000069903	30000699030198	D684-00-0479-01C	3000073402
3000069903	30000699030198	D684-00-0479-01	3000073519
3000069903	30000699030198	D684-00-0479-01U	3000073520
3000069903	30000699030198	D684-00-0479-02U	3000073521
3000069903	30000699030198	D684-00-0479-07	3000073968
3000069903	30000699030198	D684-00-0479-01U	3000074020
3000069903	30000699030198	D684-00-0479-01	3000074021
3000069903	30000699030198	D684-00-0479-02	3000074027
3000069903	30000699030198	D684-00-0479-02U	3000074028
3000069903	30000699030198	D684-00-0479-01C	3000074029
3000070172	30000701720281	D684-00-0559-01	3000072401
3000070172	30000701720281	D684-00-0559-01	3000072402
3000070172	30000701720281	D684-00-0559-01	3000072402
3000070172	30000701720284	D684-00-0559-01	3000072401
3000070172	30000701720284	D684-00-0559-01	3000072402
3000070172	30000701720284	D684-00-0559-01	3000072402
3000070172	30000701720285	D684-00-0559-01	3000072401
3000070172	30000701720285	D684-00-0559-01	3000072402
3000070172	30000701720285	D684-00-0559-01	3000072402
3000070172	30000701720288	D684-00-0559-01	3000072401
3000070172	30000701720288	D684-00-0559-01	3000072402
3000070172	30000701720288	D684-00-0559-01	3000072402
3000070172	30000701720289	D684-00-0559-01	3000072401
3000070172	30000701720289	D684-00-0559-01	3000072402
3000070172	30000701720289	D684-00-0559-01	3000072402
3000070480	30000704800130	D684-00-0480-01U	3000072916
3000070480	30000704800130	D684-00-0480-01C	3000073065
3000070480	30000704800130	D684-00-0480-02	3000073275
3000070480	30000704800130	D684-DO-0475	0
3000070480	30000704800130	D684-00-0480-01C	3000073403
3000070480	30000704800131	D684-00-0480-01U	3000072916
3000070480	30000704800131	D684-00-0480-01C	3000073065
3000070480	30000704800131	D684-00-0480-02	3000073275
3000070480	30000704800131	D684-DO-0475	0
3000070480	30000704800131	D684-00-0480-01C	3000073403
3000070480	30000704800136	D684-00-0480-01U	3000072916
3000070480	30000704800136	D684-00-0480-01C	3000073065
3000070480	30000704800136	D684-00-0480-02	3000073275
3000070480	30000704800136	D684-DO-0475	0
3000070480	30000704800136	D684-00-0480-01C	3000073403
3000070480	30000704800138	D684-00-0480-01U	3000072916
3000070480	30000704800138	D684-00-0480-01C	3000073065
3000070480	30000704800138	D684-00-0480-02	3000073275
3000070480	30000704800138	D684-DO-0475	0
3000070480	30000704800138	D684-00-0480-01C	3000073403
3000070480	30000704800139	D684-00-0480-01U	3000072916
3000070480	30000704800139	D684-00-0480-01C	3000073065
3000070480	30000704800139	D684-00-0480-02	3000073275
3000070480	30000704800139	D684-DO-0475	0
3000070480	30000704800139	D684-00-0480-01C	3000073403
3000070482	30000704820026	D684-00-0607	3000073186
3000070482	30000704820033	D684-00-0607	3000073186
3000070482	30000704820039	D684-00-0607	3000073186
3000070482	30000704820042	D684-00-0607	3000073186
3000070482	30000704820043	D684-00-0607	3000073186
3000070487	30000704870373	D684-00-0576-01	3000073278

3000070487	30000704870373	D684-00-0576-01U	3000073282
3000070487	30000704870373	D684-00-0576-01U	3000073283
3000070487	30000704870694	D684-00-0576-01	3000073278
3000070487	30000704870694	D684-00-0576-01U	3000073282
3000070487	30000704870694	D684-00-0576-01U	3000073283
3000070487	30000704870683	D684-00-0576-01	3000073278
3000070487	30000704870683	D684-00-0576-01U	3000073282
3000070487	30000704870683	D684-00-0576-01U	3000073283
3000070487	30000704870693	D684-00-0576-01	3000073278
3000070487	30000704870693	D684-00-0576-01U	3000073282
3000070487	30000704870693	D684-00-0576-01U	3000073283
3000070487	30000704870692	D684-00-0576-01	3000073278
3000070487	30000704870692	D684-00-0576-01U	3000073282
3000070487	30000704870692	D684-00-0576-01U	3000073283
3000070689	30000706890080	D684-00-0576-01U	3000073281
3000070689	30000706890081	D684-00-0576-01U	3000073281
3000070689	30000706890083	D684-00-0576-01U	3000073281
3000070689	30000706890085	D684-00-0576-01U	3000073281
3000070689	30000706890091	D684-00-0576-01U	3000073281
3000070978	30000709780104	D684-00-0559-01	3000073682
3000070978	30000709780104	D684-00-0559-01	3000073816
3000070978	30000709780104	D684-00-0559-01	3000073817
3000070978	30000709780104	D684-00-0559-01	3000074887
3000070978	30000709780124	D684-00-0559-01	3000073682
3000070978	30000709780124	D684-00-0559-01	3000073816
3000070978	30000709780124	D684-00-0559-01	3000073817
3000070978	30000709780124	D684-00-0559-01	3000074887
3000070978	30000709780109	D684-00-0559-01	3000073682
3000070978	30000709780109	D684-00-0559-01	3000073816
3000070978	30000709780109	D684-00-0559-01	3000073817
3000070978	30000709780109	D684-00-0559-01	3000074887
3000070978	30000709780108	D684-00-0559-01	3000073682
3000070978	30000709780108	D684-00-0559-01	3000073816
3000070978	30000709780108	D684-00-0559-01	3000073817
3000070978	30000709780108	D684-00-0559-01	3000074887
3000070978	30000709780116	D684-00-0559-01	3000073682
3000070978	30000709780116	D684-00-0559-01	3000073816
3000070978	30000709780116	D684-00-0559-01	3000073817
3000070978	30000709780116	D684-00-0559-01	3000074887
3000070984	30000709840137	D684-00-0607	3000073680
3000070984	30000709840137	D684-00-0607	3000073681
3000070984	30000709840137	D684-00-0607	3000074885
3000070984	30000709840132	D684-00-0607	3000073680
3000070984	30000709840132	D684-00-0607	3000073681
3000070984	30000709840132	D684-00-0607	3000074885
3000070984	30000709840128	D684-00-0607	3000073680
3000070984	30000709840128	D684-00-0607	3000073681
3000070984	30000709840128	D684-00-0607	3000074885
3000070984	30000709840131	D684-00-0607	3000073680
3000070984	30000709840131	D684-00-0607	3000073681
3000070984	30000709840131	D684-00-0607	3000074885
3000070984	30000709840144	D684-00-0607	3000073680
3000070984	30000709840144	D684-00-0607	3000073681
3000070984	30000709840144	D684-00-0607	3000074885
3000071396	30000713960208	#NV	#NV
3000071396	30000713960208	D684-00-0549-02	3000079504
3000071396	30000713960214	#NV	#NV
3000071396	30000713960214	D684-00-0549-02	3000079504

3000071396	30000713960217	#NV	#NV
3000071396	30000713960217	D684-00-0549-02	3000079504
3000071396	30000713960210	#NV	#NV
3000071396	30000713960210	D684-00-0549-02	3000079504
3000071396	30000713960207	#NV	#NV
3000071396	30000713960207	D684-00-0549-02	3000079504
3000071397	30000713970079	N/A	N/A
3000071397	30000713970076	N/A	N/A
3000071397	30000713970080	N/A	N/A
3000071397	30000713970078	N/A	N/A
3000071397	30000713970081	N/A	N/A
3000071404	30000714040003	D684-00-0292	3000079604
3000071404	30000714040005	D684-00-0292	3000079604
3000071404	30000714040007	D684-00-0292	3000079604
3000071404	30000714040004	D684-00-0292	3000079604
3000071404	30000714040002	D684-00-0292	3000079604
3000071408	30000714080039	N/A	N/A
3000071408	30000714080041	N/A	N/A
3000071408	30000714080032	N/A	N/A
3000071408	30000714080034	N/A	N/A
3000071408	30000714080029	N/A	N/A
3000071841	30000718410349	D684-00-0576-01U	3000074481
3000071841	30000718410349	D684-00-0576-01U	3000074482
3000071841	30000718410349	D684-00-0576-01	3000074484
3000071841	30000718410350	D684-00-0576-01U	3000074481
3000071841	30000718410350	D684-00-0576-01U	3000074482
3000071841	30000718410350	D684-00-0576-01	3000074484
3000071841	30000718410355	D684-00-0576-01U	3000074481
3000071841	30000718410355	D684-00-0576-01U	3000074482
3000071841	30000718410355	D684-00-0576-01	3000074484
3000071841	30000718410367	D684-00-0576-01U	3000074481
3000071841	30000718410367	D684-00-0576-01U	3000074482
3000071841	30000718410367	D684-00-0576-01	3000074484
3000071841	30000718410369	D684-00-0576-01U	3000074481
3000071841	30000718410369	D684-00-0576-01U	3000074482
3000071841	30000718410369	D684-00-0576-01	3000074484
3000071845	30000718450072	D684-00-0568-01U	3000075373
3000071845	30000718450072	D684-00-0568-01U	3000075374
3000071845	30000718450072	D684-00-0568-01	3000075678
3000071845	30000718450071	D684-00-0568-01U	3000075373
3000071845	30000718450071	D684-00-0568-01U	3000075374
3000071845	30000718450071	D684-00-0568-01	3000075678
3000071845	30000718450070	D684-00-0568-01U	3000075373
3000071845	30000718450070	D684-00-0568-01U	3000075374
3000071845	30000718450070	D684-00-0568-01	3000075678
3000071845	30000718450058	D684-00-0568-01U	3000075373
3000071845	30000718450058	D684-00-0568-01U	3000075374
3000071845	30000718450058	D684-00-0568-01	3000075678
3000071845	30000718450067	D684-00-0568-01U	3000075373
3000071845	30000718450067	D684-00-0568-01U	3000075374
3000071845	30000718450067	D684-00-0568-01	3000075678
3000072261	30000722610184	D684-00-0561-01	3000087744
3000072261	30000722610184	D684-00-0561-01	3000090528
3000072261	30000722610185	D684-00-0561-01	3000087744
3000072261	30000722610185	D684-00-0561-01	3000090528
3000072261	30000722610186	D684-00-0561-01	3000087744
3000072261	30000722610186	D684-00-0561-01	3000090528
3000072261	30000722610187	D684-00-0561-01	3000087744

3000072261	30000722610187	D684-00-0561-01	3000090528
3000072261	30000722610188	D684-00-0561-01	3000087744
3000072261	30000722610188	D684-00-0561-01	3000090528
3000072317	30000723170729	D684-00-0470-01U	3000075133
3000072317	30000723170729	D684-00-0470-01	3000075134
3000072317	30000723170729	D684-00-0470-01	3000075676
3000072317	30000723170729	D684-00-0470-01U	3000076154
3000072317	30000723170736	D684-00-0470-01U	3000075133
3000072317	30000723170736	D684-00-0470-01	3000075134
3000072317	30000723170736	D684-00-0470-01	3000075676
3000072317	30000723170736	D684-00-0470-01U	3000076154
3000072317	30000723170738	D684-00-0470-01U	3000075133
3000072317	30000723170738	D684-00-0470-01	3000075134
3000072317	30000723170738	D684-00-0470-01	3000075676
3000072317	30000723170738	D684-00-0470-01U	3000076154
3000072317	30000723170742	D684-00-0470-01U	3000075133
3000072317	30000723170742	D684-00-0470-01	3000075134
3000072317	30000723170742	D684-00-0470-01	3000075676
3000072317	30000723170742	D684-00-0470-01U	3000076154
3000072317	30000723170743	D684-00-0470-01U	3000075133
3000072317	30000723170743	D684-00-0470-01	3000075134
3000072317	30000723170743	D684-00-0470-01	3000075676
3000072317	30000723170743	D684-00-0470-01U	3000076154
3000072320	30000723200073	D684-00-0480-07	3000083859
3000072320	30000723200073	D684-00-0480-01	3000083860
3000072320	30000723200078	D684-00-0480-07	3000083859
3000072320	30000723200078	D684-00-0480-01	3000083860
3000072320	30000723200077	D684-00-0480-07	3000083859
3000072320	30000723200077	D684-00-0480-01	3000083860
3000072320	30000723200074	D684-00-0480-07	3000083859
3000072320	30000723200074	D684-00-0480-01	3000083860
3000072320	30000723200071	D684-00-0480-07	3000083859
3000072320	30000723200071	D684-00-0480-01	3000083860
3000072337	30000723370007	D684-00-0559-01	3000083328
3000072337	30000723370007	D684-00-0559-01	3000083329
3000072337	30000723370007	D684-00-0559-01	3000087746
3000072337	30000723370013	D684-00-0559-01	3000083328
3000072337	30000723370013	D684-00-0559-01	3000083329
3000072337	30000723370013	D684-00-0559-01	3000087746
3000072337	30000723370016	D684-00-0559-01	3000083328
3000072337	30000723370016	D684-00-0559-01	3000083329
3000072337	30000723370016	D684-00-0559-01	3000087746
3000072337	30000723370021	D684-00-0559-01	3000083328
3000072337	30000723370021	D684-00-0559-01	3000083329
3000072337	30000723370021	D684-00-0559-01	3000087746
3000072337	30000723370024	D684-00-0559-01	3000083328
3000072337	30000723370024	D684-00-0559-01	3000083329
3000072337	30000723370024	D684-00-0559-01	3000087746
3000072533	30000725330613	D684-00-0576-01	3000075427
3000072533	30000725330613	D684-00-0576-01U	3000075428
3000072533	30000725330613	D684-00-0576-01U	3000075429
3000072533	30000725330614	D684-00-0576-01	3000075427
3000072533	30000725330614	D684-00-0576-01U	3000075428
3000072533	30000725330614	D684-00-0576-01U	3000075429
3000072533	30000725330619	D684-00-0576-01	3000075427
3000072533	30000725330619	D684-00-0576-01U	3000075428
3000072533	30000725330619	D684-00-0576-01U	3000075429
3000072533	30000725330635	D684-00-0576-01	3000075427

3000072533	30000725330635	D684-00-0576-01U	3000075428
3000072533	30000725330635	D684-00-0576-01U	3000075429
3000072533	30000725330636	D684-00-0576-01	3000075427
3000072533	30000725330636	D684-00-0576-01U	3000075428
3000072533	30000725330636	D684-00-0576-01U	3000075429
3000072534	30000725340407	D684-00-0479-01C	3000075865
3000072534	30000725340407	D684-00-0479-01U	3000075914
3000072534	30000725340407	D684-00-0479-02	3000075916
3000072534	30000725340407	D684-00-0479-01	3000075917
3000072534	30000725340409	D684-00-0479-01C	3000075865
3000072534	30000725340409	D684-00-0479-01U	3000075914
3000072534	30000725340409	D684-00-0479-02	3000075916
3000072534	30000725340409	D684-00-0479-01	3000075917
3000072534	30000725340414	D684-00-0479-01C	3000075865
3000072534	30000725340414	D684-00-0479-01U	3000075914
3000072534	30000725340414	D684-00-0479-02	3000075916
3000072534	30000725340414	D684-00-0479-01	3000075917
3000072534	30000725340398	D684-00-0479-01C	3000075865
3000072534	30000725340398	D684-00-0479-01U	3000075914
3000072534	30000725340398	D684-00-0479-02	3000075916
3000072534	30000725340398	D684-00-0479-01	3000075917
3000072534	30000725340418	D684-00-0479-01C	3000075865
3000072534	30000725340418	D684-00-0479-01U	3000075914
3000072534	30000725340418	D684-00-0479-02	3000075916
3000072534	30000725340418	D684-00-0479-01	3000075917
3000073083	30000730830335	D684-00-0480-01U	3000075681
3000073083	30000730830335	D684-00-0480-01	3000075682
3000073083	30000730830335	D684-00-0480-02	3000075683
3000073083	30000730830335	D684-00-0480-01C	3000075684
3000073083	30000730830344	D684-00-0480-01U	3000075681
3000073083	30000730830344	D684-00-0480-01	3000075682
3000073083	30000730830344	D684-00-0480-02	3000075683
3000073083	30000730830344	D684-00-0480-01C	3000075684
3000073083	30000730830339	D684-00-0480-01U	3000075681
3000073083	30000730830339	D684-00-0480-01	3000075682
3000073083	30000730830339	D684-00-0480-02	3000075683
3000073083	30000730830339	D684-00-0480-01C	3000075684
3000073083	30000730830336	D684-00-0480-01U	3000075681
3000073083	30000730830336	D684-00-0480-01	3000075682
3000073083	30000730830336	D684-00-0480-02	3000075683
3000073083	30000730830336	D684-00-0480-01C	3000075684
3000073083	30000730830342	D684-00-0480-01U	3000075681
3000073083	30000730830342	D684-00-0480-01	3000075682
3000073083	30000730830342	D684-00-0480-02	3000075683
3000073083	30000730830342	D684-00-0480-01C	3000075684
3000073084	30000730840279	D684-00-0469-01	3000075922
3000073084	30000730840279	D684-00-0469-01U	3000076924
3000073084	30000730840279	D684-00-0469-01C	3000077203
3000073084	30000730840279	D684-00-0469-01U	3000077812
3000073084	30000730840279	D684-00-0469-01	3000077813
3000073084	30000730840269	D684-00-0469-01	3000075922
3000073084	30000730840269	D684-00-0469-01U	3000076924
3000073084	30000730840269	D684-00-0469-01C	3000077203
3000073084	30000730840269	D684-00-0469-01U	3000077812
3000073084	30000730840269	D684-00-0469-01	3000077813
3000073084	30000730840271	D684-00-0469-01	3000075922
3000073084	30000730840271	D684-00-0469-01U	3000076924
3000073084	30000730840271	D684-00-0469-01C	3000077203

[illegible]

3000073086	30000730860213	D684-00-0607	3000076710
3000073086	30000730860213	D684-00-0607	3000076711
3000073086	30000730860215	D684-00-0607	3000076710
3000073086	30000730860215	D684-00-0607	3000076711
3000073087	30000730870251	D684-00-0608	3000076242
3000073087	30000730870250	D684-00-0608	3000076242
3000073087	30000730870236	D684-00-0608	3000076242
3000073087	30000730870231	D684-00-0608	3000076242
3000073087	30000730870230	D684-00-0608	3000076242
3000073434	30000734340088	D684-00-0568-01U	3000076243
3000073434	30000734340088	D684-00-0568-01	3000076244
3000073434	30000734340084	D684-00-0568-01U	3000076243
3000073434	30000734340084	D684-00-0568-01	3000076244
3000073434	30000734340085	D684-00-0568-01U	3000076243
3000073434	30000734340085	D684-00-0568-01	3000076244
3000073434	30000734340098	D684-00-0568-01U	3000076243
3000073434	30000734340098	D684-00-0568-01	3000076244
3000073434	30000734340105	D684-00-0568-01U	3000076243
3000073434	30000734340105	D684-00-0568-01	3000076244
3000073438	30000734380514	D684-00-0576-01	3000075677
3000073438	30000734380514	D684-00-0576-01	3000075964
3000073438	30000734380514	D684-00-0576-01U	3000077086
3000073438	30000734380514	D684-00-0576-01	3000077095
3000073438	30000734380499	D684-00-0576-01	3000075677
3000073438	30000734380499	D684-00-0576-01	3000075964
3000073438	30000734380499	D684-00-0576-01U	3000077086
3000073438	30000734380499	D684-00-0576-01	3000077095
3000073438	30000734380492	D684-00-0576-01	3000075677
3000073438	30000734380492	D684-00-0576-01	3000075964
3000073438	30000734380492	D684-00-0576-01U	3000077086
3000073438	30000734380492	D684-00-0576-01	3000077095
3000073438	30000734380494	D684-00-0576-01	3000075677
3000073438	30000734380494	D684-00-0576-01	3000075964
3000073438	30000734380494	D684-00-0576-01U	3000077086
3000073438	30000734380494	D684-00-0576-01	3000077095
3000073438	30000734380501	D684-00-0576-01	3000075677
3000073438	30000734380501	D684-00-0576-01	3000075964
3000073438	30000734380501	D684-00-0576-01U	3000077086
3000073438	30000734380501	D684-00-0576-01	3000077095
3000073441	30000734410236	D684-00-0568-01U	3000077030
3000073441	30000734410236	D684-00-0568-01	3000077031
3000073441	30000734410243	D684-00-0568-01U	3000077030
3000073441	30000734410243	D684-00-0568-01	3000077031
3000073441	30000734410248	D684-00-0568-01U	3000077030
3000073441	30000734410248	D684-00-0568-01	3000077031
3000073441	30000734410252	D684-00-0568-01U	3000077030
3000073441	30000734410252	D684-00-0568-01	3000077031
3000073441	30000734410255	D684-00-0568-01U	3000077030
3000073441	30000734410255	D684-00-0568-01	3000077031
3000073510	30000735100281	D684-00-0296-01	3000075999
3000073510	30000735100281	D684-00-0296-01U	3000076000
3000073510	30000735100281	D684-00-0296-02	3000076001
3000073510	30000735100281	D684-00-0498-01	3000076131
3000073510	30000735100281	D684-00-0498-01	3000076402
3000073510	30000735100268	D684-00-0296-01	3000075999
3000073510	30000735100268	D684-00-0296-01U	3000076000
3000073510	30000735100268	D684-00-0296-02	3000076001
3000073510	30000735100268	D684-00-0498-01	3000076131

[illegible]

3000073522	30000735221024	D684-00-0479-02	3000079931
3000073522	30000735221024	D684-00-0479-01U	3000079932
3000073522	30000735221024	D684-00-0479-01C	3000080105
3000073522	30000735221025	D684-00-0479-07	3000079897
3000073522	30000735221025	D684-00-0479-01	3000079930
3000073522	30000735221025	D684-00-0479-02	3000079931
3000073522	30000735221025	D684-00-0479-01U	3000079932
3000073522	30000735221025	D684-00-0479-01C	3000080105
3000073599	30000735990362	D684-00-0480-01U	3000075131
3000073599	30000735990362	D684-00-0480-02U	3000075132
3000073599	30000735990362	D684-00-0475	3000079438
3000073599	30000735990354	D684-00-0480-01U	3000075131
3000073599	30000735990354	D684-00-0480-02U	3000075132
3000073599	30000735990354	D684-00-0475	3000079438
3000073599	30000735990363	D684-00-0480-01U	3000075131
3000073599	30000735990363	D684-00-0480-02U	3000075132
3000073599	30000735990363	D684-00-0475	3000079438
3000073599	30000735990366	D684-00-0480-01U	3000075131
3000073599	30000735990366	D684-00-0480-02U	3000075132
3000073599	30000735990366	D684-00-0475	3000079438
3000073599	30000735990345	D684-00-0480-01U	3000075131
3000073599	30000735990345	D684-00-0480-02U	3000075132
3000073599	30000735990345	D684-00-0475	3000079438
3000073600	30000736000079	D684-00-0480-01U	3000079466
3000073600	30000736000079	D684-00-0480-02U	3000079467
3000073600	30000736000079	D684-00-0480-01C	3000079569
3000073600	30000736000079	D684-00-0480-01	3000079570
3000073600	30000736000076	D684-00-0480-01U	3000079466
3000073600	30000736000076	D684-00-0480-02U	3000079467
3000073600	30000736000076	D684-00-0480-01C	3000079569
3000073600	30000736000076	D684-00-0480-01	3000079570
3000073600	30000736000080	D684-00-0480-01U	3000079466
3000073600	30000736000080	D684-00-0480-02U	3000079467
3000073600	30000736000080	D684-00-0480-01C	3000079569
3000073600	30000736000080	D684-00-0480-01	3000079570
3000073600	30000736000078	D684-00-0480-01U	3000079466
3000073600	30000736000078	D684-00-0480-02U	3000079467
3000073600	30000736000078	D684-00-0480-01C	3000079569
3000073600	30000736000078	D684-00-0480-01	3000079570
3000073600	30000736000081	D684-00-0480-01U	3000079466
3000073600	30000736000081	D684-00-0480-02U	3000079467
3000073600	30000736000081	D684-00-0480-01C	3000079569
3000073600	30000736000081	D684-00-0480-01	3000079570
3000073602	30000736020039	D684-00-0480-02U	3000079941
3000073602	30000736020039	D684-00-0480-01	3000080053
3000073602	30000736020039	D684-00-0480-01U	3000080054
3000073602	30000736020039	D684-00-0480-01C	3000080103
3000073602	30000736020041	D684-00-0480-02U	3000079941
3000073602	30000736020041	D684-00-0480-01	3000080053
3000073602	30000736020041	D684-00-0480-01U	3000080054
3000073602	30000736020041	D684-00-0480-01C	3000080103
3000073602	30000736020032	D684-00-0480-02U	3000079941
3000073602	30000736020032	D684-00-0480-01	3000080053
3000073602	30000736020032	D684-00-0480-01U	3000080054
3000073602	30000736020032	D684-00-0480-01C	3000080103
3000073602	30000736020034	D684-00-0480-02U	3000079941
3000073602	30000736020034	D684-00-0480-01	3000080053
3000073602	30000736020034	D684-00-0480-01U	3000080054

3000073602	30000736020034	D684-00-0480-01C	3000080103
3000073602	30000736020029	D684-00-0480-02U	3000079941
3000073602	30000736020029	D684-00-0480-01	3000080053
3000073602	30000736020029	D684-00-0480-01U	3000080054
3000073602	30000736020029	D684-00-0480-01C	3000080103
3000073607	30000736070272	D684-00-0560-01	3000087745
3000073607	30000736070272	D684-00-0560-01	3000090529
3000073607	30000736070272	D684-00-0560-01	3000090660
3000073607	30000736070272	D684-00-0560-01	3000092776
3000073607	30000736070268	D684-00-0560-01	3000087745
3000073607	30000736070268	D684-00-0560-01	3000090529
3000073607	30000736070268	D684-00-0560-01	3000090660
3000073607	30000736070268	D684-00-0560-01	3000092776
3000073607	30000736070269	D684-00-0560-01	3000087745
3000073607	30000736070269	D684-00-0560-01	3000090529
3000073607	30000736070269	D684-00-0560-01	3000090660
3000073607	30000736070269	D684-00-0560-01	3000092776
3000073607	30000736070270	D684-00-0560-01	3000087745
3000073607	30000736070270	D684-00-0560-01	3000090529
3000073607	30000736070270	D684-00-0560-01	3000090660
3000073607	30000736070270	D684-00-0560-01	3000092776
3000073607	30000736070267	D684-00-0560-01	3000087745
3000073607	30000736070267	D684-00-0560-01	3000090529
3000073607	30000736070267	D684-00-0560-01	3000090660
3000073607	30000736070267	D684-00-0560-01	3000092776
3000073608	30000736080263	D684-00-0295-05	3000080729
3000073608	30000736080263	D684-00-0295-02U	3000083513
3000073608	30000736080265	D684-00-0295-05	3000080729
3000073608	30000736080265	D684-00-0295-02U	3000083513
3000073608	30000736080269	D684-00-0295-05	3000080729
3000073608	30000736080269	D684-00-0295-02U	3000083513
3000073608	30000736080266	D684-00-0295-05	3000080729
3000073608	30000736080266	D684-00-0295-02U	3000083513
3000073608	30000736080270	D684-00-0295-05	3000080729
3000073608	30000736080270	D684-00-0295-02U	3000083513
3000074101	30000741010820	D684-00-0470-01U	3000076956
3000074101	30000741010820	D684-00-0470-01	3000076957
3000074101	30000741010820	D684-00-0470-01C	3000077204
3000074101	30000741010821	D684-00-0470-01U	3000076956
3000074101	30000741010821	D684-00-0470-01	3000076957
3000074101	30000741010821	D684-00-0470-01C	3000077204
3000074101	30000741010822	D684-00-0470-01U	3000076956
3000074101	30000741010822	D684-00-0470-01	3000076957
3000074101	30000741010822	D684-00-0470-01C	3000077204
3000074101	30000741010827	D684-00-0470-01U	3000076956
3000074101	30000741010827	D684-00-0470-01	3000076957
3000074101	30000741010827	D684-00-0470-01C	3000077204
3000074101	30000741010830	D684-00-0470-01U	3000076956
3000074101	30000741010830	D684-00-0470-01	3000076957
3000074101	30000741010830	D684-00-0470-01C	3000077204
3000074127	30000741270053	D684-00-0568-01	3000077765
3000074127	30000741270053	D684-00-0568-01U	3000077807
3000074127	30000741270053	D684-00-0568-01U	3000077808
3000074127	30000741270058	D684-00-0568-01	3000077765
3000074127	30000741270058	D684-00-0568-01U	3000077807
3000074127	30000741270058	D684-00-0568-01U	3000077808
3000074127	30000741270062	D684-00-0568-01	3000077765
3000074127	30000741270062	D684-00-0568-01U	3000077807

[illegible]

3000074511	30000745110179	D684-00-0561-01	3000076632
3000074511	30000745110180	D684-00-0561-01	3000076632
3000074511	30000745110183	D684-00-0561-01	3000076632
3000074511	30000745110184	D684-00-0561-01	3000076632
3000074609	30000746090179	D684-00-0559-01	3000076462
3000074609	30000746090179	D684-00-0559-01	3000076462
3000074609	30000746090180	D684-00-0559-01	3000076462
3000074609	30000746090180	D684-00-0559-01	3000076462
3000074609	30000746090183	D684-00-0559-01	3000076462
3000074609	30000746090183	D684-00-0559-01	3000076462
3000074609	30000746090184	D684-00-0559-01	3000076462
3000074609	30000746090184	D684-00-0559-01	3000076462
3000074609	30000746090185	D684-00-0559-01	3000076462
3000074609	30000746090185	D684-00-0559-01	3000076462
3000074610	30000746100137	D684-00-0478-01U	3000077341
3000074610	30000746100137	D684-00-0478-01U	3000077814
3000074610	30000746100137	D684-00-0478-01	3000078937
3000074610	30000746100137	D684-00-0478-01	3000079446
3000074610	30000746100137	D684-00-0478-01	3000080540
3000074610	30000746100143	D684-00-0478-01U	3000077341
3000074610	30000746100143	D684-00-0478-01U	3000077814
3000074610	30000746100143	D684-00-0478-01	3000078937
3000074610	30000746100143	D684-00-0478-01	3000079446
3000074610	30000746100143	D684-00-0478-01	3000080540
3000074610	30000746100145	D684-00-0478-01U	3000077341
3000074610	30000746100145	D684-00-0478-01U	3000077814
3000074610	30000746100145	D684-00-0478-01	3000078937
3000074610	30000746100145	D684-00-0478-01	3000079446
3000074610	30000746100145	D684-00-0478-01	3000080540
3000074610	30000746100151	D684-00-0478-01U	3000077341
3000074610	30000746100151	D684-00-0478-01U	3000077814
3000074610	30000746100151	D684-00-0478-01	3000078937
3000074610	30000746100151	D684-00-0478-01	3000079446
3000074610	30000746100151	D684-00-0478-01	3000080540
3000074610	30000746100153	D684-00-0478-01U	3000077341
3000074610	30000746100153	D684-00-0478-01U	3000077814
3000074610	30000746100153	D684-00-0478-01	3000078937
3000074610	30000746100153	D684-00-0478-01	3000079446
3000074610	30000746100153	D684-00-0478-01	3000080540
3000074611	30000746110103	D684-00-0295-05	3000077201
3000074611	30000746110107	D684-00-0295-05	3000077201
3000074611	30000746110110	D684-00-0295-05	3000077201
3000074611	30000746110113	D684-00-0295-05	3000077201
3000074611	30000746110117	D684-00-0295-05	3000077201
3000075054	30000750540026	D684-00-0294-01U	3000077342
3000075054	30000750540026	D684-00-0294-01	3000077343
3000075054	30000750540035	D684-00-0294-01U	3000077342
3000075054	30000750540035	D684-00-0294-01	3000077343
3000075054	30000750540039	D684-00-0294-01U	3000077342
3000075054	30000750540039	D684-00-0294-01	3000077343
3000075054	30000750540045	D684-00-0294-01U	3000077342
3000075054	30000750540045	D684-00-0294-01	3000077343
3000075054	30000750540047	D684-00-0294-01U	3000077342
3000075054	30000750540047	D684-00-0294-01	3000077343
3000075129	30000751290027	D684-00-0480-02	3000075963
3000075129	30000751290027	D684-00-0480-07	3000077205
3000075129	30000751290027	D684-00-0480-01C	3000077273
3000075129	30000751290029	D684-00-0480-02	3000075963

3000075129	30000751290029	D684-00-0480-07	3000077205
3000075129	30000751290029	D684-00-0480-01C	3000077273
3000075129	30000751290030	D684-00-0480-02	3000075963
3000075129	30000751290030	D684-00-0480-07	3000077205
3000075129	30000751290030	D684-00-0480-01C	3000077273
3000075129	30000751290033	D684-00-0480-02	3000075963
3000075129	30000751290033	D684-00-0480-07	3000077205
3000075129	30000751290033	D684-00-0480-01C	3000077273
3000075129	30000751290035	D684-00-0480-02	3000075963
3000075129	30000751290035	D684-00-0480-07	3000077205
3000075129	30000751290035	D684-00-0480-01C	3000077273
3000075135	30000751350160	D684-00-0546-01	3000077671
3000075135	30000751350149	D684-00-0546-01	3000077671
3000075135	30000751350171	D684-00-0546-01	3000077671
3000075135	30000751350157	D684-00-0546-01	3000077671
3000075135	30000751350166	D684-00-0546-01	3000077671
3000075136	30000751360136	D684-00-0480-01C	3000077519
3000075136	30000751360136	D684-00-0480-01	3000077520
3000075136	30000751360136	D684-00-0480-02	3000077521
3000075136	30000751360136	D684-00-0480-01U	3000077522
3000075136	30000751360136	D684-00-0480-01U	3000077523
3000075136	30000751360136	D684-00-0480-01U	3000077524
3000075136	30000751360138	D684-00-0480-01C	3000077519
3000075136	30000751360138	D684-00-0480-01	3000077520
3000075136	30000751360138	D684-00-0480-02	3000077521
3000075136	30000751360138	D684-00-0480-01U	3000077522
3000075136	30000751360138	D684-00-0480-01U	3000077523
3000075136	30000751360138	D684-00-0480-01U	3000077524
3000075136	30000751360142	D684-00-0480-01C	3000077519
3000075136	30000751360142	D684-00-0480-01	3000077520
3000075136	30000751360142	D684-00-0480-02	3000077521
3000075136	30000751360142	D684-00-0480-01U	3000077522
3000075136	30000751360142	D684-00-0480-01U	3000077523
3000075136	30000751360142	D684-00-0480-01U	3000077524
3000075136	30000751360143	D684-00-0480-01C	3000077519
3000075136	30000751360143	D684-00-0480-01	3000077520
3000075136	30000751360143	D684-00-0480-02	3000077521
3000075136	30000751360143	D684-00-0480-01U	3000077522
3000075136	30000751360143	D684-00-0480-01U	3000077523
3000075136	30000751360143	D684-00-0480-01U	3000077524
3000075136	30000751360144	D684-00-0480-01C	3000077519
3000075136	30000751360144	D684-00-0480-01	3000077520
3000075136	30000751360144	D684-00-0480-02	3000077521
3000075136	30000751360144	D684-00-0480-01U	3000077522
3000075136	30000751360144	D684-00-0480-01U	3000077523
3000075136	30000751360144	D684-00-0480-01U	3000077524
3000075137	30000751370029	D684-00-0608	3000077669
3000075137	30000751370030	D684-00-0608	3000077669
3000075137	30000751370035	D684-00-0608	3000077669
3000075137	30000751370038	D684-00-0608	3000077669
3000075137	30000751370045	D684-00-0608	3000077669
3000075437	30000754370430	D684-00-0479-01C	3000077560
3000075437	30000754370430	D684-00-0479-02	3000077809
3000075437	30000754370430	D684-00-0479-01U	3000077810
3000075437	30000754370430	D684-00-0479-01	3000077811
3000075437	30000754370437	D684-00-0479-01C	3000077560
3000075437	30000754370437	D684-00-0479-	

3000075437	30000754370437	D684-00-0479-01	3000077811
3000075437	30000754370431	D684-00-0479-01C	3000077560
3000075437	30000754370431	D684-00-0479-02	3000077809
3000075437	30000754370431	D684-00-0479-01U	3000077810
3000075437	30000754370431	D684-00-0479-01	3000077811
3000075437	30000754370428	D684-00-0479-01C	3000077560
3000075437	30000754370428	D684-00-0479-02	3000077809
3000075437	30000754370428	D684-00-0479-01U	3000077810
3000075437	30000754370428	D684-00-0479-01	3000077811
3000075437	30000754370432	D684-00-0479-01C	3000077560
3000075437	30000754370432	D684-00-0479-02	3000077809
3000075437	30000754370432	D684-00-0479-01U	3000077810
3000075437	30000754370432	D684-00-0479-01	3000077811
3000075438	30000754380065	D684-00-0607	3000077752
3000075438	30000754380065	D684-00-0607	3000077753
3000075438	30000754380065	D684-00-0607	3000079103
3000075438	30000754380075	D684-00-0607	3000077752
3000075438	30000754380075	D684-00-0607	3000077753
3000075438	30000754380075	D684-00-0607	3000079103
3000075438	30000754380056	D684-00-0607	3000077752
3000075438	30000754380056	D684-00-0607	3000077753
3000075438	30000754380056	D684-00-0607	3000079103
3000075438	30000754380062	D684-00-0607	3000077752
3000075438	30000754380062	D684-00-0607	3000077753
3000075438	30000754380062	D684-00-0607	3000079103
3000075438	30000754380070	D684-00-0607	3000077752
3000075438	30000754380070	D684-00-0607	3000077753
3000075438	30000754380070	D684-00-0607	3000079103
3000075439	30000754390841	D684-00-0576-01U	3000078036
3000075439	30000754390841	D684-00-0576-01U	3000078037
3000075439	30000754390841	D684-00-0576-01U	3000078038
3000075439	30000754390841	D684-00-0576-01	3000078049
3000075439	30000754390841	D684-DO-0575	0
3000075439	30000754390841	D684-DC-0575	0
3000075439	30000754390833	D684-00-0576-01U	3000078036
3000075439	30000754390833	D684-00-0576-01U	3000078037
3000075439	30000754390833	D684-00-0576-01U	3000078038
3000075439	30000754390833	D684-00-0576-01	3000078049
3000075439	30000754390833	D684-DO-0575	0
3000075439	30000754390833	D684-DC-0575	0
3000075439	30000754390824	D684-00-0576-01U	3000078036
3000075439	30000754390824	D684-00-0576-01U	3000078037
3000075439	30000754390824	D684-00-0576-01U	3000078038
3000075439	30000754390824	D684-00-0576-01	3000078049
3000075439	30000754390824	D684-DO-0575	0
3000075439	30000754390824	D684-DC-0575	0
3000075439	30000754390842	D684-00-0576-01U	3000078036
3000075439	30000754390842	D684-00-0576-01U	3000078037
3000075439	30000754390842	D684-00-0576-01U	3000078038
3000075439	30000754390842	D684-00-0576-01	3000078049
3000075439	30000754390842	D684-DO-0575	0
3000075439	30000754390842	D684-DC-0575	0
3000075439	30000754390845	D684-00-0576-01U	3000078036
3000075439	30000754390845	D684-00-0576-01U	3000078037
3000075439	30000754390845	D684-00-0576-01U	3000078038
3000075439	30000754390845	D684-00-0576-01	3000078049
3000075439	30000754390845	D684-DO-0575	0
3000075439	30000754390845	D684-DC-0575	0

[illegible]

[illegible]

[illegible]

[illegible]

3000076277	30000762770005	D684-00-0546-01	3000080310
3000076277	30000762770005	D684-00-0546-01	3000083334
3000076277	30000762770005	D684-00-0546-01	3000083335
3000076277	30000762770009	D684-00-0546-01	3000078857
3000076277	30000762770009	D684-00-0546-01	3000080306
3000076277	30000762770009	D684-00-0546-01	3000080310
3000076277	30000762770009	D684-00-0546-01	3000083334
3000076277	30000762770009	D684-00-0546-01	3000083335
3000076277	30000762770010	D684-00-0546-01	3000078857
3000076277	30000762770010	D684-00-0546-01	3000080306
3000076277	30000762770010	D684-00-0546-01	3000080310
3000076277	30000762770010	D684-00-0546-01	3000083334
3000076277	30000762770010	D684-00-0546-01	3000083335
3000076277	30000762770013	D684-00-0546-01	3000078857
3000076277	30000762770013	D684-00-0546-01	3000080306
3000076277	30000762770013	D684-00-0546-01	3000080310
3000076277	30000762770013	D684-00-0546-01	3000083334
3000076277	30000762770013	D684-00-0546-01	3000083335
3000076277	30000762770017	D684-00-0546-01	3000078857
3000076277	30000762770017	D684-00-0546-01	3000080306
3000076277	30000762770017	D684-00-0546-01	3000080310
3000076277	30000762770017	D684-00-0546-01	3000083334
3000076277	30000762770017	D684-00-0546-01	3000083335
3000076404	30000764040363	D684-00-0607	3000079104
3000076404	30000764040363	D684-00-0607	3000080873
3000076404	30000764040366	D684-00-0607	3000079104
3000076404	30000764040366	D684-00-0607	3000080873
3000076404	30000764040370	D684-00-0607	3000079104
3000076404	30000764040370	D684-00-0607	3000080873
3000076404	30000764040372	D684-00-0607	3000079104
3000076404	30000764040372	D684-00-0607	3000080873
3000076404	30000764040375	D684-00-0607	3000079104
3000076404	30000764040375	D684-00-0607	3000080873
3000076715	30000767150026	D684-00-0498-01	3000078920
3000076715	30000767150026	D684-00-0296-01	3000078921
3000076715	30000767150026	D684-00-0296-01U	3000078922
3000076715	30000767150026	D684-00-0296-02U	3000078923
3000076715	30000767150028	D684-00-0498-01	3000078920
3000076715	30000767150028	D684-00-0296-01	3000078921
3000076715	30000767150028	D684-00-0296-01U	3000078922
3000076715	30000767150028	D684-00-0296-02U	3000078923
3000076715	30000767150029	D684-00-0498-01	3000078920
3000076715	30000767150029	D684-00-0296-01	3000078921
3000076715	30000767150029	D684-00-0296-01U	3000078922
3000076715	30000767150029	D684-00-0296-02U	3000078923
3000076715	30000767150031	D684-00-0498-01	3000078920
3000076715	30000767150031	D684-00-0296-01	3000078921
3000076715	30000767150031	D684-00-0296-01U	3000078922
3000076715	30000767150031	D684-00-0296-02U	3000078923
3000076715	30000767150034	D684-00-0498-01	3000078920
3000076715	30000767150034	D684-00-0296-01	3000078921
3000076715	30000767150034	D684-00-0296-01U	3000078922
3000076715	30000767150034	D684-00-0296-02U	3000078923
3000076716	30000767160202	D684-00-0608	3000079105
3000076716	30000767160202	D684-00-0608	3000079318
3000076716	30000767160202	D684-00-0608	3000080309
3000076716	30000767160204	D684-00-0608	3000079105
3000076716	30000767160204	D684-00-0608	3000079318

[illegible]

3000077103	30000771030457	D684-00-0470-01	3000079650
3000077103	30000771030457	D684-00-0470-01U	3000080226
3000077104	30000771040204	D684-00-0480-07	3000079828
3000077104	30000771040204	D684-00-0480-02	3000079829
3000077104	30000771040207	D684-00-0480-07	3000079828
3000077104	30000771040207	D684-00-0480-02	3000079829
3000077104	30000771040213	D684-00-0480-07	3000079828
3000077104	30000771040213	D684-00-0480-02	3000079829
3000077104	30000771040215	D684-00-0480-07	3000079828
3000077104	30000771040215	D684-00-0480-02	3000079829
3000077104	30000771040223	D684-00-0480-07	3000079828
3000077104	30000771040223	D684-00-0480-02	3000079829
3000077561	30000775610063	D684-00-0294-01U	3000080542
3000077561	30000775610063	D684-00-0294-01U	3000081080
3000077561	30000775610063	D684-00-0294-01	3000084725
3000077561	30000775610065	D684-00-0294-01U	3000080542
3000077561	30000775610065	D684-00-0294-01U	3000081080
3000077561	30000775610065	D684-00-0294-01	3000084725
3000077561	30000775610070	D684-00-0294-01U	3000080542
3000077561	30000775610070	D684-00-0294-01U	3000081080
3000077561	30000775610070	D684-00-0294-01	3000084725
3000077561	30000775610072	D684-00-0294-01U	3000080542
3000077561	30000775610072	D684-00-0294-01U	3000081080
3000077561	30000775610072	D684-00-0294-01	3000084725
3000077561	30000775610074	D684-00-0294-01U	3000080542
3000077561	30000775610074	D684-00-0294-01U	3000081080
3000077561	30000775610074	D684-00-0294-01	3000084725
3000077562	30000775620270	D684-00-0295-05	3000080419
3000077562	30000775620270	D684-00-0295-01U	3000080541
3000077562	30000775620269	D684-00-0295-05	3000080419
3000077562	30000775620269	D684-00-0295-01U	3000080541
3000077562	30000775620001	D684-00-0295-05	3000080419
3000077562	30000775620001	D684-00-0295-01U	3000080541
3000077562	30000775620002	D684-00-0295-05	3000080419
3000077562	30000775620002	D684-00-0295-01U	3000080541
3000077562	30000775620003	D684-00-0295-05	3000080419
3000077562	30000775620003	D684-00-0295-01U	3000080541
3000077563	30000775630515	D684-00-0479-02	3000080102
3000077563	30000775630515	D684-00-0479-01C	3000080106
3000077563	30000775630515	D684-00-0479-01C	3000080444
3000077563	30000775630515	D684-00-0479-07	3000080464
3000077563	30000775630515	D684-00-0479-01U	3000080539
3000077563	30000775630515	D684-00-0479-01	3000080673
3000077563	30000775630515	D684-00-0479-02	3000080731
3000077563	30000775630515	D684-00-0479-01C	3000080732
3000077563	30000775630499	D684-00-0479-02	3000080102
3000077563	30000775630499	D684-00-0479-01C	3000080106
3000077563	30000775630499	D684-00-0479-01C	3000080444
3000077563	30000775630499	D684-00-0479-07	3000080464
3000077563	30000775630499	D684-00-0479-01U	3000080539
3000077563	30000775630499	D684-00-0479-01	3000080673
3000077563	30000775630499	D684-00-0479-02	3000080731
3000077563	30000775630499	D684-00-0479-01C	3000080732
3000077563	30000775630497	D684-00-0479-02	3000080102
3000077563	30000775630497	D684-00-0479-01C	3000080106
3000077563	30000775630497	D684-00-0479-01C	3000080444
3000077563	30000775630497	D684-00-0479-07	3000080464
3000077563	30000775630497	D684-00-0479-01U	3000080539

3000077563	30000775630497	D684-00-0479-01	3000080673
3000077563	30000775630497	D684-00-0479-02	3000080731
3000077563	30000775630497	D684-00-0479-01C	3000080732
3000077563	30000775630513	D684-00-0479-02	3000080102
3000077563	30000775630513	D684-00-0479-01C	3000080106
3000077563	30000775630513	D684-00-0479-01C	3000080444
3000077563	30000775630513	D684-00-0479-07	3000080464
3000077563	30000775630513	D684-00-0479-01U	3000080539
3000077563	30000775630513	D684-00-0479-01	3000080673
3000077563	30000775630513	D684-00-0479-02	3000080731
3000077563	30000775630513	D684-00-0479-01C	3000080732
3000077563	30000775630498	D684-00-0479-02	3000080102
3000077563	30000775630498	D684-00-0479-01C	3000080106
3000077563	30000775630498	D684-00-0479-01C	3000080444
3000077563	30000775630498	D684-00-0479-07	3000080464
3000077563	30000775630498	D684-00-0479-01U	3000080539
3000077563	30000775630498	D684-00-0479-01	3000080673
3000077563	30000775630498	D684-00-0479-02	3000080731
3000077563	30000775630498	D684-00-0479-01C	3000080732
3000077564	30000775640467	D684-00-0568-01	3000079679
3000077564	30000775640467	D684-00-0568-01U	3000079680
3000077564	30000775640483	D684-00-0568-01	3000079679
3000077564	30000775640483	D684-00-0568-01U	3000079680
3000077564	30000775640476	D684-00-0568-01	3000079679
3000077564	30000775640476	D684-00-0568-01U	3000079680
3000077564	30000775640471	D684-00-0568-01	3000079679
3000077564	30000775640471	D684-00-0568-01U	3000079680
3000077564	30000775640466	D684-00-0568-01	3000079679
3000077564	30000775640466	D684-00-0568-01U	3000079680
3000077655	30000776550258	D684-00-0545-01	3000083333
3000077655	30000776550258	D684-00-0545-01	3000084968
3000077655	30000776550258	D684-00-0545-01	3000087492
3000077655	30000776550258	D684-00-0545-01	3000095469
3000077655	30000776550258	D684-00-0545-01	3000096180
3000077655	30000776550257	D684-00-0545-01	3000083333
3000077655	30000776550257	D684-00-0545-01	3000084968
3000077655	30000776550257	D684-00-0545-01	3000087492
3000077655	30000776550257	D684-00-0545-01	3000095469
3000077655	30000776550257	D684-00-0545-01	3000096180
3000077655	30000776550255	D684-00-0545-01	3000083333
3000077655	30000776550255	D684-00-0545-01	3000084968
3000077655	30000776550255	D684-00-0545-01	3000087492
3000077655	30000776550255	D684-00-0545-01	3000095469
3000077655	30000776550255	D684-00-0545-01	3000096180
3000077655	30000776550260	D684-00-0545-01	3000083333
3000077655	30000776550260	D684-00-0545-01	3000084968
3000077655	30000776550260	D684-00-0545-01	3000087492
3000077655	30000776550260	D684-00-0545-01	3000095469
3000077655	30000776550260	D684-00-0545-01	3000096180
3000077655	30000776550259	D684-00-0545-01	3000083333
3000077655	30000776550259	D684-00-0545-01	3000084968
3000077655	30000776550259	D684-00-0545-01	3000087492
3000077655	30000776550259	D684-00-0545-01	3000095469
3000077655	30000776550259	D684-00-0545-01	3000096180
3000077656	30000776560089	D684-00-0480-02	3000078934
3000077656	30000776560089	D684-00-0480-02	3000080098
3000077656	30000776560089	D684-00-0480-01C	3000080104
3000077656	30000776560089	D684-00-0480-01C	3000080443

[illegible]

[illegible]

3000078605	30000786050106	D684-00-0480-01	3000081888
3000078605	30000786050106	D684-00-0480-01	3000082186
3000078605	30000786050106	D684-00-0480-01	3000082885
3000078605	30000786050106	D684-00-0480-01U	3000083094
3000078606	30000786060384	D684-00-0576-01U	3000081444
3000078606	30000786060384	D684-00-0576-01U	3000081445
3000078606	30000786060384	D684-00-0576-01	3000081446
3000078606	30000786060379	D684-00-0576-01U	3000081444
3000078606	30000786060379	D684-00-0576-01U	3000081445
3000078606	30000786060379	D684-00-0576-01	3000081446
3000078606	30000786060380	D684-00-0576-01U	3000081444
3000078606	30000786060380	D684-00-0576-01U	3000081445
3000078606	30000786060380	D684-00-0576-01	3000081446
3000078606	30000786060389	D684-00-0576-01U	3000081444
3000078606	30000786060389	D684-00-0576-01U	3000081445
3000078606	30000786060389	D684-00-0576-01	3000081446
3000078606	30000786060400	D684-00-0576-01U	3000081444
3000078606	30000786060400	D684-00-0576-01U	3000081445
3000078606	30000786060400	D684-00-0576-01	3000081446
3000078607	30000786070058	D684-00-0478-01U	3000083060
3000078607	30000786070058	D684-00-0478-01	3000084470
3000078607	30000786070058	D684-00-0478-01U	3000084471
3000078607	30000786070058	D684-00-0478-01U	3000087477
3000078607	30000786070069	D684-00-0478-01U	3000083060
3000078607	30000786070069	D684-00-0478-01	3000084470
3000078607	30000786070069	D684-00-0478-01U	3000084471
3000078607	30000786070069	D684-00-0478-01U	3000087477
3000078607	30000786070070	D684-00-0478-01U	3000083060
3000078607	30000786070070	D684-00-0478-01	3000084470
3000078607	30000786070070	D684-00-0478-01U	3000084471
3000078607	30000786070070	D684-00-0478-01U	3000087477
3000078607	30000786070052	D684-00-0478-01U	3000083060
3000078607	30000786070052	D684-00-0478-01	3000084470
3000078607	30000786070052	D684-00-0478-01U	3000084471
3000078607	30000786070052	D684-00-0478-01U	3000087477
3000078607	30000786070060	D684-00-0478-01U	3000083060
3000078607	30000786070060	D684-00-0478-01	3000084470
3000078607	30000786070060	D684-00-0478-01U	3000084471
3000078607	30000786070060	D684-00-0478-01U	3000087477
3000078789	30000787890060	D684-00-0296-01	3000080137
3000078789	30000787890061	D684-00-0296-01	3000080137
3000078789	30000787890064	D684-00-0296-01	3000080137
3000078789	30000787890067	D684-00-0296-01	3000080137
3000078789	30000787890069	D684-00-0296-01	3000080137
3000078916	30000789160079	D684-00-0480-07	3000082126
3000078916	30000789160079	D684-00-0480-02	3000082187
3000078916	30000789160088	D684-00-0480-07	3000082126
3000078916	30000789160088	D684-00-0480-02	3000082187
3000078916	30000789160082	D684-00-0480-07	3000082126
3000078916	30000789160082	D684-00-0480-02	3000082187
3000078916	30000789160083	D684-00-0480-07	3000082126
3000078916	30000789160083	D684-00-0480-02	3000082187
3000078916	30000789160098	D684-00-0480-07	3000082126
3000078916	30000789160098	D684-00-0480-02	3000082187
3000078918	30000789180084	D684-00-0296-02	3000081621
3000078918	30000789180084	D684-00-0296-01	3000081624
3000078918	30000789180084	D684-00-0498-01	3000081777
3000078918	30000789180084	D684-00-0296-01	3000082197

3000078918	30000789180080	D684-00-0296-02	3000081621
3000078918	30000789180080	D684-00-0296-01	3000081624
3000078918	30000789180080	D684-00-0498-01	3000081777
3000078918	30000789180080	D684-00-0296-01	3000082197
3000078918	30000789180085	D684-00-0296-02	3000081621
3000078918	30000789180085	D684-00-0296-01	3000081624
3000078918	30000789180085	D684-00-0498-01	3000081777
3000078918	30000789180085	D684-00-0296-01	3000082197
3000078918	30000789180087	D684-00-0296-02	3000081621
3000078918	30000789180087	D684-00-0296-01	3000081624
3000078918	30000789180087	D684-00-0498-01	3000081777
3000078918	30000789180087	D684-00-0296-01	3000082197
3000078918	30000789180086	D684-00-0296-02	3000081621
3000078918	30000789180086	D684-00-0296-01	3000081624
3000078918	30000789180086	D684-00-0498-01	3000081777
3000078918	30000789180086	D684-00-0296-01	3000082197
3000078924	30000789240001	D684-00-0607	3000081671
3000078924	30000789240002	D684-00-0607	3000081671
3000078924	30000789240013	D684-00-0607	3000081671
3000078924	30000789240003	D684-00-0607	3000081671
3000078924	30000789240015	D684-00-0607	3000081671
3000078926	30000789260038	D684-00-0480-01C	3000082659
3000078926	30000789260038	D684-00-0480-01U	3000083092
3000078926	30000789260038	D684-00-0480-01U	3000083093
3000078926	30000789260041	D684-00-0480-01C	3000082659
3000078926	30000789260041	D684-00-0480-01U	3000083092
3000078926	30000789260041	D684-00-0480-01U	3000083093
3000078926	30000789260032	D684-00-0480-01C	3000082659
3000078926	30000789260032	D684-00-0480-01U	3000083092
3000078926	30000789260032	D684-00-0480-01U	3000083093
3000078926	30000789260036	D684-00-0480-01C	3000082659
3000078926	30000789260036	D684-00-0480-01U	3000083092
3000078926	30000789260036	D684-00-0480-01U	3000083093
3000078926	30000789260035	D684-00-0480-01C	3000082659
3000078926	30000789260035	D684-00-0480-01U	3000083092
3000078926	30000789260035	D684-00-0480-01U	3000083093
3000078929	30000789290741	D684-00-0470-01	3000081778
3000078929	30000789290741	D684-00-0470-01U	3000081779
3000078929	30000789290741	D684-00-0470-01U	3000083059
3000078929	30000789290736	D684-00-0470-01	3000081778
3000078929	30000789290736	D684-00-0470-01U	3000081779
3000078929	30000789290736	D684-00-0470-01U	3000083059
3000078929	30000789290724	D684-00-0470-01	3000081778
3000078929	30000789290724	D684-00-0470-01U	3000081779
3000078929	30000789290724	D684-00-0470-01U	3000083059
3000078929	30000789290728	D684-00-0470-01	3000081778
3000078929	30000789290728	D684-00-0470-01U	3000081779
3000078929	30000789290728	D684-00-0470-01U	3000083059
3000078929	30000789290729	D684-00-0470-01	3000081778
3000078929	30000789290729	D684-00-0470-01U	3000081779
3000078929	30000789290729	D684-00-0470-01U	3000083059
3000079438	30000709830362	D684-00-0480-07	3000080463
3000079438	30000709830362	D684-00-0480-01U	3000080538
3000079438	30000709830354	D684-00-0480-07	3000080463
3000079438	30000709830354	D684-00-0480-01U	3000080538
3000079438	30000709830363	D684-00-0480-07	3000080463
3000079438	30000709830363	D684-00-0480-01U	3000080538
3000079438	30000709830366	D684-00-0480-07	3000080463

3000079438	30000709830366	D684-00-0480-01U	3000080538
3000079438	30000709830345	D684-00-0480-07	3000080463
3000079438	30000709830345	D684-00-0480-01U	3000080538
3000079458	30000794580054	D684-00-0479-08	3000082252
3000079458	30000794580054	D684-00-0479-01	3000082253
3000079458	30000794580054	D684-00-0479-01U	3000082534
3000079458	30000794580054	D684-00-0479-07	3000082658
3000079458	30000794580054	D684-00-0479-01C	3000082660
3000079458	30000794580061	D684-00-0479-08	3000082252
3000079458	30000794580061	D684-00-0479-01	3000082253
3000079458	30000794580061	D684-00-0479-01U	3000082534
3000079458	30000794580061	D684-00-0479-07	3000082658
3000079458	30000794580061	D684-00-0479-01C	3000082660
3000079458	30000794580052	D684-00-0479-08	3000082252
3000079458	30000794580052	D684-00-0479-01	3000082253
3000079458	30000794580052	D684-00-0479-01U	3000082534
3000079458	30000794580052	D684-00-0479-07	3000082658
3000079458	30000794580052	D684-00-0479-01C	3000082660
3000079458	30000794580059	D684-00-0479-08	3000082252
3000079458	30000794580059	D684-00-0479-01	3000082253
3000079458	30000794580059	D684-00-0479-01U	3000082534
3000079458	30000794580059	D684-00-0479-07	3000082658
3000079458	30000794580059	D684-00-0479-01C	3000082660
3000079458	30000794580055	D684-00-0479-08	3000082252
3000079458	30000794580055	D684-00-0479-01	3000082253
3000079458	30000794580055	D684-00-0479-01U	3000082534
3000079458	30000794580055	D684-00-0479-07	3000082658
3000079458	30000794580055	D684-00-0479-01C	3000082660
3000079459	30000794590523	D684-00-0568-01U	3000082471
3000079459	30000794590523	D684-00-0568-01	3000082472
3000079459	30000794590522	D684-00-0568-01U	3000082471
3000079459	30000794590522	D684-00-0568-01	3000082472
3000079459	30000794590532	D684-00-0568-01U	3000082471
3000079459	30000794590532	D684-00-0568-01	3000082472
3000079459	30000794590528	D684-00-0568-01U	3000082471
3000079459	30000794590528	D684-00-0568-01	3000082472
3000079459	30000794590524	D684-00-0568-01U	3000082471
3000079459	30000794590524	D684-00-0568-01	3000082472
3000079460	30000794600438	D684-00-0576-01	3000082347
3000079460	30000794600438	D684-00-0576-01U	3000082348
3000079460	30000794600438	D684-00-0576-01U	3000082349
3000079460	30000794600440	D684-00-0576-01	3000082347
3000079460	30000794600440	D684-00-0576-01U	3000082348
3000079460	30000794600440	D684-00-0576-01U	3000082349
3000079460	30000794600442	D684-00-0576-01	3000082347
3000079460	30000794600442	D684-00-0576-01U	3000082348
3000079460	30000794600442	D684-00-0576-01U	3000082349
3000079460	30000794600456	D684-00-0576-01	3000082347
3000079460	30000794600456	D684-00-0576-01U	3000082348
3000079460	30000794600456	D684-00-0576-01U	3000082349
3000079460	30000794600434	D684-00-0576-01	3000082347
3000079460	30000794600434	D684-00-0576-01U	3000082348
3000079460	30000794600434	D684-00-0576-01U	3000082349
3000079504	30000795040208	D684-00-0560-01	3000092778
3000079504	30000795040208	D684-00-0560-01	3000095472
3000079504	30000795040208	D684-00-0560-01	3000096375
3000079504	30000795040214	D684-00-0560-01	3000092778
3000079504	30000795040214	D684-00-0560-01	3000095472

3000080864	30000808640140	D684-00-0576-01U	3000084155
3000080864	30000808640142	D684-00-0576-01	3000084154
3000080864	30000808640142	D684-00-0576-01U	3000084155
3000080865	30000808650485	D684-00-0576-01U	3000083509
3000080865	30000808650485	D684-00-0576-01	3000083510
3000080865	30000808650485	D684-00-0576-01U	3000084376
3000080865	30000808650485	D684-00-0576-01U	3000084377
3000080865	30000808650492	D684-00-0576-01U	3000083509
3000080865	30000808650492	D684-00-0576-01	3000083510
3000080865	30000808650492	D684-00-0576-01U	3000084376
3000080865	30000808650492	D684-00-0576-01U	3000084377
3000080865	30000808650480	D684-00-0576-01U	3000083509
3000080865	30000808650480	D684-00-0576-01	3000083510
3000080865	30000808650480	D684-00-0576-01U	3000084376
3000080865	30000808650480	D684-00-0576-01U	3000084377
3000080865	30000808650483	D684-00-0576-01U	3000083509
3000080865	30000808650483	D684-00-0576-01	3000083510
3000080865	30000808650483	D684-00-0576-01U	3000084376
3000080865	30000808650483	D684-00-0576-01U	3000084377
3000080865	30000808650500	D684-00-0576-01U	3000083509
3000080865	30000808650500	D684-00-0576-01	3000083510
3000080865	30000808650500	D684-00-0576-01U	3000084376
3000080865	30000808650500	D684-00-0576-01U	3000084377
3000080870	30000808700087	D684-00-0480-01C	3000083600
3000080870	30000808700087	D684-00-0480-01	3000083861
3000080870	30000808700081	D684-00-0480-01C	3000083600
3000080870	30000808700081	D684-00-0480-01	3000083861
3000080870	30000808700091	D684-00-0480-01C	3000083600
3000080870	30000808700091	D684-00-0480-01	3000083861
3000080870	30000808700099	D684-00-0480-01C	3000083600
3000080870	30000808700099	D684-00-0480-01	3000083861
3000080870	30000808700094	D684-00-0480-01C	3000083600
3000080870	30000808700094	D684-00-0480-01	3000083861
3000081448	30000814480182	D684-00-0607	3000084267
3000081448	30000814480182	D684-00-0607	3000084966
3000081448	30000814480197	D684-00-0607	3000084267
3000081448	30000814480197	D684-00-0607	3000084966
3000081448	30000814480192	D684-00-0607	3000084267
3000081448	30000814480192	D684-00-0607	3000084966
3000081448	30000814480194	D684-00-0607	3000084267
3000081448	30000814480194	D684-00-0607	3000084966
3000081448	30000814480188	D684-00-0607	3000084267
3000081448	30000814480188	D684-00-0607	3000084966
3000081527	30000815270051	D684-00-0608	3000084268
3000081527	30000815270051	D684-00-0608	3000084967
3000081527	30000815270063	D684-00-0608	3000084268
3000081527	30000815270063	D684-00-0608	3000084967
3000081527	30000815270065	D684-00-0608	3000084268
3000081527	30000815270065	D684-00-0608	3000084967
3000081527	30000815270067	D684-00-0608	3000084268
3000081527	30000815270067	D684-00-0608	3000084967
3000081527	30000815270075	D684-00-0608	3000084268
3000081527	30000815270075	D684-00-0608	3000084967
3000081529	30000815290449	D684-00-0470-01	3000084071
3000081529	30000815290449	D684-00-0470-01U	3000084153
3000081529	30000815290449	D684-00-0470-01C	3000084689
3000081529	30000815290449	D684-00-0470-01	3000084788
3000081529	30000815290449	D684-00-0470-01	3000084789

[illegible]

[illegible]

3000081742	30000817420015	D684-00-0480-01	3000084565
3000081742	30000817420015	D684-00-0480-07	3000084652
3000081742	30000817420015	D684-DO-0475	0
3000081742	30000817420019	D684-00-0480-01	3000084565
3000081742	30000817420019	D684-00-0480-07	3000084652
3000081742	30000817420019	D684-DO-0475	0
3000081742	30000817420021	D684-00-0480-01	3000084565
3000081742	30000817420021	D684-00-0480-07	3000084652
3000081742	30000817420021	D684-DO-0475	0
3000081744	30000817440476	D684-00-0568-01U	3000084469
3000081744	30000817440476	D684-00-0568-01U	3000084478
3000081744	30000817440476	D684-00-0568-01U	3000085192
3000081744	30000817440476	D684-00-0568-01	3000085231
3000081744	30000817440477	D684-00-0568-01U	3000084469
3000081744	30000817440477	D684-00-0568-01U	3000084478
3000081744	30000817440477	D684-00-0568-01U	3000085192
3000081744	30000817440477	D684-00-0568-01	3000085231
3000081744	30000817440482	D684-00-0568-01U	3000084469
3000081744	30000817440482	D684-00-0568-01U	3000084478
3000081744	30000817440482	D684-00-0568-01U	3000085192
3000081744	30000817440482	D684-00-0568-01	3000085231
3000081744	30000817440490	D684-00-0568-01U	3000084469
3000081744	30000817440490	D684-00-0568-01U	3000084478
3000081744	30000817440490	D684-00-0568-01U	3000085192
3000081744	30000817440490	D684-00-0568-01	3000085231
3000081744	30000817440491	D684-00-0568-01U	3000084469
3000081744	30000817440491	D684-00-0568-01U	3000084478
3000081744	30000817440491	D684-00-0568-01U	3000085192
3000081744	30000817440491	D684-00-0568-01	3000085231
3000081750	30000817500570	D684-00-0576-01U	3000084915
3000081750	30000817500570	D684-00-0576-01U	3000084916
3000081750	30000817500570	D684-00-0576-01	3000084917
3000081750	30000817500572	D684-00-0576-01U	3000084915
3000081750	30000817500572	D684-00-0576-01U	3000084916
3000081750	30000817500572	D684-00-0576-01	3000084917
3000081750	30000817500583	D684-00-0576-01U	3000084915
3000081750	30000817500583	D684-00-0576-01U	3000084916
3000081750	30000817500583	D684-00-0576-01	3000084917
3000081750	30000817500584	D684-00-0576-01U	3000084915
3000081750	30000817500584	D684-00-0576-01U	3000084916
3000081750	30000817500584	D684-00-0576-01	3000084917
3000081750	30000817500587	D684-00-0576-01U	3000084915
3000081750	30000817500587	D684-00-0576-01U	3000084916
3000081750	30000817500587	D684-00-0576-01	3000084917
3000081751	30000817510184	D684-00-0296-01U	3000084730
3000081751	30000817510184	D684-00-0296-02U	3000084731
3000081751	30000817510184	D684-00-0296-01	3000084785
3000081751	30000817510184	D684-00-0296-02	3000084786
3000081751	30000817510184	D684-00-0498-01	3000085937
3000081751	30000817510184	D684-00-0296-01U	3000087112
3000081751	30000817510185	D684-00-0296-01U	3000084730
3000081751	30000817510185	D684-00-0296-02U	3000084731
3000081751	30000817510185	D684-00-0296-01	3000084785
3000081751	30000817510185	D684-00-0296-02	3000084786
3000081751	30000817510185	D684-00-0498-01	3000085937
3000081751	30000817510185	D684-00-0296-01U	3000087112
3000081751	30000817510187	D684-00-0296-01U	3000084730
3000081751	30000817510187	D684-00-0296-02U	3000084731

3000081751	30000817510187	D684-00-0296-01	3000084785
3000081751	30000817510187	D684-00-0296-02	3000084786
3000081751	30000817510187	D684-00-0498-01	3000085937
3000081751	30000817510187	D684-00-0296-01U	3000087112
3000081751	30000817510198	D684-00-0296-01U	3000084730
3000081751	30000817510198	D684-00-0296-02U	3000084731
3000081751	30000817510198	D684-00-0296-01	3000084785
3000081751	30000817510198	D684-00-0296-02	3000084786
3000081751	30000817510198	D684-00-0498-01	3000085937
3000081751	30000817510198	D684-00-0296-01U	3000087112
3000081751	30000817510199	D684-00-0296-01U	3000084730
3000081751	30000817510199	D684-00-0296-02U	3000084731
3000081751	30000817510199	D684-00-0296-01	3000084785
3000081751	30000817510199	D684-00-0296-02	3000084786
3000081751	30000817510199	D684-00-0498-01	3000085937
3000081751	30000817510199	D684-00-0296-01U	3000087112
3000081752	30000817520111	D684-00-0480-01U	3000084692
3000081752	30000817520111	D684-00-0480-01	3000084851
3000081752	30000817520111	D684-00-0480-01	3000085193
3000081752	30000817520112	D684-00-0480-01U	3000084692
3000081752	30000817520112	D684-00-0480-01	3000084851
3000081752	30000817520112	D684-00-0480-01	3000085193
3000081752	30000817520113	D684-00-0480-01U	3000084692
3000081752	30000817520113	D684-00-0480-01	3000084851
3000081752	30000817520113	D684-00-0480-01	3000085193
3000081752	30000817520114	D684-00-0480-01U	3000084692
3000081752	30000817520114	D684-00-0480-01	3000084851
3000081752	30000817520114	D684-00-0480-01	3000085193
3000081752	30000817520115	D684-00-0480-01U	3000084692
3000081752	30000817520115	D684-00-0480-01	3000084851
3000081752	30000817520115	D684-00-0480-01	3000085193
3000082758	30000827580517	D684-00-0608	3000084960
3000082758	30000827580517	D684-00-0608	3000084961
3000082758	30000827580518	D684-00-0608	3000084960
3000082758	30000827580518	D684-00-0608	3000084961
3000082758	30000827580519	D684-00-0608	3000084960
3000082758	30000827580519	D684-00-0608	3000084961
3000082758	30000827580520	D684-00-0608	3000084960
3000082758	30000827580520	D684-00-0608	3000084961
3000082758	30000827580521	D684-00-0608	3000084960
3000082758	30000827580521	D684-00-0608	3000084961
3000082759	30000827590002	D684-00-0607	3000084965
3000082759	30000827590009	D684-00-0607	3000084965
3000082759	30000827590010	D684-00-0607	3000084965
3000082759	30000827590013	D684-00-0607	3000084965
3000082759	30000827590023	D684-00-0607	3000084965
3000082760	30000827600234	D684-00-0568-01U	3000089101
3000082760	30000827600234	D684-00-0568-01U	3000089992
3000082760	30000827600237	D684-00-0568-01U	3000089101
3000082760	30000827600237	D684-00-0568-01U	3000089992
3000082760	30000827600238	D684-00-0568-01U	3000089101
3000082760	30000827600238	D684-00-0568-01U	3000089992
3000082760	30000827600239	D684-00-0568-01U	3000089101
3000082760	30000827600239	D684-00-0568-01U	3000089992
3000082760	30000827600243	D684-00-0568-01U	3000089101
3000082760	30000827600243	D684-00-0568-01U	3000089992
3000082761	30000827610028	D684-00-0479-01C	3000085500
3000082761	30000827610028	D684-00-0479-01	3000085756

3000083930	30000839300287	D684-00-0568-01U	3000085877
3000083930	30000839300291	D684-00-0568-01U	3000085876
3000083930	30000839300291	D684-00-0568-01U	3000085877
3000083930	30000839300292	D684-00-0568-01U	3000085876
3000083930	30000839300292	D684-00-0568-01U	3000085877
3000083930	30000839300294	D684-00-0568-01U	3000085876
3000083930	30000839300294	D684-00-0568-01U	3000085877
3000083930	30000839300299	D684-00-0568-01U	3000085876
3000083930	30000839300299	D684-00-0568-01U	3000085877
3000083935	30000839350162	D684-00-0296-01	3000086995
3000083935	30000839350162	D684-00-0296-02	3000086996
3000083935	30000839350162	D684-00-0296-01	3000087888
3000083935	30000839350162	D684-00-0296-01U	3000088376
3000083935	30000839350162	D684-00-0498-01	3000090852
3000083935	30000839350162	D684-00-0296-02	3000091755
3000083935	30000839350165	D684-00-0296-01	3000086995
3000083935	30000839350165	D684-00-0296-02	3000086996
3000083935	30000839350165	D684-00-0296-01	3000087888
3000083935	30000839350165	D684-00-0296-01U	3000088376
3000083935	30000839350165	D684-00-0498-01	3000090852
3000083935	30000839350165	D684-00-0296-02	3000091755
3000083935	30000839350166	D684-00-0296-01	3000086995
3000083935	30000839350166	D684-00-0296-02	3000086996
3000083935	30000839350166	D684-00-0296-01	3000087888
3000083935	30000839350166	D684-00-0296-01U	3000088376
3000083935	30000839350166	D684-00-0498-01	3000090852
3000083935	30000839350166	D684-00-0296-02	3000091755
3000083935	30000839350164	D684-00-0296-01	3000086995
3000083935	30000839350164	D684-00-0296-02	3000086996
3000083935	30000839350164	D684-00-0296-01	3000087888
3000083935	30000839350164	D684-00-0296-01U	3000088376
3000083935	30000839350164	D684-00-0498-01	3000090852
3000083935	30000839350164	D684-00-0296-02	3000091755
3000083935	30000839350167	D684-00-0296-01	3000086995
3000083935	30000839350167	D684-00-0296-02	3000086996
3000083935	30000839350167	D684-00-0296-01	3000087888
3000083935	30000839350167	D684-00-0296-01U	3000088376
3000083935	30000839350167	D684-00-0498-01	3000090852
3000083935	30000839350167	D684-00-0296-02	3000091755
3000083940	30000839400085	D684-00-0607	3000086799
3000083940	30000839400085	D684-00-0607	3000088548
3000083940	30000839400089	D684-00-0607	3000086799
3000083940	30000839400089	D684-00-0607	3000088548
3000083940	30000839400097	D684-00-0607	3000086799
3000083940	30000839400097	D684-00-0607	3000088548
3000083940	30000839400100	D684-00-0607	3000086799
3000083940	30000839400100	D684-00-0607	3000088548
3000083940	30000839400088	D684-00-0607	3000086799
3000083940	30000839400088	D684-00-0607	3000088548
3000083941	30000839410090	D684-00-0608	3000086800
3000083941	30000839410091	D684-00-0608	3000086800
3000083941	30000839410098	D684-00-0608	3000086800
3000083941	30000839410085	D684-00-0608	3000086800
3000083941	30000839410099	D684-00-0608	3000086800
3000083942	30000839420198	D684-00-0568-01	3000086994
3000083942	30000839420198	D684-00-0568-01U	3000087113
3000083942	30000839420198	D684-00-0568-01	3000087688
3000083942	30000839420210	D684-00-0568-01	3000086994

[illegible]

[illegible]

[illegible]

[illegible]

[illegible]

3000085803	30000858030083	D684-00-0480-01	3000088594
3000085803	30000858030083	D684-00-0480-01U	3000088620
3000085803	30000858030083	D684-00-0480-01	3000088639
3000085803	30000858030084	D684-00-0480-02U	3000088490
3000085803	30000858030084	D684-00-0480-01	3000088594
3000085803	30000858030084	D684-00-0480-01U	3000088620
3000085803	30000858030084	D684-00-0480-01	3000088639
3000085803	30000858030085	D684-00-0480-02U	3000088490
3000085803	30000858030085	D684-00-0480-01	3000088594
3000085803	30000858030085	D684-00-0480-01U	3000088620
3000085803	30000858030085	D684-00-0480-01	3000088639
3000085806	30000858060193	D684-00-0607	3000088371
3000085806	30000858060196	D684-00-0607	3000088371
3000085806	30000858060197	D684-00-0607	3000088371
3000085806	30000858060198	D684-00-0607	3000088371
3000085806	30000858060214	D684-00-0607	3000088371
3000085807	30000858070342	D684-00-0608	3000088372
3000085807	30000858070344	D684-00-0608	3000088372
3000085807	30000858070348	D684-00-0608	3000088372
3000085807	30000858070349	D684-00-0608	3000088372
3000085807	30000858070350	D684-00-0608	3000088372
3000085808	30000858080167	D684-00-0469-01U	3000088617
3000085808	30000858080167	D684-00-0469-01	3000088618
3000085808	30000858080167	D684-00-0469-01C	3000089348
3000085808	30000858080167	D684-00-0469-01	3000089991
3000085808	30000858080170	D684-00-0469-01U	3000088617
3000085808	30000858080170	D684-00-0469-01	3000088618
3000085808	30000858080170	D684-00-0469-01C	3000089348
3000085808	30000858080170	D684-00-0469-01	3000089991
3000085808	30000858080171	D684-00-0469-01U	3000088617
3000085808	30000858080171	D684-00-0469-01	3000088618
3000085808	30000858080171	D684-00-0469-01C	3000089348
3000085808	30000858080171	D684-00-0469-01	3000089991
3000085808	30000858080173	D684-00-0469-01U	3000088617
3000085808	30000858080173	D684-00-0469-01	3000088618
3000085808	30000858080173	D684-00-0469-01C	3000089348
3000085808	30000858080173	D684-00-0469-01	3000089991
3000085808	30000858080186	D684-00-0469-01U	3000088617
3000085808	30000858080186	D684-00-0469-01	3000088618
3000085808	30000858080186	D684-00-0469-01C	3000089348
3000085808	30000858080186	D684-00-0469-01	3000089991
3000085809	30000858090641	D684-00-0479-01C	3000089458
3000085809	30000858090641	D684-00-0479-01	3000089610
3000085809	30000858090641	D684-00-0479-01U	3000089732
3000085809	30000858090642	D684-00-0479-01C	3000089458
3000085809	30000858090642	D684-00-0479-01	3000089610
3000085809	30000858090642	D684-00-0479-01U	3000089732
3000085809	30000858090643	D684-00-0479-01C	3000089458
3000085809	30000858090643	D684-00-0479-01	3000089610
3000085809	30000858090643	D684-00-0479-01U	3000089732
3000085809	30000858090644	D684-00-0479-01C	3000089458
3000085809	30000858090644	D684-00-0479-01	3000089610
3000085809	30000858090644	D684-00-0479-01U	3000089732
3000085809	30000858090645	D684-00-0479-01C	3000089458
3000085809	30000858090645	D684-00-0479-01	3000089610
3000085809	30000858090645	D684-00-0479-01U	3000089732
3000086274	30000862740059	D684-00-0568-01U	3000091204
3000086274	30000862740059	D684-00-0568-01	3000091205

[illegible]

3000086277	30000862770134	D684-00-0480-01U	3000089731
3000086277	30000862770135	D684-00-0480-01C	3000089459
3000086277	30000862770135	D684-00-0480-02U	3000089725
3000086277	30000862770135	D684-00-0480-01	3000089730
3000086277	30000862770135	D684-00-0480-01U	3000089731
3000086896	30000868960605	D684-00-0480-02	3000088621
3000086896	30000868960605	D684-00-0480-01U	3000090848
3000086896	30000868960605	D684-00-0480-01C	3000090849
3000086896	30000868960606	D684-00-0480-02	3000088621
3000086896	30000868960606	D684-00-0480-01U	3000090848
3000086896	30000868960606	D684-00-0480-01C	3000090849
3000086896	30000868960607	D684-00-0480-02	3000088621
3000086896	30000868960607	D684-00-0480-01U	3000090848
3000086896	30000868960607	D684-00-0480-01C	3000090849
3000086896	30000868960608	D684-00-0480-02	3000088621
3000086896	30000868960608	D684-00-0480-01U	3000090848
3000086896	30000868960608	D684-00-0480-01C	3000090849
3000086896	30000868960609	D684-00-0480-02	3000088621
3000086896	30000868960609	D684-00-0480-01U	3000090848
3000086896	30000868960609	D684-00-0480-01C	3000090849
3000086897	30000868970548	D684-00-0576-01U	3000089988
3000086897	30000868970550	D684-00-0576-01U	3000089988
3000086897	30000868970553	D684-00-0576-01U	3000089988
3000086897	30000868970561	D684-00-0576-01U	3000089988
3000086897	30000868970562	D684-00-0576-01U	3000089988
3000087219	30000872190241	D684-00-0295-05	3000091874
3000087219	30000872190241	D684-00-0295-01U	3000092129
3000087219	30000872190241	D684-00-0295-01	3000092131
3000087219	30000872190241	D684-00-0295-01	3000092215
3000087219	30000872190242	D684-00-0295-05	3000091874
3000087219	30000872190242	D684-00-0295-01U	3000092129
3000087219	30000872190242	D684-00-0295-01	3000092131
3000087219	30000872190242	D684-00-0295-01	3000092215
3000087219	30000872190243	D684-00-0295-05	3000091874
3000087219	30000872190243	D684-00-0295-01U	3000092129
3000087219	30000872190243	D684-00-0295-01	3000092131
3000087219	30000872190243	D684-00-0295-01	3000092215
3000087219	30000872190244	D684-00-0295-05	3000091874
3000087219	30000872190244	D684-00-0295-01U	3000092129
3000087219	30000872190244	D684-00-0295-01	3000092131
3000087219	30000872190244	D684-00-0295-01	3000092215
3000087219	30000872190245	D684-00-0295-05	3000091874
3000087219	30000872190245	D684-00-0295-01U	3000092129
3000087219	30000872190245	D684-00-0295-01	3000092131
3000087219	30000872190245	D684-00-0295-01	3000092215
3000087225	30000872250647	D684-00-0470-01	3000090526
3000087225	30000872250647	D684-00-0470-01	3000090664
3000087225	30000872250647	D684-00-0470-01U	3000090665
3000087225	30000872250650	D684-00-0470-01	3000090526
3000087225	30000872250650	D684-00-0470-01	3000090664
3000087225	30000872250650	D684-00-0470-01U	3000090665
3000087225	30000872250660	D684-00-0470-01	3000090526
3000087225	30000872250660	D684-00-0470-01	3000090664
3000087225	30000872250660	D684-00-0470-01U	3000090665
3000087225	30000872250662	D684-00-0470-01	3000090526
3000087225	30000872250662	D684-00-0470-01	3000090664
3000087225	30000872250662	D684-00-0470-01U	3000090665
3000087225	30000872250664	D684-00-0470-01	3000090526

3000087225	30000872250664	D684-00-0470-01	3000090664
3000087225	30000872250664	D684-00-0470-01U	3000090665
3000087280	30000872800026	D684-00-0296-01U	3000089899
3000087280	30000872800026	D684-DO-0497	0
3000087280	30000872800026	D684-00-0296-01U	3000091754
3000087280	30000872800026	D684-00-0498-01	3000091872
3000087280	30000872800027	D684-00-0296-01U	3000089899
3000087280	30000872800027	D684-DO-0497	0
3000087280	30000872800027	D684-00-0296-01U	3000091754
3000087280	30000872800027	D684-00-0498-01	3000091872
3000087280	30000872800028	D684-00-0296-01U	3000089899
3000087280	30000872800028	D684-DO-0497	0
3000087280	30000872800028	D684-00-0296-01U	3000091754
3000087280	30000872800028	D684-00-0498-01	3000091872
3000087280	30000872800029	D684-00-0296-01U	3000089899
3000087280	30000872800029	D684-DO-0497	0
3000087280	30000872800029	D684-00-0296-01U	3000091754
3000087280	30000872800029	D684-00-0498-01	3000091872
3000087280	30000872800030	D684-00-0296-01U	3000089899
3000087280	30000872800030	D684-DO-0497	0
3000087280	30000872800030	D684-00-0296-01U	3000091754
3000087280	30000872800030	D684-00-0498-01	3000091872
3000087638	30000876380061	D684-00-0294-01U	3000090374
3000087638	30000876380061	D684-00-0294-01U	3000092556
3000087638	30000876380061	D684-00-0294-01	3000096716
3000087638	30000876380062	D684-00-0294-01U	3000090374
3000087638	30000876380062	D684-00-0294-01U	3000092556
3000087638	30000876380062	D684-00-0294-01	3000096716
3000087638	30000876380063	D684-00-0294-01U	3000090374
3000087638	30000876380063	D684-00-0294-01U	3000092556
3000087638	30000876380063	D684-00-0294-01	3000096716
3000087638	30000876380064	D684-00-0294-01U	3000090374
3000087638	30000876380064	D684-00-0294-01U	3000092556
3000087638	30000876380064	D684-00-0294-01	3000096716
3000087638	30000876380065	D684-00-0294-01U	3000090374
3000087638	30000876380065	D684-00-0294-01U	3000092556
3000087638	30000876380065	D684-00-0294-01	3000096716
3000087641	30000876410004	D684-00-0480-01	3000090853
3000087641	30000876410004	D684-00-0480-01	3000091126
3000087641	30000876410005	D684-00-0480-01	3000090853
3000087641	30000876410005	D684-00-0480-01	3000091126
3000087641	30000876410006	D684-00-0480-01	3000090853
3000087641	30000876410006	D684-00-0480-01	3000091126
3000087641	30000876410007	D684-00-0480-01	3000090853
3000087641	30000876410007	D684-00-0480-01	3000091126
3000087641	30000876410008	D684-00-0480-01	3000090853
3000087641	30000876410008	D684-00-0480-01	3000091126
3000087647	30000876470061	D684-00-0607	3000090599
3000087647	30000876470061	D684-00-0607	3000090600
3000087647	30000876470061	D684-00-0607	3000090601
3000087647	30000876470061	D684-00-0607	3000090602
3000087647	30000876470061	D684-00-0607	3000090603
3000087647	30000876470061	D684-00-0607	3000092664
3000087647	30000876470067	D684-00-0607	3000090599
3000087647	30000876470067	D684-00-0607	3000090600
3000087647	30000876470067	D684-00-0607	3000090601
3000087647	30000876470067	D684-00-0607	3000090602
3000087647	30000876470067	D684-00-0607	3000090603

3000088933	30000889330162	D684-00-0479-01	3000091698
3000088933	30000889330162	D684-00-0479-01C	3000091740
3000088933	30000889330163	D684-00-0479-02	3000091478
3000088933	30000889330163	D684-00-0479-02U	3000091479
3000088933	30000889330163	D684-00-0479-02	3000091697
3000088933	30000889330163	D684-00-0479-01	3000091698
3000088933	30000889330163	D684-00-0479-01C	3000091740
3000088933	30000889330164	D684-00-0479-02	3000091478
3000088933	30000889330164	D684-00-0479-02U	3000091479
3000088933	30000889330164	D684-00-0479-02	3000091697
3000088933	30000889330164	D684-00-0479-01	3000091698
3000088933	30000889330164	D684-00-0479-01C	3000091740
3000088933	30000889330165	D684-00-0479-02	3000091478
3000088933	30000889330165	D684-00-0479-02U	3000091479
3000088933	30000889330165	D684-00-0479-02	3000091697
3000088933	30000889330165	D684-00-0479-01	3000091698
3000088933	30000889330165	D684-00-0479-01C	3000091740
3000088979	30000889790032	D684-00-0607	3000092663
3000088979	30000889790041	D684-00-0607	3000092663
3000088979	30000889790045	D684-00-0607	3000092663
3000088979	30000889790049	D684-00-0607	3000092663
3000088979	30000889790051	D684-00-0607	3000092663
3000088981	30000889810016	D684-00-0608	3000092667
3000088981	30000889810021	D684-00-0608	3000092667
3000088981	30000889810036	D684-00-0608	3000092667
3000088981	30000889810041	D684-00-0608	3000092667
3000088981	30000889810042	D684-00-0608	3000092667
3000089607	30000896070105	D684-00-0480-01C	3000092396
3000089607	30000896070105	D684-00-0480-01	3000092503
3000089607	30000896070106	D684-00-0480-01C	3000092396
3000089607	30000896070106	D684-00-0480-01	3000092503
3000089607	30000896070107	D684-00-0480-01C	3000092396
3000089607	30000896070107	D684-00-0480-01	3000092503
3000089607	30000896070108	D684-00-0480-01C	3000092396
3000089607	30000896070108	D684-00-0480-01	3000092503
3000089607	30000896070109	D684-00-0480-01C	3000092396
3000089607	30000896070109	D684-00-0480-01	3000092503
3000089609	30000896090642	D684-00-0470-01U	3000092442
3000089609	30000896090642	D684-00-0470-01	3000092443
3000089609	30000896090647	D684-00-0470-01U	3000092442
3000089609	30000896090647	D684-00-0470-01	3000092443
3000089609	30000896090649	D684-00-0470-01U	3000092442
3000089609	30000896090649	D684-00-0470-01	3000092443
3000089609	30000896090661	D684-00-0470-01U	3000092442
3000089609	30000896090661	D684-00-0470-01	3000092443
3000089609	30000896090662	D684-00-0470-01U	3000092442
3000089609	30000896090662	D684-00-0470-01	3000092443
3000089611	30000896110235	D684-00-0296-01U	3000092134
3000089611	30000896110235	D684-00-0296-01U	3000092134
3000089611	30000896110235	D684-00-0296-02	3000092238
3000089611	30000896110235	D684-00-0296-01	3000092239
3000089611	30000896110235	D684-00-0296-01	3000092240
3000089611	30000896110235	D684-00-0296-01U	3000094449
3000089611	30000896110235	D684-00-0296-01	3000094450
3000089611	30000896110236	D684-00-0296-01U	3000092134
3000089611	30000896110236	D684-00-0296-01U	3000092134
3000089611	30000896110236	D684-00-0296-02	3000092238
3000089611	30000896110236	D684-00-0296-01	3000092239

3000089926	30000899260170	D684-00-0608	3000096268
3000089933	30000899330078	D684-00-0568-01	3000092976
3000089933	30000899330078	D684-00-0568-01U	3000093060
3000089933	30000899330089	D684-00-0568-01	3000092976
3000089933	30000899330089	D684-00-0568-01U	3000093060
3000089933	30000899330090	D684-00-0568-01	3000092976
3000089933	30000899330090	D684-00-0568-01U	3000093060
3000089933	30000899330091	D684-00-0568-01	3000092976
3000089933	30000899330091	D684-00-0568-01U	3000093060
3000089933	30000899330095	D684-00-0568-01	3000092976
3000089933	30000899330095	D684-00-0568-01U	3000093060
3000089934	30000899340045	D684-00-0480-01C	3000093113
3000089934	30000899340045	D684-00-0480-02	3000093115
3000089934	30000899340031	D684-00-0480-01C	3000093113
3000089934	30000899340031	D684-00-0480-02	3000093115
3000089934	30000899340033	D684-00-0480-01C	3000093113
3000089934	30000899340033	D684-00-0480-02	3000093115
3000089934	30000899340037	D684-00-0480-01C	3000093113
3000089934	30000899340037	D684-00-0480-02	3000093115
3000089934	30000899340036	D684-00-0480-01C	3000093113
3000089934	30000899340036	D684-00-0480-02	3000093115
3000089936	30000899360322	D684-00-0576-01U	3000093112
3000089936	30000899360322	D684-00-0576-09	3000094529
3000089936	30000899360322	D684-00-0576-09	3000094530
3000089936	30000899360325	D684-00-0576-01U	3000093112
3000089936	30000899360325	D684-00-0576-09	3000094529
3000089936	30000899360325	D684-00-0576-09	3000094530
3000089936	30000899360327	D684-00-0576-01U	3000093112
3000089936	30000899360327	D684-00-0576-09	3000094529
3000089936	30000899360327	D684-00-0576-09	3000094530
3000089936	30000899360332	D684-00-0576-01U	3000093112
3000089936	30000899360332	D684-00-0576-09	3000094529
3000089936	30000899360332	D684-00-0576-09	3000094530
3000089936	30000899360335	D684-00-0576-01U	3000093112
3000089936	30000899360335	D684-00-0576-09	3000094529
3000089936	30000899360335	D684-00-0576-09	3000094530
3000089937	30000899370005	D684-00-0295-05	3000093117
3000089937	30000899370005	D684-00-0295-05	3000093118
3000089937	30000899370005	D684-00-0295-05	3000093119
3000089937	30000899370005	D684-00-0295-05	3000094614
3000089937	30000899370006	D684-00-0295-05	3000093117
3000089937	30000899370006	D684-00-0295-05	3000093118
3000089937	30000899370006	D684-00-0295-05	3000093119
3000089937	30000899370006	D684-00-0295-05	3000094614
3000089937	30000899370011	D684-00-0295-05	3000093117
3000089937	30000899370011	D684-00-0295-05	3000093118
3000089937	30000899370011	D684-00-0295-05	3000093119
3000089937	30000899370011	D684-00-0295-05	3000094614
3000089937	30000899370013	D684-00-0295-05	3000093117
3000089937	30000899370013	D684-00-0295-05	3000093118
3000089937	30000899370013	D684-00-0295-05	3000093119
3000089937	30000899370013	D684-00-0295-05	3000094614
3000089937	30000899370015	D684-00-0295-05	3000093117
3000089937	30000899370015	D684-00-0295-05	3000093118
3000089937	30000899370015	D684-00-0295-05	3000093119
3000089937	30000899370015	D684-00-0295-05	3000094614
3000090431	30000904310190	D684-00-0568-01	3000093557
3000090431	30000904310190	D684-00-0568-01U	3000093558

3000090928	30000909280541	D684-00-0470-07	3000094141
3000090928	30000909280541	D684-00-0470-01C	3000094526
3000090928	30000909280541	D684-00-0470-09	3000094527
3000090928	30000909280541	D684-00-0470-09	3000094527
3000090928	30000909280541	D684-00-0470-01U	3000094528
3000090928	30000909280541	D684-00-0470-01U	3000094531
3000090928	30000909280541	D684-00-0470-01U	3000096053
3000090929	30000909290071	D684-00-0479-01	3000093895
3000090929	30000909290071	D684-00-0479-02	3000093943
3000090929	30000909290071	D684-00-0479-01U	3000094007
3000090929	30000909290071	D684-00-0479-10	3000094153
3000090929	30000909290072	D684-00-0479-01	3000093895
3000090929	30000909290072	D684-00-0479-02	3000093943
3000090929	30000909290072	D684-00-0479-01U	3000094007
3000090929	30000909290072	D684-00-0479-10	3000094153
3000090929	30000909290073	D684-00-0479-01	3000093895
3000090929	30000909290073	D684-00-0479-02	3000093943
3000090929	30000909290073	D684-00-0479-01U	3000094007
3000090929	30000909290073	D684-00-0479-10	3000094153
3000090929	30000909290074	D684-00-0479-01	3000093895
3000090929	30000909290074	D684-00-0479-02	3000093943
3000090929	30000909290074	D684-00-0479-01U	3000094007
3000090929	30000909290074	D684-00-0479-10	3000094153
3000090929	30000909290075	D684-00-0479-01	3000093895
3000090929	30000909290075	D684-00-0479-02	3000093943
3000090929	30000909290075	D684-00-0479-01U	3000094007
3000090929	30000909290075	D684-00-0479-10	3000094153
3000091319	30000913190011	D684-00-0480-01C	3000094005
3000091319	30000913190011	D684-00-0480-10	3000094152
3000091319	30000913190011	D684-00-0480-02	3000095291
3000091319	30000913190012	D684-00-0480-01C	3000094005
3000091319	30000913190012	D684-00-0480-10	3000094152
3000091319	30000913190012	D684-00-0480-02	3000095291
3000091319	30000913190013	D684-00-0480-01C	3000094005
3000091319	30000913190013	D684-00-0480-10	3000094152
3000091319	30000913190013	D684-00-0480-02	3000095291
3000091319	30000913190014	D684-00-0480-01C	3000094005
3000091319	30000913190014	D684-00-0480-10	3000094152
3000091319	30000913190014	D684-00-0480-02	3000095291
3000091319	30000913190015	D684-00-0480-01C	3000094005
3000091319	30000913190015	D684-00-0480-10	3000094152
3000091319	30000913190015	D684-00-0480-02	3000095291
3000091323	30000913230001	D684-00-0607	3000095467
3000091323	30000913230001	D684-00-0607	3000096264
3000091323	30000913230002	D684-00-0607	3000095467
3000091323	30000913230002	D684-00-0607	3000096264
3000091323	30000913230005	D684-00-0607	3000095467
3000091323	30000913230005	D684-00-0607	3000096264
3000091323	30000913230009	D684-00-0607	3000095467
3000091323	30000913230009	D684-00-0607	3000096264
3000091323	30000913230022	D684-00-0607	3000095467
3000091323	30000913230022	D684-00-0607	3000096264
3000091324	30000913240197	D684-00-0608	3000095468
3000091324	30000913240197	D684-00-0608	3000096267
3000091324	30000913240199	D684-00-0608	3000095468
3000091324	30000913240199	D684-00-0608	3000096267
3000091324	30000913240210	D684-00-0608	3000095468
3000091324	30000913240210	D684-00-0608	3000096267

3000091324	30000913240217	D684-00-0608	3000095468
3000091324	30000913240217	D684-00-0608	3000096267
3000091324	30000913240220	D684-00-0608	3000095468
3000091324	30000913240220	D684-00-0608	3000096267
3000091325	30000913250105	D684-00-0560-01	3000096377
3000091325	30000913250105	D684-00-0560-01	3000096378
3000091325	30000913250105	D684-00-0560-01	3000098752
3000091325	30000913250105	D684-00-0560-01	3000099588
3000091325	30000913250106	D684-00-0560-01	3000096377
3000091325	30000913250106	D684-00-0560-01	3000096378
3000091325	30000913250106	D684-00-0560-01	3000098752
3000091325	30000913250106	D684-00-0560-01	3000099588
3000091325	30000913250107	D684-00-0560-01	3000096377
3000091325	30000913250107	D684-00-0560-01	3000096378
3000091325	30000913250107	D684-00-0560-01	3000098752
3000091325	30000913250107	D684-00-0560-01	3000099588
3000091325	30000913250108	D684-00-0560-01	3000096377
3000091325	30000913250108	D684-00-0560-01	3000096378
3000091325	30000913250108	D684-00-0560-01	3000098752
3000091325	30000913250108	D684-00-0560-01	3000099588
3000091325	30000913250109	D684-00-0560-01	3000096377
3000091325	30000913250109	D684-00-0560-01	3000096378
3000091325	30000913250109	D684-00-0560-01	3000098752
3000091325	30000913250109	D684-00-0560-01	3000099588
3000091326	30000913260106	D684-00-0469-01U	3000094336
3000091326	30000913260106	D684-00-0469-01	3000094337
3000091326	30000913260106	D684-00-0469-01C	3000094338
3000091326	30000913260119	D684-00-0469-01U	3000094336
3000091326	30000913260119	D684-00-0469-01	3000094337
3000091326	30000913260119	D684-00-0469-01C	3000094338
3000091326	30000913260120	D684-00-0469-01U	3000094336
3000091326	30000913260120	D684-00-0469-01	3000094337
3000091326	30000913260120	D684-00-0469-01C	3000094338
3000091326	30000913260121	D684-00-0469-01U	3000094336
3000091326	30000913260121	D684-00-0469-01	3000094337
3000091326	30000913260121	D684-00-0469-01C	3000094338
3000091326	30000913260122	D684-00-0469-01U	3000094336
3000091326	30000913260122	D684-00-0469-01	3000094337
3000091326	30000913260122	D684-00-0469-01C	3000094338
3000091784	30000917840206	D684-00-0498-01	3000094108
3000091784	30000917840206	D684-00-0296-01	3000094451
3000091784	30000917840206	D684-00-0498-01	3000094935
3000091784	30000917840206	D684-DO-0497	0
3000091784	30000917840206	D684-DO-0497	0
3000091784	30000917840206	D684-DC-0497	0
3000091784	30000917840206	D684-00-0498-01	3000096067
3000091784	30000917840207	D684-00-0498-01	3000094108
3000091784	30000917840207	D684-00-0296-01	3000094451
3000091784	30000917840207	D684-00-0498-01	3000094935
3000091784	30000917840207	D684-DO-0497	0
3000091784	30000917840207	D684-DO-0497	0
3000091784	30000917840207	D684-DC-0497	0
3000091784	30000917840207	D684-00-0498-01	3000096067
3000091784	30000917840208	D684-00-0498-01	3000094108
3000091784	30000917840208	D684-00-0296-01	3000094451
3000091784	30000917840208	D684-00-0498-01	3000094935
3000091784	30000917840208	D684-DO-0497	0
3000091784	30000917840208	D684-DO-0497	0

3000091784	30000917840208	D684-DC-0497	0
3000091784	30000917840208	D684-00-0498-01	3000096067
3000091784	30000917840209	D684-00-0498-01	3000094108
3000091784	30000917840209	D684-00-0296-01	3000094451
3000091784	30000917840209	D684-00-0498-01	3000094935
3000091784	30000917840209	D684-DO-0497	0
3000091784	30000917840209	D684-DO-0497	0
3000091784	30000917840209	D684-DC-0497	0
3000091784	30000917840209	D684-00-0498-01	3000096067
3000091784	30000917840210	D684-00-0498-01	3000094108
3000091784	30000917840210	D684-00-0296-01	3000094451
3000091784	30000917840210	D684-00-0498-01	3000094935
3000091784	30000917840210	D684-DO-0497	0
3000091784	30000917840210	D684-DO-0497	0
3000091784	30000917840210	D684-DC-0497	0
3000091784	30000917840210	D684-00-0498-01	3000096067
3000091785	30000917850402	D684-00-0568-01U	3000094740
3000091785	30000917850407	D684-00-0568-01U	3000094740
3000091785	30000917850408	D684-00-0568-01U	3000094740
3000091785	30000917850409	D684-00-0568-01U	3000094740
3000091785	30000917850410	D684-00-0568-01U	3000094740
3000091786	30000917860580	D684-00-0479-01C	3000094658
3000091786	30000917860580	D684-00-0479-10	3000094937
3000091786	30000917860580	D684-00-0479-02U	3000094938
3000091786	30000917860580	D684-00-0479-01C	3000095998
3000091786	30000917860581	D684-00-0479-01C	3000094658
3000091786	30000917860581	D684-00-0479-10	3000094937
3000091786	30000917860581	D684-00-0479-02U	3000094938
3000091786	30000917860581	D684-00-0479-01C	3000095998
3000091786	30000917860582	D684-00-0479-01C	3000094658
3000091786	30000917860582	D684-00-0479-10	3000094937
3000091786	30000917860582	D684-00-0479-02U	3000094938
3000091786	30000917860582	D684-00-0479-01C	3000095998
3000091786	30000917860583	D684-00-0479-01C	3000094658
3000091786	30000917860583	D684-00-0479-10	3000094937
3000091786	30000917860583	D684-00-0479-02U	3000094938
3000091786	30000917860583	D684-00-0479-01C	3000095998
3000091786	30000917860585	D684-00-0479-01C	3000094658
3000091786	30000917860585	D684-00-0479-10	3000094937
3000091786	30000917860585	D684-00-0479-02U	3000094938
3000091786	30000917860585	D684-00-0479-01C	3000095998
3000091787	30000917870368	D684-00-0576-01U	3000094980
3000091787	30000917870368	D684-00-0576-01	3000095036
3000091787	30000917870371	D684-00-0576-01U	3000094980
3000091787	30000917870371	D684-00-0576-01	3000095036
3000091787	30000917870375	D684-00-0576-01U	3000094980
3000091787	30000917870375	D684-00-0576-01	3000095036
3000091787	30000917870376	D684-00-0576-01U	3000094980
3000091787	30000917870376	D684-00-0576-01	3000095036
3000091787	30000917870381	D684-00-0576-01U	3000094980
3000091787	30000917870381	D684-00-0576-01	3000095036
3000092228	30000922280252	D684-00-0295-05	3000094613
3000092228	30000922280253	D684-00-0295-05	3000094613
3000092228	30000922280260	D684-00-0295-05	3000094613
3000092228	30000922280262	D684-00-0295-05	3000094613
3000092228	30000922280263	D684-00-0295-05	3000094613
3000092230	30000922300565	D684-00-0480-01C	3000095038
3000092230	30000922300565	D684-00-0480-10	3000095152

3000092230	30000922300566	D684-00-0480-01C	3000095038
3000092230	30000922300566	D684-00-0480-10	3000095152
3000092230	30000922300567	D684-00-0480-01C	3000095038
3000092230	30000922300567	D684-00-0480-10	3000095152
3000092230	30000922300568	D684-00-0480-01C	3000095038
3000092230	30000922300568	D684-00-0480-10	3000095152
3000092230	30000922300569	D684-00-0480-01C	3000095038
3000092230	30000922300569	D684-00-0480-10	3000095152
3000092231	30000922310456	D684-00-0568-01	3000095037
3000092231	30000922310456	D684-00-0568-01U	3000095150
3000092231	30000922310456	D684-00-0568-01U	3000095151
3000092231	30000922310458	D684-00-0568-01	3000095037
3000092231	30000922310458	D684-00-0568-01U	3000095150
3000092231	30000922310458	D684-00-0568-01U	3000095151
3000092231	30000922310461	D684-00-0568-01	3000095037
3000092231	30000922310461	D684-00-0568-01U	3000095150
3000092231	30000922310461	D684-00-0568-01U	3000095151
3000092231	30000922310465	D684-00-0568-01	3000095037
3000092231	30000922310465	D684-00-0568-01U	3000095150
3000092231	30000922310465	D684-00-0568-01U	3000095151
3000092231	30000922310467	D684-00-0568-01	3000095037
3000092231	30000922310467	D684-00-0568-01U	3000095150
3000092231	30000922310467	D684-00-0568-01U	3000095151
3000092766	30000927660155	D684-00-0576-01U	3000095605
3000092766	30000927660156	D684-00-0576-01U	3000095605
3000092766	30000927660159	D684-00-0576-01U	3000095605
3000092766	30000927660165	D684-00-0576-01U	3000095605
3000092766	30000927660166	D684-00-0576-01U	3000095605
3000092767	30000927670111	D684-00-0295-05	3000095319
3000092767	30000927670112	D684-00-0295-05	3000095319
3000092767	30000927670113	D684-00-0295-05	3000095319
3000092767	30000927670114	D684-00-0295-05	3000095319
3000092767	30000927670115	D684-00-0295-05	3000095319
3000092772	30000927720031	D684-00-0498-01	3000095318
3000092772	30000927720031	D684-00-0296-01	3000095604
3000092772	30000927720031	D684-00-0296-01U	3000095692
3000092772	30000927720031	D684-00-0498-01	3000096066
3000092772	30000927720032	D684-00-0498-01	3000095318
3000092772	30000927720032	D684-00-0296-01	3000095604
3000092772	30000927720032	D684-00-0296-01U	3000095692
3000092772	30000927720032	D684-00-0498-01	3000096066
3000092772	30000927720033	D684-00-0498-01	3000095318
3000092772	30000927720033	D684-00-0296-01	3000095604
3000092772	30000927720033	D684-00-0296-01U	3000095692
3000092772	30000927720033	D684-00-0498-01	3000096066
3000092772	30000927720034	D684-00-0498-01	3000095318
3000092772	30000927720034	D684-00-0296-01	3000095604
3000092772	30000927720034	D684-00-0296-01U	3000095692
3000092772	30000927720034	D684-00-0498-01	3000096066
3000092772	30000927720035	D684-00-0498-01	3000095318
3000092772	30000927720035	D684-00-0296-01	3000095604
3000092772	30000927720035	D684-00-0296-01U	3000095692
3000092772	30000927720035	D684-00-0498-01	3000096066
3000092803	30000928030221	D684-00-0480-01U	3000095289
3000092803	30000928030221	D684-00-0480-01	3000095290
3000092803	30000928030221	D684-00-0480-02	3000095796
3000092803	30000928030222	D684-00-0480-01U	3000095289
3000092803	30000928030222	D684-00-0480-01	3000095290

3000093377	30000933770122	D684-00-0470-01U	3000096054
3000093377	30000933770122	D684-00-0470-01	3000097302
3000093378	30000933780369	D684-00-0607	3000096265
3000093378	30000933780369	D684-00-0607	3000096266
3000093378	30000933780369	D684-00-0607	3000098176
3000093378	30000933780374	D684-00-0607	3000096265
3000093378	30000933780374	D684-00-0607	3000096266
3000093378	30000933780374	D684-00-0607	3000098176
3000093378	30000933780379	D684-00-0607	3000096265
3000093378	30000933780379	D684-00-0607	3000096266
3000093378	30000933780379	D684-00-0607	3000098176
3000093378	30000933780382	D684-00-0607	3000096265
3000093378	30000933780382	D684-00-0607	3000096266
3000093378	30000933780382	D684-00-0607	3000098176
3000093378	30000933780385	D684-00-0607	3000096265
3000093378	30000933780385	D684-00-0607	3000096266
3000093378	30000933780385	D684-00-0607	3000098176
3000093885	30000938850001	D684-00-0295-01	3000095792
3000093885	30000938850001	D684-00-0295-01U	3000095793
3000093885	30000938850001	D684-00-0295-02U	3000095794
3000093885	30000938850001	D684-00-0295-05	3000095854
3000093885	30000938850001	D684-00-0295-01	3000096364
3000093885	30000938850001	D684-00-0295-09	3000096365
3000093885	30000938850001	D684-00-0295-02U	3000096366
3000093885	30000938850001	D684-00-0295-01U	3000097645
3000093885	30000938850006	D684-00-0295-01	3000095792
3000093885	30000938850006	D684-00-0295-01U	3000095793
3000093885	30000938850006	D684-00-0295-02U	3000095794
3000093885	30000938850006	D684-00-0295-05	3000095854
3000093885	30000938850006	D684-00-0295-01	3000096364
3000093885	30000938850006	D684-00-0295-09	3000096365
3000093885	30000938850006	D684-00-0295-02U	3000096366
3000093885	30000938850006	D684-00-0295-01U	3000097645
3000093885	30000938850015	D684-00-0295-01	3000095792
3000093885	30000938850015	D684-00-0295-01U	3000095793
3000093885	30000938850015	D684-00-0295-02U	3000095794
3000093885	30000938850015	D684-00-0295-05	3000095854
3000093885	30000938850015	D684-00-0295-01	3000096364
3000093885	30000938850015	D684-00-0295-09	3000096365
3000093885	30000938850015	D684-00-0295-02U	3000096366
3000093885	30000938850015	D684-00-0295-01U	3000097645
3000093885	30000938850023	D684-00-0295-01	3000095792
3000093885	30000938850023	D684-00-0295-01U	3000095793
3000093885	30000938850023	D684-00-0295-02U	3000095794
3000093885	30000938850023	D684-00-0295-05	3000095854
3000093885	30000938850023	D684-00-0295-01	3000096364
3000093885	30000938850023	D684-00-0295-09	3000096365
3000093885	30000938850023	D684-00-0295-02U	3000096366
3000093885	30000938850023	D684-00-0295-01U	3000097645
3000093885	30000938850024	D684-00-0295-01	3000095792
3000093885	30000938850024	D684-00-0295-01U	3000095793
3000093885	30000938850024	D684-00-0295-02U	3000095794
3000093885	30000938850024	D684-00-0295-05	3000095854
3000093885	30000938850024	D684-00-0295-01	3000096364
3000093885	30000938850024	D684-00-0295-09	3000096365
3000093885	30000938850024	D684-00-0295-02U	3000096366
3000093885	30000938850024	D684-00-0295-01U	3000097645
3000093886	30000938860032	D684-00-0479-01C	3000095999

[illegible]

[illegible]

[illegible]

[illegible]

[illegible]

3000095067	30000950670063	D684-00-0295-05	3000097533
3000095067	30000950670063	D684-00-0295-01	3000097534
3000095067	30000950670063	D684-00-0295-01U	3000097643
3000095067	30000950670064	D684-00-0295-05	3000097533
3000095067	30000950670064	D684-00-0295-01	3000097534
3000095067	30000950670064	D684-00-0295-01U	3000097643
3000095067	30000950670065	D684-00-0295-05	3000097533
3000095067	30000950670065	D684-00-0295-01	3000097534
3000095067	30000950670065	D684-00-0295-01U	3000097643
3000095071	30000950710002	D684-00-0480-01C	3000098249
3000095071	30000950710002	D684-00-0480-01C	3000098249
3000095071	30000950710002	D684-00-0480-02	3000098613
3000095071	30000950710002	D684-00-0480-01U	3000099419
3000095071	30000950710003	D684-00-0480-01C	3000098249
3000095071	30000950710003	D684-00-0480-01C	3000098249
3000095071	30000950710003	D684-00-0480-02	3000098613
3000095071	30000950710003	D684-00-0480-01U	3000099419
3000095072	30000950720302	D684-00-0607	3000098177
3000095170	30000951700004	D684-00-0480-10	3000097653
3000095170	30000951700004	D684-00-0480-01	3000097654
3000095170	30000951700004	D684-00-0480-01U	3000098615
3000095170	30000951700005	D684-00-0480-10	3000097653
3000095170	30000951700005	D684-00-0480-01	3000097654
3000095170	30000951700005	D684-00-0480-01U	3000098615
3000095170	30000951700002	D684-00-0480-10	3000097653
3000095170	30000951700002	D684-00-0480-01	3000097654
3000095170	30000951700002	D684-00-0480-01U	3000098615
3000095170	30000951700001	D684-00-0480-10	3000097653
3000095170	30000951700001	D684-00-0480-01	3000097654
3000095170	30000951700001	D684-00-0480-01U	3000098615
3000095170	30000951700003	D684-00-0480-10	3000097653
3000095170	30000951700003	D684-00-0480-01	3000097654
3000095170	30000951700003	D684-00-0480-01U	3000098615
3000095707	30000957070348	D684-00-0608	3000098747
3000095707	30000957070348	D684-00-0608	3000098747
3000095707	30000957070348	D684-00-0608	3000098748
3000095707	30000957070348	D684-00-0608	3000098748
3000095707	30000957070354	D684-00-0608	3000098747
3000095707	30000957070354	D684-00-0608	3000098747
3000095707	30000957070354	D684-00-0608	3000098748
3000095707	30000957070354	D684-00-0608	3000098748
3000095707	30000957070355	D684-00-0608	3000098747
3000095707	30000957070355	D684-00-0608	3000098747
3000095707	30000957070355	D684-00-0608	3000098748
3000095707	30000957070355	D684-00-0608	3000098748
3000095707	30000957070360	D684-00-0608	3000098747
3000095707	30000957070360	D684-00-0608	3000098747
3000095707	30000957070360	D684-00-0608	3000098748
3000095707	30000957070360	D684-00-0608	3000098748
3000095707	30000957070364	D684-00-0608	3000098747
3000095707	30000957070364	D684-00-0608	3000098747
3000095707	30000957070364	D684-00-0608	3000098748
3000095707	30000957070364	D684-00-0608	3000098748
3000095708	30000957081048	D684-00-0576-01	3000099219
3000095708	30000957081048	D684-00-0576-01U	3000099220
3000095708	30000957081055	D684-00-0576-01	3000099219
3000095708	30000957081055	D684-00-0576-01U	3000099220
3000095708	30000957081056	D684-00-0576-01	3000099219

3000095708	30000957081056	D684-00-0576-01U	3000099220
3000095708	30000957081058	D684-00-0576-01	3000099219
3000095708	30000957081058	D684-00-0576-01U	3000099220
3000095708	30000957081063	D684-00-0576-01	3000099219
3000095708	30000957081063	D684-00-0576-01U	3000099220
3000095709	30000957090117	D684-00-0559-01	3000098750
3000095709	30000957090118	D684-00-0559-01	3000098750
3000095709	3000095709031	D684-00-0559-01	3000098750
3000095709	30000957090133	D684-00-0559-01	3000098750
3000095709	30000957090136	D684-00-0559-01	3000098750
3000095710	30000957100081	D684-00-0560-01	3000099589
3000095710	30000957100082	D684-00-0560-01	3000099589
3000095710	30000957100083	D684-00-0560-01	3000099589
3000095710	30000957100084	D684-00-0560-01	3000099589
3000095710	30000957100085	D684-00-0560-01	3000099589
3000095711	30000957110076	D684-00-0480-01	3000099417
3000095711	30000957110076	D684-00-0480-02	3000099418
3000095711	30000957110077	D684-00-0480-01	3000099417
3000095711	30000957110077	D684-00-0480-02	3000099418
3000095711	30000957110078	D684-00-0480-01	3000099417
3000095711	30000957110078	D684-00-0480-02	3000099418
3000095711	30000957110079	D684-00-0480-01	3000099417
3000095711	30000957110079	D684-00-0480-02	3000099418
3000095711	30000957110080	D684-00-0480-01	3000099417
3000095711	30000957110080	D684-00-0480-02	3000099418
3000096068	30000960680026	D684-00-0470-01	3000098997
3000096068	30000960680026	D684-00-0470-01U	3000098998
3000096068	30000960680026	D684-00-0470-01U	3000098999
3000096068	30000960680026	D684-00-0470-09	3000099000
3000096068	30000960680028	D684-00-0470-01	3000098997
3000096068	30000960680028	D684-00-0470-01U	3000098998
3000096068	30000960680028	D684-00-0470-01U	3000098999
3000096068	30000960680028	D684-00-0470-09	3000099000
3000096068	30000960680045	D684-00-0470-01	3000098997
3000096068	30000960680045	D684-00-0470-01U	3000098998
3000096068	30000960680045	D684-00-0470-01U	3000098999
3000096068	30000960680045	D684-00-0470-09	3000099000
3000096068	30000960680047	D684-00-0470-01	3000098997
3000096068	30000960680047	D684-00-0470-01U	3000098998
3000096068	30000960680047	D684-00-0470-01U	3000098999
3000096068	30000960680047	D684-00-0470-09	3000099000
3000096068	30000960680048	D684-00-0470-01	3000098997
3000096068	30000960680048	D684-00-0470-01U	3000098998
3000096068	30000960680048	D684-00-0470-01U	3000098999
3000096068	30000960680048	D684-00-0470-09	3000099000
3000096069	30000960690002	D684-00-0479-08	3000100147
3000096069	30000960690002	D684-00-0479-01C	3000100381
3000096069	30000960690002	D684-00-0479-01U	3000100470
3000096069	30000960690002	D684-00-0479-02U	3000100471
3000096069	30000960690002	D684-00-0479-01	3000100472
3000096069	30000960690002	D684-00-0479-02	3000101842
3000096069	30000960690002	D684-00-0479-02U	3000101844
3000096069	30000960690007	D684-00-0479-08	3000100147
3000096069	30000960690007	D684-00-0479-01C	3000100381
3000096069	30000960690007	D684-00-0479-01U	3000100470
3000096069	30000960690007	D684-00-0479-02U	3000100471
3000096069	30000960690007	D684-00-0479-01	3000100472
3000096069	30000960690007	D684-00-0479-02	3000101842

3000096069	30000960690007	D684-00-0479-02U	3000101844
3000096069	30000960690009	D684-00-0479-08	3000100147
3000096069	30000960690009	D684-00-0479-01C	3000100381
3000096069	30000960690009	D684-00-0479-01U	3000100470
3000096069	30000960690009	D684-00-0479-02U	3000100471
3000096069	30000960690009	D684-00-0479-01	3000100472
3000096069	30000960690009	D684-00-0479-02	3000101842
3000096069	30000960690009	D684-00-0479-02U	3000101844
3000096069	30000960690016	D684-00-0479-08	3000100147
3000096069	30000960690016	D684-00-0479-01C	3000100381
3000096069	30000960690016	D684-00-0479-01U	3000100470
3000096069	30000960690016	D684-00-0479-02U	3000100471
3000096069	30000960690016	D684-00-0479-01	3000100472
3000096069	30000960690016	D684-00-0479-02	3000101842
3000096069	30000960690016	D684-00-0479-02U	3000101844
3000096069	30000960690017	D684-00-0479-08	3000100147
3000096069	30000960690017	D684-00-0479-01C	3000100381
3000096069	30000960690017	D684-00-0479-01U	3000100470
3000096069	30000960690017	D684-00-0479-02U	3000100471
3000096069	30000960690017	D684-00-0479-01	3000100472
3000096069	30000960690017	D684-00-0479-02	3000101842
3000096069	30000960690017	D684-00-0479-02U	3000101844
3000096070	30000960700131	D684-00-0480-01U	3000099420
3000096070	30000960700131	D684-00-0480-08	3000100000
3000096070	30000960700131	D684-00-0480-08	3000100000
3000096070	30000960700131	D684-00-0480-02U	3000100067
3000096070	30000960700135	D684-00-0480-01U	3000099420
3000096070	30000960700135	D684-00-0480-08	3000100000
3000096070	30000960700135	D684-00-0480-08	3000100000
3000096070	30000960700135	D684-00-0480-02U	3000100067
3000096070	30000960700140	D684-00-0480-01U	3000099420
3000096070	30000960700140	D684-00-0480-08	3000100000
3000096070	30000960700140	D684-00-0480-08	3000100000
3000096070	30000960700140	D684-00-0480-02U	3000100067
3000096070	30000960700146	D684-00-0480-01U	3000099420
3000096070	30000960700146	D684-00-0480-08	3000100000
3000096070	30000960700146	D684-00-0480-08	3000100000
3000096070	30000960700146	D684-00-0480-02U	3000100067
3000096070	30000960700149	D684-00-0480-01U	3000099420
3000096070	30000960700149	D684-00-0480-08	3000100000
3000096070	30000960700149	D684-00-0480-08	3000100000
3000096070	30000960700149	D684-00-0480-02U	3000100067
3000096071	30000960710107	D684-00-0568-01U	3000098713
3000096071	30000960710109	D684-00-0568-01U	3000098713
3000096071	30000960710112	D684-00-0568-01U	3000098713
3000096071	30000960710114	D684-00-0568-01U	3000098713
3000096071	30000960710125	D684-00-0568-01U	3000098713
3000096331	30000963310078	D684-00-0469-01	3000098247
3000096331	30000963310078	D684-00-0469-01U	3000098248
3000096331	30000963310086	D684-00-0469-01	3000098247
3000096331	30000963310086	D684-00-0469-01U	3000098248
3000096331	30000963310089	D684-00-0469-01	3000098247
3000096331	30000963310089	D684-00-0469-01U	3000098248
3000096331	30000963310092	D684-00-0469-01	3000098247
3000096331	30000963310092	D684-00-0469-01U	3000098248
3000096331	30000963310100	D684-00-0469-01	3000098247
3000096331	30000963310100	D684-00-0469-01U	3000098248
3000096706	30000967060030	D684-00-0295-05	3000099217

3000096706	30000967060030	D684-00-0295-02	3000099218
3000096706	30000967060030	D684-00-0295-05	3000100871
3000096706	30000967060030	D684-00-0295-01	3000100872
3000096706	30000967060030	D684-00-0295-02U	3000102115
3000096706	30000967060034	D684-00-0295-05	3000099217
3000096706	30000967060034	D684-00-0295-02	3000099218
3000096706	30000967060034	D684-00-0295-05	3000100871
3000096706	30000967060034	D684-00-0295-01	3000100872
3000096706	30000967060034	D684-00-0295-02U	3000102115
3000096706	30000967060041	D684-00-0295-05	3000099217
3000096706	30000967060041	D684-00-0295-02	3000099218
3000096706	30000967060041	D684-00-0295-05	3000100871
3000096706	30000967060041	D684-00-0295-01	3000100872
3000096706	30000967060041	D684-00-0295-02U	3000102115
3000096706	30000967060043	D684-00-0295-05	3000099217
3000096706	30000967060043	D684-00-0295-02	3000099218
3000096706	30000967060043	D684-00-0295-05	3000100871
3000096706	30000967060043	D684-00-0295-01	3000100872
3000096706	30000967060043	D684-00-0295-02U	3000102115
3000096706	30000967060045	D684-00-0295-05	3000099217
3000096706	30000967060045	D684-00-0295-02	3000099218
3000096706	30000967060045	D684-00-0295-05	3000100871
3000096706	30000967060045	D684-00-0295-01	3000100872
3000096706	30000967060045	D684-00-0295-02U	3000102115
3000096707	30000967070341	D684-00-0480-02	3000100066
3000096707	30000967070341	D684-00-0480-01U	3000100266
3000096707	30000967070341	D684-00-0480-01	3000100344
3000096707	30000967070341	D684-00-0480-02	3000100346
3000096707	30000967070355	D684-00-0480-02	3000100066
3000096707	30000967070355	D684-00-0480-01U	3000100266
3000096707	30000967070355	D684-00-0480-01	3000100344
3000096707	30000967070355	D684-00-0480-02	3000100346
3000096707	30000967070359	D684-00-0480-02	3000100066
3000096707	30000967070359	D684-00-0480-01U	3000100266
3000096707	30000967070359	D684-00-0480-01	3000100344
3000096707	30000967070359	D684-00-0480-02	3000100346
3000096707	30000967070361	D684-00-0480-02	3000100066
3000096707	30000967070361	D684-00-0480-01U	3000100266
3000096707	30000967070361	D684-00-0480-01	3000100344
3000096707	30000967070361	D684-00-0480-02	3000100346
3000096707	30000967070364	D684-00-0480-02	3000100066
3000096707	30000967070364	D684-00-0480-01U	3000100266
3000096707	30000967070364	D684-00-0480-01	3000100344
3000096707	30000967070364	D684-00-0480-02	3000100346
3000096708	30000967080273	D684-00-0546-01	3000099591
3000096708	30000967080273	D684-00-0546-01	3000112273
3000096708	30000967080280	D684-00-0546-01	3000099591
3000096708	30000967080280	D684-00-0546-01	3000112273
3000096708	30000967080286	D684-00-0546-01	3000099591
3000096708	30000967080286	D684-00-0546-01	3000112273
3000096708	30000967080293	D684-00-0546-01	3000099591
3000096708	30000967080293	D684-00-0546-01	3000112273
3000096708	30000967080294	D684-00-0546-01	3000099591
3000096708	30000967080294	D684-00-0546-01	3000112273
3000096709	30000967090295	D684-00-0568-01U	3000099522
3000096709	30000967090295	D684-00-0568-01U	3000100581
3000096709	30000967090295	D684-00-0568-01	3000100585
3000096709	30000967090301	D684-00-0568-01U	3000099522

3000096709	30000967090301	D684-00-0568-01U	3000100581
3000096709	30000967090301	D684-00-0568-01	3000100585
3000096709	30000967090305	D684-00-0568-01U	3000099522
3000096709	30000967090305	D684-00-0568-01U	3000100581
3000096709	30000967090305	D684-00-0568-01	3000100585
3000096709	30000967090307	D684-00-0568-01U	3000099522
3000096709	30000967090307	D684-00-0568-01U	3000100581
3000096709	30000967090307	D684-00-0568-01	3000100585
3000096709	30000967090308	D684-00-0568-01U	3000099522
3000096709	30000967090308	D684-00-0568-01U	3000100581
3000096709	30000967090308	D684-00-0568-01	3000100585
3000096944	30000969440106	D684-00-0480-01C	3000100378
3000096944	30000969440106	D684-00-0480-08	3000101129
3000096944	30000969440116	D684-00-0480-01C	3000100378
3000096944	30000969440116	D684-00-0480-08	3000101129
3000096944	30000969440118	D684-00-0480-01C	3000100378
3000096944	30000969440118	D684-00-0480-08	3000101129
3000096944	30000969440126	D684-00-0480-01C	3000100378
3000096944	30000969440126	D684-00-0480-08	3000101129
3000096944	30000969440128	D684-00-0480-01C	3000100378
3000096944	30000969440128	D684-00-0480-08	3000101129
3000096945	30000969450176	D684-00-0576-01U	3000099573
3000096945	30000969450176	D684-00-0576-01U	3000099927
3000096945	30000969450179	D684-00-0576-01U	3000099573
3000096945	30000969450179	D684-00-0576-01U	3000099927
3000096945	30000969450180	D684-00-0576-01U	3000099573
3000096945	30000969450180	D684-00-0576-01U	3000099927
3000096945	30000969450190	D684-00-0576-01U	3000099573
3000096945	30000969450190	D684-00-0576-01U	3000099927
3000096945	30000969450191	D684-00-0576-01U	3000099573
3000096945	30000969450191	D684-00-0576-01U	3000099927
3000096946	30000969460080	D684-00-0498-01	3000100582
3000096946	30000969460080	D684-00-0296-01U	3000100721
3000096946	30000969460080	D684-00-0296-01	3000100722
3000096946	30000969460080	D684-00-0296-01	3000105063
3000096946	30000969460081	D684-00-0498-01	3000100582
3000096946	30000969460081	D684-00-0296-01U	3000100721
3000096946	30000969460081	D684-00-0296-01	3000100722
3000096946	30000969460081	D684-00-0296-01	3000105063
3000096946	30000969460082	D684-00-0498-01	3000100582
3000096946	30000969460082	D684-00-0296-01U	3000100721
3000096946	30000969460082	D684-00-0296-01	3000100722
3000096946	30000969460082	D684-00-0296-01	3000105063
3000096946	30000969460083	D684-00-0498-01	3000100582
3000096946	30000969460083	D684-00-0296-01U	3000100721
3000096946	30000969460083	D684-00-0296-01	3000100722
3000096946	30000969460083	D684-00-0296-01	3000105063
3000096946	30000969460084	D684-00-0498-01	3000100582
3000096946	30000969460084	D684-00-0296-01U	3000100721
3000096946	30000969460084	D684-00-0296-01	3000100722
3000096946	30000969460084	D684-00-0296-01	3000105063
3000097412	30000974120282	D684-00-0469-01U	3000099928
3000097412	30000974120282	D684-00-0469-01	3000099929
3000097412	30000974120287	D684-00-0469-01U	3000099928
3000097412	30000974120287	D684-00-0469-01	3000099929
3000097412	30000974120288	D684-00-0469-01U	3000099928
3000097412	30000974120288	D684-00-0469-01	3000099929
3000097412	30000974120297	D684-00-0469-01U	3000099928

3000097412	30000974120297	D684-00-0469-01	3000099929
3000097412	30000974120298	D684-00-0469-01U	3000099928
3000097412	30000974120298	D684-00-0469-01	3000099929
3000097429	30000974290101	D684-00-0478-01	3000099930
3000097429	30000974290101	D684-00-0478-02U	3000099931
3000097429	30000974290101	D684-00-0478-01U	3000099932
3000097429	30000974290101	D684-00-0478-02	3000099933
3000097429	30000974290101	D684-00-0478-02	3000102246
3000097429	30000974290101	D684-00-0478-01U	3000103631
3000097429	30000974290102	D684-00-0478-01	3000099930
3000097429	30000974290102	D684-00-0478-02U	3000099931
3000097429	30000974290102	D684-00-0478-01U	3000099932
3000097429	30000974290102	D684-00-0478-02	3000099933
3000097429	30000974290102	D684-00-0478-02	3000102246
3000097429	30000974290102	D684-00-0478-01U	3000103631
3000097429	30000974290103	D684-00-0478-01	3000099930
3000097429	30000974290103	D684-00-0478-02U	3000099931
3000097429	30000974290103	D684-00-0478-01U	3000099932
3000097429	30000974290103	D684-00-0478-02	3000099933
3000097429	30000974290103	D684-00-0478-02	3000102246
3000097429	30000974290103	D684-00-0478-01U	3000103631
3000097429	30000974290104	D684-00-0478-01	3000099930
3000097429	30000974290104	D684-00-0478-02U	3000099931
3000097429	30000974290104	D684-00-0478-01U	3000099932
3000097429	30000974290104	D684-00-0478-02	3000099933
3000097429	30000974290104	D684-00-0478-02	3000102246
3000097429	30000974290104	D684-00-0478-01U	3000103631
3000097429	30000974290105	D684-00-0478-01	3000099930
3000097429	30000974290105	D684-00-0478-02U	3000099931
3000097429	30000974290105	D684-00-0478-01U	3000099932
3000097429	30000974290105	D684-00-0478-02	3000099933
3000097429	30000974290105	D684-00-0478-02	3000102246
3000097429	30000974290105	D684-00-0478-01U	3000103631
3000097441	30000974410101	D684-00-0561-01	3000104521
3000097441	30000974410101	D684-00-0561-01	3000115011
3000097441	30000974410110	D684-00-0561-01	3000104521
3000097441	30000974410110	D684-00-0561-01	3000115011
3000097441	30000974410111	D684-00-0561-01	3000104521
3000097441	30000974410111	D684-00-0561-01	3000115011
3000097441	30000974410112	D684-00-0561-01	3000104521
3000097441	30000974410112	D684-00-0561-01	3000115011
3000097441	30000974410114	D684-00-0561-01	3000104521
3000097441	30000974410114	D684-00-0561-01	3000115011
3000097443	30000974430046	D684-00-0559-01	3000103503
3000097443	30000974430047	D684-00-0559-01	3000103503
3000097443	30000974430048	D684-00-0559-01	3000103503
3000097443	30000974430049	D684-00-0559-01	3000103503
3000097443	30000974430050	D684-00-0559-01	3000103503
3000097445	30000974450105	D684-00-0480-07	3000100856
3000097445	30000974450105	D684-00-0480-08	3000101130
3000097445	30000974450105	D684-00-0480-02U	3000101787
3000097445	30000974450110	D684-00-0480-07	3000100856
3000097445	30000974450110	D684-00-0480-08	3000101130
3000097445	30000974450110	D684-00-0480-02U	3000101787
3000097445	30000974450114	D684-00-0480-07	3000100856
3000097445	30000974450114	D684-00-0480-08	3000101130
3000097445	30000974450114	D684-00	

3000097445	30000974450115	D684-00-0480-08	3000101130
3000097445	30000974450115	D684-00-0480-02U	3000101787
3000097445	30000974450125	D684-00-0480-07	3000100856
3000097445	30000974450125	D684-00-0480-08	3000101130
3000097445	30000974450125	D684-00-0480-02U	3000101787
3000097446	30000974460149	D684-00-0607	3000102408
3000097446	30000974460149	D684-00-0607	3000102409
3000097446	30000974460145	D684-00-0607	3000102408
3000097446	30000974460145	D684-00-0607	3000102409
3000097446	30000974460144	D684-00-0607	3000102408
3000097446	30000974460144	D684-00-0607	3000102409
3000097446	30000974460131	D684-00-0607	3000102408
3000097446	30000974460131	D684-00-0607	3000102409
3000097446	30000974460127	D684-00-0607	3000102408
3000097446	30000974460127	D684-00-0607	3000102409
3000097447	30000974470037	D684-00-0608	3000102410
3000097447	30000974470037	D684-00-0608	3000102411
3000097447	30000974470041	D684-00-0608	3000102410
3000097447	30000974470041	D684-00-0608	3000102411
3000097447	30000974470047	D684-00-0608	3000102410
3000097447	30000974470047	D684-00-0608	3000102411
3000097447	30000974470048	D684-00-0608	3000102410
3000097447	30000974470048	D684-00-0608	3000102411
3000097447	30000974470050	D684-00-0608	3000102410
3000097447	30000974470050	D684-00-0608	3000102411
3000097972	30000979720128	D684-00-0296-01U	3000102697
3000097972	30000979720128	D684-00-0296-02	3000102698
3000097972	30000979720128	D684-00-0296-01	3000102628
3000097972	30000979720128	D684-00-0296-02U	3000102699
3000097972	30000979720128	D684-00-0498-01	3000102905
3000097972	30000979720128	D684-00-0498-01	3000102905
3000097972	30000979720128	D684-00-0296-10	3000103022
3000097972	30000979720128	D684-00-0498-07	3000103039
3000097972	30000979720132	D684-00-0296-01U	3000102697
3000097972	30000979720132	D684-00-0296-02	3000102698
3000097972	30000979720132	D684-00-0296-01	3000102628
3000097972	30000979720132	D684-00-0296-02U	3000102699
3000097972	30000979720132	D684-00-0498-01	3000102905
3000097972	30000979720132	D684-00-0498-01	3000102905
3000097972	30000979720132	D684-00-0296-10	3000103022
3000097972	30000979720132	D684-00-0498-07	3000103039
3000097972	30000979720133	D684-00-0296-01U	3000102697
3000097972	30000979720133	D684-00-0296-02	3000102698
3000097972	30000979720133	D684-00-0296-01	3000102628
3000097972	30000979720133	D684-00-0296-02U	3000102699
3000097972	30000979720133	D684-00-0498-01	3000102905
3000097972	30000979720133	D684-00-0498-01	3000102905
3000097972	30000979720133	D684-00-0296-10	3000103022
3000097972	30000979720133	D684-00-0498-07	3000103039
3000097972	30000979720135	D684-00-0296-01U	3000102697
3000097972	30000979720135	D684-00-0296-02	3000102698
3000097972	30000979720135	D684-00-0296-01	3000102628
3000097972	30000979720135	D684-00-0296-02U	3000102699
3000097972	30000979720135	D684-00-0498-01	3000102905
3000097972	30000979720135	D684-00-0498-01	3000102905
3000097972	30000979720135	D684-00-0296-10	3000103022
3000097972	30000979720135	D684-00-0498-07	3000103039
3000097972	30000979720138	D684-00-0296-01U	3000102697

3000097972	30000979720138	D684-00-0296-02	3000102698
3000097972	30000979720138	D684-00-0296-01	3000102628
3000097972	30000979720138	D684-00-0296-02U	3000102699
3000097972	30000979720138	D684-00-0498-01	3000102905
3000097972	30000979720138	D684-00-0498-01	3000102905
3000097972	30000979720138	D684-00-0296-10	3000103022
3000097972	30000979720138	D684-00-0498-07	3000103039
3000097973	30000979730051	D684-00-0294-01	3000100873
3000097973	30000979730051	D684-00-0294-01U	3000102113
3000097973	30000979730051	D684-00-0294-01	3000102114
3000097973	30000979730052	D684-00-0294-01	3000100873
3000097973	30000979730052	D684-00-0294-01U	3000102113
3000097973	30000979730052	D684-00-0294-01	3000102114
3000097973	30000979730057	D684-00-0294-01	3000100873
3000097973	30000979730057	D684-00-0294-01U	3000102113
3000097973	30000979730057	D684-00-0294-01	3000102114
3000097973	30000979730061	D684-00-0294-01	3000100873
3000097973	30000979730061	D684-00-0294-01U	3000102113
3000097973	30000979730061	D684-00-0294-01	3000102114
3000097973	30000979730067	D684-00-0294-01	3000100873
3000097973	30000979730067	D684-00-0294-01U	3000102113
3000097973	30000979730067	D684-00-0294-01	3000102114
3000097974	30000979740055	D684-00-0545-01	3000103498
3000097974	30000979740055	D684-00-0545-01	3000104616
3000097974	30000979740059	D684-00-0545-01	3000103498
3000097974	30000979740059	D684-00-0545-01	3000104616
3000097974	30000979740064	D684-00-0545-01	3000103498
3000097974	30000979740064	D684-00-0545-01	3000104616
3000097974	30000979740070	D684-00-0545-01	3000103498
3000097974	30000979740070	D684-00-0545-01	3000104616
3000097974	30000979740075	D684-00-0545-01	3000103498
3000097974	30000979740075	D684-00-0545-01	3000104616
3000097975	30000979750044	D684-00-0576-01U	3000100588
3000097975	30000979750044	D684-00-0576-01U	3000100589
3000097975	30000979750050	D684-00-0576-01U	3000100588
3000097975	30000979750050	D684-00-0576-01U	3000100589
3000097975	30000979750051	D684-00-0576-01U	3000100588
3000097975	30000979750051	D684-00-0576-01U	3000100589
3000097975	30000979750055	D684-00-0576-01U	3000100588
3000097975	30000979750055	D684-00-0576-01U	3000100589
3000097975	30000979750060	D684-00-0576-01U	3000100588
3000097975	30000979750060	D684-00-0576-01U	3000100589
3000097998	30000979980005	D684-00-0479-08	3000101521
3000097998	30000979980005	D684-00-0479-07	3000101585
3000097998	30000979980005	D684-00-0479-07	3000101586
3000097998	30000979980005	D684-00-0479-01U	3000101888
3000097998	30000979980009	D684-00-0479-08	3000101521
3000097998	30000979980009	D684-00-0479-07	3000101585
3000097998	30000979980009	D684-00-0479-07	3000101586
3000097998	30000979980009	D684-00-0479-01U	3000101888
3000097998	30000979980021	D684-00-0479-08	3000101521
3000097998	30000979980021	D684-00-0479-07	3000101585
3000097998	30000979980021	D684-00-0479-07	3000101586
3000097998	30000979980021	D684-00-0479-01U	3000101888
3000097998	30000979980023	D684-00-0479-08	3000101521
3000097998	30000979980023	D684-00-0479-07	3000101585
3000097998	30000979980023	D684-00-0479-07	3000101586
3000097998	30000979980023	D684-00-0479-01U	3000101888

3000097998	30000979980024	D684-00-0479-08	3000101521
3000097998	30000979980024	D684-00-0479-07	3000101585
3000097998	30000979980024	D684-00-0479-07	3000101586
3000097998	30000979980024	D684-00-0479-01U	3000101888
3000097999	30000979990590	D684-00-0480-01	3000101784
3000097999	30000979990590	D684-00-0480-01U	3000101785
3000097999	30000979990595	D684-00-0480-01	3000101784
3000097999	30000979990595	D684-00-0480-01U	3000101785
3000097999	30000979990601	D684-00-0480-01	3000101784
3000097999	30000979990601	D684-00-0480-01U	3000101785
3000097999	30000979990602	D684-00-0480-01	3000101784
3000097999	30000979990602	D684-00-0480-01U	3000101785
3000097999	30000979990604	D684-00-0480-01	3000101784
3000097999	30000979990604	D684-00-0480-01U	3000101785
3000098601	30000986010054	D684-00-0470-01	3000101522
3000098601	30000986010054	D684-00-0470-01U	3000101523
3000098601	30000986010056	D684-00-0470-01	3000101522
3000098601	30000986010056	D684-00-0470-01U	3000101523
3000098601	30000986010069	D684-00-0470-01	3000101522
3000098601	30000986010069	D684-00-0470-01U	3000101523
3000098601	30000986010070	D684-00-0470-01	3000101522
3000098601	30000986010070	D684-00-0470-01U	3000101523
3000098601	30000986010073	D684-00-0470-01	3000101522
3000098601	30000986010073	D684-00-0470-01U	3000101523
3000098605	30000986050052	D684-00-0295-01	3000102116
3000098605	30000986050052	D684-00-0295-05	3000102117
3000098605	30000986050052	D684-00-0295-01U	3000103157
3000098605	30000986050052	D684-00-0295-02	3000104736
3000098605	30000986050056	D684-00-0295-01	3000102116
3000098605	30000986050056	D684-00-0295-05	3000102117
3000098605	30000986050056	D684-00-0295-01U	3000103157
3000098605	30000986050056	D684-00-0295-02	3000104736
3000098605	30000986050065	D684-00-0295-01	3000102116
3000098605	30000986050065	D684-00-0295-05	3000102117
3000098605	30000986050065	D684-00-0295-01U	3000103157
3000098605	30000986050065	D684-00-0295-02	3000104736
3000098605	30000986050070	D684-00-0295-01	3000102116
3000098605	30000986050070	D684-00-0295-05	3000102117
3000098605	30000986050070	D684-00-0295-01U	3000103157
3000098605	30000986050070	D684-00-0295-02	3000104736
3000098605	30000986050072	D684-00-0295-01	3000102116
3000098605	30000986050072	D684-00-0295-05	3000102117
3000098605	30000986050072	D684-00-0295-01U	3000103157
3000098605	30000986050072	D684-00-0295-02	3000104736
3000098608	30000986080115	D684-00-0480-02	3000101786
3000098608	30000986080115	D684-00-0480-01C	3000102536
3000098608	30000986080116	D684-00-0480-02	3000101786
3000098608	30000986080116	D684-00-0480-01C	3000102536
3000098608	30000986080117	D684-00-0480-02	3000101786
3000098608	30000986080117	D684-00-0480-01C	3000102536
3000098608	30000986080118	D684-00-0480-02	3000101786
3000098608	30000986080118	D684-00-0480-01C	3000102536
3000098608	30000986080119	D684-00-0480-02	3000101786
3000098608	30000986080119	D684-00-0480-01C	3000102536
3000098612	30000986120161	D684-00-0568-01U	3000102215
3000098612	30000986120165	D684-00-0568-01U	3000102215
3000098612	30000986120169	D684-00-0568-01U	3000102215
3000098612	30000986120170	D684-00-0568-01U	3000102215

3000098612	30000986120177	D684-00-0568-01U	3000102215
3000099024	30000990240152	D684-00-0479-01U	3000101932
3000099024	30000990240152	D684-00-0479-01	3000101934
3000099024	30000990240152	D684-00-0479-01C	3000102537
3000099024	30000990240152	D684-DC-0474	0
3000099024	30000990240153	D684-00-0479-01U	3000101932
3000099024	30000990240153	D684-00-0479-01	3000101934
3000099024	30000990240153	D684-00-0479-01C	3000102537
3000099024	30000990240153	D684-DC-0474	0
3000099024	30000990240159	D684-00-0479-01U	3000101932
3000099024	30000990240159	D684-00-0479-01	3000101934
3000099024	30000990240159	D684-00-0479-01C	3000102537
3000099024	30000990240159	D684-DC-0474	0
3000099024	30000990240161	D684-00-0479-01U	3000101932
3000099024	30000990240161	D684-00-0479-01	3000101934
3000099024	30000990240161	D684-00-0479-01C	3000102537
3000099024	30000990240161	D684-DC-0474	0
3000099024	30000990240162	D684-00-0479-01U	3000101932
3000099024	30000990240162	D684-00-0479-01	3000101934
3000099024	30000990240162	D684-00-0479-01C	3000102537
3000099024	30000990240162	D684-DC-0474	0
3000099048	30000990480127	D684-00-0576-01U	3000101988
3000099048	30000990480127	D684-00-0576-01	3000101990
3000099048	30000990480128	D684-00-0576-01U	3000101988
3000099048	30000990480128	D684-00-0576-01	3000101990
3000099048	30000990480134	D684-00-0576-01U	3000101988
3000099048	30000990480134	D684-00-0576-01	3000101990
3000099048	30000990480138	D684-00-0576-01U	3000101988
3000099048	30000990480138	D684-00-0576-01	3000101990
3000099048	30000990480144	D684-00-0576-01U	3000101988
3000099048	30000990480144	D684-00-0576-01	3000101990
3000099064	30000990640051	D684-00-0480-07	3000103052
3000099064	30000990640052	D684-00-0480-07	3000103052
3000099064	30000990640053	D684-00-0480-07	3000103052
3000099064	30000990640054	D684-00-0480-07	3000103052
3000099064	30000990640055	D684-00-0480-07	3000103052
3000099571	30000995710230	D684-00-0559-01	3000103504
3000099571	30000995710230	D684-00-0559-01	3000103505
3000099571	30000995710230	D684-00-0559-01	3000103506
3000099571	30000995710230	D684-00-0559-01	3000110580
3000099571	30000995710230	D684-00-0559-01	3000110581
3000099571	30000995710234	D684-00-0559-01	3000103504
3000099571	30000995710234	D684-00-0559-01	3000103505
3000099571	30000995710234	D684-00-0559-01	3000103506
3000099571	30000995710234	D684-00-0559-01	3000110580
3000099571	30000995710234	D684-00-0559-01	3000110581
3000099571	30000995710236	D684-00-0559-01	3000103504
3000099571	30000995710236	D684-00-0559-01	3000103505
3000099571	30000995710236	D684-00-0559-01	3000103506
3000099571	30000995710236	D684-00-0559-01	3000110580
3000099571	30000995710236	D684-00-0559-01	3000110581
3000099571	30000995710242	D684-00-0559-01	3000103504
3000099571	30000995710242	D684-00-0559-01	3000103505
3000099571	30000995710242	D684-00-0559-01	3000103506
3000099571	30000995710242	D684-00-0559-01	3000110580
3000099571	30000995710242	D684-00-0559-01	3000110581
3000099571	30000995710244	D684-00-0559-01	3000103504
3000099571	30000995710244	D684-00-0559-01	3000103505

3000099571	30000995710244	D684-00-0559-01	3000103506
3000099571	30000995710244	D684-00-0559-01	3000110580
3000099571	30000995710244	D684-00-0559-01	3000110581
3000099576	30000995760059	D684-00-0294-01U	3000106149
3000099576	30000995760065	D684-00-0294-01U	3000106149
3000099576	30000995760070	D684-00-0294-01U	3000106149
3000099576	30000995760071	D684-00-0294-01U	3000106149
3000099576	30000995760072	D684-00-0294-01U	3000106149
3000099577	30000995770296	D684-00-0568-01U	3000102248
3000099577	30000995770296	D684-00-0568-01U	3000102249
3000099577	30000995770296	D684-00-0568-01U	3000102249
3000099577	30000995770296	D684-00-0568-09	3000102329
3000099577	30000995770297	D684-00-0568-01U	3000102248
3000099577	30000995770297	D684-00-0568-01U	3000102249
3000099577	30000995770297	D684-00-0568-01U	3000102249
3000099577	30000995770297	D684-00-0568-09	3000102329
3000099577	30000995770298	D684-00-0568-01U	3000102248
3000099577	30000995770298	D684-00-0568-01U	3000102249
3000099577	30000995770298	D684-00-0568-01U	3000102249
3000099577	30000995770298	D684-00-0568-09	3000102329
3000099577	30000995770299	D684-00-0568-01U	3000102248
3000099577	30000995770299	D684-00-0568-01U	3000102249
3000099577	30000995770299	D684-00-0568-01U	3000102249
3000099577	30000995770299	D684-00-0568-09	3000102329
3000099577	30000995770302	D684-00-0568-01U	3000102248
3000099577	30000995770302	D684-00-0568-01U	3000102249
3000099577	30000995770302	D684-00-0568-01U	3000102249
3000099577	30000995770302	D684-00-0568-09	3000102329
3000099578	30000995780226	D684-00-0576-01U	3000102703
3000099578	30000995780226	D684-00-0576-01U	3000102704
3000099578	30000995780226	D684-00-0576-01	3000102705
3000099578	30000995780230	D684-00-0576-01U	3000102703
3000099578	30000995780230	D684-00-0576-01U	3000102704
3000099578	30000995780230	D684-00-0576-01	3000102705
3000099578	30000995780232	D684-00-0576-01U	3000102703
3000099578	30000995780232	D684-00-0576-01U	3000102704
3000099578	30000995780232	D684-00-0576-01	3000102705
3000099578	30000995780236	D684-00-0576-01U	3000102703
3000099578	30000995780236	D684-00-0576-01U	3000102704
3000099578	30000995780236	D684-00-0576-01	3000102705
3000099578	30000995780244	D684-00-0576-01U	3000102703
3000099578	30000995780244	D684-00-0576-01U	3000102704
3000099578	30000995780244	D684-00-0576-01	3000102705
3000099579	30000995790045	D684-00-0295-10	3000103156
3000099579	30000995790045	D684-00-0295-01U	3000103157
3000099579	30000995790045	D684-00-0295-05	3000103158
3000099579	30000995790045	D684-00-0295-02	3000104737
3000099579	30000995790050	D684-00-0295-10	3000103156
3000099579	30000995790050	D684-00-0295-01U	3000103157
3000099579	30000995790050	D684-00-0295-05	3000103158
3000099579	30000995790050	D684-00-0295-02	3000104737
3000099997	30000999970321	D684-00-0480-07	3000103053
3000099997	30000999970321	D684-00-0480-10	3000103055
3000099997	30000999970321	D684-00-0480-10	3000103916
3000099997	30000999970321	D684-00-0480-01C	3000103919
3000099997	30000999970321	D684-00-0480-01U	3000103920
3000099997	30000999970322	D684-00-0480-07	3000103053
3000099997	30000999970322	D684-00-0480-10	3000103055

3000099997	30000999970322	D684-00-0480-10	3000103916
3000099997	30000999970322	D684-00-0480-01C	3000103919
3000099997	30000999970322	D684-00-0480-01U	3000103920
3000099997	30000999970323	D684-00-0480-07	3000103053
3000099997	30000999970323	D684-00-0480-10	3000103055
3000099997	30000999970323	D684-00-0480-10	3000103916
3000099997	30000999970323	D684-00-0480-01C	3000103919
3000099997	30000999970323	D684-00-0480-01U	3000103920
3000099997	30000999970324	D684-00-0480-07	3000103053
3000099997	30000999970324	D684-00-0480-10	3000103055
3000099997	30000999970324	D684-00-0480-10	3000103916
3000099997	30000999970324	D684-00-0480-01C	3000103919
3000099997	30000999970324	D684-00-0480-01U	3000103920
3000099997	30000999970325	D684-00-0480-07	3000103053
3000099997	30000999970325	D684-00-0480-10	3000103055
3000099997	30000999970325	D684-00-0480-10	3000103916
3000099997	30000999970325	D684-00-0480-01C	3000103919
3000099997	30000999970325	D684-00-0480-01U	3000103920
3000100003	30001000030250	D684-00-0607	3000103469
3000100003	30001000030255	D684-00-0607	3000103469
3000100003	30001000030256	D684-00-0607	3000103469
3000100003	30001000030260	D684-00-0607	3000103469
3000100003	30001000030261	D684-00-0607	3000103469
3000100005	30001000050268	D684-00-0608	3000103470
3000100005	30001000050272	D684-00-0608	3000103470
3000100005	30001000050248	D684-00-0608	3000103470
3000100005	30001000050257	D684-00-0608	3000103470
3000100005	30001000050250	D684-00-0608	3000103470
3000100008	30001000080061	D684-00-0479-07	3000103054
3000100008	30001000080061	D684-00-0479-10	3000103155
3000100008	30001000080061	D684-00-0479-10	3000103917
3000100008	30001000080061	D684-00-0479-01C	3000103918
3000100008	30001000080062	D684-00-0479-07	3000103054
3000100008	30001000080062	D684-00-0479-10	3000103155
3000100008	30001000080062	D684-00-0479-10	3000103917
3000100008	30001000080062	D684-00-0479-01C	3000103918
3000100008	30001000080063	D684-00-0479-07	3000103054
3000100008	30001000080063	D684-00-0479-10	3000103155
3000100008	30001000080063	D684-00-0479-10	3000103917
3000100008	30001000080063	D684-00-0479-01C	3000103918
3000100008	30001000080064	D684-00-0479-07	3000103054
3000100008	30001000080064	D684-00-0479-10	3000103155
3000100008	30001000080064	D684-00-0479-10	3000103917
3000100008	30001000080064	D684-00-0479-01C	3000103918
3000100008	30001000080065	D684-00-0479-07	3000103054
3000100008	30001000080065	D684-00-0479-10	3000103155
3000100008	30001000080065	D684-00-0479-10	3000103917
3000100008	30001000080065	D684-00-0479-01C	3000103918
3000100010	30001000100081	D684-00-0469-01	3000102977
3000100010	30001000100081	D684-00-0469-01U	3000102978
3000100010	30001000100081	D684-00-0469-01C	3000103213
3000100010	30001000100081	D684-00-0469-09	3000103214
3000100010	30001000100086	D684-00-0469-01	3000102977
3000100010	30001000100086	D684-00-0469-01U	3000102978
3000100010	30001000100086	D684-00-0469-01C	3000103213
3000100010	30001000100086	D684-00-0469-09	3000103214
3000100010	30001000100096	D684-00-0469-01	3000102977

3000100010	30001000100096	D684-00-0469-01C	3000103213
3000100010	30001000100096	D684-00-0469-09	3000103214
3000100010	30001000100100	D684-00-0469-01	3000102977
3000100010	30001000100100	D684-00-0469-01U	3000102978
3000100010	30001000100100	D684-00-0469-01C	3000103213
3000100010	30001000100100	D684-00-0469-09	3000103214
3000100010	30001000100211	D684-00-0469-01	3000102977
3000100010	30001000100211	D684-00-0469-01U	3000102978
3000100010	30001000100211	D684-00-0469-01C	3000103213
3000100010	30001000100211	D684-00-0469-09	3000103214
3000100338	30001003380474	D684-00-0480-01	3000103915
3000100338	30001003380474	D684-00-0480-01	3000103915
3000100338	30001003380474	D684-00-0480-02U	3000105897
3000100338	30001003380479	D684-00-0480-01	3000103915
3000100338	30001003380479	D684-00-0480-01	3000103915
3000100338	30001003380479	D684-00-0480-02U	3000105897
3000100338	30001003380482	D684-00-0480-01	3000103915
3000100338	30001003380482	D684-00-0480-01	3000103915
3000100338	30001003380482	D684-00-0480-02U	3000105897
3000100338	30001003380490	D684-00-0480-01	3000103915
3000100338	30001003380490	D684-00-0480-01	3000103915
3000100338	30001003380490	D684-00-0480-02U	3000105897
3000100338	30001003380491	D684-00-0480-01	3000103915
3000100338	30001003380491	D684-00-0480-01	3000103915
3000100338	30001003380491	D684-00-0480-02U	3000105897
3000100339	30001003390392	D684-00-0470-01	3000103209
3000100339	30001003390392	D684-00-0470-01U	3000103210
3000100339	30001003390392	D684-00-0470-09	3000103211
3000100339	30001003390392	D684-00-0470-01C	3000103212
3000100339	30001003390394	D684-00-0470-01	3000103209
3000100339	30001003390394	D684-00-0470-01U	3000103210
3000100339	30001003390394	D684-00-0470-09	3000103211
3000100339	30001003390394	D684-00-0470-01C	3000103212
3000100339	30001003390397	D684-00-0470-01	3000103209
3000100339	30001003390397	D684-00-0470-01U	3000103210
3000100339	30001003390397	D684-00-0470-09	3000103211
3000100339	30001003390397	D684-00-0470-01C	3000103212
3000100339	30001003390399	D684-00-0470-01	3000103209
3000100339	30001003390399	D684-00-0470-01U	3000103210
3000100339	30001003390399	D684-00-0470-09	3000103211
3000100339	30001003390399	D684-00-0470-01C	3000103212
3000100339	30001003390407	D684-00-0470-01	3000103209
3000100339	30001003390407	D684-00-0470-01U	3000103210
3000100339	30001003390407	D684-00-0470-09	3000103211
3000100339	30001003390407	D684-00-0470-01C	3000103212
3000100340	30001003400091	D684-00-0296-02	3000105060
3000100340	30001003400091	D684-00-0498-01	3000105061
3000100340	30001003400091	D684-00-0296-01U	3000105062
3000100340	30001003400091	D684-00-0296-01	3000105064
3000100340	30001003400091	D684-00-0296-02	3000106965
3000100340	30001003400092	D684-00-0296-02	3000105060
3000100340	30001003400092	D684-00-0498-01	3000105061
3000100340	30001003400092	D684-00-0296-01U	3000105062
3000100340	30001003400092	D684-00-0296-01	3000105064
3000100340	30001003400092	D684-00-0296-02	3000106965
3000100340	30001003400093	D684-00-0296-02	3000105060
3000100340	30001003400093	D684-00-0498-01	3000105061
3000100340	30001003400093	D684-00-0296-01U	3000105062

3000100340	30001003400093	D684-00-0296-01	3000105064
3000100340	30001003400093	D684-00-0296-02	3000106965
3000100340	30001003400094	D684-00-0296-02	3000105060
3000100340	30001003400094	D684-00-0498-01	3000105061
3000100340	30001003400094	D684-00-0296-01U	3000105062
3000100340	30001003400094	D684-00-0296-01	3000105064
3000100340	30001003400094	D684-00-0296-02	3000106965
3000100340	30001003400095	D684-00-0296-02	3000105060
3000100340	30001003400095	D684-00-0498-01	3000105061
3000100340	30001003400095	D684-00-0296-01U	3000105062
3000100340	30001003400095	D684-00-0296-01	3000105064
3000100340	30001003400095	D684-00-0296-02	3000106965
3000101012	30001010120025	D684-00-0295-01U	3000104738
3000101012	30001010120025	D684-00-0295-05	3000104739
3000101012	30001010120025	D684-00-0295-01U	3000106751
3000101012	30001010120026	D684-00-0295-01U	3000104738
3000101012	30001010120026	D684-00-0295-05	3000104739
3000101012	30001010120026	D684-00-0295-01U	3000106751
3000101012	30001010120027	D684-00-0295-01U	3000104738
3000101012	30001010120027	D684-00-0295-05	3000104739
3000101012	30001010120027	D684-00-0295-01U	3000106751
3000101012	30001010120028	D684-00-0295-01U	3000104738
3000101012	30001010120028	D684-00-0295-05	3000104739
3000101012	30001010120028	D684-00-0295-01U	3000106751
3000101012	30001010120029	D684-00-0295-01U	3000104738
3000101012	30001010120029	D684-00-0295-05	3000104739
3000101012	30001010120029	D684-00-0295-01U	3000106751
3000101068	30001010680524	D684-00-0576-09	3000103206
3000101068	30001010680524	D684-00-0576-01	3000103207
3000101068	30001010680524	D684-00-0576-01U	3000103208
3000101068	30001010680524	D684-00-0576-01	3000104387
3000101068	30001010680524	D684-00-0576-09	3000104388
3000101068	30001010680525	D684-00-0576-09	3000103206
3000101068	30001010680525	D684-00-0576-01	3000103207
3000101068	30001010680525	D684-00-0576-01U	3000103208
3000101068	30001010680525	D684-00-0576-01	3000104387
3000101068	30001010680525	D684-00-0576-09	3000104388
3000101068	30001010680527	D684-00-0576-09	3000103206
3000101068	30001010680527	D684-00-0576-01	3000103207
3000101068	30001010680527	D684-00-0576-01U	3000103208
3000101068	30001010680527	D684-00-0576-01	3000104387
3000101068	30001010680527	D684-00-0576-09	3000104388
3000101068	30001010680535	D684-00-0576-09	3000103206
3000101068	30001010680535	D684-00-0576-01	3000103207
3000101068	30001010680535	D684-00-0576-01U	3000103208
3000101068	30001010680535	D684-00-0576-01	3000104387
3000101068	30001010680535	D684-00-0576-09	3000104388
3000101068	30001010680541	D684-00-0576-09	3000103206
3000101068	30001010680541	D684-00-0576-01	3000103207
3000101068	30001010680541	D684-00-0576-01U	3000103208
3000101068	30001010680541	D684-00-0576-01	3000104387
3000101068	30001010680541	D684-00-0576-09	3000104388
3000101069	30001010690031	D684-00-0478-01	3000103632
3000101069	30001010690031	D684-00-0478-01	3000112148
3000101069	30001010690032	D684-00-0478-01	3000103632
3000101069	30001010690032	D684-00-0478-01	3000112148
3000101069	30001010690033	D684-00-	

3000101069	30001010690034	D684-00-0478-01	3000103632
3000101069	30001010690034	D684-00-0478-01	3000112148
3000101069	30001010690035	D684-00-0478-01	3000103632
3000101069	30001010690035	D684-00-0478-01	3000112148
3000101090	30001010900011	D684-00-0480-10	3000105384
3000101090	30001010900011	D684-00-0480-01U	3000105895
3000101090	30001010900012	D684-00-0480-10	3000105384
3000101090	30001010900012	D684-00-0480-01U	3000105895
3000101090	30001010900013	D684-00-0480-10	3000105384
3000101090	30001010900013	D684-00-0480-01U	3000105895
3000101090	30001010900014	D684-00-0480-10	3000105384
3000101090	30001010900014	D684-00-0480-01U	3000105895
3000101090	30001010900015	D684-00-0480-10	3000105384
3000101090	30001010900015	D684-00-0480-01U	3000105895
3000101092	30001010920351	D684-00-0479-10	3000105386
3000101092	30001010920351	D684-00-0479-01U	3000105387
3000101092	30001010920351	D684-00-0479-01C	3000105531
3000101092	30001010920352	D684-00-0479-10	3000105386
3000101092	30001010920352	D684-00-0479-01U	3000105387
3000101092	30001010920352	D684-00-0479-01C	3000105531
3000101092	30001010920353	D684-00-0479-10	3000105386
3000101092	30001010920353	D684-00-0479-01U	3000105387
3000101092	30001010920353	D684-00-0479-01C	3000105531
3000101092	30001010920354	D684-00-0479-10	3000105386
3000101092	30001010920354	D684-00-0479-01U	3000105387
3000101092	30001010920354	D684-00-0479-01C	3000105531
3000101092	30001010920355	D684-00-0479-10	3000105386
3000101092	30001010920355	D684-00-0479-01U	3000105387
3000101092	30001010920355	D684-00-0479-01C	3000105531
3000101093	30001010930005	D684-00-0568-09	3000103404
3000101093	30001010930005	D684-00-0568-01U	3000103405
3000101093	30001010930005	D684-00-0568-01U	3000103406
3000101093	30001010930010	D684-00-0568-09	3000103404
3000101093	30001010930010	D684-00-0568-01U	3000103405
3000101093	30001010930010	D684-00-0568-01U	3000103406
3000101093	30001010930013	D684-00-0568-09	3000103404
3000101093	30001010930013	D684-00-0568-01U	3000103405
3000101093	30001010930013	D684-00-0568-01U	3000103406
3000101093	30001010930014	D684-00-0568-09	3000103404
3000101093	30001010930014	D684-00-0568-01U	3000103405
3000101093	30001010930014	D684-00-0568-01U	3000103406
3000101093	30001010930019	D684-00-0568-09	3000103404
3000101093	30001010930019	D684-00-0568-01U	3000103405
3000101093	30001010930019	D684-00-0568-01U	3000103406
3000101468	30001014680201	D684-00-0607	3000104520
3000101468	30001014680201	D684-00-0607	3000104615
3000101468	30001014680212	D684-00-0607	3000104520
3000101468	30001014680212	D684-00-0607	3000104615
3000101468	30001014680218	D684-00-0607	3000104520
3000101468	30001014680218	D684-00-0607	3000104615
3000101468	30001014680220	D684-00-0607	3000104520
3000101468	30001014680220	D684-00-0607	3000104615
3000101468	30001014680222	D684-00-0607	3000104520
3000101468	30001014680222	D684-00-0607	3000104615
3000101469	30001014690121	D684-00-0559-01	3000110582
3000101469	30001014690122	D684-00-0559-01	3000110582
3000101469	30001014690123	D684-00-0559-01	3000110582
3000101469	30001014690124	D684-00-0559-01	3000110582

3000101469	30001014690125	D684-00-0559-01	3000110582
3000101471	30001014710051	D684-00-0608	3000104519
3000101471	30001014710053	D684-00-0608	3000104519
3000101471	30001014710054	D684-00-0608	3000104519
3000101471	30001014710057	D684-00-0608	3000104519
3000101471	30001014710072	D684-00-0608	3000104519
3000101655	30001016550027	D684-00-0545-01	3000104617
3000101655	30001016550027	D684-00-0545-01	3000112275
3000101655	30001016550027	D684-00-0545-01	3000113067
3000101655	30001016550030	D684-00-0545-01	3000104617
3000101655	30001016550030	D684-00-0545-01	3000112275
3000101655	30001016550030	D684-00-0545-01	3000113067
3000101655	30001016550032	D684-00-0545-01	3000104617
3000101655	30001016550032	D684-00-0545-01	3000112275
3000101655	30001016550032	D684-00-0545-01	3000113067
3000101655	30001016550038	D684-00-0545-01	3000104617
3000101655	30001016550038	D684-00-0545-01	3000112275
3000101655	30001016550038	D684-00-0545-01	3000113067
3000101655	30001016550042	D684-00-0545-01	3000104617
3000101655	30001016550042	D684-00-0545-01	3000112275
3000101655	30001016550042	D684-00-0545-01	3000113067
3000101660	30001016600556	D684-00-0576-01U	3000104808
3000101660	30001016600556	D684-00-0576-01	3000106150
3000101660	30001016600559	D684-00-0576-01U	3000104808
3000101660	30001016600559	D684-00-0576-01	3000106150
3000101660	30001016600568	D684-00-0576-01U	3000104808
3000101660	30001016600568	D684-00-0576-01	3000106150
3000101660	30001016600569	D684-00-0576-01U	3000104808
3000101660	30001016600569	D684-00-0576-01	3000106150
3000101660	30001016600573	D684-00-0576-01U	3000104808
3000101660	30001016600573	D684-00-0576-01	3000106150
3000102059	30001020590131	D684-00-0479-01	3000105734
3000102059	30001020590131	D684-00-0479-02	3000105735
3000102059	30001020590131	D684-00-0479-01C	3000106633
3000102059	30001020590132	D684-00-0479-01	3000105734
3000102059	30001020590132	D684-00-0479-02	3000105735
3000102059	30001020590132	D684-00-0479-01C	3000106633
3000102059	30001020590133	D684-00-0479-01	3000105734
3000102059	30001020590133	D684-00-0479-02	3000105735
3000102059	30001020590133	D684-00-0479-01C	3000106633
3000102059	30001020590134	D684-00-0479-01	3000105734
3000102059	30001020590134	D684-00-0479-02	3000105735
3000102059	30001020590134	D684-00-0479-01C	3000106633
3000102059	30001020590135	D684-00-0479-01	3000105734
3000102059	30001020590135	D684-00-0479-02	3000105735
3000102059	30001020590135	D684-00-0479-01C	3000106633
3000102092	30001020920171	D684-00-0568-01	3000105260
3000102092	30001020920171	D684-00-0568-01U	3000105261
3000102092	30001020920171	D684-00-0568-01U	3000105262
3000102092	30001020920174	D684-00-0568-01	3000105260
3000102092	30001020920174	D684-00-0568-01U	3000105261
3000102092	30001020920174	D684-00-0568-01U	3000105262
3000102092	30001020920178	D684-00-0568-01	3000105260
3000102092	30001020920178	D684-00-0568-01U	3000105261
3000102092	30001020920178	D684-00-0568-01U	3000105262
3000102092	30001020920184	D684-00-0568-01	3000105260
3000102092	30001020920184	D684-00-0568-01U	3000105261
3000102092	30001020920184	D684-00-0568-01U	3000105262

3000102092	30001020920185	D684-00-0568-01	3000105260
3000102092	30001020920185	D684-00-0568-01U	3000105261
3000102092	30001020920185	D684-00-0568-01U	3000105262
3000102102	30001021020205	D684-00-0480-01	3000105898
3000102102	30001021020205	D684-00-0480-01C	3000106447
3000102102	30001021020205	D684-00-0480-10	3000106448
3000102102	30001021020205	D684-00-0480-01U	3000106449
3000102102	30001021020206	D684-00-0480-01	3000105898
3000102102	30001021020206	D684-00-0480-01C	3000106447
3000102102	30001021020206	D684-00-0480-10	3000106448
3000102102	30001021020206	D684-00-0480-01U	3000106449
3000102112	30001021120458	D684-00-0607	3000106012
3000102112	30001021120458	D684-00-0607	3000107758
3000102112	30001021120460	D684-00-0607	3000106012
3000102112	30001021120460	D684-00-0607	3000107758
3000102112	30001021120462	D684-00-0607	3000106012
3000102112	30001021120462	D684-00-0607	3000107758
3000102112	30001021120463	D684-00-0607	3000106012
3000102112	30001021120463	D684-00-0607	3000107758
3000102112	30001021120467	D684-00-0607	3000106012
3000102112	30001021120467	D684-00-0607	3000107758
3000102610	30001026100053	D684-00-0295-01	3000106750
3000102610	30001026100053	D684-00-0295-02	3000106752
3000102610	30001026100053	D684-00-0295-02U	3000106753
3000102610	30001026100053	D684-00-0295-05	3000106898
3000102610	30001026100058	D684-00-0295-01	3000106750
3000102610	30001026100058	D684-00-0295-02	3000106752
3000102610	30001026100058	D684-00-0295-02U	3000106753
3000102610	30001026100058	D684-00-0295-05	3000106898
3000102610	30001026100066	D684-00-0295-01	3000106750
3000102610	30001026100066	D684-00-0295-02	3000106752
3000102610	30001026100066	D684-00-0295-02U	3000106753
3000102610	30001026100066	D684-00-0295-05	3000106898
3000102611	30001026110275	D684-00-0608	3000106011
3000102611	30001026110277	D684-00-0608	3000106011
3000102611	30001026110278	D684-00-0608	3000106011
3000102611	30001026110279	D684-00-0608	3000106011
3000102611	30001026110280	D684-00-0608	3000106011
3000102613	30001026130081	D684-00-0480-02	3000106316
3000102613	30001026130081	D684-00-0480-02	3000106316
3000102613	30001026130081	D684-00-0480-01	3000106317
3000102613	30001026130081	D684-00-0480-01U	3000106450
3000102613	30001026130082	D684-00-0480-02	3000106316
3000102613	30001026130082	D684-00-0480-02	3000106316
3000102613	30001026130082	D684-00-0480-01	3000106317
3000102613	30001026130082	D684-00-0480-01U	3000106450
3000102613	30001026130083	D684-00-0480-02	3000106316
3000102613	30001026130083	D684-00-0480-02	3000106316
3000102613	30001026130083	D684-00-0480-01	3000106317
3000102613	30001026130083	D684-00-0480-01U	3000106450
3000102613	30001026130084	D684-00-0480-02	3000106316
3000102613	30001026130084	D684-00-0480-02	3000106316
3000102613	30001026130084	D684-00-0480-01	3000106317
3000102613	30001026130084	D684-00-0480-01U	3000106450
3000102613	30001026130085	D684-00-0480-02	3000106316
3000102613	30001026130085	D684-00-0480-02	3000106316
3000102613	30001026130085	D684-00-0480-01	3000106317
3000102613	30001026130085	D684-00-0480-01U	3000106450

[illegible]

3000103830	30001038300519	D684-00-0576-01U	3000106152
3000103830	30001038300519	D684-00-0576-09	3000108351
3000103832	30001038320854	D684-00-0470-01	3000106823
3000103832	30001038320854	D684-00-0470-09	3000106824
3000103832	30001038320854	D684-00-0470-01U	3000106825
3000103832	30001038320863	D684-00-0470-01	3000106823
3000103832	30001038320863	D684-00-0470-09	3000106824
3000103832	30001038320863	D684-00-0470-01U	3000106825
3000103832	30001038320864	D684-00-0470-01	3000106823
3000103832	30001038320864	D684-00-0470-09	3000106824
3000103832	30001038320864	D684-00-0470-01U	3000106825
3000103832	30001038320869	D684-00-0470-01	3000106823
3000103832	30001038320869	D684-00-0470-09	3000106824
3000103832	30001038320869	D684-00-0470-01U	3000106825
3000103832	30001038320870	D684-00-0470-01	3000106823
3000103832	30001038320870	D684-00-0470-09	3000106824
3000103832	30001038320870	D684-00-0470-01U	3000106825
3000103936	30001039360181	D684-00-0498-07	3000106964
3000103936	30001039360181	D684-00-0296-01U	3000107100
3000103936	30001039360181	D684-00-0296-01	3000107101
3000103936	30001039360181	D684-00-0498-01	3000107102
3000103936	30001039360182	D684-00-0498-07	3000106964
3000103936	30001039360182	D684-00-0296-01U	3000107100
3000103936	30001039360182	D684-00-0296-01	3000107101
3000103936	30001039360182	D684-00-0498-01	3000107102
3000103936	30001039360183	D684-00-0498-07	3000106964
3000103936	30001039360183	D684-00-0296-01U	3000107100
3000103936	30001039360183	D684-00-0296-01	3000107101
3000103936	30001039360183	D684-00-0498-01	3000107102
3000103936	30001039360184	D684-00-0498-07	3000106964
3000103936	30001039360184	D684-00-0296-01U	3000107100
3000103936	30001039360184	D684-00-0296-01	3000107101
3000103936	30001039360184	D684-00-0498-01	3000107102
3000103936	30001039360185	D684-00-0498-07	3000106964
3000103936	30001039360185	D684-00-0296-01U	3000107100
3000103936	30001039360185	D684-00-0296-01	3000107101
3000103936	30001039360185	D684-00-0498-01	3000107102
3000103937	30001039370351	D684-00-0480-07	3000107166
3000103937	30001039370352	D684-00-0480-07	3000107166
3000103937	30001039370353	D684-00-0480-07	3000107166
3000103937	30001039370354	D684-00-0480-07	3000107166
3000103937	30001039370355	D684-00-0480-07	3000107166
3000103938	30001039380155	D684-00-0480-07	3000107281
3000103938	30001039380155	D684-00-0480-07	3000107281
3000103938	30001039380155	D684-00-0480-07	3000107909
3000103938	30001039380155	D684-00-0480-02	3000108557
3000103938	30001039380155	D684-00-0480-02U	3000109024
3000103938	30001039380156	D684-00-0480-07	3000107281
3000103938	30001039380156	D684-00-0480-07	3000107281
3000103938	30001039380156	D684-00-0480-07	3000107909
3000103938	30001039380156	D684-00-0480-02	3000108557
3000103938	30001039380156	D684-00-0480-02U	3000109024
3000103938	30001039380157	D684-00-0480-07	3000107281
3000103938	30001039380157	D684-00-0480-07	3000107281
3000103938	30001039380157	D684-00-0480-07	3000107909
3000103938	30001039380157	D684-00-0480-02	3000108557
3000103938	30001039380157	D684-00-0480	

3000103938	30001039380158	D684-00-0480-07	3000107281
3000103938	30001039380158	D684-00-0480-07	3000107909
3000103938	30001039380158	D684-00-0480-02	3000108557
3000103938	30001039380158	D684-00-0480-02U	3000109024
3000103938	30001039380159	D684-00-0480-07	3000107281
3000103938	30001039380159	D684-00-0480-07	3000107281
3000103938	30001039380159	D684-00-0480-07	3000107909
3000103938	30001039380159	D684-00-0480-02	3000108557
3000103938	30001039380159	D684-00-0480-02U	3000109024
3000103940	30001039400104	D684-00-0469-01U	3000106754
3000103940	30001039400104	D684-00-0469-01	3000106755
3000103940	30001039400105	D684-00-0469-01U	3000106754
3000103940	30001039400105	D684-00-0469-01	3000106755
3000103940	30001039400107	D684-00-0469-01U	3000106754
3000103940	30001039400107	D684-00-0469-01	3000106755
3000103940	30001039400113	D684-00-0469-01U	3000106754
3000103940	30001039400113	D684-00-0469-01	3000106755
3000103940	30001039400118	D684-00-0469-01U	3000106754
3000103940	30001039400118	D684-00-0469-01	3000106755
3000103941	30001039410308	D684-00-0568-01	3000106962
3000103941	30001039410308	D684-00-0568-01U	3000106963
3000103941	30001039410309	D684-00-0568-01	3000106962
3000103941	30001039410309	D684-00-0568-01U	3000106963
3000103941	30001039410310	D684-00-0568-01	3000106962
3000103941	30001039410310	D684-00-0568-01U	3000106963
3000103941	30001039410315	D684-00-0568-01	3000106962
3000103941	30001039410315	D684-00-0568-01U	3000106963
3000103941	30001039410323	D684-00-0568-01	3000106962
3000103941	30001039410323	D684-00-0568-01U	3000106963
3000104652	30001046520101	D684-00-0295-05	3000107622
3000104652	30001046520101	D684-00-0295-01	3000107623
3000104652	30001046520101	D684-00-0295-01U	3000107624
3000104652	30001046520101	D684-00-0295-02U	3000112278
3000104652	30001046520102	D684-00-0295-05	3000107622
3000104652	30001046520102	D684-00-0295-01	3000107623
3000104652	30001046520102	D684-00-0295-01U	3000107624
3000104652	30001046520102	D684-00-0295-02U	3000112278
3000104652	30001046520103	D684-00-0295-05	3000107622
3000104652	30001046520103	D684-00-0295-01	3000107623
3000104652	30001046520103	D684-00-0295-01U	3000107624
3000104652	30001046520103	D684-00-0295-02U	3000112278
3000104652	30001046520104	D684-00-0295-05	3000107622
3000104652	30001046520104	D684-00-0295-01	3000107623
3000104652	30001046520104	D684-00-0295-01U	3000107624
3000104652	30001046520104	D684-00-0295-02U	3000112278
3000104652	30001046520105	D684-00-0295-05	3000107622
3000104652	30001046520105	D684-00-0295-01	3000107623
3000104652	30001046520105	D684-00-0295-01U	3000107624
3000104652	30001046520105	D684-00-0295-02U	3000112278
3000104726	30001047260031	D684-00-0607	3000107759
3000104726	30001047260032	D684-00-0607	3000107759
3000104726	30001047260034	D684-00-0607	3000107759
3000104726	30001047260036	D684-00-0607	3000107759
3000104726	30001047260037	D684-00-0607	3000107759
3000104728	30001047280375	D684-00-0608	3000107760
3000104728	30001047280375	D684-00-0608	3000110622
3000104728	30001047280380	D684-00-0608	3000107760
3000104728	30001047280380	D684-00-0608	3000110622

3000104728	30001047280385	D684-00-0608	3000107760
3000104728	30001047280385	D684-00-0608	3000110622
3000104728	30001047280391	D684-00-0608	3000107760
3000104728	30001047280391	D684-00-0608	3000110622
3000104728	30001047280393	D684-00-0608	3000107760
3000104728	30001047280393	D684-00-0608	3000110622
3000104730	30001047300131	D684-00-0296-02U	3000108432
3000104730	30001047300131	D684-DO-0497	0
3000104730	30001047300131	D684-00-0296-01	3000108962
3000104730	30001047300131	D684-00-0498-01	3000108963
3000104730	30001047300132	D684-00-0296-02U	3000108432
3000104730	30001047300132	D684-DO-0497	0
3000104730	30001047300132	D684-00-0296-01	3000108962
3000104730	30001047300132	D684-00-0498-01	3000108963
3000104730	30001047300133	D684-00-0296-02U	3000108432
3000104730	30001047300133	D684-DO-0497	0
3000104730	30001047300133	D684-00-0296-01	3000108962
3000104730	30001047300133	D684-00-0498-01	3000108963
3000104730	30001047300134	D684-00-0296-02U	3000108432
3000104730	30001047300134	D684-DO-0497	0
3000104730	30001047300134	D684-00-0296-01	3000108962
3000104730	30001047300134	D684-00-0498-01	3000108963
3000104730	30001047300135	D684-00-0296-02U	3000108432
3000104730	30001047300135	D684-DO-0497	0
3000104730	30001047300135	D684-00-0296-01	3000108962
3000104730	30001047300135	D684-00-0498-01	3000108963
3000104731	30001047310335	D684-00-0568-01U	3000107331
3000104731	30001047310335	D684-00-0568-01U	3000107332
3000104731	30001047310338	D684-00-0568-01U	3000107331
3000104731	30001047310338	D684-00-0568-01U	3000107332
3000104731	30001047310340	D684-00-0568-01U	3000107331
3000104731	30001047310340	D684-00-0568-01U	3000107332
3000104731	30001047310345	D684-00-0568-01U	3000107331
3000104731	30001047310345	D684-00-0568-01U	3000107332
3000104731	30001047310347	D684-00-0568-01U	3000107331
3000104731	30001047310347	D684-00-0568-01U	3000107332
3000104732	30001047320451	D684-00-0479-07	3000107908
3000104732	30001047320451	D684-00-0479-02	3000108048
3000104732	30001047320452	D684-00-0479-07	3000107908
3000104732	30001047320452	D684-00-0479-02	3000108048
3000104732	30001047320453	D684-00-0479-07	3000107908
3000104732	30001047320453	D684-00-0479-02	3000108048
3000104732	30001047320454	D684-00-0479-07	3000107908
3000104732	30001047320454	D684-00-0479-02	3000108048
3000104732	30001047320455	D684-00-0479-07	3000107908
3000104732	30001047320455	D684-00-0479-02	3000108048
3000104733	30001047330172	D684-00-0470-09	3000107336
3000104733	30001047330172	D684-00-0470-01U	3000107337
3000104733	30001047330172	D684-00-0470-01C	3000110211
3000104733	30001047330172	D684-00-0470-01	3000110212
3000104733	30001047330177	D684-00-0470-09	3000107336
3000104733	30001047330177	D684-00-0470-01U	3000107337
3000104733	30001047330177	D684-00-0470-01C	3000110211
3000104733	30001047330177	D684-00-0470-01	3000110212
3000104733	30001047330179	D684-00-0470-09	3000107336
3000104733	30001047330179	D684-00-0470-01U	3000107337
3000104733	30001047330179	D684-00-0470-01C	3000110211
3000104733	30001047330179	D684-00-0470-01	3000110212

3000104733	30001047330189	D684-00-0470-09	3000107336
3000104733	30001047330189	D684-00-0470-01U	3000107337
3000104733	30001047330189	D684-00-0470-01C	3000110211
3000104733	30001047330189	D684-00-0470-01	3000110212
3000104733	30001047330194	D684-00-0470-09	3000107336
3000104733	30001047330194	D684-00-0470-01U	3000107337
3000104733	30001047330194	D684-00-0470-01C	3000110211
3000104733	30001047330194	D684-00-0470-01	3000110212
3000104997	30001049970113	D684-00-0480-01C	3000109023
3000104997	30001049970113	D684-00-0480-10	3000109400
3000104997	30001049970114	D684-00-0480-01C	3000109023
3000104997	30001049970114	D684-00-0480-10	3000109400
3000104997	30001049970115	D684-00-0480-01C	3000109023
3000104997	30001049970115	D684-00-0480-10	3000109400
3000104997	30001049970116	D684-00-0480-01C	3000109023
3000104997	30001049970116	D684-00-0480-10	3000109400
3000104997	30001049970117	D684-00-0480-01C	3000109023
3000104997	30001049970117	D684-00-0480-10	3000109400
3000105002	30001050020133	D684-00-0576-01U	3000108349
3000105002	30001050020133	D684-00-0576-01	3000108350
3000105002	30001050020138	D684-00-0576-01U	3000108349
3000105002	30001050020138	D684-00-0576-01	3000108350
3000105002	30001050020145	D684-00-0576-01U	3000108349
3000105002	30001050020145	D684-00-0576-01	3000108350
3000105002	30001050020146	D684-00-0576-01U	3000108349
3000105002	30001050020146	D684-00-0576-01	3000108350
3000105002	30001050020148	D684-00-0576-01U	3000108349
3000105002	30001050020148	D684-00-0576-01	3000108350
3000105003	30001050030001	D684-00-0480-01	3000109026
3000105003	30001050030002	D684-00-0480-01	3000109026
3000105003	30001050030003	D684-00-0480-01	3000109026
3000105003	30001050030004	D684-00-0480-01	3000109026
3000105003	30001050030005	D684-00-0480-01	3000109026
3000105004	30001050040130	D684-00-0294-01	3000109816
3000105004	30001050040130	D684-00-0294-01U	3000112151
3000105004	30001050040130	D684-00-0294-02U	3000112152
3000105004	30001050040130	D684-00-0294-02	3000112153
3000105004	30001050040130	D684-00-0294-01U	3000114260
3000105004	30001050040131	D684-00-0294-01	3000109816
3000105004	30001050040131	D684-00-0294-01U	3000112151
3000105004	30001050040131	D684-00-0294-02U	3000112152
3000105004	30001050040131	D684-00-0294-02	3000112153
3000105004	30001050040131	D684-00-0294-01U	3000114260
3000105004	30001050040132	D684-00-0294-01	3000109816
3000105004	30001050040132	D684-00-0294-01U	3000112151
3000105004	30001050040132	D684-00-0294-02U	3000112152
3000105004	30001050040132	D684-00-0294-02	3000112153
3000105004	30001050040132	D684-00-0294-01U	3000114260
3000105004	30001050040133	D684-00-0294-01	3000109816
3000105004	30001050040133	D684-00-0294-01U	3000112151
3000105004	30001050040133	D684-00-0294-02U	3000112152
3000105004	30001050040133	D684-00-0294-02	3000112153
3000105004	30001050040133	D684-00-0294-01U	3000114260
3000105006	30001050060485	D684-00-0568-01U	3000108589
3000105006	30001050060485	D684-00-0568-01	3000108882
3000105006	30001050060487	D684-00-0568-01U	3000108589
3000105006	30001050060487	D684-00-0568-01	3000108882
3000105006	30001050060488	D684-00-0568-01U	3000108589

3000105006	30001050060488	D684-00-0568-01	3000108882
3000105006	30001050060498	D684-00-0568-01U	3000108589
3000105006	30001050060498	D684-00-0568-01	3000108882
3000105006	30001050060499	D684-00-0568-01U	3000108589
3000105006	30001050060499	D684-00-0568-01	3000108882
3000105506	30001055060101	D684-00-0479-01C	3000109450
3000105506	30001055060101	D684-00-0479-02U	3000109451
3000105506	30001055060101	D684-00-0479-01U	3000109452
3000105506	30001055060101	D684-00-0479-01	3000109449
3000105506	30001055060102	D684-00-0479-01C	3000109450
3000105506	30001055060102	D684-00-0479-02U	3000109451
3000105506	30001055060102	D684-00-0479-01U	3000109452
3000105506	30001055060102	D684-00-0479-01	3000109449
3000105506	30001055060103	D684-00-0479-01C	3000109450
3000105506	30001055060103	D684-00-0479-02U	3000109451
3000105506	30001055060103	D684-00-0479-01U	3000109452
3000105506	30001055060103	D684-00-0479-01	3000109449
3000105506	30001055060104	D684-00-0479-01C	3000109450
3000105506	30001055060104	D684-00-0479-02U	3000109451
3000105506	30001055060104	D684-00-0479-01U	3000109452
3000105506	30001055060104	D684-00-0479-01	3000109449
3000105506	30001055060105	D684-00-0479-01C	3000109450
3000105506	30001055060105	D684-00-0479-02U	3000109451
3000105506	30001055060105	D684-00-0479-01U	3000109452
3000105506	30001055060105	D684-00-0479-01	3000109449
3000105516	30001055160159	D684-00-0576-01U	3000108590
3000105516	30001055160159	D684-00-0576-01	3000108591
3000105516	30001055160163	D684-00-0576-01U	3000108590
3000105516	30001055160163	D684-00-0576-01	3000108591
3000105516	30001055160168	D684-00-0576-01U	3000108590
3000105516	30001055160168	D684-00-0576-01	3000108591
3000105516	30001055160171	D684-00-0576-01U	3000108590
3000105516	30001055160171	D684-00-0576-01	3000108591
3000105516	30001055160172	D684-00-0576-01U	3000108590
3000105516	30001055160172	D684-00-0576-01	3000108591
3000105530	30001055300001	D684-00-0480-01U	3000109025
3000105530	30001055300001	D684-00-0480-02	3000109401
3000105530	30001055300002	D684-00-0480-01U	3000109025
3000105530	30001055300002	D684-00-0480-02	3000109401
3000105530	30001055300003	D684-00-0480-01U	3000109025
3000105530	30001055300003	D684-00-0480-02	3000109401
3000105530	30001055300004	D684-00-0480-01U	3000109025
3000105530	30001055300004	D684-00-0480-02	3000109401
3000105530	30001055300005	D684-00-0480-01U	3000109025
3000105530	30001055300005	D684-00-0480-02	3000109401
3000106324	30001063240668	D684-00-0576-01U	3000108682
3000106324	30001063240668	D684-00-0576-01U	3000108683
3000106324	30001063240669	D684-00-0576-01U	3000108682
3000106324	30001063240669	D684-00-0576-01U	3000108683
3000106324	30001063240672	D684-00-0576-01U	3000108682
3000106324	30001063240672	D684-00-0576-01U	3000108683
3000106324	30001063240682	D684-00-0576-01U	3000108682
3000106324	30001063240682	D684-00-0576-01U	3000108683
3000106324	30001063240685	D684-00-0576-01U	3000108682
3000106324	30001063240685	D684-00-0576-01U	3000108683
3000106326	30001063260112	D684-00-0559-01	3000110583
3000106326	300010		

[illegible]

3000106387	30001063870210	D684-00-0469-01U	3000110207
3000106387	30001063870210	D684-00-0469-01	3000110208
3000106387	30001063870210	D684-00-0469-01C	3000110209
3000106387	30001063870210	D684-00-0469-09	3000110210
3000106388	30001063880190	D684-00-0607	3000110626
3000106388	30001063880190	D684-00-0607	3000113155
3000106388	30001063880193	D684-00-0607	3000110626
3000106388	30001063880193	D684-00-0607	3000113155
3000106388	30001063880204	D684-00-0607	3000110626
3000106388	30001063880204	D684-00-0607	3000113155
3000106388	30001063880206	D684-00-0607	3000110626
3000106388	30001063880206	D684-00-0607	3000113155
3000106388	30001063880210	D684-00-0607	3000110626
3000106388	30001063880210	D684-00-0607	3000113155
3000106389	30001063890091	D684-00-0480-07	3000110205
3000106389	30001063890091	D684-00-0480-01U	3000110414
3000106389	30001063890092	D684-00-0480-07	3000110205
3000106389	30001063890092	D684-00-0480-01U	3000110414
3000106389	30001063890093	D684-00-0480-07	3000110205
3000106389	30001063890093	D684-00-0480-01U	3000110414
3000106389	30001063890094	D684-00-0480-07	3000110205
3000106389	30001063890094	D684-00-0480-01U	3000110414
3000106389	30001063890095	D684-00-0480-07	3000110205
3000106389	30001063890095	D684-00-0480-01U	3000110414
3000106455	30001064550001	D684-00-0296-01U	3000109817
3000106455	30001064550001	D684-00-0296-01	3000109818
3000106455	30001064550002	D684-00-0296-01U	3000109817
3000106455	30001064550002	D684-00-0296-01	3000109818
3000106455	30001064550003	D684-00-0296-01U	3000109817
3000106455	30001064550003	D684-00-0296-01	3000109818
3000106455	30001064550004	D684-00-0296-01U	3000109817
3000106455	30001064550004	D684-00-0296-01	3000109818
3000106455	30001064550005	D684-00-0296-01U	3000109817
3000106455	30001064550005	D684-00-0296-01	3000109818
3000106669	30001066690306	D684-00-0608	3000110623
3000106669	30001066690307	D684-00-0608	3000110623
3000106669	30001066690308	D684-00-0608	3000110623
3000106669	30001066690309	D684-00-0608	3000110623
3000106669	30001066690310	D684-00-0608	3000110623
3000106670	30001066700011	D684-00-0479-07	3000110206
3000106670	30001066700011	D684-00-0479-01C	3000110736
3000106670	30001066700012	D684-00-0479-07	3000110206
3000106670	30001066700012	D684-00-0479-01C	3000110736
3000106670	30001066700013	D684-00-0479-07	3000110206
3000106670	30001066700013	D684-00-0479-01C	3000110736
3000106670	30001066700014	D684-00-0479-07	3000110206
3000106670	30001066700014	D684-00-0479-01C	3000110736
3000106670	30001066700015	D684-00-0479-07	3000110206
3000106670	30001066700015	D684-00-0479-01C	3000110736
3000106672	30001066720006	D684-00-0576-01U	3000109819
3000106672	30001066720011	D684-00-0576-01U	3000109819
3000106672	30001066720020	D684-00-0576-01U	3000109819
3000106672	30001066720022	D684-00-0576-01U	3000109819
3000106672	30001066720025	D684-00-0576-01U	3000109819
3000107115	30001071150084	D684-00-0480-02	3000110454
3000107115	30001071150084	D684-00-0480-02	3000110735
3000107115	30001071150085	D684-00-0480-02	3000110454
3000107115	30001071150085	D684-00-0480-02	3000110735

3000107115	30001071150086	D684-00-0480-02	3000110454
3000107115	30001071150086	D684-00-0480-02	3000110735
3000107115	30001071150087	D684-00-0480-02	3000110454
3000107115	30001071150087	D684-00-0480-02	3000110735
3000107115	30001071150088	D684-00-0480-02	3000110454
3000107115	30001071150088	D684-00-0480-02	3000110735
3000107117	30001071170255	D684-00-0568-01U	3000110520
3000107117	30001071170263	D684-00-0568-01U	3000110520
3000107117	30001071170265	D684-00-0568-01U	3000110520
3000107117	30001071170271	D684-00-0568-01U	3000110520
3000107117	30001071170272	D684-00-0568-01U	3000110520
3000107127	30001071270101	D684-00-0480-01C	3000110738
3000107127	30001071270101	D684-00-0480-01U	3000111358
3000107127	30001071270102	D684-00-0480-01C	3000110738
3000107127	30001071270102	D684-00-0480-01U	3000111358
3000107127	30001071270103	D684-00-0480-01C	3000110738
3000107127	30001071270103	D684-00-0480-01U	3000111358
3000107127	30001071270104	D684-00-0480-01C	3000110738
3000107127	30001071270104	D684-00-0480-01U	3000111358
3000107127	30001071270105	D684-00-0480-01C	3000110738
3000107127	30001071270105	D684-00-0480-01U	3000111358
3000107130	30001071300350	D684-00-0470-01U	3000111022
3000107130	30001071300350	D684-00-0470-01	3000111023
3000107130	30001071300350	D684-00-0470-07	3000111134
3000107130	30001071300350	D684-00-0470-09	3000111191
3000107130	30001071300358	D684-00-0470-01U	3000111022
3000107130	30001071300358	D684-00-0470-01	3000111023
3000107130	30001071300358	D684-00-0470-07	3000111134
3000107130	30001071300358	D684-00-0470-09	3000111191
3000107130	30001071300364	D684-00-0470-01U	3000111022
3000107130	30001071300364	D684-00-0470-01	3000111023
3000107130	30001071300364	D684-00-0470-07	3000111134
3000107130	30001071300364	D684-00-0470-09	3000111191
3000107130	30001071300368	D684-00-0470-01U	3000111022
3000107130	30001071300368	D684-00-0470-01	3000111023
3000107130	30001071300368	D684-00-0470-07	3000111134
3000107130	30001071300368	D684-00-0470-09	3000111191
3000107130	30001071300370	D684-00-0470-01U	3000111022
3000107130	30001071300370	D684-00-0470-01	3000111023
3000107130	30001071300370	D684-00-0470-07	3000111134
3000107130	30001071300370	D684-00-0470-09	3000111191
3000107635	30001076350005	D684-00-0480-01	3000111350
3000107635	30001076350005	D684-00-0480-01U	3000111359
3000107635	30001076350005	D684-00-0480-01	3000111973
3000107895	30001078950239	D684-00-0576-01	3000110954
3000107895	30001078950239	D684-00-0576-01U	3000110955
3000107895	30001078950244	D684-00-0576-01	3000110954
3000107895	30001078950244	D684-00-0576-01U	3000110955
3000107895	30001078950245	D684-00-0576-01	3000110954
3000107895	30001078950245	D684-00-0576-01U	3000110955
3000107895	30001078950249	D684-00-0576-01	3000110954
3000107895	30001078950249	D684-00-0576-01U	3000110955
3000107895	30001078950253	D684-00-0576-01	3000110954
3000107895	30001078950253	D684-00-0576-01U	3000110955
3000107898	30001078980125	D684-00-0469-07	3000111135
3000107898	30001078980125	D684-00-0469-01	3000111136
3000107898	30001078980125	D684-00-0469-01U	3000111137
3000107898	30001078980125	D684-00-0469-09	3000111192

3000107898	30001078980126	D684-00-0469-07	3000111135
3000107898	30001078980126	D684-00-0469-01	3000111136
3000107898	30001078980126	D684-00-0469-01U	3000111137
3000107898	30001078980126	D684-00-0469-09	3000111192
3000107898	30001078980140	D684-00-0469-07	3000111135
3000107898	30001078980140	D684-00-0469-01	3000111136
3000107898	30001078980140	D684-00-0469-01U	3000111137
3000107898	30001078980140	D684-00-0469-09	3000111192
3000108240	30001082400445	D684-00-0498-01	3000111680
3000108240	30001082400445	D684-00-0296-01U	3000111681
3000108240	30001082400446	D684-00-0498-01	3000111680
3000108240	30001082400446	D684-00-0296-01U	3000111681
3000108240	30001082400448	D684-00-0498-01	3000111680
3000108240	30001082400448	D684-00-0296-01U	3000111681
3000108240	30001082400449	D684-00-0498-01	3000111680
3000108240	30001082400449	D684-00-0296-01U	3000111681
3000108241	30001082410008	D684-00-0568-01U	3000111682
3000108241	30001082410010	D684-00-0568-01U	3000111682
3000108241	30001082410019	D684-00-0568-01U	3000111682
3000108241	30001082410020	D684-00-0568-01U	3000111682
3000108241	30001082410024	D684-00-0568-01U	3000111682
3000108352	30001083520001	D684-00-0295-01	3000112279
3000108352	30001083520001	D684-00-0295-01U	3000112280
3000108352	30001083520001	D684-00-0295-05	3000112281
3000108352	30001083520001	D684-00-0295-05	3000114254
3000108352	30001083520002	D684-00-0295-01	3000112279
3000108352	30001083520002	D684-00-0295-01U	3000112280
3000108352	30001083520002	D684-00-0295-05	3000112281
3000108352	30001083520002	D684-00-0295-05	3000114254
3000108352	30001083520003	D684-00-0295-01	3000112279
3000108352	30001083520003	D684-00-0295-01U	3000112280
3000108352	30001083520003	D684-00-0295-05	3000112281
3000108352	30001083520003	D684-00-0295-05	3000114254
3000108352	30001083520004	D684-00-0295-01	3000112279
3000108352	30001083520004	D684-00-0295-01U	3000112280
3000108352	30001083520004	D684-00-0295-05	3000112281
3000108352	30001083520004	D684-00-0295-05	3000114254
3000108352	30001083520006	D684-00-0295-01	3000112279
3000108352	30001083520006	D684-00-0295-01U	3000112280
3000108352	30001083520006	D684-00-0295-05	3000112281
3000108352	30001083520006	D684-00-0295-05	3000114254
3000108620	30001086200451	D684-00-0576-01U	3000111683
3000108620	30001086200451	D684-00-0576-01U	3000111684
3000108620	30001086200453	D684-00-0576-01U	3000111683
3000108620	30001086200453	D684-00-0576-01U	3000111684
3000108620	30001086200460	D684-00-0576-01U	3000111683
3000108620	30001086200460	D684-00-0576-01U	3000111684
3000108620	30001086200465	D684-00-0576-01U	3000111683
3000108620	30001086200465	D684-00-0576-01U	3000111684
3000108620	30001086200470	D684-00-0576-01U	3000111683
3000108620	30001086200470	D684-00-0576-01U	3000111684
3000108622	30001086220102	D684-00-0470-01	3000112113
3000108622	30001086220102	D684-00-0470-01U	3000112114
3000108622	30001086220102	D684-00-0470-09	3000112115
3000108622	30001086220108	D684-00-0470-01	3000112113
3000108622	30001086220108	D684-00-0470-01U	3000112114
3000108622	30001086220108	D684-00-0470-09	3000112115
3000108622	30001086220109	D684-00-0470-01	3000112113

3000108622	30001086220109	D684-00-0470-01U	3000112114
3000108622	30001086220109	D684-00-0470-09	3000112115
3000108622	30001086220116	D684-00-0470-01	3000112113
3000108622	30001086220116	D684-00-0470-01U	3000112114
3000108622	30001086220116	D684-00-0470-09	3000112115
3000108622	30001086220118	D684-00-0470-01	3000112113
3000108622	30001086220118	D684-00-0470-01U	3000112114
3000108622	30001086220118	D684-00-0470-09	3000112115
3000109189	30001091890060	D684-00-0296-02U	3000112154
3000109189	30001091890060	D684-00-0296-02	3000112155
3000109189	30001091890060	D684-00-0296-02	3000112156
3000109189	30001091890060	D684-00-0296-01	3000112782
3000109189	30001091890061	D684-00-0296-02U	3000112154
3000109189	30001091890061	D684-00-0296-02	3000112155
3000109189	30001091890061	D684-00-0296-02	3000112156
3000109189	30001091890061	D684-00-0296-01	3000112782
3000109189	30001091890062	D684-00-0296-02U	3000112154
3000109189	30001091890062	D684-00-0296-02	3000112155
3000109189	30001091890062	D684-00-0296-02	3000112156
3000109189	30001091890062	D684-00-0296-01	3000112782
3000109189	30001091890064	D684-00-0296-02U	3000112154
3000109189	30001091890064	D684-00-0296-02	3000112155
3000109189	30001091890064	D684-00-0296-02	3000112156
3000109189	30001091890064	D684-00-0296-01	3000112782
3000109191	30001091910277	D684-00-0607	3000112276
3000109191	30001091910277	D684-00-0607	3000113156
3000109191	30001091910279	D684-00-0607	3000112276
3000109191	30001091910279	D684-00-0607	3000113156
3000109191	30001091910285	D684-00-0607	3000112276
3000109191	30001091910285	D684-00-0607	3000113156
3000109191	30001091910292	D684-00-0607	3000112276
3000109191	30001091910292	D684-00-0607	3000113156
3000109191	30001091910297	D684-00-0607	3000112276
3000109191	30001091910297	D684-00-0607	3000113156
3000109197	30001091970850	D684-00-0608	3000112744
3000109197	30001091970850	D684-00-0608	3000113154
3000109197	30001091970857	D684-00-0608	3000112744
3000109197	30001091970857	D684-00-0608	3000113154
3000109197	30001091970859	D684-00-0608	3000112744
3000109197	30001091970859	D684-00-0608	3000113154
3000109690	30001096900001	D684-00-0479-01	3000113519
3000109690	30001096900001	D684-00-0479-01C	3000113520
3000109690	30001096900001	D684-00-0479-02	3000113622
3000109690	30001096900002	D684-00-0479-01	3000113519
3000109690	30001096900002	D684-00-0479-01C	3000113520
3000109690	30001096900002	D684-00-0479-02	3000113622
3000109690	30001096900003	D684-00-0479-01	3000113519
3000109690	30001096900003	D684-00-0479-01C	3000113520
3000109690	30001096900003	D684-00-0479-02	3000113622
3000109690	30001096900004	D684-00-0479-01	3000113519
3000109690	30001096900004	D684-00-0479-01C	3000113520
3000109690	30001096900004	D684-00-0479-02	3000113622
3000109690	30001096900005	D684-00-0479-01	3000113519
3000109690	30001096900005	D684-00-0479-01C	3000113520
3000109690	30001096900005	D684-00-0479-02	3000113622
3000109691	30001096910625	D684-00-0576-01U	3000112929
3000109691	30001096910625	D684-00-0576-01U	3000112930
3000109691	30001096910625	D684-00-0576-01	3000112931

3000109691	30001096910625	D684-00-0576-09	3000112932
3000109691	30001096910627	D684-00-0576-01U	3000112929
3000109691	30001096910627	D684-00-0576-01U	3000112930
3000109691	30001096910627	D684-00-0576-01	3000112931
3000109691	30001096910627	D684-00-0576-09	3000112932
3000109691	30001096910628	D684-00-0576-01U	3000112929
3000109691	30001096910628	D684-00-0576-01U	3000112930
3000109691	30001096910628	D684-00-0576-01	3000112931
3000109691	30001096910628	D684-00-0576-09	3000112932
3000109691	30001096910629	D684-00-0576-01U	3000112929
3000109691	30001096910629	D684-00-0576-01U	3000112930
3000109691	30001096910629	D684-00-0576-01	3000112931
3000109691	30001096910629	D684-00-0576-09	3000112932
3000109691	30001096910648	D684-00-0576-01U	3000112929
3000109691	30001096910648	D684-00-0576-01U	3000112930
3000109691	30001096910648	D684-00-0576-01	3000112931
3000109691	30001096910648	D684-00-0576-09	3000112932
3000110167	30001101670130	D684-00-0498-01	3000113152
3000110167	30001101670130	D684-00-0296-01U	3000113153
3000110167	30001101670131	D684-00-0498-01	3000113152
3000110167	30001101670131	D684-00-0296-01U	3000113153
3000110167	30001101670132	D684-00-0498-01	3000113152
3000110167	30001101670132	D684-00-0296-01U	3000113153
3000110167	30001101670133	D684-00-0498-01	3000113152
3000110167	30001101670133	D684-00-0296-01U	3000113153
3000110167	30001101670134	D684-00-0498-01	3000113152
3000110167	30001101670134	D684-00-0296-01U	3000113153
3000110617	30001106170185	D684-00-0576-01	3000113623
3000110617	30001106170185	D684-00-0576-01U	3000113624
3000110617	30001106170185	D684-00-0576-01U	3000113625
3000110617	30001106170190	D684-00-0576-01	3000113623
3000110617	30001106170190	D684-00-0576-01U	3000113624
3000110617	30001106170190	D684-00-0576-01U	3000113625
3000110617	30001106170195	D684-00-0576-01	3000113623
3000110617	30001106170195	D684-00-0576-01U	3000113624
3000110617	30001106170195	D684-00-0576-01U	3000113625
3000110617	30001106170196	D684-00-0576-01	3000113623
3000110617	30001106170196	D684-00-0576-01U	3000113624
3000110617	30001106170196	D684-00-0576-01U	3000113625
3000110834	30001108340524	D684-00-0568-01	3000113809
3000110834	30001108340524	D684-00-0568-01U	3000113810
3000110834	30001108340524	D684-00-0568-01	3000115337
3000111231	30001112310001	D684-00-0480-01C	3000115116
3000111231	30001112310001	D684-00-0480-02	3000115917
3000111231	30001112310002	D684-00-0480-01C	3000115116
3000111231	30001112310002	D684-00-0480-02	3000115917
3000111231	30001112310003	D684-00-0480-01C	3000115116
3000111231	30001112310003	D684-00-0480-02	3000115917
3000111231	30001112310004	D684-00-0480-01C	3000115116
3000111231	30001112310004	D684-00-0480-02	3000115917
3000111231	30001112310005	D684-00-0480-01C	3000115116
3000111231	30001112310005	D684-00-0480-02	3000115917
3000111839	30001118390051	D684-00-0479-01C	3000115115
3000111839	30001118390051	D684-00-0479-02	3000115918
3000111839	30001118390051	D684-00-0479-01	3000116780
3000111839	30001118390051	D684-00-0479-01	3000116780
3000111839	30001118390052	D684-00-0479-01C	3000115115
3000111839	30001118390052	D684-00-0479-02	3000115918

3000111839	30001118390052	D684-00-0479-01	3000116780
3000111839	30001118390052	D684-00-0479-01	3000116780
3000111839	30001118390053	D684-00-0479-01C	3000115115
3000111839	30001118390053	D684-00-0479-02	3000115918
3000111839	30001118390053	D684-00-0479-01	3000116780
3000111839	30001118390053	D684-00-0479-01	3000116780
3000111839	30001118390055	D684-00-0479-01C	3000115115
3000111839	30001118390055	D684-00-0479-02	3000115918
3000111839	30001118390055	D684-00-0479-01	3000116780
3000111839	30001118390055	D684-00-0479-01	3000116780
3000112417	30001124170005	D684-00-0480-01U	3000116581
3000112417	30001124170005	D684-00-0480-10	3000116582
3000112417	30001124170005	D684-00-0480-01C	3000116665
3000112417	30001124170007	D684-00-0480-01U	3000116581
3000112417	30001124170007	D684-00-0480-10	3000116582
3000112417	30001124170007	D684-00-0480-01C	3000116665
3000113628	30001136280459	D684-00-0576-01U	3000117020
3000113628	30001136280463	D684-00-0576-01U	3000117020
3000113628	30001136280470	D684-00-0576-01U	3000117020
3000113628	30001136280471	D684-00-0576-01U	3000117020
3000113628	30001136280472	D684-00-0576-01U	3000117020