

December-2022

URGENT – FIELD SAFETY NOTICE

Type of Action	Recall
Teleflex Reference	EIF-000523
Batch/Lot Number	Refer to Appendix 2
Commercial Name	Product Code
Slick Set [®] Cuffed Endotracheal Tube and Stylet Set	170060; 170070; 170075; 170080; 170085; 170090
Flexi-Set Cuffed Endotracheal Tube and Stylet Set, oral/nasal	504550; 504560; 504580; 504585; 504590
Endotracheal Tube oral/nasal Magill, High Volume, Low Pressure Cuff	112080-000065; 112080-000080
Endotracheal Tube oral/nasal Murphy Eye, High Volume, Low Pressure Cuff	112082-000050; 112082-000055; 112082-000060; 112082-000065; 112082-000070; 112082-000075; 112082-000080; 112082-000085; 112082050; 112082055; 112082060; 112082065; 112082070; 112082075; 112082080; 112082085; 112082095
Flexiset Super Safety Clear Tracheal Tube, oral/nasal with Cuff and Insertion Aid	112481-000055; 112481-000060; 112481-000070; 112481-000075; 112481-000080; 112481-000085; 112481-000090
RUSCHELIT [®] Safety Clear Tracheal Tube, oral/nasal, Magill	112480-000050; 112480-000055; 112480-000060; 112480-000065; 112480-000070; 112480-000075; 112480-000080; 112480-000085; 112480-000090; 112480-000095
RUSCHELIT [®] Super Safety Clear Tracheal Tube, oral/nasal, Murphy	112482065LA; 112482070LA; 112482075LA; 112482080LA; 112482085LA; 112482095LA; 112482100LA
RUSCHELIT [®] Safety Clear Tracheal Tube, oral/nasal, Murphy	112482-000050; 112482-000055; 112482-000060; 112482-000065; 112482-000070; 112482-000075; 112482-000080; 112482-000085; 112482-000090; 112482-000095

Dear Customer,

Details of affected devices

Teleflex Medical Europe Limited has initiated a voluntary Field Safety Corrective Action (“FSCA”) for the above listed products; refer to Appendix 2 for Unique Device Identifier (“UDI”) and lot number information.

Description of the problem & immediate actions required

This voluntary FSCA for the above-listed products has been initiated due to reports of pilot balloon non-inflation or cuff non-deflation.

In the event of pilot balloon non-inflation, there is a potential risk of injury to the patient as the pilot balloon no longer reflects the condition of the cuff, which may lead to over-inflation of the cuff. The risk persists in situ, therefore clinicians should evaluate the risk benefit to individual patients of extubation and reintubation.

In the event of cuff non-deflation, there is the potential risk of injury if the device is removed from the patient with the cuff still inflated. If extubation and removal of the endotracheal tube is planned

and the cuff cannot be deflated via the Bepak valve, then the inflation line should be cut in order to facilitate passive deflation of the cuff.

No deaths or long-term patient injuries have been reported at this time.

Our records indicate you have received products that are subject to this FSCA.

Depending on your device location please adhere to the following Action list:

Device location	Action List Number
Medical facilities (hospitals, medical staff, etc.)	1
Distributors	2

Action list number 1 – Medical facilities

1. We request that you immediately check your inventory for product within the scope of this FSCA. Users should cease use and distribution of affected product and immediately quarantine the affected product.
2. If you have impacted product, mark the applicable checkbox on the Acknowledgement Form (Appendix 1) and contact Teleflex Customer Service by calling the phone number provided below. Teleflex Customer Service will issue a Return Goods Authorisation (RGA) number to you. Write the (RGA) number into the respective field in the Acknowledgement Form and promptly return this form to the e-mail address below.
3. If you do not have impacted product, mark the applicable checkbox on the Acknowledgement Form (Appendix 1) and return the form to Teleflex at the contact details provided.
4. Teleflex (or your local dealer) will issue a credit note upon receipt of the returned affected product.

Action list number 2 – Distributors

1. Provide this field safety notice to all customers who have received impacted product. Each of your customers is then required to complete the Acknowledgement Form and return it to you.
2. We request that you immediately check your inventory for impacted product. Cease use and distribution of impacted product and immediately quarantine the affected product. You may then return all product in scope.
3. As a distributor, you are then required to confirm to Teleflex that you have completed the field activity outlined in actions 1 and 2 of this Action List Number 2. Upon completion of your actions, please forward the completed Acknowledgement Form to the e-mail address below.
Important - Please ensure you only list batch numbers in scope of this Field Safety Notice when completing this form.
4. Please be aware that all European Economic Area/Switzerland, United Kingdom (EEA/CH/UK) and Turkey (TR) Competent Authorities in which Teleflex distribute directly will be notified by Teleflex.
5. If you have further distributed product outside of your country, please notify Teleflex Customer Service by return e-mail to the e-mail address below.
6. If you are a distributor and/or have a reporting responsibility within or outside the EEA/CH/UK/TR area, please notify your local Competent Authority of this action. Please forward the notification and all communication with your local competent authority to Teleflex.

Adverse reactions or quality problems experienced with the use of this product should be reported to Teleflex Customer Service at the contact information below.

Transmission of this Field Safety Notice

This notice should be passed on to all persons who need to be aware within your organisation or to any organisation where the potentially affected devices have been transferred.

Please consider end users, clinicians, risk managers, supply chain/distribution centres, etc., in the circulation of this notice. Please maintain awareness of this notice until all required actions have been completed in your organisation.

Contact reference person

Should you require any further information or support concerning this issue, please contact:

Customer Service:

Contact: Nicole Morawiec

Telephone: +41 (0) 31 818 40 90

FAX: +41 (0) 31 818 40 93

Email: info.ch@teleflex.com

Teleflex is committed to providing high quality, safe and effective products. We regret any inconvenience this action may cause your operations. If you have any other questions, please contact your local Teleflex sales representative or Teleflex Customer Service.

For and on behalf of Teleflex,

Padraig Hegarty

Padraig Hegarty VP, Global QA (Manufacturing)

FIELD SAFETY CORRECTIVE ACTION
ACKNOWLEDGEMENT FORM

PRODUCT FIELD ACTION BY TELEFLEX – IMMEDIATE ATTENTION REQUIRED

Ref. EIF-000523

RETURN COMPLETED FORM IMMEDIATELY TO:

E-mail: info.ch@teleflex.com

<input type="checkbox"/> We confirm receipt of this FSN and completion of the required actions contained therein. We further confirm that our inventory does NOT include products impacted by this Field Action.	<input type="checkbox"/> We confirm receipt of this FSN and completion of the required actions contained therein. We further confirm our inventory DOES include products impacted by this Field Action. The use and further distribution of the affected products is stopped. All products are put on hold and the amount below will be returned. Return Goods Authorisation No _____
---	--

Complete this Acknowledgement Form and return the completed form immediately using the contact information above.

Product code	Lot/batch number	Quantity returning
Important - Please ensure you only list batch numbers in scope of this Field Safety Notice when completing this form.		
<ul style="list-style-type: none"> • Include a copy of the completed Acknowledgement Form in the returns package with the returned units • Ensure the RGA number is clearly visible on the returns package • Please label returns as "Field Safety Returns" 		
Note: Non-FSCA product returns should be processed per standard product return processes.		

INSTITUTION NAME (E.G., NAME OF HOSPITAL, HEALTH CARE ORGANISATION)	
INSTITUTION ADDRESS	PHONE/FAX/E-MAIL
FORM COMPLETED BY	STAMP
PRINT NAME: _____	
SIGNATURE: _____	
DATE	

Appendix 2: Product in scope of EIF-000523

Product Code	Batch	UDI Information (Unit Level)
112482065LA	KME21E0467	(01)14026704629142(17)260428(10)KME21E0467
112482070LA	KME21E0450	(01)14026704629159(17)260428(10)KME21E0450
112482075LA	KME21E0372	(01)14026704629166(17)260428(10)KME21E0372
112482075LA	KME21E0378	(01)14026704629166(17)260428(10)KME21E0378
112482075LA	KME21E0451	(01)14026704629166(17)260428(10)KME21E0451
112482075LA	KME21E0489	(01)14026704629166(17)260428(10)KME21E0489
112482080LA	KME21E0381	(01)14026704629173(17)260428(10)KME21E0381
112482085LA	KME21F0101	(01)14026704629180(17)260528(10)KME21F0101
112482085LA	KME21E0779	(01)14026704629180(17)260428(10)KME21E0779
112482085LA	KME21E0956	(01)14026704629180(17)260428(10)KME21E0956
112482085LA	KME21F0520	(01)14026704629180(17)260528(10)KME21F0520
112482085LA	KME21F0637	(01)14026704629180(17)260528(10)KME21F0637
112482095LA	KME21E0742	(01)14026704629203(17)260428(10)KME21E0742
112482100LA	KME21D1667	(01)14026704629210(17)260328(10)KME21D1667
112080-000065	KME21D2250	(01)14026704195937(17)260328(10)KME21D2250
112080-000080	KME21E1336	(01)14026704195968(17)260428(10)KME21E1336
112082-000050	KME21F0726	(01)14026704196019(17)260528(10)KME21F0726
112082-000055	KME21E1231	(01)14026704196026(17)260428(10)KME21E1231
112082-000055	KME21F0945	(01)14026704196026(17)260528(10)KME21F0945
112082-000060	KME21E1395	(01)14026704196033(17)260428(10)KME21E1395
112082-000065	KME21F0383	(01)14026704196040(17)260528(10)KME21F0383
112082-000065	KME21F0833	(01)14026704196040(17)260528(10)KME21F0833
112082-000065	KME21F0931	(01)14026704196040(17)260528(10)KME21F0931
112082-000070	KME21E0582	(01)14026704196057(17)260428(10)KME21E0582
112082-000070	KME21E1267	(01)14026704196057(17)260428(10)KME21E1267
112082-000070	KME21E1350	(01)14026704196057(17)260428(10)KME21E1350
112082-000070	KME21F0435	(01)14026704196057(17)260528(10)KME21F0435
112082-000070	KME21E0683	(01)14026704196057(17)260428(10)KME21E0683
112082-000070	KME21E1183	(01)14026704196057(17)260428(10)KME21E1183
112082-000070	KME21E1266	(01)14026704196057(17)260428(10)KME21E1266
112082-000070	KME21E1289	(01)14026704196057(17)260428(10)KME21E1289
112082-000070	KME21F0008	(01)14026704196057(17)260528(10)KME21F0008
112082-000070	KME21F0836	(01)14026704196057(17)260528(10)KME21F0836
112082-000075	KME21D1576	(01)14026704196064(17)260328(10)KME21D1576
112082-000075	KME21E1185	(01)14026704196064(17)260428(10)KME21E1185
112082-000075	KME21E0686	(01)14026704196064(17)260428(10)KME21E0686
112082-000075	KME21F0438	(01)14026704196064(17)260528(10)KME21F0438
112082-000075	KME21D2439	(01)14026704196064(17)260328(10)KME21D2439
112082-000075	KME21E0521	(01)14026704196064(17)260428(10)KME21E0521
112082-000075	KME21E0606	(01)14026704196064(17)260428(10)KME21E0606
112082-000075	KME21F0319	(01)14026704196064(17)260528(10)KME21F0319
112082-000080	KME21D1578	(01)14026704196071(17)260328(10)KME21D1578
112082-000080	KME21F0651	(01)14026704196071(17)260528(10)KME21F0651
112082-000080	KME21E1453	(01)14026704196071(17)260428(10)KME21E1453
112082-000085	KME21E0900	(01)14026704196088(17)260428(10)KME21E0900
112082-000085	KME21F0646	(01)14026704196088(17)260528(10)KME21F0646
112480-000050	KME21D1824	(01)14026704196200(17)260328(10)KME21D1824
112480-000050	KME21F0153	(01)14026704196200(17)260528(10)KME21F0153
112480-000050	KME21F0241	(01)14026704196200(17)260528(10)KME21F0241
112480-000050	KME21F0354	(01)14026704196200(17)260528(10)KME21F0354
112480-000050	KME21F0542	(01)14026704196200(17)260528(10)KME21F0542
112480-000050	KME21F0425	(01)14026704196200(17)260528(10)KME21F0425

Appendix 2: Product in scope of EIF-000523

Product Code	Batch	UDI Information (Unit Level)
112480-000050	KME21F0714	(01)14026704196200(17)260528(10)KME21F0714
112480-000055	KME21F0437	(01)14026704196217(17)260528(10)KME21F0437
112480-000055	KME21D2253	(01)14026704196217(17)260328(10)KME21D2253
112480-000055	KME21D2511	(01)14026704196217(17)260328(10)KME21D2511
112480-000055	KME21D2669	(01)14026704196217(17)260328(10)KME21D2669
112480-000055	KME21D2826	(01)14026704196217(17)260328(10)KME21D2826
112480-000055	KME21E0798	(01)14026704196217(17)260428(10)KME21E0798
112480-000055	KME21D2901	(01)14026704196217(17)260328(10)KME21D2901
112480-000055	KME21E0585	(01)14026704196217(17)260428(10)KME21E0585
112480-000055	KME21F0436	(01)14026704196217(17)260528(10)KME21F0436
112480-000060	KME21D1703	(01)14026704196224(17)260328(10)KME21D1703
112480-000060	KME21E1148	(01)14026704196224(17)260428(10)KME21E1148
112480-000060	KME21D1796	(01)14026704196224(17)260328(10)KME21D1796
112480-000060	KME21D1823	(01)14026704196224(17)260328(10)KME21D1823
112480-000060	KME21D1958	(01)14026704196224(17)260328(10)KME21D1958
112480-000060	KME21D2255	(01)14026704196224(17)260328(10)KME21D2255
112480-000060	KME21E0803	(01)14026704196224(17)260428(10)KME21E0803
112480-000060	KME21E0865	(01)14026704196224(17)260428(10)KME21E0865
112480-000060	KME21E1222	(01)14026704196224(17)260428(10)KME21E1222
112480-000060	KME21E1100	(01)14026704196224(17)260428(10)KME21E1100
112480-000060	KME21F0297	(01)14026704196224(17)260528(10)KME21F0297
112480-000065	KME21D2714	(01)14026704196231(17)260328(10)KME21D2714
112480-000065	KME21F0519	(01)14026704196231(17)260528(10)KME21F0519
112480-000065	KME21F0780	(01)14026704196231(17)260528(10)KME21F0780
112480-000065	KME21F0009	(01)14026704196231(17)260528(10)KME21F0009
112480-000065	KME21D2776	(01)14026704196231(17)260328(10)KME21D2776
112480-000065	KME21D2902	(01)14026704196231(17)260328(10)KME21D2902
112480-000065	KME21E0802	(01)14026704196231(17)260428(10)KME21E0802
112480-000065	KME21E1181	(01)14026704196231(17)260428(10)KME21E1181
112480-000065	KME21E1347	(01)14026704196231(17)260428(10)KME21E1347
112480-000065	KME21E1408	(01)14026704196231(17)260428(10)KME21E1408
112480-000065	KME21F0010	(01)14026704196231(17)260528(10)KME21F0010
112480-000065	KME21E1499	(01)14026704196231(17)260428(10)KME21E1499
112480-000065	KME21F0427	(01)14026704196231(17)260528(10)KME21F0427
112480-000065	KME21F0690	(01)14026704196231(17)260528(10)KME21F0690
112480-000065	KME21F0830	(01)14026704196231(17)260528(10)KME21F0830
112480-000070	KME21D2010	(01)14026704196248(17)260328(10)KME21D2010
112480-000070	KME21F0116	(01)14026704196248(17)260528(10)KME21F0116
112480-000070	KME21E1368	(01)14026704196248(17)260428(10)KME21E1368
112480-000070	KME21D1574	(01)14026704196248(17)260328(10)KME21D1574
112480-000070	KME21D1701	(01)14026704196248(17)260328(10)KME21D1701
112480-000070	KME21D1797	(01)14026704196248(17)260328(10)KME21D1797
112480-000070	KME21D1826	(01)14026704196248(17)260328(10)KME21D1826
112480-000070	KME21D1879	(01)14026704196248(17)260328(10)KME21D1879
112480-000070	KME21D2019	(01)14026704196248(17)260328(10)KME21D2019
112480-000070	KME21D2258	(01)14026704196248(17)260328(10)KME21D2258
112480-000070	KME21D2508	(01)14026704196248(17)260328(10)KME21D2508
112480-000070	KME21D2538	(01)14026704196248(17)260328(10)KME21D2538
112480-000070	KME21D2644	(01)14026704196248(17)260328(10)KME21D2644
112480-000070	KME21D2596	(01)14026704196248(17)260328(10)KME21D2596
112480-000070	KME21D2715	(01)14026704196248(17)260328(10)KME21D2715
112480-000070	KME21D2753	(01)14026704196248(17)260328(10)KME21D2753

Appendix 2: Product in scope of EIF-000523

Product Code	Batch	UDI Information (Unit Level)
112480-000070	KME21D2827	(01)14026704196248(17)260328(10)KME21D2827
112480-000070	KME21D2904	(01)14026704196248(17)260328(10)KME21D2904
112480-000070	KME21E0019	(01)14026704196248(17)260428(10)KME21E0019
112480-000070	KME21E0499	(01)14026704196248(17)260428(10)KME21E0499
112480-000070	KME21E0586	(01)14026704196248(17)260428(10)KME21E0586
112480-000070	KME21E1500	(01)14026704196248(17)260428(10)KME21E1500
112480-000070	KME21F0180	(01)14026704196248(17)260528(10)KME21F0180
112480-000070	KME21F0278	(01)14026704196248(17)260528(10)KME21F0278
112480-000070	KME21F0413	(01)14026704196248(17)260528(10)KME21F0413
112480-000070	KME21F0544	(01)14026704196248(17)260528(10)KME21F0544
112480-000070	KME21F0713	(01)14026704196248(17)260528(10)KME21F0713
112480-000075	KME21D1827	(01)14026704196255(17)260328(10)KME21D1827
112480-000075	KME21E1282	(01)14026704196255(17)260428(10)KME21E1282
112480-000075	KME21F0183	(01)14026704196255(17)260528(10)KME21F0183
112480-000075	KME21F0777	(01)14026704196255(17)260528(10)KME21F0777
112480-000075	KME21D2700	(01)14026704196255(17)260328(10)KME21D2700
112480-000075	KME21D2828	(01)14026704196255(17)260328(10)KME21D2828
112480-000075	KME21E0331	(01)14026704196255(17)260428(10)KME21E0331
112480-000075	KME21E0744	(01)14026704196255(17)260428(10)KME21E0744
112480-000075	KME21E0903	(01)14026704196255(17)260428(10)KME21E0903
112480-000075	KME21E0928	(01)14026704196255(17)260428(10)KME21E0928
112480-000075	KME21E1014	(01)14026704196255(17)260428(10)KME21E1014
112480-000075	KME21E1015	(01)14026704196255(17)260428(10)KME21E1015
112480-000075	KME21E1047	(01)14026704196255(17)260428(10)KME21E1047
112480-000075	KME21E1108	(01)14026704196255(17)260428(10)KME21E1108
112480-000075	KME21E1146	(01)14026704196255(17)260428(10)KME21E1146
112480-000075	KME21E1168	(01)14026704196255(17)260428(10)KME21E1168
112480-000075	KME21E1190	(01)14026704196255(17)260428(10)KME21E1190
112480-000075	KME21E1192	(01)14026704196255(17)260428(10)KME21E1192
112480-000075	KME21E1209	(01)14026704196255(17)260428(10)KME21E1209
112480-000075	KME21E1263	(01)14026704196255(17)260428(10)KME21E1263
112480-000075	KME21E1268	(01)14026704196255(17)260428(10)KME21E1268
112480-000075	KME21E1348	(01)14026704196255(17)260428(10)KME21E1348
112480-000075	KME21E1394	(01)14026704196255(17)260428(10)KME21E1394
112480-000075	KME21F0100	(01)14026704196255(17)260528(10)KME21F0100
112480-000075	KME21E1501	(01)14026704196255(17)260428(10)KME21E1501
112480-000075	KME21F0012	(01)14026704196255(17)260528(10)KME21F0012
112480-000075	KME21F0237	(01)14026704196255(17)260528(10)KME21F0237
112480-000075	KME21F0296	(01)14026704196255(17)260528(10)KME21F0296
112480-000075	KME21F0433	(01)14026704196255(17)260528(10)KME21F0433
112480-000075	KME21F0533	(01)14026704196255(17)260528(10)KME21F0533
112480-000075	KME21F0712	(01)14026704196255(17)260528(10)KME21F0712
112480-000075	KME21F0875	(01)14026704196255(17)260528(10)KME21F0875
112480-000075	KME21F0941	(01)14026704196255(17)260528(10)KME21F0941
112480-000075	KME21F0970	(01)14026704196255(17)260528(10)KME21F0970
112480-000075	KME21F1005	(01)14026704196255(17)260528(10)KME21F1005
112480-000075	KME21F1066	(01)14026704196255(17)260528(10)KME21F1066
112480-000075	KME21F1108	(01)14026704196255(17)260528(10)KME21F1108
112480-000075	KME21F1141	(01)14026704196255(17)260528(10)KME21F1141
112480-000080	KME21D1880	(01)14026704196262(17)260328(10)KME21D1880
112480-000080	KME21F1139	(01)14026704196262(17)260528(10)KME21F1139
112480-000080	KME21D2256	(01)14026704196262(17)260328(10)KME21D2256

Appendix 2: Product in scope of EIF-000523

Product Code	Batch	UDI Information (Unit Level)
112480-000080	KME21D2788	(01)14026704196262(17)260328(10)KME21D2788
112480-000080	KME21D2789	(01)14026704196262(17)260328(10)KME21D2789
112480-000080	KME21D2906	(01)14026704196262(17)260328(10)KME21D2906
112480-000080	KME21D2938	(01)14026704196262(17)260328(10)KME21D2938
112480-000080	KME21E0020	(01)14026704196262(17)260428(10)KME21E0020
112480-000080	KME21E0500	(01)14026704196262(17)260428(10)KME21E0500
112480-000080	KME21F0636	(01)14026704196262(17)260528(10)KME21F0636
112480-000080	KME21E0468	(01)14026704196262(17)260428(10)KME21E0468
112480-000080	KME21E0592	(01)14026704196262(17)260428(10)KME21E0592
112480-000080	KME21E1223	(01)14026704196262(17)260428(10)KME21E1223
112480-000080	KME21E1346	(01)14026704196262(17)260428(10)KME21E1346
112480-000080	KME21E1371	(01)14026704196262(17)260428(10)KME21E1371
112480-000080	KME21E1476	(01)14026704196262(17)260428(10)KME21E1476
112480-000080	KME21F0151	(01)14026704196262(17)260528(10)KME21F0151
112480-000080	KME21F0011	(01)14026704196262(17)260528(10)KME21F0011
112480-000080	KME21F0186	(01)14026704196262(17)260528(10)KME21F0186
112480-000080	KME21F0276	(01)14026704196262(17)260528(10)KME21F0276
112480-000080	KME21F0333	(01)14026704196262(17)260528(10)KME21F0333
112480-000080	KME21F0386	(01)14026704196262(17)260528(10)KME21F0386
112480-000080	KME21F0419	(01)14026704196262(17)260528(10)KME21F0419
112480-000080	KME21F0683	(01)14026704196262(17)260528(10)KME21F0683
112480-000080	KME21F0705	(01)14026704196262(17)260528(10)KME21F0705
112480-000080	KME21F0973	(01)14026704196262(17)260528(10)KME21F0973
112480-000080	KME21F1074	(01)14026704196262(17)260528(10)KME21F1074
112480-000080	KME21F1159	(01)14026704196262(17)260528(10)KME21F1159
112480-000085	KME21E1189	(01)14026704196279(17)260428(10)KME21E1189
112480-000085	KME21E1443	(01)14026704196279(17)260428(10)KME21E1443
112480-000085	KME21F0235	(01)14026704196279(17)260528(10)KME21F0235
112480-000090	KME21D1982	(01)14026704196286(17)260328(10)KME21D1982
112480-000095	KME21E1409	(01)14026704196293(17)260428(10)KME21E1409
112480-000095	KME21E1134	(01)14026704196293(17)260428(10)KME21E1134
112481-000055	KME21D1840	(01)14026704196330(17)260328(10)KME21D1840
112481-000055	KME21F0812	(01)14026704196330(17)260528(10)KME21F0812
112481-000060	KME21E1386	(01)14026704196347(17)260428(10)KME21E1386
112481-000070	KME21E0866	(01)14026704196361(17)260428(10)KME21E0866
112481-000070	KME21E1094	(01)14026704196361(17)260428(10)KME21E1094
112481-000070	KME21F0793	(01)14026704196361(17)260528(10)KME21F0793
112481-000070	KME21D1468	(01)14026704196361(17)260328(10)KME21D1468
112481-000070	KME21D1562	(01)14026704196361(17)260328(10)KME21D1562
112481-000075	KME21D1146	(01)14026704196378(17)260328(10)KME21D1146
112481-000075	KME21D2825	(01)14026704196378(17)260328(10)KME21D2825
112481-000075	KME21E1092	(01)14026704196378(17)260428(10)KME21E1092
112481-000075	KME21F0794	(01)14026704196378(17)260528(10)KME21F0794
112481-000075	KME21D1686	(01)14026704196378(17)260328(10)KME21D1686
112481-000080	KME21F1073	(01)14026704196385(17)260528(10)KME21F1073
112481-000080	KME21D1153	(01)14026704196385(17)260328(10)KME21D1153
112481-000080	KME21D1296	(01)14026704196385(17)260328(10)KME21D1296
112481-000080	KME21D1560	(01)14026704196385(17)260328(10)KME21D1560
112481-000080	KME21F0895	(01)14026704196385(17)260528(10)KME21F0895
112481-000080	KME21F0995	(01)14026704196385(17)260528(10)KME21F0995
112481-000085	KME21E1119	(01)14026704196392(17)260428(10)KME21E1119
112481-000085	KME21D0707	(01)14026704196392(17)260328(10)KME21D0707

Appendix 2: Product in scope of EIF-000523

Product Code	Batch	UDI Information (Unit Level)
112481-000085	KME21D1685	(01)14026704196392(17)260328(10)KME21D1685
112481-000090	KME21E1115	(01)14026704196408(17)260428(10)KME21E1115
112481-000090	KME21E1116	(01)14026704196408(17)260428(10)KME21E1116
112482-000050	KME21E0688	(01)14026704196460(17)260428(10)KME21E0688
112482-000050	KME21F0720	(01)14026704196460(17)260528(10)KME21F0720
112482-000050	KME21D1704	(01)14026704196460(17)260328(10)KME21D1704
112482-000050	KME21D2260	(01)14026704196460(17)260328(10)KME21D2260
112482-000050	KME21D2534	(01)14026704196460(17)260328(10)KME21D2534
112482-000050	KME21D2499	(01)14026704196460(17)260328(10)KME21D2499
112482-000050	KME21D2701	(01)14026704196460(17)260328(10)KME21D2701
112482-000050	KME21D2562	(01)14026704196460(17)260328(10)KME21D2562
112482-000050	KME21D2790	(01)14026704196460(17)260328(10)KME21D2790
112482-000050	KME21D2842	(01)14026704196460(17)260328(10)KME21D2842
112482-000050	KME21E0382	(01)14026704196460(17)260428(10)KME21E0382
112482-000050	KME21E0465	(01)14026704196460(17)260428(10)KME21E0465
112482-000050	KME21E0801	(01)14026704196460(17)260428(10)KME21E0801
112482-000050	KME21F0236	(01)14026704196460(17)260528(10)KME21F0236
112482-000050	KME21E1502	(01)14026704196460(17)260428(10)KME21E1502
112482-000050	KME21F0355	(01)14026704196460(17)260528(10)KME21F0355
112482-000050	KME21F0420	(01)14026704196460(17)260528(10)KME21F0420
112482-000050	KME21F0518	(01)14026704196460(17)260528(10)KME21F0518
112482-000050	KME21F0536	(01)14026704196460(17)260528(10)KME21F0536
112482-000050	KME21F0593	(01)14026704196460(17)260528(10)KME21F0593
112482-000050	KME21F1004	(01)14026704196460(17)260528(10)KME21F1004
112482-000055	KME21F0718	(01)14026704196477(17)260528(10)KME21F0718
112482-000055	KME21F0832	(01)14026704196477(17)260528(10)KME21F0832
112482-000055	KME21D1587	(01)14026704196477(17)260328(10)KME21D1587
112482-000055	KME21D1966	(01)14026704196477(17)260328(10)KME21D1966
112482-000055	KME21D2261	(01)14026704196477(17)260328(10)KME21D2261
112482-000055	KME21D2385	(01)14026704196477(17)260328(10)KME21D2385
112482-000055	KME21D2417	(01)14026704196477(17)260328(10)KME21D2417
112482-000055	KME21E0219	(01)14026704196477(17)260428(10)KME21E0219
112482-000055	KME21E0593	(01)14026704196477(17)260428(10)KME21E0593
112482-000055	KME21E0637	(01)14026704196477(17)260428(10)KME21E0637
112482-000055	KME21E0689	(01)14026704196477(17)260428(10)KME21E0689
112482-000055	KME21E0800	(01)14026704196477(17)260428(10)KME21E0800
112482-000055	KME21E0882	(01)14026704196477(17)260428(10)KME21E0882
112482-000055	KME21F1068	(01)14026704196477(17)260528(10)KME21F1068
112482-000055	KME21F1069	(01)14026704196477(17)260528(10)KME21F1069
112482-000060	KME21D2898	(01)14026704196484(17)260328(10)KME21D2898
112482-000060	KME21E0682	(01)14026704196484(17)260428(10)KME21E0682
112482-000060	KME21F0543	(01)14026704196484(17)260528(10)KME21F0543
112482-000060	KME21F0806	(01)14026704196484(17)260528(10)KME21F0806
112482-000060	KME21D2531	(01)14026704196484(17)260328(10)KME21D2531
112482-000060	KME21D2647	(01)14026704196484(17)260328(10)KME21D2647
112482-000060	KME21D2751	(01)14026704196484(17)260328(10)KME21D2751
112482-000060	KME21D2829	(01)14026704196484(17)260328(10)KME21D2829
112482-000060	KME21D2927	(01)14026704196484(17)260328(10)KME21D2927
112482-000060	KME21E0018	(01)14026704196484(17)260428(10)KME21E0018
112482-000060	KME21E0220	(01)14026704196484(17)260428(10)KME21E0220
112482-000060	KME21E0549	(01)14026704196484(17)260428(10)KME21E0549
112482-000060	KME21E0327	(01)14026704196484(17)260428(10)KME21E0327

Appendix 2: Product in scope of EIF-000523

Product Code	Batch	UDI Information (Unit Level)
112482-000060	KME21E0444	(01)14026704196484(17)260428(10)KME21E0444
112482-000060	KME21E0550	(01)14026704196484(17)260428(10)KME21E0550
112482-000060	KME21E0466	(01)14026704196484(17)260428(10)KME21E0466
112482-000060	KME21E0931	(01)14026704196484(17)260428(10)KME21E0931
112482-000060	KME21E1149	(01)14026704196484(17)260428(10)KME21E1149
112482-000060	KME21F0473	(01)14026704196484(17)260528(10)KME21F0473
112482-000060	KME21F0592	(01)14026704196484(17)260528(10)KME21F0592
112482-000060	KME21F0649	(01)14026704196484(17)260528(10)KME21F0649
112482-000060	KME21F0807	(01)14026704196484(17)260528(10)KME21F0807
112482-000060	KME21F0929	(01)14026704196484(17)260528(10)KME21F0929
112482-000060	KME21F0979	(01)14026704196484(17)260528(10)KME21F0979
112482-000060	KME21F1075	(01)14026704196484(17)260528(10)KME21F1075
112482-000060	KME21F1113	(01)14026704196484(17)260528(10)KME21F1113
112482-000065	KME21D2900	(01)14026704196491(17)260328(10)KME21D2900
112482-000065	KME21F0552	(01)14026704196491(17)260528(10)KME21F0552
112482-000065	KME21F0729	(01)14026704196491(17)260528(10)KME21F0729
112482-000065	KME21D1573	(01)14026704196491(17)260328(10)KME21D1573
112482-000065	KME21D1650	(01)14026704196491(17)260328(10)KME21D1650
112482-000065	KME21D1793	(01)14026704196491(17)260328(10)KME21D1793
112482-000065	KME21D1944	(01)14026704196491(17)260328(10)KME21D1944
112482-000065	KME21D1981	(01)14026704196491(17)260328(10)KME21D1981
112482-000065	KME21D2384	(01)14026704196491(17)260328(10)KME21D2384
112482-000065	KME21D2418	(01)14026704196491(17)260328(10)KME21D2418
112482-000065	KME21D2526	(01)14026704196491(17)260328(10)KME21D2526
112482-000065	KME21D2611	(01)14026704196491(17)260328(10)KME21D2611
112482-000065	KME21D2598	(01)14026704196491(17)260328(10)KME21D2598
112482-000065	KME21D2883	(01)14026704196491(17)260328(10)KME21D2883
112482-000065	KME21E0115	(01)14026704196491(17)260428(10)KME21E0115
112482-000065	KME21E0221	(01)14026704196491(17)260428(10)KME21E0221
112482-000065	KME21E0520	(01)14026704196491(17)260428(10)KME21E0520
112482-000065	KME21E0594	(01)14026704196491(17)260428(10)KME21E0594
112482-000065	KME21E0551	(01)14026704196491(17)260428(10)KME21E0551
112482-000065	KME21E0684	(01)14026704196491(17)260428(10)KME21E0684
112482-000065	KME21E0799	(01)14026704196491(17)260428(10)KME21E0799
112482-000065	KME21E0806	(01)14026704196491(17)260428(10)KME21E0806
112482-000065	KME21E1225	(01)14026704196491(17)260428(10)KME21E1225
112482-000065	KME21F0635	(01)14026704196491(17)260528(10)KME21F0635
112482-000065	KME21F0584	(01)14026704196491(17)260528(10)KME21F0584
112482-000065	KME21F1002	(01)14026704196491(17)260528(10)KME21F1002
112482-000070	KME21D2699	(01)14026704196507(17)260328(10)KME21D2699
112482-000070	KME21E1370	(01)14026704196507(17)260428(10)KME21E1370
112482-000070	KME21F0316	(01)14026704196507(17)260528(10)KME21F0316
112482-000070	KME21E1140	(01)14026704196507(17)260428(10)KME21E1140
112482-000070	KME21F1077	(01)14026704196507(17)260528(10)KME21F1077
112482-000070	KME21D1567	(01)14026704196507(17)260328(10)KME21D1567
112482-000070	KME21D1651	(01)14026704196507(17)260328(10)KME21D1651
112482-000070	KME21D1702	(01)14026704196507(17)260328(10)KME21D1702
112482-000070	KME21D1799	(01)14026704196507(17)260328(10)KME21D1799
112482-000070	KME21D1945	(01)14026704196507(17)260328(10)KME21D1945
112482-000070	KME21D1952	(01)14026704196507(17)260328(10)KME21D1952
112482-000070	KME21D2009	(01)14026704196507(17)260328(10)KME21D2009
112482-000070	KME21D2016	(01)14026704196507(17)260328(10)KME21D2016

Appendix 2: Product in scope of EIF-000523

Product Code	Batch	UDI Information (Unit Level)
112482-000070	KME21D2263	(01)14026704196507(17)260328(10)KME21D2263
112482-000070	KME21D2640	(01)14026704196507(17)260328(10)KME21D2640
112482-000070	KME21D2383	(01)14026704196507(17)260328(10)KME21D2383
112482-000070	KME21D2528	(01)14026704196507(17)260328(10)KME21D2528
112482-000070	KME21D2497	(01)14026704196507(17)260328(10)KME21D2497
112482-000070	KME21D2535	(01)14026704196507(17)260328(10)KME21D2535
112482-000070	KME21D2612	(01)14026704196507(17)260328(10)KME21D2612
112482-000070	KME21D2599	(01)14026704196507(17)260328(10)KME21D2599
112482-000070	KME21D2750	(01)14026704196507(17)260328(10)KME21D2750
112482-000070	KME21D2843	(01)14026704196507(17)260328(10)KME21D2843
112482-000070	KME21D2798	(01)14026704196507(17)260328(10)KME21D2798
112482-000070	KME21D2800	(01)14026704196507(17)260328(10)KME21D2800
112482-000070	KME21D2865	(01)14026704196507(17)260328(10)KME21D2865
112482-000070	KME21D2899	(01)14026704196507(17)260328(10)KME21D2899
112482-000070	KME21E0015	(01)14026704196507(17)260428(10)KME21E0015
112482-000070	KME21E0282	(01)14026704196507(17)260428(10)KME21E0282
112482-000070	KME21E0328	(01)14026704196507(17)260428(10)KME21E0328
112482-000070	KME21E0530	(01)14026704196507(17)260428(10)KME21E0530
112482-000070	KME21E0470	(01)14026704196507(17)260428(10)KME21E0470
112482-000070	KME21E0552	(01)14026704196507(17)260428(10)KME21E0552
112482-000070	KME21E0563	(01)14026704196507(17)260428(10)KME21E0563
112482-000070	KME21E1040	(01)14026704196507(17)260428(10)KME21E1040
112482-000070	KME21E0749	(01)14026704196507(17)260428(10)KME21E0749
112482-000070	KME21E0807	(01)14026704196507(17)260428(10)KME21E0807
112482-000070	KME21E0850	(01)14026704196507(17)260428(10)KME21E0850
112482-000070	KME21E0863	(01)14026704196507(17)260428(10)KME21E0863
112482-000070	KME21E0878	(01)14026704196507(17)260428(10)KME21E0878
112482-000070	KME21E0913	(01)14026704196507(17)260428(10)KME21E0913
112482-000070	KME21E0915	(01)14026704196507(17)260428(10)KME21E0915
112482-000070	KME21E0982	(01)14026704196507(17)260428(10)KME21E0982
112482-000070	KME21E1006	(01)14026704196507(17)260428(10)KME21E1006
112482-000070	KME21E1007	(01)14026704196507(17)260428(10)KME21E1007
112482-000070	KME21E1086	(01)14026704196507(17)260428(10)KME21E1086
112482-000070	KME21E1087	(01)14026704196507(17)260428(10)KME21E1087
112482-000070	KME21E1101	(01)14026704196507(17)260428(10)KME21E1101
112482-000070	KME21E1188	(01)14026704196507(17)260428(10)KME21E1188
112482-000070	KME21E1166	(01)14026704196507(17)260428(10)KME21E1166
112482-000070	KME21E1198	(01)14026704196507(17)260428(10)KME21E1198
112482-000070	KME21E1208	(01)14026704196507(17)260428(10)KME21E1208
112482-000070	KME21E1317	(01)14026704196507(17)260428(10)KME21E1317
112482-000070	KME21E1334	(01)14026704196507(17)260428(10)KME21E1334
112482-000070	KME21E1345	(01)14026704196507(17)260428(10)KME21E1345
112482-000070	KME21E1369	(01)14026704196507(17)260428(10)KME21E1369
112482-000070	KME21E1441	(01)14026704196507(17)260428(10)KME21E1441
112482-000070	KME21E1477	(01)14026704196507(17)260428(10)KME21E1477
112482-000070	KME21F0126	(01)14026704196507(17)260528(10)KME21F0126
112482-000070	KME21F0014	(01)14026704196507(17)260528(10)KME21F0014
112482-000070	KME21F0138	(01)14026704196507(17)260528(10)KME21F0138
112482-000070	KME21F0148	(01)14026704196507(17)260528(10)KME21F0148
112482-000070	KME21F0181	(01)14026704196507(17)260528(10)KME21F0181
112482-000070	KME21F0214	(01)14026704196507(17)260528(10)KME21F0214
112482-000070	KME21F0234	(01)14026704196507(17)260528(10)KME21F0234

Appendix 2: Product in scope of EIF-000523

Product Code	Batch	UDI Information (Unit Level)
112482-000070	KME21F0308	(01)14026704196507(17)260528(10)KME21F0308
112482-000070	KME21F0942	(01)14026704196507(17)260528(10)KME21F0942
112482-000070	KME21F0418	(01)14026704196507(17)260528(10)KME21F0418
112482-000070	KME21F0515	(01)14026704196507(17)260528(10)KME21F0515
112482-000070	KME21F0640	(01)14026704196507(17)260528(10)KME21F0640
112482-000070	KME21F0688	(01)14026704196507(17)260528(10)KME21F0688
112482-000070	KME21F0779	(01)14026704196507(17)260528(10)KME21F0779
112482-000070	KME21F0831	(01)14026704196507(17)260528(10)KME21F0831
112482-000070	KME21F0856	(01)14026704196507(17)260528(10)KME21F0856
112482-000070	KME21F0873	(01)14026704196507(17)260528(10)KME21F0873
112482-000070	KME21F0971	(01)14026704196507(17)260528(10)KME21F0971
112482-000070	KME21F1107	(01)14026704196507(17)260528(10)KME21F1107
112482-000070	KME21F1142	(01)14026704196507(17)260528(10)KME21F1142
112482-000070	KME21F1157	(01)14026704196507(17)260528(10)KME21F1157
112482-000075	KME21D2667	(01)14026704196514(17)260328(10)KME21D2667
112482-000075	KME21E0491	(01)14026704196514(17)260428(10)KME21E0491
112482-000075	KME21F0697	(01)14026704196514(17)260528(10)KME21F0697
112482-000075	KME21F0531	(01)14026704196514(17)260528(10)KME21F0531
112482-000075	KME21D2519	(01)14026704196514(17)260328(10)KME21D2519
112482-000075	KME21D1579	(01)14026704196514(17)260328(10)KME21D1579
112482-000075	KME21D1648	(01)14026704196514(17)260328(10)KME21D1648
112482-000075	KME21D1795	(01)14026704196514(17)260328(10)KME21D1795
112482-000075	KME21D1881	(01)14026704196514(17)260328(10)KME21D1881
112482-000075	KME21D1891	(01)14026704196514(17)260328(10)KME21D1891
112482-000075	KME21D1960	(01)14026704196514(17)260328(10)KME21D1960
112482-000075	KME21D2013	(01)14026704196514(17)260328(10)KME21D2013
112482-000075	KME21D2017	(01)14026704196514(17)260328(10)KME21D2017
112482-000075	KME21D2264	(01)14026704196514(17)260328(10)KME21D2264
112482-000075	KME21D2404	(01)14026704196514(17)260328(10)KME21D2404
112482-000075	KME21D2420	(01)14026704196514(17)260328(10)KME21D2420
112482-000075	KME21D2440	(01)14026704196514(17)260328(10)KME21D2440
112482-000075	KME21D2503	(01)14026704196514(17)260328(10)KME21D2503
112482-000075	KME21D2518	(01)14026704196514(17)260328(10)KME21D2518
112482-000075	KME21D2642	(01)14026704196514(17)260328(10)KME21D2642
112482-000075	KME21D2564	(01)14026704196514(17)260328(10)KME21D2564
112482-000075	KME21D2600	(01)14026704196514(17)260328(10)KME21D2600
112482-000075	KME21D2794	(01)14026704196514(17)260328(10)KME21D2794
112482-000075	KME21D2926	(01)14026704196514(17)260328(10)KME21D2926
112482-000075	KME21E0017	(01)14026704196514(17)260428(10)KME21E0017
112482-000075	KME21E0222	(01)14026704196514(17)260428(10)KME21E0222
112482-000075	KME21E0263	(01)14026704196514(17)260428(10)KME21E0263
112482-000075	KME21E0325	(01)14026704196514(17)260428(10)KME21E0325
112482-000075	KME21E0490	(01)14026704196514(17)260428(10)KME21E0490
112482-000075	KME21E0895	(01)14026704196514(17)260428(10)KME21E0895
112482-000075	KME21E0662	(01)14026704196514(17)260428(10)KME21E0662
112482-000075	KME21E1284	(01)14026704196514(17)260428(10)KME21E1284
112482-000075	KME21E1449	(01)14026704196514(17)260428(10)KME21E1449
112482-000075	KME21F0103	(01)14026704196514(17)260528(10)KME21F0103
112482-000075	KME21F0179	(01)14026704196514(17)260528(10)KME21F0179
112482-000075	KME21F0215	(01)14026704196514(17)260528(10)KME21F0215
112482-000075	KME21F0300	(01)14026704196514(17)260528(10)KME21F0300
112482-000075	KME21F0356	(01)14026704196514(17)260528(10)KME21F0356

Appendix 2: Product in scope of EIF-000523

Product Code	Batch	UDI Information (Unit Level)
112482-000075	KME21F0382	(01)14026704196514(17)260528(10)KME21F0382
112482-000075	KME21F0469	(01)14026704196514(17)260528(10)KME21F0469
112482-000075	KME21F0472	(01)14026704196514(17)260528(10)KME21F0472
112482-000075	KME21F0517	(01)14026704196514(17)260528(10)KME21F0517
112482-000075	KME21F0648	(01)14026704196514(17)260528(10)KME21F0648
112482-000075	KME21F0587	(01)14026704196514(17)260528(10)KME21F0587
112482-000075	KME21F0585	(01)14026704196514(17)260528(10)KME21F0585
112482-000075	KME21F0588	(01)14026704196514(17)260528(10)KME21F0588
112482-000075	KME21F0684	(01)14026704196514(17)260528(10)KME21F0684
112482-000075	KME21F0685	(01)14026704196514(17)260528(10)KME21F0685
112482-000075	KME21F0858	(01)14026704196514(17)260528(10)KME21F0858
112482-000075	KME21F0857	(01)14026704196514(17)260528(10)KME21F0857
112482-000075	KME21F0871	(01)14026704196514(17)260528(10)KME21F0871
112482-000075	KME21F0918	(01)14026704196514(17)260528(10)KME21F0918
112482-000075	KME21F0972	(01)14026704196514(17)260528(10)KME21F0972
112482-000075	KME21F1001	(01)14026704196514(17)260528(10)KME21F1001
112482-000075	KME21F1031	(01)14026704196514(17)260528(10)KME21F1031
112482-000075	KME21F1057	(01)14026704196514(17)260528(10)KME21F1057
112482-000075	KME21F1112	(01)14026704196514(17)260528(10)KME21F1112
112482-000080	KME21D1825	(01)14026704196521(17)260328(10)KME21D1825
112482-000080	KME21E0016	(01)14026704196521(17)260428(10)KME21E0016
112482-000080	KME21F0514	(01)14026704196521(17)260528(10)KME21F0514
112482-000080	KME21F0776	(01)14026704196521(17)260528(10)KME21F0776
112482-000080	KME21F1030	(01)14026704196521(17)260528(10)KME21F1030
112482-000080	KME21D1882	(01)14026704196521(17)260328(10)KME21D1882
112482-000080	KME21D1663	(01)14026704196521(17)260328(10)KME21D1663
112482-000080	KME21D2441	(01)14026704196521(17)260328(10)KME21D2441
112482-000080	KME21D2506	(01)14026704196521(17)260328(10)KME21D2506
112482-000080	KME21D2498	(01)14026704196521(17)260328(10)KME21D2498
112482-000080	KME21D2520	(01)14026704196521(17)260328(10)KME21D2520
112482-000080	KME21D2646	(01)14026704196521(17)260328(10)KME21D2646
112482-000080	KME21D2566	(01)14026704196521(17)260328(10)KME21D2566
112482-000080	KME21D2717	(01)14026704196521(17)260328(10)KME21D2717
112482-000080	KME21D2791	(01)14026704196521(17)260328(10)KME21D2791
112482-000080	KME21D2844	(01)14026704196521(17)260328(10)KME21D2844
112482-000080	KME21D2925	(01)14026704196521(17)260328(10)KME21D2925
112482-000080	KME21E0501	(01)14026704196521(17)260428(10)KME21E0501
112482-000080	KME21E1038	(01)14026704196521(17)260428(10)KME21E1038
112482-000080	KME21E0638	(01)14026704196521(17)260428(10)KME21E0638
112482-000080	KME21E0748	(01)14026704196521(17)260428(10)KME21E0748
112482-000080	KME21E0883	(01)14026704196521(17)260428(10)KME21E0883
112482-000080	KME21E0935	(01)14026704196521(17)260428(10)KME21E0935
112482-000080	KME21E1210	(01)14026704196521(17)260428(10)KME21E1210
112482-000080	KME21E1265	(01)14026704196521(17)260428(10)KME21E1265
112482-000080	KME21F0650	(01)14026704196521(17)260528(10)KME21F0650
112482-000080	KME21E1269	(01)14026704196521(17)260428(10)KME21E1269
112482-000080	KME21E1318	(01)14026704196521(17)260428(10)KME21E1318
112482-000080	KME21E1387	(01)14026704196521(17)260428(10)KME21E1387
112482-000080	KME21E1442	(01)14026704196521(17)260428(10)KME21E1442
112482-000080	KME21E1454	(01)14026704196521(17)260428(10)KME21E1454
112482-000080	KME21F0110	(01)14026704196521(17)260528(10)KME21F0110
112482-000080	KME21F0007	(01)14026704196521(17)260528(10)KME21F0007

Appendix 2: Product in scope of EIF-000523

Product Code	Batch	UDI Information (Unit Level)
112482-000080	KME21F0139	(01)14026704196521(17)260528(10)KME21F0139
112482-000080	KME21F0149	(01)14026704196521(17)260528(10)KME21F0149
112482-000080	KME21F0349	(01)14026704196521(17)260528(10)KME21F0349
112482-000080	KME21F0381	(01)14026704196521(17)260528(10)KME21F0381
112482-000080	KME21F0476	(01)14026704196521(17)260528(10)KME21F0476
112482-000080	KME21F0502	(01)14026704196521(17)260528(10)KME21F0502
112482-000080	KME21F0532	(01)14026704196521(17)260528(10)KME21F0532
112482-000080	KME21F0687	(01)14026704196521(17)260528(10)KME21F0687
112482-000080	KME21F0717	(01)14026704196521(17)260528(10)KME21F0717
112482-000080	KME21F0774	(01)14026704196521(17)260528(10)KME21F0774
112482-000080	KME21F0804	(01)14026704196521(17)260528(10)KME21F0804
112482-000080	KME21F0826	(01)14026704196521(17)260528(10)KME21F0826
112482-000080	KME21F0880	(01)14026704196521(17)260528(10)KME21F0880
112482-000080	KME21F0967	(01)14026704196521(17)260528(10)KME21F0967
112482-000080	KME21F1110	(01)14026704196521(17)260528(10)KME21F1110
112482-000080	KME21F0968	(01)14026704196521(17)260528(10)KME21F0968
112482-000080	KME21F1072	(01)14026704196521(17)260528(10)KME21F1072
112482-000080	KME21F1158	(01)14026704196521(17)260528(10)KME21F1158
112482-000085	KME21D1883	(01)14026704196538(17)260328(10)KME21D1883
112482-000085	KME21E1319	(01)14026704196538(17)260428(10)KME21E1319
112482-000085	KME21D1577	(01)14026704196538(17)260328(10)KME21D1577
112482-000085	KME21D1800	(01)14026704196538(17)260328(10)KME21D1800
112482-000085	KME21D1889	(01)14026704196538(17)260328(10)KME21D1889
112482-000085	KME21E0916	(01)14026704196538(17)260428(10)KME21E0916
112482-000085	KME21E0957	(01)14026704196538(17)260428(10)KME21E0957
112482-000085	KME21E1039	(01)14026704196538(17)260428(10)KME21E1039
112482-000085	KME21E1206	(01)14026704196538(17)260428(10)KME21E1206
112482-000085	KME21E1388	(01)14026704196538(17)260428(10)KME21E1388
112482-000085	KME21F0638	(01)14026704196538(17)260528(10)KME21F0638
112482-000085	KME21F0837	(01)14026704196538(17)260528(10)KME21F0837
112482-000085	KME21F1003	(01)14026704196538(17)260528(10)KME21F1003
112482-000085	KME21F1032	(01)14026704196538(17)260528(10)KME21F1032
112482-000085	KME21F1114	(01)14026704196538(17)260528(10)KME21F1114
112482-000090	KME21E1475	(01)14026704196545(17)260428(10)KME21E1475
112482-000090	KME21E1041	(01)14026704196545(17)260428(10)KME21E1041
112482-000090	KME21E1107	(01)14026704196545(17)260428(10)KME21E1107
112482-000090	KME21E1135	(01)14026704196545(17)260428(10)KME21E1135
112482-000090	KME21F0154	(01)14026704196545(17)260528(10)KME21F0154
112482-000090	KME21F0721	(01)14026704196545(17)260528(10)KME21F0721
112482-000095	KME21E0519	(01)14026704196552(17)260428(10)KME21E0519
112482-000095	KME21E1109	(01)14026704196552(17)260428(10)KME21E1109
112482-000095	KME21E1450	(01)14026704196552(17)260428(10)KME21E1450
112482-000095	KME21F0724	(01)14026704196552(17)260528(10)KME21F0724
112082050	KME21F0727	(01)14026704341389(17)260528(10)KME21F0727
112082050	KME21F0321	(01)14026704341389(17)260528(10)KME21F0321
112082055	KME21E0864	(01)14026704341396(17)260428(10)KME21E0864
112082055	KME21F0946	(01)14026704341396(17)260528(10)KME21F0946
112082055	KME21E1230	(01)14026704341396(17)260428(10)KME21E1230
112082060	KME21E0660	(01)14026704341402(17)260428(10)KME21E0660
112082060	KME21E1391	(01)14026704341402(17)260428(10)KME21E1391
112082065	KME21E0607	(01)14026704341419(17)260428(10)KME21E0607
112082065	KME21E0685	(01)14026704341419(17)260428(10)KME21E0685

Appendix 2: Product in scope of EIF-000523

Product Code	Batch	UDI Information (Unit Level)
112082065	KME21F0384	(01)14026704341419(17)260528(10)KME21F0384
112082065	KME21F0516	(01)14026704341419(17)260528(10)KME21F0516
112082065	KME21F0694	(01)14026704341419(17)260528(10)KME21F0694
112082065	KME21F0834	(01)14026704341419(17)260528(10)KME21F0834
112082065	KME21F0930	(01)14026704341419(17)260528(10)KME21F0930
112082070	KME21D2935	(01)14026704341426(17)260328(10)KME21D2935
112082070	KME21E0002	(01)14026704341426(17)260428(10)KME21E0002
112082070	KME21E0363	(01)14026704341426(17)260428(10)KME21E0363
112082070	KME21E1390	(01)14026704341426(17)260428(10)KME21E1390
112082070	KME21D1580	(01)14026704341426(17)260328(10)KME21D1580
112082070	KME21D1652	(01)14026704341426(17)260328(10)KME21D1652
112082070	KME21D1706	(01)14026704341426(17)260328(10)KME21D1706
112082070	KME21D1878	(01)14026704341426(17)260328(10)KME21D1878
112082070	KME21D2020	(01)14026704341426(17)260328(10)KME21D2020
112082070	KME21D2929	(01)14026704341426(17)260328(10)KME21D2929
112082070	KME21E0001	(01)14026704341426(17)260428(10)KME21E0001
112082070	KME21F0102	(01)14026704341426(17)260528(10)KME21F0102
112082070	KME21F0140	(01)14026704341426(17)260528(10)KME21F0140
112082070	KME21F0239	(01)14026704341426(17)260528(10)KME21F0239
112082070	KME21F0434	(01)14026704341426(17)260528(10)KME21F0434
112082075	KME21E0288	(01)14026704341433(17)260428(10)KME21E0288
112082075	KME21E1150	(01)14026704341433(17)260428(10)KME21E1150
112082075	KME21E1178	(01)14026704341433(17)260428(10)KME21E1178
112082075	KME21E1186	(01)14026704341433(17)260428(10)KME21E1186
112082075	KME21F0557	(01)14026704341433(17)260528(10)KME21F0557
112082075	KME21D1705	(01)14026704341433(17)260328(10)KME21D1705
112082075	KME21D1780	(01)14026704341433(17)260328(10)KME21D1780
112082075	KME21D1923	(01)14026704341433(17)260328(10)KME21D1923
112082075	KME21D1956	(01)14026704341433(17)260328(10)KME21D1956
112082075	KME21D2014	(01)14026704341433(17)260328(10)KME21D2014
112082075	KME21D2018	(01)14026704341433(17)260328(10)KME21D2018
112082075	KME21D2262	(01)14026704341433(17)260328(10)KME21D2262
112082075	KME21E0021	(01)14026704341433(17)260428(10)KME21E0021
112082075	KME21E0498	(01)14026704341433(17)260428(10)KME21E0498
112082075	KME21E0469	(01)14026704341433(17)260428(10)KME21E0469
112082075	KME21F0351	(01)14026704341433(17)260528(10)KME21F0351
112082075	KME21E1264	(01)14026704341433(17)260428(10)KME21E1264
112082075	KME21F0184	(01)14026704341433(17)260528(10)KME21F0184
112082075	KME21F0385	(01)14026704341433(17)260528(10)KME21F0385
112082075	KME21F0594	(01)14026704341433(17)260528(10)KME21F0594
112082075	KME21F0725	(01)14026704341433(17)260528(10)KME21F0725
112082080	KME21D2012	(01)14026704341440(17)260328(10)KME21D2012
112082080	KME21F0981	(01)14026704341440(17)260528(10)KME21F0981
112082080	KME21D1653	(01)14026704341440(17)260328(10)KME21D1653
112082080	KME21D1946	(01)14026704341440(17)260328(10)KME21D1946
112082080	KME21D2011	(01)14026704341440(17)260328(10)KME21D2011
112082080	KME21E0584	(01)14026704341440(17)260428(10)KME21E0584
112082080	KME21E0687	(01)14026704341440(17)260428(10)KME21E0687
112082080	KME21E1452	(01)14026704341440(17)260428(10)KME21E1452
112082080	KME21E0770	(01)14026704341440(17)260428(10)KME21E0770
112082080	KME21F1006	(01)14026704341440(17)260528(10)KME21F1006
112082085	KME21E0902	(01)14026704341457(17)260428(10)KME21E0902

Appendix 2: Product in scope of EIF-000523

Product Code	Batch	UDI Information (Unit Level)
112082085	KME21E1291	(01)14026704341457(17)260428(10)KME21E1291
112082085	KME21F0741	(01)14026704341457(17)260528(10)KME21F0741
112082095	KME21E1393	(01)14026704341471(17)260428(10)KME21E1393
170060	KME21E1411	(01)04026704319558(17)260428(10)KME21E1411
170070	KME21D1430	(01)04026704319572(17)260328(10)KME21D1430
170070	KME21E1290	(01)04026704319572(17)260428(10)KME21E1290
170075	KME21D1785	(01)04026704319589(17)260328(10)KME21D1785
170080	KME21D1204	(01)04026704319596(17)260328(10)KME21D1204
170080	KME21E1287	(01)04026704319596(17)260428(10)KME21E1287
170085	KME21E1296	(01)4026704319602(17)260428(10)KME21E1296
170090	KME21E1232	(01)4026704319619(17)260428(10)KME21E1232
504550	KME21D1964	(01)4026704547579(17)260328(10)KME21D1964
504560	KME21E1412	(01)4026704547593(17)260428(10)KME21E1412
504580	KME21D1784	(01)4026704547630(17)260328(10)KME21D1784
504580	KME21E1286	(01)4026704547630(17)260428(10)KME21E1286
504585	KME21E1295	(01)4026704547647(17)260428(10)KME21E1295
504590	KME21E1114	(01)4026704547654(17)260428(10)KME21E1114

-End of Document-