



THUASNE®

IMPORTANT INFORMATION PRODUCT RECALL

**Globe-Trotter
Globe-Trotter+**

You trusted us by choosing our **Globe-Trotter and Globe-Trotter+** products and we thank you.


We would like to inform you of a **risk of non-compliance** with our

Globe-Trotter and Globe-Trotter+ crutches

Until now, our products were recommended for a maximum weight of 150 kg. **As a precautionary measure we have decided to reduce this maximum weight to 130 kg.**

CONSEQUENTLY:

If your weight is greater than 130 kg, please no longer use our products and bring them back to your supplier.



If your weight is less than 130 kg, you can continue to use our products. However, please contact your supplier who will provide you with an addendum to the instructions for use as well as the sticker to affix to the aluminium tube.

**W2016021010
Globe-Trotter**



**W2016022010
Globe-Trotter**



**W2017022006
Globe-Trotter+**



*Please accept our apologies for the inconvenience caused by this operation. More than ever, we are committed to ensuring patient safety.
We remain at your disposal for any questions.*



Thuasne
120, rue Marius Auphan
92300 Levallois-Perret (France)
www.thuasne.com

©Thuasne - 0000000 - (2024-05)



Réserve pour labels certif. papier
FSC ou PEFC ou
À CHANGER avant envoi de l'EN, cf BDC